



Disposable Shoe Covers  
**Soft Line™**

DISCOVER BATRIK  
**SHOOZ**COVERS™



Disposable Shoe Covers  
**Coated Line™**



Disposable Shoe Covers  
**Tough Line™**



**Grippez Shooz™**  
Disposable Shoe Covers  
Rhino Grey



**Grippez Shooz™**  
Disposable Shoe Covers  
Gecko Navy

# ABOUT US

The **SURGMED Group of Companies** initially began its operations in 1992 as a medical device distribution company known as Advanced Surgi-Pharm Inc. Quickly becoming recognized as an industry leader serving the OBGYN, L&D, NICU, and Operating Room.

SURGMED made strategic acquisitions and greatly expanded into the direct manufacturing of medical devices under the brands of **BATRIK** and **IMPERIAL SURGICAL**. Today The SURGMED Group's North American manufactured medical devices are being sold in over 85 countries.

2165 46th Ave., Lachine, QC, Canada, H8T 2P1  
Phone: +1.514.631.7988 | Toll Free: 800.661.5432 | Fax.: +1.514.631.9083  
Email: [info@surgmed.com](mailto:info@surgmed.com) | Website: [www.batrik.com](http://www.batrik.com)



# SHOOZCOVERS™ | SOFT LINE™


Batrik ShoozCovers™ Soft Line™ are manufactured using highest quality 35 g/m<sup>2</sup> fabric supported with 100% latex-free USA manufactured elastics for performance and durability you can rely on. Our anti-skid tread reduces the risk of slipping and our smart design delivers easy on-and-off convenience and great fit.

## SHOOZ-REPELLENT™ TECHNOLOGY:

The Soft Line™ integrates Shooz-Repellent™ technology which creates a protective barrier that repels fluids, delivering basic incidental fluid contact repellent. Shoe covers are not waterproof or water-resistant.

## ANTI-SKID TREAD:

Our anti-skid tread reduces the risk of slipping with triple ultrasonic welds that increases tear strength and delivers easy on-and-off convenience.



MADE IN CANADA 



## OUR SOFT LINE™ IS PERFECT FOR:



MEDICAL  
ESTABLISHMENTS



PHARMACEUTICAL  
LABS



GMP  
SITES



ELECTRONIC  
MANUFACTURING

LOLLY BLUE | **SOFT LINE™**

WHY YOU'LL LOVE THE SOFT LINE™:



- Batrik ShoozCovers™ Soft Line™ are manufactured using highest quality 35 g/m² fabric supported with 100% latex-free USA manufactured elastics.
- This line is particle count tested and complies with IEST-RP-CC003.4.
- Other colors available upon request.

UN XL

AVAILABLE IN TWO VERSATILE SIZES

BUBBLE GUM PINK | **SOFT LINE™**



SUPPORT A CAUSE, WEAR PINK:

- A portion of our Bubble Gum Pink shoe cover sales will be donated to charity to help further breast cancer research. Join us in making a meaningful contribution by choosing Pink ShoozCovers™!

UN

AVAILABLE IN UNIVERSAL

PRODUCT SELECTION:

ITEM #	DESCRIPTION	COLOR	SIZE	PACKAGING	CASE DIMENSIONS
SC-1210-PB	ShoozCovers™ - Soft Line™ <i>polypropylene non-woven fabric</i>	Lolly Blue	Universal Up to shoe size 12 ~ 42 cm x 16.5 cm / 16.5" x 6.5"	Case of 1,000 (10 bags of 100 units)	20"x 16"x 16"
SC-1216-PDB	ShoozCovers™ - Soft Line™ <i>polypropylene non-woven fabric</i>	Lolly Blue	Universal Shoe size 12 ~ 42 cm x 16.5 cm / 16.5" x 6.5"	Case of 600 (6 Dispenser Boxes of 100 units)	16.25"x 8.875"x 18.75"
SC-1216-PB	ShoozCovers™ - Soft Line™ <i>polypropylene non-woven fabric</i>	Lolly Blue	Universal Shoe size 12 ~ 42 cm x 16.5 cm / 16.5" x 6.5"	Case of 600 (6 bags of 100 units)	19.5"x 19.5"x 14.25"
SC-1226-PDB	ShoozCovers™ - Soft Line™ <i>polypropylene non-woven fabric</i>	Lolly Blue	XL Shoe size 12-15 ~ 47 cm x 19 cm / 18.5" x 7.5"	Case of 600 (6 Dispenser Boxes of 100 units)	17.25"x 9.375"x 21.75"
SC-1220-PB	ShoozCovers™ - Soft Line™ <i>polypropylene non-woven fabric</i>	Lolly Blue	XL Shoe size 12-15 ~ 47 cm x 19 cm / 18.5" x 7.5"	Case of 1,000 (10 bags of 100 units)	22"x 16"x 16"
SC-1616-PDB	ShoozCovers™ - Soft Line™ <i>polypropylene non-woven fabric</i>	Bubble Gum Pink	Universal Shoe size 12-15 ~ 42 cm x 16.5 cm / 16.5" x 6.5"	Case of 600 (6 Dispenser Boxes of 100 units)	16.25"x 8.875"x 18.75"
SC-1616-PB	ShoozCovers™ - Soft Line™ <i>polypropylene non-woven fabric</i>	Bubble Gum Pink	Universal Shoe size 12-15 ~ 42 cm x 16.6 cm / 16.5" x 6.5"	Case of 600 (6 bags of 100 units)	19.5"x 19.5"x 14.25"

ALTERNATIVE PACKAGING: Add "-PDB" to end of code for a dispenser box instead of bags.

<sup>1</sup>Many suppliers claim their products are latex free, however this often only relates to the fabric and not the elastics used. We source our ShoozCovers™ elastics from the USA to ensure that our components conform fully as a 100% latex free product.




# SHOOZCOVERS™ | COATED LINE™

## STAY DRY:

Our Batrik ShoozCovers™ Coated Line™ is made from a ~55 g/m<sup>2</sup> blend of polypropylene and polyethylene, providing high-quality water resistance.<sup>1</sup> Ideal for environments where maximum protection against fluids is required.

## TRIPLE WELDED TECHNOLOGY:

These shoe covers are reinforced with **triple ultrasonic welds**, a technique that joins materials together without the need for adhesives or stitching. This technique **increase tear strength and reduce linting**, for long lasting, resilient covers.



MADE IN CANADA 



## OUR COATED LINE™ IS PERFECT FOR:



Medical  
Establishments



Pharmaceutical  
Labs



GMP Sites



Electronic  
Manufacturing

<sup>1</sup>Our ShoozCovers™ Coated Line™ is resistant against all non petroleum based liquids.

MARSHMALLOW WHITE | COATED LINE™

WHY YOU'LL LOVE THE COATED LINE™:



- **Superior protection using our unique blend of fabrics** which protects against not only water but all non-petroleum-based liquids, making them the ideal choice for environments where spills and leaks are a concern.
- **Elastics around the ankle create a strong and secure fit** which is especially important when working in hazardous environments where a misplaced shoe cover could lead to contamination or injury.
- **Smart design** with a large opening, providing a better fit and easy on-and-off convenience.
- **Other colors available upon request.**



AVAILABLE IN TWO VERSATILE SIZES

PRODUCT SELECTION:

ITEM #	DESCRIPTION	COLOR	SIZE	PACKAGING
SC-2112-NPB	ShoozCovers™ Coated Line™ <i>Polypropylene + polyethylene non-woven fabric, waterproof</i>	Marshmallow White	Universal Up to shoe size 12 ~ 42 cm x 17.5 cm / 16.5" x 6.88"	Case of 300 (2 bags of 150 units)
SC-2112-NPDB	ShoozCovers™ Coated Line™ <i>Polypropylene + polyethylene non-woven fabric, waterproof</i>	Marshmallow White	Universal Up to shoe size 12 ~ 42 cm x 17.5 cm / 16.5" x 6.88"	Case of 300 (2 bags of 150 units)
SC-2122-NPB	ShoozCovers™ Coated Line™ <i>Polypropylene + polyethylene non-woven fabric, waterproof</i>	Marshmallow White	XL Shoe size 12-15 ~ 47 cm x 19 cm / 18.5" x 6.88"	Case of 240 (2 bags of 120 units)
SC-2122-NPDB	ShoozCovers™ Coated Line™ <i>Polypropylene + polyethylene non-woven fabric, waterproof</i>	Marshmallow White	XL Shoe size 12-15 ~ 47 cm x 19 cm / 18.5" x 6.88"	Case of 240 (2 bags of 120 units)

ALTERNATIVE PACKAGING: Add "-PDB" to end of code for a dispenser box instead of bags

# SHOOZCOVERS™ | TOUGH LINE™

## EXTRA-LARGE FOR EXTRA-TOUGH WORK:


Our Batrik ShoozCovers™ Tough Line™ comes in an **extra-large** size that fits over a **men's size 15 steel-toed boot**, delivering optimal full coverage and protection while you work.<sup>1</sup>

## ENHANCED DURABILITY:

These shoe covers are formulated with our **thickest fabric yet!** Made of **~55 g/m² polypropylene fabric**, these highly durable covers provide added protection, increase tear strength and reduce linting.



SHOE SIZES 12-15 (45-49.5 EU)  
~ 47 CM X 19 CM / 18.5" X 7.5"



MADE IN CANADA 



## OUR TOUGH LINE™ IS PERFECT FOR:



HOMES OR  
COMMERCIAL  
BUILDINGS



SHIELDING  
AGAINST MUD  
& GREASE



MEDICAL  
ESTABLISHMENTS



PHARMACEUTICAL  
LABS

<sup>1</sup> SHOOZCOVERS™ ARE MADE LONGER (42 & 47CM) FOR ACCURATE SHOE SIZE SELECTION.

BLACK LIQUORICE | TOUGH LINE™

WHY YOU’LL LOVE THE TOUGH LINE™:

- Triple ultrasonic welds that are long lasting to reduce waste
- Batrik Tough Line™ are manufactured using highest quality ~ 55 g/m² fabric supported with 100% latex-free USA manufactured elastics to contract the shoe cover around your foot and deliver a strong and secure fit².
- Shooz-Repellent™ technology delivers basic incidental fluid contact repellency³
- Other colors available upon request



XL

AVAILABLE IN EXTRA LARGE

ANY TIME, ANYWHERE:



RESIDENTIAL



COMMERCIAL



INDUSTRIAL

PRODUCT SELECTION:

ITEM #	DESCRIPTION	COLOR	SIZE	PACKAGING
SC-1922-PB	ShoozCovers™ - Tough Line™ ~55 g/m² polypropylene non-woven fabric	Black Liquorice	XL Shoe size 12-15 ~ 47 cm x 19 cm / 18.5" x 7.5"	Case of 200 (2 bags of 100 units)

ALTERNATIVE PACKAGING: Add "-PDB" to end of code for a dispenser box instead of bags

² Many suppliers claim their products are latex free, however this often only relates to the fabric and not the elastics used. We source our ShoozCovers™ elastics from the USA to ensure that our components conform fully as a 100% latex free product. ³ Shooz-Repellent™ technology provides basic incidental fluid contact repellent. It is not waterproof nor water resistant.



# GRIPPEEZSHOOZ™



## RELIABLE, DURABLE, AND TRULY ANTI-SLIP SHOE COVERS

**ShoozCovers™ Grippeez** are the ultimate solution for professionals navigating slippery, high-risk environments. Designed for **superior grip** and **waterproof protection**, these premium shoe covers enhance safety, hygiene, and performance. Stay confident, stay protected—every step of the way.

MADE IN CANADA 🇨🇦



## PERFECT FOR:



AEROSPACE



FOOD  
PROCESSING  
INDUSTRY



OIL, GAS AND  
MINING SITES



CONSTRUCTION &  
RENOVATION SITES

<sup>1</sup>Industrial strength elastics

## GRIPPEEZSHOOZ™ ARE A PERFECT FIT FOR EVERYONE:

- **Fast, Reliable Fit:** Stretch material ensures a secure, easy-on fit.
- Comes in an **extra-large size** that fits over a **men's size 15 steel-toed boot**.
- **Enhanced Safety:** Anti-slip rubberized coating provides reliable grip on all surfaces, helping to prevent sliding in critical areas.
- Available in two colors: **Gecko Navy** and **Rhino Grey**.

### RHINO GREY



### GECKO NAVY



## PRODUCT SELECTION:

ITEM #	DESCRIPTION	COLOR	SIZE	PACKAGING	CASE DIMENSIONS
<b>SC-29312-NPDB</b>	Grippez ShoozCovers™ Polypropylene + polyethylene non-woven fabric, waterproof	Rhino Grey	Universal Up to shoe size 12 ~ 42 cm x 17.5 cm / 16.5" x 6.88"	Case of 300 (2 dispenser boxes of 150 units)	16.62" x 14.875" x 8.75"
<b>SC-29312-NPB</b>	Grippez ShoozCovers™ Polypropylene + polyethylene non-woven fabric, waterproof	Rhino Grey	Universal Shoe size 12 ~ 42 cm x 17.5 cm / 16.5" x 6.88"	Case of 300 (2 bags of 150 units)	16" x 13.375" x 10"
<b>SC-29322-NPDB</b>	Grippez ShoozCovers™ Polypropylene + polyethylene non-woven fabric, waterproof	Rhino Grey	XL Up to Size 12-15 ~ 47 cm x 17.5 cm / 18.5" x 6.88"	Case of 240 (2 dispenser boxes of 120 units)	16.62" x 14.875" x 8.75"
<b>SC-29322-NPB</b>	Grippez ShoozCovers™ Polypropylene + polyethylene non-woven fabric, waterproof	Rhino Grey	XL Up to Size 12-15 ~ 47 cm x 17.5 cm / 18.5" x 6.88"	Case of 240 (2 bags of 120 units)	16" x 13.375" x 10"
<b>SC-22212-NPDB</b>	Grippez ShoozCovers™ Polypropylene + polyethylene non-woven fabric, waterproof	Gecko Navy	Universal Up to shoe size 12 ~ 42 cm x 17.5 cm / 16.5" x 6.88"	Case of 300 (2 dispenser boxes of 150 units)	16.62" x 14.875" x 8.75"
<b>SC-22212-NPB</b>	Grippez ShoozCovers™ Polypropylene + polyethylene non-woven fabric, waterproof	Gecko Navy	Universal Up to shoe size 12 ~ 42 cm x 17.5 cm / 16.5" x 6.88"	Case of 300 (2 bags of 150 units)	16" x 13.375" x 10"
<b>SC-22222-NPDB</b>	Grippez ShoozCovers™ Polypropylene + polyethylene non-woven fabric, waterproof	Gecko Navy	XL Up to Size 12-15 ~ 47 cm x 17.5 cm / 18.5" x 6.88"	Case of 240 (2 dispenser boxes of 120 units)	16.62" x 14.875" x 8.75"
<b>SC-22222-NPB</b>	Grippez ShoozCovers™ Polypropylene + polyethylene non-woven fabric, waterproof	Gecko Navy	XL Up to Size 12-15 ~ 47 cm x 17.5 cm / 18.5" x 6.88"	Case of 240 (2 bags of 120 units)	16" x 13.375" x 10"



2165 46th Ave., Lachine, QC, Canada, H8T 2P1  
Phone: +1.514.631.7988 | Toll Free: 800.661.5432 | Fax.: +1.514.631.9083  
Email: [info@surgmed.com](mailto:info@surgmed.com) | Website: [www.batrik.com](http://www.batrik.com)