

BATRIK

TAG-IT™

Batrik Tag-It™ is our new line of information tags developed to help medical professionals identify and document a wide range of instructions. Whether you need to identify instruments, trays or carts, both before or after sterilization, Tag-It™ has got you covered. These durable and tear-resistant tags are perfect for indicating special handling instructions.

WHY YOU'LL LOVE THEM:

- The paper, tie wraps, and smudge-proof ink are designed to withstand high temperatures of up to 273.2°F (134°C).
- The synthetic paper is durable and will not wrinkle or absorb water. Made from recyclable plastic promoting environmental sustainability.
- Robust and easy to use tie wraps designed to securely attach to a wide range of trays, instruments, and other medical devices.
- Available in over 8 languages.


KEY PRODUCT FEATURES:

- **High Quality Permanent Marker**
- **Steam or Plasma indicator Dots** provide confirmation of processing
- **Eco-friendly, Waterproof & Tearproof Tags**
- **Large Quantities:** Each bundle of Tag-It™ comes with 100 versatile tags, to quickly identify instruments that require attention.
- **Available in multiple linguistic choices**



TI-3503

All tags available with indicator dots
Provides visual confirmation that the tag was processed



PRODUCT SELECTION:

ITEM #	COLOR	DESCRIPTION	SIZE	PACKAGING
TI-3502	Blue	VENDOR	3" x 5" / 8 cm x 13cm	100 / PKG
TI-3503	Green	LOANER	3" x 5" / 8 cm x 13cm	100 / PKG
TI-3504	Red	PRIORITY	3" x 5" / 8 cm x 13cm	100 / PKG
TI-3505	Yellow	RE-RUN	3" x 5" / 8 cm x 13cm	100 / PKG
TI-3506	Pink	MISSING	3" x 5" / 8 cm x 13cm	100 / PKG
TI-3507	Orange	BROKEN	3" x 5" / 8 cm x 13cm	100 / PKG
TI-3509	Black	BLANK	3" x 5" / 8 cm x 13cm	100 / PKG
TI-35B1	White	BLANK	3" x 5" / 8 cm x 13cm	100 / PKG
TI-35B4	Red	BROKEN	3" x 5" / 8 cm x 13cm	100 / PKG
TI-35C22	Light Blue	VENDOR CHECK IN + OUT	3" x 5" / 8 cm x 13cm	100 / PKG
TI-35L3	Green	LOANER CHECK IN + OUT	3" x 5" / 8 cm x 13cm	100 / PKG
TI-35T2	Blue	QUICK TURNOVER	3" x 5" / 8 cm x 13cm	100 / PKG
TI-PM09	Black	FINE PERMANENT MARKERS	1.5" Ø x 5.5" / 3.81cm Ø x 12.7 cm	36 / Bag

All information tags are single-use. **PACKAGING:** 1 pack of 100 tags, 100 tie wraps, and 1 permanent marker.

PROCESSING INDICATOR DOTS: Add suffix "ST" to end of code for steam dots.

PROCESSING INDICATOR DOTS: Add suffix "PL" to end of code for plasma dots.

PRODUCT SELECTION:

VENDOR / BLUE

ITEM #	SIZE
TI-3502	3" x 5" (8 cm x 13cm)



VENDOR

Company Name _____


Tray Name _____

Date of Surgery _____

WWW.BATRIK.COM TAG-IT™ TI-3502

BLANK / WHITE


ITEM #	SIZE
TI-35B1	3" x 5" (8 cm x 13cm)



WWW.BATRIK.COM TAG-IT™ TI-35B1

LOANER / GREEN

ITEM #	SIZE
TI-3503	3" x 5" (8 cm x 13cm)



LOANER

Loaned From _____

Loaned To _____

By _____ Date _____

WWW.BATRIK.COM TAG-IT™ TI-3503

BROKEN / RED

ITEM #	SIZE
TI-35B4	3" x 5" (8 cm x 13cm)



BROKEN DO NOT USE

DATE: _____ PHONE: _____

DEPARTMENT/CONTACT: _____

ITEM: _____

CRITICAL: (☠) Yes No

INSTRUMENT SET: _____

ASSEMBLED BY: _____

WWW.BATRIK.COM TAG-IT™ TI-35B4

PRIORITY / RED

ITEM #	SIZE
TI-3504	3" x 5" (8 cm x 13cm)



PRIORITY

Room # _____

Time Needed _____

Initials _____

WWW.BATRIK.COM TAG-IT™ TI-3504

VENDOR IN + OUT / LIGHT BLUE



VENDOR
Check In

Tray Received From: _____

Tray Name: _____

Received By: _____

Date: _____ Surgery Date: _____

☐ CPD Notified Prior to Receipt of Trays
☐ Received >24 Hours Before Case
☐ Inventory List Available
☐ Inventory Check Complete
☐ All Instruments in Good Condition
☐ Containers in Good Condition
☐ All Trays Labeled (Patient & Surgeon)
☐ All Trays <25 pounds
☐ Written Instructions For Use Available
☐ Disassembled/Decontaminated

WWW.BATRIK.COM TAG-IT™ TI-35C22

VENDOR
Check Out

☐ Disassembled/Decontaminated
☐ Inventory Check Complete
☐ All Items Complete and Accounted For

Returned By: _____

Received By: _____

Date: _____

Notes: _____

WWW.BATRIK.COM TAG-IT™ TI-35C22

ITEM #
TI-35C22
SIZE
5" x 3" (13 cm x 8 cm)

RE-RUN / YELLOW

ITEM #	SIZE
TI-3505	3" x 5" (8 cm x 13cm)



RE-RUN

Room # _____

WWW.BATRIK.COM TAG-IT™ TI-3505

MISSING / PINK


ITEM #	SIZE
TI-3506	3" x 5" (8 cm x 13cm)



MISSING

WWW.BATRIK.COM TAG-IT™ TI-3506

LOANER IN + OUT / GREEN



LOANER
Check In

Tray Received From: _____

Tray Name: _____

Received By: _____

Date: _____ Surgery Date: _____

☐ CPD Notified Prior to Receipt of Trays
☐ Received >24 Hours Before Case
☐ Inventory List Available
☐ Inventory Check Complete
☐ All Instruments in Good Condition
☐ Containers in Good Condition
☐ All Trays Labeled (Patient & Surgeon)
☐ All Trays <25 pounds
☐ Written Instructions For Use Available
☐ Disassembled/Decontaminated

WWW.BATRIK.COM TAG-IT™ TI-35L3

LOANER
Check Out

☐ Disassembled/Decontaminated
☐ Inventory Check Complete
☐ All Items Complete and Accounted For

Returned By: _____

Received By: _____

Date: _____


Notes: _____

WWW.BATRIK.COM TAG-IT™ TI-35L3

ITEM #
TI-35L3
SIZE
5" x 3" (13 cm x 8 cm)

BROKEN / ORANGE

ITEM #	SIZE
TI-3507	3" x 5" (8 cm x 13cm)



BROKEN
DO NOT USE

DATE: _____ PHONE: _____

DEPARTMENT/CONTACT: _____

ITEM: _____

CRITICAL: (☠) Yes No


WWW.BATRIK.COM TAG-IT™ TI-3507

BLANK / BLACK

ITEM #	SIZE
TI-3509	3" x 5" (8 cm x 13cm)



WWW.BATRIK.COM TAG-IT™ TI-35B1



QUICK TURNOVER

Tech Initials: _____

Room Number: _____

Item Number: _____

Time Needed: _____

☐ Decontam/LLD/LD/HLD
☐ Immediate Use Sterilization
☐ Terminal Sterilization in SPD

WWW.BATRIK.COM TAG-IT™ TI-35T2

QUICK TURNOVER/ BLUE

ITEM #	SIZE
TI-35T2	3" x 5" (8 cm x 13cm)

All information tags are single-use. PACKAGING: 1 pack of 100 tags, 100 tie wraps, and 1 permanent marker.