

BATRIK

TAG-IT™

Batrik Tag-It™ is our new line of information tags developed to help medical professionals identify and document a wide range of instructions. Whether you need to identify instruments, trays or carts, both before or after sterilization, Tag-It™ has got you covered. These durable and tear-resistant tags are perfect for indicating special handling instructions.

WHY YOU'LL LOVE THEM:

- The paper, tie wraps, and smudge-proof ink are designed to withstand high temperatures of up to 273.2°F (134°C) and are compatible with Gas plasma sterilization (STERRAD®100 NX Sterilization System).
- The synthetic paper is durable and will not wrinkle or absorb water. Made from recyclable plastic promoting environmental sustainability.
- Robust and easy to use tie wraps designed to securely attach to a wide range of trays, instruments, and other medical devices.
- Available in over 8 languages.


KEY PRODUCT FEATURES:

- **High Quality Permanent Marker**
- **Steam or Plasma indicator Dots** provide confirmation of processing
- **Eco-friendly, Waterproof & Tearproof Tags**
- **Large Quantities:** Each bundle of Tag-It™ comes with 100 versatile tags, to quickly identify instruments that require attention.
- **Available in multiple linguistic choices**



TI-3503

All tags available with indicator dots
Provides visual confirmation that the tag was processed



PRODUCT SELECTION:

ITEM #	COLOR	DESCRIPTION	SIZE	PACKAGING
TI-3502	Blue	VENDOR	3" x 5" / 8 cm x 13cm	100 / PKG
TI-3503	Green	LOANER	3" x 5" / 8 cm x 13cm	100 / PKG
TI-3504	Red	PRIORITY	3" x 5" / 8 cm x 13cm	100 / PKG
TI-3505	Yellow	RE-RUN	3" x 5" / 8 cm x 13cm	100 / PKG
TI-3506	Pink	MISSING	3" x 5" / 8 cm x 13cm	100 / PKG
TI-3507	Orange	BROKEN	3" x 5" / 8 cm x 13cm	100 / PKG
TI-3509	Black	BLANK	3" x 5" / 8 cm x 13cm	100 / PKG
TI-35B1	White	BLANK	3" x 5" / 8 cm x 13cm	100 / PKG
TI-35B4	Red	BROKEN	3" x 5" / 8 cm x 13cm	100 / PKG
TI-35C22	Light Blue	VENDOR CHECK IN + OUT	3" x 5" / 8 cm x 13cm	100 / PKG
TI-35L3	Green	LOANER CHECK IN + OUT	3" x 5" / 8 cm x 13cm	100 / PKG
TI-35T2	Blue	QUICK TURNOVER	3" x 5" / 8 cm x 13cm	100 / PKG
TI-PM09	Black	FINE PERMANENT MARKERS	1.5" Ø x 5.5" / 3.81cm Ø x 12.7 cm	36 / Bag

All information tags are single-use. **PACKAGING: 1 pack of 100 tags, 100 tie wraps, and 1 permanent marker.**

PROCESSING INDICATOR DOTS: Add suffix "ST" to end of code for steam dots.

PROCESSING INDICATOR DOTS: Add suffix "PL" to end of code for plasma dots.

PRODUCT SELECTION:

<p>VENDOR / BLUE</p> <p>ITEM # SIZE</p> <hr/> <p>TI-3502 3" x 5" (8 cm x 13cm)</p>			<p>BLANK / WHITE</p> <p>ITEM # SIZE</p> <hr/> <p>TI-35B1 3" x 5" (8 cm x 13cm)</p>
<p>LOANER / GREEN</p> <p>ITEM # SIZE</p> <hr/> <p>TI-3503 3" x 5" (8 cm x 13cm)</p>			<p>BROKEN / RED</p> <p>ITEM # SIZE</p> <hr/> <p>TI-35B4 3" x 5" (8 cm x 13cm)</p>
<p>PRIORITY / RED</p> <p>ITEM # SIZE</p> <hr/> <p>TI-3504 3" x 5" (8 cm x 13cm)</p>			<p>VENDOR IN + OUT / LIGHT BLUE</p> <p>ITEM #</p> <hr/> <p>TI-35C22</p> <p>SIZE</p> <hr/> <p>5" x 3" (13 cm x 8 cm)</p>
<p>RE-RUN / YELLOW</p> <p>ITEM # SIZE</p> <hr/> <p>TI-3505 3" x 5" (8 cm x 13cm)</p>			<p>LOANER IN + OUT / GREEN</p> <p>ITEM #</p> <hr/> <p>TI-35L3</p> <p>SIZE</p> <hr/> <p>5" x 3" (13 cm x 8 cm)</p>
<p>MISSING / PINK</p> <p>ITEM # SIZE</p> <hr/> <p>TI-3506 3" x 5" (8 cm x 13cm)</p>			<p>BROKEN / ORANGE</p> <p>ITEM # SIZE</p> <hr/> <p>TI-3507 3" x 5" (8 cm x 13cm)</p>
<p>BLANK / BLACK</p> <p>ITEM # SIZE</p> <hr/> <p>TI-3509 3" x 5" (8 cm x 13cm)</p>			<p>QUICK TURNOVER / BLUE</p> <p>ITEM # SIZE</p> <hr/> <p>TI-35T2 3" x 5" (8 cm x 13cm)</p>

All information tags are single-use. PACKAGING: 1 pack of 100 tags, 100 tie wraps, and 1 permanent marker.