# When Simply Broken doesn't cut it -



TAG IT, TELL IT,

OWN IT!





56 01 60

24 25


80 23 TÜE


> 04 05 20


F

JUL NUL

MAR APR MAY











#### NEW! BATRIK



# WHEN SIMPLY BROKEN DOESN'T CUT IT—TAG IT, TELL IT, OWN IT!

Tag-It<sup>™</sup> Notes offer the freedom to **identify and document instructions for instruments, trays, and carts** before and after **sterilization**. Tag-It<sup>™</sup> Notes are **tear-resistant, withstand extreme chemicals, high temperature** and ensures **reliable communication**.

#### WHY YOU'LL LOVE THEM:

- Tag-It™ Notes gives you the freedom to clearly document important information without the clutter.
- Robust and easy-to-use tie wraps designed to securely attach to instruments (134°C & chemical tested).
- Includes a fast-drying, heat-resistant permanent marker for writing your important notes.
- Choose from other customizable cards for personalized needs or pre-made cards for fast, convenient use.
- When you need more than a "repair" or "sharpen" seal to identify an issue. Instead, Tag it, Tell it, Own it.



#### WATERPROOF

The synthetic paper withstands spills and liquid exposure without any damage.

#### **ECO-FRIENDLY**

Made from sustainable materials.

#### **HEAT RESISTANT**

Can resist high temperatures of up to 273.2°F (134°C). Compatible with Gas Plasma Sterilization (STERAD®100 NX Sterilization).

#### **DURABLE & TEAR-PROOF**

Handles rigorous cleaning, tough conditions, and resists tearing.

#### **PRODUCT SELECTION**

ITEM #	COLOR		SIZE	PACKAGING
TI-NOTE-R0	Yellow	NOTES	3" x 5" / 8 cm x 13cm	100 / PKG

All Tag-It™ are single-use. PACKAGING: 1 pack of 100 tags, 1 permanent marker, and 100 ties.





Tag-It™ Cards are an effortless way to add information with either pre-printed or blank formats. These **durable** tags streamline your workflow, **withstand tough conditions**, and feature **(optional) indicator dots** for quick verification. They help keep **instruments**, **trays**, and **carts** organized, ensuring efficiency before and after sterilization/decontamination.

#### **KEY PRODUCT FEATURES:**

- Special Permanent Marker: The ink will not run or appear blurry after exposure to chemicals or high temperature (max. 260°C/500°F) sterilization. (Regular pens or markers can be used only if the tags are not going to be exposed to any sterilization/decontamination).
- **Tie Wraps:** Assist in securing the tags to almost any instrument or container. (*Tie wraps have been tested for chemical and steam sterilization*).
- Our Tag-It™: Will not shrink, absorb water/chemicals, or lose its intended integrity after chemical or steam sterilization.



(Optional)
All Tags Available
with Indicator Dots
Provides visual confirmation that
the tag was processed.

CLICK BELOW TO WATCH OUR VIDEO



#### PRODUCT SELECTION:

ITEM #	COLOR	DESCRIPTION	SIZE	PACKAGING
TI-3502	Blue	VENDOR	3" x 5" / 8 cm x 13cm	100 / PKG
TI-3503	Green	LOANER	3" x 5" / 8 cm x 13cm	100 / PKG
TI-3504	Red	PRIORITY	3" x 5" / 8 cm x 13cm	100 / PKG
TI-3505	Yellow	RE-RUN	3" x 5" / 8 cm x 13cm	100 / PKG
TI-3506	Pink	MISSING	3" x 5" / 8 cm x 13cm	100 / PKG
TI-3507	Orange	BROKEN	3" x 5" / 8 cm x 13cm	100 / PKG
TI-3509	Black	BLANK	3" x 5" / 8 cm x 13cm	100 / PKG
TI-35B1	White	BLANK	3" x 5" / 8 cm x 13cm	100 / PKG
TI-35B4	Red	BROKEN	3" x 5" / 8 cm x 13cm	100 / PKG
TI-35C22	Light Blue	VENDOR CHECK IN + OUT	3" x 5" / 8 cm x 13cm	100 / PKG
TI-35L3	Green	LOANER CHECK IN + OUT	3" x 5" / 8 cm x 13cm	100 / PKG
TI-35T2	Blue	QUICK TURNOVER	3" x 5" / 8 cm x 13cm	100 / PKG
TI-PM09	Black	FINE PERMANENT MARKERS	1.5" Ø x 5.5" / 3.81cm Ø x 12.7 cm	36 / Bag

PACKAGING: 1 pack of 100 tags, 100 tie wraps, and 1 permanent marker.

OPTIONAL PROCESSING INDICATOR DOTS: Add suffix "ST" to end of code for steam dots, "PS' for plasma dots.

Available in other languages or customize your own.



### PRODUCT SELECTION: (customization also available)












#### **LOANER / GREEN**

ITEM #	SIZE
TI-3503	3" x 5" (8 cm x 13cm)





ITEM #	SIZE
TI-35B1	3" x 5" (8 cm x 13cm)

#### **PRIORITY / RED**

ITEM #	SIZE
TI-3504	3" x 5" (8 cm x 13cm)





DATE:	_ PHONE:	
DEPARTMENT/CONTACT:		ITEM
CRITICAL: (%) Yes No		TI-35E
ASSEMBLED BY:	TI-3584	

# **BROKEN / RED**

ITEM #	SIZE
TI-35B4	3" x 5" (8 cm x 13cm)

#### **RE-RUN / YELLOW**

ITEM #	SIZE
TI-3505	3" x 5" (8 cm x 13cm)



#### **VENDOR IN + OUT / LIGHT BLUE**



ITEM #		
TI-35C22		
SIZE		
5" x 3"		
(13 cm x 8 cm)		

#### **MISSING / PINK**

ITEM #	SIZE
TI-3506	3" x 5" (8 cm x 13cm)



#### **LOANER IN + OUT / GREEN**



ITEM #			
TI-35L3			
SIZE			
5" x 3"			
(13 cm x 8 cm)			


#### **BROKEN / ORANGE**

ITEM #	SIZE
TI-3507	3" x 5" (8 cm x 13cm)



#### **BLANK / BLACK**

ITEM #	SIZE
TI-3509	3" x 5" (8 cm x 13cm)



# **QUICK TURNOVER**

Time Needed:	Decontam/LLD/ILD/HLD	
	☐ Immediate Use Sterilization ☐ Terminal Sterilization in SPD	
WWW.BATRIK.COM	TAG-IT <sup>TM</sup>	TI-35T

#### **QUICK TURNOVER/ BLUE**

ITEM #	SIZE
TI-35T2	3" x 5" (8 cm x 13cm)