



COMPLETE PRODUCT CATALOGUE

SERVING CANADIAN HEALTHCARE WORKERS SINCE 1992

ADVANCED SURGI-PHARM INC.

Founded in 1992, we are the distribution arm of Surgmed™ in Canada, representing leading US, international and domestic manufacturers for a wide range of specialty products with unique benefits for applications in the operating room, women's health, labour & delivery, ICU, and sterile processing areas. Equipped with a direct sales team, we provide effective nationwide coverage as the Canadian distributor of Batrik, Vygon, Rocket Medical, Cooper Surgical, SunMed, and other industry leading products including a wide-ranging portfolio of top-rated PPE products.



SURGMED

GROUP 

Proudly Canadian, **The Surgmed™ Group of Companies** began operations in 1992 as a single distribution company known as **Advanced Surgi-Pharm Inc.** The company quickly became known as an industry leader and specialty distributor of value-added surgical products in the Canadian hospital market. Five years later, **Surgmed expanded upon its success with the launch of Batrik Medical Manufacturing Inc.** in 1997 and the subsequent acquisition of **Imperial Surgical Ltd.** one year later in 1998.

Together, our three unique companies serve the Canadian and global healthcare markets with high quality and competitively priced products, all backed by world-class customer care.



Founded in 1992, Advanced Surgi-Pharm Inc. is the distribution arm of Surgmed™ in Canada, representing leading US, international and domestic manufacturers for a wide range of specialty products with unique benefits for applications in the Operating Room, Women's Health, Labor & Delivery, ICU, and Sterile Processing areas. Equipped with a direct sales team, ASP provides effective nationwide coverage as the Canadian distributor of Batrik, Vygon, Rocket Medical, Cooper Surgical, SunMed, and other industry leading products including a wide-ranging portfolio of top-rated PPE products.

Established in 1997 and sold in 85 countries globally, Batrik Medical Manufacturing Inc. specializes in the manufacturing of endoscopy accessories, enzymatic solutions, infection control, central sterilization, and surgical devices. Batrik's in-house controlled environment production facility adheres to strictly enforced quality and regulatory control systems. They are registered as a device manufacturing establishment with the United States FDA and are in full compliance with CE and Health Canada Regulations. All Batrik products are rigorously tested, liability insured, and meet/exceed high international quality standards.

Established in 1935, Imperial Surgical Ltd. is one of Canada's oldest hospital capital equipment manufacturers, well-known for their stainless steel healthcare equipment for Sterile Processing, Operating Room, Examination and Treatment, Patient Handling, and Mortuary/ Pathology areas. Their growing portfolio of products also has applications in Infection Control & Prevention, an important area of focus for their future.

**SURGICAL PRODUCTS
&
OPERATING ROOM**

- 06—Suture Clamp Tags
- 07—Clamp Jaw Covers
- 08—Vascular Silicone Ties
- 09—G.O.L.F.F.™ Anti-Fog Solution
- 10—SKINS™ Surgical Multi-Tip Skin Markers
- 11—PandaPatch™ Eye Safety Tape **NEW!**
- 13—Magnetic Mat
- 14—Pin Guards
- 15—Cautery Tip Cleaner

- 16—Safe-Soak™ SPAC & Floor Mats
- 17—Blue Cloud-Gel™ Anti-Fatigue Mats
- 18—Blue-Cloud Nimbus™ Anti-Fatigue Mats
- 19—Stackable Step Stools
- 20—Tourny™+ Tourniquets
- 22—Surgi-Kit® Instrument Holders
- 23—Lap-Kit® Instrument Holders
- 24—Warm Gel® Instant Heat Packs
- 25—K-Wires, Steinmann Pins & Sterilization Racks

**ENDOSCOPY
PRODUCTS**

- 28—Endoscopy BO2™ Bite Blocks
- 29—Endoscopy Bitestrong™ Bite Blocks
- 30—Couvrette™ Valve Covers
- 31—Daybrush™ Channel Brushes

- 32—ARC™ Channel Brushes
- 35—Valve / Control Brushes
- 36—V-Brushes

CLEANING BRUSHES

- 37—ARC™ Wire Style Brushes
- 41—Double ARC™ Wire Brushes
Asymmetric / Symmetric
- 44—In-Brush™ Antimicrobial Twisted Wire Brushes
- 46—BrightBrush - Color Beaded Twisted Wire Brush **NEW!**
- 47—Stainless Steel ARC™ Micro Brushes
- 48—Instrument Cleaning Flat Brushes

- 49—Precision Brush
- 50—Trifecta™ Brush **NEW!**
- 51—All Purpose Reusable Brushes
- 52—Acetabular Reamer Brush
- 53—Orthopedic Rasps
- 54—Bottle & Drain Clesning Brushes
- 55—Autoclave Cleaning Brush Kits

**INFECTION CONTROL
&
STERILE PROCESSING**

- 56—Silicone Pin Mat
- 57—Tag-It™ **NEW!**
- 62—Corner Pockets™ & Soft Pockets™
- 63—Lila's Instrument Staine Remover
- 64—Instrument Protectors
- 67—PrAmLi Power³ Lube™ Instrument Lubricant
- 68—S.T.I.C.S.
- 70—ID Tape
- 71—Tape-B-Gone™ Tape Remover
- 72—Padlocks Security Tags
- 76—Scrub Pad
- 77—Your First ✓™ Enzymatic & Non-Enzymatic

- Cleaners
- 78—Your First ✓™ Enviro-Sachet™
- 79—Your First ✓™ Foaming Spray
- 80—Your First ✓™ F1rst Wipes
- 81—Your First ✓™ Dry Sponges
- 82—Your First ✓™ RightSwipe™ Enzymatic Sponges
& Bedside Sponge Kits **NEW!**
- 83—Your First ✓™ RightSwipe™ Non-Enzymatic
Detergent Bedside Sponge Kits **NEW!**
- 84—Your First ✓™ Transport Bags For Endoscope
And Medical Instruments

**PERSONAL
PROTECTIVE
EQUIPMENT**

FACE SHIELDS

- 85—SAVEFACE™ Draped Shields

SHOE COVERS

- 87—SHOOZCOVERS™ SOFT LINE™
- 89—SHOOZCOVERS™ GRIPPEEZ™
- 91—SHOOZCOVERS™ COATED LINE™
- 93—SHOOZCOVERS™ TOUGH LINE™

REUSABLE EQUIPMENT

- 95—Antimicrobial IV Poles With Cuverro Shield™
- 97—Plexiglass Protection Panels
- 98—Stainless Steel Hands-Free Door Opener



| | | |
|---------------------------------------|--|---|
| NEONATAL | ENTERAL NUTRITION 99—Naso-Gastric Feeding Tubes 100—Dual-Flow Gastric Feeding Tubes 101—Gastro-Duodenal Feeding Tubes 102—NutriSafe2® Feeding System | THERMOREGULATION 117—Neohelp™ Heat Loss Prevention Suit 121—Warm Gel® Infant Heel Warmer 122—TransWarmer® Transport Mattress |
| | VASCULAR ACCESS 106—Premicath® PICC 1Fr Catheters 108—Nutriline® PICC 2Fr Catheters 109—Nutriline Twinflo® PICC 2Fr Dual Lumen Catheters 110—Introducers & Microflash 111—Grip-Lok® Catheter Securement Device 112—Umbilical Catheters 115—Exchange Transfusion Tray 116—Neoventa - Goldtrace® | RESPIRATORY SUPPORT 123—Surfcath™ Surfactant Administration Catheter 125—Endotracheal Tubes 126—Mucus Extractors 127—Nasal Ventilation Tubes 128—Safe Spot Resuscitator/Ventilator 130—Inca® Infant Nasal CPAP Assembly 131—Neo-Fit™ Neonatal Endotracheal Tube Grip & Neo-Prep™ Skin Barrier |
| OBSTETRICS/ WOMEN'S HEALTH | 132—Os Finder™ Cervical Os and Canal Dilator 133—Tampa Infusion Catheter 134—H/S Catheters 135—H/S Elliptosphere™ Catheter Set and Tray 136—Comfort Belly Bands | 137—Laminaria Cervical Dilator |
| IV THERAPY | 138—LeaderCath Arterial Catheters 139—Lectrospiral® Coiled Extension Tubes 140—Lectrocath Extension Lines 141—Biocap Needle-free Drug Aspiration | 142—IV Accessories 143—Infu-Surg® Pressure Infusion Bag 144—Infu-Surg® Clarity Pressure Infusion Bag |
| DRAINAGE | 145—Jet-Vac® Bottles | 148—ASEPT® Pleural & Peritoneal Drainage System |
| RESPIRATORY TRACT | 150—Bronchitrac® "L" Set 151—COMFIT™ Endotracheal Tube Holder | 152—Esophageal Balloon Catheter Set |
| EMERGENCY | 153—Tum-E-Vac® Gastric Lavage Kit | 154—Lavage Tubing |
| WOUND IRRIGATION | 155—Irrijet® Splatter Guard | 157—Combiport® + Combiguard II™ |
| PATIENT HANDLING | 158—Smooth Mover Patient Shifters | 159—Adjustable Patient Positioning Devices |

SUTURE CLAMP TAGS

Non-toxic, radiopaque, and latex free, our Suture Clamp Tags are designed to protect and safely hold sutures without any damage.

Additionally, its foam base, affixed to a sterile area with the clamps, allows for easy insertion.

- Packed in a convenient foam base with easy peel-off adhesive backing, even while wearing gloves
- Available in sterile and non-sterile
- Available in 2, 3, or 5 pairs
- Pre-counted



OUR "EASY PEEL CORNER" IS A FEATURE UNIQUE TO BATRIK

PRODUCT SELECTION

| ITEM # | STERILE | ITEM# | NON-STERILE | COLOR | DIMENSIONS | PACKAGING |
|---------|---------|----------|-------------|----------|---|---------------------|
| 11-1002 | | 911-1002 | | Yellow | Standard ø (0.06" x 0.5" / 0.15 cm x 1.27 cm) | 2 pairs / foam base |
| 11-1003 | | 911-1003 | | Yellow | Standard ø (0.06" x 0.5" L / 0.15 cm x 1.27 cm) | 3 pairs / foam base |
| 11-1005 | | 911-1005 | | Yellow | Standard ø (0.06" x 0.5" L / 0.15 cm x 1.27 cm) | 5 pairs / foam base |
| 11-1025 | | 911-1025 | | Assorted | Standard ø (0.06" x 0.5" L / 0.15 cm x 1.27 cm) | 5 pairs / foam base |
| 11-1045 | | 911-1045 | | White | Standard ø (0.06" x 0.5" L / 0.15 cm x 1.27 cm) | 5 pairs / foam base |
| 11-1065 | | 911-1065 | | Blue | Standard ø (0.06" x 0.5" L / 0.15 cm x 1.27 cm) | 5 pairs / foam base |
| 11-1085 | | 911-1085 | | Red | Standard ø (0.06" x 0.5" L / 0.15 cm x 1.27 cm) | 5 pairs / foam base |
| 11-3002 | | 911-3002 | | Pink | Large ø (0.09" x 0.75" L / 0.23 cm x 1.91 cm) | 2 pairs / foam base |
| 11-3003 | | 911-3003 | | Pink | Large ø (0.09" x 0.75" L / 0.23 cm x 1.91 cm) | 3 pairs / foam base |
| 11-3005 | | 911-3005 | | Pink | Large ø (0.09" x 0.75" L / 0.23 cm x 1.91 cm) | 5 pairs / foam base |

STERILE PACKAGING: 1 Foam Base / Pouch – 5 Pouches / Box

NON-STERILE PACKAGING: Sold per foam base

CLAMP JAW COVERS

Manufactured from medical grade silicone, our Clamp Jaw Covers reduce trauma to vessels requiring occlusion, particularly when a hemostatic clamp is utilized. The thread-like indents allow for gentle, constant and non-slip tension.

GOOD TO KNOW: OUR CLAMP JAW COVERS CONVERT ALMOST ANY INSTRUMENT TO ATRAUMATIC



- Gamma sterilized
- Available in 3 colors and 3 widths to fit most clamp sizes
- Radiopaque and latex free

PRODUCT SELECTION

| ITEM # STERILE | ITEM # NON-STERILE | DIMENSIONS | COLOR |
|----------------|--------------------|-------------------------------|-------|
| 12-1000 | 912-1000 | 0.1" X 5" / 0.25 cm X 12.7 cm | Red |
| 12-1001 | 912-1001 | 0.2" X 5" / 0.51 cm X 12.7 cm | White |
| 12-1002 | 912-1002 | 0.3" X 5" / 0.76 cm X 12.7 cm | Blue |

STERILE PACKAGING: 2 clamps per sterile pouch - 10 pouches per box

NON-STERILE PACKAGING: 20 clamps per bag (also available 100 units / bag)


VASCULAR SILICONE TIES

Our Vascular Silicone Ties permit safe and smooth retraction, occlusion and identification of vessels during surgical procedures.

Unlike clamps and retractors which can pinch the vessels flat, or even crack them, the Batrik elliptical shape allows for the gentle handling of arteries, veins, nerves, tendons, spermatic cords, ureters, or whenever delicate retraction / occlusion is required.

- Manufactured of medical grade silicone
- Do not twist on themselves when unpacked
- Easy to separate, eliminating the need for cardboard inserts
- Sterile, radiopaque, and latex free

AVAILABLE IN 4 SIZES AND 4 COLORS



PRODUCT SELECTION

| ITEM # | DIMENSIONS | COLOR |
|----------------|-------------------------------|------------------|
| 10-1000 | 0.075" X 16" (0.2 cm x 40 cm) | Mini-Red |
| 10-1001 | 0.075" X 16" (0.2 cm x 40 cm) | Mini-Blue |
| 10-1002 | 0.075" X 16" (0.2 cm x 40 cm) | Mini-Yellow |
| 10-1003 | 0.075" X 16" (0.2 cm x 40 cm) | Mini-White |
| 10-1010 | 0.120" x 16" (0.3 cm x 40 cm) | Maxi-Red |
| 10-1011 | 0.120" x 16" (0.3 cm x 40 cm) | Maxi-Blue |
| 10-1012 | 0.120" x 16" (0.3 cm x 40 cm) | Maxi-Yellow |
| 10-1013 | 0.120" x 16" (0.3 cm x 40 cm) | Maxi-White |
| 10-1101 | 0.04" x 16" (0.1 cm x 40 cm) | Extra-Mini-Blue |
| 10-1111 | 0.200" x 22" (0.5 cm x 56 cm) | Extra-Maxi-Blue |
| 10-1113 | 0.200" x 22" (0.5 cm x 56 cm) | Extra-Maxi-White |

| ITEM # | DIMENSIONS | COLOR |
|----------------|-------------------------------|-------------|
| 10-3200 | 0.075" x 32" (0.2 cm x 81 cm) | Mini-Red |
| 10-3201 | 0.075" x 32" (0.2 cm x 81 cm) | Mini-Blue |
| 10-3202 | 0.075" x 32" (0.2 cm x 81 cm) | Mini-Yellow |
| 10-3203 | 0.075" x 32" (0.2 cm x 81 cm) | Mini-White |
| 10-3210 | 0.120" x 32" (0.3 cm x 81 cm) | Maxi-Red |
| 10-3211 | 0.120" x 32" (0.3 cm x 81 cm) | Maxi-Blue |
| 10-3212 | 0.120" x 32" (0.3 cm x 81 cm) | Maxi-Yellow |
| 10-3213 | 0.120" x 32" (0.3 cm x 81 cm) | Maxi-White |

PACKAGING: 2 ties per sterile pouch – 10 pouches per box

G.O.L.F.F.™ ANTI-FOG SOLUTION

Our G.O.L.F.F.™ Anti-Fog Solution solves the fogging problem often encountered during Laparoscopic surgical procedures.

- Prevents the occurrence of fog on lenses
- Ensures clearer vision and avoids unnecessary interruptions during surgery
- Set includes an easy-squeeze 6 mL bottle and an adhesive radiopaque sponge for secure placement
- OEM (available)
- Sterile
- Latex free



PRODUCT SELECTION

| ITEM # | DESCRIPTION | PACKAGING | CASE DIMENSIONS |
|-------------------|---|-----------|-----------------|
| G.O.L.F.F. | Anti-Fog Solution & Adhesive Sponge set | 20 / Box | 20" x 16" x 18" |

SKINS™ SURGICAL MULTI-TIP SKIN MARKERS

Our SKINS™ Surgical Multi-Tip Skin Markers are designed with a beveled tip to enable the creation of fine or bold tracing lines by holding the pen at different angles.

- Available with or without 6 labels and a flexible 6" calibrated (in/mm) ruler
- Available in non-sterile for pre-operative use, or sterile for use in the operating room
- Surgical prep resistant
- Made of lasting gentian violet ink that is non-toxic and non-tattooing
- Latex free




PRODUCT SELECTION

| ITEM # | TYPE | DESCRIPTION | PACKAGING |
|----------------|-------------|---|--------------|
| S-3000 | Sterile | Multi-Tip fine marker with ruler and label set | 50 / Box |
| S-3001 | Sterile | Multi-Tip fine marker | 50 / Box |
| S-30005 | Sterile | Multi-Tip fine marker with ruler and label set (Volume Packaging) | 500 / Box |
| S-30015 | Sterile | Multi-Tip fine marker (Volume Packaging) | 500 / Box |
| S-93000 | Non-Sterile | Multi-Tip fine marker with ruler and label set | 1000 / Box |
| S-93001 | Non-Sterile | Multi-Tip fine marker | 1000 / Box |
| S-94000 | Non-Sterile | Ruler | 100 / Bag |
| S-94001 | Non-Sterile | Ruler (Bulk) | 20,000 / Box |
| S-95000 | Non-Sterile | Set of 3 labels | 100 / Bag |

SAVE by ordering our SKINS™ MARKERS in volume packaging. Cost effective and reduces waste.

NEW!




Introducing PandaPatch™: Innovative Adhesive Eye Protection for Safer Patient Care

PandaPatch™ offers an effective alternative to traditional tape for protecting patients' eyes during surgery. These adhesive eye protectors are designed to reduce the risk of corneal injuries and ophthalmic complications for patients under general anesthesia.

With PandaPatch™, you ensure superior patient outcome by preventing avoidable eye damage associated with the use of standard surgical tape.

KEY FEATURES:

- **Ergonomic Design:** Contours comfortably around the eyes.
- **Post-Surgery Friendly:** Particularly effective for children, the elderly, and individuals with sensitive or delicate skin.
- **Latex-Free & Hypoallergenic Adhesive:** Made from medical-grade material, gentle on the skin and designed to minimize irritation.
- **Non-Adhesive Colored Tabs:** Allowing quick and hassle free application and removal, facilitating effortless monitoring of pupil dilation during surgery.



70-04ADU
Adult Sensitive Eye Protector

HERE'S WHY IT MATTERS:

- During prolonged surgery under general anesthesia, **up to 59%¹** of patients may not fully close their eyes, risking corneal abrasions.
- PandaPatch™ gentle adhesive eye protectors ensure eyelids remain safely and comfortably closed, minimizing the risks typically associated with the use of standard surgical tape.



REFERENCES

1. Anson, Jonathan, M.D. "Perioperative Corneal Abrasions: Etiology, Prevention, and Management." Clinical Updates, Penn State University College of Medicine, 2021.




PANDAPATCH™ RANGE:

Three adhesive eye protectors to accommodate different age groups and skin sensitivities. Providing optimal protection for every patient.

- Adult Extra Sensitive Adhesive (Blue Tabs)
- Adult Adhesive (Red Tabs)
- Pediatric Adhesive (Red Tabs)



70-04ADU
Adult Adhesive Eye Protector




70-04PED
Pediatric Adhesive Eye Protector

PRODUCT SELECTION

| ITEM # | DESCRIPTION | DIMENSIONS | PACKAGING |
|-----------------|--|---------------|-----------|
| 70-02ADU | ADULT EXTRA SENSITIVE ADHESIVE EYE PROTECTORS (BLUE) | 2.15" x 1.35" | 50 / Box |
| 70-04ADU | ADULT ADHESIVE EYE PROTECTORS (RED) | 2.15" x 1.35" | 50 / Box |
| 70-04PED | PEDIATRIC ADHESIVE EYE PROTECTORS (RED) | 1.85" x 0.93" | 50 / Box |

MAGNETIC MAT

Our Magnetic Mat provides a sterile, neutral field for the safe placement and hands-free transfer of instruments. It is latex free, reusable, and autoclavable (up to 134° C).



PRODUCT SELECTION


| ITEM # | DESCRIPTION | DIMENSIONS | PACKAGING |
|-----------------|---------------------------------|---------------------------|-----------|
| MAG-02-1 | Magnetic Mat | 12" x 15" / 30 cm x 38 cm | 1 Unit |
| MAG-02-5 | Magnetic Mat (Volume Packaging) | 12" x 15" / 30 cm x 38 cm | 5 Units |

*Independent laboratory tested to validate pre-vacuum steam sterilization efficacy for reprocessing of medical devices in a healthcare setting

BATRIK PIN GUARDS



- Safely cover the sharp ends of your K-Wire and Steinmann Pins
- Color coded by wire size which eliminates the risk of selecting of the wrong size
- A removable sleeve is provided to avoid injury and facilitate application (sterile only)
- Non-sterile model is available with or without plastic sleeves
- Packaging designed to reduce waste
- Latex free



PRODUCT SELECTION

| ITEM # STERILE PIN GUARDS + SLEEVE | ITEM # NON-STERILE PIN GUARDS + SLEEVE | ITEM # NON-STERILE PIN GUARDS (NO SLEEVE) | COLOR | DIMENSIONS (WIRE SIZE) |
|--|--|---|------------|---------------------------|
| 14-1028 | 914-1028 - SLV | 914-1028 | Yellow | 0.028" / 0.71 mm |
| 14-1035 | 914-1035 - SLV | 914-1035 | Light Blue | 0.035" / 0.89 mm |
| 14-1045 | 914-1045 - SLV | 914-1045 | White | 0.045" / 1.14 mm |
| 14-1054 | 914-1054 - SLV | 914-1054 | Dark Red | 0.054" / 1.40 mm |
| 14-1062 | 914-1062 - SLV | 914-1062 | Green | 0.062" / 1.57 mm |
| 14-1078 | 914-1078 - SLV | 914-1078 | Pink | 0.078" / 1.98 mm |
| 14-1093 | 914-1093 - SLV | 914-1093 | Grey | 0.093" / 2.36 mm |
| 14-1109 | 914-1109 - SLV | 914-1109 | Beige | 0.109" / 2.77 mm |
| 14-1125 | 914-1125 - SLV | 914-1125 | Purple | 0.125" / 3.17 mm |
| 14-1141 | 914-1141- SLV | 914-1141 | Brown | 0.141" / 3.58 mm |
| 14-1156 | 914-1156 - SLV | 914-1156 | Dark Blue | 0.156" / 3.96 mm |
| 14-1188 | 914-1188 - SLV | 914-1188 | Red | 0.188" / 4.77 mm |

STERILE PACKAGING: 80 units/box (2 units per sterile pouch - 40 pouches per box)

NON-STERILE PACKAGING: 100 units/box

CAUTERY TIP CLEANER

Our Cautery Tip Cleaner assists in removing material from cautery tips for both mono-polar and bi-polar electrodes.

- Foam pad with integrated abrasive material
- Adheres to disposable and non-disposable drapes
- Radiopaque
- Available exclusively in Non-Sterile
- Latex free



PRODUCT SELECTION

| ITEM # | DESCRIPTION | DIMENSIONS | PACKAGING |
|----------|------------------------------------|-------------------------------|------------|
| 60-10000 | Cautery Tip Cleaner Non-Sterile | 1.8" x 1.8" / 4.5 cm x 4.5 cm | 2,200 /Box |

SAFE-SOAK™ SPAC

Our Safe-Soak™ SPAC floor suction disk quickly removes fluid waste from the operating room floor through the application of suction.



- Highly durable: it won't be damaged if stepped on
- Easily maneuverable by foot around operating room floor
- Latex free, disposable, and non-sterile
- The bottom opening permits fluid collection through a highly absorbent material that comes into contact with its surface
- The hose is securely affixed on its side to reduce risk of catching equipment
- Its highly absorbent base material allows for quick fluid wicking and removal when suction is applied

BATRIK SAFE-SOAK™ ABSORBENT FLOOR MATS

Our Safe-Soak™ Absorbent Floor Mats are designed to absorb a variety of fluids from the operating room floor quickly and efficiently, ensuring the floor beneath remains slip-free, safe, and dry during surgical procedures.

- Highly absorbent
- Facilitates quick and easy clean up
- Latex free, single-use, disposable, and non-sterile
- Conforms to uneven surfaces
- Skid-resistant backing (blue absorbent mat only)



PRODUCT SELECTION

| ITEM # | DESCRIPTION | ABSORBENCY | SIZE | PACKAGING |
|-------------------|--|-------------|----------------------------|-----------|
| SS-20-1000 | Safe-Soak™ / White extra absorbent mat | 4.5 L / mat | 28" x 48" (71 cm x 122 cm) | 20 / Box |
| SS-20-2000 | Safe-Soak™ / Blue absorbent mat with skid-resistant backer | 2.5 L / mat | 28" x 40" (71 cm x 102 cm) | 20 / Box |
| SS-20-3000 | Safe-Soak™ / Yellow absorbent mat | 3 L / mat | 30" x 52" (76 cm x 132 cm) | 15 / Box |
| SS-21-1000 | Safe-Soak™ SPAC | N/A | 8.625" (22 cm diameter) | 12 / Box |

BLUE CLOUD-GEL™ ANTI-FATIGUE MATS


Our Blue Cloud-Gel™ Anti-Fatigue Mats have been scientifically engineered for use in the O.R. using patented gel technology and premium foam to provide the ultimate comfort and support for surgeons, nurses, and technicians particularly when movement is restricted for prolonged durations.

- Proprietary 3/4" Blue Cloud-Gel™ foam technology will not "bottom-out" like conventional foam mats
- Conforms to and supports feet promoting efficient circulation, reducing physical discomfort and fatigue
- Promotes proper posture while diminishing back, knee, and foot discomfort
- Latex free, water resistant, and nonabsorbent
- Reusable, easy to clean, and Cart Washer compatible (contact us for more details)
- Excellent friction properties when damp or dry surgical booties are worn
- Completely sealed with an RF weld for excellent protection against moisture absorption
- US Standards: CALL 117 and NFPA 260 Compliant
- 1 Year Warranty

DID YOU KNOW? STUDIES HAVE DEMONSTRATED THAT ANTI-FATIGUE MATS IMPROVE PERFORMANCE, ENERGY, AND POST-OPERATIVE DISCOMFORT.

PRODUCT SELECTION

| ITEM # | DIMENSIONS | PACKAGING |
|-----------------|---|-----------|
| 800-2401 | 18" X 24" / 46 cm x 61 cm | 1 / Box |
| 800-2405 | 18" X 24" / 46 cm x 61 cm | 5 / Box |
| 800-3201 | 20" X 32" / 51 cm x 81 cm | 1 / Box |
| 800-3203 | 20" X 32" / 51 cm x 81 cm | 3 / Box |
| 800-4801 | 20" X 48" / 51 cm x 122 cm | 1 / Box |
| 800-4802 | 20" X 48" / 51 cm x 122 cm | 2 / Box |
| 800-7201 | 30" X 72" / 76 cm x 183 cm (Special Order) | 1 / Box |




BLUE-CLOUD NIMBUS™ ANTI-FATIGUE MATS

Our Blue-Cloud Nimbus™ Anti-Fatigue Mats are designed for general hospital use. Made with proprietary high-density “energy-return” polyurethane materials, they provide unmatched comfort, relief, and ergonomic support.

- Our Energy-Return properties will help reduce leg, foot, lower back discomfort, and work-related injuries
- Longer lasting support as 3/4” proprietary foam blend will not bottom out like most other foam based anti-fatigue mats
- Antimicrobial properties blended into the mat provides added protection
- Highly water resistant surface will not absorb water or dirt
- Lightweight and easy to clean
- Textured design on surface and high-traction base ensures greater safety. (US Standard NFSI 101-C)
- Beveled edging reduces tripping hazards and anti-slip backing enhances safety
- Durable design is ideal for commercial applications, yet still ergonomically and aesthetically pleasing for home or office environments
- Flammability/Burn resistant per CFR 16 Parts 1630/1631
- Made with eco-friendly materials, is non-toxic, and latex free
- 5 Year Warranty (non-industrial usage), 4 Year Warranty (industrial usage)

GOOD TO KNOW: IF WITHIN 30 DAYS YOU ARE NOT 100% SATISFIED, CONTACT YOUR DEALER FOR A REFUND (EXCLUDES TRANSPORT, HANDLING, AND TAXES)



Available in Taupe or Black

PRODUCT SELECTION

| ITEM # | DIMENSIONS INCHES / CM | PACKAGING |
|--------------------|----------------------------|-----------|
| BLACK COLOR | | |
| 809-3201 | 20" X 32" / 51 cm x 81 cm | 1 / Box |
| 809-3204 | 20" X 32" / 51 cm x 81 cm | 4 / Box |
| 809-4801 | 20" X 48" / 51 cm x 122 cm | 1 / Box |
| 809-4804 | 20" X 48" / 51 cm x 122 cm | 4 / Box |
| 809-3601 | 24" X 36" / 61 cm x 91 cm | 1 / Box |
| 809-7201 | 20" X 72" / 51 cm x 183 cm | 1 / Box |
| 809-6001 | 36" X 60" / 91 cm x 152 cm | 1 / Box |
| TAUPE COLOR | | |
| 810-3201 | 20" X 32" / 51 cm x 81 cm | 1 / Box |
| 810-3204 | 20" X 32" / 51 cm x 81 cm | 4 / Box |
| 810-4801 | 20" X 48" / 51 cm x 122 cm | 1 / Box |
| 810-4804 | 20" X 48" / 51 cm x 122 cm | 4 / Box |
| 810-3601 | 24" X 36" / 61 cm x 91 cm | 1 / Box |
| 810-7201 | 20" X 72" / 51 cm x 183 cm | 1 / Box |
| 810-6001 | 36" X 60" / 91 cm x 152 cm | 1 / Box |

STACKABLE STEP STOOLS

Made of lightweight cast aluminium and protected with an antimicrobial powder coating, our patented Stackable Step Stools have been designed with infection prevention as a prime focus.



PRODUCT SELECTION

| ITEM # | DESCRIPTION | COLOR | CASE DIMENSIONS |
|-------------|----------------------------------|-----------|-----------------|
| OR-36363 | Stackable Step Stool with mat | Black | 19" x 14" x 6" |
| OR-36363-02 | Stackable Step Stool with mat | Dark Blue | 19" x 14" x 6" |
| OR-36363-03 | Stackable Step Stool with mat | Green | 19" x 14" x 6" |
| OR-36363-07 | Stackable Step Stool with mat | Orange | 19" x 14" x 6" |
| OR-36364 | Stackable Step Stool without mat | N/A | 19" x 14" x 6" |

PACKAGING: Sold by Unit

REPLACEMENT PARTS

| ITEM # | DESCRIPTION | PACKAGING | CASE DIMENSIONS |
|-------------|------------------------------|-----------|-----------------------------------|
| OR-36363-1A | Replacement mat for OR-36363 | 1 Unit | 15" x 18" x 5" (Ziplock) |
| OR-36363-1B | Replacement Foot Pegs | 1 Unit | Ziplock Based on Quantity ordered |

OPTIONAL ACCESSORIES

| ITEM # | DESCRIPTION | PACKAGING | CASE DIMENSIONS |
|----------|--|-----------|-----------------|
| OR-3693 | Transport cart | 1 Unit | 34" x 17" x 7" |
| 809-1705 | Reusable Anti-Fatigue Mat with adhesive strips to secure it in place | 5 / Box | 12" x 12" x 4" |
| 899-1710 | Single-Use Anti-Fatigue Mat with adhesive strips to secure it in place | 10 / Box | 12" x 12" x 8" |

PRODUCT FEATURES

- Weight: 8.3 lbs (3.8 kg)
- Dimensions: 17.4" L x 13" W x 4.7" H (44.2 cm L x 33 cm W x 11.9 cm H)
- Transport Cart Available
- Weight Capacity: 300 lb (136 kg)
- Warranty: 7 years (step stool only)

ANTIMICROBIAL COATING:

Prevents the spread of hazardous microorganisms and suppresses the growth of mold, bacteria, yeast, and mildew

ANTI-SLIP SURFACE:

Provides enhanced stability.

RUBBER-FREE CORNERS:

Avoids fluid collection.

WASHABLE CART:

Model OR-36364 (without mat) can be sterilized.

HYGIENIC DESIGN:


No welds or large hollow openings to trap debris or fluids



OR-3693

TOURNY™ + TOURNIQUETS

Our Tourniquets are comfortable and easy to use on all patients, including the elderly, undernourished, children, and newborns.



AVAILABLE COLORS



- Available in different sizes including compact sizes, which are convenient for phlebotomy trays
- Innovative latex free material has a more comfortable texture for patients, puts less stress on skin, and safeguards hospital staff and patients from latex allergens

SINGLE-USE BENEFITS

- Reduced risk of cross-contamination (patient to patient or patient to healthcare provider)
- Reduced risk of spread of nosocomial infections
- May improve your facility's compliance with OSHA's 1910-1030 (d)(1) "Blood-borne Pathogens Standard"



FOR YOUR CONVENIENCE:
25 UNITS PER ROLL
COME IN A DISPENSER BOX
FOR EASY ACCESS



PRODUCT SELECTION


| ITEM # | DESCRIPTION | COLOR | DIMENSIONS | PACKAGING |
|-------------------|--------------------------------|---|-------------------------------|-----------------------|
| K-002 | TOURNY™+ Non-Sterile | Blue | 1" x 18" / 2.5 cm x 45 cm | 1 / Bag |
| K-004 | TOURNY™+ Non-Sterile | Blue | 1" x 18" / 2.5 cm x 45 cm | 100 / Box |
| K-009 | TOURNY™+ Non-Sterile | Blue | 0.200" x 22" / 0.5 cm x 56 cm | 20 / Box |
| K-010 | TOURNY™+ Non-Sterile | White | 0.200" x 22" / 0.5 cm x 56 cm | 20 / Box |
| K-400 | TOURNY™+ Non-Sterile | Blue | 1" x 18" / 2.5 cm x 45 cm | 16 rolls 400 / Box |
| K-400PB | TOURNY™+ Non-Sterile (Pebbled) | Light Blue | 1" x 18" / 2.5 cm x 45 cm | 16 rolls 400 / Box |
| K-100 | TOURNY™+ Non-Sterile | Assorted (3 rolls of each) Blue, Orange, Green, and Pink | 1" x 18" / 2.5 cm x 45 cm | 12 rolls 300 / Box |
| K-254 | TOURNY™+ Non-Sterile (Pebbled) | Light Blue | 1" x 18" / 2.5 cm x 45 cm | 250 / Box |
| K-2-18RL | TOURNY™+ Non-Sterile | Blue | 1" x 18" / 2.5 cm x 45 cm | 1 Roll |
| K-2-18PBRL | TOURNY™+ Non-Sterile (Pebbled) | Light Blue | 1" x 18" / 2.5 cm x 45 cm | 1 Roll |
| K-3-18RL | TOURNY™+ Non-Sterile | Green | 1" x 18" / 2.5 cm x 45 cm | 1 Roll |
| K-6-18RL | TOURNY™+ Non-Sterile | Pink | 1" x 18" / 2.5 cm x 45 cm | 1 Roll |
| K-7-18RL | TOURNY™+ Non-Sterile | Orange | 1" x 18" / 2.5 cm x 45 cm | 1 Roll |
| K-2-29RL | Long TOURNY™+ Non-Sterile | Blue | 1" x 29.5" / 2.5 cm x 75 cm | 1 Roll |

SURGI-KIT®

LATEX FREE | STERILE | SINGLE USE | DISPOSABLE

Ethox

- Disposable pockets for sterile delivery to OR
- Keeps laparoscopic instruments close to operating field
- Maintains sterile operating field while keeping instruments within arm's reach
- Reduces time required for locating instruments
- Eliminates cleaning procedures
- Multi-pocket sleeves to accommodate a variety of instrument sizes
- Available with peel-n-stick tape to adhere to operating drape




8 Pocket
SK-100
SK-100N




4 Pocket
SK-150
SK-150N




2 Pocket
SK-200
SK-200N




8 Pocket + Tape
SK-300
SK-300N



4 Pocket + Tape
SK-350
SK-350N



2 Pocket + Tape
SK-400
SK-400N



10 Pocket + Tape
SK-650
SK-650N

STERILE

| ITEM# | DESCRIPTION | SIZE | PK |
|--------|-------------------|-------------|-----|
| SK-100 | 8 pockets | 19cm x 23cm | 50 |
| SK-150 | 4 pockets | 19cm x 11cm | 100 |
| SK-200 | 2 pockets | 28cm x 18cm | 50 |
| SK-300 | 8 pockets + tape | 19cm x 23cm | 50 |
| SK-350 | 4 pockets + tape | 19cm x 11cm | 75 |
| SK-400 | 2 pockets + tape | 28cm x 18cm | 50 |
| SK-650 | 10 pockets + tape | 5cm x 19cm | 50 |

NON-STERILE

| ITEM# | DESCRIPTION | SIZE | PK |
|---------|-------------------|-------------|-----|
| SK-100N | 8 pockets | 19cm x 23cm | 50 |
| SK-150N | 4 pockets | 19cm x 11cm | 100 |
| SK-200N | 2 pockets | 28cm x 18cm | 50 |
| SK-300N | 8 pockets + tape | 19cm x 23cm | 50 |
| SK-350N | 4 pockets + tape | 19cm x 11cm | 75 |
| SK-400N | 2 pockets + tape | 28cm x 18cm | 50 |
| SK-650N | 10 pockets + tape | 5cm x 19cm | 50 |

* Non-sterile kits available only to kit packers.

Instrument Holder


LAP-KIT®

LATEX FREE | STERILE | SINGLE USE | DISPOSABLE


Ethox




- Disposable pockets for sterile delivery to OR
- Keeps laparoscopic instruments close to operating field
- Maintains sterile operating field while keeping instruments within arm's reach
- Reduces time required for locating instruments
- Eliminates cleaning procedures
- Multi-pocket sleeves to accommodate a variety of instrument sizes
- Available with peel-n-stick tape to adhere to operating drape



2 Pockets
LK-700
LK-700N



4 Pockets
LK-800



3 Pockets
LK-900
LK-900N

STERILE

| ITEM | DESCRIPTION | SIZE | PK |
|--------|------------------|---|----|
| LK-700 | 2 Pockets + Tape | 6.5 cm x 30.5 cm | 25 |
| LK-800 | 4 Pockets + Tape | 3.2 cm x 30.5 cm | 25 |
| LK-900 | 3 Pockets + Tape | 3.2 cm x 30.5 cm (2), 6.5 cm x 30.5 cm | 25 |

NON-STERILE

| ITEM | DESCRIPTION | SIZE | PK |
|---------|------------------|---|----|
| LK-700N | 2 Pockets + Tape | 6.5 cm x 30.5 cm | 25 |
| LK-900N | 3 Pockets + Tape | 3.2 cm x 30.5 cm (2), 6.5 cm x 30.5 cm | 25 |

* Non-sterile kits available only to kit packers.

Warm Gel®

Get relief where you need it when you need it



Features and Benefits

- Pliable WarmGel® conforms to any body part
- Easy to activate
- Soft non-woven fabric for patient comfort
- Consistent temperature
- Available in two convenient sizes

WarmGel® is a versatile and reliable source for instant heat

Many applications in:

- Maternity
- IV Therapy
- Labor & Delivery
- Rehabilitation

Temperature Range:

WarmGel® has been designed to achieve 112°F (44°C) when activated at room temperature of 75 °F (24 °C)

| Item # | Description |
|--------|---|
| 20419 | Medium 6" x 6" Adult WarmGel® Instant Heat Pack - 40/Case |
| 20420 | Large 7" x 9" Adult WarmGel® Instant Heat Pack - 24/Case |

K-Wires & Steinmann Pins

Packaging: 6 or 10 units per plastic tube / 1 tube per package

K-Wires


| | DOUBLE TROCAR | SINGLE TROCAR | | | | |
|--------------------|-----------------------|-----------------------|--------------------|--------------------|----------------------|----------------------|
| | Item # | Item # | Length (in) | Length (cm) | Diameter (in) | Diameter (mm) |
| SMOOTH | SU78.2000 | - | 4" | 10.1 cm | .031" | 0.8 mm |
| | SU78.2010 | - | 4" | 10.1 cm | .039" | 1.0 mm |
| | SU78.2020 | - | 4" | 10.1 cm | .047" | 1.2 mm |
| | SU78.2040 | - | 4" | 10.1 cm | .054" | 1.4 mm |
| | SU78.2030 | - | 4" | 10.1 cm | .062" | 1.6 mm |
| | SU78.1210 | - | 6" | 15.2 cm | .031" | 0.8 mm |
| | SU78.1220 | SU78.2810 | 6" | 15.2 cm | .039" | 1.0 mm |
| | SU78.1230 | SU78.2820 | 6" | 15.2 cm | .047" | 1.2 mm |
| | SU78.1240 | - | 6" | 15.2 cm | .054" | 1.4 mm |
| | SU78.1250 | SU78.2850 | 6" | 15.2 cm | .062" | 1.6 mm |
| | SU78.2105 | SU78.2500 | 9" | 22.8 cm | .031" | 0.8 mm |
| | SU78.2110 | SU78.2510 | 9" | 22.8 cm | .039" | 1.0 mm |
| | SU78.2120 | SU78.2520 | 9" | 22.8 cm | .047" | 1.2 mm |
| | SU78.2140 | - | 9" | 22.8 cm | .054" | 1.4 mm |
| | SU78.2130 | SU78.2530 | 9" | 22.8 cm | .062" | 1.6 mm |
| | FULL THREAD | SU78.4030 | - | 9" | 22.8 cm | .062" |
| SMOOTH | DOUBLE DIAMOND | SINGLE DIAMOND | | | | |
| | Item # | Item # | Length (in) | Length (cm) | Diameter (in) | Diameter (mm) |
| | - | SU78.2810 | 6" | 15.2 cm | .039" | 1.0 mm |
| | - | SU78.2820 | 6" | 15.2 cm | .047" | 1.2 mm |
| | - | SU78.2850 | 6" | 15.2 cm | .062" | 1.6 mm |
| | SU78.3110 | SU78.3410 | 9" | 22.8 cm | .039" | 1.0 mm |
| | SU78.3120 | SU78.3420 | 9" | 22.8 cm | .047" | 1.2 mm |
| SU78.3130 | SU78.3430 | 9" | 22.8 cm | .062" | 1.6 mm | |
| FULL THREAD | - | SU78.4070 | 9" | 22.8 cm | .062" | 1.6 mm |

K-Wires & Steinmann Pins

Steinmann Pins

| DOUBLE TROCAR | | SINGLE TROCAR | | | | | |
|---------------|----------------|---------------|----------------|---------------|---------------|---------------|--|
| Item # | Item # | Length (in) | Length (cm) | Diameter (in) | Diameter (mm) | | |
| - | SU78.6525 | 6" | 15.2 cm | .098" | 2.5 mm | | |
| - | SU78.6429 | 6" | 15.2 cm | .110" | 2.8 mm | | |
| SMOOTH | SU78.5500 | SU78.6100 | 9" | 22.8 cm | .078" | 2.0 mm | |
| | SU78.5530 | SU78.6130 | 9" | 22.8 cm | .098" | 2.5 mm | |
| | SU78.5560 | SU78.6160 | 9" | 22.8 cm | .110" | 2.8 mm | |
| | SU78.5590 | SU78.6190 | 9" | 22.8 cm | .118" | 3.0 mm | |
| | SU78.5620 | SU78.6220 | 9" | 22.8 cm | .141" | 3.5 mm | |
| | SU78.5650 | SU78.6250 | 9" | 22.8 cm | .157" | 4.0 mm | |
| | SU78.5680 | SU78.6280 | 9" | 22.8 cm | .177" | 4.5 mm | |
| | SU78.5690 | SU78.6290 | 9" | 22.8 cm | .197" | 5.0 mm | |
| FULL THREAD | SU78.8500 | - | 9" | 22.8 cm | .078" | 2.0 mm | |
| | SU78.8530 | - | 9" | 22.8 cm | .098" | 2.5 mm | |
| | SU78.8560 | - | 9" | 22.8 cm | .110" | 2.8 mm | |
| | SU78.8590 | - | 9" | 22.8 cm | .118" | 3.0 mm | |
| | SU78.8620 | - | 9" | 22.8 cm | .141" | 3.5 mm | |
| | SU78.8650 | - | 9" | 22.8 cm | .157" | 4.0 mm | |
| SMOOTH | DOUBLE DIAMOND | | SINGLE DIAMOND | | | | |
| | Item # | Item # | Length (in) | Length (cm) | Diameter (in) | Diameter (mm) | |
| | SU78.7000 | SU78.7780 | 9" | 22.8 cm | .078" | 2.0 mm | |
| | SU78.7060 | SU78.7660 | 9" | 22.8 cm | .110" | 2.8 mm | |
| | - | SU78.7690 | 9" | 22.8 cm | .118" | 3.0 mm | |
| | - | SU78.7720 | 9" | 22.8 cm | .141" | 3.5 mm | |
| | - | SU78.7782 | 9" | 22.8 cm | .177" | 4.5 mm | |
| FULL THREAD | SU78.8300 | - | 9" | 22.8 cm | .078" | 2.0 mm | |
| | - | SU78.8030 | 9" | 22.8 cm | .098" | 2.5 mm | |
| | - | SU78.8120 | 9" | 22.8 cm | .141" | 3.5 mm | |
| 3 SHANK END | - | SU78.7630 | 9" | 22.8 cm | .110" | 2.8 mm | |

K-Wires, Steinmann Pins & Sterilization Racks



Sterilization Racks

Sterilization racks are designed for use during washing, sterilizing, storing/transporting, and dispensing K-wires, Steinmann Pins, or Twist Drills. The racks are lightweight but made to withstand repeated use. They fold closed for sterilization and easily open to convert to a table top stand for use in the O.R. Dimensions are clearly marked on the rack to assist with quick identification.



K-Wire and Pin Rack

| Number | Details | Dimensions | Quantity / Package |
|--------------|---|---|--------------------|
| SU98.5409-RO | Anodized aluminium, holds 9" K-wires and pins Holds 6 of each of the following diameters: 0.028" (0.7 mm) 0.045" (1.1 mm) 0.062" (1.6 mm) 0.093" (2.4 mm) 0.125" (3.2 mm) 0.156" (4.0 mm) 0.035" (0.9 mm) 0.054" (1.4 mm) 0.078" (2.0 mm) 0.109" (2.8 mm) 0.140" (3.5 mm) 0.177" (4.5 mm) | 12.5" x 5.25" x 1.5" (32 x 13 x 3.8) | 1 |



K-Wire Rack

| Number | Details | Dimensions | Quantity / Package |
|--------------|---|---|--------------------|
| SU98.5404-RO | Anodized aluminium, holds 4 - 6" K-wires and pins Holds 6 of each of the following diameters: 0.028" (0.7 mm) 0.045" (1.1 mm) 0.062" (1.6 mm) 0.035" (0.9 mm) 0.054" (1.4 mm) | 12.5" x 5.25" x 1.5" (32 x 13 x 3.8) | 1 |


Ver.2021-07-21

ENDOSCOPY BO2™ BITE BLOCKS

WITH NASAL OXYGEN DELIVERY

Our BO2™ Bite Blocks significantly enhance patient safety by providing supplemental oxygenation during gastrointestinal endoscopy, preventing hypoxia in sedated patients, and ensuring superior patient comfort.

- Saves time and money by eliminating the need for separate nasal prongs and oxygen masks
- Designed to be free of sharp edges to ensure patient comfort and safety
- Available exclusively in Single-Use
- Optional retention strap is available to ensure a secure fit
- The Bite Blocks and optional retention straps are all latex free, lead free and DEHP free



PRODUCT SELECTION

| ITEM # | TYPE | DESCRIPTION | PACKAGING | CASE DIMENSIONS |
|--------------|------------|-----------------------------|-----------|-----------------|
| BO2 | Single-Use | 18mm / 54FR, includes strap | 100 / Box | 10" x 10" x 10" |
| BO2-N | Single-Use | 18mm / 54FR, without strap | 100 / Box | 10" x 10" x 10" |

ENDOSCOPY BITESTRONG™ BITE BLOCKS

ADULT & PEDIATRIC BITE BLOCKS

Our Bitestrong™ Bite Blocks provide a safe passage through the patient's mouth during TEE, protecting the patient's teeth and endoscope.

Our molding process ensures a smooth finish with no sharp edges (flashing) that can often irritate a patient's mouth. All Bite Blocks are individually wrapped in sealed, easy to open PE bags and are latex free, lead free, and DEHP free. Additional straps are available, if required

SINGLE USE

- Available without dental retention

REUSABLE

- Designed for repeated usage and will withstand commonly practiced chemical and steam sterilization cycles (134°C / 275°F).
- Available with or without an adjustable retention strap

PRODUCT SELECTION

| ITEM # | TYPE | DESCRIPTION | PACKAGING | CASE DIMENSIONS |
|-----------------|---|---|-----------|------------------|
| B1 | Single-Use Pediatric | 16 mm / 48FR, includes strap | 100 / Box | 10" x 10" x 10" |
| B1-25 | Single-Use Pediatric | 16 mm / 48FR, includes strap | 25 / Box | 10" x 7" x 7" |
| B1-N | Single-Use Pediatric | 16 mm / 48FR, without strap | 100 / Box | 10" x 10" x 10" |
| B3 | Single-Use | 20 mm / 60FR, includes strap | 100 / Box | 10" x 10" x 10" |
| B3-N | Single-Use | 20 mm / 60FR, without strap | 100 / Box | 10" x 10" x 10" |
| B3-SA | Single-Use | 20 mm / 60FR, with strap attached to one side | 100 / Box | 10" x 10" x 10" |
| B3-NDR | Single-Use | 20 mm / 60FR, no dental retention rim | 100 / Box | 10" x 10" x 10" |
| B3-NDR-N | Single-Use | Same as B3-N, no dental retention rim, without strap. | 100 / Box | 10" x 10" x 10" |
| B4 | Single-Use | 19 mm / 60FR | 100 / Box | 10" x 10" x 10" |
| B4T | Single-Use | 19 mm / 60FR, includes strap | 100 / Box | 10" x 10" x 10" |
| B4T-N | Single-Use | 19 mm / 60FR, without strap | 100 / Box | 10" x 10" x 10" |
| B6 | Reusable | 19 mm / 60FR | 100 / Box | 10" x 10" x 10" |
| B6T | Reusable | 19 mm / 60FR, includes strap | 100 / Box | 10" x 10" x 10" |
| B6T-N | Reusable | 19 mm / 60FR, without strap | 100 / Box | 10" x 10" x 10" |
| B6-25 | Reusable | 19 mm / 60FR | 25 / Box | 10" x 7" x 7" |
| B6T-25 | Reusable | 19 mm / 60FR, includes strap | 25 / Box | 10" x 7" x 7" |
| B6T-N-25 | Reusable | 19 mm / 60FR, without strap | 25 / Box | 10" x 7" x 7" |
| BS | Single-Use | Bite-Block straps 21.7" x 0.79" / 55 cm x 2 cm | 100 / Box | 10" x 13" x 1.5" |



BS



B1



B3



B4



B6

COUVRETTE™ VALVE COVERS

Our Couvrette™ Valve Covers eliminate the need for manual cleaning and reprocessing of reusable valves. Our single use Valve Covers effectively reduce cross contamination between patients and also reduce fluid leakage from the parts during procedures.

- Easily and securely fit to all valves
- Reduce the potential for errors
- Single-Use
- Non-Sterile
- Latex free

TIP: REUSABLE VALVE COVERS ARE NOTORIOUSLY DIFFICULT TO CLEAN

SAVE TIME AND REDUCE RISK WITH COUVRETTE™



16-0002

NOW AVAILABLE!



16-0022



16-0044



PRODUCT SELECTION

| ITEM # | DESCRIPTION | COLOR | PACKAGING |
|---------|--|------------|------------------------------------|
| 16-0002 | Biopsy Valve Cover - Olympus® | Blue | 100 / Box |
| 16-0022 | Biopsy Valve Cover - Olympus® & Fujinon® | Light Blue | 100 / Box NOW AVAILABLE! |
| 16-0044 | Biopsy Valve Cover - Pentax® | Pink | 100 / Box NOW AVAILABLE! |


NOTE: Contact us to verify compatibility with your scope

DAYBRUSH™ CHANNEL BRUSHES

REUSABLE

Our Reusable DayBrush™ Channel Brushes can be used for cleaning endoscopes for an entire day. In fact, it's been proven to effectively clean up to 10 scopes before the need to discard.*

- Soft tip ensures easy manipulation within the channel with no risk of damage
- The back-end is sealed and softened to prevent fluid from getting trapped inside the tube. The softening process also prevents damage if the cleaning brush is used in reverse.
- Can be steam sterilized.
- Bristle volume is 20mm in length – nearly double most brushes on the market. A large brush provides more surface contact during the cleaning process.
- Significant per use savings when compared to other reusable or disposable channel cleaning brushes.



() Discard after cleaning 10 scopes or sooner if any brush degradation is noticed. After ten applications, the brush may have the appearance that it can continue to clean effectively, but we discourage prolonged use due to risk of cross contamination. Actual number of uses will vary based on handling and processing.*

PRODUCT SELECTION

5 MM CYLINDRICAL BRUSHES – SHEATH Ø 1.75 MM

| ITEM # | COLOUR | LENGTH | BRUSH LENGTH | BRUSH (Ø) | SHEATH (Ø) | BRUSH ENDS |
|--------------------|--------|--------------|---------------|--------------|-----------------|------------|
| 30-18050 | White | 70" / 180 cm | 0.79" / 20 mm | 0.20" / 5 mm | 0.07" / 1.75 mm | Single |
| 30-18050-T | Yellow | 70" / 180 cm | 0.79" / 20 mm | 0.20" / 5 mm | 0.07" / 1.75 mm | Single |
| 30-18050-2T | Yellow | 70" / 180 cm | 0.79" / 20 mm | 0.20" / 5 mm | 0.07" / 1.75 mm | Double |
| 30-23050 | White | 90" / 230 cm | 0.79" / 20 mm | 0.20" / 5 mm | 0.07" / 1.75 mm | Single |
| 30-23050-T | Yellow | 90" / 230 cm | 0.79" / 20 mm | 0.20" / 5 mm | 0.07" / 1.75 mm | Single |
| 30-23050-2T | Yellow | 90" / 230 cm | 0.79" / 20 mm | 0.20" / 5 mm | 0.07" / 1.75 mm | Double |

PACKAGING : 50 / BOX

OPTIONAL KIT ADD-ONS

ADD THE SUFFIX:

"F" to receive a kit with our Flat Sponge (Packaging: Box of 40). Applicable on certain models only, call for details.

"K" to receive a kit with the single use Valve / Control Cleaning Brush 5/11 Double Brushes (item # 27-57787). Applicable on all models.

"FK" to receive both the sponge and valve cleaning brushes.

ARC™ CHANNEL BRUSHES

SINGLE - USE

BATRIK offers one of the largest selection of Endoscopy Channel Cleaning Brushes in the industry. Our brushes are available in a variety of lengths, materials, and styles to accommodate the widest range of endoscopic devices used in the healthcare industry today.

No matter what your requirements, you'll be able to find the brush that suits your specific needs. If you don't see the brush you require, simply call us. We have over 50 different brushes available and can custom fabricate to your precise specifications.

A fan tip protects the biopsy channel from damage, vs. metal-ball tips that can scratch or completely dislodge in the channel

At least 25%-75% longer brush-heads than most competitors

Latex free

Higher bristle-density and longer brush provides increased surface contact for more effective cleaning

The sheath is attached to the brush with proprietary welding (no glue) and guaranteed not to separate

No handle at the end of the brush so that it passes completely through the channel for superior cleaning

Available in 2 different materials:
Polypropylene: The industry standard offering excellent flexibility and low cost.
Delrin: More rigid for easier insertion and handling, saving time and resources . Please note: all small-diameter brushes are offered exclusively with Delrin).

OPTIONAL KIT ADD-ONS



ADD THE SUFFIX:

"F" to receive a kit with our Flat Sponge (Packaging: Box of 40). Applicable on certain models only, call for details.

"K" to receive a kit with the single use Valve / Control Cleaning Brush 5/11 Double Brushes (item # 27-57787). Applicable on all models.

"FK" to receive both the sponge and valve cleaning brushes.

SERIES 311 & 333



0.94" / 24 mm

PRODUCT SELECTION

5 MM CYLINDRICAL BRUSHES – SHEATH Ø 1.85 MM - PUNCHY PINK POLYPROPYLENE

| ITEM # | LENGTH | BRUSH LENGTH | BRUSH (Ø) | SHEATH (Ø) | BRUSH ENDS |
|--------------------|--------------|---------------|--------------|------------------|------------|
| 311-10050 | 39" / 100 cm | 0.94" / 24 mm | 0.20" / 5 mm | 0.073" / 1.85 mm | Single |
| 311-18050 | 70" / 180 cm | 0.94" / 24 mm | 0.20" / 5 mm | 0.073" / 1.85 mm | Single |
| 311-23050 | 90" / 230 cm | 0.94" / 24 mm | 0.20" / 5 mm | 0.073" / 1.85 mm | Single |
| 311-23050-2 | 90" / 230 cm | 0.94" / 24 mm | 0.20" / 5 mm | 0.073" / 1.85 mm | Double |

PACKAGING: 50 / BOX

5/6 MM TAPERED BRUSHES – SHEATH Ø 1.85 MM - PUNCHY PINK POLYPROPYLENE

| ITEM # | LENGTH | BRUSH LENGTH | BRUSH (Ø) | SHEATH (Ø) | BRUSH ENDS |
|--------------------|--------------|---------------|------------------------|------------------|------------|
| 333-23056 | 90" / 230 cm | 0.94" / 24 mm | 0.20", 0.24" / 5, 6 mm | 0.073" / 1.85 mm | Single |
| 333-23056-2 | 90" / 230 cm | 0.94" / 24 mm | 0.20", 0.24" / 5, 6 mm | 0.073" / 1.85 mm | Double |

PACKAGING: 50 / Box

SERIES 300 & 322



0.79" / 20 mm

PRODUCT SELECTION

5 MM CYLINDRICAL BRUSHES – SHEATH Ø 1.75 MM - METALLIC GREY DELRIN

| ITEM # | LENGTH | BRUSH LENGTH | BRUSH (Ø) | SHEATH (Ø) | BRUSH ENDS |
|--------------------|--------------|---------------|--------------|-----------------|------------|
| 300-18050 | 70" / 180 cm | 0.79" / 20 mm | 0.20" / 5 mm | 0.07" / 1.75 mm | Single |
| 300-23050 | 90" / 230 cm | 0.79" / 20 mm | 0.20" / 5 mm | 0.07" / 1.75 mm | Single |
| 300-23050-2 | 90" / 230 cm | 0.79" / 20 mm | 0.20" / 5 mm | 0.07" / 1.75 mm | Double |


PACKAGING: 50 / Box

5/6 MM TAPERED BRUSHES – SHEATH Ø 1.75 MM - METALLIC GREY DELRIN

| ITEM # | LENGTH | BRUSH LENGTH | BRUSH (Ø) | SHEATH (Ø) | BRUSH ENDS |
|--------------------|--------------|---------------|------------------------|-----------------|------------|
| 322-23056 | 90" / 230 cm | 0.79" / 20 mm | 0.20", 0.24" / 5, 6 mm | 0.07" / 1.75 mm | Single |
| 322-23056-2 | 90" / 230 cm | 0.79" / 20 mm | 0.20", 0.24" / 5, 6 mm | 0.07" / 1.75 mm | Double |

PACKAGING: 50 / Box

SERIES 344




PRODUCT SELECTION

3 MM CYLINDRICAL BRUSHES – SHEATH Ø 1.55 MM - METALLIC PURPLE DELRIN

| ITEM # | LENGTH | BRUSH LENGTH | BRUSH (Ø) | SHEATH (Ø) | BRUSH ENDS |
|------------------|--------------|---------------|--------------|-----------------|------------|
| 344-10030 | 39" / 100 cm | 0.79" / 20 mm | 0.12" / 3 mm | 0.06" / 1.55 mm | Single |
| 344-18030 | 70" / 180 cm | 0.79" / 20 mm | 0.12" / 3 mm | 0.06" / 1.55 mm | Single |
| 344-23030 | 90" / 230 cm | 0.79" / 20 mm | 0.12" / 3 mm | 0.06" / 1.55 mm | Single |

PACKAGING: 50 / Box

SERIES 35



PRODUCT SELECTION

2.5 MM CYLINDRICAL BRUSHES – SHEATH Ø 1.00 MM - YELLOW DELRIN

| ITEM # | LENGTH | BRUSH LENGTH | BRUSH (Ø) | SHEATH (Ø) | BRUSH ENDS |
|-------------------|----------------|---------------|-----------------|------------------|------------|
| 35-12025 | 47.2" / 120 cm | 0.59" / 15 mm | 0.098" / 2.5 mm | 0.039" / 1.00 mm | Single |
| 35-23025 | 90" / 230 cm | 0.59" / 15 mm | 0.098" / 2.5 mm | 0.039" / 1.00 mm | Single |
| 35-23025-2 | 90" / 230 cm | 0.59" / 15 mm | 0.098" / 2.5 mm | 0.039" / 1.00 mm | Double |


PACKAGING: 50 / Box

2 MM CYLINDRICAL BRUSHES – SHEATH Ø 1.00 MM - YELLOW DELRIN

| ITEM # | LENGTH | BRUSH LENGTH | BRUSH (Ø) | SHEATH (Ø) | BRUSH ENDS |
|-------------------|----------------|---------------|-----------------|------------------|------------|
| 35-12020 | 47.2" / 120 cm | 0.59" / 15 mm | 0.079" / 2.0 mm | 0.039" / 1.00 mm | Single |
| 35-18020 | 70.8" / 180 cm | 0.59" / 15 mm | 0.079" / 2.0 mm | 0.039" / 1.00 mm | Single |
| 35-23020 | 90" / 230 cm | 0.59" / 15 mm | 0.079" / 2.0 mm | 0.039" / 1.00 mm | Single |
| 35-23020-2 | 90" / 230 cm | 0.59" / 15 mm | 0.079" / 2.0 mm | 0.039" / 1.00 mm | Double |

PACKAGING: 50 / Box

VALVE / CONTROL BRUSHES




PRODUCT SELECTION

| ITEM # | DESCRIPTION | LENGTH | BRUSH 1 LENGTH | BRUSH 2 LENGTH | BRUSH 1 (Ø) | BRUSH 2 (Ø) |
|--|--|-----------------|----------------|----------------|---------------|--------------|
| 27-05197 | Double Trumpet 5/7 White Color | 5.00" / 13 cm | 0.20" / 5 mm | 1.00" / 25 mm | 0.20" / 5 mm | 0.28" / 7 mm |
| PACKAGING: Bag of 3 (single-use recommended). Volume Packaging: Bag of 50 use ITEM # 27-05197-50 | | | | | | |
| 27-57787 | Double. Universal 5/11 Yellow Color | 6.00" / 15 cm | 1.18" / 30 mm | 0.79" / 20 mm | 0.43" / 11 mm | 0.20" / 5 mm |
| PACKAGING: Bag of 3 (single-use recommended). Volume Packaging: Bag of 50 use ITEM # 27-57787-50 | | | | | | |
| 40-18511 | Double Universal 5/11 Extra Long Blue Color | 8.07" / 20.5 cm | 1.57" / 40 mm | 0.79" / 20 mm | 0.43" / 11 mm | 0.20" / 5 mm |
| PACKAGING: Bag of 3 (reusable and can be steam sterilized). Volume Packaging: Bag of 50 use ITEM # 40-18511-50 | | | | | | |
| 36-10010 | Single end large Red Color | 4.3" / 11 cm | 1.57" / 40 mm | N/A | 0.39" / 10 mm | N/A |
| PACKAGING: Bag of 3 (single-use recommended). Volume Packaging: Bag of 50 use ITEM # 36-10010-50 | | | | | | |
| 27-15311 | Single end large Blue Color | 4" / 10.3 cm | 1.57" / 40 mm | N/A | 0.39" / 10 mm | N/A |
| PACKAGING: Bag of 3 (single-use recommended). Volume Packaging: Bag of 50 use ITEM # 27-15311-50 | | | | | | |

V-BRUSHES

EXTERNAL CLEANING BRUSHES

Our External Cleaning Brushes for endoscopes / optical instruments provide excellent cleaning capabilities. With their unique U-shape bristles and innovative design, our brushes are the ultimate tool for cleaning the exterior of your endoscopic instruments. Free from sharp or metallic ends that can scratch or cause damage to valuable equipment, both single-use and reusable models are latex free and feature ergonomic handles.



PRODUCT SELECTION

| ITEM # | DESCRIPTION | BRUSH LENGTH / (Ø) | HANDLE LENGTH / (Ø) | TOTAL LENGTH | PACKAGING |
|-----------------|------------------------------------|--------------------|---------------------|--------------|-----------|
| 36-20005 | Single-Use External Cleaning Brush | 70 mm / 25 mm | 80 mm / 3 mm | 150 mm | 50 / Box |
| 36-20010 | Reusable External Cleaning Brush | 90 mm / 25 mm | 120 mm / 12 mm | 210 mm | 3 / Bag |

ARC™ WIRE STYLE BRUSHES

SINGLE-USE AND REUSABLE

You've got options with our ARC™ Wire Brushes. Our single-use and reusable brushes are available in both flexible and rigid varieties and come in a multitude of sizes. ARC™ Wire Brushes are guaranteed to last longer than twisted wire and preserve the integrity of your valuable instruments.

BRUSH SELECTION TIPS:

INSTRUMENT SPECIFICATIONS

What is the maximum length of the instrument(s) to be cleaned?

Add at least 2" (5cm) to that length, so you have an adequate length to work with.

What is the inner diameter of the instrument?

Add at most 0.125" (3mm) to determine the diameter of the bristle to ensure effective cleaning. The bristle diameter has to slightly exceed the inner tube diameter. Bristles that are too large will bend too much and not clean effectively.

BRUSH SELECTION

Simply find the section in our catalogue that best corresponds to the overall length and diameter of your instrument. As a rule, if your present brush is 8" (20.32cm) in length and we only offer 10" (25.4cm) length with your required diameter, that is entirely acceptable.


In many cases, our ARC™ Wire Brushes will have more bristles on the brush. As an example, many 5mm (0.19") diameter brushes are only 1" (2.5cm) in length. At Batrik, we strived to provide the greatest amount of bristles to ensure the most effective cleaning.

CUSTOMIZABLE WIRE BRUSHES


A unique aspect of single-use ARC™ Wire Brushes is that they can be customized to any length:

- Up to 2000mm for Wire Brushes on plastic semi-rigid sheaths with brushes from 5 to 20mm in diameter.
- Flexible plastic sheaths with brushes from 2-5mm in diameter.

RIGID



FLEXIBLE



12" / 30 CM SINGLE-END WIRE BRUSH

| | ITEM # SINGLE-USE | ITEM # REUSABLE | COLOR | TOTAL LENGTH | BRUSH (Ø) | BRUSH LENGTH | SHEATH (Ø) |
|-----------------------------------|----------------------|--------------------|--------|-------------------|-----------------|----------------|-----------------|
| FLEXIBLE BRUSH & SHEATH | 41-31502-S | 41-31502 | White | 12.401" / 31.5 cm | 0.078" / 2 mm | 0.59" / 15 mm | 0.39" / 1 mm |
| | 41-32503-S | 41-32503 | Yellow | 12.795" / 32.5 cm | 0.118" / 3 mm | 0.98" / 25 mm | 0.39" / 1 mm |
| | 411-32504-S | 411-32504 | Violet | 12.795" / 32.5 cm | 0.157" / 4 mm | 0.98" / 25 mm | 0.61" / 1.55 mm |
| | 411-32505-S | 411-32505 | Grey | 12.795" / 32.5 cm | 0.19" / 5 mm | 0.98" / 25 mm | 0.68" / 1.75 mm |
| RIGID BRUSH & SHEATH | 41-35005-S | 41-35005 | White | 13.779" / 35 cm | 0.19" / 5 mm | 1.968" / 50 mm | 0.118" / 3 mm |
| | 41-35006-S | 41-35006 | Green | 13.779" / 35 cm | 0.23" / 6 mm | 1.968" / 50 mm | 0.118" / 3 mm |
| | 41-35008-S | 41-35008 | Blue | 13.779" / 35 cm | 0.31" / 8 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| | 41-35010-S | 41-35010 | Red | 13.779" / 35 cm | 0.39" / 10 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| | 41-35012-S | 41-35012 | Yellow | 13.779" / 35 cm | 0.49" / 12.5 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| | 41-35015-S | 41-35015 | Green | 13.779" / 35 cm | 0.59" / 15 mm | 1.968" / 50 mm | 0.228" / 5.8 mm |
| | 41-35020-S | 41-35020 | Blue | 13.779" / 35 cm | 0.78" / 20 mm | 1.968" / 50 mm | 0.228" / 5.8 mm |

PACKAGING Single-Use: Bag of 20 / Reusable: Bag of 3

18" / 45 CM SINGLE-END WIRE BRUSH

| | ITEM # SINGLE-USE | ITEM # REUSABLE | COLOR | TOTAL LENGTH | BRUSH (Ø) | BRUSH LENGTH | SHEATH (Ø) |
|-----------------------------------|----------------------|--------------------|--------|-------------------|-----------------|----------------|-----------------|
| FLEXIBLE BRUSH & SHEATH | 41-46502-S | 41-46502 | White | 18.307" / 46.5 cm | 0.078" / 2 mm | 0.59" / 15 mm | 0.39" / 1 mm |
| | 41-47503-S | 41-47503 | Yellow | 18.7" / 47.5 cm | 0.118" / 3 mm | 0.98" / 25 mm | 0.39" / 1 mm |
| | 411-47504-S | 411-47504 | Violet | 18.7" / 47.5 cm | 0.157" / 4 mm | 0.98" / 25 mm | 0.61" / 1.55 mm |
| | 411-47505-S | 411-47505 | Grey | 18.7" / 47.5 cm | 0.19" / 5 mm | 0.98" / 25 mm | 0.68" / 1.75 mm |
| RIGID BRUSH & SHEATH | 41-50005-S | 41-50005 | White | 19.685" / 50 cm | 0.19" / 5 mm | 1.968" / 50 mm | 0.118" / 3 mm |
| | 41-50006-S | 41-50006 | Green | 19.685" / 50 cm | 0.23" / 6 mm | 1.968" / 50 mm | 0.118" / 3 mm |
| | 41-50008-S | 41-50008 | Blue | 19.685" / 50 cm | 0.31" / 8 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| | 41-50010-S | 41-50010 | Red | 19.685" / 50 cm | 0.39" / 10 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| | 41-50012-S | 41-50012 | Yellow | 19.685" / 50 cm | 0.49" / 12.5 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| | 41-50015-S | 41-50015 | Green | 19.685" / 50 cm | 0.59" / 15 mm | 1.968" / 50 mm | 0.228" / 5.8 mm |
| | 41-50020-S | 41-50020 | Blue | 19.685" / 50 cm | 0.78" / 20 mm | 1.968" / 50 mm | 0.228" / 5.8 mm |

PACKAGING: Single -Use: Bag of 20 / Reusable: Bag of 3

24" / 60 CM SINGLE-END WIRE BRUSH


| | ITEM # SINGLE-USE | ITEM # REUSABLE | COLOR | TOTAL LENGTH | BRUSH (Ø) | BRUSH LENGTH | SHEATH (Ø) |
|-----------------------------------|----------------------|--------------------|--------|-------------------|-----------------|----------------|-----------------|
| FLEXIBLE BRUSH & SHEATH | 41-61502-S | 41-61502 | White | 24.212" / 61.5 cm | 0.078" / 2 mm | 0.59" / 15 mm | 0.39" / 1mm |
| | 41-62503-S | 41-62503 | Yellow | 24.6" / 62.5 cm | 0.118" / 3 mm | 0.98" / 25 mm | 0.39" / 1mm |
| | 411-62504-S | 411-62504 | Violet | 24.6" / 62.5 cm | 0.157" / 4 mm | 0.98" / 25 mm | 0.61" / 1.55mm |
| | 411-62505-S | 411-62505 | Grey | 24.6" / 62.5 cm | 0.19" / 5 mm | 0.98" / 25 mm | 0.68" / 1.75mm |
| RIGID BRUSH & SHEATH | 41-65005-S | 41-65005 | White | 25.59" / 65 cm | 0.19" / 5 mm | 1.968" / 50 mm | 0.118" / 3mm |
| | 41-65006-S | 41-65006 | Green | 25.59" / 65 cm | 0.23" / 6 mm | 1.968" / 50 mm | 0.118" / 3mm |
| | 41-65008-S | 41-65008 | Blue | 25.59" / 65 cm | 0.31" / 8 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| | 41-65010-S | 41-65010 | Red | 25.59" / 65 cm | 0.39" / 10 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| | 41-65012-S | 41-65012 | Yellow | 25.59" / 65 cm | 0.49" / 12.5 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| | 41-65015-S | 41-65015 | Green | 25.59" / 65 cm | 0.59" / 15 mm | 1.968" / 50 mm | 0.228" / 5.8 mm |
| | 41-65020-S | 41-65020 | Blue | 25.59" / 65 cm | 0.78" / 20 mm | 1.968" / 50 mm | 0.228" / 5.8 mm |

PACKAGING: Single-Use: Bag of 20 / Reusable: Bag of 3


30" / 75 CM SINGLE-END WIRE BRUSH

| | ITEM # SINGLE-USE | ITEM # REUSABLE | COLOR | TOTAL LENGTH | BRUSH (Ø) | BRUSH LENGTH | SHEATH (Ø) |
|-----------------------------------|----------------------|--------------------|--------|-------------------|-----------------|----------------|-----------------|
| FLEXIBLE BRUSH & SHEATH | 41-76502-S | 41-76502 | White | 30.118" / 76.5 cm | 0.078" / 2 mm | 0.59" / 15 mm | 0.39" / 1 mm |
| | 41-77503-S | 41-77503 | Yellow | 30.511" / 77.5 cm | 0.118" / 3 mm | 0.98" / 25 mm | 0.39" / 1 mm |
| | 411-77504-S | 411-77504 | Violet | 30.511" / 77.5 cm | 0.157" / 4 mm | 0.98" / 25 mm | 0.61" / 1.55 mm |
| | 411-77505-S | 411-77505 | Grey | 30.511" / 77.5 cm | 0.19" / 5 mm | 0.98" / 25 mm | 0.68" / 1.75 mm |
| RIGID BRUSH & SHEATH | 41-80005-S | 41-80005 | White | 31.496" / 80 cm | 0.19" / 5 mm | 1.968" / 50 mm | 0.118" / 3 mm |
| | 41-80006-S | 41-80006 | Green | 31.496" / 80 cm | 0.23" / 6 mm | 1.968" / 50 mm | 0.118" / 3 mm |
| | 41-80008-S | 41-80008 | Blue | 31.496" / 80 cm | 0.31" / 8 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| | 41-80010-S | 41-80010 | Red | 31.496" / 80 cm | 0.39" / 10 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| | 41-80012-S | 41-80012 | Yellow | 31.496" / 80 cm | 0.49" / 12.5 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| | 41-80015-S | 41-80015 | Green | 31.496" / 80 cm | 0.59" / 15 mm | 1.968" / 50 mm | 0.228" / 5.8 mm |
| | 41-80020-S | 41-80020 | Blue | 31.496" / 80 cm | 0.78" / 20 mm | 1.968" / 50 mm | 0.228" / 5.8 mm |

PACKAGING: Single-Use: Bag of 20 / Reusable: Bag of 3



FLEXIBLE
Flexible brush diameter (Ø) ranging from 2 mm to 5 mm



RIGID
Rigid brush diameter (Ø) ranging from 5 mm to 20 mm

BRUSH SELECTION SYSTEM (Actual Size)

| 2mm | 3 mm | 4 mm | 5 mm | 6 mm | 8 mm | 10 mm | 12 mm | 15 mm | 20 mm |
|-------|--------|--------|--------|--------|--------|--------|--------|--------|--------|
| 0.078 | 0.118" | 0.157" | 0.196" | 0.236" | 0.314" | 0.394" | 0.472" | 0.591" | 0.787" |
| | | | | | | | | | |

REUSABLE DOUBLE ARC™ WIRE BRUSHES

ASYMMETRIC

14" / 35 CM ASYMMETRIC DOUBLE-END WIRE BRUSHES

| ITEM # | COLOR | TOTAL LENGTH | BRUSH 1 (Ø) | BRUSH 2 (Ø) | BRUSH LENGTH | SHEATH (Ø) |
|----------|--------|-----------------|---------------|-----------------|----------------|-----------------|
| 40-35056 | White | 13.779" / 35 cm | 0.19" / 5 mm | 0.23" / 6 mm | 1.968" / 50 mm | 0.118" / 3 mm |
| 40-35810 | Blue | 13.779" / 35 cm | 0.31" / 8 mm | 0.39" / 10 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| 40-35112 | Yellow | 13.779" / 35 cm | 0.39" / 10 mm | 0.49" / 12.5 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| 40-35152 | Green | 13.779" / 35 cm | 0.59" / 15 mm | 0.78" / 20 mm | 1.968" / 50 mm | 0.228" / 5.8 mm |

PACKAGING: Bag of 3

20" / 50 CM ASYMMETRIC DOUBLE-END WIRE BRUSHES

| ITEM # | COLOR | TOTAL LENGTH | BRUSH 1 (Ø) | BRUSH 2 (Ø) | BRUSH LENGTH | SHEATH (Ø) |
|----------|--------|-----------------|---------------|-----------------|----------------|-----------------|
| 40-50056 | White | 19.685" / 50 cm | 0.19" / 5 mm | 0.23" / 6 mm | 1.968" / 50 mm | 0.118" / 3 mm |
| 40-50810 | Blue | 19.685" / 50 cm | 0.31" / 8 mm | 0.39" / 10 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| 40-50112 | Yellow | 19.685" / 50 cm | 0.39" / 10 mm | 0.49" / 12.5 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| 40-50152 | Green | 19.685" / 50 cm | 0.59" / 15 mm | 0.78" / 20 mm | 1.968" / 50 mm | 0.228" / 5.8 mm |

PACKAGING: Bag of 3

26" / 65 CM ASYMMETRIC DOUBLE-END WIRE BRUSHES

| ITEM # | COLOR | TOTAL LENGTH | BRUSH 1 (Ø) | BRUSH 2 (Ø) | BRUSH LENGTH | SHEATH (Ø) |
|----------|--------|----------------|---------------|-----------------|----------------|-----------------|
| 40-65056 | White | 25.59" / 65 cm | 0.19" / 5 mm | 0.23" / 6 mm | 1.968" / 50 mm | 0.118" / 3 mm |
| 40-65810 | Blue | 25.59" / 65 cm | 0.31" / 8 mm | 0.39" / 10 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| 40-65112 | Yellow | 25.59" / 65 cm | 0.39" / 10 mm | 0.49" / 12.5 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| 40-65152 | Green | 25.59" / 65 cm | 0.59" / 15 mm | 0.78" / 20 mm | 1.968" / 50 mm | 0.228" / 5.8 mm |

PACKAGING: Bag of 3

31" / 80 CM ASYMMETRIC DOUBLE-END WIRE BRUSHES

| ITEM # | COLOR | TOTAL LENGTH | BRUSH 1 (Ø) | BRUSH 2 (Ø) | BRUSH LENGTH | SHEATH (Ø) |
|----------|--------|-----------------|---------------|-----------------|----------------|-----------------|
| 40-80056 | White | 31.496" / 80 cm | 0.19" / 5 mm | 0.23" / 6 mm | 1.968" / 50 mm | 0.118" / 3 mm |
| 40-80810 | Blue | 31.496" / 80 cm | 0.31" / 8 mm | 0.39" / 10 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| 40-80112 | Yellow | 31.496" / 80 cm | 0.39" / 10 mm | 0.49" / 12.5 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| 40-80152 | Green | 31.496" / 80 cm | 0.59" / 15 mm | 0.78" / 20 mm | 1.968" / 50 mm | 0.228" / 5.8 mm |

PACKAGING: Bag of 3

REUSABLE DOUBLE ARC™ WIRE BRUSHES SYMMETRIC

14" / 35 CM SYMMETRIC DOUBLE-END WIRE BRUSHES

| ITEM # | COLOR | TOTAL LENGTH | BRUSH 1 (Ø) | BRUSH LENGTH | SHEATH (Ø) |
|----------|--------|-----------------|-----------------|----------------|-----------------|
| 40-35005 | White | 13.779" / 35 cm | 0.19" / 5 mm | 1.968" / 50 mm | 0.118" / 3 mm |
| 40-35006 | Green | 13.779" / 35 cm | 0.23" / 6 mm | 1.968" / 50 mm | 0.118" / 3 mm |
| 40-35008 | Blue | 13.779" / 35 cm | 0.31" / 8 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| 40-35010 | Red | 13.779" / 35 cm | 0.39" / 10 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| 40-35012 | Yellow | 13.779" / 35 cm | 0.49" / 12.5 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| 40-35015 | Green | 13.779" / 35 cm | 0.59" / 15 mm | 1.968" / 50 mm | 0.228" / 5.8 mm |
| 40-35020 | Blue | 13.779" / 35 cm | 0.78" / 20 mm | 1.968" / 50 mm | 0.228" / 5.8 mm |

PACKAGING: Bag of 3

20" / 50 CM SYMMETRIC DOUBLE-END WIRE BRUSHES

| ITEM # | COLOR | TOTAL LENGTH | BRUSH 1 (Ø) | BRUSH LENGTH | SHEATH (Ø) |
|----------|--------|-----------------|-----------------|----------------|-----------------|
| 40-50005 | White | 19.685" / 50 cm | 0.19" / 5 mm | 1.968" / 50 mm | 0.118" / 3 mm |
| 40-50006 | Green | 19.685" / 50 cm | 0.23" / 6 mm | 1.968" / 50 mm | 0.118" / 3 mm |
| 40-50008 | Blue | 19.685" / 50 cm | 0.31" / 8 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| 40-50010 | Red | 19.685" / 50 cm | 0.39" / 10 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| 40-50012 | Yellow | 19.685" / 50 cm | 0.49" / 12.5 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| 40-50015 | Green | 19.685" / 50 cm | 0.59" / 15 mm | 1.968" / 50 mm | 0.228" / 5.8 mm |
| 40-50020 | Blue | 19.685" / 50 cm | 0.78" / 20 mm | 1.968" / 50 mm | 0.228" / 5.8 mm |

PACKAGING: Bag of 3

26" / 65 CM SYMMETRIC DOUBLE-END WIRE BRUSHES

| ITEM # | COLOR | TOTAL LENGTH | BRUSH 1 (Ø) | BRUSH LENGTH | SHEATH (Ø) |
|----------|--------|----------------|-----------------|----------------|-----------------|
| 40-65005 | White | 25.59" / 65 cm | 0.19" / 5 mm | 1.968" / 50 mm | 0.118" / 3 mm |
| 40-65006 | Green | 25.59" / 65 cm | 0.23" / 6 mm | 1.968" / 50 mm | 0.118" / 3 mm |
| 40-65008 | Blue | 25.59" / 65 cm | 0.31" / 8 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| 40-65010 | Red | 25.59" / 65 cm | 0.39" / 10 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| 40-65012 | Yellow | 25.59" / 65 cm | 0.49" / 12.5 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| 40-65015 | Green | 25.59" / 65 cm | 0.59" / 15 mm | 1.968" / 50 mm | 0.228" / 5.8 mm |
| 40-65020 | Blue | 25.59" / 65 cm | 0.78" / 20 mm | 1.968" / 50 mm | 0.228" / 5.8 mm |

PACKAGING: Bag of 3

31" / 80 CM SYMMETRIC DOUBLE-END WIRE BRUSHES

| ITEM # | COLOR | TOTAL LENGTH | BRUSH 1 (Ø) | BRUSH LENGTH | SHEATH (Ø) |
|----------|--------|-----------------|-----------------|----------------|-----------------|
| 40-80005 | White | 31.496" / 80 cm | 0.19" / 5 mm | 1.968" / 50 mm | 0.118" / 3 mm |
| 40-80006 | Green | 31.496" / 80 cm | 0.23" / 6 mm | 1.968" / 50 mm | 0.118" / 3 mm |
| 40-80008 | Blue | 31.496" / 80 cm | 0.31" / 8 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| 40-80010 | Red | 31.496" / 80 cm | 0.39" / 10 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| 40-65112 | Yellow | 31.496" / 80 cm | 0.49" / 12.5 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| 40-80015 | Green | 31.496" / 80 cm | 0.59" / 15 mm | 1.968" / 50 mm | 0.228" / 5.8 mm |
| 40-80020 | Blue | 31.496" / 80 cm | 0.78" / 20 mm | 1.968" / 50 mm | 0.228" / 5.8 mm |

PACKAGING: Bag of 3

ASYMMETRIC BRUSHES

5 MM

6 MM



8 MM

10 MM



10 MM

12.5 MM



15 MM

20 MM




REUSABLE IN-BRUSH™ ANTIMICROBIAL TWISTED WIRE BRUSHES

Our In-Brush™ Antimicrobial Reusable Twisted Wire Brushes are designed for use with Frazier/Baron suction tubes. The twisted wires have a fine/medium diameter yet are rigid enough to ensure excellent cleaning capabilities.

Our wide assortment of Nylon Twisted Wire Brushes are ideal for large bore, uniquely shaped, and endoscopic instruments.

GOOD TO KNOW: OUR TWISTED WIRE BRUSHES FEATURE A PROTECTIVE BLUE TIP FOR ADDED SAFETY.



10" / 250 MM FRAZIER /BARON SUCTION TUBE BRUSHES

| ITEM # | SIZE | BRUSH LENGTH | BRUSH (Ø) |
|----------|-------|----------------|------------------|
| 21-10039 | 3 FR | 0.5" / 12.7 mm | 0.039" / 1 mm |
| 21-10067 | 5 FR | 0.5" / 12.7 mm | 0.065" / 1.67 mm |
| 20-10059 | 5 FR | 1" / 25.4 mm | 0.059" / 1.50 mm |
| 20-10079 | 6 FR | 1" / 25.4 mm | 0.079" / 2 mm |
| 20-10090 | 7 FR | 1" / 25.4 mm | 0.090" / 2.3 mm |
| 20-10106 | 8 FR | 1" / 25.4 mm | 0.106" / 2.7 mm |
| 20-10118 | 9 FR | 1" / 25.4 mm | 0.118" / 3 mm |
| 20-10129 | 10 FR | 1.5" / 38 mm | 0.129" / 3.3 mm |
| 20-10146 | 11 FR | 1.5" / 38 mm | 0.0146" / 3.7 mm |
| 20-10157 | 12 FR | 1.5" / 38 mm | 0.157" / 4 mm |

PACKAGING: Bag of 3.
Also available in bags of 50, add the suffix "-50".

12" / 305 MM NYLON TWISTED WIRE BRUSHES

| ITEM # | BRUSH LENGTH | BRUSH (Ø) |
|----------|---------------|----------------|
| 22-12078 | 2.28" / 58 mm | 0.078" / 2 mm |
| 22-12118 | 2.28" / 58 mm | 0.118" / 3 mm |
| 22-12157 | 3" / 76 mm | 0.157" / 4 mm |
| 22-12197 | 3" / 76 mm | 0.197" / 5 mm |
| 22-12236 | 3" / 76 mm | 0.236" / 6 mm |
| 22-12393 | 3" / 76 mm | 0.394" / 10 mm |

PACKAGING: Bag of 3.
Also available in bags of 50, add the suffix "- 50".

16" / 406 MM NYLON TWISTED WIRE BRUSHES

| ITEM # | BRUSH LENGTH | BRUSH (Ø) |
|----------|----------------|----------------|
| 23-16197 | 2.5" / 63.5 mm | 0.197" / 5 mm |
| 23-16236 | 2.5" / 63.5 mm | 0.236" / 6 mm |
| 23-16276 | 3" / 76 mm | 0.276" / 7 mm |
| 23-16394 | 3" / 76 mm | 0.394" / 10 mm |
| 23-16472 | 3" / 76 mm | 0.472" / 12 mm |
| 23-16591 | 3" / 76 mm | 0.591" / 15 mm |
| 23-16787 | 3" / 76 mm | 0.787" / 20 mm |

PACKAGING: Bag of 3.
Also available in bags of 50, add the suffix "-50".

18" / 457 MM NYLON TWISTED WIRE BRUSHES

| ITEM # | BRUSH LENGTH | BRUSH (Ø) |
|----------|---------------|----------------|
| 19-18010 | 0.47" / 12 mm | 0.039" / 1 mm |
| 19-18015 | 0.5" / 12 mm | 0.06" / 1.6 mm |

PACKAGING: Bag of 3.
Also available in bags of 50, add the suffix "-50".

24" / 610 MM NYLON TWISTED WIRE BRUSHES

| ITEM # | BRUSH LENGTH | BRUSH (Ø) |
|----------|-----------------|-----------------|
| 24-24098 | 2.5" / 63.5 mm | 0.098" / 2.5 mm |
| 24-24197 | 2.5" / 63.5 mm | 0.197" / 5 mm |
| 24-24276 | 3.25" / 82.6 mm | 0.276" / 7 mm |
| 24-24394 | 3.25" / 82.6 mm | 0.394" / 10 mm |
| 24-24472 | 3.25" / 82.6 mm | 0.472" / 12 mm |

PACKAGING: Bag of 3.
Also available in bags of 50, add the suffix "-50".

TRACHEOSTOMY BRUSHES

| ITEM # | DESCRIPTION | TOTAL LENGTH | BRUSH LENGTH | BRUSH (Ø) |
|----------|----------------------|--------------|---------------|----------------|
| 25-08500 | Large 6-12.3 3 / 4" | 8" / 203 mm | 3.75" / 95 mm | 0.5" / 12.7 mm |
| 25-37550 | Small 000-5.1 7 / 8" | 4" / 93 mm | 1.88" / 48 mm | 0.5" / 12.7 mm |

PACKAGING: Bag of 3. Also available in bags of 50, add the suffix "-50".



BRONCHO BRUSHES


| ITEM # | DESCRIPTION | BRUSH LENGTH | BRUSH (Ø) INCHES / MM |
|----------|-------------|--------------|-----------------------|
| 26-06175 | Small | 6" / 152 mm | 1.75" / 44.5 mm |
| 26-08325 | Large | 8" / 203 mm | 3.25" / 82.6 mm |

PACKAGING: Bag of 3.
Also available in bags of 50, add the suffix "-50".

BRIGHTBRUSH

THE RIGHT BRUSH FOR EVERY JOB

Tackle even the toughest cleaning jobs with ease! Our Color-Coded Twisted Wire Brushes are designed for precision, convenience, and superior performance—ensuring spotless results every time.



Color-Coded Bead

Each brush features a color-coded bead that indicates the bristle diameter and size, making it easy to select the right tool quickly and efficiently.



Hanging Loop

The brush includes a loop at the end for convenient hanging.

PRODUCT SELECTION

| ITEM # | WIRE (Ø) (mm) | BEAD (Ø) (mm) | BEAD COLOR | DIMENSIONS Length x Brush (Ø) x Brush Length |
|---|------------------|------------------|------------|---|
| 12" Color-coded Twisted Wire Brushes | | | | |
| 22-12078B | 0.7 | 4 | Red | 12" x .078" (2M)x 58mm |
| 22-12118B | 0.7 | 4 | White | 12" x .118" (3M)x 58mm |
| 22-12157B | 0.8 | 4 | Blue | 12" x .157" (4M)x 76mm |
| 22-12197B | 0.8 | 4 | Green | 12" x .197" (5M)x 76mm |
| 22-12236B | 0.8 | 4 | Black | 12" x .236" (6M)x 76mm |
| 16" Color-coded Twisted Wire Brushes | | | | |
| 23-16197B | 1.0 | 4 | Green | 16" X .197" (5M)x 63.5mm |
| 23-16276B | 1.0 | 4 | Yellow | 16" X .276" (7M)x 76mm |
| 23-16394B | 1.0 | 4 | Pink | 16" X .394" (10M)x 76mm |
| 18" Color-coded Twisted Wire Brushes | | | | |
| 19-18078B | 0.6 | 4 | Red | 18" x .078" (2M)x 63.5mm |
| 19-18118B | 0.8 | 4 | White | 18" x .118" (3M)x 63.5mm |
| 19-18157B | 0.8 | 4 | Blue | 18" x .157" (4M)x 63.5mm |
| 19-18197B | 1.0 | 4 | Green | 18" x .197" (5M)x 63.5mm |
| 19-18236B | 1.0 | 4 | Black | 18" x .236" (6M)x 76mm |
| 24" Color-coded Twisted Wire Brushes | | | | |
| 24-24098B | 0.8 | 4 | Orange | 24" X .098" (2.5M)x 63.5mm |
| 24-24197B | 0.8 | 4 | Green | 24" X .197" (5M)x 63.5mm |
| 24-24276B | 1.0 | 4 | Yellow | 24" X .276" (7M)x 82.6mm |
| 24-24394B | 1.0 | 4 | Pink | 24" X .394" (10M)x 82.6mm |

PACKAGING: Bag of 3. Also available in bags of 50, add the suffix "-50".

STAINLESS STEEL ARC™ MICRO BRUSHES

Our Stainless Steel ARC™ Micro Brushes are some of the smallest brushes in the industry. They are designed to clean channels with an inner diameter ranging from 0.5 mm to 1.40 mm.

- Flexible stainless steel sheath will not bend like the regular metal twisted wire brushes
- Polyamide bristles prevents or slows bacterial growth between sterilizations
- Reusable*
- Available in two lengths

TIP: USING THE PULL-THROUGH
TECHNIQUE IS HIGHLY
RECOMMENDED WITH THESE
BRUSHES. SIMPLY INSERT THE
BLUNT, NON-BRUSH END FIRST.



PRODUCT SELECTION

| ITEM # | TOTAL LENGTH | BRUSH SIZE (Ø) | BRUSH LENGTH | SHEATH (Ø) | PACKAGING |
|-----------------|--------------|----------------|--------------|-----------------|-----------|
| 35-03015 | 12" / 30 cm | 0.06" / 1.5 mm | 0.6" / 15 mm | 0.02" / 0.40 mm | 6 / Tube |
| 35-04515 | 18" / 45 cm | 0.06" / 1.5 mm | 0.6" / 15 mm | 0.02" / 0.40 mm | 6 / Tube |

INSTRUMENT CLEANING FLAT BRUSHES

Our Instrument Cleaning Flat Brushes are specifically designed for hospital use in the pre-clean and preparation of instruments.

- Long ergonomic handle permits ease-of-use
- Extra long length of brush head (60 [±5] mm in Single-Use and 75 [±5] mm in Reusable)
- Transparent brush fibers allow for quick verification of its cleanliness
- The plastic handle will not release any toxic components
- Available in 6 models (4 reusable/autoclavable and 2 single-use)



FLEXIBLE

This brush is designed for large instruments that require a more gentle cleaning application such as optics, plastics, or hard to reach spot.



EXTRA FLEXIBLE

The very long, flexible bristles eliminate any stress during sanitation, and makes this the ideal brush for fragile instruments as well as both pliable and rigid endoscopes. The length of the bristles assist in cleaning areas that would otherwise be difficult to reach, such as the accessory channel or tip deflection of pliable endoscopes.



RIGID

The rigidity of this brush, along with the length of the bristles and handle, improves and simplifies instrument cleaning.



EXTRA RIGID

Instrument safe polyamide fibers offer steel-like rigidity and allows deep cleaning of soiled instruments, particularly in grooved or uneven surfaces.



50-40919



50-40921



50-40920



50-40918



50-40922



50-40917

PRODUCT SELECTION

| ITEM # | BRUSH TYPE | DESCRIPTION | COLOR | BRUSH LENGTH | BRUSH SIZE | PACKAGING |
|----------|----------------|----------------------------------|------------|--------------|-----------------------|-----------|
| 50-40917 | Extra Rigid | Reusable Autoclavable Flat Brush | Blue | 220 mm | 75 mm x 20 mm x 10 mm | Bag of 3 |
| 50-40918 | Rigid | Reusable Autoclavable Flat Brush | Green | 220 mm | 75 mm x 20 mm x 15 mm | Bag of 3 |
| 50-40919 | Flexible | Reusable Autoclavable Flat Brush | Yellow | 220 mm | 75 mm x 20 mm x 15 mm | Bag of 3 |
| 50-40920 | Extra Flexible | Reusable Autoclavable Flat Brush | Red | 220 mm | 75 mm x 20 mm x 25 mm | Bag of 3 |
| 50-40921 | Flexible | Single-Use Flat Brush | Orange | 220 mm | 60 mm x 20 mm x 15 mm | Bag of 20 |
| 50-40922 | Rigid | Single-Use Flat Brush | Light Blue | 220 mm | 60 mm x 20 mm x 15 mm | Bag of 20 |

Reusable models are also available in bags of 50. Add the suffix "- 50".

REUSABLE PRECISION BRUSH

INSTRUMENT CLEANING

Our Reusable Precision Instrument Cleaning Brush is the ideal brush for the quick removal of difficult debris. It cleans externally as well as deep into instrument jaws. Although the firm and durable bristles feel like stainless steel, there is absolutely no risk of instrument damage or personal injury.

- Bristles have an extra firm density and superior durability
- Polyamide (TYNEX®) bristles
- Steamable thermoplastic ergonomic handle
- Steam / chemical sterilizable
- Latex free



GOOD TO KNOW: ELIMINATES THE NEED FOR THE STAINLESS STEEL TOOTHBRUSH-STYLE BRUSH, WHICH CAN DAMAGE VALUABLE INSTRUMENTS!

PRODUCT SELECTION

| ITEM # | DESCRIPTION | TOTAL LENGTH |
|----------|--|----------------|
| 50-40800 | Reusable Precision External Cleaning Brush | 6.29" / 160 mm |

PACKAGING: Bag of 3

NEW!

TRIFECTA™ BRUSH

INSTRUMENT CLEANING

Introducing the **Trifecta™ Brush**, a versatile **3-in-1 tool** designed to meet all your instrument cleaning needs. This reusable instrument brush is available in two styles: one with **all-nylon** bristles for gentle cleaning and another with nylon bristles on two sides and **stainless steel bristles** on the third, ideal for tackling tougher buildup and stubborn cleaning tasks. This innovative design ensures effective, thorough maintenance for a wide range of instruments.

TRIFECTA™ BRUSH - NYLON

1 NYLON BRISTLES

Gentle, transparent nylon bristles provide an effective yet safe cleaning solution. This side allows you to verify cleanliness easily, without risking damage to delicate surfaces. Ideal for a light, non-invasive clean.

2 PRECISION TIP

Reach those hard-to-clean spots with the precision tip. This small, firm-bristle tip is designed for tight, hard-to-reach areas, such as instrument jaws and box joints. It ensures a deep clean where it's most needed, without the risk of damage or injury to delicate parts.

3 FIRM NYLON BRISTLES

For tougher residues, the firm nylon side delivers the cleaning power you need. With bristles that feel similar to stainless steel but are safe for your instruments, this side easily handles stubborn dirt and buildup.



TRIFECTA™ BRUSH - NYLON & STAINLESS STEEL



29-23012

28-23012

PRODUCT SELECTION

| ITEM # | COLOR | DESCRIPTION | LENGTH TOTAL | BRISTLE LENGTH | BRISTLES 1 Ø | BRISTLES 2 Ø | BRISTLES 3 Ø |
|----------|-------|--------------------------------|----------------|---|-----------------|-----------------|-----------------|
| 28-23012 | Blue | Triple Nylon Brush | 8.66" / 220 mm | 1.20" x 1.20" x 0.47" / 30.5 x 30.5 x 12.0 mm | 0.47" / 12.0 mm | 0.47" / 12.0 mm | 0.315" / 8.0 mm |
| 29-23012 | Blue | Triple Nylon / Stainless Steel | 8.66" / 220 mm | 1.20" x 1.20" x 0.47" / 30.5 x 30.5 x 12.0 mm | 0.47" / 12.0 mm | 0.47" / 12.0 mm | 0.47" / 12.0 mm |

PACKAGING: Bag of 3 - For a bag of 50, add "-50" following the item code.

ALL PURPOSE REUSABLE BRUSHES

Plastic Handle Brushes

We offer a wide selection of brushes for instrument cleaning applications. These brushes are ideal for large bore, uniquely shaped, and endoscopic instruments.

In-Brush™ Plastic Handle

Available in regular toothbrush style (nylon or stainless) and double-ended brushes.



PRODUCT SELECTION

| ITEM # | DESCRIPTION | COLOR | BRUSH LENGTH | BRISTLE SIZE | PACKAGING |
|-----------------|--|-------|--------------|---------------------------|-----------|
| 28-72550 | Plastic Single End Nylon Reusable Handle Brush | BLACK | 184mm | 38mm x 12.7mm | 3 |
| 29-72550 | Plastic Single End Stainless Steel Reusable Handle Brush | BLACK | 184mm | 38mm x 12.7mm | 3 |
| 29-70050 | Plastic Single End Brass Reusable Handle Brush | BLACK | 184mm | 38mm x 12.7mm | 3 |
| 28-69394 | Plastic Double End Nylon Reusable Handle Brush | BLACK | 175mm | 35mm x 10mm 25mm x 5mm | 3 |
| 29-69394 | Plastic Double End Stainless Steel Reusable Handle Brush | BLACK | 175mm | 35mm x 10mm 25mm x 5mm | 3 |
| 50-40911 | Hand Scrub Nylon Brush Reusable autoclavable | WHITE | 104mm | 95mm x 28mm | 3 |

Reusable models are also available in bags of 50. Add the suffix "- 50"

NEW!
COMING SOON

Engineered for precision and durability, this reusable dual-ended is designed to tackle even the most complex instruments. Its ultra-fine tip reaches into the smallest, hardest-to-access areas—delivering exceptional hygiene with every use. Built from long-lasting materials, this brush not only enhances cleaning performance but also supports hospital sustainability goals by reducing single-use waste and associated costs.



50-40517 EXTRA-RIGID **50-40518** RIGID **50-40520** EXTRA-FLEXIBLE

PRODUCT SELECTION

| ITEM # | DESCRIPTION | COLOR | BRUSH LENGTH | BRISTLE SIZE 1 | BRISTLE SIZE 2 | PACKAGING |
|-----------------|---|-------|--------------|-------------------|-------------------|-----------|
| 50-40517 | Reusable Extra Rigid Flat Brush Dual Ended | BLUE | 183mm | 38mm x 9mm x 10mm | 28mm x 3mm x 8mm | 3 |
| 50-40518 | Reusable Rigid Flat Brush Dual Ended | GREEN | 183mm | 38mm x 9mm x 15mm | 28mm x 3mm x 13mm | 3 |
| 50-40520 | Reusable Extra Flexible Flat Brush Dual Ended | RED | 183mm | 38mm x 9mm x 25mm | 28mm x 3mm x 23mm | 3 |

Reusable models are also available in bags of 50. Add the suffix "- 50"

ACETABULAR REAMER BRUSH

Our double-ended cleaning brush is manufactured exclusively for cleaning acetabular bone reamers. The “Pig Tail” rounded style circular head brush is designed to effectively clean inside the reamer, while the rigid straight brush on the reverse side is used for cleaning individual cutting holes. Special fusion of the brush heads to the hand piece along with its high quality antimicrobial medical grade bristles ensures a durable, reusable, and highly effective cleaning tool.



STRAIGHT BRUSH



ROUND BRUSH



PRODUCT SELECTION

| ITEM # | DESCRIPTION | HANDLE LENGTH | STRAIGHT BRUSH | ROUND BRUSH HEAD |
|------------------|----------------------------------|-----------------------------------|-----------------------------------|-------------------------------|
| 50-ACETAB | Acetabular Reamer Cleaning Brush | 4.72" / 120 mm ø 0.47" / 12 mm | 1.18" / 30 mm ø 0.3125" / 8 mm | 2" / 50 mm ø 1.75" / 45 mm |


PACKAGING: Bag of 3

NOTE: Specifications subject to small variances

REUSABLE ORTHOPEDIC RASPS CLEANING BRUSH

Our Reusable Orthopedic Rasps Cleaning Brush is used to clean orthopedic rasps and other similar instruments used in orthopedic surgery.

- Long lasting and can withstand many sterilization cycles
- Ergonomic steamable thermoplastic handle
- Antibacterial polyamide bristles
- Reusable
- Latex free



PRODUCT DESCRIPTION

| ITEM # | TOTAL LENGTH | HANDLE LENGTH | BRUSH SIZE | PACKAGING |
|---------|--------------|---------------|---|-----------|
| 50-RASP | 9" / 23 cm | 5.5" / 14 cm | 0.78" x 2.7" x 3.5" 2 cm x 7 cm x 9 cm | 3 / Bag |

BOTTLE & DRAIN CLEANING BRUSHES


Designed specifically for hospital applications, our Bottle & Drain Cleaning Brushes feature long and flexible bristles to maximize the effectiveness of large-surface cleaning.

- High quality polyamide clear bristles
- Easy to identify debris
- Autoclavable
- Latex free

REUSABLE JAR/BOTTLE BRUSH



STERILIZER DRAIN BRUSH



Ergonomic Handle

Blue Tip for added safety

Heavy Duty Stainless Wire

PRODUCT SELECTION

| ITEM # | DESCRIPTION | HANDLE LENGTH | BRUSH LENGTH | BRUSH (Ø) | TOTAL LENGTH |
|-----------------|---------------------------------|---------------|---------------|--------------|-----------------|
| 50-10295 | Reusable Jar / Bottle Brush | 120 mm ± 2 mm | 120 mm ± 5 mm | 7 mm ± 2 mm | 13.19" / 355 mm |
| 50-10224 | Sterilizer Drain Cleaning Brush | N/A | 4" / 10.1 cm | 2" / 5.08 cm | 24" / 609 mm |

PACKAGING: Bag of 3

AUTOCLAVE CLEANING BRUSH KITS



ROUND BRUSH

- Features 3 holes for adjustable pole placement
- Adjustable brush placement (parallel or perpendicular) and longer brush length provides easy access to hard-to-reach corners
- Made with FDA approved polyester bristles

HAND BRUSH

- Ergonomic handle for better hand comfort
- Longer, FDA approved polyester bristles
- Easy to clean (no cavity)

PRODUCT SELECTION

REGULAR AUTOCLAVE CLEANING BRUSH KIT

| ITEM # (KIT) | DESCRIPTION | ITEM # (UNIT) | DESCRIPTION | DIMENSIONS |
|--------------|--------------------------------------|---------------|-------------------------------------|--|
| 0-CLAVEREG | Autoclave REGULAR Brush Cleaning Kit | 50-35900 | 1 Hand Brush | Length: 3.5" x 8" / 8.89 cm x 20.32 cm |
| | | 50-55800 | 1 Rounded Corner Brush (Head Brush) | Head size: 5.5 x 10" / 14 x 25.4 cm |
| | | 50-EXTPOLE | 1 Anodized Aluminum Extension Pole | 2-4 ft / 60-122 cm |

SUPERIOR AUTOCLAVE CLEANING BRUSH KIT

| ITEM # (KIT) | DESCRIPTION | ITEM # (UNIT) | DESCRIPTION | DIMENSIONS | CASE DIMENSIONS |
|--------------|---------------------------------------|---------------|-------------------------------------|---------------------------------------|---------------------------------|
| 0-CLAVESUP | Autoclave SUPERIOR Brush Cleaning Kit | 50-35900 | 1 Hand Brush | Length: 3.5" x 8" / 11.5 cm x 20.32cm | |
| | | 50-55800 | 1 Rounded Corner Brush (Head Brush) | Head size: 5.5 x 10" / 14 x 25.4 cm | |
| | | 50-EXTPOLE | 1 Anodized Aluminum Extension Pole | 2-4 ft / 60-122 cm | |
| | | SF-240 | 2 SAVEFACE™ Face Shields | 9" / 23 cm | CASE OF 240 / 26.5" x 21" x 15" |

REPLACEMENT COMPONENTS



HAND BRUSH
50-35900

Packaging: 1



BRUSH HEAD
50-55800

Packaging: 1



EXTENSION POLE
50-EXTPOLE

Packaging: 1



SAVEFACE™
FACE SHIELDS SF-240



Packaging: 100 / Box




SILICONE PIN MAT

Our Silicone Pin Mat cushions, protects, and organizes delicate surgical instruments for sterilization. The elevated nubs provide greater air-flow ventilation and drying, significantly reducing residual moisture collection.

- Reusable
- Autoclavable
- Latex free



FRONT SIDE
Pins keep instruments securely in place



BACK SIDE
Numerous holes and elevated nubs provide superior venting and reduce water collection

PRODUCT SELECTION

| ITEM # | DIMENSIONS | PACKAGING |
|----------------|---|-----------|
| 13-2001 | 13" x 18" (33 cm x 45.72 cm) 11.97 mm high from base of mat to end of pin 15.53 mm high from base of elevated nub to end of pin | 1 unit |
| 13-2005 | 13" x 18" (33 cm x 45.72 cm) 11.97 mm high from base of mat to end of pin 15.53 mm high from base of elevated nub to end of pin | 5 units |

GOOD TO KNOW:
CUSTOM SIZES
ARE AVAILABLE
UPON REQUEST



DISTRIBUTED BY:
ADVANCED SURGI-PHARM INC.
www.advancedsurgipharm.com
Tel:1-800-661-5432


TAG-IT™ | NOTES

WHEN SIMPLY BROKEN DOESN'T CUT IT—TAG IT, TELL IT, OWN IT!

Tag-It™ Notes offer the freedom to **identify and document instructions for instruments, trays, and carts** before and after **sterilization**. Tag-It™ Notes are **tear-resistant, withstand extreme chemicals, high temperature** and ensures **reliable communication**.

WHY YOU'LL LOVE THEM:

- Tag-It™ Notes gives you the **freedom** to clearly document important information **without the clutter**.
- Robust and easy-to-use **tie wraps** designed to **securely attach** to instruments (134°C & chemical tested).
- Includes a **fast-drying, heat-resistant permanent marker** for writing your important notes.
- Choose from other **customizable cards** for personalized needs or **pre-made cards** for fast, convenient use.
- When you need more than a "repair" or "sharpen" seal to identify an issue. Instead, **Tag it, Tell it, Own it.**




WATERPROOF

The synthetic paper withstands spills and liquid exposure without any damage.




ECO-FRIENDLY

Made from sustainable materials.



HEAT RESISTANT

Can resist high temperatures of up to 273.2°F (134°C). Compatible with Gas Plasma Sterilization (STERAD®100 NX Sterilization).



DURABLE & TEAR-PROOF

Handles rigorous cleaning, tough conditions, and resists tearing.

PRODUCT SELECTION

| ITEM # | COLOR | SIZE | PACKAGING |
|------------|--------|-----------------------|-----------|
| TI-NOTE-R0 | Yellow | 3" x 5" / 8 cm x 13cm | 100 / PKG |

All Tag-It™ are single-use. **PACKAGING: 1 pack of 100 tags, 1 permanent marker, and 100 ties.**



Tag-It™ Cards are an effortless way to add information with either pre-printed or blank formats. These **durable** tags streamline your workow, **withstand tough conditions**, and feature **(optional) indicator dots** for quick verification. They help keep **instruments, trays, and carts** organized, ensuring efficiency before and after sterilization/decontamination.

KEY PRODUCT FEATURES:

- **Special Permanent Marker:** The ink will not run or appear blurry after exposure to chemicals or high temperature (max. 260°C/500°F) sterilization. *(Regular pens or markers can be used only if the tags are not going to be exposed to any sterilization/decontamination).*
- **Tie Wraps:** Assist in securing the tags to almost any instrument or container. *(Tie wraps have been tested for chemical and steam sterilization).*
- **Our Tag-It™:** Will not shrink, absorb water/chemicals, or lose its intended integrity after chemical or steam sterilization.



(Optional)
All Tags Available with Indicator Dots
 Provides visual confirmation that the tag was processed.

[CLICK BELOW TO WATCH OUR VIDEO](#)



PRODUCT SELECTION:

| ITEM # | COLOR | DESCRIPTION | SIZE | PACKAGING |
|----------|------------|------------------------|------------------------------------|-----------|
| TI-3502 | Blue | VENDOR | 3" x 5" / 8 cm x 13cm | 100 / PKG |
| TI-3503 | Green | LOANER | 3" x 5" / 8 cm x 13cm | 100 / PKG |
| TI-3504 | Red | PRIORITY | 3" x 5" / 8 cm x 13cm | 100 / PKG |
| TI-3505 | Yellow | RE-RUN | 3" x 5" / 8 cm x 13cm | 100 / PKG |
| TI-3506 | Pink | MISSING | 3" x 5" / 8 cm x 13cm | 100 / PKG |
| TI-3507 | Orange | BROKEN | 3" x 5" / 8 cm x 13cm | 100 / PKG |
| TI-3509 | Black | BLANK | 3" x 5" / 8 cm x 13cm | 100 / PKG |
| TI-35B1 | White | BLANK | 3" x 5" / 8 cm x 13cm | 100 / PKG |
| TI-35B4 | Red | BROKEN | 3" x 5" / 8 cm x 13cm | 100 / PKG |
| TI-35C22 | Light Blue | VENDOR CHECK IN + OUT | 3" x 5" / 8 cm x 13cm | 100 / PKG |
| TI-35L3 | Green | LOANER CHECK IN + OUT | 3" x 5" / 8 cm x 13cm | 100 / PKG |
| TI-35T2 | Blue | QUICK TURNOVER | 3" x 5" / 8 cm x 13cm | 100 / PKG |
| TI-PM09 | Black | FINE PERMANENT MARKERS | 1.5" Ø x 5.5" / 3.81cm Ø x 12.7 cm | 36 / Bag |

PACKAGING: 1 pack of 100 tags, 100 tie wraps, and 1 permanent marker.

OPTIONAL | PROCESSING INDICATOR DOTS: Add suffix "ST" to end of code for steam dots, "PS" for plasma dots.

Available in other languages or customize your own.

PRODUCT SELECTION: (customization also available)

VENDOR / BLUE

| ITEM # | SIZE |
|---------|--------------------------|
| TI-3502 | 3" x 5" (8 cm x 13cm) |

NOTES / YELLOW

| ITEM # | SIZE |
|------------|--------------------------|
| TI-NOTE-RO | 3" x 5" (8 cm x 13cm) |

LOANER / GREEN

| ITEM # | SIZE |
|---------|--------------------------|
| TI-3503 | 3" x 5" (8 cm x 13cm) |

BLANK / WHITE

| ITEM # | SIZE |
|---------|--------------------------|
| TI-35B1 | 3" x 5" (8 cm x 13cm) |

PRIORITY / RED

| ITEM # | SIZE |
|---------|--------------------------|
| TI-3504 | 3" x 5" (8 cm x 13cm) |

BROKEN / RED

| ITEM # | SIZE |
|---------|--------------------------|
| TI-35B4 | 3" x 5" (8 cm x 13cm) |

RE-RUN / YELLOW

| ITEM # | SIZE |
|---------|--------------------------|
| TI-3505 | 3" x 5" (8 cm x 13cm) |

VENDOR IN + OUT / LIGHT BLUE

| ITEM # | SIZE |
|----------|---------------------------|
| TI-35C22 | 5" x 3" (13 cm x 8 cm) |

MISSING / PINK

| ITEM # | SIZE |
|---------|--------------------------|
| TI-3506 | 3" x 5" (8 cm x 13cm) |

LOANER IN + OUT / GREEN

| ITEM # | SIZE |
|---------|---------------------------|
| TI-35L3 | 5" x 3" (13 cm x 8 cm) |

BROKEN / ORANGE

| ITEM # | SIZE |
|---------|--------------------------|
| TI-3507 | 3" x 5" (8 cm x 13cm) |

BLANK / BLACK

| ITEM # | SIZE |
|---------|--------------------------|
| TI-3509 | 3" x 5" (8 cm x 13cm) |

QUICK TURNOVER / BLUE

| ITEM # | SIZE |
|---------|--------------------------|
| TI-35T2 | 3" x 5" (8 cm x 13cm) |

All information tags are single-use. PACKAGING: 1 pack of 100 tags, 100 tie wraps, and 1 permanent marker.

TAG-IT™ | SCOPE

Batrik Scope Tags are tear-resistant, double-sided tags that will prove indispensable throughout the reprocessing journey. Our tags make it easy to record important information, including processing dates, with a user-friendly feature that allows you to select the month, day, and date. Crafted from waterproof synthetic paper, these tags feature a self-looping design, making it effortless to attach them to the scope without the need for additional fastening materials. Elevate your reprocessing protocols with Batrik!

KEY FEATURES

- The paper, tie wraps, and smudge-proof ink are designed to withstand high temperatures of up to 273.2°F (134°C) and are compatible with Gas plasma sterilization (STERRAD®100 NX Sterilization System).
- Tags include the month, day, and date selection on one side, making it convenient to mark the processing date accurately



COLOR-CODED AND CUSTOMIZABLE

- Available in seven colors
- The tags facilitate color coding for easy identification

SELF LOOPING

- The self-looping design allows for easy attachment to the scope
- Eliminating the need for additional fastening materials


SINGLE-USE AND LATEX FREE

- Batrik Scope Tags are designed for single use and are latex-free
- Promoting hygiene and reducing the risk of cross-contamination

PRODUCT SELECTION


| ITEM # | COLOR | SIZE | PACKAGING |
|---------|--------|--------------------------------|-----------|
| TI-ST05 | Yellow | 11.8" x 1.55" / 33.8 cm x 4 cm | 100 / PKG |
| TI-ST02 | Blue | 11.8" x 1.55" / 33.8 cm x 4 cm | 100 / PKG |
| TI-ST03 | Green | 11.8" x 1.55" / 33.8 cm x 4 cm | 100 / PKG |
| TI-ST07 | Orange | 11.8" x 1.55" / 33.8 cm x 4 cm | 100 / PKG |
| TI-ST06 | Purple | 11.8" x 1.55" / 33.8 cm x 4 cm | 100 / PKG |
| TI-ST04 | Red | 11.8" x 1.55" / 33.8 cm x 4 cm | 100 / PKG |
| TI-ST01 | White | 11.8" x 1.55" / 33.8 cm x 4 cm | 100 / PKG |

All information tags are single-use. PACKAGING: 1 pack of 100 tags and 1 permanent marker.




TAG-IT™ | SCOPE REPROCESSING

Batrik Scope Reprocessing Tags are durable and can withstand handling during reprocessing **without tearing**, ensuring that information remains intact. These tags are vital in recording important information such as **date to reprocess on, scope model and serial number, sterilizer type, meeting sterilizer parameter conditions, alcohol flush, and hang date**. These tags have been meticulously designed to elevate **operational efficiency** and contribute to the delivery of the highest standard of patient care.



TI-RST



TI-MPB

*Suitable for various information tagging purposes

LOGGING IMPORTANT INFORMATION

- Easy for staff to indicate when the scope was processed
- Key scope, technician, and cycle-related information easily reported and identified.

SMART DESIGN

- The self-looping design allows for easy attachment to the scope
- No requirement for extra fastening materials
- Secure fastening as the loop is tightened around the instrument
- Features an integrated barcode area




PRODUCT SELECTION

| ITEM # | DESCRIPTION | DIMENSIONS | PACKAGING |
|--------|-------------------------|---------------------------|-----------|
| TI-RST | REPROCESSING SCOPE TAG | 3" x 10.5" / 8 cm x 27 cm | 100 / PKG |
| TI-MPB | MULTI-PURPOSE BLANK TAG | 3" x 10.5" / 8 cm x 27 cm | 100 / PKG |


All information tags are single-use. **PACKAGING: 1 pack of 100 tags and 1 permanent marker.**

CORNER POCKETS™

With three models to choose from, our Corner Pockets™ portfolio is a truly comprehensive line of products designed to prevent tearing associated with sharp corners or feet on surgical trays.



18-CP02



18-CP07




18-CP06

| COLOUR | BLUE | ORANGE | PINK |
|-----------------------|---------------------------------|--|--|
| USAGE | SINGLE-USE | REUSABLE <i>withstands repeated sterilization up to 134°C</i> | REUSABLE <i>withstands repeated sterilization up to 134°C</i> |
| DIMENSIONS | 2.8" x 4" (7.3 cm x 10.3 cm) | 2.8" x 4" (7.3 cm x 10.3 cm) | 3.4" x 4.7" (8.5 cm x 12 cm) |
| LATEX-FREE | ✓ | ✓ | ✓ |
| NO MOISTURE RETENTION | ✓ | ✓ | ✓ |
| FABRICATION | Flexible, silicone-like plastic | Flexible, silicone-like plastic | Our most flexible plastic for easiest placement and most secure protection |

SOFT POCKETS™

Our Soft Pockets™ Foam Corner Protectors are made from high density foam to effectively cushion and protect the corners of surgical trays to prevent tearing. They are especially useful with heavier instrument sets.



PRODUCT SELECTION

| ITEM # | DESCRIPTION | PACKAGING | CASE DIMENSIONS |
|------------|---|---------------------------------|------------------------------|
| 18-CP02** | Corner Pockets™ Single-Use, Blue | 240 / Bag | 80 / Box (24 Boxes / Pallet) |
| 18-EVAL2** | Corner Pockets™ Single-Use, Blue | 24 / Bag | - |
| 18-CP07 | Corner Pockets™ Reusable, Orange | 80 / Bag | - |
| 18-CORNER | Corner Pockets™ Reusable, Orange | 320 / Box (4 Bags of 80 pcs) | - |
| 18-EVAL7 | Corner Pockets™ Reusable, Orange | 24 / Bag | - |
| 18-CP06 | Corner Pockets™ Reusable, Pink | 80 / Bag | - |
| 18-CPF-02 | Soft Pockets™ Single-Use Foam Corner Protectors | 80 / Box (24 Boxes / Pallet) | - |

** Although the single use Corner Pockets™ appear similar to the reusable models they are blended differently and have only been validated as single-use. Re-using them may compromise sterility.

LILA'S INSTRUMENT STAIN REMOVER

Lila's special blend of naturally sourced powders provides outstanding removal of stubborn stains.

Convenient-to-use formula, you can achieve spotless results every time!

It can be used on surgical/dental instruments, autoclaves, ultrasonic cleaners, stainless steel countertops, trays, bowls and more.

Specifically designed for use on stainless steel, chrome, nickel, brass, copper, and silver.



LSR-8501

EASE TO USE AS 1,2,3...

1. Dampen a cloth or grab a Lila's Low Linting Dura-Towel and dab into the jar.
2. Gently rub the area that is stained.¹
3. Rinse thoroughly for a flawless shine.²

A little Goes a long Way!

Even great to tackle stubborn stains on steel pots and pans in your kitchen with ease.

Lila's - If you have used Spectra-Scrub or Integra's instrument stain remover then you will love Lila's!



60-46003



17-DT1213C

LILA'S SCRUB PAD

- Say goodbye to extra stubborn stains with Lila's Scrub Pad without leaving any metal residue.
- Simply add a small amount of Lila's Stain Remover with one of our Scrub Pad's for very large areas.
- The gentle yet effective scrub is the perfect choice for Ophthalmic Instruments.¹

LILA'S LOW-LINTING DURA-TOWEL

- Lila's Dura-Towel has long lasting strength and has a minimal lint composition.
- These single-use towels are clean room certified and are specifically designed for instrument wiping.
- Its ideal dimensions and unique attributes guarantee optimal performance, setting it apart from ordinary paper towels.

PRODUCT SELECTION


| ITEM # | DESCRIPTION | DIMENSIONS | PACKAGING |
|-------------------|--|--|---------------------------|
| LSR-8501 | LILA'S INSTRUMENT STAIN REMOVER 3 oz (85g) plastic jar with a convenient screw-on lid | Ø 2.8" / 70 mm | 1 Jar / Bag |
| LSR-8512 | LILA'S INSTRUMENT STAIN REMOVER 3 oz (85g) plastic jar with a convenient screw-on lid | Ø 2.8" / 70 mm | 12 Jars / Box |
| 17-DT1213C | LILA'S LOW LINTING DURA-TOWEL 1080 towels | 12" x 13" / 31 cm x 33 cm Folded: 3" x 3.25" / 7.5 cm x 8.25 cm | 90 / Bag 12 Bags / Box |
| 60-46003 | LILA'S SCRUB PAD | 4" x 6" / 10 cm x 15 cm | 3 / Bag |
| LSR-BUNDLE | LILA'S STAIN REMOVER + LILA'S DURA-TOWEL 3 oz (85g) plastic jar + 90 towels | Ø 2.8" / 70 MM 12" x 13" / 31 cm x 33 cm | 1 Jar / Bag 90 / Bag |

1 : Although Lila's is completely safe, it is not recommended for use on ophthalmic instruments. Use Lila's Scrub Pad on its own.
2: The stain remover should not be mixed with bleach and is for external use only to remove instrument stains.

INSTRUMENT PROTECTORS

Our Instrument Protectors protect the tips and edges of instruments to prevent damage during handling and sterilization. Available in a wide array of color coded sizes for almost every application, including osteotome, retractor, double, and perforated.

- Steam and gas permeable
- Perforated protectors reduce vacuum seal and moisture build-up
- Validated for pre-vacuum sterilization*
- Multiple color-coded sizes
- Latex Free
- Non-Sterile



*Independent laboratory tested to validate pre-vacuum steam sterilization efficacy in a healthcare setting

PRODUCT SELECTION

REGULAR & PERFORATED / CLEAR


RC-07

PC-07

| ITEM # REGULAR | ITEM # PERFORATED | COLOR | DIMENSIONS |
|-------------------|----------------------|-------|---|
| RC-01 | PC-01 | Clear | 0.062" x 0.75" (1.57 mm x 19 mm) |
| RC-02 | PC-02 | Clear | 0.075" x 0.75" (2 mm x 19 mm) |
| RC-03 | PC-03 | Clear | 0.109" x 0.75" (2.77 mm x 19 mm) |
| RC-04 | PC-04 | Clear | 0.125" x 1" (3.18 mm x 25.4 mm) |
| RC-05 | PC-05 | Clear | 0.187" x 1" (4.75 mm x 25.4 mm) |
| RC-06 | PC-06 | Clear | 0.406" x 0.75" (10.3 mm x 19 mm) |
| RC-07 | PC-07 | Clear | 0.062" x 0.375" x 1" (1.57 mm x 9.53 mm x 25.4 mm) |
| RC-08 | PC-08 | Clear | 0.078" x 0.625" x 1" (2 mm x 15.875 mm x 25.4 mm) |
| RC-09 | PC-09 | Clear | 0.125" x 1" x 1" (3.18 mm x 25.4 mm x 25.4 mm) |
| RC-10 | PC-10 | Clear | Assorted |

PACKAGING: Bag of 100

REGULAR & PERFORATED / TINTED


RT-07

PT-07

| | | | |
|--------------|--------------|-----------------|---|
| RT-01 | PT-01 | Tinted White | 0.062" x 0.75" (1.57 mm x 19 mm) |
| RT-02 | PT-02 | Tinted Blue | 0.075" x 0.75" (2 mm x 19 mm) |
| RT-03 | PT-03 | Tinted Green | 0.109" x 0.75" (2.77 mm x 19 mm) |
| RT-04 | PT-04 | Tinted Red | 0.125" x 1" (3.18 mm x 25.4 mm) |
| RT-05 | PT-05 | Tinted Yellow | 0.187" x 1" (4.75 mm x 25.4 mm) |
| RT-06 | PT-06 | Tinted Maroon | 0.406" x 0.75" (10.3 mm x 19 mm) |
| RT-07 | PT-07 | Tinted Orange | 0.062" x 0.375" x 1" (1.57 mm x 9.53 mm x 25.4 mm) |
| RT-08 | PT-08 | Tinted Brown | 0.078" x 0.625" x 1" (2 mm x 15.875 mm x 25.4 mm) |
| RT-09 | PT-09 | Tinted Black | 0.125" x 1" x 1" (3.18 mm x 25.4 mm x 25.4 mm) |
| RT-10 | PT-10 | Tinted Assorted | Assorted |

PACKAGING: Bag of 100

REGULAR & PERFORATED / SOLID


RS-07

PS-07

| | | | |
|--------------|--------------|----------------|---|
| RS-01 | PS-01 | Solid White | 0.062" x 0.75" (1.57 mm x 19 mm) |
| RS-02 | PS-02 | Solid Blue | 0.075" x 0.75" (2 mm x 19 mm) |
| RS-03 | PS-03 | Solid Green | 0.109" x 0.75" (2.77 mm x 19 mm) |
| RS-04 | PS-04 | Solid Red | 0.125" x 1" (3.18 mm x 25.4 mm) |
| RS-05 | PS-05 | Solid Yellow | 0.187" x 1" (4.75 mm x 25.4 mm) |
| RS-06 | PS-06 | Solid Maroon | 0.406" x 0.75" (10.3 mm x 19 mm) |
| RS-07 | PS-07 | Solid Orange | 0.062" x 0.375" x 1" (1.57 mm x 9.53 mm x 25.4 mm) |
| RS-08 | PS-08 | Solid Brown | 0.078" x 0.625" x 1" (2 mm x 15.875 mm x 25.4 mm) |
| RS-09 | PS-09 | Solid Black | 0.125" x 1" x 1" (3.18 mm x 25.4 mm x 25.4 mm) |
| RS-10 | PS-10 | Solid Assorted | Assorted |

PACKAGING: Bag of 100

| ITEM # | DIMENSIONS |
|--------------|-------------------------------------|
| OP-11 | 0.25" x 0.75" (6.35 mm x 19.05 mm) |
| OP-12 | 0.375" x 1" (9.53 mm x 25.4 mm) |
| OP-13 | 0.5" x 1" (12.7 mm x 25.4 mm) |
| OP-14 | 0.75" x 1" (19.05 mm x 25.4 mm) |
| OP-15 | 1" x 1" (25.4 mm x 25.4 mm) |
| OP-16 | 1.25" x 1.25" (31.75 mm x 31.75 mm) |
| OP-17 | 1.5" x 1.5" (38.1 mm x 38.1 mm) |
| OP-18 | 1.75" x 1.25" (44.45 mm x 31.75 mm) |
| OP-20 | Assorted |

PACKAGING: Bag of 50

| ITEM # | DIMENSIONS | PACKAGING |
|---------------|--------------|-----------|
| RGS-01 | 1.5" (38 mm) | 55 / Bag |
| RGS-02 | 2" (51 mm) | 45 / Bag |
| RGS-03 | 3" (76.2 mm) | 40 / Bag |


| ITEM # | COLOR | SIZE |
|---------------|-------|--------|
| RGM-01 | Pink | Small |
| RGM-02 | Clear | Small |
| RGM-03 | Green | Medium |
| RGM-04 | Clear | Medium |

PACKAGING: Bag of 50

| ITEM # | COLOR | DIMENSIONS |
|--------------|----------|---|
| DG-01 | Clear | 0.0625" x 0.75" (1.59 mm x 19.05 mm) |
| DG-02 | Blue | 0.078" x 0.75" (1.98 mm x 19.05 mm) |
| DG-03 | Green | 0.109" x 0.75" (2.78 mm x 19.05 mm) |
| DG-04 | Red | 0.125" x 1" (3.18 mm x 25.4 mm) |
| DG-05 | Yellow | 0.187" x 1" (4.76 mm x 25.4 mm) |
| DG-07 | Orange | 0.062" x 0.375" x 1" (1.59 mm x 9.53 mm x 25.4 mm) |
| DG-10 | Assorted | Assorted |

PACKAGING: Bag of 50

OSTEOTOME PROTECTORS




SINGLE-TIP RETRACTOR PROTECTORS



WEITLANER RETRACTOR PROTECTORS



PERFORATED DOUBLE PROTECTORS



PrAmLi Power³ LUBE™ Instrument Lubricant

Specially formulated to lubricate surgical instruments, hinges, box locks, and general equipment during autoclaving. PrAmLi Power³ LUBE™ is uniquely formulated with an all natural corrosion inhibitor.

- Forms a protective coating which can extend the instrument's, useful life
- Bacteriostatic
- Neutral pH 6.5 - 7.5
- Prevents corrosion

TIP:
FOR MILK BATH SOLUTION, DILUTE 1 PART SOLUTION IN 6 PARTS WATER



PRODUCT SELECTION

| ITEM # | DILUTION LUBRICANT / WATER | PACKAGING | CASE DIMENSIONS |
|------------|---|-----------------------------------|---------------------|
| SS-21-6001 | 1 to 4 ml / L (0.125 to 0.5 oz / Gallon) | 1 x 3.8 L (1 Gallon) | 6" x 6" x 13" |
| SS-21-6004 | 1 to 4 ml / L (0.125 to 0.5 oz / Gallon) | 4 x 3.8 L (4 Gallons) | 13" x 13" x 14" |
| SS-21-6025 | 1 to 4 ml / L (0.125 to 0.5 oz / Gallon) | 2 x 9.5 L (2 x 2.5 Gallons) | 14.5" x 9.25" x 15" |
| SS-21-6032 | No dilution required | 1 x 946 ml (32 oz) Spray bottle | 10" x 4" |
| SS-21-6012 | No dilution required | 12 X 946 ml (32 oz) Spray Bottles | 18" x 12" x 11" |


S.T.I.C.S.

Our Simple Tracking Identification Color System, or S.T.I.C.S., is an easy to use, latex free, tool that allows you to quickly identify your instruments by department and/or procedure. The tape adheres securely and every tape sheet comes with a protective cover to ensure it lasts as long as you need.



Each 8.5" x 11" (216 mm x 279 mm) sheet contains the following adhesive strips:

- 1" (25.4 mm) x 220
- 1.5" (38.1 mm) x 44
- 2" (50.8 mm) x 44
- Total tape per sheet 374" (9500 mm)



PRODUCT SELECTION

SOLID COLORS

| ITEM # | COLOR | |
|------------|-------------|--|
| T-850-2-01 | White | |
| T-850-2-02 | Royal Blue | |
| T-850-2-03 | Kelly Green | |
| T-850-2-04 | Red | |
| T-850-2-05 | Yellow | |
| T-850-2-06 | Purple | |
| T-850-2-07 | Orange | |
| T-850-2-08 | Brown | |
| T-850-2-09 | Black | |
| T-850-2-10 | Gray | |
| T-850-2-11 | Lavender | |
| T-850-2-12 | Beige | |
| T-850-2-13 | Light Blue | |
| T-850-2-14 | Lime Green | |
| T-850-2-15 | Pink | |
| T-850-2-17 | Plum | |
| T-850-2-18 | Mint Green | |

STARS

| ITEM # | COLOR | |
|------------|-------------------|--|
| T-850-6-01 | Stars Fuchsia | |
| T-850-6-02 | Stars Medium Blue | |
| T-850-6-03 | Stars Green | |
| T-850-6-04 | Stars Violet | |

TRIANGLES

| ITEM # | COLOR | |
|-------------|------------------|--|
| T-850-10-01 | Triangle Fuchsia | |
| T-850-10-04 | Triangle Violet | |

WHITE STRIPES

| ITEM # | COLOR | |
|------------|-------------------|--|
| T-850-4-02 | White/Royal Blue | |
| T-850-4-03 | White/Kelly Green | |
| T-850-4-04 | White/Red | |
| T-850-4-05 | White/Yellow | |
| T-850-4-07 | White/Orange | |
| T-850-4-08 | White/Brown | |
| T-850-4-11 | White/Lavender | |
| T-850-4-13 | White/Light Blue | |
| T-850-4-14 | White/Lime Green | |
| T-850-4-15 | White/Pink | |
| T-850-4-17 | White/Plum | |

BLACK STRIPES

| ITEM # | COLOR | |
|------------|-------------------|--|
| T-850-3-01 | Black/White | |
| T-850-3-03 | Black/Kelly Green | |
| T-850-3-04 | Black/Red | |
| T-850-3-05 | Black/Yellow | |
| T-850-3-07 | Black/Orange | |
| T-850-3-10 | Black/Gray | |
| T-850-3-11 | Black/Lavender | |
| T-850-3-12 | Black/Beige | |
| T-850-3-13 | Black/Light Blue | |
| T-850-3-14 | Black/Lime Green | |
| T-850-3-15 | Black/Pink | |
| T-850-3-17 | Black/Plum | |

DOTS

| ITEM # | COLOR | |
|------------|--------------|--|
| T-850-1-01 | Dots Fuchsia | |
| T-850-1-03 | Dots Green | |

FOUR-COLOR

| ITEM # | COLOR | |
|------------|--------------------------|--|
| T-850-2-31 | Red/Yellow/Purple/Orange | |
| T-850-2-32 | Blue/White/Black/Green | |

TRI-COLOR

| ITEM # | COLOR | |
|------------|-----------------|--|
| T-850-5-51 | Red/White/Blue | |
| T-850-5-52 | Red/White/Black | |
| T-850-5-54 | Red/Black/Green | |

STRIPES

| ITEM # | COLOR | |
|------------|--------------|--|
| T-850-7-01 | Blue/Green | |
| T-850-7-02 | Blue/Red | |
| T-850-7-03 | Yellow/Blue | |
| T-850-7-04 | Yellow/Red | |
| T-850-7-05 | Yellow/Green | |

SPECIALTY

| ITEM # | COLOR | |
|------------|------------------------|--|
| T-850-8-01 | Red Hearts/White | |
| T-850-8-02 | Smiley Yellow/Black | |
| T-850-8-03 | Alphabet Fuchsia/White | |
| T-850-8-04 | White Hearts/Red | |
| T-850-9-03 | Diagonal Blue/Black | |
| T-850-9-04 | Diagonal Red/Black | |

NUMBERS


| ITEM # | COLOR | |
|-------------|-----------------|--|
| T-850-11-01 | Numbers Fuchsia | |

BATRIK ID TAPE

A simple solution for marking surgical instruments with the aid of color-coding, our ID Tape is durable and will not chip, flake, or peel. It can be used in autoclaves, withstanding repeated exposure to sterilization. Independent validation tests demonstrated the efficiency of sterilization of marked instruments, even underneath the BATRIK ID Tape.



AVAILABLE WIDTHS



Quickly identifies instruments by department or procedure



Easy to apply



Withstands repeated exposure to sterilization



Has more per tape roll (300" / 7620 mm)



Available in 7 colors and 2 widths



Conveniently packaged in a clear box



Non-Sterile

PRODUCT SELECTION

| ITEM # | ITEM# | COLOR |
|-----------------|----------------|--------|
| 0.125" / 318 MM | 0.25" / 635 MM | |
| T-125-01 | T-250-01 | White |
| T-125-02 | T-250-02 | Blue |
| T-125-03 | T-250-03 | Green |
| T-125-04 | T-250-04 | Red |
| T-125-05 | T-250-05 | Yellow |
| T-125-07 | T-250-07 | Orange |
| T-125-09 | T-250-09 | Black |

TAPE-B-GONE™ TAPE REMOVER

Tape-B-Gone™ is a uniquely formulated and ready-to-use proprietary formula that easily removes adhesives from medical equipment and instruments. Unlike many other brands, Tape-B-Gone™ does not contain highly dangerous or caustic chemicals, is not flammable or combustible, and is safe for use on metal, glass, plastic, and rubber without any special equipment or ventilation.



PRODUCT SELECTION

| ITEM # | DESCRIPTION | PACKAGING |
|---------------------|-------------------------------|-----------|
| T-08RMV | 8 oz (237 ml) bottle | 1 Unit |
| T-32RMV | 8 oz (237 ml) bottle | 4 / Pack |
| T-SAMP-RMV | 0.5 oz (15 ml) bottle | 1 Unit |
| T-96RMVSC-R0 | 8OZ (237 ML) BOTTLE SCREWCAP* | 12 / Box |

*Samples will be available in Batrik by April 21st

Ver.2021-04-22

DISTRIBUTED BY:
 ADVANCED SURGI-PHARM INC.
 www.advancedsurgipharm.com
 Tel.1-800-661-5432



PADLOCKS SECURITY TAGS

Our Padlock Security Tags for instrument containers and general hospital applications are specifically designed for the healthcare market. Our material selection, fabrication process, and engineering ensure consistent quality and high temperature compliance with the various steam sterilization temperatures in worldwide hospital applications.

Available with optional reactive ink indicator dots that provide visual sterility indication. Light blue dots turn brown/dark brown when steam sterilized.



Made of polypropylene



Latex and DEHP Free



Will not melt or soften with any established steam sterilization protocol




Fully open design (45°) facilitates insertion into instrument containers



Easy twist breaking design



Large base for easier handling



Available in different size




Custom Colors and OEM available upon request

AVAILABLE IN THREE FORMATS


700-SERIES



800-SERIES



900-SERIES



PRINTED & CUSTOM NUMBERED PADLOCKS

Available in 700, 800 & 900 SERIES

PRODUCT SELECTION

| COLOR | PRINTED WORD |
|-------------------------|-----------------|
| PRINTED PADLOCKS | |
| WHITE | CLEAN |
| BLUE | STERILE |
| RED | REPAIR |
| YELLOW | SHARPEN |
| YELLOW | LOANER |
| BLACK | DISCARD |
| CUSTOM NUMBERED | |
| ALL COLORS | up to 8 numbers |



CUSTOM NUMBERED
Custom printed padlocks (up to 8 characters)



PRINTED
6 different printed labels available

NOTE: Printed & Custom Numbered Padlocks are available in 700, 800, and 900 Series and can be ordered in 100, 200, and 500 unit quantities. Please contact us for ordering details.

INDICATOR DOTS

REGULAR


700-SERIES PADLOCKS

The 700 series padlocks are perfect for "Aesculap" containers and those with a short, large opening.

PRODUCT SELECTION

REGULAR PADLOCKS

| COLOR | 100 UNITS | 200 UNITS | 500 UNITS |
|-------------------|-----------|-----------|-----------|
| WHITE | 15-701 | 15-711 | 15-710 |
| BLUE | 15-702 | 15-712 | 15-720 |
| GREEN | 15-703 | 15-713 | 15-730 |
| RED | 15-704 | 15-714 | 15-740 |
| YELLOW | 15-705 | 15-715 | 15-750 |
| PINK | 15-706 | 15-716 | 15-760 |
| ORANGE | 15-707 | 15-717 | 15-770 |
| BLACK | 15-709 | 15-719 | 15-790 |
| LIGHT BLUE | 15-722 | 15-732 | 15-742 |

WITH INDICATORS DOTS

| COLOR | 100 UNITS | 200 UNITS | 500 UNITS |
|-------------------|------------|------------|------------|
| WHITE | 15-701-IND | 15-711-IND | 15-710-IND |
| BLUE | 15-702-IND | 15-712-IND | 15-720-IND |
| GREEN | 15-703-IND | 15-713-IND | 15-730-IND |
| RED | 15-704-IND | 15-714-IND | 15-740-IND |
| YELLOW | 15-705-IND | 15-715-IND | 15-750-IND |
| PINK | 15-706-IND | 15-716-IND | 15-760-IND |
| ORANGE | 15-707-IND | 15-717-IND | 15-770-IND |
| BLACK | 15-709-IND | 15-719-IND | 15-790-IND |
| LIGHT BLUE | 15-722-IND | 15-732-IND | 15-742-IND |

Double locking hooks can be inserted in a large guide hole



The grooves on the padlock facilitate handling



Ver.2021-04-22

DISTRIBUTED BY:
 ADVANCED SURGI-PHARM INC.
 www.advancedsurgipharma.com
 Tel.1-800-661-5432

INDICATOR DOTS



REGULAR



800-SERIES PADLOCKS

The 800 series padlocks have a longer arm (32 mm) and have a single, thinner hook than the 700 series. They are perfect for almost any type of container, even those with long and thin channels.

PRODUCT SELECTION

REGULAR PADLOCKS

| COLOR | 100 UNITS | 200 UNITS | 500 UNITS |
|------------|-----------|-----------|-----------|
| WHITE | 15-801 | 15-811 | 15-810 |
| BLUE | 15-802 | 15-812 | 15-820 |
| GREEN | 15-803 | 15-813 | 15-830 |
| RED | 15-804 | 15-814 | 15-840 |
| YELLOW | 15-805 | 15-815 | 15-850 |
| PINK | 15-806 | 15-816 | 15-860 |
| ORANGE | 15-807 | 15-817 | 15-870 |
| BLACK | 15-809 | 15-819 | 15-890 |
| LIGHT BLUE | 15-822 | 15-832 | 15-842 |

WITH INDICATORS DOTS

| COLOR | 100 UNITS | 200 UNITS | 500 UNITS |
|------------|------------|------------|------------|
| WHITE | 15-801-IND | 15-811-IND | 15-810-IND |
| BLUE | 15-802-IND | 15-812-IND | 15-820-IND |
| GREEN | 15-803-IND | 15-813-IND | 15-830-IND |
| RED | 15-804-IND | 15-814-IND | 15-840-IND |
| YELLOW | 15-805-IND | 15-815-IND | 15-850-IND |
| PINK | 15-806-IND | 15-816-IND | 15-860-IND |
| ORANGE | 15-807-IND | 15-817-IND | 15-870-IND |
| BLACK | 15-809-IND | 15-819-IND | 15-890-IND |
| LIGHT BLUE | 15-822-IND | 15-832-IND | 15-842-IND |

Single locking hook can be inserted in a large guide hole



The locking system allows for **easy closure**



**INDICATOR
DOTS**

REGULAR


900-SERIES PADLOCKS

The 900 series padlocks are identical to the 800 series, but with a larger rectangular base which allows for the placement of a 20 mm x 30 mm tracking label.

PRODUCT SELECTION

REGULAR PADLOCKS

| COLOR | 100 UNITS | 200 UNITS | 500 UNITS |
|-------------------|-----------|-----------|-----------|
| WHITE | 15-901 | 15-911 | 15-910 |
| BLUE | 15-902 | 15-912 | 15-920 |
| GREEN | 15-903 | 15-913 | 15-930 |
| RED | 15-904 | 15-914 | 15-940 |
| YELLOW | 15-905 | 15-915 | 15-950 |
| PINK | 15-906 | 15-916 | 15-960 |
| ORANGE | 15-907 | 15-917 | 15-970 |
| BLACK | 15-909 | 15-919 | 15-990 |
| LIGHT BLUE | 15-922 | 15-932 | 15-942 |

WITH INDICATORS DOTS

| COLOR | 100 UNITS | 200 UNITS | 500 UNITS |
|-------------------|------------|------------|------------|
| WHITE | 15-901-IND | 15-911-IND | 15-910-IND |
| BLUE | 15-902-IND | 15-912-IND | 15-920-IND |
| GREEN | 15-903-IND | 15-913-IND | 15-930-IND |
| RED | 15-904-IND | 15-914-IND | 15-940-IND |
| YELLOW | 15-905-IND | 15-915-IND | 15-950-IND |
| PINK | 15-906-IND | 15-916-IND | 15-960-IND |
| ORANGE | 15-907-IND | 15-917-IND | 15-970-IND |
| BLACK | 15-909-IND | 15-919-IND | 15-990-IND |
| LIGHT BLUE | 15-922-IND | 15-932-IND | 15-942-IND |

Single locking hook
for use with smaller
diameters



Ideal for adding labels/
indicators indicators or
markings to the large
rectangular base




Ver.2021-04-22

SCRUB PAD

Our Scrub Pad is ideal for pre-cleaning your surgical instruments without leaving any metal residue that often occurs when using steel wool or metal sponges. The abrasiveness is perfectly balanced to provide you with efficient cleaning without damaging your instruments. Its fine fibers allow the user to polish, clean, and debur all types of material surfaces.

- Abrasive Material: Aluminum Oxide
- Grit: 400
- Grade: Very Fine
- Style: Flexible
- Pad material: Non-Woven
- Color: Maroon
- Type: Hand Pad



PRODUCT SELECTION

| ITEM # | DESCRIPTION | DIMENSIONS | PACKAGING |
|----------|-------------|-------------------------|--|
| 60-46050 | Scrub Pad | 4" x 6" / 10 cm x 15 cm | 150 / Box (50 bags of 3 units each) |
| 60-46003 | Scrub Pad | 4" x 6" / 10 cm x 15 cm | 3 / Bag |

Ver.2021-04-22

Your First ™

ENZYMATIC & NON-ENZYMATIC CLEANERS

Our Enzymatic and Non-Enzymatic Cleaners are specially formulated to remove organic material from medical equipment and instrumentation.

Our **Enzymatic Cleaner** formats are formulated with PrAmLi Power3 Enzymes™, our proprietary unique blend of enzymes, and a surfactant to efficiently remove all proteins, blood, fatty tissue, and more. They also remove stains from surgical scrubs, linens, and medical clothing.

Our **Non-Enzymatic Cleaner** formats are formulated with a free-rinsing, neutral pH detergent blend for effective cleaning.



Both formats are powerful solutions that break down all bioburden for effective surface and instrument pre-cleaning*

PRODUCT SELECTION

| ENZYMATIC CLEANER | NON-ENZYMATIC CLEANER | DILUTION DETERGENT / WATER | PACKAGING | CASE DIMENSIONS |
|--------------------------|-----------------------|--|-------------------|-----------------|
| LOW FOAM | | | | |
| SS-21-1001 | SS-721-1001 | 4 to 8 ml / L (0.5 to 1 oz / Gal) | 1 x 3.8 L (1 Gal) | 6x6x13" |
| SS-21-1004 | SS-721-1004 | 4 to 8 ml / L (0.5 to 1 oz / Gal) | 4 x 3.8 L (4 Gal) | 13x13x14" |
| MILD FOAM | | | | |
| SS-21-2001 | SS-721-2001 | 4 to 8 ml / L (0.5 to 1 oz / Gal) | 1 x 3.8 L (1 Gal) | 6x6x13" |
| SS-21-2004 | SS-721-2004 | 4 to 8 ml / L (0.5 to 1 oz / Gal) | 4 x 3.8 L (4 Gal) | 13x13x14" |
| SUPER CONCENTRATE | | | | |
| SS-21-3001 | SS-721-3001 | 1 to 2 ml / L (0.125 to 0.25 oz / Gal) | 1 x 3.8 L (1 Gal) | 6x6x13" |
| SS-21-3004 | SS-721-3004 | 1 to 2 ml / L (0.125 to 0.25 oz / Gal) | 4 x 3.8 L (4 Gal) | 13x13x14" |

GOOD TO KNOW:

OUR LOW FOAM & SUPER CONCENTRATE DETERGENTS (ENZYMATIC & NON-ENZYMATIC) MAY BE USED IN WASHERS / ULTRASONIC CLEANERS.

WE RECOMMEND CONFIRMING SUGGESTED DILUTION RATIO AND WATER TEMPERATURE WITH THE MANUFACTURER OF YOUR WASHER / ULTRASONIC CLEANER AS MACHINES / SETTINGS VARY.

NOTE: Our Low Foam & Super Concentrate Detergent formulas are available in a 2.5 Gallon format based on a minimum order quantity. Contact us for more information.

*Caution: We recommend consulting your hospital pre-cleaning/cleaning protocols for various equipment/instrument cleaning detergent preferences.

Ver. 2022-05-31

Your First ™

ENVIRO-SACHET™

ENZYMATIC GEL & NON-ENZYMATIC CLEANER

Your First ™ Enviro-Sachet™, our unique cleaning solutions now available in Enzymatic & Non-Enzymatic formats.



- Our Enzymatic Gel format is formulated with the same triple action **PrAmLi Power³ Enzymes™** that you know and trust.
- Our Non-Enzymatic Cleaner format is formulated with a free-rinsing, neutral pH detergent blend.

Both formats are powerful solutions that break down all bioburden for effective surface and instrument pre-cleaning.

PRODUCT SELECTION

| ENZYMATIC GEL | NON-ENZYMATIC CLEANER | DESCRIPTION | PACKAGING | CASE DIMENSIONS |
|---------------|-----------------------|--|----------------------------|-------------------|
| 17-166 | 717-166 | 1 Dry Tubular Sponge & Enviro-Sachet™ set | 25 / Box 4 Boxes / Case | 21" x 14" x 13.5" |
| 17-266 | 717-266 | 1 Dry Flat Sponge & Enviro-Sachet™ set | 50 / Box 4 Boxes / Case | 21" x 14" x 13.5" |
| 17-177 | 717-177 | 1 Dry Tubular Sponge, Enviro-Sachet™ & Stainless Steel Bowl kit | 1 Kit | 7" x 6" x 3" |
| 17-277 | 717-277 | 1 Dry Flat Sponge, Enviro-Sachet™ & Stainless Steel Bowl kit | 1 Kit | 7" x 6" x 3" |
| 17-EZSAMP | 717-EZSAMP | SAMPLE KIT INCLUDES: 3 Dry Flat Sponges with Enviro-Sachet™ sets 1 Dry Tubular Sponge with Enviro-Sachet™ set 1 Reusable Stainless Steel Bowl | 1 Kit | 7" x 6" x 3" |
| H-EZ500 | | Reusable Stainless Steel Bowl 4.3" x 5.3" x 3" / 11 cm x 13.5 cm x 7.5 cm | 1 Unit | 7" x 6" x 3" |

Your First ™

FOAMING SPRAY

READY-TO-USE ENZYME FOAMING SPRAY

Our Foaming Spray is ideal for instrument pre-cleaning as the thick foam spray enzymatically breaks down all bioburden to allow for effective cleaning and sterilization.

Formulated with **PrAmLi Power³ Enzymes™**, a unique blend of three powerful enzymes along with a surfactant to provide the ultimate in cleaning efficacy



INSTRUCTIONS



PRODUCT SELECTION

| ITEM # | PACKAGING | CASE DIMENSIONS |
|-------------------|-----------------------------------|-----------------|
| SS-21-4001 | 1 x 946 ml (32 oz) Spray Bottle | 10" x 4" |
| SS-21-4012 | 12 x 946 ml (32 oz) Spray Bottles | 18" x 12" x 11" |

Important: In addition to our enzymatic formulation, our spray bottle contains as much as 25% more solution when compared to the competition.

DISTRIBUTED BY:
 ADVANCED SURGI-PHARM INC.
www.advancedsurgipharma.com
 Tel.1-800-661-5432



Your First ™

F1RST WIPES

Our cleanroom wipes are designed to do the first wipe-down of an endoscope after a procedure. They are made from hydroentangled non-woven material composed of 45% polyester and 55% cellulose. These wipes are manufactured and packaged in a cleanroom environment.



- Can be used for a cleanroom environment of ISO Class 5 through 8
- Production process guarantees extremely low particle generation and chemical extractables
- Compatible with virtually any solvent or solution for cleaning and wiping (such as **PrAmLi Power³ Enzymes™**)
- Double-bag packaging
- Available in two versatile sizes

PRODUCT SELECTION

| ITEM # | DESCRIPTION | DIMENSIONS | PACKAGING | CASE DIMENSIONS |
|-----------------|---------------------------|------------------------------|-----------------------------|-----------------|
| 17-A0909 | Non-woven Cleanroom Wipes | 9" x 9" (23 cm x 23 cm) | 300 / Bag 14 Bags / Case | 18.5x10x25.25" |
| 17-A1212 | Non-woven Cleanroom Wipes | 12" x 12" (31 cm x 31 cm) | 150 / Bag 10 Bags / Case | 13x13x19" |

Your First ™ DRY SPONGES

Our sponges are specifically fabricated for use with all medical instrumentation. Unlike most competitor sponges, our sponges are white to better evidence the effectiveness of the cleaning process. Dark colored sponges easily hide contaminants, but with our white sponges you can visually evaluate the pre-clean efficacy.



Available in two formats, Flat and Tubular:

Eco-friendly

100% medical grade high-density foam

Compatible with enzymatic and disinfectant solutions

Will not flake or scratch expensive equipment

TIP: OUR TUBULAR SPONGE CAN ALSO BE USED AS A BUMPER TO PROTECT YOUR SCOPES

High absorbency retains more solution than conventional sponges or paper wipes

OUR WHITE SPONGES MAY APPEAR SLIGHTLY YELLOWED. THIS IS ONLY DUE TO UV EXPOSURE, IS NOT A REFLECTION OF ITS AGE.

PRODUCT SELECTION

| ITEM # | DESCRIPTION | DIMENSIONS | PACKAGING |
|--------|---|---|------------------------------|
| 17-100 | DRY TUBULAR SPONGE White Dispenser Box | 2.0" ø x 4.25" (2.7 cm ø x 10.8 cm) | 25 / Box (4 Boxes / Case) |
| 17-101 | DRY TUBULAR SPONGE Bulk | 2.0" ø x 4.25" (2.7 cm ø x 10.8 cm) | 560 / Case |
| 17-200 | DRY FLAT SPONGE White Dispenser Box | 4" x 5" x 0.5" (10.2 cm x 12.7 cm x 1.3 cm) | 50 / Box (4 Boxes / Case) |
| 17-201 | DRY FLAT SPONGE Bulk | 4" x 5" x 0.5" (10.2 cm x 12.7 cm x 1.3 cm) | 1,000 / Case |

DISTRIBUTED BY:
ADVANCED SURGI-PHARM INC.
www.advancedsurgipharma.com
Tel.1-800-661-5432



ENZYMATIC BEDSIDE SPONGE KIT

RightSwipe™

Introducing RightSwipe™, our newly launched Enzymatic Bedside Sponge Kits now feature sponges treated with an enhanced formula for a more effective pre-clean or surface clean solution. Our improved three-enzyme blend ensures thorough and safe cleaning and a simplified application process, making RightSwipe™ the perfect match for your endoscope care needs. The tubular sponge design is compatible with diverse endoscope models, providing versatility alongside the superior benefits of our enzymatic solution.



FLAT SPONGE



TUBULAR SPONGE



FREESTANDING POUCH

ENVIRO-SACHET™
NON-ENZYMATIC CLEANER

DISPENSER BOX

PRODUCT SELECTION

| ITEM # | DESCRIPTION | PACKAGING | CASE DIMENSIONS |
|--------|--|---------------------------|-------------------|
| 17-111 | PRE-TREATED TUBULAR SPONGE Add 100-500 ml of water to activate detergent solution. Comes in plastic bag ¹ | 25 / Box (4 Boxes / Case) | 21" x 14" x 13.5" |
| 17-211 | PRE-TREATED FLAT SPONGE Add 100-500 ml of water to activate detergent solution. Comes in plastic bag ¹ | 50 / Box (4 Boxes / Case) | 21" X 14" X 13.5" |
| 17-122 | PRE-TREATED TUBULAR SPONGE Comes in a freestanding pouch to use as a basin. Simply add 500 ml of water to the pouch. | 15 / Box (4 Boxes / Case) | 21" X 14" X 13.5" |
| 17-222 | PRE-TREATED FLAT SPONGE Comes in a freestanding pouch to use as a basin. Simply add 500 ml of water to the pouch. | 25 / Box (4 Boxes / Case) | 21" X 14" X 13.5" |
| 17-133 | READY TO USE TUBULAR SPONGE Comes pre-treated in 500 ml of an enzymatic solution in freestanding pouch. | 24 / Case | 14" X 12" X 9" |
| 17-233 | READY TO USE FLAT SPONGE Comes pre-treated in 500 ml of an enzymatic solution in freestanding pouch. | 24 / Case | 14" X 12" X 9" |
| 17-144 | READY TO USE DUAL PURPOSE TUBULAR SPONGE Comes pre-treated in 100 ml of an enzymatic solution in freestanding pouch. | 12 / Box (4 Boxes / Case) | 21" X 14" X 13.5" |
| 17-244 | READY TO USE DUAL PURPOSE FLAT SPONGE Comes pre-treated in 100 ml of an enzymatic solution in freestanding pouch. | 12 / Box (4 Boxes / Case) | 21" X 14" X 13.5" |
| 17-155 | READY TO USE TUBULAR SPONGE Comes pre-treated in 220 ml of an enzymatic solution in freestanding pouch. | 24 / Case | 14" x 12" x 9" |
| 17-255 | READY TO USE FLAT SPONGE Comes pre-treated in 220 ml of an enzymatic solution in freestanding pouch. | 24 / Case | 14" X 12" X 9" |
| 17-277 | ENZYMATIC GEL Comes with 1 Dry Flat Sponge, Tubular Sponge, Enviro-Sachet™ & Stainless Steel Bowl kit | 1 Kit | 1 Kit |

¹Our stainless steel bowl can be sold separately and used as a basin for easy and convenient sponge activation.

Your First ™

NON-ENZYMATIC BEDSIDE SPONGE KIT

RightSwipe™

Presenting our enhanced Non-Enzymatic Bedside Sponge Kits. Each kit includes sponges pre-treated with a free-rinsing, neutral pH detergent blend from RightSwipe™, ensuring an effective pre-cleaning and surface cleaning solution for medical instruments. This kit stands as Your First™ step in the medical instrument cleaning process, showcasing the advanced benefits for superior instrument care. Our enhanced bedside kits will provide you with a streamlined workflow.



FLAT SPONGE



TUBULAR SPONGE



FREESTANDING POUCH


 ENVIRO-SACHET™
 NON-ENZYMATIC CLEANER


DISPENSER BOX

PRODUCT SELECTION

| ITEM # | DESCRIPTION | PACKAGING | CASE DIMENSIONS |
|----------------|---|---------------------------|-------------------|
| 717-122 | PRE-TREATED TUBULAR SPONGE Comes in a freestanding pouch to use as a basin. Simply add 500 ml of water to the pouch. ¹ | 15 / Box (4 BOXES / CASE) | 21" X 14" X 13.5" |
| 717-222 | PRE-TREATED FLAT SPONGE Comes in a freestanding pouch to use as a basin. Simply add 500 ml of water to the pouch. ¹ | 25 / Box (4 Boxes / Case) | 21" X 14" X 13.5" |
| 717-133 | READY TO USE TUBULAR SPONGE Comes pre-treated in 500 ml of non-enzymatic solution in freestanding pouch. | 24 / Case | 14" X 12" X 9" |
| 717-233 | READY TO USE FLAT SPONGE Comes pre-treated in 500 ml of non-enzymatic solution in freestanding pouch. | 24 / Case | 14" X 12" X 9" |
| 717-144 | READY TO USE DUAL PURPOSE TUBULAR SPONGE Comes pre-treated in 100 ml of non-enzymatic solution in freestanding pouch. | 12 / Box (4 Boxes / Case) | 21" X 14" X 13.5" |
| 717-244 | READY TO USE DUAL PURPOSE FLAT SPONGE Comes pre-treated in 100 ml of non-enzymatic solution in freestanding pouch. | 12 / Box (4 Boxes / Case) | 21" X 14" X 13.5" |
| 717-155 | READY TO USE TUBULAR SPONGE Comes pre-treated in 220 ml of non-enzymatic solution in freestanding pouch. | 24 / Case | 14" x 12" x 9" |
| 717-255 | READY TO USE FLAT SPONGE Comes pre-treated in 220 ml of non-enzymatic solution in freestanding pouch. | 24 / Case | 14" X 12" X 9" |
| 717-277 | NON-ENZYMATIC GEL Comes with 1 Dry Flat Sponge, Tubular Sponge, Enviro-Sachet™ & Stainless Steel Bowl kit | 1 Kit | 1 Kit |

¹Our stainless steel bowl can be sold separately and used as a basin for easy and convenient sponge activation.



Your First ™ TRANSPORT BAGS FOR ENDOSCOPE AND MEDICAL INSTRUMENTS

Thoughtfully designed to eliminate touch-points during endoscope transport from point-of-care to the reprocessing area, and back, focused on safety of patients and medical staff.


- Made of recyclable LDPE material known to be quite flexible, highly resilient, and tough
- Latex and DEHP-free
- Excellent resistance to dilute and concentrated acids, alcohols, bases and esters*
- Good resistance to aldehydes, ketones, and vegetable-based oils*
- Can withstand high temperatures of up to 90 °C (194°F)
- Tamper-evident seal for added safety, also aids in keeping the instrument moist to avoid drying of bio-debris
- Clearly marked bags with color-coded visuals facilitates immediate identification of clean vs. used instruments
- Users can write important information directly on the bag in spaces marked with visual icons
- Available in two convenient sizes

GOOD TO KNOW: USING **Your First ™** SEALABLE TRANSPORT BAGS REDUCES USER ERROR, PREVENTS CROSS-CONTAMINATION, AND PROMOTES PATIENT AND MEDICAL STAFF SAFETY.

CLEAN INSTRUMENTS BAG



SOILED INSTRUMENTS BAG



MADE IN CANADA

PRODUCT SELECTION

| ITEM # | DESCRIPTION | DIMENSIONS | PACKAGING | CASE DIMENSIONS |
|-----------------|------------------------|---------------------------------|------------|-----------------|
| 17-HAP24 | Clean Instruments Bag | 24" x 22" / 61.0 cm x 56.0 cm | 200 / Case | 18" x 14" x 5" |
| 17-SAD24 | Soiled Instruments Bag | 24" x 22" / 61.0 cm x 56.0 cm | 200 / Case | 18" x 14" x 5" |
| 17-HAP30 | Clean Instruments Bag | 30.5" x 26" / 77.0 cm x 66.0 cm | 200 / Case | 21" x 19" x 7" |
| 17-SAD30 | Soiled Instruments Bag | 30.5" x 26" / 77.0 cm x 66.0 cm | 200 / Case | 21" x 19" x 7" |

*Reference on file



DISTRIBUTED BY:
ADVANCED SURGI-PHARM INC.
www.advancedsurgipharm.com
Tel.1-800-661-5432

SAVEFACE™ DRAPED SHIELDS

REDUCE PERSONAL INJURY
SUPERIOR INSTRUMENT DECONTAMINATION PROTECTION

SAVEFACE DRAPED-FACE SHIELD

In the Sterile Processing Department (SPD)/Central Sterile Supply Department (CSSD), personal safety is paramount.

Our Draped Face Shield provides essential protection during high-risk instrument decontamination.

EFFECTIVE PROTECTION IS NON-NEGOTIABLE

SAVEFACE Draped dynamically reduces harmful exposure to caustic chemicals as well as the risk of cross-contamination. Effective face protection is non-negotiable.

PPE CATEGORY 2 CERTIFIED

SAVEFACE DRAPED is tested to protect. Wear it with confidence.



PATENT PENDING



INCLUDES CLOSING SYSTEM TO REDUCE RISK OF FLUID CONTACT FROM BELOW THE SHIELD

SAVEFACE: KEY FEATURES

ANSI/AAMI PB70 LEVEL 3 COMPLIANT DRAPE:

- Provides a high level of fluid penetration protection.
- Tuck into clothing or pinch together.

EXTRA-LONG 9"/23CM WRAP-AROUND SHIELD:

- Ensures full coverage, reducing the risk of exposure from different angles.

COMFORTABLE:

- The soft, self-shaping foam headpiece molds comfortably to your forehead, allowing for extended wear without discomfort.

OPTICAL CLARITY WITH ANTI-FOG AND ANTI-GLARE

DON'T MISS THE OPPORTUNITY TO EXPLORE OUR BROADER RANGE OF PROTECTIVE SHIELDS.

Visit our website www.surgmed.com or consult our detailed literature to find the shield that best fits your needs.

VERY IMPORTANT

IS YOUR PRESENT FACE SHIELD CATEGORY 2 CERTIFIED?

Batrik face shields are independently certified to comply with the PPE Category 2 standard. Many facilities are now requesting compliance to this standard.

SHIELD LENGTH
9" / 23 CM

DRAPE LENGTH
6" / 15 CM



SF-240



SF-8240



PRODUCT SELECTION


| ITEM # | DESCRIPTION | DIMENSIONS | TOTAL LENGTH | PACKAGING | CASE DIMENSIONS |
|---------|---|---|--------------|-----------------------------------|-------------------|
| SF-240 | BATRIK SAVEFACE™ no drape | Shield Length: 9" / 23cm Shield width: 13' / 33cm | 9" (23cm) | 240 / Case (8 Boxes of 30 Pcs) | 26.5" x 21" x 15" |
| SF-8240 | BATRIK SAVEFACE™ with ANSI/AAMI PB70 Level 3 compliant protection drape | Shield Length: 9" / 23cm Shield width: 13' / 33cm Shield width: 6" / 15cm | 58 cm / 23" | 240 / Case (8 Boxes of 30 Pcs) | 26.5" x 21" x 15" |

SHOOZCOVERS™ | SOFT LINE™

Batrik ShoozCovers™ Soft Line™ is manufactured in Canada from high quality **35 g/m² polypropylene fabric and 100% latex free elastics** for performance and durability you can rely on. Our anti-skid tread reduces the risk of slipping and our smart design delivers easy on-and-off convenience and great fit.



MADE IN CANADA



SHOOZ-REPELLENT™ TECHNOLOGY:

The **Soft Line™** integrates Shooz-Repellent™ technology which creates a protective barrier that repels fluids, delivering basic incidental fluid contact repellent. Shoe covers are not waterproof or water-resistant.



ANTI-SKID THREAD:

Our anti-skid tread reduces the risk of slipping with **triple ultrasonic welds** that increases tear strength and delivers easy on-and-off convenience.

OUR SOFT LINE™ IS PERFECT FOR:



MEDICAL ESTABLISHMENTS



PHARMACEUTICAL LABS



GMP SITES



ELECTRONIC MANUFACTURING

DISTRIBUTED BY:
ADVANCED SURGI-PHARM INC.
www.advancedsurgipharm.com
Tel.1-800-661-5432




WHY YOU'LL LOVE THE SOFT LINE™:


- Both shoe covers are manufactured in Canada with top-quality 35 g/m² polypropylene fabric and 100% latex-free USA-made elastics!
- This line is particle count tested and complies with IEST-RP-CC003.4.
- Other colors available upon request.

SUPPORT A CAUSE, WEAR PINK:

- A portion of our Bubble Gum Pink shoe cover sales will be donated to charity to help further breast cancer research. Join us in making a meaningful contribution by choosing Pink ShoozCovers™!



AVAILABLE IN TWO VERSATILE SIZES



AVAILABLE IN UNIVERSAL

PRODUCT SELECTION:

| ITEM# | DESCRIPTION | COLOR | SIZE | PACKAGING | CASE DIMENSIONS |
|-------------|--|-----------------|---|--|--------------------------|
| SC-1210-PB | ShoozCovers™ - Soft Line™ <i>polypropylene non-woven fabric</i> | Lolly Blue | Universal Up to shoe size 12 ~ 42 cm x 16.5 cm / 16.5" x 6.5" | Case of 1,000 (10 bags of 100 units) | 20" x 16" x 16" |
| SC-1216-PDB | ShoozCovers™ - Soft Line™ <i>polypropylene non-woven fabric</i> <i>*Covers are available with or without skid-resistant tread.</i> | Lolly Blue | Universal Up to shoe size 12 ~ 42 cm x 16.5 cm / 16.5" x 6.5" | Case of 600 (6 Dispenser Boxes of 100 units) | 16.25" x 8.875" x 18.75" |
| SC-1216-PB | ShoozCovers™ - Soft Line™ <i>polypropylene non-woven fabric</i> <i>*Covers are available with or without skid-resistant tread.</i> | Lolly Blue | Universal Up to shoe size 12 ~ 42 cm x 16.5 cm / 16.5" x 6.5" | Case of 600 (6 bags of 100 units) | 19.5" x 19.5" x 14.25" |
| SC-1226-PDB | ShoozCovers™ - Soft Line™ <i>polypropylene non-woven fabric</i> <i>*Covers are available with or without skid-resistant tread.</i> | Lolly Blue | XL Shoe size 12-15 ~ 47 cm x 19 cm / 18.5" x 7.5" | Case of 600 (6 Dispenser Boxes of 100 units) | 17.25" x 9.375" x 21.75" |
| SC-1220-PB | ShoozCovers™ - Soft Line™ <i>polypropylene non-woven fabric</i> | Lolly Blue | XL Shoe size 12-15 ~ 47 cm x 19 cm / 18.5" x 7.5" | Case of 1,000 (10 bags of 100 units) | 22" x 16" x 16" |
| SC-1616-PDB | ShoozCovers™ - Soft Line™ <i>polypropylene non-woven fabric</i> | Bubble Gum Pink | Universal Up to shoe size 12 ~ 42 cm x 16.5 cm / 16.5" x 6.5" | Case of 600 (6 Dispenser Boxes of 100 units) | 16.25" x 8.875" x 18.75" |
| SC-1616-PB | ShoozCovers™ - Soft Line™ <i>polypropylene non-woven fabric</i> | Bubble Gum Pink | Universal Up to shoe size 12 ~ 42 cm x 16.5 cm / 16.5" x 6.5" | Case of 600 (6 bags of 100 units) | 19.5" x 19.5" x 14.25" |

*Many suppliers claim their products are latex free, however this often only relates to the fabric and not the elastics used. We source our ShoozCovers™ elastics from the USA to ensure that our components conform fully as a 100% latex free product.

PALLET INFORMATION:

| PRODUCT | CASES PER PALLET* | PALLET DIMENSIONS | PALLET WEIGHT | SIZE | PACKAGING UOM |
|----------|-------------------|-------------------|---------------|------|---------------------------------------|
| Softline | 40 | 40 x 48 x 48 in | 420 lbs | U | 6 Dispenser Box / Case (100 un. each) |

*Quantity of Cases per Pallet may vary

GrippezShooz™

EXTRA-LARGE FOR EXTRA-TOUGH WORK:

GrippezShooz™ provide exceptional protection and **enhanced traction** for professionals working in the most demanding, slippery environments. Featuring waterproof protection, our shoe covers ensure safety and hygiene without compromising on quality. They are the ideal solution for tough jobs.



MADE IN CANADA



SURGICAL PRODUCTS & OPERATING ROOM

PERFECT FOR:



AEROSPACE



FOOD PROCESSING
INDUSTRY



OIL, GAS AND
MINING SITES



CONSTRUCTION AND
RENOVATION SITES

DISTRIBUTED BY:
ADVANCED SURGI-PHARM INC.
www.advancedsurgipharm.com
Tel.1-800-661-5432



GrippezShooz™ ARE A PERFECT FIT FOR EVERYONE:

- Fast, Reliable Fit: Stretch material ensures a secure, easy-on fit. Coating provides protection against slippery surfaces.
- Comes in an extra-large size that fits over a men's size 15 steel-toed boot.
- Enhanced Safety: Anti-slip rubberized coating provides reliable grip on all surfaces, helping to prevent sliding in critical areas.
- Available in two colors: Gecko Navy and Rhino Grey.

RHINO GREY



GECKO NAVY



PRODUCT SELECTION:

| ITEM# | DESCRIPTION | COLOR | SIZE | PACKAGING | CASE DIMENSIONS |
|---------------|---|------------|---|---|--------------------------|
| SC-29312-NPDB | GrippezShooz™ <i>Polypropylene + polyethylene non-woven fabric, waterproof</i> | Rhino Grey | Universal Up to shoe size 12 ~ 42 cm x 17.5 cm / 16.5" x 6.9" | Case of 300 (2 dispenser boxes of 150 units) | 16.62" x 14.875" x 8.75" |
| SC-29312-NPB | GrippezShooz™ <i>Polypropylene + polyethylene non-woven fabric, waterproof</i> | Rhino Grey | Universal Up to shoe size 12 ~ 42 cm x 17.5 cm / 16.5" x 6.9" | Case of 300 (2 bags of 150 units) | 16" x 13.375" x 10" |
| SC-29322-NPDB | GrippezShooz™ <i>Polypropylene + polyethylene non-woven fabric, waterproof</i> | Rhino Grey | XL Shoe size 12-15 ~ 47 cm x 17.5 cm / 18.5" x 6.9" | Case of 240 (2 dispenser boxes of 120 units) | 16.62" x 14.875" x 8.75" |
| SC-29322-NPB | GrippezShooz™ <i>Polypropylene + polyethylene non-woven fabric, waterproof</i> | Rhino Grey | XL Shoe size 12-15 ~ 47 cm x 17.5 cm / 18.5" x 6.9" | Case of 240 (2 bags of 120 units) | 16" x 13.375" x 10" |
| SC-22212-NPDB | GrippezShooz™ <i>Polypropylene + polyethylene non-woven fabric, waterproof</i> | Gecko Navy | Universal Up to shoe size 12 ~ 42 cm x 17.5 cm / 16.5" x 6.9" | Case of 300 (2 dispenser boxes of 150 units) | 16.62" x 14.875" x 8.75" |
| SC-22212-NPB | GrippezShooz™ <i>Polypropylene + polyethylene non-woven fabric, waterproof</i> | Gecko Navy | Universal Up to shoe size 12 ~ 42 cm x 17.5 cm / 16.5" x 6.9" | Case of 300 (2 bags of 150 units) | 16" x 13.375" x 10" |
| SC-22222-NPDB | GrippezShooz™ <i>Polypropylene + polyethylene non-woven fabric, waterproof</i> | Gecko Navy | XL Shoe size 12-15 ~ 47 cm x 17.5 cm / 18.5" x 6.9" | Case of 240 (2 dispenser boxes of 120 units) | 16.62" x 14.875" x 8.75" |
| SC-22222-NPB | GrippezShooz™ <i>Polypropylene + polyethylene non-woven fabric, waterproof</i> | Gecko Navy | XL Shoe size 12-15 ~ 47 cm x 17.5 cm / 18.5" x 6.9" | Case of 240 (2 bags of 120 units) | 16" x 13.375" x 10" |

PALLET INFORMATION:

| PRODUCT | CASES PER PALLET* | PALLET DIMENSIONS | PALLET WEIGHT | SIZE | PACKAGING UOM |
|---------------|-------------------|-------------------|---------------|------|---------------------------------------|
| GRIPPEEZSHOOZ | 50 | 40 x 48 x 88 in | 446-479 lbs | U | 2 Dispenser Box / Case (150 un. each) |

*Quantity of Cases per Pallet may vary

SHOOZCOVERS™ | COATED LINE™


STAY DRY:

Our **Batrik ShoozCovers™ Coated Line™** is made from a **~55 g/m² blend of polypropylene and polyethylene**, providing high-quality water resistance.¹ Ideal for environments where maximum protection against fluids is required.



SURGICAL PRODUCTS & OPERATING ROOMS

MADE IN CANADA 



TRIPLE WELDED TECHNOLOGY:

These shoe covers are reinforced with triple ultrasonic welds, a technique that joins materials together without the need for adhesives or stitching. This technique increase tear strength and reduce linting, for long lasting, resilient covers.

OUR COATED LINE™ IS PERFECT FOR:



MEDICAL
ESTABLISHMENTS



PHARMACEUTICAL
LABS



GMP SITES



ELECTRONIC
MANUFACTURING

¹OUR SHOOZCOVERS™ COATED LINE™ IS RESISTANT AGAINST ALL NON PETROLEUM BASED LIQUIDS.

MARSHMALLOW WHITE | COATED LINE™



WHY YOU'LL LOVE THE COATED Line™:

- Superior protection using our unique blend of fabrics which protects against not only water but all non-petroleum-based liquids, making them the ideal choice for environments where spills and leaks are a concern
- Elastics around the ankle create a strong and secure fit which is especially important when working in hazardous environments where a misplaced shoe cover could lead to contamination or injury
- Smart design with a large opening, providing a better fit and easy on-and-off convenience
- Other colors available upon request



UN - XL

AVAILABLE IN TWO VERSATILE SIZES

PRODUCT SELECTION:

| ITEM# | DESCRIPTION | COLOR | SIZE | PACKAGING | CASE DIMENSIONS |
|---------------------|---|-------------------|---|---|--------------------------|
| SC-2112-NPB | ShoozCovers™ · Coated Line™ <i>Polypropylene + polyethylene non-woven fabric, waterproof</i> | Marshmallow White | Universal Up to shoe size 12 ~ 42 cm x 17.5 cm / 16.5" x 6.9" | Case of 300 (2 bags of 150 units) | 16" x 13.375" x 10" |
| SC-2112-NPDB | ShoozCovers™ · Coated Line™ <i>Polypropylene + polyethylene non-woven fabric, waterproof</i> | Marshmallow White | Universal Up to shoe size 12 ~ 42 cm x 17.5 cm / 16.5" x 6.9" | 150 / Dispenser Box 2 Dispenser Boxes / Case | 16.62" x 14.875" x 8.75" |
| SC-2122-NPB | ShoozCovers™ · Coated Line™ <i>Polypropylene + polyethylene non-woven fabric, waterproof</i> | Marshmallow White | XL Shoe size 12-15 ~ 47 cm x 17.5 cm / 18.5" x 6.9" | Case of 240 (2 bags of 120 units) | 16" x 13.375" x 10" |
| SC-2122-NPDB | ShoozCovers™ · Coated Line™ <i>Polypropylene + polyethylene non-woven fabric, waterproof</i> | Marshmallow White | XL Shoe size 12-15 ~ 47 cm x 17.5 cm / 18.5" x 6.9" | 120 / Dispenser Box 2 Dispenser Boxes / Case | 16.62" x 14.875" x 8.75" |

PALLET INFORMATION:

| PRODUCT | CASES PER PALLET* | PALLET DIMENSIONS | PALLET WEIGHT | SIZE | PACKAGING UOM |
|---------------|-------------------|-------------------|---------------|------|---------------------------------------|
| Coated | 72 | 40 x 48 x 90 in | 481 lbs | XL | 2 Dispenser Box / Case (120 un. each) |

*Quantity of Cases per Pallet may vary

SHOOZCOVERS™ | TOUGH LINE™



MADE IN CANADA 



EXTRA-LARGE FOR EXTRA-TOUGH WORK:

SHOE SIZES 12-15 (45-49.5 EU)
~ 47 CM X 19 CM / 18.5" X 7.5"

Our **Batrik ShoozCovers™ Tough Line™** comes in an extra-large size that fits over a **men's size 15 steel-toed boot**, delivering optimal full coverage and protection while you work.¹



ENHANCED DURABILITY:

These shoe covers are formulated with our **thickest fabric yet!** Made of **~55 g/m² polypropylene fabric**, these highly durable covers provide added protection, increase tear strength and reduce linting.

OUR TOUGH LINE™ IS PERFECT FOR:



HOMES OR
COMMERCIAL
BUILDINGS



SHIELDING AGAINST
MUD & GREASE



MEDICAL
ESTABLISHMENTS



PHARMACEUTICAL
LABS

¹ SHOOZCOVERSTM ARE MADE LONGER (42 & 47CM) FOR ACCURATE SHOE SIZE SELECTION.



WHY YOU'LL LOVE THE TOUGH LINE™:

- Triple ultrasonic welds that are long lasting to reduce waste
- High quality made in the USA latex-free elastics are triple reinforced to contract the shoe cover around your foot and deliver a strong and secure fit²
- Shooz-Repellent™ technology delivers basic incidental fluid contact repellency³
- Other colors available upon request



XL

AVAILABLE IN EXTRA LARGE

ANY TIME, ANYWHERE:



RESIDENTIAL



COMMERCIAL



INDUSTRIAL

PRODUCT SELECTION:

| ITEM# | DESCRIPTION | COLOR | SIZE | PACKAGING | CASE DIMENSIONS |
|-------------------|---|--------------------|---|--------------------------------------|---------------------|
| SC-1922-PB | ShoozCovers™ - Tough Line™ <i>Polypropylene non-woven fabric</i> | Black Liquorice | XL Shoe size 12-15 ~ 47 cm x 19 cm / 18.5" x 7.5" | Case of 200 (2 bags of 100 units) | 16" x 13.375" x 10" |

ALTERNATIVE PACKAGING: Add "-PDB" to end of code for a dispenser box instead of bags.

PALLET INFORMATION:

| | CASES PER PALLET* | PALLET DIMENSIONS | PALLET WEIGHT | SIZE | PACKAGING UOM |
|------------------|-------------------|-------------------|---------------|------|---------------------------------------|
| Toughline | 72 | 40 x 48 x 90 in | 481 lbs | XL | 2 Dispenser Box / Case (120 un. each) |

*Quantity of Cases per Pallet may vary

² Many suppliers claim their products are latex free, however this often only relates to the fabric and not the elastics used. We source our ShoozCovers™ elastics from the USA to ensure that our components conform fully as a 100% latex free product. ³ Shooz-Repellent™ technology provides basic incidental fluid contact repellent. It is not waterproof nor water resistant.

ANTIMICROBIAL IV POLES WITH CUVERRO SHIELD™

Frequently touched surfaces like IV poles have long been identified as reservoirs for the spread of pathogenic microbes.¹ These can easily contaminate the hands and equipment of healthcare professionals who in turn can transmit these pathogens to patients during routine care leading to Healthcare-associated infections.

Our Antimicrobial IV Poles protected with CuVerro Shield™ by Aereus Technologies* offers you a solution to this problem.


CuVerro Shield™ is a unique coating of naturally antimicrobial copper that is applied to our stainless steel poles, ram horns, and u-hooks to deliver outstanding protection to your highest touch zones.

The result? CuVerro® kills more than 99.9% of bacteria** within two hours and continues to kill 99% of bacteria** even after repeated contamination.

A US government funded study conducted by researchers at the US National Institutes of Health and the Centers for Disease Control and Protection found that the HCoV-19 virus remained viable for up to 2 to 3 days on plastic and stainless steel surfaces versus up to 4 hours on copper.²



OR-36770-CUV



Viability of SARS-CoV-1 and SARS-CoV-2 in Aerosols and on Various Surfaces.²

* CuVerro Shield™ by Aereus Technologies is a Health Canada Pest Management Regulatory Agency (PMRA) registered process.

**Laboratory testing shows that, when cleaned regularly, CuVerro surfaces kill greater than 99.9% of the following bacteria within 2 hours of exposure: Methicillin-Resistant Staphylococcus aureus, Staphylococcus aureus, Enterobacter aerogenes, Pseudomonas aeruginosa, E.coli 0157:H7, and Vancomycin-Resistant Enterococcus faecalis (VRE). The use of CuVerro® bactericidal copper products is a supplement to and not a substitute for standard infection control practices; users must continue to follow all current infection control practices, including those practices related to cleaning and disinfection of environmental surfaces. This surface has been shown to reduce microbial contamination, but it does not necessarily prevent cross contamination.

¹Oie S, Hosokawa I, Kamiya A. 2002. Contamination of room door handles by methicillin sensitive/methicillin-resistant Staphylococcus aureus. J. Hosp. Infect. 2002;51:140-143.

²Doremalen NV, Bushmaker T, Morris DH, et al. Aerosol and surface stability of HCoV-19 (SARS-CoV-2) compared to SARS-CoV-1. N Engl J Med 2020.



BASES & CASTERS

| ITEM # | DESCRIPTION |
|-------------------------|---|
| OR-36770-CUV | CuVERRO SHIELD™ IV POLE, CENTRE COLLAR, AND KNOB with 23" cream aluminium base (V-40016) and 2" casters |
| OR-36760-25-CUV | CuVERRO SHIELD™ IV POLE, CENTRE COLLAR, AND KNOB with 25" stainless steel base (V-40013-25) and 2" casters |
| OR-36760-253-CUV | with 25" stainless steel base (V-40013-25) and 3" casters |
| OR-36771-25-CUV | CuVERRO SHIELD™ IV POLE, CENTRE COLLAR, AND KNOB with 25" stainless steel base (V-40013-25), 2" casters, and weight kit of 12.5 lbs. (V-50004) |
| OR-36771-253-CUV | with 25" stainless steel base (V-40013-25), 3" casters, and weight kit of 12.5 lbs. (V-50004) |
| OR-36772-253-CUV | with 25" steel base (V-40014-25-ST) and 3" casters |

SMOOTH-ROLLING
NON-MARKING
DURABLE CASTERS




OPTIONAL WEIGHT KIT
AVAILABLE FOR INCREASED
STABILITY

HOOKS & RAM HORNS

| ITEM # | DESCRIPTION |
|---------------------|--|
| OR-36499-CUV | CuVERRO SHIELD™ STAINLESS STEEL RAM HORNS, AND ENDCAP Set of 2 |
| OR-36500-CUV | CuVERRO SHIELD™ STAINLESS STEEL RAM HORNS, AND ENDCAP Set of 4 |
| OR-36501-CUV | CuVERRO SHIELD™ STAINLESS STEEL U-HOOKS, AND ENDCAP Set of 2 |
| OR-36502-CUV | CuVERRO SHIELD™ STAINLESS STEEL U-HOOKS, AND ENDCAP Set of 4 |
| OR-36504-CUV | CuVERRO SHIELD™ STAINLESS STEEL RAM HORNS (2), STAINLESS STEEL U-HOOK (2), AND ENDCAP Set of 4 |
| V-50002-CUV | CuVERRO SHIELD™ STAINLESS STEEL URINARY DRAIN HOOK AND, COLLAR WITH KNOB Set of 1 |

CUSTOM HOOK DESIGNS
AVAILABLE TO SUIT YOUR
SPECIFIC NEEDS



HEAVY-DUTY RAM HORNS,
U-HOOKS OR A COMBINATION
OF BOTH AVAILABLE

POLES & HANDLES

| ITEM # | DESCRIPTION |
|--------------------|--|
| V-CUV-50003 | CuVERRO SHIELD™ ROUND STEEL HANDLE |
| V-CUV-50006 | CuVERRO SHIELD™ STAINLESS STEEL TRIANGULAR HANDLE |
| V-CUV-50007 | CuVERRO SHIELD™ SEMI-ROUND STEEL HANDLE |

THE LARGE PRONG KNOB FOR
HEIGHT ADJUSTMENT IS EASY TO
TURN EVEN WITH WET HANDS



GOOD TO KNOW: In the unlikely event your CuVerro coated IV Pole should become deeply scratched, a study conducted by the University of Waterloo demonstrated that antimicrobial activity of the Aereus coating is maintained even when completely removed up to a width of between 3 and 5 mm, under the test conditions used in this study.³

³ Tanvir S, Anderson WA. 2020. Antimicrobial Activity of "Deeply Scratched" Aereus-Coated Surfaces. Department of Medical Engineering. University of Waterloo. 2020.




PURCHASE DIRECTLY FROM OUR SISTER COMPANY
IMPERIAL SURGICAL LTD.
www.imperialsurgical.com
Tel:1-800-661-5432

PLEXIGLASS PROTECTION PANELS

Our Plexiglass Protection Panels are made of 100% acrylic clear plastic to safeguard face-to-face interaction. They are available in a variety of ready-made standard sizes plus we can create custom sizes/designs to accommodate your specific needs. Our simple and ergonomic design provides the necessary protection for your patients, customers, and community without an imposing presence.

FEATURES

- Multiple dimensions and opening sizes available
- 100% acrylic clear plastic
- Easy to install and move as needed
- Easy to disinfect



PPE - REUSABLE EQUIPMENT

PRODUCT SELECTION


| ITEM # | DIMENSIONS | ACRYLIC THICKNESS | OPENING SIZE |
|----------------|-------------------|-------------------|--------------|
| PG-2331 | 23.5" W x 31.5" H | 3/16" | 12" W x 4" H |
| PG-2336 | 23.5" W x 36" H | 1/4" | 8" W x 2" H |
| PG-2432 | 24" W x 32" H | 3/16" | 12" W x 4" H |
| PG-2436 | 24" W x 36" H | 3/16" | 12" W x 4" H |
| PG-3531 | 35.5" W x 31.5" H | 3/16" | 12" W x 4" H |
| PG-3632 | 36" W x 32" H | 3/16" | 5" W x 10" H |
| PG-3636 | 36" W x 36" H | 3/16" | 5" W x 10" H |
| PG-3735 | 37.5" W x 35.5" H | 1/4" | 8" W x 2" H |
| PG-4831 | 48" W x 31.5" H | 3/16" | 12" W x 4" H |

STAINLESS STEEL HANDS-FREE DOOR OPENER

The Imperial Surgical Ltd. Hands-Free Door Opener is a 100% stainless steel, foot operated door opener that works on any latch-less door.

Easy to install and easy to clean, the no-touch opener helps minimize the spread of germs by reducing the need to touch door handles where bacteria and viruses may live.

The hands-free door opener comes with an eye catching sticker to put on the door to encourage people to use it. We also provide the hardware so you have everything you need for installation.



MADE IN CANADA 

PRODUCT SELECTION

| ITEM # | DIMENSIONS | PACKAGING |
|----------|--|-----------|
| 11656-R0 | 5" Length x 1.5" Height x 0.14" Thickness (127 mm Length x 38 mm Height x 3.5 mm Thickness) | 1 Unit |



Ver.2021-04-22

naso-gastric

Feeding Tubes

Feeding tubes in DEHP-free PVC (for short-term use) and polyurethane (for medium-term use) with flexible female mount and cap.

- Smoothly rounded, closed distal tip and opposed lateral eyes
- Radiopaque
- Available in various lengths: Centimeter graduations (5cm to 25cm)



Standard tube

PVC tubes

| Size | 04 Fr | 05 Fr | 06 Fr | 07 Fr | 08 Fr | 09 Fr | 10 Fr | 12 Fr | Units/Case |
|--------|--------|--------|--------|--------|--------|--------|--------|--------|------------|
| 40 cm | 310.04 | 310.05 | 310.06 | 310.07 | - | - | - | - | 50 |
| 50 cm | - | - | - | - | 310.08 | 310.09 | 310.10 | 310.12 | 50 |
| 75 cm | 311.04 | 311.05 | 311.06 | - | 311.08 | - | 311.10 | - | 50 |
| 125 cm | 312.04 | - | 312.06 | 312.07 | 312.08 | 312.09 | 312.10 | 312.12 | 25 |


PUR tubes

| Size | 04 Fr | 05 Fr | 06 Fr | 08 Fr | Units/Case |
|-------|---------|---------|---------|---------|------------|
| 40 cm | 1310.04 | 1310.05 | 1310.06 | - | 25 |
| 50 cm | - | - | - | 1310.08 | 25 |
| 75 cm | - | 1311.05 | 1311.06 | 1311.08 | 25 |

dual-flowgastric

Feeding Tubes

The dual-flow gastric tubes are used to prevent the accumulation of gastric secretions and air, thus allowing decompression of the stomach by continuous aspiration without trauma to the gastric mucosa.



Dual-flow tube

Dual-flow gastric tubes

| Product Code | | Size | Length | Number of eyes | Units/Case |
|--------------|---------|-------|--------|----------------|------------|
| PVC | PUR | | | | |
| 340.06 | 1340.06 | 6 Fr | 60 cm | 6 | 10 |
| 340.08 | 1340.08 | 8 Fr | 60 cm | 6 | 10 |
| 340.10 | 1340.10 | 10 Fr | 90 cm | 10 | 10 |
| 340.12 | 1340.12 | 12 Fr | 90 cm | 10 | 10 |
| 340.14 | - | 14 Fr | 120 cm | 10 | 10 |
| 340.16 | 1340.16 | 16 Fr | 120 cm | 10 | 10 |
| 340.18 | 1340.18 | 18 Fr | 120 cm | 10 | 10 |

gastro-duodenal

Feeding Tubes

125 cm tubes for placement in the duodenum.

Weighted tubes: Silicone (for both medium and long-term use) tube with weighted metal tip to facilitate placement

Levin type tubes: DEHP-free PVC (for short-term use) tubes for gastric decompression or gavage feeding.



Silicone tube



Levin type tubes with standard connector

Weighted tubes (silicone)

| Product Code | Size | Units/Case |
|--------------|------|------------|
| 2395.05 | 5 Fr | 25 |
| 2395.06 | 6 Fr | 25 |
| 2395.09 | 9 Fr | 25 |

Levin type tubes with standard connector (PVC)

| Product Code | Size | Units/Case |
|--------------|-------|------------|
| 391.10 | 10 Fr | 50 |
| 391.12 | 12 Fr | 50 |
| 391.14 | 14 Fr | 50 |
| 391.16 | 16 Fr | 50 |

nutriSafe2®

A complete range dedicated to safe enteral nutrition. Designed to be fail-safe, the nutriSafe2® system incorporates a unique connection that locks securely so it won't slip apart and yet cannot be connected to a standard 6% tapered luer fitting.

Safe - The risk of misconnection errors is significantly reduced for patient safety.

Secure - Accidental disconnection is nearly impossible, ensuring precious nutrients are received.

Simple - It does not change the technique of the users and therefore does not require any special training.

nutriSafe2® Feeding tubes

nutriSafe2® feeding tubes are available in a range of sizes, lengths and materials to suit patient needs. Numbered graduations every centimeter aid placement and monitoring, and the rounded tip provides atraumatic insertion.

PVC (DEHP-free): for use up to 5 days.

Polyurethane (PUR): for use up to 30 days.

Silicone: for use up to 30 days.

The nutriSafe2® connection allows the use of these tubes only with other nutriSafe2® devices..

PVC tubes



| Size | 04 Fr ● | 05 Fr ● | 06 Fr ● | 08 Fr ● | 10 Fr ● | Units/Case |
|--------|---------|---------|---------|---------|---------|------------|
| 40 cm | 364.042 | 364.052 | 364.062 | 364.082 | - | 25 |
| 50 cm | 361.042 | 361.052 | 361.062 | 361.082 | 361.102 | 50 |
| 75 cm | 363.042 | 363.052 | 363.062 | 363.082 | 363.102 | 25 |
| 125 cm | 362.042 | 362.052 | 362.062 | 362.082 | 362.102 | 25 |

PUR tubes



| Size | 04 Fr ● | 05 Fr ● | 06 Fr ● | 08 Fr ● | 10 Fr ● | Units/Case |
|--------|----------|----------|----------|----------|----------|------------|
| 40 cm | 1364.042 | 1364.052 | 1364.062 | 1364.082 | - | 25 |
| 50 cm | 1361.042 | 1361.052 | 1361.062 | 1361.082 | 1361.102 | 50 |
| 75 cm | 1363.042 | 1363.052 | 1363.062 | 1363.082 | 1363.102 | 25 |
| 125 cm | 1362.042 | 1362.052 | 1362.062 | 1362.082 | 1362.102 | 25 |

Silicone tubes



| Size | 04 Fr ● | 05 Fr ● | 06 Fr ● | 08 Fr ● | 10 Fr ● | Units/Case |
|-------------------|-----------|-----------|-----------|-----------|----------|------------|
| 50 cm | 2331.042G | 2331.052G | 2331.062G | 2331.082G | 2331.102 | 10 |
| 125 cm | 2332.042G | - | 2332.062G | 2332.082G | 2332.102 | 10 |
| 125 cm (weighted) | - | - | 2395.062G | 2395.082G | - | 10 |

nutrisafe² Syringes

A complete range of syringes, from 1 ml to 60ml for administration of enteral nutrition and oral medication.

Amber colored syringes allow preservation of light sensitive drugs. They are available with caps.

These syringes only connect with other items in the nutrisafe²® range.



Transparent syringes



Amber syringes

C-GON - Transparent syringes

| Product Code | Volume | Graduations | Connection | Units/Case |
|--------------|--------|-------------|-----------------|------------|
| 1015.012 | 1 ml | 0.01 ml | Centered / slip | 100 |
| 1015.022 | 2.5 ml | 0.1 ml | Centered / slip | 100 |
| 1015.052 | 5 ml | 0.2 ml | Centered / slip | 100 |
| 1015.102 | 10 ml | 0.2 ml | Centered / lock | 100 |
| 1015.202 | 20 ml | 1 ml | Centered / lock | 50 |
| 1015.352 | 35 ml | 1 ml | Centered / lock | 50 |
| 1015.602 | 60 ml | 1 ml | Offset / lock | 25 |

Monoject - Transparent syringes

| Product Code | Volume | Graduations | Connection | Units/Case |
|--------------|--------|-------------|-----------------|------------|
| 1015.062M | 6 ml | 0.1 ml | Centered / slip | 100 |
| 1015.202M | 20 ml | 1 ml | Centered / lock | 50 |
| 1015.352M | 35 ml | 1 ml | Offset / lock | 25 |
| 1015.602M | 60 ml | 1 ml | Offset / lock | 25 |

Amber syringes (bulk non-sterile)

| Product Code | Volume | Graduations | Connection | Units/Case |
|--------------|--------|-------------|-----------------|------------|
| 1022.012 | 1 ml | 0.01 ml | Centered / slip | 50 |
| 1022.022 | 2.5 ml | 0.1 ml | Centered / slip | 50 |
| 1022.052 | 5 ml | 0.2 ml | Centered / lock | 50 |
| 1022.102 | 10 ml | 1 ml | Centered / lock | 50 |

nutrisafe² Extension lines

Extension lines to allow for positioning of the syringe driver away from the incubator and to move handling away from the newborn.

Available in standard bore for breast milk or in large bore for fortified milk.



| Product Code | Int. ø | Length | Clamp | Units/Case |
|--------------|--------|--------|-------|------------|
| 368.032-50 | 1.5 mm | 30 cm | No | 50 |
| 368.152-50 | 1.5 mm | 150 cm | No | 50 |
| 369.152 | 2.8 mm | 150 cm | No | 25 |

nutrisafe² Caps

A range of caps to secure the contents of vials containing oral / enteral medications. These caps allow access to the medication only from a nutrisafe² syringe, to avoid risk of accidental injection through IV systems.

Caps for syringes and tubes prevent contamination of the feed or drugs between feeding sessions.

Needleless vial access caps

| Product Code | Connector | Int. ø | Units/Case |
|--------------|-----------|--------|------------|
| 821.022 | male | 20 mm | 25 |



Needleless vial access caps

Caps for syringes and tubes

| Product Code | Description | Units/Case |
|--------------|----------------------|------------|
| 828.01 | Cap for syringe (x1) | 100 |
| 828.08 | Cap for syringe (x8) | 100 |



Syringe cap

Vented caps for syringes - for gravity feeding

| Product Code | Description | Units/Case |
|--------------|-------------|------------|
| 818.22 | 20 ml | 50 |
| 818.62 | 60 ml | 50 |



Vented cap

nutrisafe2® sampling straws



nutrisafe2® drawing-up straws allow the filling of nutrisafe2® syringes with feed or medication from all sizes of containers and bottles.

| Product Code | Size | Bore | Units/Case |
|--------------|-------|---------------|------------|
| 817.052 | 5 cm | Standard bore | 50 |
| 817.152 | 15 cm | Standard bore | 50 |
| 817.153 | 15 cm | Large bore | 50 |
| 817.202 | 20 cm | Standard bore | 50 |
| 817.203 | 20 cm | Large bore | 50 |
| 817.302 | 30 cm | Large bore | 50 |
| 817.502 | 50 cm | Standard bore | 20 |

nutrisafe2® accessories

Connectors, stopcocks and accessories to facilitate daily use of the nutrisafe2® range.

Connectors

| Product Code | Description | Units/Case |
|--------------|---|------------|
| 368.02 | Adaptor for feeding sets | 100 |
| 368.12 | Taper connector | 30 |
| 368.32 | Extension tube with connector for suction | 25 |
| 368.42R | Y connector with oral port | 100 |



Adaptor for feeding sets



Y connector



Taper Connector



Extension Tube with connector

Stopcock

| Product Code | Description | Units/Case |
|--------------|----------------|------------|
| 5802.01 | 3-way stopcock | 50 |



Stopcock

Adapters

| Product Code | Description | Units/Case |
|--------------|----------------------|------------|
| 828.88 | Oral tip adapter | 1000 |
| 828.99 | Medicine cap adapter | 1000 |



Oral tip adapter





Medicine cap adapter

Premicath®

A 1Fr PICC line designed specially for newborns under 1kg

Infants born prematurely weighing under 1kg often need the nutritional and therapeutic support that can only be delivered into their tiny vessels via venous access. At just 1Fr, Premicath® has been specially designed for the smallest and most fragile babies to deliver optimum accuracy during placement as well as minimizing trauma for patients.

Available with choice of introducers including splitting needle and microflash®



Premicath®

| Product Code | Catheter | | | Introducer | | Stylet | Units/Box |
|---------------|----------|----------------|----------------------------|-----------------|------|--------|-----------|
| | Length | Priming Volume | Flow rate (Pressure 1 bar) | Type | Size | | |
| 1261.20G | 20 cm | 0.09 ml | 0.7 ml / min | SplittingNeedle | 24 G | No | 10 |
| 1261.203G | 20 cm | 0.09 ml | 0.7 ml / min | SplittingNeedle | 24 G | Yes | 10 |
| 1261.205G | 20 cm | 0.09 ml | 0.7 ml / min | Neocath Split | 22 G | Yes | 10 |
| 1261.208G | 20 cm | 0.09 ml | 0.7 ml / min | Microflash | 20 G | Yes | 10 |
| 1261.306G | 30 cm | 0.11 ml | 0.6 ml / min | No Introducer | - | Yes | 10 |
| Higher Barium | | | | | | | |
| 1261.203A | 20 cm | 0.09 ml | 0.7 ml / min | SplittingNeedle | 24 G | Yes | 10 |
| 1261.306A | 30 cm | 0.11 ml | 0.6 ml / min | No Introducer | - | Yes | 10 |

Kits & Accessories

| Product Code | Description | Units/Box |
|--------------|---------------------------------|-----------|
| 7371.19 | 20 G Microflash Introducer | 10 |
| 1261.20N | 24 G Splitting Needle | 20 |
| AMS-850TC | V-Cutter Catheter Trimming Tool | 20 |



V-Cutter



Splitting needle




Microflash



nutri^{line}[®]

A PICC line with a unique & easy-to-use peelable needle available in 2 Fr

nutri^{line}'s one-piece catheter provides clinicians with a high degree of safety. Available in a choice of lengths to ensure accurate tip placement for most IV access sites.



| Product code | Length | Size | Priming volume | Flow rate | Introducer | Stylet | Units/Case |
|--------------|--------|------|----------------|--------------|-----------------------|--------|------------|
| 1252.30G | 30 cm | 2 Fr | 0.12 ml | 5.0 ml / min | Splitting needle 20 G | No | 10 |
| 1252.31G | 30 cm | 2 Fr | 0.12 ml | 5.0 ml / min | No Introducer | No | 10 |
| 1252.35G | 30 cm | 2 Fr | 0.12 ml | 5.0 ml / min | Microflash 20 G | No | 10 |
| 1252.030G | 30 cm | 2 Fr | 0.12 ml | 5.0 ml / min | Microflash 20 G | Yes | 10 |
| 1252.031 | 30 cm | 2 Fr | 0.12 ml | 5.0 ml / min | No Introducer | Yes | 10 |




DISTRIBUTED BY:
ADVANCED SURGI-PHARM INC.
www.advancedsurgipharm.com
Tel:1-800-661-5432

nutri^{line}[®]

TWINFLO

A 2Fr dual lumen catheter with peelable needle

nutri^{line}[®] twinflo has been designed specifically for babies who require multiple infusions. Its dual lumen tip design offers enhanced fluid management and avoids incompatible drugs mixing.



| Product Code | Catheter | | | Introducer | | Stylet | Units / Case |
|--------------|----------|---------------------|---------------------------|---------------|------|--------|--------------|
| | Length | Priming Case Volume | Flow rate (pressure 1bar) | Type | Size | | |
| I252.232G | 30 cm | 2 x 0.2 ml | 2 x 1.45 ml / min | No introducer | - | No | 10 |
| I252.232M | 30 cm | 2 x 0.2 ml | 2 x 1.45 ml / min | Microflash | 20 G | No | 10 |
| I252.230G | 30 cm | 2 x 0.2 ml | 2 x 1.45 ml / min | Microflash | 20 G | Yes | 10 |

Kits & Accessories

| Product Code | Description | Units/Box |
|--------------|---------------------------------|-----------|
| 7371.19 | 20 G Microflash Introducer | 10 |
| I262.30N | 20 G Splitting Needle | 20 |
| AMS-850TC | V-Cutter Catheter Trimming Tool | 20 |





introducers

For PICC catheters

Stand alone introducers for neonatal PICC catheters: splitting needle, microflash® peelable cannula and IV cannula.



IV cannula

IV cannula

| Product Code | Size | For catheter | Units/Case |
|--------------|------|--------------|------------|
| CMS-207 | 24 G | 1 Fr | 10 |



Splitting needle

Splitting needle

| Product Code | Size | For catheter | Units/Case |
|--------------|------|--------------|------------|
| 1261.20N | 24 G | 1 Fr | 20 |
| 1262.30N | 20G | 2 Fr | 20 |




microflash™

Peelable Introducer/Peelable Safety Introducer*

Peelable Introducer for the insertion of PICC & Midline catheters, with notched needle and clear cannula sheath for early detection of blood flashback

| Product Code | Size | For catheter | Units/Case |
|--------------|------|--------------|------------|
| 7371.19 | 20 G | 1 Fr / 2 Fr | 10 |

- A** Top hole needle shaft permits blood flashback between needle & cannula for immediate indication of vessel access
- B** Short, light weight safety housing for balance and ease of handling
- C** Curved hub for secure ergonomic fingertip placement
- D** Tri-bevel needle tip for ease of skin penetration
- E** Smooth needle to cannula transition reduces drag through the tissue
- F** Transparent cannula permits immediate blood visualization upon vessel penetration
- G** Funnel cannula opening to facilitate easy catheter entry
- H** Radiopaque stripe on cannula for patient safety



*The microflash® Introducer is intended for use as a split cannula introducer used to puncture peripheral veins.




DISTRIBUTED BY:
ADVANCED SURGI-PHARM INC.
www.advancedsurgipharm.com
Tel:1-800-661-5432

Griplok

Catheter Securement Device

Our Vygon Grip-Lok® is especially efficient with Vygon Neonatal Catheters. It is hypoallergenic and latex free, which means low risk of irritation and allergic reaction. Its zinc oxide adhesive allows for gentle, yet efficient skin fixation. The catheter fixation is secure to prevent migration and its suture-less fixation increases ease of use.

- Hypoallergenic
- Zinc Oxide Adhesive
- Specially designed hub insert
- Suture-less fixation
- Sterile
- Latex free



| Product Code | Description | Units/Case |
|--------------|-------------------------------------|------------|
| 5804.08 | Grip-Lok Catheter Securement Device | 20 |

umbilical catheters

Single-lumen umbilical catheters

A single lumen umbilical catheter used for venous or arterial access.

Umbilical venous catheters are used for:

- Parenteral nutrition & fluid administration
- Drugs administration
- Venous blood sampling
- Transfusion of blood or blood products
- Exchange transfusion

Umbilical arterial catheters are used for:

- Arterial blood sampling
- Arterial pressure measurement
- Blood pH and blood gas analysis
- Fluids and drugs administration



Three-way tap

- Colour-coded caps for venous / arterial identification


Luer-lock connector

- For efficient connection



Numerical markings every cm

- Make it easy to check the length of catheter inserted



| Product Code | Size | Length | Priming vol | Flow rate | Material | Units/Box |
|--------------|--------|--------|-------------|---------------|----------|-----------|
| 1270.02 | 2.5 Fr | 30 cm | 0.12 ml | >3 ml / min | PUR | 8 |
| 1270.03 | 3.5 Fr | 40 cm | 0.27 ml | >10 ml / min | PUR | 8 |
| 1270.04 | 4 Fr | 40 cm | 0.29 ml | >10 ml / min | PUR | 20 |
| 1270.05 | 5 Fr | 40 cm | 0.40 ml | >23 ml / min | PUR | 20 |
| 1270.08 | 8 Fr | 40 cm | 0.81 ml | >100 ml / min | PUR | 20 |
| 270.04 | 4 Fr | 37 cm | 0.27 ml | >16 ml / min | PVC | 15 |
| 270.05 | 5 Fr | 37 cm | 0.38 ml | >31 ml / min | PVC | 15 |
| 270.06 | 6 Fr | 37 cm | 0.51 ml | >52 ml / min | PVC | 15 |
| 270.08 | 8 Fr | 37 cm | 0.75 ml | >106 ml / min | PVC | 15 |

Double lumen umbilical catheter

A double lumen umbilical catheter for short term catheterisation through the umbilical vein when multiple access is needed.

- Infusion of incompatible drugs
- Reduces the need for additional access
- Infusion of medication and nutrition at the same time



Colour-coded hubs

- For easy identification and safety



Double lumen

- Clearly marked lines with separate exit points prevent incompatible drugs mixing

Neonatal clamps

- Small line clamp for line management

Numerical markings every cm

- Make it easy to check the length of catheter inserted



X-ray opaque PUR catheter

- X-ray opaque for accurate tip location without additional contrast medium

| | Product Code | Size | Exit Points | Length | Priming volume | Flow rate | Units/Box |
|-----|--------------|------|---------------------------|--------|----------------|----------------|-----------|
| PUR | 1272.14 | 4Fr | One lateral One distal | 20 cm | 2 x 0.19 ml | 11 & 12 ml/min | 10 |
| | 1274.14 | 4Fr | | 40 cm | 2 x 0.27 ml | 6 & 7 ml/min | 10 |
| | 1274.17 | 5Fr | | 40 cm | 2 x 0.26 ml | 7 & 9 ml/min | 10 |

Accessories for umbilical catheters

Lockable three-way stopcock with identification caps (red and blue) to distinguish between venous and arterial access.

- User feels digitally the position chosen
- Safety Luer-lock prevents accidental disconnections
- Lipid-resistant



| Product Code | Description | Units/Case |
|--------------|---|------------|
| 876.20 | Lockable three-way stopcock with identification caps (red and blue) | 50 |

exchange transfusion tray

Set containing the accessories necessary to perform removal and replacement of a patient's blood.



| Product Code | Contents | Units/Case |
|------------------|---|------------|
| 275.001 | 1x four-way stopcock equipped with 1 injection site (for injections of additional medication) | 1 |
| | 1x extension tube for evacuation of the discarded blood | |
| | 1x 5Fr exchange transfusion catheter (PVC - XRO) 37 cm long, with 2 sides holes and centimetric markings | |
| | 1x 7Fr exchange transfusion catheter (PVC - XRO) 37 cm long, with 2 sides holes and centimetric markings | |
| | 2x 20 ml Luer-lock syringes | |
| | 1x 10 ml Luer-lock syringe | |
| | 1x 16x0.5 mm (25G) Luer hypodermic needle | |
| | 1x graduated plastic container (2000 ml) | |
| | 1x transfusion set, 163 cm long, 18.8 ml dead volume | |
| | 1x 15 cm ruler for venous pressure measurement | |
| | 3x 50 x 50 mm gauze dressings | |
| | 1x 50 x 60 cm fenestrated drape | |
| | 1x air-venting needle | |
| 1x control sheet | | |

goldtrace

The Gold-plated Fetal Scalp Electrode

A message from the heart

Goldtrace® Fetal Spiral Electrode (FSE) is used when fetal heart rate monitoring during labor requires continuous and reliable internal registration. Goldtrace FSE is made of materials whose properties ensure the best-possible-quality signal for CTG.

Goldtrace FSE comes in a flexible protective guide tube that keeps the spiral tip inside the tube during application, protecting the user, mother and baby. The flexibility of the guide tube makes it possible to form according to desired anatomical curvature.

Specifications

- Gold-plated spiral electrode
- Single use only
- Includes disposable applicator
- Not made with natural rubber latex
- Sterile product

Article number

N-CNS000004 Goldtrace Fetal Spiral Electrode in boxes of 75 single packed pcs

Materials

| | |
|-----------------|---------------------------------------|
| Spiral needle | 316 Stainless Steel with 10-25µm Gold |
| Reference plate | 316 Stainless Steel with 10-25µm Gold |
| Hub | PE |
| Connector | PE and brass, gold-plated |
| Wires | Copper with PVC jacket |
| Drive tube | PE |
| Guide tube | PE |
| Drive tube grip | PE |



FECG Legplate Adapter for Philips

The Connection Between Mother and Child

The FECG Legplate Adapter for Philips provides a signal path between the Fetal Spiral Electrode (FSE), Skin Electrode, and Fetal Monitor when continuously monitoring fetal heart rate during labor and delivery procedures.



Use a FECG Legplate that is compatible with your Fetal Monitor and connect only compatible FSE to the FECG Legplate.

The FECG Legplate Adapter for Philips is compatible with:

- Philips Avalon Monitor

Article Number

N-ACC300120-R0

FECG Legplate Adapter for Philips

neohelp™

A sterile heat loss prevention suit to deliver life-saving thermal care

Immediately after a premature baby is delivered by c-section, the priority is to place them in a warm and sterile environment. neohelp™ provides this essential protection so other essential Golden Hour care and treatment, including delayed cord clamping, can take place.^(1,2)

Integrated adjustable hood

Decreases heat loss through radiation:

- More efficient than a stockinette cap allowing air to pass through the material
- The toggle allows you to adjust to the baby's head

Double layer of soft clear polyethylene

Decreases heat loss through convection and evaporation:

- Creates a warm and humid environment, mimicking the incubator effect
- The thin inner allows for excellent skin contact
- Creates a barrier against drafts
- Allows passage of radiant heat from an additional warming device (if used)⁽³⁾
- Transparent material allows vital observations⁽⁴⁾

Pre-shaped foam support


Decreases heat loss through conduction:

- Helps to maintain an open airway by raising the shoulders
- Stabilises the baby's position
- Provides thermal care during transportation
- Provides comfort


Central VELCRO® opening

Provides optimum seal to ensure heat conservation:


- Quick and easy to place around the baby
- Allows full access to the baby's body
- Designed for easy placement of monitoring equipment, IV, umbilical catheters and carrying out Golden Hour care




How to use neohelp™




1. Unfold neohelp™ completely.
2. Place it on the resuscitation table, under the heat lamp (if applicable).
3. Fully open the Velcro and the hood to facilitate the baby's reception.



4. Do not dry the baby at birth.
5. Place the baby immediately after birth in neohelp™
6. Align the shoulders at the upper edges of the foam.
7. Close the Velcro tightly along its entire length.



8. Adjust the hood to the head of the baby.
Do not cover the airways due to the risk of suffocation.
9. Other heating appliances (e.g. heating) will need to be adjusted accordingly due to risk of hyperthermia.



- 10 The baby should be kept wrapped until temperature stabilization.
- 11 Temperature of the baby should be monitored continuously or consecutively (ideally every 5 minutes).
- 12 Do not use more than 24 hours.

Why is thermal care important?

Thermal care is vital for a preterm infant because they may have unbalanced skin-surface to weight ratio, very little or no capacity to generate heat (brown adipose tissue), inadequate stores of subcutaneous (insulating) fat and immature epidermal barrier.^(4,5)

At this vital time, neohelp™ prevents heat loss through its double layer of soft, clear polyethylene, integrated adjustable hood and VELCRO® seal.

For every 1°C decrease:

Risk of sepsis increases by 11%
Risk of death increases by 28%⁽³⁾

In the first 10-20 minutes,

without any protection,
temperature can fall by 2-4°C⁽⁶⁾

Why is delayed cord clamping important?

Provided the baby can be kept warm and does not need immediate resuscitation, the Resuscitation Council UK (RCUK) recommends delayed cord clamping (DCC) for at least 60 seconds whilst breathing is established.⁽⁷⁾

Delayed cord clamping has been shown to reduce the relative risk of:

- Intraventricular haemorrhage by 41% (RR 0.59, 95% CI 0.41 to 0.85)⁽⁸⁾
- Necrotising enterocolitis by 38% (RR 0.62, 95% CI 0.43 to 0.90)⁽⁸⁾


In addition, this procedure increases circulation of blood volume after birth and an improvement in cardiovascular stability, reducing the need for a blood transfusion.

Neonatal hypothermia: a worldwide issue

Hypothermia is an **important factor** in **morbidity** and **mortality** of all birthweights and gestational ages, and **particularly for vulnerable preterm infants**⁽⁹⁾.


Incidence of hypothermia on admission in the NICU from the delivery room is ⁽¹⁰⁾:

- ≥ 56% for infants < 750g
- ≥ 25% for infants ≤ 2500g



Consequences of neonatal hypothermia

Consequences of neonatal hypothermia ^(9,4)



Neonatal heat loss prevention suit

neohelp™ is a sterile suit to swaddle the baby immediately after birth (before resuscitation).

- Polyethylene
- Occlusive
- Transparent

neohelp™ prevents heat loss due to convection, conduction, radiation and evaporation.

| Product Code | Baby weight | Description | Dimensions | Units/Box |
|--------------|-------------|-----------------|---------------------|-----------|
| 37.09.14 | < 1 kg | neohelp™ Small | L. 38 cm x W. 30 cm | 10 |
| 37.09.15 | 1 – 2.5 kg | neohelp™ Medium | L. 44 cm x W. 38 cm | 10 |
| 37.09.16 | > 2.5 kg | neohelp™ Large | L. 50 cm x W. 38 cm | 10 |

International recommendations

“The transparency of bags makes it easier for caregivers to observe and manage the infant with minimal disruption of the wrap.”⁽⁴⁾

“The transport incubator used to limit heat loss can be cumbersome and difficult to obtain. It may be «replaced» by a stockinette cap and a transparent polyethylene bag wrapping whilst the baby is still wet. This greatly reduces the risk of hypothermia.”⁽¹¹⁾

“Meta-analysis of [...] studies found that plastic wraps (polyurethane or polyethylene bag) were statistically significantly more effective than routine care in reducing heat losses in infants aged < 28 weeks' of gestation. Stockinette caps were not effective in reducing heat loss in infants”.⁽¹²⁾

References

1. Rabe H, Diaz-Rossello JL, Duley L, Dowswell T. Effect of timing of umbilical cord clamping and other strategies to influence placental transfusion at preterm birth on maternal and infant outcomes. *Cochrane Database Syst Rev* 2012;8 CD003248.
2. Backes CH, Rivera BK, Haque U, et al. Placental transfusion strategies in very preterm neonates: a systematic review and metaanalysis. *ObstetGynecol*2014; 124:47-56
3. W.W.M.V. Mori R. Interventions to prevent hypothermia at birth in preterm and/or low birth weight infant. RHL.
4. T. Cordaro and al. Hypothermia and occlusive skin wrap in the low birth weight premature infant. *NAINR*. 2012;12(2):78-85.
5. B. Mathew and al. Vinyl Bags prevent hypothermia At Birth in Preterm Infants. 2006.
6. World Health Organisation (WHO). Thermal protection of the newborn: practical guide. 1997.
7. Resuscitation-Council-UK-(RCUK). Newborn Life Support – National Resuscitation Guidelines. 2015.
8. Rabe H, Diaz-Rossello JL, Duley L, Dowswell T. Early cord clamping versus delayed cord clamping or cord milking for preterm babies. *Cochrane Database Syst Rev* 2012;15 CD003248.
9. The Cochrane Collaboration. Interventions to prevent hypothermia at birth in preterm and/or low birthweight infants (review), 2010.
10. DR Bhatt, R. White and al. Transitional hypothermia in preterm newborns. *Journal of Perinatology*, 2007.
11. J.F. Diependaele and A. Fily. Management of a newborn baby, 51st congress of French Society of Anesthesia and Intensive care, 2009.
12. Wariki WMV and Mori R. Interventions to prevent hypothermia at birth in preterm and/or lowbirth-weight infants: RHL commentary (last revised: 1 June 2010). The WHO Reproductive Health Library; Geneva: World Health Organization.

Warm Gel® Infant Heel Warmer

Uniform heat to reduce heel sticks

WarmGel® Infant Heel Warmer is a gel-filled, disposable heel warmer designed to provide gentle heat to increase blood flow and enhance infant blood sample collection.

WarmGel® Technology

Effective

- Conforms to infant's shape, keeping warmth where needed
- Even heat distribution increases blood flow

Safe

- Non-toxic food-grade gel won't harm infant's skin or eyes
- Single-patient use minimizes risk of infection

Reliable

- Reaches maximum temperature in less than 60 seconds, maintains optimal 104° F for 10 minutes¹
- Trigger disc minimizes auto-activation during shipment



WarmGel® Versus Liquid

- Conforms to infant's shape
- Provides even heat distribution
- Resists migration, keeping warmth where needed
- Competitively priced compared to liquid-filled products

| Item # | Description |
|--------|---------------------------------------|
| 24401 | Liquid Heel Warmer (100/Box) |
| 20418 | WarmGel® Infant Heel Warmer (100/Box) |

¹"Neonatal Thermoregulation: Guidelines for Practice," National Association of Neonatal Nurses, 1997.



TransWarmer®

Transport Mattress

Consistent, Safe Warming for Newborns



The TransWarmer® Infant Transport Mattress provides safe and reliable heat for up to two hours when cold stress is a concern. This thermostable, nontoxic gel mattress nests the infant in soothing warmth, while minimizing vibration during transport. Based on a 24 °C/75 °F start temperature, the mattress reaches a peak temperature of 40 °C/104 °F in less than 60 seconds - quickly improving the infant's admission temperature and reducing the risk of hypothermia.



Features and Benefits

- Comfortable, gel-filled mattress cradles baby in therapeutic warmth
- Gently conforms to the patient's shape for better nesting
- Minimizes the risk of hypothermia in cesarean deliveries
- Better heat distribution increases blood flow
- Disposable heat pack minimizes risk of cross contamination and infection
- Patented disc mechanism is easy to activate
- TransWarmer is nontoxic, latex- and PVC-free
- Use in labour and delivery, NICU, nursery, OR, ER, and transport units
- 10" wide x 16" tall

| Item # | Description |
|--------|--|
| 20421 | TransWarmer® Infant Transport Mattress (6/Box) |



DISTRIBUTED BY:
ADVANCED SURGI-PHARM INC.
www.advancedsurgipharm.com
Tel: 1-800-661-5432

surfcath™

Respiratory Distress Syndrome (RDS)



RDS is a pulmonary disorder resulting from a surfactant deficiency which commonly occurs in infants whose lungs have not yet fully developed.



In Europe, RDS is observed for about 90% of babies born at 24 weeks of gestation and for 80% of babies born at 28 weeks of gestation.¹

International recommendations¹

“Preterm infants should be managed without mechanical ventilation where possible”

“CPAP with early rescue surfactant is considered optimal management for babies with RDS”

“LISA is the preferred mode of surfactant administration for spontaneously breathing babies on CPAP, provided that clinicians are experienced with this technique”

European Consensus Guidelines on the Management of RDS - 2019


LISA method: Less Invasive Surfactant Administration

LISA method consists of a surfactant administration through a thin catheter inserted with Magill forceps through the vocal cords while maintaining a non-invasive ventilation.

Clinical Evidence²

Objective: A systematic review of 6 randomized controlled trials, enrolling a total of 895 preterm infants, comparing LISA method with surfactant delivery methods using an endotracheal tube.

Meta-analysis of the clinical outcomes :




Conclusion: “LISA technique for surfactant delivery results in a lesser need for mechanical ventilation in infants with RDS, reduction in the composite outcome of death or BPD at 36 weeks, and BPD₃₆ among survivors”

References

- David G. Sweet et al., European Consensus Guidelines on the Management of Respiratory Distress Syndrome – 2019 Update
- Aldana-Aguirre JC, Pinto M, Featherstone RM, et al. Arch Dis Child Fetal Neonatal Ed 2017;102:F17– F23

A catheter for surfactant administration using the LISA technique

For babies needing treatment for respiratory distress syndrome (RDS), our new surfcath™ uses the LISA (Less Invasive Surfactant Administration) technique to place the catheter. It was specifically designed to improve manoeuvrability during placement whilst also eliminating the need for the commonly used Magill forceps.



2cm soft distal black tip

- 2cm mark shows when surfcath™ is in place
- Soft tip minimises risk of tracheal lesions and prevents kinking
- Pre-curved to follow the airway anatomy and eases the passage between the vocal cords

Transparent bendable thermoplastic material

- Transparent material allows visual check on the delivery of surfactant
- Thermosensitive material allows you to curve surfcath™ prior to use meaning no need for Magill forceps
- Semi-rigid material also allows high manoeuvrability, helping to follow the airway anatomy


Less invasive 6Fr gauge

- No obstruction of the airways allowing spontaneous breathing
- Low dead space (0.2mL)



20cm length with cm markings

- Longer length eases manipulation of the surfactant syringe away from the patient's head
- Centimetre markings indicate the inserted length and check that surfcath™ stays in place



Scan to watch our video
on how to use surfcath™

Please note: The link will take you to an external website

| Code | Size | Length | Priming vol | Units/Box |
|----------|------|--------|-------------|-----------|
| 5590.106 | 6 Fr | 20 cm | 0.2 ml | 10 |

endotracheal tubes

Single and double lumen tubes in PVC (DEHP-free) for emergency oral or nasal intubation and routine anesthesia. Double lumen tubes have secondary lumen for surfactant administration or airway monitoring.

- Radiopaque
- 15mm adapter
- Markings every half centimeter

Single lumen tubes are also available as a soft tube designed to minimize tracheal damage during long term intubation.



Plain endotracheal tubes

| Product Code | | Ext. ø | Int. ø | Length | Units/Case |
|--------------|-----------|--------|--------|--------|------------|
| Standard | Soft tube | | | | |
| 516.20 | 520.20 | 3.4 mm | 2.0 mm | 165 mm | 25 |
| 516.25 | 520.25 | 4.1 mm | 2.5 mm | 165 mm | 25 |
| 516.30 | 520.30 | 4.6 mm | 3.0 mm | 165 mm | 25 |
| 516.35 | 520.35 | 5.2 mm | 3.5 mm | 165 mm | 25 |
| 516.40 | 520.40 | 5.7 mm | 4.0 mm | 230 mm | 20 |
| 516.45 | 520.45 | 6.2 mm | 4.5 mm | 230 mm | 20 |

Endotracheal tubes with secondary lumen

| Product Code | Tube | | | Secondary lumen | | Units/Case |
|--------------|--------|--------|--------|-----------------|-------------|------------|
| | Ext. ø | Int. ø | Length | Lumen diam. | Flow rate | |
| 5516.20 | 3.4 mm | 2.0 mm | 165 mm | 0.40 mm | 2.15 ml/min | 20 |
| 5516.25 | 4.1 mm | 2.5 mm | 165 mm | 0.55 mm | 3.35 ml/min | 20 |
| 5516.30 | 4.6 mm | 3.0 mm | 165 mm | 0.65 mm | 5.00 ml/min | 20 |
| 5516.35 | 5.2 mm | 3.5 mm | 165 mm | 0.80 mm | 15.0 ml/min | 20 |
| 5516.40 | 5.7 mm | 4.0 mm | 230 mm | 0.95 mm | 35.0 ml/min | 20 |

mucus extractors

Mucaid allows for quick and safe suctioning and sampling of the upper respiratory tract in the newborn using mechanical / wall suction. This avoids the dangers of self-contamination inherent in oral suctioning (De Lee method).

Included in the range are Mucaid with a secure threaded cap for safe transportation of sample for testing and Mucaid with a filter to prevent any mucus or amniotic fluid from passing through into the line.

Meconium aspirators allows for suctioning below the vocal cords to clear meconium from the lungs and airways.

- Finger tip control
- Endotracheal tube and endotracheal tube hub
- Rounded, soft tip
- Centimeters every half centimeter



Mucaid



Meconium aspirator

Mucaid: Mucus extractor

| Product Code | Volume | Tubing to suction source length | Catheter | Units/Case |
|--------------|--------|---------------------------------|---------------|------------|
| 542.06 | 25 ml | 24 cm | 6 Fr - 33 cm | 25 |
| 542.08 | 25 ml | 24 cm | 8 Fr - 33 cm | 25 |
| 542.10 | 25 ml | 24 cm | 10 Fr - 33 cm | 25 |

Mucaid with filter

| Product Code | Volume | Tubing to suction source length | Catheter | Units/Case |
|--------------|--------|---------------------------------|---------------|------------|
| 80.542.10B | 25 ml | 24 cm | 10 Fr - 33 cm | 25 |

Meconium aspirator

| Product Code | Int. Ø | Length | Units/Case |
|---------------|--------|--------|------------|
| 80.516.30D-25 | 3.0 mm | 165 mm | 25 |
| 80.516.35D-25 | 3.5 mm | 165 mm | 25 |

nasal ventilation

Single ventilation tube

Single silicone nasal tubes designed to reduce the stress related to ventilation and to improve the comfort of premature babies and newborns.

- Open tip
- 15mm adapter



Single ventilation tube

Double tube in silicone for nasal neonatal ventilation

Length of tube on 2 branches: 6 cm, except code 2596.30: 7 cm.

Open and rounded distal tips and centimetric markings.

The proximal end has an ISO standard 15 mm adaptor. This tube is specially designed to avoid tracheal intubation which is very traumatic, and to improve the comfort of premature babies and neonates.

4 sizes are available: Ext. Ø 2.7 mm - 3.4 mm - 4.0 mm - 5.0 mm



Double tube in silicone for nasal neonatal ventilation


| Product Code | Int. ø | Ext. ø | Diameter | Length | Units/Case |
|--------------------|--------|--------|----------|--------|------------|
| Single Tube | | | | | |
| 2595.15 | 1.5 mm | 2.7 mm | 8 Fr | 9 cm | 50 |
| 2595.20 | 2.0 mm | 3.4 mm | 10 Fr | 9 cm | 50 |
| 2595.25 | 2.5 mm | 4.0 mm | 12 Fr | 9 cm | 50 |
| 2595.30 | 3.0 mm | 5.0 mm | 14 Fr | 9 cm | 50 |
| Double Tube | | | | | |
| 2596.15 | 1.5 mm | 2.7 mm | 8 Fr | 6 cm | 50 |
| 2596.20 | 2.0 mm | 3.4 mm | 10 Fr | 6 cm | 50 |
| 2596.25 | 2.5 mm | 4.0 mm | 12 Fr | 6 cm | 50 |
| 2596.30 | 3.0 mm | 5.0 mm | 15 Fr | 7 cm | 50 |

Manual Resuscitator / Ventilator

SAFE SPOT

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE

Ventlab



SafeSpot™
SS3200-Series
Shown with integrated pressure manometer, pop-off valve, inflatable mask and inflatable O₂ reservoir

Safe Spot™

- Single patient use resuscitation/ventilation system designed specifically for infant patients
- Bag features ultra-thin 'spots' to provide an enhanced tactile feel for lung compliance and increased patient safety
- Integrated pressure manometer promotes proper delivery (<20 cm H₂O) to guard against over-pressurization of babies' fragile lungs; the tri-color fields alert to the risk level of delivered pressure
- Accurate pressure delivery helps reduce barotrauma or pneumothoraces
- Pop-off valve (40 cm H₂O and 25 cm H₂O available) enables clinicians to use the pop off level appropriate for a given patient population creating a greater margin of safety for the baby
- Optional Vent Mask in various sizes available
- Select from a variety of components to assemble a custom resuscitation bag to meet your unique patient needs

| ITEM# | SIZE | VOLUME | STROKE VOLUME | BODY MASS | PK |
|----------|--------|--------|---------------|------------------|----|
| SS3200MB | Infant | 250 mL | 150-170 mL | ≤ 10 kg (22 lbs) | 10 |

CUSTOMIZE YOUR MANUAL RESUSCITATOR

Choose from the following components to create a custom resuscitation/ventilation bag to best meet your needs:

SAFETY COMPONENTS



INTEGRATED PRESSURE MANOMETER

Ventlab's innovative color-coded dial clearly denotes pressure delivery:

Green - Target pressure level <math><20\text{ cm H}_2\text{O}</math>

Yellow - Potentially dangerous pressure level

Red - Dangerous pressure level



POP-OFF VALVES

40 cm H₂O & 25 cm H₂O

Allows clinician to use an appropriate pop-off level



FILTERS

A bacterial/viral filter provides effective protection for the healthcare provider during expiration



CO₂ INDICATOR

Choose from STAT-Check II or MaxCap for effective assessment of CO₂ exchange



PEEP VALVES

Adjustable PEEP valves;
5-20 cm/H₂O

SafeSpot™, SS3200-Series

Shown with integrated pressure manometer, pop-off valve, inflatable mask and corrugated tubing



SafeSpot™ bags come standard with an inflatable Vent Mask

OXYGEN RESERVOIRS



INFLATABLE O₂ BAG

Closed system oxygen reservoir delivers higher FDO₂ and provides visual indication of proper airflow



CORRUGATED TUBING

For use with secondary oxygen supply; available in 24" and 48" lengths

CONVENIENCE COMPONENTS



FACE MASKS

Inflatable and non-inflatable face masks available



FLEX NECK

Removable, flexible joint allows corrugated tubing to be joined with the patient valve; also reduces the torque on ET tube

INCA®

Infant Nasal CPAP Assembly

INCA® Infant Nasal CPAP Assembly is the most respected Continuous Positive Airway Pressure (CPAP) device among clinicians. INCA affords all the benefits of CPAP/Bubble CPAP for newborns in your care.

Features enhanced precision fit for neonates:

- Five anatomically proportional nasal cannula sizes for neonates
- Clear, extra-soft tubing provides custom shaping and maximal stability
- Chin strap is available to help prevent air leakage around the mouth
- Available in Complete and Replacement Sets



INCA COMPLETE SET CA CE

Consists of a nasal cannula, two lengths of Stay-Flex tubing, two tubing adapters, one in-line "T" pressure adapter, a knit cap with hook/loop strap, two foam blocks and a hook/loop chin strap. One INCA Sizing Gauge is included in every box of five sets.




| Item # | Description |
|---------|----------------------------------|
| 44-2707 | Set with 7.5 Fr Cannula (5/Box) |
| 44-2709 | Set with 9 Fr Cannula (5/Box) |
| 44-2710 | Set with 10.5 Fr Cannula (5/Box) |
| 44-2712 | Set with 12 Fr Cannula (5/Box) |
| 44-2715 | Set with 15 Fr Cannula (5/Box) |

INCA ACCESSORY PACK CE

Includes five knit caps with hook/loop strap, five chin straps with ten hook fasteners (not attached) and ten foam blocks.

1 Pack = 5 Hats and 10 Foam Blocks



| Item # | Description |
|---------|--|
| 44-2001 | X-Small Pack (5 XS Caps, 10 Foam Blocks) |
| 44-2002 | Small Pack (5 S Caps, 10 Foam Blocks) |
| 44-2003 | Medium Pack (5 M Caps, 10 Foam Blocks) |
| 44-2004 | Large Pack (5 L Caps, 10 Foam Blocks) |

INCA REPLACEMENT SET CA CE

Includes the nasal cannula, two lengths of Stay-Flex tubing, two tubing adapters and one in-line pressure "T" adapter.



| Item # | Description |
|---------|----------------------------------|
| 44-0707 | Set with 7.5 Fr Cannula (5/Box) |
| 44-0709 | Set with 9 Fr Cannula (5/Box) |
| 44-0710 | Set with 10.5 Fr Cannula (5/Box) |
| 44-0712 | Set with 12 Fr Cannula (5/Box) |
| 44-0715 | Set with 15 Fr Cannula (5/Box) |

INCA NASAL SIZING GAUGE CA CE

For appropriate nasal prong fitting, use the INCA Sizing Gauge. Wipe gauge with sterile 70 percent isopropyl alcohol prior to each use. Hold the gauge up to the nasal area to determine appropriate cannula size.



| Item # | Description |
|---------|----------------------------------|
| 025-001 | INCA Nasal Sizing Gauge (1/Pack) |

NEO-fit™

Neonatal Endotracheal Tube Grip CA CE Facilitates Endotracheal Tube Management

The NEO-fit provides quick, secure endotracheal tube stabilization. Designed to firmly grip a Neonatal endotracheal tube, NEO-fit utilizes metal clips to mechanically grip the soft material of the tube. The device is positioned on baby's face using a hydrocolloid adhesive that is hypoallergenic and Hook / Loop strap that allows for rapid release and readjustment of the tube. MRI Conditional.



| Item # | Description |
|---------|---|
| 42-2540 | NEO-Fit Neonatal Endotracheal Tube Grip Set (20/Box) Each Set includes: (1) plastic tube grip, with hook/loop strap, (2) foam foot pads, (2) foam "lollipop" tapes. Non-sterile. |

NEO-prep™

Skin Barrier CE

A must-have product for your CPAP patients. Protect their tender tissues during this critical therapy. Neo-prep Skin Barriers are precut and non-sterile.



| Item # | Description |
|---------|---|
| 44-8000 | Precut Hypoallergenic Skin Barrier, non-sterile, small (10/Box) |
| 44-8001 | Precut Hypoallergenic Skin Barrier, non-sterile, large (10/Box) |

Os Finder™

Cervical Os and Canal Dilator

The Os Finder naturally seeks and gently opens the cervical Os. Made of an extremely smooth and malleable Teflon-like material, this instrument is ideal for use in intracervical dilation or uterine sounding, or for intrauterine procedures.

DEVICE DESCRIPTION

The sterile, single use Os Finder™ naturally seeks and gently opens the cervical Os. It is made of extremely smooth and malleable teflon-like material. The dilator is 21.5 cm long with a taper from 2 to 4 mm. CAUTION: Federal restricts this device to sale by or on the order of physician.

CONTRAINDICATIONS

This instrument should not be used under any of the following conditions:

- Pregnancy or suspicion of pregnancy
- Acute or active inflammation of the cervix, endometrium, fallopian tube, ovary or peritoneum
- In the case of cervicitis, endometritis, tubo-ovarian inflammatory disease or pelvic inflammatory disease
- Invasive cancer that is visible on examination

PRECAUTIONS

Do not use the Os Finder if the package or instrument are damaged.

INDICATIONS FOR USE

This instrument is ideal for use in intracervical dilation, uterine sounding, or for intrauterine procedures.

DIRECTIONS FOR USE

The Os Finder is used to find the external Os by simply pushing the narrow end gently against the cervix until it finds its way. Further insertion of the flexible tip will follow the natural path of the canal until it reaches and goes through the internal Os. Dispose of in accordance with all applicable Federal, State and local Medical/Hazardous waste practices.



| Item # | Description |
|--------|--------------------|
| 1176 | Os Finder (50/Box) |



Ver.2021-04-22

Tampa Infusion Catheter

For sonohysterography - the simple, cost-effective choice for ultrasound examination of the uterus.

- Unique introducer sheath (33 cm long) facilitates insertion into cervical Os with minimal trauma
- 7 cm marker provides optimal placement of catheter tip in uterus
- Balloonless configuration is ideal for postmenopausal patients with stenotic cervical canals, as well as for many nulliparous patients
- Small catheter size allows easy placement next to ultrasonography probe



| Item # | Description |
|---------|--|
| 61-2005 | Tampa Catheter (20/Box) 5 French, clear, flexible catheter with introducer sheath, 33 cm long, sterile |



Ver.2021-04-22

DISTRIBUTED BY:
ADVANCED SURGI-PHARM INC.
www.advancedsurgipharma.com
Tel.1-800-661-5432



H/S Catheters


For Enhanced Visualization of Uterine Pathology

H/S CATHETER SET CA CE

These easy-to-use devices provide enhanced visualization. Designed for use in sonohysterography and hysterosalpingography, the catheter can be used with oil- or saline-based media depending on the procedure requirements.

Features and Benefits

- Available with integrated stylet for easy insertion
- Unique placement sheath eliminates the need for a tenaculum
- Comes in both 5 and 7 French
- Ideal for office use
- Latex-Free construction



Pre-Saline Infusion

Post-Saline Infusion

| Item # | Description |
|---------|--|
| 61-5005 | H/S Catheter, Kraton, 5 Fr (10/Box) |
| 61-5007 | H/S Catheter, Kraton, 7 Fr (10/Box) |
| 61-0020 | H/S Stylet Catheter Stiffener (10/Box) |
| 61-3005 | H/S Catheter with integrated Stylet 5FR (10/Box) |
| 61-3007 | H/S Catheter with integrated Stylet 7FR (10/Box) |

Ver.2021-04-22



DISTRIBUTED BY:
ADVANCED SURGI-PHARM INC.
www.advancedsurgipharm.com
Tel:1-800-661-5432

H/S Elliptosphere™ Catheter Set and Tray

The Safe and Simple Answer for Both
Sonohysterography and Hysterosalpingography



The patented Elliptosphere H/S catheter offers several advantages over the traditional balloon catheters including the ability to gradually inflate the balloon for placement in any size cervical canal, the ability to see the entire lower uterine segment during an SIS procedure, and greater patient comfort.

The Elliptosphere balloon contour and variable size permit secure placement in the cervical canal or uterine cavity for optimal visualization. The balloon's unique design prevents retrograde leakage of distension media, which facilitates peak uterine distention. A small catheter size enables placement next to ultrasound probe for minimal patient discomfort. The H/S Elliptosphere Catheter can be used with oil- or saline-based medium. Available in two versions – with and without an integrated stylet – the H/S Elliptosphere features latex-free construction.



| Item # | Description |
|----------|--|
| 61-4005 | H/S Elliptosphere Catheter, 5 Fr (10/Box) |
| 61-4005S | H/S Elliptosphere Catheter Set with Integrated Stylet (10/Box) |



Comfort Belly Bands

Abdominal Fetal Monitoring Band

A lighter, cooler abdominal fetal monitoring band, Belly Bands will be appreciated by your patients as they hold vital tocos and transducers firmly in place. Belly Bands are available in a wide range of widths to comfortably cover any patient abdomen.



Features and Benefits

- Breathable material helps wick away moisture from the mother's skin
- Lighter and cooler to ensure patient comfort
- Flexible precut lengths cover the belly and secure transducer placement
- Available in 30-foot rolls or individually wrapped packages of precut 14" lengths
- Choose from 10-inch to 18-inch widths
- Economical, easy-to-use product

| Item # | Description |
|--------|--|
| 13010 | 10" Width Belly Bands on Roll (10" x 10 yds) (1 Roll) |
| 13014 | 14" Width Belly Bands on Roll (14" x 10 yds) (1 Roll) |
| 13015 | 14" Width Belly Bands - Precut (14" x 14") (25/Box) |




LAMINARIA

MedGyn's Laminaria is all-natural cervical disposable and sterile dilator made from a sea-grown plant containing no synthetic materials. Laminaria provides gentle, gradual dilation.

Features & Benefits

MedGyn's disposable and sterile Laminaria is all-natural cervical dilator made from a sea-grown plant containing no synthetic materials. Safely used for over a century, Laminaria provides gentle, gradual dilation. Once inserted into cervical canal, Laminaria absorbs liquid, swelling several times its initial dried diameter, gently dilating the cervix.



McMahon Laminaria Inserter

MedGyn Laminaria Inserter

Indication for use:

Laminaria should be considered whenever cervical dilation and/or softening of the cervix is desired. Some examples are: cervical stenosis, related to dysmenorrhea, considered a possible cause of infertility, resulting from cauterization or conization, placement and removal of intrauterine devices, induction of labor, radium placement, curettage and drainage of uterine cavity, endometrial biopsy, and suction cannula aspiration.

Standard, approx. length 60mm

021003 Laminaria 3mm Small

leadercath

Arterial LeaderCath (PE)

Although it can be inserted in a variety of arteries, the most common sites for the placement of this catheter are the radial and the femoral artery.

The Seldinger technique makes arterial puncture easier particularly in hypotensive patients with a weak peripheral pulse and when arteries are atherosclerosed or narrow.

The product is supplied in a rigid blister pack containing:

- 1 introducer needle
- 1 straight guidewire
- 1 transparent and radiopaque catheter in polyethylene with winged hub



| Product Code | Catheter | | | | | | Puncture Needle | | | | Guidewire | | Unit/Box |
|--------------|-----------|-----------|-----------|---------|------|------------------|-----------------|--------|---------|-----------|-----------|-----------|----------|
| | Int. Ø mm | Ext. Ø mm | Length cm | Gauge G | Ø Fr | Flow rate ml/min | Int. Ø mm | Ext. Ø | Gauge G | Length mm | Ø mm | Length mm | |
| 115.090 | 0.6 | 0.9 | 8 | 20 | 3 | 24 | 0.6 | 0.9 | 20 | 38 | 0.53 | 20 | 20 |
| 115.11 | 0.8 | 1.2 | 10 | 18 | 4 | 24 | 0.8 | 1.1 | 19 | 54 | 0.71 | 30 | 20 |

lectrospiral®

Coiled extension tube


Lectrospiral® is a flexible coiled extension tube that provides patients with freedom of movement.

Flexible

- Coiled tubing expands and contracts providing patients with greater freedom of movement than straight extension lines.
- Ideal to use with pediatric, home care and oncology patients.

Safe

- Polyethylene material allows for compatibility with commonly infused drugs and reduces adsorption, absorption or permeation that may occur with PVC.
- Extension tubing minimizes potential of trauma at IV site.



| Product Code | Length | Inner Diameter | Outer Diameter | Priming Volume | Units/Case |
|--------------|--------|----------------|----------------|----------------|------------|
| 1155.60 | 100 cm | 1.0 mm | 2.0 mm | 1.1 ml | 15 |
| 1155.65 | 150 cm | 1.0 mm | 2.0 mm | 1.5 ml | 15 |
| 1155.70 | 200 cm | 1.0 mm | 2.0 mm | 2.0 ml | 15 |
| 1155.80 | 300 cm | 1.0 mm | 2.0 mm | 2.8 ml | 15 |
| 1155.90 | 400 cm | 1.0 mm | 2.0 mm | 3.7 ml | 15 |
| 1159.60 | 100 cm | 1.5 mm | 2.5 mm | 1.9 ml | 15 |
| 1159.65 | 150 cm | 1.5 mm | 2.5 mm | 2.7 ml | 15 |
| 1159.70 | 200 cm | 1.5 mm | 2.5 mm | 3.6 ml | 15 |
| 1159.80 | 300 cm | 1.5 mm | 2.5 mm | 5.2 ml | 15 |
| 1159.90 | 400 cm | 1.5 mm | 2.5 mm | 6.9 ml | 15 |


Ver.2021-05-20

lectrocath

A range of polyethylene extension lines, suitable for arterial pressure monitoring and venous medication delivery.

The advantages include:

- a Polyethylene material which allows for compatibility with commonly used drugs and reduces adsorption, absorption or permeation that may occur with PVC
- Latex-Free
- DEHP-Free
- Pyrogen-Free



| Feature | Benefit |
|--------------------|---------------------------------------|
| PE tubing | Drug-compatible |
| Colour-coded clips | For arterial or venous identification |

| Product Code | ID - OD | Length | Capacity | Units/Box |
|--------------|--------------|--------|----------|-----------|
| 1155.01 | 1.0 - 2.0 mm | 11 cm | 0.4 ml | 50 |
| 1155.03 | 1.0 - 2.0 mm | 30 cm | 0.6 ml | 50 |
| 1155.05 | 1.0 - 2.0 mm | 50 cm | 0.7 ml | 50 |
| 1155.10 | 1.0 - 2.0 mm | 100 cm | 1.2 ml | 50 |
| 1155.15 | 1.0 - 2.0 mm | 150 cm | 1.6 ml | 50 |
| 1155.20 | 1.0 - 2.0 mm | 200 cm | 2.0 ml | 50 |
| 1159.01 | 1.5 - 2.5 mm | 11 cm | 0.5 ml | 50 |
| 1159.03 | 1.5 - 2.5 mm | 30 cm | 0.8 ml | 50 |
| 1159.05 | 1.5 - 2.5 mm | 50 cm | 1.1 ml | 50 |
| 1159.10 | 1.5 - 2.5 mm | 100 cm | 2.0 ml | 50 |
| 1159.15 | 1.5 - 2.5 mm | 150 cm | 2.8 ml | 50 |
| 1159.20 | 1.5 - 2.5 mm | 200 cm | 3.6 ml | 50 |

biocap

Needle-free drug aspiration; for multiple vial access

- Contains Latex: No
- Contains DEHP: No
- Contains biological or animal-based product: No
- Non-pyrogenic: Yes

| Product Code | Cap Ø | Units/Box |
|--------------|-------|-----------|
| 821.21 | 14 mm | 50 |
| 821.22 | 20 mm | 50 |



bionector

Bionector is a multi-purpose, closed, needle-free connection device.

When Bionector is in place, you can infuse, inject sample and change the I.V. tubing without the I.V. circuit to the atmosphere. Bionector's membrane automatically opens the fluid pathway only when a male Luer is disconnected, the membrane automatically seals the fluid pathway.

- Bionector preserves constantly a closed system
- 7 days or 360 connections
- Straight fluid pathway
- Neutral displacement of fluid
- Minimal dead space
- Can be used with power injectors (CT-scan).
Pressure resistance: 350 psi
- MRI conditional
- Contains Latex: No
- Contains DEHP: No
- Contains biological or animal-based product: No
- Non-pyrogenic: Yes




| Product Code | Flow rate | Priming volume | Equivalent Gauge G | Number of accesses | Max. duration of use | Length | Diameter | Units/Box |
|--------------|------------|----------------|--------------------|--------------------|----------------------|---------|----------|-----------|
| 896.038 | 105 ml/min | 0.03 ml | - | 360 | 7 days | 25.7 mm | 10 mm | 50 |

Ver.2021-04-22

IV Accessories

Ramp of stopcocks (1 way)

| Product Code | Description | Extension tube | | | Lateral accesses | Units/Box |
|--------------|---------------------------|----------------|--------|--------|------------------|-----------|
| | | Ext. Ø | Int. Ø | Length | number | |
| 873.24 | Ramp of stopcocks (1 way) | 4.0mm | 2.5 mm | 200 cm | 4 | 15 |



Ramp of Stopcocks (1 way): female Luer-lock/female Luer-lock

Stopcock

| Product Code | Description | Units/Box |
|--------------|---------------------------|-----------|
| 872.10 | Lockable one-way stopcock | 50 |



Lockable one-way stopcock

Connectors

| Product Code | Description | Units/Box |
|--------------|--|-----------|
| 892.00 | Female Luer-lock to female Luer-lock connector | 30 |
| 893.00 | Male Luer Lock to Male Luer Lock Connector | 30 |



Female Luer Lock Connector



Male Luer-lock to male Luer-lock connector

Extension tubes (single lumen)

| Product Code | Description | Extension tube | | | Units/Box |
|--------------|---|----------------|--------|--------|-----------|
| | | Ext. Ø | Int. Ø | Length | |
| 220.050 | Lectroflex extension tube | 4.0 mm | 2.5 mm | 50 cm | 25 |
| 221.100 | Male Luer-lock to male Luer-lock extension tube | 4.0 mm | 2.5 mm | 100 cm | 25 |



Lectroflex



Male Luer-lock to male Luer-lock extension tube

Rapid Infusion

INFU-SURG® PRESSURE INFUSION BAG

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE

Ethox



Three-way stopcock ensures precise control of pressure while the large oval bulb allows for quick bladder inflation, all with ONE hand



Solution bag hook enables easy bag changing without having to remove the Infu-Surg from the IV pole; keeps fluid bag in proper position while in use



360° view of the color-coded gauge facilitates accurate pressure monitoring while the safety pressure relief valve prevents over-inflation

- Ideal for intravenous and intra-arterial infusion treatment including A-line pressure monitoring
- Pressure relief valve prevents over inflation (330 mm Hg pressure relief)
- Large, oval shaped bulb allows for quick and easy inflation of bladder
- Single-handed inflation and deflation design makes it easy to use and requires minimal training
- Dedicated fluid bag hook enables efficient loading and unloading of the system without having to remove the infusion bag from the IV pole
- Reinforced hanging loop allows for quick and easy hanging on IV pole
- Color-coded gauge makes for accurate pressure monitoring (150-300 mm Hg)
- Three-way stopcock ensures precise control of pressure
- Piston-style gauge allows for 360° visibility

| ITEM# | DESCRIPTION | VOLUME | PK |
|-------|----------------------|---------|----|
| 4005H | Infu-Surg® with Hook | 500 mL | 5 |
| 4010H | Infu-Surg® with Hook | 1000 mL | 5 |
| 4030 | Infu-Surg® | 3000 mL | 5 |

Ver.2021-04-22

DISTRIBUTED BY:
 ADVANCED SURGI-PHARM INC.
 www.advancedsurgipharm.com
 Tel.1-800-661-5432



Rapid Infusion

INFU-SURG® CLARITY PRESSURE INFUSION BAG

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE

Ethox



Clear panel enables full visibility of solution bag to monitor fluid level and easily scan patient barcode



- Transparent panel for clear visibility of solution bag
- Ideal for intravenous and intra-arterial infusion treatment including A-line pressure monitoring
- Dedicated hook enables quick change of solution bag while providing proper fluid bag positioning in pocket
- Pressure relief valve prevents over-inflation (330 mm Hg)
- Color-coded gauge provides accurate pressure monitoring (150-300 mm Hg)
- Single-handed inflation and deflation with oval-shaped bulb makes it easy to use and requires minimal training
- Welded panel and hook construction are more durable than sewn designs to withstand repeated use
- Reinforced hanging loop allows for quick and easy hanging on IV pole
- Three-way stopcock ensures precise control of pressure
- Piston-style gauge allows for 360° visibility



Three-way stopcock ensures precise control of pressure while the large oval bulb allows for quick bladder inflation, all with ONE hand



Solution bag hook enables easy bag changing without having to remove the Infu-Surg from the IV pole; keeps fluid bag in proper position while in use



360° view of the color-coded gauge facilitates accurate pressure monitoring while the safety pressure relief valve prevents over-inflation

| ITEM | DESCRIPTION | VOLUME | PK |
|------|--------------------|---------|----|
| 4050 | Infu-Surg® Clarity | 500 mL | 5 |
| 4100 | Infu-Surg® Clarity | 1000 mL | 25 |

Ver.2021-04-22



DISTRIBUTED BY:
ADVANCED SURGI-PHARM INC.
www.advancedsurgipharm.com
Tel.1-800-661-5432

Jet-Vac Bottles

Needle-Free Luer Locking Connection

- Removes need for spikes or needles to access the bottle
- Sample fluids aseptically

Shatter-Proof Plastic Design

- Reduces risk for hospital personnel and patients by avoiding glass breaks
- Reduces hazardous waste disposal costs compared to glass containers

Easy to Read Vacuum Indicator

- Provides visual confirmation of the vacuum inside the bottle
- Indicates when a change of the bottle is necessary



JET-VAC 600

- Ideal for thoracentesis procedures
- Compact size is optimal for sample transfer between hospital departments
- Patent pending**



JET-VAC 1000

- Ideal for paracentesis procedures
- Excellent substitute for glass bottles
- Patent pending**



JET-VAC 1500

- Ideal for large volume paracentesis procedures
- Large size increases efficiency of procedure and reduces waste
- Patented design**

Ver.2023-02-02

**U.S. Design Patent No. D818119; Canadian Design Patent No. 172153

DISTRIBUTED BY:
ADVANCED SURGI-PHARM INC.
www.advancedsurgipharm.com
Tel.1-800-661-5432



Ordering Information

Jet-Vac Bottles

| Item No. | Item | Description | Case Qty |
|----------|---------------|--|----------|
| M7090 | Jet-Vac 1000 | 1,000 mL Evacuated Drainage Bottle in Valve-Only Configuration | 10/Case |
| M7092 | Jet-Vac 600 | 600 mL Evacuated Drainage Bottle in Valve-Only Configuration | 10/Case |
| M7093 | Jet-Vac 1500 | 1,500 mL Evacuated Drainage Bottle in Valve-Only Configuration | 10/Case |
| M7094 | Jet-Vac 1000D | 1,000 mL Evacuated Drainage Bottle with Drain Line | 10/Case |
| M7095 | Jet-Vac 1500D | 1,500 mL Evacuated Drainage Bottle with Drain Line | 10/Case |
| M7096 | Jet-Vac 600D | 600 mL Evacuated Drainage Bottle with Drain Line | 10/Case |

Accessories

| Item No. | Description | Box Qty |
|----------|--|---------|
| M7005 | Drainage Line Set (Includes Drainage Line with End Caps, and 5-in-1 Adapter) | 10/Box |
| M7098 | Jet-Vac Large Bore Drainage Line | 5/Box |
| M7099 | Jet-Vac Drainage Line | 10/Box |





Jet-Vac Centesis Catheters, Jet-Vac Procedure Trays and Jet-Vac Bottles provide a complete drainage system for physicians, nurses, and hospitals. The Jet-Vac drainage bottle is a disposable vacuum device used to aspirate, remove, or sample body fluids.

Jet-Vac Centesis Catheters

Large circular side holes designed to increase fluid drainage

Tapered catheter tip

- Provides a smooth transition
- Helps minimize patient discomfort

HDPE material

- Reduces catheter collapsing
- Promotes a smooth catheter insertion

Clear needle hub with catheter lock


- Increases fluid visibility
- Ergonomic design provides a secure grip
- Prevents needle/catheter separation

Plastic vent plug

- Allows for air to exit
- Provides added safety for physician

Echogenic needle

- Ensures needle visibility during ultrasound placement



Ordering Information

Jet-Vac Centesis Catheters


| Item No. | Catheter Size | Length | Needle Gauge | Box Qty |
|----------|---------------|--------|--------------|---------|
| M7511 | 4F | 7 cm | 20G | 5/Box |
| M7512 | 4F | 10 cm | 20G | 5/Box |
| M7513 | 5F | 7 cm | 19G | 5/Box |
| M7514 | 5F | 10 cm | 19G | 5/Box |
| M7515 | 5F | 15 cm | 19G | 5/Box |

ASEPT®

Pleural & Peritoneal Drainage System

ASEPT Catheter and Insertion Kits

The ASEPT Pleural and Peritoneal Drainage Systems are designed for long-term, intermittent drainage of symptomatic, recurrent fluid accumulation.



- | | | | |
|---|--|--------------------------------|------------------------------|
| 1. ASEPT Catheter | 6. 16G x 1.0" (25.4mm) Needle | 12. Introducer Needle | 18. Tunneler |
| 2. ASEPT Drainage Line | 7. 25G x 1.0" (25.4mm) Needle | 13. Tearaway Introducer | 19. Forceps |
| 3. 5-in-1 Drain Line Adapter | 8. 22G x 1.5" (38.1mm) Needle | 14. 8F (2.67mm) Dilator | 20. Safety Scalpel |
| 4. Chloraprep® 10.5mL Vial | 9. Tegaderm® Self Adhesive Dressing | 15. 12F (4.0mm) Dilator | 21. Fenestrated Drape |
| 5. Synthetic Suture, 2-0 Straight and Synthetic Suture, 3-0 Curved | 10. Gauze Pads, 4"x4" | 16. 10mL Syringe | |
| | 11. Foam Catheter Pad | 17. J-Tip Guidewire | |

The primary component of the ASEPT Drainage system is the catheter, which has an aseptic, self-sealing valve attached that provides a secure and reliable connection. The valve can be connected to an ASEPT vacuum drainage bottle, drainage bag, or a stand-alone drainage line to connect to other fluid collection equipment.

Features:

- Secure connection
- Cap-less and Aseptic design
- Needle-free and swabable
- Easy to replace
- Self-sealing

ASEPT Drainage Bottles and Bags

The ASEPT Drainage kit configuration includes a 600mL or 1000mL vacuum bottle set, a 2000mL gravity drainage bag, and the procedure pack with all of the dressing change items needed to complete an at-home drainage procedure. A standalone drain line is available to connect to the ASEPT drainage catheter for drainage with other collection equipment during in-patient procedures.



Bottle Features:

- Convenient vacuum release clamp
- Green vacuum indicator
- Luer-lock connector



Bag Features:

- Large 2000mL volume
- Fluid release valve for ease of bag emptying

ASEPT®

ASEPT Accessories

PFM Medical, Inc. is proud to offer a full line of accessories ensuring that the physician or patient will have everything or she needs to complete the drainage procedure. Our accessory portfolio contains the following items:

ASEPT Replacement Valve



ASEPT Drainage Line Set



ASEPT Procedure Pack



Ordering Information

ASEPT Pleural/Peritoneal Insertion Kits

| REF | Description | Box Qty (ea) |
|-----------|-----------------------------------|--------------|
| M7001-CAN | Pleural Drainage Insertion Kit | 1 |
| M7008-CAN | Peritoneal Drainage Insertion Kit | 1 |

ASEPT Drainage Bottles & Bags

| | | |
|-------|--|----|
| M7002 | 600 mL Drainage Kit* | 10 |
| M7080 | 600 mL Drainage Starter Kit** | 5 |
| M7072 | 600 mL Bottle Only | 10 |
| M7000 | 1000 mL Drainage Kit* | 10 |
| M7081 | 1000 mL Drainage Starter Kit** | 5 |
| M7070 | 1000 mL Bottle Only | 10 |
| M7060 | 2000 mL Drainage Bag Starter Kit** (Peritoneal Use Only) | 10 |
| M7058 | 2000 mL Drainage Bag (Peritoneal Use Only) | 10 |

ASEPT Accessories

| | | |
|-------|--|----|
| M7040 | Procedure Pack (Traditional Dressing Change) | 10 |
| M7004 | Catheter Replacement Valve | 5 |
| M7005 | Drainage Line Set (includes 5-in-1 adapter) | 10 |

*ASEPT Drainage Kit includes bottle and traditional dressing change procedure pack.

** ASEPT Drainage Starter Kit recommended for patient discharge. Includes bottles & traditional dressing change procedure packs.

Indications for Use

ASEPT Pleural Indications

The ASEPT Pleural Drainage System is intended for long term, intermittent drainage of symptomatic, recurrent, pleural effusions, including malignant pleural effusions and other pleural effusions that do not respond to treatment of the underlying disease.

ASEPT Peritoneal Indications

The ASEPT Peritoneal Drainage System is indicated for periodic drainage of recurrent and symptomatic malignant ascites. The catheter is intended for long term access of the peritoneal cavity in order to relieve symptoms such as dyspnea.

Bronchitrac® “L” Set

There is no need to twist the catheter or patient's body. Simply advance the Bronchitrac® “L” through a properly placed endotracheal or tracheostomy tube. The dual-plane, patented curvature of the catheter ensures consistent entry into the left main-stem bronchus.

The time savings and reduced patient trauma provided by the Bronchitrac® “L” make it the obvious choice when suctioning the left-bronchial airway.

- Ideal for suctioning of post-cardiac surgical, neuro and critical care patients
- Can be used for deep-lobe suctioning
- Completely latex-free



| Item # | Description |
|---------|--|
| 41-7814 | 21" long, 14 Fr Bronchial Suction Catheter, Slide-cover Connector and (2) Vinyl Gloves, Sterile (20/Box) |

TRAKE-fit™ Collar Set

Tracheostomy Tube Holder For Maximum Patient Comfort

- Available in single- or double-collar sets
- Installs in seconds, remains secure
- Maximum comfort with cotton-lined neckband
- No adhesive, no narrow tape
- Completely adjustable
- Narrow tabs available for pediatric and small adult use
- Third-party reimbursable



| Item # | Description |
|---------|---|
| 43-1850 | Double Collar Set (25/Box) Each Set includes: (2) neckbands, (2) wide hook/loop tabs |
| 43-1851 | Single Collar Set (25/Box) Each Set includes: (1) neckband, (2) wide hook/loop tabs |
| 43-1862 | Double Collar/Narrow Tab Set (25/Box) Each Set includes: (2) neckbands, (2) narrow hook/loop tabs |



COMFIT™ Endotracheal Tube Holder

The Tapeless Way to Secure an Endotracheal Tube



COMFIT™ is a reliable, tapeless solution used to secure endotracheal tubes and better manage patient airways. Physicians appreciate COMFIT'S compliance with Advanced Cardiac Life Support (ACLS) Guidelines that recommend endotracheal tube stabilization devices as a substitute for tape. For patients, COMFIT™ provides a comfortable fit without the skin irritation that tape can cause.

Features and Benefits

- Significantly decreases the incidence of self-extubation¹
- Wide cotton-lined neckband minimizes skin irritation, providing maximum patient comfort
- No tape on the tube or the patient
- Completely adjustable
- Minimal plastic loop around the endotracheal tube allows access to the oral cavity
- A low-cost device that requires infrequent changing
- Latex-free product
- Complies with American Heart Association's ACLS Guidelines recommending the use of commercially available endotracheal tube holders instead of tape
- Plastic double securing loop, ½- x 19-inch neckband



1 McAndrew J, RR. Reduction of self extubation using Comfit™ tracheal tube holder. Respiratory Care. November 1994; 39 (11).

| Item # | Description |
|---------|---|
| 42-1820 | Tapeless Endotracheal Tube Holder (25/Box) Each set includes: (1) plastic tube grip, with hook/loop strap, foam foot pads, (2) foam "lollipop" tapes |




DISTRIBUTED BY:
ADVANCED SURGI-PHARM INC.
www.advancedsurgipharma.com
Tel.1-800-661-5432



Esophageal Balloon Catheter Set

For Intrathoracic Pressure Monitoring

- Ease of Insertion - The catheter comes with a new, easy-to-remove PTFE-coated stylet to aid the insertion procedure
- Comfort - The 5 French catheter is soft and pliable, offering maximum patient comfort over extended time periods
- Ease of Positioning - The new proximal Y-connector allows for precise balloon positioning with the stylet in place
- Depth Markings - The catheter is radio opaque with numerical depth markings to facilitate positioning
- Ease of Connection - The catheter has a standard female Luer connector
- Several Applications - Ideal for monitoring ventilator and sleep apnea patients and calculating lung compliance



| Item # | Description |
|---------|---|
| 47-9005 | 5 Fr Catheter, latex-free Balloon, Stylet, 36" Pressure Monitoring Extension, 3-way Stopcock (5/Box) |



Gastric Lavage

TUM-E-VAC® GASTRIC LAVAGE KIT

LATEX FREE | SINGLE USE | DISPOSABLE

Ethox



- Pre-assembled kits are quick and easy to use, reducing the time for lavage procedures
- Promotes quick, efficient lavage of stomach contents
- Closed system helps prevent exposure to stomach contents and charcoal solutions
- Available with activated charcoal suspension (with or without Sorbitol)
- Twist-Lok™ expressing key aids in charcoal administration
- Easy-to-fill top port
- Available with Universal Solution Spike (for use with pre-filled solution containers)
- No clean up or disinfecting procedures necessary
- All kits come standard with Bite Block

WITHOUT CHARCOAL

| ITEM# | DESCRIPTION | PK |
|-------|--|----|
| 2054 | Tum-E-Vac + 3L Solution Bag + 32Fr Lavage Tubing | 10 |
| 2055 | Tum-E-Vac + 3L Solution Bag | 10 |
| 2075 | Tum-E-Vac + Universal Solution Spike | 10 |

WITH ACTIVATED CHARCOAL


| ITEM# | DESCRIPTION | PK |
|-------|--|----|
| 2057 | Tum-E-Vac + 32Fr Lavage Tubing + (2) 25g Activated Charcoal | 10 |
| 2058 | Tum-E-Vac + 32Fr Lavage Tubing + (2) 25g Activated Charcoal with Sorbitol | 10 |
| 2059 | Tum-E-Vac + 32Fr Lavage Tubing + (1) 25g Activated Charcoal with Sorbitol + (1) 25g Activated Charcoal | 10 |
| 2081 | Tum-E-Vac + Universal Solution Spike + (2) 25g Activated Charcoal with Sorbitol | 10 |

LAVAGE TUBING

LATEX FREE | SINGLE USE | DISPOSABLE

Ethox

Radiopaque strip confirms proper tube placement



Reference marks at 18", 20" and 22" from distal end




- Open bullet tip design with multiple eyelets to create free-flowing movement with low resistance
- Soft clear tube with smooth tip
- Suction connector on proximal end
- Radiopaque stripe
- All tubing is 48" long with reference marks at 18", 20" and 22" from distal end
- Individually packaged

| ITEM# | SIZE | DESCRIPTION | PK |
|-------|------|-------------------|----|
| 2124 | 24Fr | 48" Lavage Tubing | 10 |
| 2128 | 28Fr | 48" Lavage Tubing | 10 |
| 2132 | 32Fr | 48" Lavage Tubing | 10 |
| 2136 | 36Fr | 48" Lavage Tubing | 10 |

Irrijet®

Splatter-Gard™

Simple, one-hand, thumb-controlled pressure provides an appropriately forceful stream (8 psi) directly to the desired area (or a gentler stream if needed). The soft flexible tip has an inside diameter equivalent to an 18-19 gauge needle.



- Splatter-Gard™ prevents contamination from body secretions – new concave design provides maximum splash-back protection
- Soft, flexible plastic tip eliminates needle-stick for user, minimizes trauma to patient
- Irrigant is delivered with force equivalent up to 8 psi
- Latex-free
- Sterile


| Item # | Description |
|---------|------------------------------|
| 55-0004 | Irrijet GP (5/Box) |
| 55-0504 | Multi-Use Irrijet (5/Box) |
| 55-2004 | Single-Use Irrijet Z (5/Box) |
| 55-2104 | Single-Use Irrijet (5/Box) |
| 55-2134 | Single-Use Irrijet (5/Box) |
| 55-6504 | Irrijet DS (5/Box) |
| 59-5014 | Infuset GP (5/Box) |



Irrijet GP

55-0004


Spring-loaded syringe connected to a double check-valve, through a side arm and 6-ft. length of tubing to a spike for entry into a bag of irrigating solution. Double check-valve prevents retrograde flow. Alternate movement of the syringe plunger delivers irrigating solution as needed. Syringe Capacity: 12cc.



Multi-Use Irrijet

55-0504


Same basic assembly as 55-0004, plus one package of five individually sealed jet tips.



Single-Use Irrijet Z

55-2004

Same basic assembly as 55-0004, plus one Zerowet Splashshield®.




Single-Use Irrijet

55-2104

Same basic assembly as 55-0004, plus includes permanent flexible plastic jet tip and Splatter-Gard™.

55-2134

Same basic assembly as 55-2104, but with 35 cc capacity syringe.



Wound Irrigation

COMBIPORT® + COMBIGUARD II™

LATEX FREE | STERILE | SINGLE PATIENT USE | DISPOSABLE

Ethox



The Combiguard II™ attaches to a syringe (or Combiport) to assist in wound cleansing

Fluid flows through angled port, directed away from cleansed area and toward the outlet ports

Working together, a saline bag uses a Combiport® attached to a Combiguard II™ to effectively cleanse the wound without splash-back on the caregiver



COMBIPORT® Wound Irrigation Device

COMBIGUARD II™ Irrigation Splash Guard

COMBIPORT® WOUND IRRIGATION DEVICE

- Attaches to a standard saline bag to administer a high volume of fluid in a short period of time
- Irrigation helps remove foreign bodies, reduce pathogen numbers and moisten wound tissue
- Patented spike attaches to standard 500 mL or 1000 mL saline bag
- For use only with Combiguard Irrigation Splash Guard (not sold separately)
- Removable cap for easy stop/start during cleansing
- Easy to assemble and use

COMBIGUARD II™ IRRIGATION SPLASH GUARD

- Egg-shaped cup covers wound site to prevent irrigation splash-back, helping to protect clinician
- Dual outlet ports enable contaminated fluids to safely escape
- Fluid delivery inlet delivers a stream equivalent to an 18 g needle
- Angled delivery inlet directs fluid away from cleansed area toward outlet ports
- Rounded edges provide a smooth surface against the skin
- Fits on 30-40 cc luer-tipped syringe
- Can be combined with Combiport Wound Irrigation Device

STERILE

| ITEM# | DESCRIPTION | PK |
|-------|--|----|
| 10222 | Combiguard II™ Irrigation Splash Guard | 50 |
| 10100 | Combiguard II™ + Combiport® | 50 |

NON-STERILE

| ITEM# | DESCRIPTION | PK |
|---------|--|-----|
| 10222NS | Combiguard II™ Irrigation Splash Guard Non-Sterile | 500 |

SMOOTH MOVER



PRODUCT DESCRIPTION

An economical ultra-smooth polyethylene patented board especially designed to furnish expert patient transference with minimal trauma to the patient. The Smooth Mover is capable of safely moving extremely large patients with minimal staffing requirements. It is radiotranslucent; it can be left in place during x-ray procedures. Is lightweight yet durable. Perfect for daily patient moving needs: EMS or hospital, nursing home or clinic applications. Easy on ...easy off. Yes, even one person can move a large patient!

Optional Accessory:

Pull Straps... loops straps slip through hand-holes enabling easy board pick-up on hard surfaces tables (x-ray, etc.) and minimize reach required when transferring patients from bed to stretcher to table during routine moving.

PRODUCT ORDERING INFORMATION

Hospital Shifter Model SM3


Is widely used for patient shifting, bed to stretcher to table, etc..
 Dimensions: 3/16" x 22" x 72"
 Weight: 9 lbs.

Ambulance Patient Shifter Model SM40

Is used as an aid in moving patients from vehicle to stretcher or backboard.
 Dimensions: 3/16" x 18" x 48"
 Weight: 5½ lbs.

ADJUSTABLE PATIENT POSITIONING DEVICES

Patented, inflatable, adjustable OR Comfort™ Patient Positioning Devices provide patient stability and comfort, as well as optimal surgical positioning. Devices can be adjusted throughout a procedure without having to manually lift or move the patient.



SHOULDER-FLOAT™ ADJUSTABLE AXILLARY SUPPORT DEVICE

- Relieves patient's shoulder discomfort during surgical procedures in the lateral decubitus position. Minimizes pressure on shoulder and head to help reduce post surgical pain.
- Useful during total hip replacement, thoracotomy, shoulder arthroscopy and kidney surgery. It can benefit mildly sedated patients undergoing epidural anesthesia and patients undergoing lengthy procedures.

| Shoulder-Float Device | | | |
|-----------------------|-----------------------|--------|-----------------------|
| PART # | PACKAGING DESCRIPTION | PART # | PACKAGING DESCRIPTION |
| SF05 | 5 per BX | SF10 | 10 per BX |



PELVIC-TILT™ ADJUSTABLE DISPLACEMENT DEVICE

- Designed to lift a patient's hip between 20-40 degrees during surgical procedures. Includes a non-inflatable flap to place underneath the patient to prevent a lateral shifting of the device.
- Useful in obstetrics to avoid compression of the vena cava and/or aorta during cesarean sections and deliveries in the supine position. Practical during kidney surgery and surgery on the side of the leg or ankle.

| Pelvic-Tilt Device | | | |
|--------------------|-----------------------|--------|-----------------------|
| PART # | PACKAGING DESCRIPTION | PART # | PACKAGING DESCRIPTION |
| PT05 | 5 per BX | PT10 | 10 per BX |



DELGADO™ INFLATABLE POST

- Slides over the stabilizer post of the operating room table to provide an adjustable "pillow" between the post and the patient.
- Maximizes the stability of the patient and distributes the pressure, minimizing the likelihood of tissue damage or discomfort after surgery.

| Delgado™ Inflatable Post | | | |
|--------------------------|-----------------------|--------|-----------------------|
| PART # | PACKAGING DESCRIPTION | PART # | PACKAGING DESCRIPTION |
| DG05 | 5 per BX | DG10 | 10 per BX |

REQUEST YOUR COPY OF THE COMPLETE IMPERIAL SURGICAL CATALOGUE TODAY!



THE IMPERIAL SURGICAL PORTFOLIO INCLUDES:

- Warming Cabinets
- IV Poles
- Carts & Cabinets
- Tables
- Maternity Bassinets
- Stools & Chairs
- Imperial Ware, Stands & Receptacles
- Specialized Equipment
- Sinks & Sink Systems



Learn more and view the complete digital catalogue at

www.imperialsurgical.com



Ver. 2025-10-14



2165 46th Ave., Lachine, Quebec H8T 2P1
Phone: +1.514.631.7988 | Toll Free: 800.661.5432 | Fax: +1.514.631.9083
Email: info@surgmed.com | Website: www.advancedsurgipharm.com