



COMPLETE PRODUCT CATALOGUE

SERVING CANADIAN HEALTHCARE WORKERS SINCE 1992

ADVANCED SURGI-PHARM INC.

Founded in 1992, we are the distribution arm of Surgmed™ in Canada, representing leading US, international and domestic manufacturers for a wide range of specialty products with unique benefits for applications in the operating room, women's health, labour & delivery, ICU, and sterile processing areas. Equipped with a direct sales team, we provide effective nationwide coverage as the Canadian distributor of Batrik, Vygon, Rocket Medical, Cooper Surgical, SunMed, and other industry leading products including a wide-ranging portfolio of top-rated PPE products.



SURGMED

GROUP 

Proudly Canadian, **The Surgmed™ Group of Companies** began operations in 1992 as a single distribution company known as **Advanced Surgi-Pharm Inc.** The company quickly became known as an industry leader and specialty distributor of value-added surgical products in the Canadian hospital market. Five years later, **Surgmed expanded upon its success with the launch of Batrik Medical Manufacturing Inc.** in 1997 and the subsequent acquisition of **Imperial Surgical Ltd.** one year later in 1998.

Together, our three unique companies serve the Canadian and global healthcare markets with high quality and competitively priced products, all backed by world-class customer care.



Founded in 1992, Advanced Surgi-Pharm Inc. is the distribution arm of Surgmed™ in Canada, representing leading US, international and domestic manufacturers for a wide range of specialty products with unique benefits for applications in the Operating Room, Women's Health, Labor & Delivery, ICU, and Sterile Processing areas. Equipped with a direct sales team, ASP provides effective nationwide coverage as the Canadian distributor of Batrik, Vygon, Rocket Medical, Cooper Surgical, SunMed, and other industry leading products including a wide-ranging portfolio of top-rated PPE products.

Established in 1997 and sold in 85 countries globally, Batrik Medical Manufacturing Inc. specializes in the manufacturing of endoscopy accessories, enzymatic solutions, infection control, central sterilization, and surgical devices. Batrik's in-house controlled environment production facility adheres to strictly enforced quality and regulatory control systems. They are registered as a device manufacturing establishment with the United States FDA and are in full compliance with CE and Health Canada Regulations. All Batrik products are rigorously tested, liability insured, and meet/exceed high international quality standards.

Established in 1935, Imperial Surgical Ltd. is one of Canada's oldest hospital capital equipment manufacturers, well-known for their stainless steel healthcare equipment for Sterile Processing, Operating Room, Examination and Treatment, Patient Handling, and Mortuary/ Pathology areas. Their growing portfolio of products also has applications in Infection Control & Prevention, an important area of focus for their future.

**SURGICAL PRODUCTS
&
OPERATING ROOM**

- 06—Suture Clamp Tags
- 07—Clamp Jaw Covers
- 08—Vascular Silicone Ties
- 09—G.O.L.F.F.™ Anti-Fog Solution
- 10—SKINS™ Surgical Multi-Tip Skin Markers
- 11—PandaPatch™ Eye Safety Tape **NEW!**
- 13—Magnetic Mat
- 14—Pin Guards
- 15—Cautery Tip Cleaner

- 16—Safe-Soak™ SPAC & Floor Mats
- 17—Blue Cloud-Gel™ Anti-Fatigue Mats
- 18—Blue-Cloud Nimbus™ Anti-Fatigue Mats
- 19—Stackable Step Stools
- 20—Tourny™+ Tourniquets
- 22—Surgi-Kit® Instrument Holders
- 23—Lap-Kit® Instrument Holders
- 24—Warm Gel® Instant Heat Packs
- 25—K-Wires, Steinmann Pins & Sterilization Racks

**ENDOSCOPY
PRODUCTS**

- 28—Endoscopy BO2™ Bite Blocks
- 29—Endoscopy Bitestrong™ Bite Blocks
- 30—Couvrette™ Valve Covers
- 31—Daybrush™ Channel Brushes

- 32—ARC™ Channel Brushes
- 35—Valve / Control Brushes
- 36—V-Brushes

CLEANING BRUSHES

- 37—ARC™ Wire Style Brushes
- 41—Double ARC™ Wire Brushes
Asymmetric / Symmetric
- 44—In-Brush™ Antimicrobial Twisted Wire Brushes
- 46—BrightBrush - Color Beaded Twisted Wire Brush **NEW!**
- 47—Stainless Steel ARC™ Micro Brushes
- 48—Instrument Cleaning Flat Brushes

- 49—Precision Brush
- 50—Trifecta™ Brush
- 51—All Purpose Reusable Brushes
- 52—Acetabular Reamer Brush
- 53—Orthopedic Rasps
- 54—Bottle & Drain Clesning Brushes
- 55—Autoclave Cleaning Brush Kits

**INFECTION CONTROL
&
STERILE PROCESSING**

- 56—Silicone Pin Mat
- 57—Tag-It™
- 62—Corner Pockets™ & Soft Pockets™
- 63—Lila's Instrument Staine Remover
- 64—Instrument Protectors
- 67—PrAmLi Power³ Lube™ Instrument Lubricant
- 68—S.T.I.C.S.
- 70—ID Tape
- 71—Tape-B-Gone™ Tape Remover
- 72—Padlocks Security Tags
- 76—Scrub Pad
- 77—Your First ™ Enzymatic & Non-Enzymatic

- Cleaners
- 78—Your First ™ Enviro-Sachet™
 - 79—Your First ™ Foaming Spray
 - 80—Your First ™ F1rst Wipes
 - 81—Your First ™ Dry Sponges
 - 82—Your First ™ RightSwipe™ Enzymatic Sponges
& Bedside Sponge Kits
 - 83—Your First ™ RightSwipe™ Non-Enzymatic
Detergent Bedside Sponge Kits
 - 84—Your First ™ Transport Bags For Endoscope
And Medical Instruments

**PERSONAL
PROTECTIVE
EQUIPMENT**

FACE SHIELDS

- 85—SAVEFACE™ Draped Shields

SHOE COVERS

- 87—SHOOZCOVERS™ SOFT LINE™
- 89—SHOOZCOVERS™ GRIPPEEZ™
- 91—SHOOZCOVERS™ COATED LINE™
- 93—SHOOZCOVERS™ TOUGH LINE™

REUSABLE EQUIPMENT

- 95—Antimicrobial IV Poles With Cuverro Shield™
- 97—Plexiglass Protection Panels
- 98—Stainless Steel Hands-Free Door Opener

NEONATAL

ENTERAL NUTRITION

- 99—Naso-Gastric Feeding Tubes
- 100—Dual-Flow Gastric Feeding Tubes
- 101—Gastro-Duodenal Feeding Tubes
- 102—NutriSafe2® Feeding System
- 106—Nutrifit

VASCULAR ACCESS

- 113—Premicath® PICC 1Fr Catheters
- 115—Nutriline® PICC 2Fr Catheters
- 116—Nutriline Twinflo®
PICC 2Fr Dual Lumen Catheters
- 117—Introducers & Microflash
- 118—Grip-Lok® Catheter Securement Device
- 119—Umbilical Catheters
- 123—Umbilical Placement Pack
- 125—Exchange Transfusion Tray

- 125—PICC Insertion Trays
- 127—Neoventa - Goldtrace®

THERMOREGULATION

- 128—Neohelp™ Heat Loss Prevention Suit
- 130—Warm Gel® Infant Heel Warmer
- 131—TransWarmer® Transport Mattress

RESPIRATORY SUPPORT

- 132—Surfcath™ Surfactant Administration Catheter
- 134—Endotracheal Tubes
- 135—Mucus Extractors
- 136—Nasal Ventilation Tubes
- 137—Safe Spot Resuscitator/Ventilator
- 139—Inca® Infant Nasal CPAP Assembly
- 140—Neo-Fit™ Neonatal Endotracheal Tube Grip & Neo-Prep™ Skin Barrier

**OBSTETRICS/
WOMEN'S HEALTH**

- 141—Os Finder™ Cervical Os and Canal Dilator
- 142—Tampa Infusion Catheter
- 143—H/S Catheters
- 144—H/S Elliptosphere™ Catheter Set and Tray
- 145—Comfort Belly Bands

- 146—Laminaria Cervical Dilator

IV THERAPY

- 147—LeaderCath Arterial Catheters
- 148—Lectrospiral® Coiled Extension Tubes
- 149—Lectrocath Extension Lines
- 150—Biocap Needle-free Drug Aspiration

- 151—IV Accessories
- 152—Infu-Surg® Pressure Infusion Bag
- 153—Infu-Surg® Clarity Pressure Infusion Bag

DRAINAGE

- 154—ASEPT® Catheter and Insertion Kits

- 157—Jet-Vac® Bottles

RESPIRATORY TRACT

- 162—Bronchitrac® "L" Set
- 163—COMFIT™ Endotracheal Tube Holder

- 164—Esophageal Balloon Catheter Set

EMERGENCY

- 165—Tum-E-Vac® Gastric Lavage Kit

- 166—Lavage Tubing

WOUND IRRIGATION

- 167—Irrijet® Splatter Guard

- 169—Combiport® + Combiguard II™

PATIENT HANDLING

- 170—Smooth Mover Patient Shifters



- 171—Adjustable Patient Positioning Devices

SUTURE CLAMP TAGS

Non-toxic, radiopaque, and latex free, our **Suture Clamp Tags** are designed to protect and safely hold sutures without any damage.

Additionally, its foam base, affixed to a sterile area with the clamps, allows for easy insertion.

- Packed in a convenient foam base with easy peel-off adhesive backing, even while wearing gloves
- Available in sterile and non-sterile
- Available in 2, 3, or 5 pairs
- Pre-counted



OUR "EASY PEEL CORNER" IS A FEATURE UNIQUE TO Batrik

PRODUCT SELECTION

ITEM # STERILE	ITEM# NON-STERILE	COLOR	DIMENSIONS	PACKAGING
11-1002	911-1002	Yellow	Standard \varnothing (0.06" x 0.5" / 0.15 cm x 1.27 cm)	2 pairs / foam base
11-1003	911-1003	Yellow	Standard \varnothing (0.06" x 0.5" L / 0.15 cm x 1.27 cm)	3 pairs / foam base
11-1005	911-1005	Yellow	Standard \varnothing (0.06" x 0.5" L / 0.15 cm x 1.27 cm)	5 pairs / foam base
11-1025	911-1025	Assorted	Standard \varnothing (0.06" x 0.5" L / 0.15 cm x 1.27 cm)	5 pairs / foam base
11-1045	911-1045	White	Standard \varnothing (0.06" x 0.5" L / 0.15 cm x 1.27 cm)	5 pairs / foam base
11-1065	911-1065	Blue	Standard \varnothing (0.06" x 0.5" L / 0.15 cm x 1.27 cm)	5 pairs / foam base
11-1085	911-1085	Red	Standard \varnothing (0.06" x 0.5" L / 0.15 cm x 1.27 cm)	5 pairs / foam base
11-3002	911-3002	Pink	Large \varnothing (0.09" x 0.75" L / 0.23 cm x 1.91 cm)	2 pairs / foam base
11-3003	911-3003	Pink	Large \varnothing (0.09" x 0.75" L / 0.23 cm x 1.91 cm)	3 pairs / foam base
11-3005	911-3005	Pink	Large \varnothing (0.09" x 0.75" L / 0.23 cm x 1.91 cm)	5 pairs / foam base

STERILE PACKAGING: 1 Foam Base / Pouch – 5 Pouches / Box

NON-STERILE PACKAGING: Sold per foam base

CLAMP JAW COVERS

Manufactured from medical grade silicone, our **Clamp Jaw Covers** reduce trauma to vessels requiring occlusion, particularly when a hemostatic clamp is utilized. The thread-like indents allow for gentle, constant and non-slip tension.



GOOD TO KNOW OUR CLAMP JAW COVERS CONVERT ALMOST ANY INSTRUMENT TO ATRAUMATIC



- Gamma sterilized
- Available in 3 colors and 3 widths to fit most clamp sizes
- Radiopaque and latex free

PRODUCT SELECTION

ITEM #	STERILE	ITEM #	NON-STERILE	DIMENSIONS	COLOR
12-1000		912-1000		0.1" X 5" / 0.25 cm X 12.7 cm	Red
12-1001		912-1001		0.2" X 5" / 0.51 cm X 12.7 cm	White
12-1002		912-1002		0.3" X 5" / 0.76 cm X 12.7 cm	Blue

STERILE PACKAGING: 2 clamps per sterile pouch - 10 pouches per box

NON-STERILE PACKAGING: 20 clamps per bag (also available 100 units / bag)


VASCULAR SILICONE TIES

Our **Vascular Silicone Ties** permit safe and smooth retraction, occlusion and identification of vessels during surgical procedures.

Unlike clamps and retractors which can pinch the vessels flat, or even crack them, the **BATRIK** elliptical shape allows for the gentle handling of arteries, veins, nerves, tendons, spermatic cords, ureters, or whenever delicate retraction / occlusion is required.

- Manufactured of medical grade silicone
- Do not twist on themselves when unpacked
- Easy to separate, eliminating the need for cardboard inserts
- Sterile, radiopaque, and latex free

AVAILABLE IN 4
SIZES AND 4 COLORS



PRODUCT SELECTION

ITEM #	DIMENSIONS	COLOR	ITEM #	DIMENSIONS	COLOR
10-1000	0.075" X 16" (0.2 cm x 40 cm)	Mini-Red	10-3200	0.075" x 32" (0.2 cm x 81 cm)	Mini-Red
10-1001	0.075" X 16" (0.2 cm x 40 cm)	Mini-Blue	10-3201	0.075" x 32" (0.2 cm x 81 cm)	Mini-Blue
10-1002	0.075" X 16" (0.2 cm x 40 cm)	Mini-Yellow	10-3202	0.075" x 32" (0.2 cm x 81 cm)	Mini-Yellow
10-1003	0.075" X 16" (0.2 cm x 40 cm)	Mini-White	10-3203	0.075" x 32" (0.2 cm x 81 cm)	Mini-White
10-1010	0.120" x 16" (0.3 cm x 40 cm)	Maxi-Red	10-3210	0.120" x 32" (0.3 cm x 81 cm)	Maxi-Red
10-1011	0.120" x 16" (0.3 cm x 40 cm)	Maxi-Blue	10-3211	0.120" x 32" (0.3 cm x 81 cm)	Maxi-Blue
10-1012	0.120" x 16" (0.3 cm x 40 cm)	Maxi-Yellow	10-3212	0.120" x 32" (0.3 cm x 81 cm)	Maxi-Yellow
10-1013	0.120" x 16" (0.3 cm x 40 cm)	Maxi-White	10-3213	0.120" x 32" (0.3 cm x 81 cm)	Maxi-White
10-1101	0.04" x 16" (0.1 cm x 40 cm)	Extra-Mini-Blue			
10-1111	0.200" x 22" (0.5 cm x 56 cm)	Extra-Maxi-Blue			
10-1113	0.200" x 22" (0.5 cm x 56 cm)	Extra-Maxi-White			

PACKAGING: 2 ties per sterile pouch – 10 pouches per box

G.O.L.F.F.™ ANTI-FOG SOLUTION

Our G.O.L.F.F.™ Anti-Fog Solution solves the fogging problem often encountered during Laparoscopic surgical procedures.

- Prevents the occurrence of fog on lenses
- Ensures clearer vision and avoids unnecessary interruptions during surgery
- Set includes an easy-squeeze 6 mL bottle and an adhesive radiopaque sponge for secure placement
- OEM (available)
- Sterile
- Latex free



PRODUCT SELECTION

ITEM #	DESCRIPTION	PACKAGING	CASE DIMENSIONS
G.O.L.F.F.	Anti-Fog Solution & Adhesive Sponge set	20 / Box	20" x 16" x 18"

SKINS™ SURGICAL MULTI-TIP SKIN MARKERS

Our SKINS™ Surgical Multi-Tip Skin Markers are designed with a beveled tip to enable the creation of fine or bold tracing lines by holding the pen at different angles.

- Available with or without 6 labels and a flexible 6" calibrated (in/mm) ruler
- Available in non-sterile for pre-operative use, or sterile for use in the operating room
- Surgical prep resistant
- Made of lasting gentian violet ink that is non-toxic and non-tattooing
- Latex free





PRODUCT SELECTION

ITEM #	TYPE	DESCRIPTION	PACKAGING
S-3000	Sterile	Multi-Tip fine marker with ruler and label set	50 / Box
S-3001	Sterile	Multi-Tip fine marker	50 / Box
S-30005	Sterile	Multi-Tip fine marker with ruler and label set (Volume Packaging)	500 / Box
S-30015	Sterile	Multi-Tip fine marker (Volume Packaging)	500 / Box
S-93000	Non-Sterile	Multi-Tip fine marker with ruler and label set	1000 / Box
S-93001	Non-Sterile	Multi-Tip fine marker	1000 / Box
S-94000	Non-Sterile	Ruler	100 / Bag
S-94001	Non-Sterile	Ruler (Bulk)	20,000 / Box
S-95000	Non-Sterile	Set of 3 labels	100 / Bag

SAVE by ordering our SKINS™ MARKERS in volume packaging. Cost effective and reduces waste.

NEW!



Introducing PandaPatch™: Innovative Adhesive Eye Protection for Safer Patient Care

PandaPatch™ offers an effective alternative to traditional tape for protecting patients' eyes during surgery. These adhesive eye protectors are designed to reduce the risk of corneal injuries and ophthalmic complications for patients under general anesthesia.

With PandaPatch™, you ensure superior patient outcome by preventing avoidable eye damage associated with the use of standard surgical tape.

KEY FEATURES:

- **Ergonomic Design:** Contours comfortably around the eyes.
- **Post-Surgery Friendly:** Particularly effective for children, the elderly, and individuals with sensitive or delicate skin.
- **Latex-Free & Hypoallergenic Adhesive:** Made from medical-grade material, gentle on the skin and designed to minimize irritation.
- **Non-Adhesive Colored Tabs:** Allowing quick and hassle free application and removal, facilitating effortless monitoring of pupil dilation during surgery.



70-04ADU
Adult Sensitive Eye Protector

HERE'S WHY IT MATTERS:

- During prolonged surgery under general anesthesia, **up to 59%¹** of patients may not fully close their eyes, risking corneal abrasions.
- **PandaPatch™** gentle adhesive eye protectors ensure eyelids remain safely and comfortably closed, minimizing the risks typically associated with the use of standard surgical tape.



REFERENCES

1. Anson, Jonathan, M.D. "Perioperative Corneal Abrasions: Etiology, Prevention, and Management." Clinical Updates, Penn State University College of Medicine, 2021.




PANDAPATCH™ RANGE:

Three adhesive eye protectors to accommodate different age groups and skin sensitivities. Providing optimal protection for every patient.

- Adult Extra Sensitive Adhesive (Blue Tabs)
- Adult Adhesive (Red Tabs)
- Pediatric Adhesive (Red Tabs)



70-04ADU
Adult Adhesive Eye Protector




70-04PED
Pediatric Adhesive Eye Protector

PRODUCT SELECTION

ITEM #	DESCRIPTION	DIMENSIONS	PACKAGING
70-02ADU	ADULT EXTRA SENSITIVE ADHESIVE EYE PROTECTORS (BLUE)	2.15" x 1.35"	50 / Box
70-04ADU	ADULT ADHESIVE EYE PROTECTORS (RED)	2.15" x 1.35"	50 / Box
70-04PED	PEDIATRIC ADHESIVE EYE PROTECTORS (RED)	1.85" x 0.93"	50 / Box

MAGNETIC MAT

Our **Magnetic Mat** provides a sterile, neutral field for the safe placement and hands-free transfer of instruments. It is latex free, reusable, and autoclavable (up to 134° C).



PRODUCT SELECTION


ITEM #	DESCRIPTION	DIMENSIONS	PACKAGING
MAG-02-1	Magnetic Mat	12" x 15" / 30 cm x 38 cm	1 Unit
MAG-02-5	Magnetic Mat (Volume Packaging)	12" x 15" / 30 cm x 38 cm	5 Units

*Independent laboratory tested to validate pre-vacuum steam sterilization efficacy for reprocessing of medical devices in a healthcare setting

PIN GUARDS



- Safely cover the sharp ends of your K-Wire and Steinmann Pins
- Color coded by wire size which eliminates the risk of selecting of the wrong size
- A removable sleeve is provided to avoid injury and facilitate application (sterile only)
- Non-sterile model is available with or without plastic sleeves
- Packaging designed to reduce waste
- Latex free



PRODUCT SELECTION

ITEM #STERILE PIN GUARDS + SLEEVE	ITEM # NON-STERILE PIN GUARDS + SLEEVE	ITEM # NON-STERILE PIN GUARDS (NO SLEEVE)	COLOR	DIMENSIONS (WIRE SIZE)
14-1028	914-1028 - SLV	914-1028	Yellow	0.028" / 0.71 mm
14-1035	914-1035 - SLV	914-1035	Light Blue	0.035" / 0.89 mm
14-1045	914-1045 - SLV	914-1045	White	0.045" / 1.14 mm
14-1054	914-1054 - SLV	914-1054	Dark Red	0.054" / 1.40 mm
14-1062	914-1062 - SLV	914-1062	Green	0.062" / 1.57 mm
14-1078	914-1078 - SLV	914-1078	Pink	0.078" / 1.98 mm
14-1093	914-1093 - SLV	914-1093	Grey	0.093" / 2.36 mm
14-1109	914-1109 - SLV	914-1109	Beige	0.109" / 2.77 mm
14-1125	914-1125 - SLV	914-1125	Purple	0.125" / 3.17 mm
14-1141	914-1141- SLV	914-1141	Brown	0.141" / 3.58 mm
14-1156	914-1156 - SLV	914-1156	Dark Blue	0.156" / 3.96 mm
14-1188	914-1188 - SLV	914-1188	Red	0.188" / 4.77 mm

STERILE PACKAGING: 80 units/box (2 units per sterile pouch - 40 pouches per box)

NON-STERILE PACKAGING: 100 units/box

CAUTERY TIP CLEANER

Our **Cautery Tip Cleaner** assists in removing material from cautery tips for both mono-polar and bi-polar electrodes.

- Foam pad with integrated abrasive material
- Adheres to disposable and non-disposable drapes
- Radiopaque
- Available exclusively in Non-Sterile
- Latex free



PRODUCT SELECTION

ITEM #	DESCRIPTION	DIMENSIONS	PACKAGING
60-10000	Cautery Tip Cleaner Non-Sterile	1.8" x 1.8" / 4.5 cm x 4.5 cm	2,200 /Box

SAFE-SOAK™ SPAC

Our Safe-Soak™ SPAC floor suction disk quickly removes fluid waste from the operating room floor through the application of suction.



- Highly durable: it won't be damaged if stepped on
- Easily maneuverable by foot around operating room floor
- Latex free, disposable, and non-sterile
- The bottom opening permits fluid collection through a highly absorbent material that comes into contact with its surface
- The hose is securely affixed on its side to reduce risk of catching equipment
- Its highly absorbent base material allows for quick fluid wicking and removal when suction is applied

BATRIK SAFE-SOAK™ ABSORBENT FLOOR MATS

Our Safe-Soak™ Absorbent Floor Mats are designed to absorb a variety of fluids from the operating room floor quickly and efficiently, ensuring the floor beneath remains slip-free, safe, and dry during surgical procedures.

- Highly absorbent
- Facilitates quick and easy clean up
- Latex free, single-use, disposable, and non-sterile
- Conforms to uneven surfaces
- Skid-resistant backing (blue absorbent mat only)



PRODUCT SELECTION

ITEM #	DESCRIPTION	ABSORBENCY	SIZE	PACKAGING
SS-20-1000	Safe-Soak™ / White extra absorbent mat	4.5 L / mat	28" x 48" (71 cm x 122 cm)	20 / Box
SS-20-2000	Safe-Soak™ / Blue absorbent mat with skid-resistant backer	2.5 L / mat	28" x 40" (71 cm x 102 cm)	20 / Box
SS-20-3000	Safe-Soak™ / Yellow absorbent mat	3 L / mat	30" x 52" (76 cm x 132 cm)	15 / Box
SS-21-1000	Safe-Soak™ SPAC	N/A	8.625" (22 cm diameter)	12 / Box

BLUE CLOUD-GEL™ ANTI-FATIGUE MATS


Our **Blue Cloud-Gel™ Anti-Fatigue Mats** have been scientifically engineered for use in the O.R. using patented gel technology and premium foam to provide the ultimate comfort and support for surgeons, nurses, and technicians particularly when movement is restricted for prolonged durations.

- Proprietary 3/4" thickness **Blue Cloud-Gel™** foam technology will not "bottom-out" like conventional foam mats
- Conforms to and supports feet promoting efficient circulation, reducing physical discomfort and fatigue
- Promotes proper posture while diminishing back, knee, and foot discomfort
- Latex free, water resistant, and nonabsorbent
- Reusable, easy to clean, and Cart Washer compatible (contact us for more details)
- Excellent friction properties when damp or dry surgical booties are worn
- Completely sealed with an RF weld for excellent protection against moisture absorption
- US Standards: CALL 117 and NFPA 260 Compliant
- 1 Year Warranty

ⓘ DID YOU KNOW: STUDIES HAVE DEMONSTRATED THAT ANTI-FATIGUE MATS IMPROVE PERFORMANCE, ENERGY, AND POST-OPERATIVE DISCOMFORT.

PRODUCT SELECTION


ITEM #	DIMENSIONS	PACKAGING
800-2401	18" X 24" / 46 cm x 61 cm	1 / Box
800-2405	18" X 24" / 46 cm x 61 cm	5 / Box
800-3201	20" X 32" / 51 cm x 81 cm	1 / Box
800-3203	20" X 32" / 51 cm x 81 cm	3 / Box
800-4801	20" X 48" / 51 cm x 122 cm	1 / Box
800-4802	20" X 48" / 51 cm x 122 cm	2 / Box
800-7201	30" X 72" / 76 cm x 183 cm (Special Order)	1 / Box



18" X 24"



20" X 32"



20" X 48"

BLUE-CLOUD NIMBUS™ ANTI-FATIGUE MATS


Our **Blue-Cloud Nimbus™ Anti-Fatigue Mats** are designed for general hospital use. Made with proprietary high-density “energy-return” polyurethane materials, they provide unmatched comfort, relief, and ergonomic support.

Our Energy-Return properties will help reduce leg, foot, lower back discomfort, and work-related injuries

- Longer lasting support as 3/4” thickness proprietary foam blend will not bottom out like most other foam based anti- fatigue mats
- Antimicrobial properties blended into the mat provides added protection
- Highly water resistant surface will not absorb water or dirt
- Lightweight and easy to clean
- Textured design on surface and high-traction base ensures greater safety. (US Standard NFSI 101-C)
- Beveled edging reduces tripping hazards and anti-slip backing enhances safety
- Durable design is ideal for commercial applications, yet still ergonomically and aesthetically pleasing for home or office environments
- Flammability/Burn resistant per CFR 16 Parts 1630/1631
- Made with eco-friendly materials, is non-toxic, and latex free
- 5 Year Warranty (non-industrial usage), 4 Year Warranty (industrial usage)

GOOD TO KNOW: IF WITHIN 30 DAYS YOU ARE NOT 100% SATISFIED, CONTACT YOUR DEALER FOR A REFUND

(EXCLUDES TRANSPORT, HANDLING, AND TAXES)




PRODUCT SELECTION

ITEM #	DIMENSIONS	PACKAGING
BLACK COLOR		
809-3201	20" X 32" / 51 cm x 81 cm	1 / Box
809-3204	20" X 32" / 51 cm x 81 cm	4 / Box
809-4801	20" X 48" / 51 cm x 122 cm	1 / Box
809-4804	20" X 48" / 51 cm x 122 cm	4 / Box
809-3601	24" X 36" / 61 cm x 91 cm	1 / Box
809-7201	20" X 72" / 51 cm x 183 cm	1 / Box
809-6001	36" X 60" / 91 cm x 152 cm	1 / Box
TAUPE COLOR		
810-3201	20" X 32" / 51 cm x 81 cm	1 / Box
810-3204	20" X 32" / 51 cm x 81 cm	4 / Box
810-4801	20" X 48" / 51 cm x 122 cm	1 / Box
810-4804	20" X 48" / 51 cm x 122 cm	4 / Box
810-3601	24" X 36" / 61 cm x 91 cm	1 / Box
810-7201	20" X 72" / 51 cm x 183 cm	4 / Box
810-6001	36" X 60" / 91 cm x 152 cm	1 / Box

STACKABLE STEP STOOLS

Made of lightweight cast aluminium and protected with an antimicrobial powder coating, our patented **Stackable Step Stools** have been designed with infection prevention as a prime focus.



PRODUCT FEATURES

- Weight: 8.3 lbs (3.8 kg)
- Dimensions: 17.4" L x 13" W x 4.7" H (44.2 cm L x 33 cm W x 11.9 cm H)
- Transport Cart Available
- Weight Capacity: 300 lb (136 kg)
- Warranty: 7 years (step stool only)

ANTIMICROBIAL COATING:

Prevents the spread of hazardous microorganisms and suppresses the growth of mold, bacteria, yeast, and mildew

ANTI-SLIP SURFACE:

Provides enhanced stability.

RUBBER-FREE CORNERS:

Avoids fluid collection.

WASHABLE CART:

Model OR-36364 (without mat) can be sterilized.

HYGIENIC DESIGN:

No welds or large hollow openings to trap debris or fluids

PRODUCT SELECTION

ITEM #	DESCRIPTION	COLOR	CASE DIMENSIONS
OR-36363	Stackable Step Stool with mat	Black	19" x 14" x 6"
OR-36363-02	Stackable Step Stool with mat	Dark Blue	19" x 14" x 6"
OR-36363-03	Stackable Step Stool with mat	Green	19" x 14" x 6"
OR-36363-07	Stackable Step Stool with mat	Orange	19" x 14" x 6"
OR-36364	Stackable Step Stool without mat	N/A	19" x 14" x 6"

PACKAGING: Sold by Unit

REPLACEMENT PARTS

ITEM #	DESCRIPTION	PACKAGING	CASE DIMENSIONS
OR-36363-1A	Replacement mat for OR-36363	1 Unit	15" x 18" x 5" (Ziplock)
OR-36363-1B	Replacement Foot Pegs	1 Unit	Ziplock Based on Quantity ordered

OPTIONAL ACCESSORIES


ITEM #	DESCRIPTION	PACKAGING	CASE DIMENSIONS
OR-3693	Transport cart	1 Unit	34" x 17" x 7"
809-1705	Reusable Anti-Fatigue Mat with adhesive strips to secure it in place	5 / Box	12" x 12" x 4"
899-1710	Single-Use Anti-Fatigue Mat with adhesive strips to secure it in place	10 / Box	12" x 12" x 8"



OR-3693

TOURNY™ + TOURNIQUETS

Our **Tourniquets** are comfortable and easy to use on all patients, including the elderly, undernourished, children, and newborns.



AVAILABLE COLORS



- Available in different sizes including compact sizes, which are convenient for phlebotomy trays
- Innovative latex free material has a more comfortable texture for patients, puts less stress on skin, and safeguards hospital staff and patients from latex allergens

SINGLE-USE BENEFITS

- Reduced risk of cross-contamination (patient to patient or patient to healthcare provider)
- Reduced risk of spread of nosocomial infections
- May improve your facility's compliance with OSHA's 1910-1030 (d)(1) "Blood-borne Pathogens Standard"



FOR YOUR CONVENIENCE:
25 UNITS PER ROLL
COME IN A DISPENSER BOX
FOR EASY ACCESS



PRODUCT SELECTION


ITEM #	DESCRIPTION	COLOR	DIMENSIONS	PACKAGING
K-002	TOURNY™+ Non-Sterile	Blue	1" x 18" / 2.5 cm x 45 cm	1 / Bag
K-004	TOURNY™+ Non-Sterile	Blue	1" x 18" / 2.5 cm x 45 cm	100 / Box
K-009	TOURNY™+ Non-Sterile	Blue	0.200" x 22" / 0.5 cm x 56 cm	20 / Box
K-010	TOURNY™+ Non-Sterile	White	0.200" x 22" / 0.5 cm x 56 cm	20 / Box
K-400	TOURNY™+ Non-Sterile	Blue	1" x 18" / 2.5 cm x 45 cm	16 rolls 400 / Box
K-400PB	TOURNY™+ Non-Sterile (Pebbled)	Light Blue	1" x 18" / 2.5 cm x 45 cm	16 rolls 400 / Box
K-100	TOURNY™+ Non-Sterile	Assorted (3 rolls of each) Blue, Orange, Green, and Pink	1" x 18" / 2.5 cm x 45 cm	12 rolls 300/ Box
K-254	TOURNY™+ Non-Sterile (Pebbled)	Light Blue	1" x 18" / 2.5 cm x 45 cm	250 / Box
K-2-18RL	TOURNY™+ Non-Sterile	Blue	1" x 18" / 2.5 cm x 45 cm	1 Roll
K-2-18PBRL	TOURNY™+ Non-Sterile (Pebbled)	Light Blue	1" x 18" / 2.5 cm x 45 cm	1 Roll
K-3-18RL	TOURNY™+ Non-Sterile	Green	1" x 18" / 2.5 cm x 45 cm	1 Roll
K-6-18RL	TOURNY™+ Non-Sterile	Pink	1" x 18" / 2.5 cm x 45 cm	1 Roll
K-7-18RL	TOURNY™+ Non-Sterile	Orange	1" x 18" / 2.5 cm x 45 cm	1 Roll
K-2-29RL	Long TOURNY™+ Non-Sterile	Blue	1" x 29.5" / 2.5 cm x 75 cm	1 Roll

SURGI-KIT®

LATEX FREE | STERILE | SINGLE USE | DISPOSABLE

Ethox


- Disposable pockets for sterile delivery to OR
- Keeps laparoscopic instruments close to operating field
- Maintains sterile operating field while keeping instruments within arm's reach
- Reduces time required for locating instruments
- Eliminates cleaning procedures
- Multi-pocket sleeves to accommodate a variety of instrument sizes
- Available with peel-n-stick tape to adhere to operating drape




8 Pocket
SK-100
SK-100N




4 Pocket
SK-150
SK-150N




2 Pocket
SK-200
SK-200N



8 Pocket + Tape
SK-300
SK-300N



4 Pocket + Tape
SK-350
SK-350N



2 Pocket + Tape
SK-400
SK-400N



10 Pocket + Tape
SK-650
SK-650N

STERILE

ITEM#	DESCRIPTION	SIZE	PK
SK-100	8 pockets	19cm x 23cm	50
SK-150	4 pockets	19cm x 11cm	100
SK-200	2 pockets	28cm x 18cm	50
SK-300	8 pockets + tape	19cm x 23cm	50
SK-350	4 pockets + tape	19cm x 11cm	75
SK-400	2 pockets + tape	28cm x 18cm	50
SK-650	10 pockets + tape	5cm x 19cm	50

NON-STERILE

ITEM#	DESCRIPTION	SIZE	PK
SK-100N	8 pockets	19cm x 23cm	50
SK-150N	4 pockets	19cm x 11cm	100
SK-200N	2 pockets	28cm x 18cm	50
SK-300N	8 pockets + tape	19cm x 23cm	50
SK-350N	4 pockets + tape	19cm x 11cm	75
SK-400N	2 pockets + tape	28cm x 18cm	50
SK-650N	10 pockets + tape	5cm x 19cm	50

* Non-sterile kits available only to kit packers.

Instrument Holder


LAP-KIT®

LATEX FREE | STERILE | SINGLE USE | DISPOSABLE


Ethox




- Disposable pockets for sterile delivery to OR
- Keeps laparoscopic instruments close to operating field
- Maintains sterile operating field while keeping instruments within arm's reach
- Reduces time required for locating instruments
- Eliminates cleaning procedures
- Multi-pocket sleeves to accommodate a variety of instrument sizes
- Available with peel-n-stick tape to adhere to operating drape



2 Pockets
LK-700
LK-700N



4 Pockets
LK-800



3 Pockets
LK-900
LK-900N

STERILE

ITEM	DESCRIPTION	SIZE	PK
LK-700	2 Pockets + Tape	6.5 cm x 30.5 cm	25
LK-800	4 Pockets + Tape	3.2 cm x 30.5 cm	25
LK-900	3 Pockets + Tape	3.2 cm x 30.5 cm (2), 6.5 cm x 30.5 cm	25

NON-STERILE

ITEM	DESCRIPTION	SIZE	PK
LK-700N	2 Pockets + Tape	6.5 cm x 30.5 cm	25
LK-900N	3 Pockets + Tape	3.2 cm x 30.5 cm (2), 6.5 cm x 30.5 cm	25

* Non-sterile kits available only to kit packers.

Warm Gel®

Get relief where you need it when you need it



Features and Benefits

- Pliable WarmGel® conforms to any body part
- Easy to activate
- Soft non-woven fabric for patient comfort
- Consistent temperature
- Available in two convenient sizes

WarmGel® is a versatile and reliable source for instant heat

Many applications in:

- Maternity
- IV Therapy
- Labor & Delivery
- Rehabilitation

Temperature Range:

WarmGel® has been designed to achieve 112°F (44°C) when activated at room temperature of 75 °F (24 °C)

Item #	Description
20419	Medium 6" x 6" Adult WarmGel® Instant Heat Pack - 40/Case
20420	Large 7" x 9" Adult WarmGel® Instant Heat Pack - 24/Case

K-Wires & Steinmann Pins

Packaging: 6 or 10 units per plastic tube / 1 tube per package

K-Wires

	DOUBLE TROCAR	SINGLE TROCAR				
	Item #	Item #	Length (in)	Length (cm)	Diameter (in)	Diameter (mm)
SMOOTH	SU78.2000	-	4"	10.1 cm	.031"	0.8 mm
	SU78.2010	-	4"	10.1 cm	.039"	1.0 mm
	SU78.2020	-	4"	10.1 cm	.047"	1.2 mm
	SU78.2040	-	4"	10.1 cm	.054"	1.4 mm
	SU78.2030	-	4"	10.1 cm	.062"	1.6 mm
	SU78.1210	-	6"	15.2 cm	.031"	0.8 mm
	SU78.1220	SU78.2810	6"	15.2 cm	.039"	1.0 mm
	SU78.1230	SU78.2820	6"	15.2 cm	.047"	1.2 mm
	SU78.1240	-	6"	15.2 cm	.054"	1.4 mm
	SU78.1250	SU78.2850	6"	15.2 cm	.062"	1.6 mm
	SU78.2105	SU78.2500	9"	22.8 cm	.031"	0.8 mm
	SU78.2110	SU78.2510	9"	22.8 cm	.039"	1.0 mm
	SU78.2120	SU78.2520	9"	22.8 cm	.047"	1.2 mm
	SU78.2140	-	9"	22.8 cm	.054"	1.4 mm
	SU78.2130	SU78.2530	9"	22.8 cm	.062"	1.6 mm
FULL THREAD	SU78.4030	-	9"	22.8 cm	.062"	1.6 mm
SMOOTH	DOUBLE DIAMOND	SINGLE DIAMOND				
	Item #	Item #	Length (in)	Length (cm)	Diameter (in)	Diameter (mm)
	-	SU78.2810	6"	15.2 cm	.039"	1.0 mm
	-	SU78.2820	6"	15.2 cm	.047"	1.2 mm
	-	SU78.2850	6"	15.2 cm	.062"	1.6 mm
	SU78.3110	SU78.3410	9"	22.8 cm	.039"	1.0 mm
	SU78.3120	SU78.3420	9"	22.8 cm	.047"	1.2 mm
SU78.3130	SU78.3430	9"	22.8 cm	.062"	1.6 mm	
FULL THREAD	-	SU78.4070	9"	22.8 cm	.062"	1.6 mm


K-Wires & Steinmann Pins

Steinmann Pins

DOUBLE TROCAR		SINGLE TROCAR					
Item #	Item #	Length (in)	Length (cm)	Diameter (in)	Diameter (mm)		
-	SU78.6525	6"	15.2 cm	.098"	2.5 mm		
-	SU78.6429	6"	15.2 cm	.110"	2.8 mm		
SMOOTH	SU78.5500	SU78.6100	9"	22.8 cm	.078"	2.0 mm	
	SU78.5530	SU78.6130	9"	22.8 cm	.098"	2.5 mm	
	SU78.5560	SU78.6160	9"	22.8 cm	.110"	2.8 mm	
	SU78.5590	SU78.6190	9"	22.8 cm	.118"	3.0 mm	
	SU78.5620	SU78.6220	9"	22.8 cm	.141"	3.5 mm	
	SU78.5650	SU78.6250	9"	22.8 cm	.157"	4.0 mm	
	SU78.5680	SU78.6280	9"	22.8 cm	.177"	4.5 mm	
	SU78.5690	SU78.6290	9"	22.8 cm	.197"	5.0 mm	
FULL THREAD	SU78.8500	-	9"	22.8 cm	.078"	2.0 mm	
	SU78.8530	-	9"	22.8 cm	.098"	2.5 mm	
	SU78.8560	-	9"	22.8 cm	.110"	2.8 mm	
	SU78.8590	-	9"	22.8 cm	.118"	3.0 mm	
	SU78.8620	-	9"	22.8 cm	.141"	3.5 mm	
	SU78.8650	-	9"	22.8 cm	.157"	4.0 mm	
SMOOTH	DOUBLE DIAMOND	SINGLE DIAMOND					
	Item #	Item #	Length (in)	Length (cm)	Diameter (in)	Diameter (mm)	
	SU78.7000	SU78.7780	9"	22.8 cm	.078"	2.0 mm	
	SU78.7060	SU78.7660	9"	22.8 cm	.110"	2.8 mm	
	-	SU78.7690	9"	22.8 cm	.118"	3.0 mm	
	-	SU78.7720	9"	22.8 cm	.141"	3.5 mm	
	-	SU78.7782	9"	22.8 cm	.177"	4.5 mm	
FULL THREAD	SU78.8300	-	9"	22.8 cm	.078"	2.0 mm	
	-	SU78.8030	9"	22.8 cm	.098"	2.5 mm	
	-	SU78.8120	9"	22.8 cm	.141"	3.5 mm	
3 SHANK END	-	SU78.7630	9"	22.8 cm	.110"	2.8 mm	

Ver.2021-04-22

K-Wires, Steinmann Pins & Sterilization Racks



Sterilization Racks

Sterilization racks are designed for use during washing, sterilizing, storing/transporting, and dispensing K-wires, Steinmann Pins, or Twist Drills. The racks are lightweight but made to withstand repeated use. They fold closed for sterilization and easily open to convert to a table top stand for use in the O.R. Dimensions are clearly marked on the rack to assist with quick identification.



K-Wire and Pin Rack

Number	Details	Dimensions	Quantity / Package
SU98.5409-RO	Anodized aluminium, holds 9" K-wires and pins Holds 6 of each of the following diameters: 0.028" (0.7 mm) 0.045" (1.1 mm) 0.062" (1.6 mm) 0.093" (2.4 mm) 0.125" (3.2 mm) 0.156" (4.0 mm) 0.035" (0.9 mm) 0.054" (1.4 mm) 0.078" (2.0 mm) 0.109" (2.8 mm) 0.140" (3.5 mm) 0.177" (4.5 mm)	12.5" x 5.25" x 1.5" (32 x 13 x 3.8)	1



K-Wire Rack

Number	Details	Dimensions	Quantity / Package
SU98.5404-RO	Anodized aluminium, holds 4 - 6" K-wires and pins Holds 6 of each of the following diameters: 0.028" (0.7 mm) 0.045" (1.1 mm) 0.062" (1.6 mm) 0.035" (0.9 mm) 0.054" (1.4 mm)	12.5" x 5.25" x 1.5" (32 x 13 x 3.8)	1

Ver.2021-07-21

ENDOSCOPY BO2™ BITE BLOCKS


WITH NASAL OXYGEN DELIVERY

Our BO2™ Bite Blocks significantly enhance patient safety by providing supplemental oxygenation during gastrointestinal endoscopy, preventing hypoxia in sedated patients, and ensuring superior patient comfort.

- Saves time and money by eliminating the need for separate nasal prongs and oxygen masks
- Designed to be free of sharp edges to ensure patient comfort and safety
- Available exclusively in Single-Use
- Optional retention strap is available to ensure a secure fit
- The Bite Blocks and optional retention straps are all latex free, lead free and DEHP free

BO2
(FRONT SIDE)

BO2
(BACK SIDE)



PRODUCT SELECTION

ITEM #	TYPE	DESCRIPTION	PACKAGING	CASE DIMENSION
BO2	Single-Use	18mm / 54FR, includes strap	100 / Box	10" x 10" x 10"
BO2-N	Single-Use	18mm / 54FR, without strap	100 / Box	10" x 10" x 10"

ENDOSCOPY BITESTRONG™ BITE BLOCKS

Our Bitestrong™ Bite Blocks provide a safe passage through the patient's mouth during TEE, protecting the patient's teeth and endoscope.


Our molding process ensures a smooth finish with no sharp edges (flashing) that can often irritate a patient's mouth. All Bite Blocks are individually wrapped in sealed, easy to open PE bags and are latex free, lead free, and DEHP free. Additional straps are available, if required

SINGLE USE

Available without dental retention

REUSABLE

- Designed for repeated usage and will withstand commonly practiced chemical and steam sterilization cycles (134°C / 275°F).
- Available with or without an adjustable retention strap



PRODUCT SELECTION

ITEM #	TYPE	DESCRIPTION	PACKAGING	CASE DIMENSIONS
B1	Single-Use Pediatric	16 mm / 48FR, includes strap	100 / Box	10" x 10" x 10"
B1-25	Single-Use Pediatric	16 mm / 48FR, includes strap	25 / Box	10" x 7" x 7"
B1-N	Single-Use Pediatric	16 mm / 48FR, without strap	100 / Box	10" x 10" x 10"
B3	Single-Use	20 mm / 60FR, includes strap	100 / Box	10" x 10" x 10"
B3-N	Single-Use	20 mm / 60FR, without strap	100 / Box	10" x 10" x 10"
B3-SA	Single-Use	20 mm / 60FR, with strap attached to one side	100 / Box	10" x 10" x 10"
B3-NDR	Single-Use	20 mm / 60FR, no dental retention rim	100 / Box	10" x 10" x 10"
B3-NDR-N	Single-Use	Same as B3-N, no dental retention rim, without strap.	100 / Box	10" x 10" x 10"
B4	Single-Use	19 mm / 60FR	100 / Box	10" x 10" x 10"
B4T	Single-Use	19 mm / 60FR, includes strap	100 / Box	10" x 10" x 10"
B4T-N	Single-Use	19 mm / 60FR, without strap	100 / Box	10" x 10" x 10"
B6	Reusable	19 mm / 60FR	100 / Box	10" x 10" x 10"
B6T	Reusable	19 mm / 60FR, includes strap	100 / Box	10" x 10" x 10"
B6T-N	Reusable	19 mm / 60FR, without strap	100 / Box	10" x 10" x 10"
B6-25	Reusable	19 mm / 60FR	25 / Box	10" x 7" x 7"
B6T-25	Reusable	19 mm / 60FR, includes strap	25 / Box	10" x 7" x 7"
B6T-N-25	Reusable	19 mm / 60FR, without strap	25 / Box	10" x 7" x 7"
BS	Single-Use	Bite-Block straps 21.7" x 0.79" / 55 cm x 2 cm	100 / Box	10" x 13" x 1.5"

COUVRETTE™ VALVE COVERS

Our Couvrette™ Valve Covers eliminate the need for manual cleaning and reprocessing of reusable valves. Our single use Valve Covers effectively reduce cross contamination between patients and also reduce fluid leakage from the parts during procedures.

- Easily and securely fit to all valves
- Reduce the potential for errors
- Single-Use
- Non-Sterile
- Latex free

TIP: REUSABLE VALVE COVERS ARE NOTORIOUSLY DIFFICULT TO CLEAN

SAVE TIME AND REDUCE RISK WITH COUVRETTETM



16-0002



16-0022



16-0044



PRODUCT SELECTION

ITEM #	DESCRIPTION	COLOR	PACKAGING
16-0002	Biopsy Valve Cover - Olympus®	Blue	100 / Box
16-0022	Biopsy Valve Cover - Olympus® & Fujinon®	Light Blue	25 / Box
16-0044	Biopsy Valve Cover - Pentax®	Pink	100 / Box


NOTE: Contact us to verify compatibility with your scope

DAYBRUSH™ CHANNEL BRUSHES

REUSABLE

Our Reusable **DAYBRUSH™ Channel Brushes** can be used for cleaning endoscopes for an entire day. In fact, it's been proven to effectively clean up to 10 scopes before the need to discard.*

- Soft tip ensures easy manipulation within the channel with no risk of damage
- The back-end is sealed and softened to prevent fluid from getting trapped inside the tube. The softening process also prevents damage if the cleaning brush is used in reverse.
- Can be steam sterilized.
- Bristle volume is 20 mm in length – nearly double most brushes on the market. A large brush provides more surface contact during the cleaning process.
- Significant per use savings when compared to other reusable or disposable channel cleaning brushes.



20 mm / 0.79"

(*Discard after cleaning 10 scopes or sooner if any brush degradation is noticed. After ten applications, the brush may have the appearance that it can continue to clean effectively, but we discourage prolonged use due to risk of cross contamination. Actual number of uses will vary based on handling and processing.

PRODUCT SELECTION

5 MM CYLINDRICAL BRUSHES – SHEATH Ø 1.75 MM

ITEM #	COLOUR	LENGTH	BRUSH LENGTH	BRUSH (Ø)	SHEATH (Ø)	BRUSH ENDS	PACKAGING
30-18050	White	70" / 180 cm	0.79" / 20 mm	0.20" / 5 mm	0.07" / 1.75 mm	Single	Box of 50
30-18050-T	Yellow	70" / 180 cm	0.79" / 20 mm	0.20" / 5 mm	0.07" / 1.75 mm	Single	Box of 50
30-18050-2T	Yellow	70" / 180 cm	0.79" / 20 mm	0.20" / 5 mm	0.07" / 1.75 mm	Double	Box of 50
30-23050	White	90" / 230 cm	0.79" / 20 mm	0.20" / 5 mm	0.07" / 1.75 mm	Single	Box of 50
30-23050-T	Yellow	90" / 230 cm	0.79" / 20 mm	0.20" / 5 mm	0.07" / 1.75 mm	Single	Box of 50
30-23050-2T	Yellow	90" / 230 cm	0.79" / 20 mm	0.20" / 5 mm	0.07" / 1.75 mm	Double	Box of 50

OPTIONAL KIT ADD-ONS

ADD THE SUFFIX:

"F" to receive a kit with our Flat Sponge. Applicable on certain models only, call for details. (Box of 40)

"K" to receive a kit with the single use Valve / Control Cleaning Brush 5/11 Double Brushes (item # 27-57787). Applicable on all models. (Box of 50)


"FK" to receive both the sponge and valve cleaning brushes. (Box of 40)

ARC™ CHANNEL BRUSHES

SINGLE - USE

BATRIK offers one of the largest selection of Endoscopy Channel Cleaning Brushes in the industry. Our brushes are available in a variety of lengths, materials, and styles to accommodate the widest range of endoscopic devices used in the healthcare industry today.

No matter what your requirements, you'll be able to find the brush that suits your specific needs. If you don't see the brush you require, simply call us. We have over 50 different brushes available and can custom fabricate to your precise specifications.



OPTIONAL KIT ADD-ONS

ADD THE SUFFIX:

"F" to receive a kit with our Flat Sponge. Applicable on certain models only, call for details. (Box of 40)

"K" to receive a kit with the single use Valve / Control Cleaning Brush 5/11 Double Brushes (item # 27-57787). Applicable on all models. (Box of 50)

"FK" to receive both the sponge and valve cleaning brushes. (Box of 40)

SERIES 333 & 311



0.94" / 24 mm

PRODUCT SELECTION

5 MM CYLINDRICAL BRUSHES – SHEATH Ø 1.85 MM - PUNCHY PINK POLYPROPYLENE

ITEM #	LENGTH	BRUSH LENGTH	BRUSH (Ø)	SHEATH (Ø)	BRUSH ENDS	PACKAGING
311-23050	90" / 230 cm	0.94" / 24 mm	0.20" / 5 mm	0.073" / 1.85 mm	Single	Box of 50
311-23050-2	90" / 230 cm	0.94" / 24 mm	0.20" / 5 mm	0.073" / 1.85 mm	Double	Box of 50
311-18050	70" / 180 cm	0.94" / 24 mm	0.20" / 5 mm	0.073" / 1.85 mm	Single	Box of 50
311-10050	39" / 100 cm	0.94" / 24 mm	0.20" / 5 mm	0.073" / 1.85 mm	Single	Box of 50

5/6 MM TAPERED BRUSHES – SHEATH Ø 1.85 MM - PUNCHY PINK POLYPROPYLENE

ITEM #	LENGTH	BRUSH LENGTH	BRUSH (Ø)	SHEATH (Ø)	BRUSH ENDS	PACKAGING
333-23056	90" / 230 cm	0.79" / 20 mm	0.20", 0.24" / 5, 6 mm	0.073" / 1.85 mm	Single	Box of 50
333-23056-2	90" / 230 cm	0.79" / 20 mm	0.20", 0.24" / 5, 6 mm	0.073" / 1.85 mm	Double	Box of 50

SERIES 300 & 322



0.79" / 20 mm

PRODUCT SELECTION

5 MM CYLINDRICAL BRUSHES – SHEATH Ø 1.75 MM - METALLIC GREY DELRIN

ITEM #	LENGTH	BRUSH LENGTH	BRUSH (Ø)	SHEATH (Ø)	BRUSH ENDS	PACKAGING
300-18050	70" / 180 cm	0.79" / 20 mm	0.20" / 5 mm	0.070" / 1.75 mm	Single	Box of 50
300-23050	90" / 230 cm	0.79" / 20 mm	0.20" / 5 mm	0.070" / 1.75 mm	Single	Box of 50
300-23050-2	90" / 230 cm	0.79" / 20 mm	0.20" / 5 mm	0.070" / 1.75 mm	Double	Box of 50

5/6 MM TAPERED BRUSHES – SHEATH Ø 1.75 MM - METALLIC GREY DELRIN

ITEM #	LENGTH	BRUSH LENGTH	BRUSH (Ø)	SHEATH (Ø)	BRUSH ENDS	PACKAGING
322-23056	90" / 230 cm	0.79" / 20 mm	0.20", 0.24" / 5, 6 mm	0.070" / 1.75 mm	Single	Box of 50
322-23056-2	90" / 230 cm	0.79" / 20 mm	0.20", 0.24" / 5, 6 mm	0.070" / 1.75 mm	Single	Box of 50

SERIES 344




PRODUCT SELECTION

3 MM CYLINDRICAL BRUSHES – SHEATH Ø 1.55 MM - METALLIC PURPLE DELRIN

ITEM #	LENGTH	BRUSH LENGTH	BRUSH (Ø)	SHEATH (Ø)	BRUSH ENDS	PACKAGING
344-23030	90" / 230 cm	0.79" / 20 mm	0.12" / 3 mm	0.060" / 1.55 mm	Single	Box of 50
344-18030	70" / 180 cm	0.79" / 20 mm	0.12" / 3 mm	0.060" / 1.55 mm	Single	Box of 50
344-10030	39" / 100 cm	0.79" / 20 mm	0.12" / 3 mm	0.060" / 1.55 mm	Single	Box of 50

SERIES 35



PRODUCT SELECTION

2.5 MM CYLINDRICAL BRUSHES – SHEATH Ø 1.00 MM - YELLOW DELRIN - WITH TIP LIMITED AVAILABILITY

ITEM #	LENGTH	BRUSH LENGTH	BRUSH (Ø)	SHEATH (Ø)	BRUSH ENDS	PACKAGING
35-23025	90" / 230 cm	0.59" / 15 mm	0.098" / 2.5 mm	0.039" / 1.00 mm	Single	Box of 50
35-12025	47.2" / 120 cm	0.59" / 15 mm	0.098" / 2.5 mm	0.039" / 1.00 mm	Single	Box of 50


2 MM CYLINDRICAL BRUSHES – SHEATH Ø 1.00 MM - YELLOW DELRIN - WITHOUT TIP

ITEM #	LENGTH	BRUSH LENGTH	BRUSH (Ø)	SHEATH (Ø)	BRUSH ENDS	PACKAGING
35-23020	90" / 230 cm	0.59" / 15 mm	0.079" / 2.0 mm	0.039" / 1.00 mm	Single	Box of 50
35-12020	47.2" / 120 cm	0.59" / 15 mm	0.079" / 2.0 mm	0.039" / 1.00 mm	Single	Box of 50
35-18020	70.8" / 180 cm	0.59" / 15 mm	0.079" / 2.0 mm	0.039" / 1.00 mm	Single	Box of 50

VALVE / CONTROL BRUSHES

27-05197

27-57787

40-18511

36-10010

27-15311


PRODUCT SELECTION


ITEM #	DESCRIPTION	LENGTH	BRUSH 1 LENGTH	BRUSH 2 LENGTH	BRUSH 1 (Ø)	BRUSH 2 (Ø)
27-05197	Double Trumpet 5/7 White Color	5.00" / 13 cm	0.20" / 5 mm	1.00" / 25 mm	0.20" / 5 mm	0.28" / 7 mm
PACKAGING: Bag of 3 (single-use recommended). Volume Packaging: Bag of 50 use ITEM # 27-05197-50						
27-57787	Double, Universal 5/11 Yellow Color	6.00" / 15 cm	1.18" / 30 mm	0.79" / 20 mm	0.43" / 11 mm	0.20" / 5 mm
PACKAGING: Bag of 3 (single-use recommended). Volume Packaging: Bag of 50 use ITEM # 27-57787-50						
40-18511	Double Universal 5/11 Extra Long Blue Color	8.07" / 20.5 cm	1.57" / 40 mm	0.79" / 20 mm	0.43" / 11 mm	0.20" / 5 mm
PACKAGING: Bag of 3 (reusable and can be steam sterilized). Volume Packaging: Bag of 50 use ITEM # 40-18511-50						
36-10010	Single end large Red Color	4.3" / 11 cm	1.57" / 40 mm	N/A	0.39" / 10 mm	N/A
PACKAGING: Bag of 3 (single-use recommended). Volume Packaging: Bag of 50 use ITEM # 36-10010-50						
27-15311	Single end large Blue Color	4" / 10.3 cm	1.57" / 40 mm	N/A	0.39" / 10 mm	N/A

PACKAGING: Bag of 3 (single-use recommended). Volume Packaging: Bag of 50 use ITEM # 27-15311-50

V-BRUSHES

EXTERNAL CLEANING BRUSHES

Our External Cleaning Brushes for endoscopes / optical instruments provide excellent cleaning capabilities. With their unique U-shape bristles and innovative design, our brushes are the ultimate tool for cleaning the exterior of your endoscopic instruments. Free from sharp or metallic ends that can scratch or cause damage to valuable equipment, both single-use and reusable models are latex free and feature ergonomic handles.



PRODUCT SELECTION

ITEM #	DESCRIPTION	BRUSH LENGTH	HANDLE LENGTH / (Ø)	TOTAL LENGTH	PACKAGING
36-20005	Single-Use External Cleaning Brush	70 mm / 25 mm	80 mm / 3 mm	150 mm	50 / Box
36-20010	Reusable External Cleaning Brush	90 mm / 25 mm	120 mm / 12 mm	210 mm	3 / Bag

ARC™ WIRE STYLE BRUSHES

SINGLE-USE AND REUSABLE

You've got options with our ARC™ Wire Brushes. Our single-use and reusable brushes are available in both flexible and rigid varieties and come in a multitude of sizes. ARC™ Wire Brushes are guaranteed to last longer than twisted wire and preserve the integrity of your valuable instruments.

BRUSH SELECTION TIPS:

INSTRUMENT SPECIFICATIONS

What is the maximum length of the instrument(s) to be cleaned?

Add at least 2" (5 cm) to that length, so you have an adequate length to work with.

What is the inner diameter of the instrument?

Add at most 0.125" (3 mm) to determine the diameter of the bristle to ensure effective cleaning. The bristle diameter has to slightly exceed the inner tube diameter. Bristles that are too large will bend too much and not clean effectively.

BRUSH SELECTION

Simply find the section in our catalogue that best corresponds to the overall length and diameter of your instrument. As a rule, if your present brush is 8" (20.32 cm) in length and we only offer 10" (25.4 cm) length with your required diameter, that is entirely acceptable.

In many cases, our ARC™ Wire Brushes will have more bristles on the brush. As an example, many 5 mm (0.19") diameter brushes are only 1" (2.5 cm) in length. At Batrik, we strived to provide the greatest amount of bristles to ensure the most effective cleaning.

Too small

Bristles will not make sufficient contact and surfaces may go untouched.



Too large

Applies unnecessary pressure on bristles and may cause brush to get stuck in device.



Just right


Bristles make sufficient contact with the instrument.



RIGID



FLEXIBLE



CUSTOMIZABLE WIRE BRUSHES

A unique aspect of single-use ARC™ Wire Brushes is that they can be customized to any length:

- Up to 2000 mm for Wire Brushes on plastic semi-rigid sheaths with brushes from 5 to 20 mm in diameter.
- Flexible plastic sheaths with brushes from 2-5 mm in diameter.

12" / 30 CM SINGLE-END WIRE BRUSH

	ITEM # SINGLE-USE	ITEM # REUSABLE	COLOR	TOTAL LENGTH	BRUSH (Ø)	BRUSH LENGTH	SHEATH (Ø)
FLEXIBLE BRUSH & SHEATH	41-31502-S	41-31502	White	12.401" / 31.5 cm	0.078" / 2 mm	0.59" / 15 mm	0.39" / 1 mm
	41-32503-S	41-32503	Yellow	12.795" / 32.5 cm	0.118" / 3 mm	0.98" / 25 mm	0.39" / 1 mm
	411-32504-S	411-32504	Violet	12.795" / 32.5 cm	0.157" / 4 mm	0.98" / 25 mm	0.61" / 1.55 mm
	411-32505-S	411-32505	Grey	12.795" / 32.5 cm	0.19" / 5 mm	0.98" / 25 mm	0.68" / 1.75 mm
RIGID BRUSH & SHEATH	41-35005-S	41-35005	White	13.779" / 35 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
	41-35006-S	41-35006	Green	13.779" / 35 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
	41-35008-S	41-35008	Blue	13.779" / 35 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
	41-35010-S	41-35010	Red	13.779" / 35 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
	41-35012-S	41-35012	Yellow	13.779" / 35 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
	41-35015-S	41-35015	Green	13.779" / 35 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
	41-35020-S	41-35020	Blue	13.779" / 35 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Single-Use: Bag of 20 / Reusable: Bag of 3

18" / 45 CM SINGLE-END WIRE BRUSH

	ITEM # SINGLE-USE	ITEM # REUSABLE	COLOR	TOTAL LENGTH	BRUSH (Ø)	BRUSH LENGTH	SHEATH (Ø)
FLEXIBLE BRUSH & SHEATH	41-46502-S	41-46502	White	18.307" / 46.5 cm	0.078" / 2 mm	0.59" / 15 mm	0.39" / 1 mm
	41-47503-S	41-47503	Yellow	18.7" / 47.5 cm	0.118" / 3 mm	0.98" / 25 mm	0.39" / 1 mm
	411-47504-S	411-47504	Violet	18.7" / 47.5 cm	0.157" / 4 mm	0.98" / 25 mm	0.61" / 1.55 mm
	411-47505-S	411-47505	Grey	18.7" / 47.5 cm	0.19" / 5 mm	0.98" / 25 mm	0.68" / 1.75 mm
RIGID BRUSH & SHEATH	41-50005-S	41-50005	White	19.685" / 50 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
	41-50006-S	41-50006	Green	19.685" / 50 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
	41-50008-S	41-50008	Blue	19.685" / 50 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
	41-50010-S	41-50010	Red	19.685" / 50 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
	41-50012-S	41-50012	Yellow	19.685" / 50 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
	41-50015-S	41-50015	Green	19.685" / 50 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
	41-50020-S	41-50020	Blue	19.685" / 50 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Single-Use: Bag of 20 / Reusable: Bag of 3

24" / 60 CM SINGLE-END WIRE BRUSH


	ITEM # SINGLE-USE	ITEM # REUSABLE	COLOR	TOTAL LENGTH	BRUSH (Ø)	BRUSH LENGTH	SHEATH (Ø)
FLEXIBLE BRUSH & SHEATH	41-61502-S	41-61502	White	24.212" / 61.5 cm	0.078" / 2 mm	0.59" / 15 mm	0.39" / 1 mm
	41-62503-S	41-62503	Yellow	24.6" / 62.5 cm	0.118" / 3 mm	0.98" / 25 mm	0.39" / 1 mm
	411-62504-S	411-62504	Violet	24.6" / 62.5 cm	0.157" / 4 mm	0.98" / 25 mm	0.61" / 1.55 mm
	411-62505-S	411-62505	Grey	24.6" / 62.5 cm	0.19" / 5 mm	0.98" / 25 mm	0.68" / 1.75 mm
RIGID BRUSH & SHEATH	41-65005-S	41-65005	White	25.59" / 65 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
	41-65006-S	41-65006	Green	25.59" / 65 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
	41-65008-S	41-65008	Blue	25.59" / 65 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
	41-65010-S	41-65010	Red	25.59" / 65 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
	41-65012-S	41-65012	Yellow	25.59" / 65 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
	41-65015-S	41-65015	Green	25.59" / 65 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
	41-65020-S	41-65020	Blue	25.59" / 65 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Single-Use: Bag of 20 / Reusable: Bag of 3

30" / 75 CM SINGLE-END WIRE BRUSH

	ITEM # SINGLE-USE	ITEM # REUSABLE	COLOR	TOTAL LENGTH	BRUSH (Ø)	BRUSH LENGTH	SHEATH (Ø)
FLEXIBLE BRUSH & SHEATH	41-76502-S	41-76502	White	30.118" / 76.5 cm	0.078" / 2 mm	0.59" / 15 mm	0.39" / 1 mm
	41-77503-S	41-77503	Yellow	30.511" / 77.5 cm	0.118" / 3 mm	0.98" / 25 mm	0.39" / 1 mm
	411-77504-S	411-77504	Violet	30.511" / 77.5 cm	0.157" / 4 mm	0.98" / 25 mm	0.61" / 1.55 mm
	411-77505-S	411-77505	Grey	30.511" / 77.5 cm	0.19" / 5 mm	0.98" / 25 mm	0.68" / 1.75 mm
RIGID BRUSH & SHEATH	41-80005-S	41-80005	White	31.496" / 80 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
	41-80006-S	41-80006	Green	31.496" / 80 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
	41-80008-S	41-80008	Blue	31.496" / 80 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
	41-80010-S	41-80010	Red	31.496" / 80 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
	41-80012-S	41-80012	Yellow	31.496" / 80 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
	41-80015-S	41-80015	Green	31.496" / 80 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
	41-80020-S	41-80020	Blue	31.496" / 80 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Single-Use: Bag of 20 / Reusable: Bag of 3



FLEXIBLE 2 MM

FLEXIBLE 3 MM


FLEXIBLE 4 MM

FLEXIBLE 5 MM



FLEXIBLE

Flexible brush diameter (Ø) ranging from 2 mm to 5 mm



RIGID 20 MM

RIGID 15 MM

RIGID 12 MM

RIGID 10 MM

RIGID 8 MM

RIGID 6 MM

RIGID 5 MM

RIGID

Rigid brush diameter (Ø) ranging from 5 mm to 20 mm

BRUSH SELECTION SYSTEM (Actual Size)

2 mm 0.078"	3 mm 0.118"	4 mm 0.157"	5 mm 0.196"	6 mm 0.236"	8 mm 0.314"	10 mm 0.394"	12 mm 0.472"	15 mm 0.591"	20 mm 0.787"

REUSABLE DOUBLE ARC™ WIRE BRUSHES

MIXED SIZE BRUSHES

14" / 35 CM MIXED SIZE DOUBLE-END WIRE BRUSHES

ITEM #	COLOR	TOTAL LENGTH	BRUSH 1 (Ø)	BRUSH 2 (Ø)	BRUSH LENGTH	SHEATH (Ø)
40-35056	White	13.779" / 35 cm	0.19" / 5 mm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
40-35810	Blue	13.779" / 35 cm	0.31" / 8 mm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
40-35112	Yellow	13.779" / 35 cm	0.39" / 10 mm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
40-35152	Green	13.779" / 35 cm	0.59" / 15 mm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Bag of 3

20" / 50 CM MIXED SIZE DOUBLE-END WIRE BRUSHES

ITEM #	COLOR	TOTAL LENGTH	BRUSH 1 (Ø)	BRUSH 2 (Ø)	BRUSH LENGTH	SHEATH (Ø)
40-50056	White	19.685" / 50 cm	0.19" / 5 mm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
40-50810	Blue	19.685" / 50 cm	0.31" / 8 mm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
40-50112	Yellow	19.685" / 50 cm	0.39" / 10 mm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
40-50152	Green	19.685" / 50 cm	0.59" / 15 mm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Bag of 3

26" / 65 CM MIXED SIZE DOUBLE-END WIRE BRUSHES

ITEM #	COLOR	TOTAL LENGTH	BRUSH 1 (Ø)	BRUSH 2 (Ø)	BRUSH LENGTH	SHEATH (Ø)
40-65056	White	25.59" / 65 cm	0.19" / 5 mm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
40-65810	Blue	25.59" / 65 cm	0.31" / 8 mm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
40-65112	Yellow	25.59" / 65 cm	0.39" / 10 mm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
40-65152	Green	25.59" / 65 cm	0.59" / 15 mm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Bag of 3

31" / 80 CM MIXED SIZE DOUBLE-END WIRE BRUSHES

ITEM #	COLOR	TOTAL LENGTH	BRUSH 1 (Ø)	BRUSH 2 (Ø)	BRUSH LENGTH	SHEATH (Ø)
40-80056	White	31.496" / 80 cm	0.19" / 5 mm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
40-80810	Blue	31.496" / 80 cm	0.31" / 8 mm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
40-80112	Yellow	31.496" / 80 cm	0.39" / 10 mm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
40-80152	Green	31.496" / 80 cm	0.59" / 15 mm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Bag of 3

REUSABLE DOUBLE ARC™ WIRE BRUSHES

SAME SIZE BRUSHES

14" / 35 CM SAME SIZE DOUBLE-END WIRE BRUSHES

ITEM #	COLOR	TOTAL LENGTH	BRUSH 1 (Ø)	BRUSH LENGTH	SHEATH (Ø)
40-35005	White	13.779" / 35 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
40-35006	Green	13.779" / 35 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
40-35008	Blue	13.779" / 35 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
40-35010	Red	13.779" / 35 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
40-35012	Yellow	13.779" / 35 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
40-35015	Green	13.779" / 35 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
40-35020	Blue	13.779" / 35 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Bag of 3

20" / 50 CM SAME SIZE DOUBLE-END WIRE BRUSHES

ITEM #	COLOR	TOTAL LENGTH	BRUSH 1 (Ø)	BRUSH LENGTH	SHEATH (Ø)
40-50005	White	19.685" / 50 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
40-50006	Green	19.685" / 50 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
40-50008	Blue	19.685" / 50 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
40-50010	Red	19.685" / 50 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
40-50012	Yellow	19.685" / 50 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
40-50015	Green	19.685" / 50 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
40-50020	Blue	19.685" / 50 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Bag of 3

26" / 65 CM SAME SIZE DOUBLE-END WIRE BRUSHES

ITEM #	COLOR	TOTAL LENGTH	BRUSH 1 (Ø)	BRUSH LENGTH	SHEATH (Ø)
40-65005	White	25.59" / 65 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
40-65006	Green	25.59" / 65 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
40-65008	Blue	25.59" / 65 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
40-65010	Red	25.59" / 65 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
40-65012	Yellow	25.59" / 65 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
40-65015	Green	25.59" / 65 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
40-65020	Blue	25.59" / 65 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Bag of 3

31" / 80 CM SAME SIZE DOUBLE-END WIRE BRUSHES

ITEM #	COLOR	TOTAL LENGTH	BRUSH 1 (Ø)	BRUSH LENGTH	SHEATH (Ø)
40-80005	White	31.496" / 80 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
40-80006	Green	31.496" / 80 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
40-80008	Blue	31.496" / 80 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
40-80010	Red	31.496" / 80 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
40-65112	Yellow	31.496" / 80 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
40-80015	Green	31.496" / 80 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
40-80020	Blue	31.496" / 80 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Bag of 3

MIXED SIZE BRUSHES

5 MM

6 MM



8 MM

10 MM



10 MM

12.5 MM



15 MM

20 MM




IN-BRUSH™ ANTIMICROBIAL TWISTED WIRE BRUSHES

Our In-Brush™ Antimicrobial Reusable Twisted Wire Brushes are designed for use with Frazier/Baron suction tubes. The twisted wires have a fine/medium diameter yet are rigid enough to ensure excellent cleaning capabilities.

Our wide assortment of nylon twisted wire brushes are ideal for large bore, uniquely shaped, and endoscopic instruments.

GOOD TO KNOW: OUR TWISTED WIRE BRUSHES FEATURE A PROTECTIVE BLUE TIP FOR ADDED SAFETY.



10" / 250 MM FRAZIER /BARON SUCTION TUBE BRUSHES

ITEM #	SIZE	BRUSH LENGTH	BRUSH (Ø)
21-10039	3 FR	0.5" / 12.7 mm	0.039" / 1 mm
21-10067	5 FR	0.5" / 12.7 mm	0.065" / 1.67 mm
20-10059	5 FR	1" / 25.4 mm	0.059" / 1.50 mm
20-10079	6 FR	1" / 25.4 mm	0.079" / 2 mm
20-10090	7 FR	1" / 25.4 mm	0.090" / 2.3 mm
20-10106	8 FR	1" / 25.4 mm	0.106" / 2.7 mm
20-10118	9 FR	1" / 25.4 mm	0.118" / 3 mm
20-10129	10 FR	1.5" / 38 mm	0.129" / 3.3 mm
20-10146	11 FR	1.5" / 38 mm	0.146" / 3.7 mm
20-10157	12 FR	1.5" / 38 mm	0.157" / 4 mm

PACKAGING: Bag of 3.
Also available in bags of 50, add the suffix "-50".

12" / 305 MM NYLON TWISTED WIRE BRUSHES

ITEM #	BRUSH LENGTH	BRUSH (Ø)
22-12078	2.28" / 58 mm	0.078" / 2 mm
22-12118	2.28" / 58 mm	0.118" / 3 mm
22-12157	3" / 76 mm	0.157" / 4 mm
22-12197	3" / 76 mm	0.197" / 5 mm
22-12236	3" / 76 mm	0.236" / 6 mm
22-12393	3" / 76 mm	0.394" / 10 mm

PACKAGING: Bag of 3.
Also available in bags of 50, add the suffix "- 50".

16" / 406 MM NYLON TWISTED WIRE BRUSHES

ITEM #	BRUSH LENGTH	BRUSH (Ø)
23-16197	2.5" / 63.5 mm	0.197" / 5 mm
23-16236	2.5" / 63.5 mm	0.236" / 6 mm
23-16276	3" / 76 mm	0.276" / 7 mm
23-16394	3" / 76 mm	0.394" / 10 mm
23-16472	3" / 76 mm	0.472" / 12 mm
23-16591	3" / 76 mm	0.591" / 15 mm
23-16787	3" / 76 mm	0.787" / 20 mm

PACKAGING: Bag of 3.
Also available in bags of 50, add the suffix "-50".

18" / 457 MM NYLON TWISTED WIRE BRUSHES

ITEM #	BRUSH LENGTH	BRUSH (Ø)
19-18010	0.47" / 12 mm	0.039" / 1 mm
19-18015	0.5" / 12 mm	0.06" / 1.6 mm

PACKAGING: Bag of 3.
Also available in bags of 50, add the suffix "-50".

24" / 610 MM NYLON TWISTED WIRE BRUSHES

ITEM #	BRUSH LENGTH	BRUSH (Ø)
24-24098	2.5" / 63.5 mm	0.098" / 2.5 mm
24-24197	2.5" / 63.5 mm	0.197" / 5 mm
24-24276	3.25" / 82.6 mm	0.276" / 7 mm
24-24394	3.25" / 82.6 mm	0.394" / 10 mm
24-24472	3.25" / 82.6 mm	0.472" / 12 mm

PACKAGING: Bag of 3.
Also available in bags of 50, add the suffix "-50".

TRACHEOSTOMY BRUSHES

ITEM #	DESCRIPTION	TOTAL LENGTH	BRUSH LENGTH	BRUSH (Ø)
25-08500	Large 6-12.3 3 / 4"	8" / 203 mm	3.75" / 95 mm	0.5" / 12.7 mm
25-37550	Small 000-5.1 7 / 8"	4" / 93 mm	1.88" / 48 mm	0.5" / 12.7 mm

PACKAGING: Bag of 3. Also available in bags of 50, add the suffix "-50".



BRONCHO BRUSHES


ITEM #	DESCRIPTION	BRUSH LENGTH	BRUSH (Ø) INCHES / MM
26-06175	Small	6" / 152 mm	1.75" / 44.5 mm
26-08325	Large	8" / 203 mm	3.25" / 82.6 mm

PACKAGING: Bag of 3.
Also available in bags of 50, add the suffix "-50".

BRIGHTBRUSH^{NEW!}

THE RIGHT BRUSH FOR EVERY JOB

Tackle even the toughest cleaning jobs with ease! Our Color-Coded Twisted Wire Brushes are designed for precision, convenience, and superior performance—ensuring spotless results every time.



Color-Coded Bead

Each brush features a color-coded bead that indicates the bristle diameter and size, making it easy to select the right tool quickly and efficiently.



Hanging Loop

The brush includes a loop at the end for convenient hanging.

PRODUCT SELECTION

ITEM #	WIRE (Ø) (mm)	BEAD (Ø) (mm)	BEAD COLOR	DIMENSIONS Length x Brush (Ø) x Brush Length
12" COLOR-CODED TWISTED WIRE BRUSHES				
22-12078B	0.7 mm	4 mm	Red	12" x .078" (2M)x 58 mm
22-12118B	0.7 mm	4 mm	White	12" x .118" (3M)x 58 mm
22-12157B	0.8 mm	4 mm	Blue	12" x .157" (4M)x 76 mm
22-12197B	0.8 mm	4 mm	Green	12" x .197" (5M)x 76 mm
22-12236B	0.8 mm	4 mm	Black	12" x .236" (6M)x 76 mm
16" COLOR-CODED TWISTED WIRE BRUSHES				
23-16197B	1.0 mm	4 mm	Green	16" X .197" (5M)x 63.5 mm
23-16276B	1.0 mm	4 mm	Yellow	16" X .276" (7M)x 76 mm
23-16394B	1.0 mm	4 mm	Pink	16" X .394" (10M)x 76 mm
18" COLOR-CODED TWISTED WIRE BRUSHES				
19-18078B	0.6 mm	4 mm	Red	18" x .078" (2M)x 63.5 mm
19-18118B	0.8 mm	4 mm	White	18" x .118" (3M)x 63.5 mm
19-18157B	0.8 mm	4 mm	Blue	18" x .157" (4M)x 63.5 mm
19-18197B	1.0 mm	4 mm	Green	18" x .197" (5M)x 63.5 mm
19-18236B	1.0 mm	4 mm	Black	18" x .236" (6M)x 76 mm
24" COLOR-CODED TWISTED WIRE BRUSHES				
24-24098B	0.8 mm	4 mm	Orange	24" X .098" (2.5M)x 63.5 mm
24-24197B	0.8 mm	4 mm	Green	24" X .197" (5M)x 63.5 mm
24-24276B	1.0 mm	4 mm	Yellow	24" X .276" (7M)x 82.6 mm
24-24394B	1.0 mm	4 mm	Pink	24" X .394" (10M)x 82.6 mm

PACKAGING: Bag of 3. Also available in bags of 50, add the suffix "-50".

STAINLESS STEEL ARC™ MICRO BRUSHES

Our **Stainless Steel ARC™** Micro Brushes are some of the smallest brushes in the industry. They are designed to clean channels with an inner diameter ranging from 0.5 mm to 1.40 mm.

- Flexible stainless steel sheath will not bend like the regular metal twisted wire brushes
- Polyamide bristles prevents or slows bacterial growth between sterilizations
- Reusable*
- Available in two lengths

TIP: USING THE PULL-THROUGH TECHNIQUE IS HIGHLY RECOMMENDED WITH THESE BRUSHES. SIMPLY INSERT THE BLUNT, NON-BRUSH END FIRST.



PRODUCT SELECTION

ITEM #	TOTAL LENGTH	BRUSH SIZE (Ø)	BRUSH LENGTH	SHEATH (Ø)	PACKAGING
35-03015	12" / 30 cm	0.06" / 1.5 mm	0.6" / 15 mm	0.02" / 0.40 mm	6 / Tube
35-04515	18" / 45 cm	0.06" / 1.5 mm	0.6" / 15 mm	0.02" / 0.40 mm	6 / Tube

*Although all materials used in the fabrication of these brushes can be re-processed repeatedly, due to the delicate nature of these brushes re-use is conditional on their handling and application. Re-use it at the user's discretion.

INSTRUMENT CLEANING FLAT BRUSHES

Our Instrument Cleaning Flat Brushes are specifically designed for hospital use in the pre-clean and preparation of instruments.

- Long ergonomic handle permits ease-of-use
- Extra long length of brush head (60 [±5] mm in Single-Use and 75 [±5] mm in Reusable)
- Transparent brush fibers allow for quick verification of its cleanliness
- The plastic handle will not release any toxic components
- Available in 6 models (4 reusable/autoclavable and 2 single-use)



FLEXIBLE

This brush is designed for large instruments that require a more gentle cleaning application such as optics, plastics, or hard to reach spot.



50-40919



50-40921



EXTRA FLEXIBLE

The very long, flexible bristles eliminate any stress during sanitation, and makes this the ideal brush for fragile instruments as well as both pliable and rigid endoscopes. The length of the bristles assist in cleaning areas that would otherwise be difficult to reach, such as the accessory channel or tip deflection of pliable endoscopes.



50-40920



RIGID

The rigidity of this brush, along with the length of the bristles and handle, improves and simplifies instrument cleaning.



50-40918



50-40922



EXTRA RIGID

Instrument safe polyamide fibers offer steel-like rigidity and allows deep cleaning of soiled instruments, particularly in grooved or uneven surfaces.



50-40917

PRODUCT SELECTION

ITEM #	BRUSH TYPE	DESCRIPTION	COLOR	BRUSH LENGTH	BRUSH SIZE	PACKAGING
50-40917	Extra Rigid	Reusable Autoclavable Flat Brush	Blue	220 mm	75 mm x 20 mm x 10 mm	Bag of 3
50-40918	Rigid	Reusable Autoclavable Flat Brush	Green	220 mm	75 mm x 20 mm x 15 mm	Bag of 3
50-40919	Flexible	Reusable Autoclavable Flat Brush	Yellow	220 mm	75 mm x 20 mm x 15 mm	Bag of 3
50-40920	Extra Flexible	Reusable Autoclavable Flat Brush	Red	220 mm	75 mm x 20 mm x 25 mm	Bag of 3
50-40921	Flexible	Single-Use Flat Brush	Orange	220 mm	60 mm x 20 mm x 15 mm	Bag of 20
50-40922	Rigid	Single-Use Flat Brush	Light Blue	220 mm	60 mm x 20 mm x 15 mm	Bag of 20


Reusable models are also available in bags of 50. Add the suffix "- 50".

REUSABLE PRECISION BRUSH

INSTRUMENT CLEANING

Our Reusable Precision Instrument Cleaning Brush is the ideal brush for the quick removal of difficult debris. It cleans externally as well as deep into instrument jaws. Although the firm and durable bristles feel like stainless steel, there is absolutely no risk of instrument damage or personal injury.

- Bristles have an extra firm density and superior durability
- Polyamide (TYNEX®) bristles
- Steamable thermoplastic ergonomic handle
- Steam / chemical sterilizable
- Latex free



GOOD TO KNOW: ELIMINATES THE NEED FOR THE STAINLESS STEEL TOOTHBRUSH-STYLE BRUSH, WHICH CAN DAMAGE VALUABLE INSTRUMENTS!

PRODUCT SELECTION

ITEM #	DESCRIPTION	TOTAL LENGTH
50-40800	Reusable Precision External Cleaning Brush	6.29" / 160 mm

PACKAGING: Bag of 3

TRIFECTA™ BRUSH

INSTRUMENT CLEANING

Introducing the **Trifecta™** Brush, a versatile 3-in-1 tool designed to meet all your instrument cleaning needs. This reusable instrument brush is available in two styles: one with all-nylon bristles for gentle cleaning and another with nylon bristles on two sides and stainless steel bristles on the third, ideal for tackling tougher buildup and stubborn cleaning tasks. This innovative design ensures effective, thorough maintenance for a wide range of instruments.

TRIFECTA™ BRUSH - NYLON

1 NYLON BRISTLES


Gentle, transparent nylon bristles provide an effective yet safe cleaning solution. This side allows you to verify cleanliness easily, without risking damage to delicate surfaces. Ideal for a light, non-invasive clean.

2 PRECISION TIP

Reach those hard-to-clean spots with the precision tip. This small, firm-bristle tip is designed for tight, hard-to-reach areas, such as instrument jaws and box joints. It ensures a deep clean where it's most needed, without the risk of damage or injury to delicate parts.

3 FIRM NYLON BRISTLES

For tougher residues, the firm nylon side delivers the cleaning power you need. With bristles that feel similar to stainless steel but are safe for your instruments, this side easily handles stubborn dirt and buildup.



28-23012

TRIFECTA™ BRUSH - NYLON & STAINLESS STEEL



29-23012

→ STAINLESS STEEL BRISTLES ARE MORE DURABLE AND EFFECTIVE FOR CLEANING STUBBORN BUILDUP.

PRODUCT SELECTION

ITEM #	COLOR	DESCRIPTION	LENGTH TOTAL	BRISTLE LENGTH	BRISTLES 1 Ø	BRISTLES 2 Ø	BRISTLES 3 Ø
28-23012	Blue	Triple Nylon Brush	8.66" / 220 mm	1.20" x 1.20" x 0.47" / 30.5 x 30.5x 12.0 mm	0.47" / 12.0 mm	0.47" / 12.0 mm	0.315" / 8.0 mm
29-23012	Blue	Triple Nylon / Stainless Steel	8.66" / 220 mm	1.20" x 1.20" x 0.47" / 30.5 x 30.5 x 12.0 mm	0.47" / 12.0 mm	0.47" / 12.0 mm	0.47" / 12.0 mm

PACKAGING: Bag of 3 - For a bag of 50, add "-50" following the item code.

ALL PURPOSE REUSABLE BRUSHES

Plastic Handle Brushes

We offer a wide selection of brushes for instrument cleaning applications. These brushes are ideal for large bore, uniquely shaped, and endoscopic instruments.

In-Brush™ Plastic Handle

Available in regular toothbrush style (nylon or stainless) and double-ended brushes.



PRODUCT SELECTION

ITEM #	DESCRIPTION	COLOR	BRUSH LENGTH	BRISTLE SIZE	PACKAGING
28-72550	Plastic Single End Nylon Reusable Handle Brush	Black	184 mm	38 mm x 12.7 mm	Bag of 3
29-72550	Plastic Single End Stainless Steel Reusable Handle Brush	Black	184 mm	38 mm x 12.7 mm	Bag of 3
29-70050	Plastic Single End Brass Reusable Handle Brush	Black	184 mm	38 mm x 12.7 mm	Bag of 3
28-69394	Plastic Double End Nylon Reusable Handle Brush	Black	175 mm	35 mm x 10 mm 25 mm x 5 mm	Bag of 3
29-69394	Plastic Double End Stainless Steel Reusable Handle Brush	Black	175 mm	35 mm x 10 mm 25 mm x 5 mm	Bag of 3
50-40911	Hand Scrub Nylon Brush Reusable autoclavable	White	104 mm	95 mm x 28 mm	Bag of 3

Reusable models are also available in bags of 50. Add the suffix "- 50"

NEW!

Engineered for precision and durability, this reusable dual-ended is designed to tackle even the most complex instruments. Its ultra-fine tip reaches into the smallest, hardest-to-access areas—delivering exceptional hygiene with every use. Built from long-lasting materials, this brush not only enhances cleaning performance but also supports hospital sustainability goals by reducing single-use waste and associated costs.



PRODUCT SELECTION

ITEM #	DESCRIPTION	COLOR	BRUSH LENGTH	BRISTLE SIZE 1	BRISTLE SIZE 2	PACKAGING
50-40517	Reusable Extra Rigid Flat Brush Dual Ended	Blue	183 mm	38 mm x 9 mm x 10 mm	28 mm x 3 mm x 8 mm	Bag of 3
50-40518	Reusable Rigid Flat Brush Dual Ended	Green	183 mm	38 mm x 9 mm x 15 mm	28 mm x 3 mm x 13 mm	Bag of 3
50-40520	Reusable Extra Flexible Flat Brush Dual Ended	Red	183 mm	38 mm x 9 mm x 25 mm	28 mm x 3 mm x 23 mm	Bag of 3

Reusable models are also available in bags of 50. Add the suffix "- 50"

ACETABULAR REAMER BRUSH

Our double-ended cleaning brush is manufactured exclusively for cleaning acetabular bone reamers. The “Pig Tail” rounded style circular head brush is designed to effectively clean inside the reamer, while the rigid straight brush on the reverse side is used for cleaning individual cutting holes. Special fusion of the brush heads to the hand piece along with its high quality antimicrobial medical grade bristles ensures a durable, reusable, and highly effective cleaning tool.



STRAIGHT BRUSH



ROUND BRUSH



PRODUCT SELECTION


ITEM #	DESCRIPTION	HANDLE LENGTH	STRAIGHT BRUSH	ROUND BRUSH HEAD	PACKAGING
50-ACETAB	Acetabular Reamer Cleaning Brush	4.72" / 120 mm Ø 0.47" / 12 mm	1.18" / 30 mm Ø 0.3125" / 8 mm	2" / 50 mm Ø 1.75" / 45 mm	Bag of 3

NOTE: Specifications subject to small variances

REUSABLE ORTHOPEDIC RASPS CLEANING BRUSH

Our Reusable Orthopedic Rasps Cleaning Brush is used to clean orthopedic rasps and other similar instruments used in orthopedic surgery.

- Long lasting and can withstand many sterilization cycles
- Ergonomic steamable thermoplastic handle
- Antibacterial polyamide bristles
- Reusable
- Latex free



PRODUCT SELECTION

ITEM #	TOTAL LENGTH	HANDLE LENGTH	BRUSH SIZE	PACKAGING
50-RASP	9" / 23 cm	5.5" / 14 cm	0.78" x 2.7" x 3.5" 2 cm x 7 cm x 9 cm	Bag of 3


DISTRIBUTED BY:
ADVANCED SURGI-PHARM INC.
www.advancedsurgipharma.com
Tel.1-800-661-5432



BOTTLE & DRAIN CLEANING BRUSHES

Designed specifically for hospital applications, our Reusable Instrument Cleaning Brushes feature long and flexible bristles to maximize the effectiveness of large-surface cleaning.

- High quality polyamide clear bristles
- Easy to identify debris
- Autoclavable
- Latex free



PRODUCT SELECTION

BLUE DELRIN PLASTIC

ITEM #	DESCRIPTION	HANDLE LENGTH	BRUSH LENGTH	BRUSH (Ø)	TOTAL LENGTH	PACKAGING
50-10295	Reusable Jar / Bottle Brush	120 mm ± 2 mm	120 mm ± 5 mm	70 mm ± 2 mm	13.19" / 355 mm	Bag of 3
50-10224	Sterilizer Drain Cleaning Brush	N/A	4" / 10.1 cm	2" / 5.08 cm	24" / 609 mm	Bag of 3

AUTOCLAVE CLEANING BRUSH KITS



ROUND BRUSH

- Features 3 holes for adjustable pole placement
- Adjustable brush placement (parallel or perpendicular) and longer brush length provides easy access to hard-to-reach corners
- Made with FDA approved polyester bristles

HAND BRUSH

- Ergonomic handle for better hand comfort
- Longer, FDA approved polyester bristles
- Easy to clean (no cavity)

PRODUCT SELECTION

REGULAR AUTOCLAVE CLEANING BRUSH KIT

ITEM # (KIT)	DESCRIPTION	ITEM # (UNIT)	DESCRIPTION	DIMENSIONS
50-CLAVEREG	Autoclave REGULAR Brush Cleaning Kit	50-35900	1 Hand Brush	Length: 3.5" x 8.3" / 9 cm x 21 cm
		50-55800	1 Rounded Corner Brush (Head Brush)	Head size: 5.5 x 9.8" / 14 x 25 cm
		50-EXTPOLE	1 Anodized Aluminum Extension Pole	2-4 ft / 68-120 cm

SUPERIOR AUTOCLAVE CLEANING BRUSH KIT

ITEM # (KIT)	DESCRIPTION	ITEM # (UNIT)	DESCRIPTION	DIMENSIONS
50-CLAVESUP	Autoclave SUPERIOR Brush Cleaning Kit	50-35900	1 Hand Brush	Length: 3.5" x 8.3" / 9 cm x 21 cm
		50-55800	1 Rounded Corner Brush (Head Brush)	Head size: 5.5 x 9.8" / 14 x 25 cm
		50-EXTPOLE	1 Anodized Aluminum Extension Pole	2-4 ft / 68-120 cm
		SF-100	2 SAVEFACE™ Face Shields	9" / 23 cm

REPLACEMENT COMPONENTS



HAND BRUSH
50-35900

Packaging: 1



BRUSH HEAD
50-55800

Packaging: 1



EXTENSION POLE
50-EXTPOLE

Packaging: 1





SAVEFACE™
FACE SHIELDS SF-100

Packaging: 100 / Box


SILICONE PIN MAT

Our **Silicone Pin Mat** cushions, protects, and organizes delicate surgical instruments for sterilization. The elevated nubs provide greater air-flow ventilation and drying, significantly reducing residual moisture collection.

- Reusable
- Autoclavable
- Latex free



FRONT SIDE
Pins keep instruments securely in place



BACK SIDE
Numerous holes and elevated nubs provide superior venting and reduce water collection

PRODUCT SELECTION

ITEM #	DIMENSIONS	PACKAGING
13-2001	13" x 18" (33 cm x 45.72 cm) 11.97 mm high from base of mat to end of pin 15.53 mm high from base of elevated nub to end of pin	1 unit
13-2005	13" x 18" (33 cm x 45.72 cm) 11.97 mm high from base of mat to end of pin 15.53 mm high from base of elevated nub to end of pin	5 units

GOOD TO KNOW:

CUSTOM SIZES ARE AVAILABLE UPON REQUEST

TAG-IT™ | NOTES

WHEN SIMPLY BROKEN DOESN'T CUT IT—TAG IT, TELL IT, OWN IT!

Tag-It™ Notes offer the freedom to **identify and document instructions for instruments, trays, and carts** before and after **sterilization**. Tag-It™ Notes are **tear-resistant, withstand extreme chemicals, high temperature** and ensures **reliable communication**.

WHY YOU'LL LOVE THEM:

- Tag-It™ Notes gives you the **freedom** to clearly document important information **without the clutter**.
- Robust and easy-to-use **tie wraps** designed to **securely attach** to instruments (134°C & chemical tested).
- Includes a **fast-drying, heat-resistant permanent marker** for writing your important notes.
- Choose from other **customizable cards** for personalized needs or **pre-made cards** for fast, convenient use.
- When you need more than a "repair" or "sharpen" seal to identify an issue. Instead, **Tag it, Tell it, Own it.**

WATERPROOF
The synthetic paper withstands spills and liquid exposure without any damage.

ECO-FRIENDLY
Made from sustainable materials.

HEAT RESISTANT
Can resist high temperatures of up to 273.2°F (134°C). Compatible with Gas Plasma Sterilization (STERAD®100 NX Sterilization).

DURABLE & TEAR-PROOF
Handles rigorous cleaning, tough conditions, and resists tearing.

PRODUCT SELECTION

ITEM #	COLOR	SIZE	PACKAGING
TI-NOTE-RO	Yellow	3" x 5" / 8 cm x 13cm	100 / PKG


All Tag-It™ are single-use. **PACKAGING: 1 pack of 100 tags, 1 permanent marker, and 100 ties.**



Tag-It™ Cards are an effortless way to add information with either pre-printed or blank formats. These **durable** tags streamline your work ow, **withstand tough conditions**, and feature **(optional) indicator dots** for quick verification. They help keep **instruments, trays, and carts** organized, ensuring efficiency before and after sterilization/decontamination.

KEY PRODUCT FEATURES:

- **Special Permanent Marker:** The ink will not run or appear blurry after exposure to chemicals or high temperature (max. 260C/500°F) sterilization. *(Regular pens or markers can be used only if the tags are not going to be exposed to any sterilization/decontamination).*
- **Tie Wraps:** Assist in securing the tags to almost any instrument or container. *(Tie wraps have been tested for chemical and steam sterilization).*
- **Our Tag-It™:** Will not shrink, absorb water/chemicals, or lose its intended integrity after chemical or steam sterilization.



(Optional)
All Tags Available with Indicator Dots
 Provides visual confirmation that the tag was processed.

SCAN QR TO WATCH OUR VIDEO



PRODUCT SELECTION:

ITEM #	COLOR	DESCRIPTION	SIZE	PACKAGING
TI-3502	Blue	VENDOR	3" x 5" / 8 cm x 13cm	100 / PKG
TI-3503	Green	LOANER	3" x 5" / 8 cm x 13cm	100 / PKG
TI-3504	Red	PRIORITY	3" x 5" / 8 cm x 13cm	100 / PKG
TI-3505	Yellow	RE-RUN	3" x 5" / 8 cm x 13cm	100 / PKG
TI-3506	Pink	MISSING	3" x 5" / 8 cm x 13cm	100 / PKG
TI-3507	Orange	BROKEN	3" x 5" / 8 cm x 13cm	100 / PKG
TI-3509	Black	BLANK	3" x 5" / 8 cm x 13cm	100 / PKG
TI-35B1	White	BLANK	3" x 5" / 8 cm x 13cm	100 / PKG
TI-35B4	Red	BROKEN	3" x 5" / 8 cm x 13cm	100 / PKG
TI-35C22	Light Blue	VENDOR CHECK IN + OUT	3" x 5" / 8 cm x 13cm	100 / PKG
TI-35L3	Green	LOANER CHECK IN + OUT	3" x 5" / 8 cm x 13cm	100 / PKG
TI-35T2	Blue	QUICK TURNOVER	3" x 5" / 8 cm x 13cm	100 / PKG
TI-PM09	Black	FINE PERMANENT MARKERS	1.5" Ø x 5.5" / 3.81cm Ø x 12.7 cm	36 / Bag

PACKAGING: 1 pack of 100 tags, 100 tie wraps, and 1 permanent marker.

OPTIONAL PROCESSING INDICATOR DOTS: Add suffix "ST" to end of code for steam dots, "PS" for plasma dots.

Available in other languages or customize your own.

PRODUCT SELECTION: (customization also available)


VENDOR / BLUE

ITEM #	SIZE
TI-3502	3" x 5" (8 cm x 13cm)




NOTES / YELLOW

ITEM #	SIZE
TI-NOTE-RO	3" x 5" (8 cm x 13cm)



LOANER / GREEN

ITEM #	SIZE
TI-3503	3" x 5" (8 cm x 13cm)



BLANK / WHITE

ITEM #	SIZE
TI-35B1	3" x 5" (8 cm x 13cm)



PRIORITY / RED

ITEM #	SIZE
TI-3504	3" x 5" (8 cm x 13cm)



BROKEN / RED

ITEM #	SIZE
TI-35B4	3" x 5" (8 cm x 13cm)



RE-RUN / YELLOW

ITEM #	SIZE
TI-3505	3" x 5" (8 cm x 13cm)



VENDOR IN + OUT / LIGHT BLUE

ITEM #	SIZE
TI-35C22	5" x 3" (13 cm x 8 cm)




MISSING / PINK

ITEM #	SIZE
TI-3506	3" x 5" (8 cm x 13cm)



LOANER IN + OUT / GREEN

ITEM #	SIZE
TI-35L3	5" x 3" (13 cm x 8 cm)




BROKEN/ ORANGE

ITEM #	SIZE
TI-3507	3" x 5" (8 cm x 13cm)



BLANK / BLACK

ITEM #	SIZE
TI-3509	3" x 5" (8 cm x 13cm)



QUICK TURNOVER / BLUE

ITEM #	SIZE
TI-35T2	3" x 5" (8 cm x 13cm)



All information tags are single-use. PACKAGING: 1 pack of 100 tags, 100 tie wraps, and 1 permanent marker.

TAG-IT™ | SCOPE

Batrik Scope Tags are tear-resistant, double-sided tags that will prove indispensable throughout the reprocessing journey. Our tags make it easy to **record important information**, including processing dates, with a user-friendly feature that allows you to select the month, day, and date. Crafted from **waterproof synthetic paper**, these tags feature a **self-looping design**, making it effortless to attach them to the scope without the need for additional fastening materials. Elevate your reprocessing protocols with Batrik!

KEY FEATURES

- The paper, tie wraps, and smudge-proof ink are designed to withstand high temperatures of up to 273.2°F (134°C) and are compatible with Gas plasma sterilization (STERRAD®100 NX Sterilization System).
- Tags include the month, day, and date selection on one side, making it convenient to mark the processing date accurately



COLOR-CODED AND CUSTOMIZABLE

- Available in seven colors
- The tags facilitate color coding for easy identification

SELF LOOPING

- The self-looping design allows for easy attachment to the scope
- Eliminating the need for additional fastening materials


SINGLE-USE AND LATEX FREE

- Batrik Scope Tags are designed for single use and are latex-free
- Promoting hygiene and reducing the risk of cross-contamination

PRODUCT SELECTION


ITEM #	COLOR	SIZE	PACKAGING
TI-ST05	Yellow	11.8" x 1.55" / 33.8 cm x 4 cm	100 / PKG
TI-ST02	Blue	11.8" x 1.55" / 33.8 cm x 4 cm	100 / PKG
TI-ST03	Green	11.8" x 1.55" / 33.8 cm x 4 cm	100 / PKG
TI-ST07	Orange	11.8" x 1.55" / 33.8 cm x 4 cm	100 / PKG
TI-ST06	Purple	11.8" x 1.55" / 33.8 cm x 4 cm	100 / PKG
TI-ST04	Red	11.8" x 1.55" / 33.8 cm x 4 cm	100 / PKG
TI-ST01	White	11.8" x 1.55" / 33.8 cm x 4 cm	100 / PKG

All information tags are single-use. PACKAGING: 1 pack of 100 tags and 1 permanent marker.




TAG-IT™ | SCOPE REPROCESSING

Batrik Scope Reprocessing Tags are durable and can withstand handling during reprocessing **without tearing**, ensuring that information remains intact. These tags are vital in recording important information such as **date to reprocess on, scope model and serial number, sterilizer type, meeting sterilizer parameter conditions, alcohol flush, and hang date**. These tags have been meticulously designed to elevate **operational efficiency** and contribute to the delivery of the highest standard of patient care.



TI-RST



TI-MPB

*Suitable for various information tagging purposes

LOGGING IMPORTANT INFORMATION

- Easy for staff to indicate when the scope was processed
- Key scope, technician, and cycle-related information easily reported and identified.

SMART DESIGN

- The self-looping design allows for easy attachment to the scope
- No requirement for extra fastening materials
- Secure fastening as the loop is tightened around the instrument
- Features an integrated barcode area



PRODUCT SELECTION

ITEM #	DESCRIPTION	DIMENSIONS	PACKAGING
TI-RST	REPROCESSING SCOPE TAG	3" x 10.5" / 8 cm x 27 cm	100 / PKG
TI-MPB	MULTI-PURPOSE BLANK TAG	3" x 10.5" / 8 cm x 27 cm	100 / PKG

All information tags are single-use. PACKAGING: 1 pack of 100 tags and 1 permanent marker.

CORNER POCKETS™

With three models to choose from, our **Corner Pockets™** portfolio is a truly comprehensive line of products designed to prevent tearing associated with sharp corners or feet on surgical trays.



18-CP02



18-CP07




18-CP06

COLOUR	BLUE	ORANGE	PINK
USAGE	SINGLE-USE	REUSABLE <small>withstands repeated sterilization up to 134°C</small>	REUSABLE <small>withstands repeated sterilization up to 134°C</small>
DIMENSIONS	2.8" x 4" (7.3 cm x 10.3 cm)	2.8" x 4" (7.3 cm x 10.3 cm)	3.4" x 4.7" (8.5 cm x 12 cm)
LATEX-FREE	✓	✓	✓
NO MOISTURE RETENTION	✓	✓	✓
FABRICATION	Flexible, silicone-like plastic	Flexible, silicone-like plastic	Our most flexible plastic for easiest placement and most secure protection

SOFT POCKETS™

Our **Soft Pockets™** Foam Corner Protectors are made from high density foam to effectively cushion and protect the corners of surgical trays to prevent tearing. They are especially useful with heavier instrument sets.



PRODUCT SELECTION

ITEM #	DESCRIPTION	PACKAGING	CASE DIMENSIONS
18-CP02**	Corner Pockets™ Single-Use, Blue	240 / Bag	80 / Box (24 Boxes / Pallet)
18-EVAL2**	Corner Pockets™ Single-Use, Blue	24 / Bag	-
18-CP07	Corner Pockets™ Reusable, Orange	80 / Bag	-
18-CORNER	Corner Pockets™ Reusable, Orange	320 / Box (4 Bags of 80 pcs)	-
18-EVAL7	Corner Pockets™ Reusable, Orange	24 / Bag	-
18-CP06	Corner Pockets™ Reusable, Pink	80 / Bag	-
18-CPF-02	Soft Pockets™ Single-Use Foam Corner Protectors	80 / Box (24 Boxes / Pallet)	-

** Although the single use Corner Pockets™ appear similar to the reusable models they are blended differently and have only been validated as single-use. Re-using them may compromise sterility.

LILA'S INSTRUMENT STAIN REMOVER

Lila's special blend of naturally sourced powders provides outstanding removal of stubborn stains.

Convenient-to-use formula, you can achieve spotless results every time!

It can be used on surgical/dental instruments, autoclaves, ultrasonic cleaners, stainless steel countertops, trays, bowls and more.

Specifically designed for use on stainless steel, chrome, nickel, brass, copper, and silver.



LSR-8501


EASE TO USE AS 1,2,3...

1. Dampen a cloth or grab a Lila's Low Linting Dura-Towel and dab into the jar.
2. Gently rub the area that is stained.¹
3. Rinse thoroughly for a flawless shine.²

A little Goes a long Way!

Even great to tackle stubborn stains on steel pots and pans in your kitchen with ease.

Lila's - If you have used Spectra-Scrub or Integra's instrument stain remover then you will love Lila's!



60-46003

LILA'S SCRUB PAD

- Say goodbye to extra stubborn stains with Lila's Scrub Pad without leaving any metal residue.
- Simply add a small amount of Lila's Stain Remover with one of our Scrub Pad's for very large areas.
- The gentle yet effective scrub is the perfect choice for Ophthalmic Instruments.¹



17-DT1213C

LILA'S LOW-LINTING DURA-TOWEL

- Lila's Dura-Towel has long lasting strength and has a minimal lint composition.
- These single-use towels are clean room certified and are specifically designed for instrument wiping.
- Its ideal dimensions and unique attributes guarantee optimal performance, setting it apart from ordinary paper towels.

PRODUCT SELECTION

ITEM #	DESCRIPTION	DIMENSIONS	PACKAGING
LSR-8501	LILA'S INSTRUMENT STAIN REMOVER 3 oz (85g) plastic jar with a convenient screw-on lid	Ø 2.8" / 70 mm	1 Jar / Bag
LSR-8512	LILA'S INSTRUMENT STAIN REMOVER 3 oz (85g) plastic jar with a convenient screw-on lid	Ø 2.8" / 70 mm	12 Jars / Box
17-DT1213C	LILA'S LOW LINTING DURA-TOWEL 1080 towels	12" x 13" / 31 cm x 33 cm Folded: 3" x 3.25" / 7.5 cm x 8.25 cm	90 / Bag 12 Bags / Box
60-46003	LILA'S SCRUB PAD	4" x 6" / 10 cm x 15 cm	3 / Bag
LSR-BUNDLE	LILA'S STAIN REMOVER + LILA'S DURA-TOWEL 3 oz (85g) plastic jar + 90 towels	Ø 2.8" / 70 MM 12" x 13" / 31 cm x 33 cm	1 Jar / Bag 90 / Bag


1 : Although Lila's is completely safe, it is not recommended for use on ophthalmic instruments. Use Lila's Scrub Pad on its own.
2: The stain remover should not be mixed with bleach and is for external use only to remove instrument stains.

INSTRUMENT PROTECTORS


Our **Instrument Protectors** protect the tips and edges of instruments to prevent damage during handling and sterilization. Available in a wide array of color coded sizes for almost every application, including osteotome, retractor, double, and perforated.

- Steam and gas permeable
- Perforated protectors reduce vacuum seal and moisture build-up
- Validated for pre-vacuum sterilization*
- Multiple color-coded sizes
- Latex Free
- Non-Sterile

REGULAR / CLEAR



PERFORATED / CLEAR



REGULAR / TINTED



PERFORATED / TINTED



REGULAR / SOLID COLOR



PERFORATED / SOLID COLOR



*Independent laboratory tested to validate pre-vacuum steam sterilization efficacy in a healthcare setting

PRODUCT SELECTION

REGULAR & PERFORATED / CLEAR



RC-07



PC-07

ITEM # REGULAR	ITEM # PERFORATED	COLOR	DIMENSIONS
RC-01	PC-01	Clear	0.062" x 0.75" (1.57 mm x 19 mm)
RC-02	PC-02	Clear	0.075" x 0.75" (2 mm x 19 mm)
RC-03	PC-03	Clear	0.109" x 0.75" (2.77 mm x 19 mm)
RC-04	PC-04	Clear	0.125" x 1" (3.18 mm x 25.4 mm)
RC-05	PC-05	Clear	0.187" x 1" (4.75 mm x 25.4 mm)
RC-06	PC-06	Clear	0.406" x 0.75" (10.3 mm x 19 mm)
RC-07	PC-07	Clear	0.062" x 0.375" x 1" (1.57 mm x 9.53 mm x 25.4 mm)
RC-08	PC-08	Clear	0.078" x 0.625" x 1" (2 mm x 15.875 mm x 25.4 mm)
RC-09	PC-09	Clear	0.125" x 1" x 1" (3.18 mm x 25.4 mm x 25.4 mm)
RC-10	PC-10	Clear	Assorted

PACKAGING: Bag of 100

REGULAR & PERFORATED / TINTED



RT-07



PT-07

RT-01	PT-01	Tinted White	0.062" x 0.75" (1.57 mm x 19 mm)
RT-02	PT-02	Tinted Blue	0.075" x 0.75" (2 mm x 19 mm)
RT-03	PT-03	Tinted Green	0.109" x 0.75" (2.77 mm x 19 mm)
RT-04	PT-04	Tinted Red	0.125" x 1" (3.18 mm x 25.4 mm)
RT-05	PT-05	Tinted Yellow	0.187" x 1" (4.75 mm x 25.4 mm)
RT-06	PT-06	Tinted Maroon	0.406" x 0.75" (10.3 mm x 19 mm)
RT-07	PT-07	Tinted Orange	0.062" x 0.375" x 1" (1.57 mm x 9.53 mm x 25.4 mm)
RT-08	PT-08	Tinted Brown	0.078" x 0.625" x 1" (2 mm x 15.875 mm x 25.4 mm)
RT-09	PT-09	Tinted Black	0.125" x 1" x 1" (3.18 mm x 25.4 mm x 25.4 mm)
RT-10	PT-10	Tinted Assorted	Assorted

PACKAGING: Bag of 100

REGULAR & PERFORATED / SOLID



RS-07



PS-07

RS-01	PS-01	Solid White	0.062" x 0.75" (1.57 mm x 19 mm)
RS-02	PS-02	Solid Blue	0.075" x 0.75" (2 mm x 19 mm)
RS-03	PS-03	Solid Green	0.109" x 0.75" (2.77 mm x 19 mm)
RS-04	PS-04	Solid Red	0.125" x 1" (3.18 mm x 25.4 mm)
RS-05	PS-05	Solid Yellow	0.187" x 1" (4.75 mm x 25.4 mm)
RS-06	PS-06	Solid Maroon	0.406" x 0.75" (10.3 mm x 19 mm)
RS-07	PS-07	Solid Orange	0.062" x 0.375" x 1" (1.57 mm x 9.53 mm x 25.4 mm)
RS-08	PS-08	Solid Brown	0.078" x 0.625" x 1" (2 mm x 15.875 mm x 25.4 mm)
RS-09	PS-09	Solid Black	0.125" x 1" x 1" (3.18 mm x 25.4 mm x 25.4 mm)
RS-10	PS-10	Solid Assorted	Assorted

PACKAGING: Bag of 100

ITEM #	DIMENSIONS
OP-11	0.25" x 0.75" (6.35 mm x 19.05 mm)
OP-12	0.375" x 1" (9.53 mm x 25.4 mm)
OP-13	0.5" x 1" (12.7 mm x 25.4 mm)
OP-14	0.75" x 1" (19.05 mm x 25.4 mm)
OP-15	1" x 1" (25.4 mm x 25.4 mm)
OP-16	1.25" x 1.25" (31.75 mm x 31.75 mm)
OP-17	1.5" x 1.5" (38.1 mm x 38.1 mm)
OP-18	1.75" x 1.25" (44.45 mm x 31.75 mm)
OP-20	Assorted

PACKAGING: Bag of 50

ITEM #	DIMENSIONS	PACKAGING
RGS-01	1.5" (38 mm)	55 / Bag
RGS-02	2" (51 mm)	45 / Bag
RGS-03	3" (76.2 mm)	40 / Bag


ITEM #	COLOR	SIZE
RGM-01	Pink	Small
RGM-02	Clear	Small
RGM-03	Green	Medium
RGM-04	Clear	Medium

PACKAGING: Bag of 50

ITEM #	COLOR	DIMENSIONS
DG-01	Clear	0.0625" x 0.75" (1.59 mm x 19.05 mm)
DG-02	Blue	0.078" x 0.75" (1.98 mm x 19.05 mm)
DG-03	Green	0.109" x 0.75" (2.78 mm x 19.05 mm)
DG-04	Red	0.125" x 1" (3.18 mm x 25.4 mm)
DG-05	Yellow	0.187" x 1" (4.76 mm x 25.4 mm)
DG-07	Orange	0.062" x 0.375" x 1" (1.59 mm x 9.53 mm x 25.4 mm)
DG-10	Assorted	Assorted

PACKAGING: Bag of 50


OSTEOTOME PROTECTORS




SINGLE-TIP RETRACTOR PROTECTORS



WEITLANER RETRACTOR PROTECTORS



PERFORATED DOUBLE PROTECTORS



PrAmLi Power³ LUBE™ Instrument Lubricant

Specially formulated to lubricate surgical instruments, hinges, box locks, and general equipment during autoclaving. PrAmLi Power³ LUBE™ is uniquely formulated with an all natural corrosion inhibitor.

- Forms a protective coating which can extend the instrument's, useful life
- Bacteriostatic
- Neutral pH 6.5 - 7.5
- Prevents corrosion

TIP:
FOR MILK BATH SOLUTION,
DILUTE 1 PART SOLUTION
IN 6 PARTS WATER

Formulated with
PrAmLi Power³ Enzymes™

ITEM #	DIMENSIONS
OP-11	0.25" x 0.75" (6.35 mm x 19.05 mm)
OP-12	0.375" x 1" (9.53 mm x 25.4 mm)
OP-13	0.5" x 1" (12.7 mm x 25.4 mm)
OP-14	0.75" x 1" (19.05 mm x 25.4 mm)
OP-15	1" x 1" (25.4 mm x 25.4 mm)
OP-16	1.25" x 1.25" (31.75 mm x 31.75 mm)
OP-17	1.5" x 1.5" (38.1 mm x 38.1 mm)
OP-18	1.75" x 1.25" (44.45 mm x 31.75 mm)
OP-20	Assorted



PACKAGING: Bag of 50

PRODUCT SELECTION

ITEM #	DILUTION LUBRICANT / WATER	PACKAGING	CASE DIMENSIONS
SS-21-6001	1 to 4 ml / L (0.125 to 0.5 oz / Gallon)	1 x 3.8 L (1 Gallon)	6" x 6" x 13"
SS-21-6004	1 to 4 ml / L (0.125 to 0.5 oz / Gallon)	4 x 3.8 L (4 Gallons)	13" x 13" x 14"
SS-21-6025	1 to 4 ml / L (0.125 to 0.5 oz / Gallon)	2 x 9.5 L (2 x 2.5 Gallons)	14.5" x 9.25" x 15"
SS-21-6032	No dilution required	1 x 946 ml (32 oz) Spray bottle	10" x 4"
SS-21-6012	No dilution required	12 X 946 ml (32 oz) Spray Bottles	18" x 12" x 11"


S.T.I.C.S.

Our Simple Tracking Identification Color System, or **S.T.I.C.S.**, is an easy to use, latex free, tool that allows you to quickly identify your instruments by department and/or procedure. The tape adheres securely and every tape sheet comes with a protective cover to ensure it lasts as long as you need.



Each 8.5" x 11" (216 mm x 279 mm) sheet contains the following adhesive strips:

- 1" (25.4 mm) x 220
- 1.5" (38.1 mm) x 44
- 2" (50.8 mm) x 44
- Total tape per sheet 374" (9500 mm)



PRODUCT SELECTION

SOLID COLORS

ITEM #	COLOR	
T-850-2-01	White	
T-850-2-02	Royal Blue	
T-850-2-03	Kelly Green	
T-850-2-04	Red	
T-850-2-05	Yellow	
T-850-2-06	Purple	
T-850-2-07	Orange	
T-850-2-08	Brown	
T-850-2-09	Black	
T-850-2-10	Gray	
T-850-2-11	Lavender	
T-850-2-12	Beige	
T-850-2-13	Light Blue	
T-850-2-14	Lime Green	
T-850-2-15	Pink	
T-850-2-17	Plum	
T-850-2-18	Mint Green	

STARS

ITEM #	COLOR	
T-850-6-01	Stars Fuchsia	
T-850-6-02	Stars Medium Blue	
T-850-6-03	Stars Green	
T-850-6-04	Stars Violet	

TRIANGLES

ITEM #	COLOR	
T-850-10-01	Triangle Fuchsia	
T-850-10-04	Triangle Violet	

WHITE STRIPES

ITEM #	COLOR	
T-850-4-02	White/Royal Blue	
T-850-4-03	White/Kelly Green	
T-850-4-04	White/Red	
T-850-4-05	White/Yellow	
T-850-4-07	White/Orange	
T-850-4-08	White/Brown	
T-850-4-11	White/Lavender	
T-850-4-13	White/Light Blue	
T-850-4-14	White/Lime Green	
T-850-4-15	White/Pink	
T-850-4-17	White/Plum	

BLACK STRIPES

ITEM #	COLOR	
T-850-3-01	Black/White	
T-850-3-03	Black/Kelly Green	
T-850-3-04	Black/Red	
T-850-3-05	Black/Yellow	
T-850-3-07	Black/Orange	
T-850-3-10	Black/Gray	
T-850-3-11	Black/Lavender	
T-850-3-12	Black/Beige	
T-850-3-13	Black/Light Blue	
T-850-3-14	Black/Lime Green	
T-850-3-15	Black/Pink	
T-850-3-17	Black/Plum	

DOTS

ITEM #	COLOR	
T-850-1-01	Dots Fuchsia	
T-850-1-03	Dots Green	

FOUR-COLOR

ITEM #	COLOR	
T-850-2-31	Red/Yellow/ Purple/Orange	
T-850-2-32	Blue/White/ Black/Green	

TRI-COLOR

ITEM #	COLOR	
T-850-5-51	Red/White/Blue	
T-850-5-52	Red/White/Black	
T-850-5-54	Red/Black/Green	

STRIPES

ITEM #	COLOR	
T-850-7-01	Blue/Green	
T-850-7-02	Blue/Red	
T-850-7-03	Yellow/Blue	
T-850-7-04	Yellow/Red	
T-850-7-05	Yellow/Green	

SPECIALTY

ITEM #	COLOR	
T-850-8-01	Red Hearts/White	
T-850-8-02	Smiley Yellow /Black	
T-850-8-03	Alphabet Fuchsia /White	
T-850-8-04	White Hearts/Red	
T-850-9-03	Diagonal Blue /Black	
T-850-9-04	Diagonal Red /Black	

NUMBERS


ITEM #	COLOR	
T-850-11-01	Numbers Fuchsia	

ID TAPE










A simple solution for marking surgical instruments with the aid of color-coding, our ID Tape is durable and will not chip, flake, or peel. It can be used in autoclaves, withstanding repeated exposure to sterilization. Independent validation tests demonstrated the efficiency of sterilization of marked instruments, even underneath the Batrik ID Tape.

AVAILABLE WIDTHS



PRODUCT SELECTION

ITEM # 0.125" / 318 MM	ITEM# 0.25" / 635 MM	COLOR
T-125-01	T-250-01	White
T-125-02	T-250-02	Blue
T-125-03	T-250-03	Green
T-125-04	T-250-04	Red
T-125-05	T-250-05	Yellow
T-125-07	T-250-07	Orange
T-125-09	T-250-09	Black

-  Quickly identifies instruments by department or procedure
-  Easy to apply
-  Withstands repeated exposure to sterilization
-  Has more per tape roll (300" / 7620 mm)
-  Available in 7 colors and 2 widths
-  Conveniently packaged in a clear box
-  Non-Sterile

TAPE-B-GONE™ TAPE REMOVER

Tape-B-Gone™ is a uniquely formulated and ready-to-use proprietary formula that easily removes adhesives from medical equipment and instruments. Unlike many other brands, Tape-B-Gone™ does not contain highly dangerous or caustic chemicals, is not flammable or combustible, and is safe for use on metal, glass, plastic, and rubber without any special equipment or ventilation.



Easy to use Spray

Apply liberally to any label or adhesive. Allow to soak for 3-5 minutes. Lift corner and peel tape off slowly. Wipe clean with a damp cloth and rinse off any remaining formula

Not flammable or combustible

We always recommend incorporating best hygiene and safety practices when handling any chemicals

Trade secret formula with a neutral pH of 7.0 - 8.0 surfactants

Safe for metal, glass, plastic and rubber

Safe for land transport (DOT), sea transport (IMDG/IMO), and air transport (IATA)

PRODUCT SELECTION

ITEM #	DESCRIPTION	PACKAGING
T-08RMV	8 oz (237 ml) bottle	1 Unit
T-32RMV	8 oz (237 ml) bottle	4 / Pack
T-SAMP-RMV	0.5 oz (15 ml) bottle	1 Unit
T-96RMVSC-R0	8oz (237 ml) bottle Screwcap*	12 / Box


PADLOCKS SECURITY TAGS

Our **Padlock Security Tags** for instrument containers and general hospital applications are specifically designed for the healthcare market. Our material selection, fabrication process, and engineering ensure consistent quality and high temperature compliance with the various steam sterilization temperatures in worldwide hospital applications.


Available with optional reactive ink indicator dots that provide visual sterility indication. Light blue dots turn brown/dark brown when steam sterilized.

AVAILABLE IN THREE FORMATS


700-SERIES



800-SERIES



900-SERIES



Made of polypropylene



Latex and DEHP Free



Will not melt or soften with any established steam sterilization protocol



Fully open design (45°) facilitates insertion into instrument containers



Easy twist breaking design



Large base for easier handling



Available in different size



Custom Colors and OEM available upon request

PRINTED & CUSTOM NUMBERED PADLOCKS

Available in 700, 800 & 900 SERIES

PRODUCT SELECTION

COLOR	PRINTED WORD
PRINTED PADLOCKS	
WHITE	CLEAN
BLUE	STERILE
RED	REPAIR
YELLOW	SHARPEN
YELLOW	LOANER
BLACK	DISCARD
CUSTOM NUMBERED	
ALL COLORS	up to 8 numbers



CUSTOM NUMBERED
Custom printed padlocks (up to 8 characters)



PRINTED
6 different printed labels available

NOTE: Printed & Custom Numbered Padlocks are available in 700, 800, and 900 Series and can be ordered in 100, 200, and 500 unit quantities. Please contact us for ordering details.

INDICATOR DOTS

REGULAR


700-SERIES PADLOCKS

The 700 series padlocks are perfect for "Aesculap" containers and those with a short, large opening.

PRODUCT SELECTION

REGULAR PADLOCKS

COLOR	100 UNITS	200 UNITS	500 UNITS
WHITE	15-701	15-711	15-710
BLUE	15-702	15-712	15-720
GREEN	15-703	15-713	15-730
RED	15-704	15-714	15-740
YELLOW	15-705	15-715	15-750
PINK	15-706	15-716	15-760
ORANGE	15-707	15-717	15-770
BLACK	15-709	15-719	15-790
LIGHT BLUE	15-722	15-732	15-742

WITH INDICATORS DOTS

COLOR	100 UNITS	200 UNITS	500 UNITS
WHITE	15-701-IND	15-711-IND	15-710-IND
BLUE	15-702-IND	15-712-IND	15-720-IND
GREEN	15-703-IND	15-713-IND	15-730-IND
RED	15-704-IND	15-714-IND	15-740-IND
YELLOW	15-705-IND	15-715-IND	15-750-IND
PINK	15-706-IND	15-716-IND	15-760-IND
ORANGE	15-707-IND	15-717-IND	15-770-IND
BLACK	15-709-IND	15-719-IND	15-790-IND
LIGHT BLUE	15-722-IND	15-732-IND	15-742-IND

Double locking hooks can be inserted in a large guide hole



The **grooves** on the padlock facilitate handling



**INDICATOR
DOTS**

REGULAR


800-SERIES PADLOCKS

The 800 series padlocks have a longer arm (32 mm) and have a single, thinner hook than the 700 series. They are perfect for almost any type of container, even those with long and thin channels.

PRODUCT SELECTION

REGULAR PADLOCKS

COLOR	100 UNITS	200 UNITS	500 UNITS
WHITE	15-801	15-811	15-810
BLUE	15-802	15-812	15-820
GREEN	15-803	15-813	15-830
RED	15-804	15-814	15-840
YELLOW	15-805	15-815	15-850
PINK	15-806	15-816	15-860
ORANGE	15-807	15-817	15-870
BLACK	15-809	15-819	15-890
LIGHT BLUE	15-822	15-832	15-842

WITH INDICATORS DOTS

COLOR	100 UNITS	200 UNITS	500 UNITS
WHITE	15-801-IND	15-811-IND	15-810-IND
BLUE	15-802-IND	15-812-IND	15-820-IND
GREEN	15-803-IND	15-813-IND	15-830-IND
RED	15-804-IND	15-814-IND	15-840-IND
YELLOW	15-805-IND	15-815-IND	15-850-IND
PINK	15-806-IND	15-816-IND	15-860-IND
ORANGE	15-807-IND	15-817-IND	15-870-IND
BLACK	15-809-IND	15-819-IND	15-890-IND
LIGHT BLUE	15-822-IND	15-832-IND	15-842-IND

Single locking hook
can be inserted
in a large guide hole



The locking system
allows for **easy closure**



**INDICATOR
DOTS**

REGULAR


900-SERIES PADLOCKS

The 900 series padlocks are identical to the 800 series, but with a larger rectangular base which allows for the placement of a 20 mm x 30 mm tracking label.

PRODUCT SELECTION

REGULAR PADLOCKS

COLOR	100 UNITS	200 UNITS	500 UNITS
WHITE	15-901	15-911	15-910
BLUE	15-902	15-912	15-920
GREEN	15-903	15-913	15-930
RED	15-904	15-914	15-940
YELLOW	15-905	15-915	15-950
PINK	15-906	15-916	15-960
ORANGE	15-907	15-917	15-970
BLACK	15-909	15-919	15-990
LIGHT BLUE	15-922	15-932	15-942

WITH INDICATORS DOTS

COLOR	100 UNITS	200 UNITS	500 UNITS
WHITE	15-901-IND	15-911-IND	15-910-IND
BLUE	15-902-IND	15-912-IND	15-920-IND
GREEN	15-903-IND	15-913-IND	15-930-IND
RED	15-904-IND	15-914-IND	15-940-IND
YELLOW	15-905-IND	15-915-IND	15-950-IND
PINK	15-906-IND	15-916-IND	15-960-IND
ORANGE	15-907-IND	15-917-IND	15-970-IND
BLACK	15-909-IND	15-919-IND	15-990-IND
LIGHT BLUE	15-922-IND	15-932-IND	15-942-IND

Single locking hook
for use with smaller
diameters




Ideal for adding labels/
indicators indicators or
markings to the large
rectangular base



SCRUB PAD

Our **Scrub Pad** is ideal for pre-cleaning your surgical instruments without leaving any metal residue that often occurs when using steel wool or metal sponges. The abrasiveness is perfectly balanced to provide you with efficient cleaning without damaging your instruments. Its fine fibers allow the user to polish, clean, and debur all types of material surfaces.

- Abrasive Material: Aluminum Oxide
- Grit: 400
- Grade: Very Fine
- Style: Flexible
- Pad material: Non-Woven
- Color: Maroon
- Type: Hand Pad



PRODUCT SELECTION

ITEM #	DESCRIPTION	DIMENSIONS	PACKAGING
60-46050	Scrub Pad	4" x 6" / 10 cm x 15 cm	150 / Box (50 bags of 3 units each)
60-46003	Scrub Pad	4" x 6" / 10 cm x 15 cm	3 / Bag

Your First ™

ENZYMATIC & NON-ENZYMATIC CLEANERS

Our Enzymatic and Non-Enzymatic Cleaners are specially formulated to remove organic material from medical equipment and instrumentation.

Our **Enzymatic Cleaner** formats are formulated with **PrAmLi Power3 Enzymes™**, our proprietary unique blend of enzymes, and a surfactant to efficiently remove all proteins, blood, fatty tissue, and more. They also remove stains from surgical scrubs, linens, and medical clothing.

Our **Non-Enzymatic Cleaner** formats are formulated with a free-rinsing, neutral pH detergent blend for effective cleaning.



Both formats are powerful solutions that break down all bioburden for effective surface and instrument pre-cleaning*

PRODUCT SELECTION

ENZYMATIC CLEANER	NON-ENZYMATIC CLEANER	DILUTION DETERGENT / WATER	PACKAGING	CASE DIMENSIONS
LOW FOAM				
SS-21-1001	SS-721-1001	4 to 8 ml / L (0.5 to 1 oz / Gal)	1 x 3.8 L (1 Gal)	6" x 6" x 13"
SS-21-1004	SS-721-1004	4 to 8 ml / L (0.5 to 1 oz / Gal)	4 x 3.8 L (4 Gal)	13" x 13" x 14"
MILD FOAM				
SS-21-2001	SS-721-2001	4 to 8 ml / L (0.5 to 1 oz / Gal)	1 x 3.8 L (1 Gal)	6" x 6" x 13"
SS-21-2004	SS-721-2004	4 to 8 ml / L (0.5 to 1 oz / Gal)	4 x 3.8 L (4 Gal)	13" x 13" x 14"
SUPER CONCENTRATE				
SS-21-3001	SS-721-3001	1 to 2 ml / L (0.125 to 0.25 oz / Gal)	1 x 3.8 L (1 Gal)	6" x 6" x 13"
SS-21-3004	SS-721-3004	1 to 2 ml / L (0.125 to 0.25 oz / Gal)	4 x 3.8 L (4 Gal)	13" x 13" x 14"

GOOD TO KNOW:

OUR LOW FOAM & SUPER CONCENTRATE DETERGENTS (ENZYMATIC & NON-ENZYMATIC) MAY BE USED IN WASHERS / ULTRASONIC CLEANERS.

WE RECOMMEND CONFIRMING SUGGESTED DILUTION RATIO AND WATER TEMPERATURE WITH THE MANUFACTURER OF YOUR WASHER / ULTRASONIC CLEANER AS MACHINES / SETTINGS VARY.

NOTE: Our Low Foam & Super Concentrate Detergent formulas are available in a 2.5 Gallon format based on a minimum order quantity. Contact us for more information.

*Caution: We recommend consulting your hospital pre-cleaning/cleaning protocols for various equipment/instrument cleaning detergent preferences.



Your First ™

ENVIRO-SACHET™

ENZYMATIC GEL & NON-ENZYMATIC CLEANER

Your First ™ Enviro-Sachet™, our unique cleaning solutions now available in Enzymatic & Non-Enzymatic formats.



- Our Enzymatic Gel format is formulated with the same triple action **PrAmLi Power³ Enzymes™** that you know and trust.
- Our Non-Enzymatic Cleaner format is formulated with a free-rinsing, neutral pH detergent blend.

Both formats are powerful solutions that break down all bioburden for effective surface and instrument pre-cleaning.

PRODUCT SELECTION

ENZYMATIC GEL	NON-ENZYMATIC CLEANER	DESCRIPTION	PACKAGING	CASE DIMENSIONS
17-166	717-166	1 Dry Tubular Sponge & Enviro-Sachet™ set	25 / Box 4 Boxes / Case	21" x 14" x 13.5"
17-266	717-266	1 Dry Flat Sponge & Enviro-Sachet™ set	50 / Box 4 Boxes / Case	21" x 14" x 13.5"
17-177	717-177	1 Dry Tubular Sponge, Enviro-Sachet™ & Stainless Steel Bowl kit	1 Kit	7" x 6" x 3"
17-277	717-277	1 Dry Flat Sponge, Enviro-Sachet™ & Stainless Steel Bowl kit	1 Kit	7" x 6" x 3"
17-EZSAMP	717-EZSAMP	SAMPLE KIT INCLUDES: 3 Dry Flat Sponges with Enviro-Sachet™ sets 1 Dry Tubular Sponge with Enviro-Sachet™ set 1 Reusable Stainless Steel Bowl	1 Kit	7" x 6" x 3"
H-EZ500		Reusable Stainless Steel Bowl 4.3" x 5.3" x 3" / 11 cm x 13.5 cm x 7.5 cm	1 Unit	7" x 6" x 3"



DISTRIBUTED BY:
ADVANCED SURGI-PHARM INC.
www.advancedsurgipharm.com
Tel:1-800-661-5432

Your First ™

FOAMING SPRAY

READY-TO-USE ENZYME FOAMING SPRAY

Our **Foaming Spray** is ideal for instrument pre-cleaning as the thick foam spray enzymatically breaks down all bioburden to allow for effective cleaning and sterilization.

Formulated with **PrAmLi Power³ Enzymes™**, a unique blend of three powerful enzymes along with a surfactant to provide the ultimate in cleaning efficacy.




INSTRUCTIONS



Shake well



Spray immediately after procedure



Keep wet until decontamination process is started

PRODUCT SELECTION

ITEM #	PACKAGING	CASE DIMENSIONS
SS-21-4001	1 x 946 ml (32 oz) Spray Bottle	10" x 4"
SS-21-4012	12 x 946 ml (32 oz) Spray Bottles	18" x 12" x 11"

Important: In addition to our enzymatic formulation, our spray bottle contains as much as 25% more solution when compared to the competition.

Your First ™

F1RST WIPES

Our cleanroom wipes are designed to do the first wipe-down of an endoscope after a procedure. They are made from hydroentangled non-woven material composed of 45% polyester and 55% cellulose. These wipes are manufactured and packaged in a cleanroom environment.



- Can be used for a cleanroom environment of ISO Class 5 through 8
- Production process guarantees extremely low particle generation and chemical extractables
- Compatible with virtually any solvent or solution for cleaning and wiping (such as **PrAmLi Power³ Enzymes™**)
- Double-bag packaging
- Available in two versatile sizes

PRODUCT SELECTION

ITEM #	DESCRIPTION	DIMENSIONS	PACKAGING	CASE DIMENSIONS
17-A0909	Non-woven Cleanroom Wipes	9" x 9" (23 cm x 23 cm)	300 / Bag 14 Bags / Case	18.5" x 10" x 25.25"
17-A1212	Non-woven Cleanroom Wipes	12" x 12" (31 cm x 31 cm)	150 / Bag 10 Bags / Case	13" x 13" x 19"

Your First ™

DRY SPONGES

Our sponges are specifically fabricated for use with all medical instrumentation. Unlike most competitor sponges, our sponges are white to better evidence the effectiveness of the cleaning process. Dark colored sponges easily hide contaminants, but with our white sponges you can visually evaluate the pre-clean efficacy.



Available in two formats, Flat and Tubular:

Will not flake or scratch expensive equipment

Eco-friendly

100% medical grade high-density foam

Compatible with enzymatic and disinfectant solutions

TIP: OUR TUBULAR SPONGE CAN ALSO BE USED AS A BUMPER TO PROTECT YOUR SCOPES

High absorbency retains more solution than conventional sponges or paper wipes

PRODUCT SELECTION


ITEM #	DESCRIPTION	DIMENSIONS	PACKAGING
17-100	DRY TUBULAR SPONGE White Dispenser Box	2.0" ø x 4.25" (2.7 cm ø x 10.8 cm)	25 / Box (4 Boxes / Case)
17-101	DRY TUBULAR SPONGE Bulk	2.0" ø x 4.25" (2.7 cm ø x 10.8 cm)	560 / Case
17-200	DRY FLAT SPONGE White Dispenser Box	4" x 5" x 0.5" (10.2 cm x 12.7 cm x 1.3 cm)	50 / Box (4 Boxes / Case)
17-201	DRY FLAT SPONGE Bulk	4" x 5" x 0.5" (10.2 cm x 12.7 cm x 1.3 cm)	1,000 / Case

Your First ™

ENZYMATIC BEDSIDE SPONGE KIT

RightSwipe™

Introducing **RightSwipe™**, our newly launched Enzymatic Bedside Sponge Kits now feature sponges treated with an enhanced formula for a more effective pre-clean or surface clean solution. Our improved three-enzyme blend ensures thorough and safe cleaning and a simplified application process, making **RightSwipe™** the perfect match for your endoscope care needs. The tubular sponge design is compatible with diverse endoscope models, providing versatility alongside the superior benefits of our enzymatic solution.



FLAT SPONGE



TUBULAR SPONGE



FREESTANDING POUCH



DISPENSER BOX

PRODUCT SELECTION

ITEM #	DESCRIPTION	PACKAGING	CASE DIMENSIONS
17-111	PRE-TREATED TUBULAR SPONGE Add 100-500 ml of water to activate detergent solution. Comes in plastic bag ¹	25 / Box (4 Boxes / Case)	21" x 14" x 13.5"
17-211	PRE-TREATED FLAT SPONGE Add 100-500 ml of water to activate detergent solution. Comes in plastic bag ¹	50 / Box (4 Boxes / Case)	21" X 14" X 13.5"
17-122	PRE-TREATED TUBULAR SPONGE Comes in a freestanding pouch to use as a basin. Simply add 500 ml of water to the pouch.	15 / Box (4 Boxes / Case)	21" X 14" X 13.5"
17-222	PRE-TREATED FLAT SPONGE Comes in a freestanding pouch to use as a basin. Simply add 500 ml of water to the pouch.	25 / Box (4 Boxes / Case)	21" X 14" X 13.5"
17-133	READY TO USE TUBULAR SPONGE Comes pre-treated in 500 ml of an enzymatic solution in freestanding pouch.	24 / Case	14" X 12" X 9"
17-233	READY TO USE FLAT SPONGE Comes pre-treated in 500 ml of an enzymatic solution in freestanding pouch.	24 / Case	14" X 12" X 9"
17-144	READY TO USE DUAL PURPOSE TUBULAR SPONGE Comes pre-treated in 100 ml of an enzymatic solution in freestanding pouch.	12 / Box (4 Boxes / Case)	21" X 14" X 13.5"
17-244	READY TO USE DUAL PURPOSE FLAT SPONGE Comes pre-treated in 100 ml of an enzymatic solution in freestanding pouch.	12 / Box (4 Boxes / Case)	21" X 14" X 13.5"
17-155	READY TO USE TUBULAR SPONGE Comes pre-treated in 220 ml of an enzymatic solution in freestanding pouch.	24 / Case	14" x 12" x 9"
17-255	READY TO USE FLAT SPONGE Comes pre-treated in 220 ml of an enzymatic solution in freestanding pouch.	24 / Case	14" X 12" X 9"

¹Our stainless steel bowl can be sold separately and used as a basin for easy and convenient sponge activation.

Your First ™

NON-ENZYMATIC BEDSIDE SPONGE KIT

RightSwipe™

Presenting our enhanced Non-Enzymatic Bedside Sponge Kits. Each kit includes sponges pre-treated with a free-rinsing, neutral pH detergent blend from **RightSwipe™**, ensuring an effective pre-cleaning and surface cleaning solution for medical instruments. This kit stands as Your First™ step in the medical instrument cleaning process, showcasing the advanced benefits for superior instrument care. Our enhanced bedside kits will provide you with a streamlined workflow.



FLAT SPONGE



TUBULAR SPONGE



FREESTANDING POUCH



DISPENSER BOX

PRODUCT SELECTION

ITEM #	DESCRIPTION	PACKAGING	CASE DIMENSIONS
717-122	PRE-TREATED TUBULAR SPONGE Comes in a freestanding pouch to use as a basin. Simply add 500 ml of water to the pouch. ¹	15 / Box (4 Boxes / Case)	21" X 14" X 13.5"
717-222	PRE-TREATED FLAT SPONGE Comes in a freestanding pouch to use as a basin. Simply add 500 ml of water to the pouch. ¹	25 / Box (4 Boxes / Case)	21" X 14" X 13.5"
717-133	READY TO USE TUBULAR SPONGE Comes pre-treated in 500 ml of non-enzymatic solution in freestanding pouch.	24 / Case	14" X 12" X 9"
717-233	READY TO USE FLAT SPONGE Comes pre-treated in 500 ml of non-enzymatic solution in freestanding pouch.	24 / Case	14" X 12" X 9"
717-144	READY TO USE DUAL PURPOSE TUBULAR SPONGE Comes pre-treated in 100 ml of non-enzymatic solution in freestanding pouch.	12 / Box (4 Boxes / Case)	21" X 14" X 13.5"
717-244	READY TO USE DUAL PURPOSE FLAT SPONGE Comes pre-treated in 100 ml of non-enzymatic solution in freestanding pouch.	12 / Box (4 Boxes / Case)	21" X 14" X 13.5"
717-155	READY TO USE TUBULAR SPONGE Comes pre-treated in 220 ml of non-enzymatic solution in freestanding pouch.	24 / Case	14" x 12" x 9"
717-255	READY TO USE FLAT SPONGE Comes pre-treated in 220 ml of non-enzymatic solution in freestanding pouch.	24 / Case	14" X 12" X 9"

¹Our stainless steel bowl can be sold separately and used as a basin for easy and convenient sponge activation.



Your First ™

TRANSPORT BAGS FOR ENDOSCOPE AND MEDICAL INSTRUMENTS

Thoughtfully designed to eliminate touch-points during endoscope transport from point-of-care to the reprocessing area, and back, focused on safety of patients and medical staff.

- Made of recyclable LDPE material known to be quite flexible, highly resilient, and tough
- Latex and DEHP-free
- Excellent resistance to dilute and concentrated acids, alcohols, bases and esters*
- Good resistance to aldehydes, ketones, and vegetable-based oils*
- Can withstand high temperatures of up to 90 °C (194°F)
- Tamper-evident seal for added safety, also aids in keeping the instrument moist to avoid drying of bio-debris
- Clearly marked bags with color-coded visuals facilitates immediate identification of clean vs. used instruments
- Users can write important information directly on the bag in spaces marked with visual icons
- Available in two convenient sizes

GOOD TO KNOW USING Your First ™ SEALABLE TRANSPORT BAGS REDUCES USER ERROR, PREVENTS CROSS-CONTAMINATION, AND PROMOTES PATIENT AND MEDICAL STAFF SAFETY.

CLEAN INSTRUMENTS BAG



SOILED INSTRUMENTS BAG



MADE IN CANADA

PRODUCT SELECTION

ITEM #	DESCRIPTION	DIMENSIONS	PACKAGING	CASE DIMENSIONS
17-HAP24	Clean Instruments Bag	24" x 22" / 61.0 cm x 56.0 cm	200 / Case	18" x 14" x 5"
17-SAD24	Soiled Instruments Bag	24" x 22" / 61.0 cm x 56.0 cm	200 / Case	18" x 14" x 5"
17-HAP30	Clean Instruments Bag	30.5" x 26" / 77.0 cm x 66.0 cm	200 / Case	21" x 19" x 7"
17-SAD30	Soiled Instruments Bag	30.5" x 26" / 77.0 cm x 66.0 cm	200 / Case	21" x 19" x 7"

*Reference on file



DISTRIBUTED BY:
ADVANCED SURGI-PHARM INC.
www.advancedsurgipharm.com
Tel: 1-800-661-5432

SAVEFACE™ DRAPED SHIELDS

REDUCE PERSONAL INJURY
SUPERIOR INSTRUMENT DECONTAMINATION PROTECTION

SAVEFACE DRAPED-FACE SHIELD

In the Sterile Processing Department (SPD) Central Sterile Supply Department (CSSD), personal safety is paramount.

Our Draped Face Shield provides essential protection during high-risk instrument decontamination.

EFFECTIVE PROTECTION IS NON-NEGOTIABLE

SAVEFACE Draped dynamically reduces harmful exposure to caustic chemicals as well as the risk of cross-contamination. Effective face protection is non-negotiable.

PPE CATEGORY 2 CERTIFIED

SAVEFACE DRAPED is tested to protect. Wear it with confidence.



PATENT PENDING



INCLUDES CLOSING SYSTEM TO REDUCE RISK OF FLUID CONTACT FROM BELOW THE SHIELD

SAVEFACE: KEY FEATURES

ANSI/AAMI PB70 LEVEL 3 COMPLIANT DRAPE:

- Provides a high level of fluid penetration protection.
- Tuck into clothing or pinch together.

EXTRA-LONG 9"/23 CM WRAP-AROUND SHIELD:

- Ensures full coverage, reducing the risk of exposure from different angles.

COMFORTABLE:

- The soft, self-shaping foam headpiece molds comfortably to your forehead, allowing for extended wear without discomfort.

OPTICAL CLARITY WITH ANTI-FOG AND ANTI-GLARE

DON'T MISS THE OPPORTUNITY TO EXPLORE OUR BROADER RANGE OF PROTECTIVE SHIELDS.

Visit our website www.surgmed.com or consult our detailed literature to find the shield that best fits your needs.

VERY IMPORTANT

IS YOUR PRESENT FACE SHIELD CATEGORY 2 CERTIFIED?

Batrik face shields are independently certified to comply with the PPE Category 2 standard. Many facilities are now requesting compliance to this standard.

SHIELD LENGTH
9" / 23 CM

DRAPE LENGTH
6" / 15 CM



SF-8240

SF-240



PRODUCT SELECTION

ITEM #	DESCRIPTION	DIMENSIONS	TOTAL LENGTH	PACKAGING	CASE DIMENSIONS
SF-240	BATRIK SAVEFACE™ no drape	Shield Length: 9" / 23 cm Shield width: 13' / 33 cm	9" (23 cm)	240 / Case (8 Boxes of 30 Pcs)	26.5" x 21" x 15"
SF-8240	BATRIK SAVEFACE™ with ANSI/AAMI PB70 Level 3 compliant protection drape	Shield Length: 9" / 23 cm Shield width: 13' / 33 cm Shield width: 6" / 15 cm	23" (58 cm)	240 / Case (8 Boxes of 30 Pcs)	26.5" x 21" x 15"


SHOOZCOVERS™ | SOFT LINE™

Batrik ShoozCovers™ Soft Line™ is manufactured in Canada from high quality **35 g/m² polypropylene fabric and 100% latex free elastics** for performance and durability you can rely on. Our anti-skid tread reduces the risk of slipping and our smart design delivers easy on-and-off convenience and great fit.

MADE IN CANADA



SURGICAL PRODUCTS & OPERATING ROOM



SHOOZ-REPELLENT™ TECHNOLOGY:

The **Soft Line™** integrates Shooz-Repellent™ technology which creates a protective barrier that repels fluids, delivering basic incidental fluid contact repellent. Shoe covers are not waterproof or water-resistant.



ANTI-SKID TREAD:

Our anti-skid tread reduces the risk of slipping with **triple ultrasonic welds** that increases tear strength and delivers easy on-and-off convenience.

OUR SOFT LINE™ IS PERFECT FOR:



MEDICAL ESTABLISHMENTS



PHARMACEUTICAL LABS



GMP SITES



ELECTRONIC MANUFACTURING

DISTRIBUTED BY:
ADVANCED SURGI-PHARM INC.
www.advancedsurgipharm.com
Tel.1-800-661-5432



WHY YOU'LL LOVE THE SOFT LINE™:

- Both shoe covers are manufactured in Canada with top-quality 35 g/m² polypropylene fabric and 100% latex-free USA-made elastics¹.
- This line is particle count tested and complies with IEST-RP-CC003.4.
- Other colors available upon request.

SUPPORT A CAUSE, WEAR PINK:

- A portion of our Bubble Gum Pink shoe cover sales will be donated to charity to help further breast cancer research. Join us in making a meaningful contribution by choosing Pink ShoozCovers™!



AVAILABLE IN TWO VERSATILE SIZES



AVAILABLE IN UNIVERSAL

PRODUCT SELECTION:

ITEM#	DESCRIPTION	COLOR	SIZE	PACKAGING	CASE DIMENSIONS
SC-1210-PB	ShoozCovers™ - Soft Line™ <i>polypropylene non-woven fabric</i>	Lolly Blue	Universal Up to shoe size 12 ~ 42 cm x 16.5 cm / 16.5" x 6.5"	Case of 1,000 (10 bags of 100 units)	20" x 16" x 16"
SC-1216-PDB	ShoozCovers™ - Soft Line™ <i>polypropylene non-woven fabric</i> <i>*Covers are available with or without skid-resistant tread.</i>	Lolly Blue	Universal Up to shoe size 12 ~ 42 cm x 16.5 cm / 16.5" x 6.5"	Case of 600 (6 Dispenser Boxes of 100 units)	16.25" x 8.875" x 18.75"
SC-1216-PB	ShoozCovers™ - Soft Line™ <i>polypropylene non-woven fabric</i> <i>*Covers are available with or without skid-resistant tread.</i>	Lolly Blue	Universal Up to shoe size 12 ~ 42 cm x 16.5 cm / 16.5" x 6.5"	Case of 600 (6 bags of 100 units)	19.5" x 19.5" x 14.25"
SC-1226-PDB	ShoozCovers™ - Soft Line™ <i>polypropylene non-woven fabric</i> <i>*Covers are available with or without skid-resistant tread.</i>	Lolly Blue	XL Shoe size 12-15 ~ 47 cm x 19 cm / 18.5" x 7.5"	Case of 600 (6 Dispenser Boxes of 100 units)	17.25" x 9.375" x 21.75"
SC-1220-PB	ShoozCovers™ - Soft Line™ <i>polypropylene non-woven fabric</i>	Lolly Blue	XL Shoe size 12-15 ~ 47 cm x 19 cm / 18.5" x 7.5"	Case of 1,000 (10 bags of 100 units)	22" x 16" x 16"
SC-1616-PDB	ShoozCovers™ - Soft Line™ <i>polypropylene non-woven fabric</i>	Bubble Gum Pink	Universal Up to shoe size 12 ~ 42 cm x 16.5 cm / 16.5" x 6.5"	Case of 600 (6 Dispenser Boxes of 100 units)	16.25" x 8.875" x 18.75"
SC-1616-PB	ShoozCovers™ - Soft Line™ <i>polypropylene non-woven fabric</i>	Bubble Gum Pink	Universal Up to shoe size 12 ~ 42 cm x 16.5 cm / 16.5" x 6.5"	Case of 600 (6 bags of 100 units)	19.5" x 19.5" x 14.25"

¹Many suppliers claim their products are latex free, however this often only relates to the fabric and not the elastics used. We source our ShoozCovers™ elastics from the USA to ensure that our components conform fully as a 100% latex free product.

PALLET INFORMATION:

PRODUCT	CASES PER PALLET*	PALLET DIMENSIONS	PALLET WEIGHT	SIZE	PACKAGING UOM
Soft Line™	40	40 x 48 x 48 in	420 lbs	U	6 Dispenser Box / Case (100 un. each)

*Quantity of Cases per Pallet may vary



GrippezShooz™

EXTRA-LARGE FOR EXTRA-TOUGH WORK:

GrippezShooz™ provide exceptional protection and **enhanced traction** for professionals working in the most demanding, slippery environments. Featuring waterproof protection, our shoe covers ensure safety and hygiene without compromising on quality. They are the ideal solution for tough jobs.



MADE IN CANADA



PERFECT FOR:



AEROSPACE



FOOD PROCESSING
INDUSTRY



OIL, GAS AND
MINING SITES



CONSTRUCTION AND
RENOVATION SITES

¹ INDUSTRIAL STRENGTH ELASTICS

GrippezShooz™ ARE A PERFECT FIT FOR EVERYONE:

- Fast, Reliable Fit: Stretch material ensures a secure, easy-on fit. Coating provides protection against slippery surfaces.
- Comes in an extra-large size that fits over a men's size 15 steel-toed boot.
- Enhanced Safety: Anti-slip rubberized coating provides reliable grip on all surfaces, helping to prevent sliding in critical areas.
- Available in two colors: Gecko Navy and Rhino Grey.

RHINO GREY



GECKO NAVY



PRODUCT SELECTION:

ITEM#	DESCRIPTION	COLOR	SIZE	PACKAGING	CASE DIMENSIONS
SC-29312-NPDB	GrippezShooz™ Polypropylene + polyethylene non-woven fabric, waterproof	Rhino Grey	Universal Up to shoe size 12 ~ 42 cm x 17.5 cm / 16.5" x 6.9"	Case of 300 (2 dispenser boxes of 150 units)	16.62" x 14.875" x 8.75"
SC-29312-NPB	GrippezShooz™ Polypropylene + polyethylene non-woven fabric, waterproof	Rhino Grey	Universal Up to shoe size 12 ~ 42 cm x 17.5 cm / 16.5" x 6.9"	Case of 300 (2 bags of 150 units)	16" x 13.375" x 10"
SC-29322-NPDB	GrippezShooz™ Polypropylene + polyethylene non-woven fabric, waterproof	Rhino Grey	XL Shoe size 12-15 ~ 47 cm x 17.5 cm / 18.5" x 6.9"	Case of 240 (2 dispenser boxes of 120 units)	16.62" x 14.875" x 8.75"
SC-29322-NPB	GrippezShooz™ Polypropylene + polyethylene non-woven fabric, waterproof	Rhino Grey	XL Shoe size 12-15 ~ 47 cm x 17.5 cm / 18.5" x 6.9"	Case of 240 (2 bags of 120 units)	16" x 13.375" x 10"
SC-22212-NPDB	GrippezShooz™ Polypropylene + polyethylene non-woven fabric, waterproof	Gecko Navy	Universal Up to shoe size 12 ~ 42 cm x 17.5 cm / 16.5" x 6.9"	Case of 300 (2 dispenser boxes of 150 units)	16.62" x 14.875" x 8.75"
SC-22212-NPB	GrippezShooz™ Polypropylene + polyethylene non-woven fabric, waterproof	Gecko Navy	Universal Up to shoe size 12 ~ 42 cm x 17.5 cm / 16.5" x 6.9"	Case of 300 (2 bags of 150 units)	16" x 13.375" x 10"
SC-22222-NPDB	GrippezShooz™ Polypropylene + polyethylene non-woven fabric, waterproof	Gecko Navy	XL Shoe size 12-15 ~ 47 cm x 17.5 cm / 18.5" x 6.9"	Case of 240 (2 dispenser boxes of 120 units)	16.62" x 14.875" x 8.75"
SC-22222-NPB	GrippezShooz™ Polypropylene + polyethylene non-woven fabric, waterproof	Gecko Navy	XL Shoe size 12-15 ~ 47 cm x 17.5 cm / 18.5" x 6.9"	Case of 240 (2 bags of 120 units)	16" x 13.375" x 10"

PALLET INFORMATION:

PRODUCT	CASES PER PALLET*	PALLET DIMENSIONS	PALLET WEIGHT	SIZE	PACKAGING UOM
GrippezShooz™	50	40 x 48 x 88 in	446-479 lbs	U	2 Dispenser Box / Case (150 un. each)

*Quantity of Cases per Pallet may vary


SHOOZCOVERS™ | COATED LINE™

STAY DRY:

Our **Batrik ShoozCovers™ Coated Line™** is made from a **~55 g/m² blend of polypropylene and polyethylene**, providing high-quality water resistance.¹ Ideal for environments where maximum protection against fluids is required.



MADE IN CANADA



TRIPLE WELDED TECHNOLOGY:

These shoe covers are reinforced with triple ultrasonic welds, a technique that joins materials together without the need for adhesives or stitching. This technique increase tear strength and reduce linting, for long lasting, resilient covers.

OUR COATED LINE™ IS PERFECT FOR:



MEDICAL ESTABLISHMENTS



PHARMACEUTICAL LABS



GMP SITES



ELECTRONIC MANUFACTURING

¹OUR SHOOZCOVERS™ COATED LINE™ IS RESISTANT AGAINST ALL NON PETROLEUM BASED LIQUIDS.

MARSHMALLOW WHITE | COATED LINE™



WHY YOU'LL LOVE THE COATED LINE™:

- Superior protection using our unique blend of fabrics which protects against not only water but all non-petroleum-based liquids, making them the ideal choice for environments where spills and leaks are a concern
- Elastics around the ankle create a strong and secure fit which is especially important when working in hazardous environments where a misplaced shoe cover could lead to contamination or injury
- Smart design with a large opening, providing a better fit and easy on-and-off convenience
- Other colors available upon request


UN - XL

AVAILABLE IN TWO VERSATILE SIZES

PRODUCT SELECTION:

ITEM#	DESCRIPTION	COLOR	SIZE	PACKAGING	CASE DIMENSIONS
SC-2112-NPB	ShoozCovers™ - Coated Line™ <i>Polypropylene + polyethylene non-woven fabric, waterproof</i>	Marshmallow White	Universal Up to shoe size 12 ~ 42 cm x 17.5 cm / 16.5" x 6.9"	Case of 300 (2 bags of 150 units)	16" x 13.375" x 10"
SC-2112-NPDB	ShoozCovers™ - Coated Line™ <i>Polypropylene + polyethylene non-woven fabric, waterproof</i>	Marshmallow White	Universal Up to shoe size 12 ~ 42 cm x 17.5 cm / 16.5" x 6.9"	150 / Dispenser Box 2 Dispenser Boxes / Case	16.62" x 14.875" x 8.75"
SC-2122-NPB	ShoozCovers™ - Coated Line™ <i>Polypropylene + polyethylene non-woven fabric, waterproof</i>	Marshmallow White	XL Shoe size 12-15 ~ 47 cm x 17.5 cm / 18.5" x 6.9"	Case of 240 (2 bags of 120 units)	16" x 13.375" x 10"
SC-2122-NPDB	ShoozCovers™ - Coated Line™ <i>Polypropylene + polyethylene non-woven fabric, waterproof</i>	Marshmallow White	XL Shoe size 12-15 ~ 47 cm x 17.5 cm / 18.5" x 6.9"	120 / Dispenser Box 2 Dispenser Boxes / Case	16.62" x 14.875" x 8.75"

PALLET INFORMATION:

PRODUCT	CASES PER PALLET*	PALLET DIMENSIONS	PALLET WEIGHT	SIZE	PACKAGING UOM
Coated Line™	72	40 x 48 x 90 in	481 lbs	XL	2 Dispenser Box / Case (120 un. each)

*Quantity of Cases per Pallet may vary

SHOOZCOVERS™ | TOUGH LINE™



MADE IN CANADA 



EXTRA-LARGE FOR EXTRA-TOUGH WORK:

SHOE SIZES 12-15 (45-49.5 EU)
~ 47 CM X 19 CM / 18.5" X 7.5"

Our **Batrik ShoozCovers™ Tough Line™** comes in an extra-large size that fits over a **men's size 15 steel-toed boot**, delivering optimal full coverage and protection while you work.¹



ENHANCED DURABILITY:

These shoe covers are formulated with our **thickest fabric yet!** Made of **~55 g/m² polypropylene fabric**, these highly durable covers provide added protection, increase tear strength and reduce linting.

OUR TOUGH LINE™ IS PERFECT FOR:



HOMES OR
COMMERCIAL
BUILDINGS



SHIELDING AGAINST
MUD & GREASE



MEDICAL
ESTABLISHMENTS




PHARMACEUTICAL
LABS

¹ SHOOZCOVERSTM ARE MADE LONGER (42 & 47CM) FOR ACCURATE SHOE SIZE SELECTION.



WHY YOU'LL LOVE THE TOUGH LINE™:

- Triple ultrasonic welds that are long lasting to reduce waste
- High quality made in the USA latex-free elastics are triple reinforced to contract the shoe cover around your foot and deliver a strong and secure fit²
- Shooz-Repellent™ technology delivers basic incidental fluid contact repellency³
- Other colors available upon request



XL

AVAILABLE IN EXTRA LARGE

ANY TIME, ANYWHERE:



RESIDENTIAL



COMMERCIAL



INDUSTRIAL

PRODUCT SELECTION:

ITEM#	DESCRIPTION	COLOR	SIZE	PACKAGING	CASE DIMENSIONS
SC-1922-PB	ShoozCovers™ - Tough Line™ <i>Polypropylene non-woven fabric</i>	Black Liquorice	XL Shoe size 12-15 ~ 47 cm x 19 cm / 18.5" x 7.5"	Case of 200 (2 bags of 100 units)	16" x 13.375" x 10"

ALTERNATIVE PACKAGING: Add "-PDB" to end of code for a dispenser box instead of bags.

PALLET INFORMATION:

PRODUCT	CASES PER PALLET*	PALLET DIMENSIONS	PALLET WEIGHT	SIZE	PACKAGING UOM
Tough Line™	72	40 x 48 x 90 in	481 lbs	XL	2 Dispenser Box / Case (120 un. each)

*Quantity of Cases per Pallet may vary

² Many suppliers claim their products are latex free, however this often only relates to the fabric and not the elastics used. We source our ShoozCovers™ elastics from the USA to ensure that our components conform fully as a 100% latex free product. ³ Shooz-Repellent™ technology provides basic incidental fluid contact repellent. It is not waterproof nor water resistant.

ANTIMICROBIAL IV POLES WITH CUVERRO SHIELD™

Frequently touched surfaces like IV poles have long been identified as reservoirs for the spread of pathogenic microbes.¹ These can easily contaminate the hands and equipment of healthcare professionals who in turn can transmit these pathogens to patients during routine care leading to Healthcare-associated infections.

Our Antimicrobial IV Poles protected with CuVerro Shield™ by Aereus Technologies* offers you a solution to this problem.


CuVerro Shield™ is a unique coating of naturally antimicrobial copper that is applied to our stainless steel poles, ram horns, and u-hooks to deliver outstanding protection to your highest touch zones.

The result? CuVerro® kills more than 99.9% of bacteria** within two hours and continues to kill 99% of bacteria** even after repeated contamination.

A US government funded study conducted by researchers at the US National Institutes of Health and the Centers for Disease Control and Protection found that the HCoV-19 virus remained viable for up to 2 to 3 days on plastic and stainless steel surfaces versus up to 4 hours on copper.²



OR-36770-CUV



Viability of SARS-CoV-1 and SARS-CoV-2 in Aerosols and on Various Surfaces.²

* CuVerro Shield™ by Aereus Technologies is a Health Canada Pest Management Regulatory Agency (PMRA) registered process.

**Laboratory testing shows that, when cleaned regularly, CuVerro surfaces kill greater than 99.9% of the following bacteria within 2 hours of exposure: Methicillin-Resistant Staphylococcus aureus, Staphylococcus aureus, Enterobacter aerogenes, Pseudomonas aeruginosa, E.coli O157:H7, and Vancomycin-Resistant Enterococcus faecalis (VRE). The use of CuVerro® bactericidal copper products is a supplement to and not a substitute for standard infection control practices; users must continue to follow all current infection control practices, including those practices related to cleaning and disinfection of environmental surfaces. This surface has been shown to reduce microbial contamination, but it does not necessarily prevent cross contamination.

¹Oie S, Hosokawa I, Kamiya A. 2002. Contamination of room door handles by methicillin sensitive/methicillin-resistant Staphylococcus aureus. J. Hosp. Infect. 2002;51:140-143.

²Doremalen NV, Bushmaker T, Morris DH, et al. Aerosol and surface stability of HCoV-19 (SARS-CoV-2) compared to SARS-CoV-1. N Engl J Med 2020.



BASES & CASTERS

ITEM #	DESCRIPTION
OR-36770-CUV	CuVERRO SHIELD™ IV POLE, CENTRE COLLAR, AND KNOB with 23" cream aluminium base (V-40016) and 2" casters
OR-36760-25-CUV	CuVERRO SHIELD™ IV POLE, CENTRE COLLAR, AND KNOB with 25" stainless steel base (V-40013-25) and 2" casters
OR-36760-253-CUV	with 25" stainless steel base (V-40013-25) and 3" casters
OR-36771-25-CUV	CuVERRO SHIELD™ IV POLE, CENTRE COLLAR, AND KNOB with 25" stainless steel base (V-40013-25), 2" casters, and weight kit of 12.5 lbs. (V-50004)
OR-36771-253-CUV	with 25" stainless steel base (V-40013-25), 3" casters, and weight kit of 12.5 lbs. (V-50004)
OR-36772-253-CUV	with 25" steel base (V-40014-25-ST) and 3" casters

SMOOTH-ROLLING
NON-MARKING
DURABLE CASTERS




OPTIONAL WEIGHT KIT
AVAILABLE FOR INCREASED
STABILITY

HOOKS & RAM HORNS

ITEM #	DESCRIPTION
OR-36499-CUV	CuVERRO SHIELD™ STAINLESS STEEL RAM HORNS, AND ENDCAP Set of 2
OR-36500-CUV	CuVERRO SHIELD™ STAINLESS STEEL RAM HORNS, AND ENDCAP Set of 4
OR-36501-CUV	CuVERRO SHIELD™ STAINLESS STEEL U-HOOKS, AND ENDCAP Set of 2
OR-36502-CUV	CuVERRO SHIELD™ STAINLESS STEEL U-HOOKS, AND ENDCAP Set of 4
OR-36504-CUV	CuVERRO SHIELD™ STAINLESS STEEL RAM HORNS (2), STAINLESS STEEL U-HOOK (2), AND ENDCAP Set of 4
V-50002-CUV	CuVERRO SHIELD™ STAINLESS STEEL URINARY DRAIN HOOK AND, COLLAR WITH KNOB Set of 1

CUSTOM HOOK DESIGNS
AVAILABLE TO SUIT YOUR
SPECIFIC NEEDS




HEAVY-DUTY RAM HORNS,
U-HOOKS OR A COMBINATION
OF BOTH AVAILABLE

POLES & HANDLES

ITEM #	DESCRIPTION
V-CUV-50003	CuVERRO SHIELD™ ROUND STEEL HANDLE
V-CUV-50006	CuVERRO SHIELD™ STAINLESS STEEL TRIANGULAR HANDLE
V-CUV-50007	CuVERRO SHIELD™ SEMI-ROUND STEEL HANDLE

THE LARGE PRONG KNOB FOR
HEIGHT ADJUSTMENT IS EASY TO
TURN EVEN WITH WET HANDS



GOOD TO KNOW: In the unlikely event your CuVerro coated IV Pole should become deeply scratched, a study conducted by the University of Waterloo demonstrated that antimicrobial activity of the Aereus coating is maintained even when completely removed up to a width of between 3 and 5 mm, under the test conditions used in this study.³

³ Tanvir S, Anderson WA. 2020. Antimicrobial Activity of "Deeply Scratched" Aereus-Coated Surfaces. Department of Medical Engineering. University of Waterloo. 2020.




PURCHASE DIRECTLY FROM OUR SISTER COMPANY
IMPERIAL SURGICAL LTD.
www.imperialsurgical.com
Tel:1-800-661-5432

PLEXIGLASS PROTECTION PANELS

Our Plexiglass Protection Panels are made of 100% acrylic clear plastic to safeguard face-to-face interaction. They are available in a variety of ready-made standard sizes plus we can create custom sizes/designs to accommodate your specific needs. Our simple and ergonomic design provides the necessary protection for your patients, customers, and community without an imposing presence.

FEATURES

- Multiple dimensions and opening sizes available
- 100% acrylic clear plastic
- Easy to install and move as needed
- Easy to disinfect



PPE - REUSABLE EQUIPMENT

PRODUCT SELECTION


ITEM #	DIMENSIONS	ACRYLIC THICKNESS	OPENING SIZE
PG-2331	23.5" W x 31.5" H	3/16"	12" W x 4" H
PG-2336	23.5" W x 36" H	1/4"	8" W x 2" H
PG-2432	24" W x 32" H	3/16"	12" W x 4" H
PG-2436	24" W x 36" H	3/16"	12" W x 4" H
PG-3531	35.5" W x 31.5" H	3/16"	12" W x 4" H
PG-3632	36" W x 32" H	3/16"	5" W x 10" H
PG-3636	36" W x 36" H	3/16"	5" W x 10" H
PG-3735	37.5" W x 35.5" H	1/4"	8" W x 2" H
PG-4831	48" W x 31.5" H	3/16"	12" W x 4" H

STAINLESS STEEL HANDS-FREE DOOR OPENER

The Imperial Surgical Ltd. Hands-Free Door Opener is a 100% stainless steel, foot operated door opener that works on any latch-less door.

Easy to install and easy to clean, the no-touch opener helps minimize the spread of germs by reducing the need to touch door handles where bacteria and viruses may live.

The hands-free door opener comes with an eye catching sticker to put on the door to encourage people to use it. We also provide the hardware so you have everything you need for installation.



MADE IN CANADA 

PRODUCT SELECTION

ITEM #	DIMENSIONS	PACKAGING
11656-R0	5" Length x 1.5" Height x 0.14" Thickness (127 mm Length x 38 mm Height x 3.5 mm Thickness)	1 Unit




Ver.2021-04-22

naso-gastric

Feeding Tubes

Feeding tubes in DEHP-free PVC (for short-term use) and polyurethane (for medium-term use) with flexible female mount and cap.

- Smoothly rounded, closed distal tip and opposed lateral eyes
- Radiopaque
- Available in various lengths: Centimeter graduations (5cm to 25cm)



Standard tube

PVC tubes

Size	04 Fr	05 Fr	06 Fr	07 Fr	08 Fr	09 Fr	10 Fr	12 Fr	Units/Case
40 cm	310.04	310.05	310.06	310.07	-	-	-	-	50
50 cm	-	-	-	-	310.08	310.09	310.10	310.12	50
75 cm	311.04	311.05	311.06	-	311.08	-	311.10	-	50
125 cm	312.04	-	312.06	312.07	312.08	312.09	312.10	312.12	25

PUR tubes

Size	04 Fr	05 Fr	06 Fr	08 Fr	Units/Case
40 cm	1310.04	1310.05	1310.06	-	25
50 cm	-	-	-	1310.08	25
75 cm	-	1311.05	1311.06	1311.08	25

dual-flowgastric

Feeding Tubes

The dual-flow gastric tubes are used to prevent the accumulation of gastric secretions and air, thus allowing decompression of the stomach by continuous aspiration without trauma to the gastric mucosa.



Dual-flow tube

Dual-flow gastric tubes

Product Code		Size	Length	Number of eyes	Units/Case
PVC	PUR				
340.06	1340.06	6 Fr	60 cm	6	10
340.08	1340.08	8 Fr	60 cm	6	10
340.10	1340.10	10 Fr	90 cm	10	10
340.12	1340.12	12 Fr	90 cm	10	10
340.14	-	14 Fr	120 cm	10	10
340.16	1340.16	16 Fr	120 cm	10	10
340.18	1340.18	18 Fr	120 cm	10	10

gastro-duodenal

Feeding Tubes

125 cm tubes for placement in the duodenum.

Weighted tubes: Silicone (for both medium and long-term use) tube with weighted metal tip to facilitate placement

Levin type tubes: DEHP-free PVC (for short-term use) tubes for gastric decompression or gavage feeding.



Silicone tube



Levin type tubes with standard connector

Weighted tubes (silicone)

Product Code	Size	Units/Case
2395.05	5 Fr	25
2395.06	6 Fr	25
2395.09	9 Fr	25

Levin type tubes with standard connector (PVC)

Product Code	Size	Units/Case
391.10	10 Fr	50
391.12	12 Fr	50
391.14	14 Fr	50
391.16	16 Fr	50

nutriSafe2®

A complete range dedicated to safe enteral nutrition. Designed to be fail-safe, the nutriSafe2® system incorporates a unique connection that locks securely so it won't slip apart and yet cannot be connected to a standard 6% tapered luer fitting.

Safe - The risk of misconnection errors is significantly reduced for patient safety.

Secure - Accidental disconnection is nearly impossible, ensuring precious nutrients are received.

Simple - It does not change the technique of the users and therefore does not require any special training.

nutriSafe2® Feeding tubes

nutriSafe2® feeding tubes are available in a range of sizes, lengths and materials to suit patient needs. Numbered graduations every centimeter aid placement and monitoring, and the rounded tip provides atraumatic insertion.


PVC (DEHP-free): for use up to 5 days.

Polyurethane (PUR): for use up to 30 days.

Silicone: for use up to 30 days.

The nutriSafe2® connection allows the use of these tubes only with other nutriSafe2® devices..

PVC tubes




Size	04 Fr ●	05 Fr ●	06 Fr ●	08 Fr ●	10 Fr ●	Units/Case
40 cm	364.042	364.052	364.062	364.082	-	25
50 cm	361.042	361.052	361.062	361.082	361.102	50
75 cm	363.042	363.052	363.062	363.082	363.102	25
125 cm	362.042	362.052	362.062	362.082	362.102	25

PUR tubes



Size	04 Fr ●	05 Fr ●	06 Fr ●	08 Fr ●	10 Fr ●	Units/Case
40 cm	1364.042	1364.052	1364.062	1364.082	-	25
50 cm	1361.042	1361.052	1361.062	1361.082	1361.102	50
75 cm	1363.042	1363.052	1363.062	1363.082	1363.102	25
125 cm	1362.042	1362.052	1362.062	1362.082	1362.102	25

Silicone tubes



Size	04 Fr ●	05 Fr ●	06 Fr ●	08 Fr ●	10 Fr ●	Units/Case
50 cm	2331.042G	2331.052G	2331.062G	2331.082G	2331.102	10
125 cm	2332.042G	-	2332.062G	2332.082G	2332.102	10
125 cm (weighted)	-	-	2395.062G	2395.082G	-	10

nutrisafe² Syringes

A complete range of syringes, from 1 ml to 60ml for administration of enteral nutrition and oral medication.

Amber colored syringes allow preservation of light sensitive drugs. They are available with caps.

These syringes only connect with other items in the nutrisafe²® range.



Transparent syringes



Amber syringes

C-GON - Transparent syringes

Product Code	Volume	Graduations	Connection	Units/Case
1015.012	1 ml	0.01 ml	Centered / slip	100
1015.022	2.5 ml	0.1 ml	Centered / slip	100
1015.052	5 ml	0.2 ml	Centered / slip	100
1015.102	10 ml	0.2 ml	Centered / lock	100
1015.202	20 ml	1 ml	Centered / lock	50
1015.352	35 ml	1 ml	Centered / lock	50
1015.602	60 ml	1 ml	Offset / lock	25

Monoject - Transparent syringes

Product Code	Volume	Graduations	Connection	Units/Case
1015.062M	6 ml	0.1 ml	Centered / slip	100
1015.202M	20 ml	1 ml	Centered / lock	50
1015.352M	35 ml	1 ml	Offset / lock	25
1015.602M	60 ml	1 ml	Offset / lock	25

Amber syringes (bulk non-sterile)

Product Code	Volume	Graduations	Connection	Units/Case
1022.012	1 ml	0.01 ml	Centered / slip	50
1022.022	2.5 ml	0.1 ml	Centered / slip	50
1022.052	5 ml	0.2 ml	Centered / lock	50
1022.102	10 ml	1 ml	Centered / lock	50

nutrisafe² Extension lines

Extension lines to allow for positioning of the syringe driver away from the incubator and to move handling away from the newborn.

Available in standard bore for breast milk or in large bore for fortified milk.



Product Code	Int. ø	Length	Clamp	Units/Case
368.032-50	1.5 mm	30 cm	No	50
368.152-50	1.5 mm	150 cm	No	50
369.152	2.8 mm	150 cm	No	25

nutrisafe² Caps

A range of caps to secure the contents of vials containing oral / enteral medications. These caps allow access to the medication only from a nutrisafe² syringe, to avoid risk of accidental injection through IV systems.

Caps for syringes and tubes prevent contamination of the feed or drugs between feeding sessions.

Needleless vial access caps

Product Code	Connector	Int. ø	Units/Case
821.022	male	20 mm	25



Needleless vial access caps

Caps for syringes and tubes

Product Code	Description	Units/Case
828.01	Cap for syringe (x1)	100
828.08	Cap for syringe (x8)	100



Syringe cap

Vented caps for syringes - for gravity feeding

Product Code	Description	Units/Case
818.22	20 ml	50
818.62	60 ml	50



Vented cap

nutrisafe2® sampling straws



nutrisafe2® drawing-up straws allow the filling of nutrisafe2® syringes with feed or medication from all sizes of containers and bottles.

Product Code	Size	Bore	Units/Case
817.052	5 cm	Standard bore	50
817.152	15 cm	Standard bore	50
817.153	15 cm	Large bore	50
817.202	20 cm	Standard bore	50
817.203	20 cm	Large bore	50
817.302	30 cm	Large bore	50
817.502	50 cm	Standard bore	20

nutrisafe2® accessories

Connectors, stopcocks and accessories to facilitate daily use of the nutrisafe2® range.

Connectors

Product Code	Description	Units/Case
368.02	Adaptor for feeding sets	100
368.12	Taper connector	30
368.32	Extension tube with connector for suction	25
368.42R	Y connector with oral port	100



Adaptor for feeding sets



Y connector



Taper Connector



Extension Tube with connector

Stopcock

Product Code	Description	Units/Case
5802.01	3-way stopcock	50



Stopcock

Adapters

Product Code	Description	Units/Case
828.88	Oral tip adapter	1000
828.99	Medicine cap adapter	1000



Oral tip adapter



Medicine cap adapter



nutrifit


The perfect fit for pediatrics and adults

Enteral nutrition is recommended when oral food intake proves impossible or insufficient to meet the patient's nutritional needs. Enteral nutrition contributes to:

- ✓ Preserving normal physiology¹, digestion, and nutrient absorption
- ↘ Reducing:
 - Length of stay²
 - Infectious complications³
 - Costs⁴
- ⊕ The reinforcement of the immune system⁵



nutrifit
The ENFit™ brand of Vygon



DISTRIBUTED BY:
ADVANCED SURGI-PHARM INC.
www.advancedsurgipharm.com
Tel:1-800-661-5432

PVC FEEDING TUBES

Adapted for children and adults, the PVC nutritif feeding tubes are intended to be used for short-term enteral feeding (up to 5 days) and feature the following characteristics:

- Made of PVC
- Closed distal tip with 2 lateral eyes (4Fr to 12Fr) and 4 lateral eyes (14Fr to 21Fr)
- Small cap with color code to quickly identify the size of the tube
- Numerical markings every centimeter:
 - from 5cm to 35cm for codes 364, 361 & 363
 - from 5cm to 70cm for codes 366, 362
- Radiopaque
- Transparent to better visualize the administered solution.

Code	PVC tube					Color code	Quantity
	Length cm	Fr	Ext. Ø mm	Int. Ø mm	Prim. vol. ml		Units/case
364.047	40	4	1.5	0.8	0.27	Purple	25 / case
364.057	40	5	1.7	1.0	0.39	Grey	25 / case
364.067	40	6	2.0	1.2	0.54	Light green	25 / case
364.087	40	8	2.5	1.5	0.81	Blue	25 / case
361.047	50	4	1.5	0.8	0.32	Purple	50 / case
361.057	50	5	1.7	1.0	0.47	Grey	50 / case
361.067	50	6	2.0	1.2	0.65	Light green	50 / case
361.087	50	8	2.5	1.5	0.99	Blue	50 / case
361.107	50	10	3.3	2.0	1.7	Black	50 / case
363.047	75	4	1.5	0.8	0.45	Purple	25 / case
363.057	75	5	1.7	1.0	0.67	Grey	25 / case
363.067	75	6	2.0	1.2	0.94	Light green	25 / case
363.087	75	8	2.5	1.5	1.4	Blue	25 / case
363.107	75	10	3.3	2.0	2.5	Black	25 / case
363.127	75	12	4.0	2.5	3.9	White	25 / case
366.087	90	8	2.5	1.5	1.7	Blue	25 / case
366.107	90	10	3.3	2.0	3.0	Black	25 / case
366.127	90	12	4.0	2.5	4.6	White	25 / case
366.147	90	14	4.5	3.0	6.6	Dark green	25 / case
366.167	90	16	5.0	3.5	8.9	Orange	25 / case
366.187	90	18	6.0	4.0	11.7	Red	25 / case
362.047	125	4	1.5	0.8	0.7	Purple	25 / case
362.057	125	5	1.7	1.0	1.06	Grey	25 / case
362.067	125	6	2.0	1.2	1.5	Light green	25 / case
362.087	125	8	2.5	1.5	2.3	Blue	25 / case
362.107	125	10	3.3	2.0	4.0	Black	25 / case
362.127	125	12	4.0	2.5	6.3	White	25 / case
362.147	125	14	4.5	3.0	9.0	Dark green	25 / case
362.167	125	16	5.0	3.5	12.3	Orange	25 / case
362.187	125	18	6.0	4.0	16.0	Red	25 / case

POLYURETHANE FEEDING TUBES WITHOUT STYLET

Adapted for children and adults, the PUR nutrit**fit** feeding tubes are intended to be used for medium-term (1 to 4 weeks) and feature the following characteristics:

- Made of PUR
- Closed distal tip with 2 lateral eyes (4Fr to 12Fr) and 4 lateral eyes (14Fr to 18Fr)
- Small cap with color code to quickly identify the size
- Numerical markings every centimeter:
 - from 5cm to 35cm for codes 1364, 1361 & 1363
 - from 5cm to 70cm for codes 1366, 1362 & 1365
- Radiopaque line from 4Fr to 6Fr, and 3/4 radiopaque tubing from 8Fr
- Transparent to better visualize the administered solution.




Code	PUR tube without stylet					Color code	Quantity
	Length cm	Fr	Ext. Ø mm	Int. Ø mm	Prim. vol. ml		Units/case
1364.047	40	4	1.5	0.8	0.27	Purple	25 / case
1364.057	40	5	1.7	1.0	0.39	Grey	25 / case
1364.067	40	6	2.0	1.2	0.54	Light green	25 / case
1364.087	40	8	2.5	1.5	0.81	Blue	25 / case
1361.047	50	4	1.5	0.8	0.32	Purple	50 / case
1361.057	50	5	1.7	1.0	0.47	Grey	50 / case
1361.067	50	6	2.0	1.2	0.65	Light green	50 / case
1361.087	50	8	2.5	1.5	0.99	Blue	50 / case
1361.107	50	10	3.3	2.0	1.7	Black	50 / case
1363.047	75	4	1.5	0.8	0.45	Purple	25 / case
1363.057	75	5	1.7	1.0	0.67	Grey	25 / case
1363.067	75	6	2.0	1.2	0.94	Light green	25 / case
1363.087	75	8	2.5	1.5	1.4	Blue	25 / case
1363.107	75	10	3.3	2.0	2.5	Black	25 / case
1363.127	75	12	4.0	2.5	3.9	White	25 / case
1366.087	90	8	2.5	1.5	1.7	Blue	25 / case
1366.107	90	10	3.3	2.0	3.0	Black	25 / case
1366.127	90	12	4.0	2.5	4.6	White	25 / case
1362.047	125	4	1.5	0.8	0.7	Purple	25 / case
1362.057	125	5	1.7	1.0	1.06	Grey	25 / case
1362.067	125	6	2.0	1.2	1.5	Light green	25 / case
1362.087	125	8	2.5	1.5	2.3	Blue	25 / case
1362.107	125	10	3.3	2.0	4.0	Black	25 / case
1362.127	125	12	4.0	2.5	6.3	White	25 / case
1362.147	125	14	4.5	3.0	9.0	Dark green	45 / case
1362.167	125	16	5.0	3.5	12.3	Orange	45 / case

POLYURETHANE FEEDING TUBES WITH STYLET

Adapted for children and adults, the PUR nutritif feeding tubes are intended to be used for medium-term enteral feeding (1 to 4 weeks) and feature the following characteristics:

- Made of PUR
- Stylet protected by a ball end to facilitate the placement
- Distal tip with 2 lateral eyes (6Fr to 12Fr) and 4 lateral eyes (14Fr to 18Fr)
- Small cap with color code to quickly identify the size
- Numerical markings every centimeter
 - from 5cm to 35cm for code 1363
 - from 5cm to 70cm for codes 1366 & 1362
- Radiopaque.



Code	PUR tube with stylet						Color code	Quantity
	Length cm	Fr	Ext. Ø mm	Int. Ø mm	Distal tip	Prim. vol. ml		Units/case
1366.0671	90	6	2.0	1.2	Closed	1.1	Light green	25 / case
1366.0871	90	8	2.5	1.7	Closed	2.2	Blue	25 / case
1366.1071	90	10	3.3	2.5	Closed	4.6	Black	25 / case
1366.1271	90	12	4.0	3.2	Closed	7.5	White	25 / case
1366.1471	90	14	4.5	3.5	Closed	9.0	Dark green	25 / case
1366.1671	90	16	5.3	4.0	Closed	11.7	Orange	20 / case
1366.1871	90	18	6.0	4.5	Closed	14.8	Red	20 / case
1366.0872	90	8	2.5	1.7	Open	2.2	Blue	25 / case
1366.1072	90	10	3.3	2.5	Open	4.6	Black	25 / case
1366.1272	90	12	4.0	3.2	Open	7.5	White	25 / case
1362.0671	125	6	2.0	1.2	Closed	1.5	Light green	25 / case
1362.0871	125	8	2.5	1.7	Closed	2.9	Blue	25 / case
1362.1071	125	10	3.3	2.5	Closed	6.3	Black	25 / case
1362.1271	125	12	4.0	3.2	Closed	10.3	White	25 / case
1362.1471	125	14	4.5	3.5	Closed	12.3	Dark green	25 / case
1362.1671	125	16	5.3	4.0	Closed	16.0	Orange	20 / case
1362.1871	125	18	6.0	4.5	Closed	20.3	Red	20 / case
1362.0872	125	8	2.5	1.7	Open	2.9	Blue	25 / case
1362.1072	125	10	3.3	2.5	Open	6.3	Black	25 / case

SYRINGES

LOW DOSE TIP SYRINGES

The Low Dose Tip (LDT) **nutrifit** syringes incorporate a purple-coloured plunger and "Enteral" print, enabling instant identification for enteral use. They are sterile and individually packed.



Code	Transparent LDT syringes			Cap	Quantity
	Volume	Graduation	ENFit® connection		Units/case
1015.018	1 ml	0.01 ml	Centered	No	100 / case
1015.028	2.5 ml	0.1 ml	Centered	No	100 / case
1015.058	5 ml	0.2 ml	Centered	No	100 / case

SYRINGES

TRANSPARENT SYRINGES

The transparent **nutrifit** syringes incorporate a purple-coloured plunger and "Enteral" print, enabling instant identification for enteral use. They are sterile and individually packed.



Code	Transparent syringes			Cap	Quantity
	Volume	Graduation	ENFit® connection		Units/case
1015.107	10 ml	0.5 ml	Offset	No	100 / case
1015.207	20 ml	0.5 ml	Offset	No	80 / case
1015.607	60 ml	1 ml	Offset	No	60 / case
1015.1007	100 ml	2 ml	Centered	No	50 / case

SYRINGES ACCESSORIES

MALE CAPS FOR **nutrifit** SYRINGES

These caps are specifically designed for **nutrifit** syringes. Purple-coloured, it enables instant identification for enteral use. Compatible with Low Dose Tip syringes.



Code	Packaging	Units/case
828.017	Individually packed	200 / case
828.087	Blister containing 8 caps	200 / case

VENTED CAPS

Vented cap are used for **nutrifit** syringes (20 and 60ml) to enable gravity feeding and protecting the solution of enteral feed from contamination.

Code 818.07
Quantity Case: 30




ACCESSORIES

EXTENSION TUBES

STANDARD EXTENSION TUBES

These nutri**fit** extension tubes keep the feeding tube away from the syringe pump or from the feeding pump to improve the patient's comfort.

- Purple-coloured PVC extension tube
- Attached cap to close the tube and avoid contamination depending on the reference.




Code	Tube					Quantity
	Length cm	Ext. Ø mm	Int. Ø mm	Priming volume ml	Cap	Units/case
368.097	90	2.5	1.5	1.38	No	200 / case
368.107	100	2.5	1.5	1.9	No	45 / case
368.157	150	2.5	1.5	2.8	No	45 / case
368.1572	150	2.5	1.5	2.3	No	200 / case

EXTENSION TUBES WITH CLAMP

This nutri**fit** device gives extended connections to pumps and allows greater freedom of movement for ambulatory patients:

- Purple - coloured PVC extension tube
- With Roberts clamp
- Attached cap to close the tube and avoid contamination.



Code	Tube				Quantity
	Length cm	Ext. Ø mm	Int. Ø mm	Priming volume ml	Units/case
368.1071	100	2.5	1.5	1.9	45 / case
368.1571	150	2.5	1.5	2.8	45 / case

SAMPLING DEVICES

DRAWING-UP STRAWS

The nutri**fit** drawing-up straw is a device connected to a nutri**fit** syringe for medication sampling.



Code	Drawing-up straw				Quantity
	Length cm	Ext. Ø mm	Int. Ø mm	Priming volume ml	Units/case
817.057	5	2.0	1.2	0.13	100 / case
817.0572	5	2.0	1.0	0.10	200 / case
817.157	15	3.3	2.0	0.57	50 / case
817.207	20	4.1	2.8	1.3	200 / case
817.208	20	7.4	4.9	3.75	200 / case

BOTTLE CAPS FOR SAMPLING

These nutri**fit** caps secure the medicine bottle and allow the sampling with a nutri**fit** syringe.



Code	Bottle caps for sampling	Quantity
	Ø mm	Units/case
820.1297	12.9	30 / case
820.2047	20.4	30 / case
820.2777	27.7	30 / case

OTHER ACCESSORIES

CHRISTMAS TREE ADAPTERS

These nutri**fit** adapters help customize your sets and help manage the ENFit® transition period.





Code	Product	Quantity Units/case
0VEACC406	Male ENFit® / Christmas tree adaptor (without cap, non-sterile)	50 / case
0VEACC416	Male ENFit® / Christmas tree adaptor (with cap, non-sterile)	50 / case

Premicath®

A 1Fr PICC line designed specially for newborns under 1kg

Infants born prematurely weighing under 1kg often need the nutritional and therapeutic support that can only be delivered into their tiny vessels via venous access. At just 1Fr, Premicath® has been specially designed for the smallest and most fragile babies to deliver optimum accuracy during placement as well as minimising trauma for patients.

Available with choice of introducers including splitting needle and microflash®



Premicath®

Product Code	Catheter			Introducer		Stylet	Units/Box
	Length	Priming Volume	Flow rate (Pressure 1 bar)	Type	Size		
I261.20G	20 cm	0.09 ml	0.7 ml / min	Splitting Needle	24 G	No	10
I261.203G	20 cm	0.09 ml	0.7 ml / min	Splitting Needle	24 G	Yes	10
I261.205G	20 cm	0.09 ml	0.7 ml / min	Neocath Split	22 G	Yes	10
I261.208G	20 cm	0.09 ml	0.7 ml / min	Microflash	20 G	Yes	10
I261.306G	30 cm	0.11 ml	0.6 ml / min	No Introducer	-	Yes	10
Higher Barium							
I261.203A	20 cm	0.09 ml	0.7 ml / min	Splitting Needle	24 G	Yes	10
I261.306A	30 cm	0.11 ml	0.6 ml / min	No Introducer	-	Yes	10

Kits & Accessories

Product Code	Description	Units/Box
I261.20N	24 G Splitting Needle	20
AMS-850TC	V-Cutter Catheter Trimming Tool	20



V-Cutter




Splitting needle



nutri^{line}®

A PICC line with a unique & easy-to-use peelable needle available in 2 Fr

nutri^{line}'s one-piece catheter provides clinicians with a high degree of safety. Available in a choice of lengths to ensure accurate tip placement for most IV access sites.



Product code	Length	Size	Priming volume	Flow rate	Introducer	Stylet	Units/Case
I252.30G	30 cm	2 Fr	0.12 ml	5.0 ml / min	Splitting needle 20 G	No	10
I252.31G	30 cm	2 Fr	0.12 ml	5.0 ml / min	No Introducer	No	10
I252.35G	30 cm	2 Fr	0.12 ml	5.0 ml / min	Microflash 20 G	No	10
I252.030G	30 cm	2 Fr	0.12 ml	5.0 ml / min	Microflash 20 G	Yes	10
I252.031	30 cm	2 Fr	0.12 ml	5.0 ml / min	No Introducer	Yes	10




DISTRIBUTED BY:
ADVANCED SURGI-PHARM INC.
www.advancedsurgipharm.com
Tel:1-800-661-5432

nutriline®

TWINFLO

A 2Fr dual lumen catheter with peelable needle

nutriline® twinflo has been designed specifically for babies who require multiple infusions. Its dual lumen tip design offers enhanced fluid management and avoids incompatible drugs mixing.



Product Code	Catheter			Introducer		Stylet	Units / Case
	Length	Priming Case Volume	Flow rate (pressure 1bar)	Type	Size		
1252.232G	30 cm	2 x 0.2 ml	2 x 1.45 ml / min	No introducer	-	No	10
1252.232M	30 cm	2 x 0.2 ml	2 x 1.45 ml / min	Microflash	20 G	No	10
1252.230G	30 cm	2 x 0.2 ml	2 x 1.45 ml / min	Microflash	20 G	Yes	10

Kits & Accessories

Product Code	Description	Units/Box
7371.19	20 G Microflash Introducer	10
1261.20N	24 G Splitting Needle	20
AMS-850TC	V-Cutter Catheter Trimming Tool	20





introducers

For PICC catheters

Stand alone introducers for neonatal PICC catheters: splitting needle, microflash® peelable cannula and IV cannula.



IV cannula



Splitting needle

IV cannula

Product Code	Size	For catheter	Units/Case
CMS-207	24 G	1 Fr	10

Splitting needle

Product Code	Size	For catheter	Units/Case
1261.20N	24 G	1 Fr	20
1262.30N	20G	2 Fr	20




microflash™

Peelable Introducer/Peelable Safety Introducer*

Peelable Introducer for the insertion of PICC & Midline catheters, with notched needle and clear cannula sheath for early detection of blood flashback

- A** Top hole needle shaft permits blood flashback between needle & cannula for immediate indication of vessel access
- B** Short, light weight safety housing for balance and ease of handling
- C** Curved hub for secure ergonomic fingertip placement
- D** Tri-bevel needle tip for ease of skin penetration
- E** Smooth needle to cannula transition reduces drag through the tissue
- F** Transparent cannula permits immediate blood visualization upon vessel penetration
- G** Funnel cannula opening to facilitate easy catheter entry
- H** Radiopaque stripe on cannula for patient safety

Product Code	Size	For catheter	Units/Case
7371.19	20 G	1 Fr / 2 Fr	10



*The microflash® Introducer is intended for use as a split cannula introducer used to puncture peripheral veins.




DISTRIBUTED BY:
ADVANCED SURGI-PHARM INC.
www.advancedsurgipharm.com
Tel:1-800-661-5432

Griplok

Catheter Securement Device

Our Vygon Grip-Lok® is especially efficient with Vygon Neonatal Catheters. It is hypoallergenic and latex free, which means low risk of irritation and allergic reaction. Its zinc oxide adhesive allows for gentle, yet efficient skin fixation. The catheter fixation is secure to prevent migration and its suture-less fixation increases ease of use.

- Hypoallergenic
- Zinc Oxide Adhesive
- Specially designed hub insert
- Suture-less fixation
- Sterile
- Latex free



Product Code	Description	Units/Case
5804.08	Grip-Lok Catheter Securement Device	20

umbilical catheters

Single-lumen umbilical catheters

A single lumen umbilical catheter used for venous or arterial access.

Umbilical venous catheters are used for:

- Parenteral nutrition & fluid administration
- Drugs administration
- Venous blood sampling
- Transfusion of blood or blood products
- Exchange transfusion

Umbilical arterial catheters are used for:

- Arterial blood sampling
- Arterial pressure measurement
- Blood pH and blood gas analysis
- Fluids and drugs administration **Umbilical**



Three-way tap


- Colour-coded caps for venous / arterial identification

Luer-lock connector

- For efficient connection

Numerical markings every cm

- Make it easy to check the length of catheter inserted



Product Code	Size	Length	Priming vol	Flow rate	Material	Units/Box
1270.02	2.5 Fr	30 cm	0.12 ml	>3 ml / min	PUR	8
1270.03	3.5 Fr	40 cm	0.27 ml	>10 ml/min	PUR	8
1270.04	4 Fr	40 cm	0.29 ml	>10 ml/min	PUR	20
1270.05	5 Fr	40 cm	0.40 ml	>23 ml/min	PUR	20
1270.08	8 Fr	40 cm	0.81 ml	>100 ml/min	PUR	20
270.04	4 Fr	37 cm	0.27 ml	>16 ml/min	PVC	15
270.05	5 Fr	37 cm	0.38 ml	>31 ml/min	PVC	15
270.06	6 Fr	37 cm	0.51 ml	>52 ml/min	PVC	15
270.08	8 Fr	37 cm	0.75 ml	>106 ml/min	PVC	15

Double lumen umbilical catheter

A double lumen umbilical catheter for short term catheterisation through the umbilical vein when multiple access is needed.

- Infusion of incompatible drugs
- Reduces the need for additional access
- Infusion of medication and nutrition at the same time



Colour-coded hubs

- For easy identification and safety



Double lumen

- Clearly marked lines with separate exit points prevent incompatible drugs mixing



Neonatal clamps

- Small line clamp for line management

Numerical markings every cm

- Make it easy to check the length of catheter inserted



X-ray opaque PUR catheter

- X-ray opaque for accurate tip location without additional contrast medium

	Product Code	Size	Exit Points	Length	Priming volume	Flow rate	Units/Box
PUR	1272.14	4Fr	One lateral One distal	20 cm	2 x 0.19 ml	11 & 12 ml/min	10
	1274.14	4Fr		40 cm	2 x 0.27 ml	6 & 7 ml/min	10
	1274.17	5Fr		40 cm	2 x 0.26 ml	7 & 9 ml/min	10

Accessories for umbilical catheters

Lockable three-way stopcock with identification caps (red and blue) to distinguish between venous and arterial access.


- User feels digitally the position chosen
- Safety Luer-lock prevents accidental disconnections
- Lipid-resistant



Product Code	Description	Units/Case
876.20	Lockable three-way stopcock with identification caps (red and blue)	50

Triple lumen umbilical catheter

A triple lumen umbilical catheter to administer TPN and for the injection of drugs; particularly those with incompatible features or when one drug needs to be injected continually.



	Product Code	Size	Exit Points	Length	Priming volume	Flow rate	Units/Box
PUR	V02127312	4.5Fr	One lateral One distal One proximal	20 cm	2 x 0.26 m	2, 2.1 & 7.4 ml/min	10

umbilical placement pack

Our umbilical placement pack contains everything required to safely & efficiently place an umbilical catheter



The benefits of our umbilical packs



Cost saving

A complete pack costs less than buying the individual components.



Delivers better patient care

Pack enables Aseptic Non Touch Technique for catheter placement to minimise infections. Also allows traceability of pack components used.



Reinforces best practice

Pack includes a heat sealed drape (to maximise the sterile barrier) and a safety hypodermic needle (to minimise the risk of accidental needlestick injuries).



Standardises procedures

Consistent clinical practice assured with identical pack contents.



More environmentally friendly

Less waste with a pack compared to multiple packaging of individual components.

Each pack contains

	18G hypodermic safety needle		20G hypodermic safety needle		Scalpel
	Drapes without adhesive		40x40cm fenestrated transparent drape with "Easy-peel"		Drape towels
	1x pair of suture scissors		10x swabs 10x10cm, 8-ply		10x swabs 5x5cm, 8-ply
	1x transparent gallipot 60mL		1x red gallipot 60mL		Outer wrap
	Instructions for use		Umbilical tape		Measuring guide
	Mosquito haemostat curved		Mosquito haemostat straight		Iris forceps full curved
	Iris forceps half curve		Iris forceps straight		Iris forceps straight with teeth
	Needle holder		Measuring tape		1mL syringe
	3mL syringe		5mL syringe		Vessel dilator probe
	3.0 silk suture with curved cutting needle		Tray 20x15x4cm		2x fixation strips

exchange transfusion tray

Set containing the accessories necessary to perform removal and replacement of a patient's blood.



Product Code	Contents	Units/Case
275.001	1x four-way stopcock equipped with 1 injection site (for injections of additional medication)	1
	1x extension tube for evacuation of the discarded blood	
	1x 5Fr exchange transfusion catheter (PVC - XRO) 37 cm long, with 2 sides holes and centimetric markings	
	1x 7Fr exchange transfusion catheter (PVC - XRO) 37 cm long, with 2 sides holes and centimetric markings	
	2x 20 ml Luer-lock syringes	
	1x 10 ml Luer-lock syringe	
	1x 16x0.5 mm (25G) Luer hypodermic needle	
	1x graduated plastic container (2000 ml)	
	1x transfusion set, 163 cm long, 18.8 ml dead volume	
	1x 15 cm ruler for venous pressure measurement	
	3x 50 x 50 mm gauze dressings	
	1x 50 x 60 cm fenestrated drape	
	1x air-venting needle	
	1x control sheet	

PICC insertion trays

Kits containing the instrumentation and accessories needed to perform a PICC insertion.

- All components in one pack
- Single use components
- Instruments dedicated to neonates



PICC insertion tray

Product Code	Contents	Packaging
AMS-8198-1	1x Autowrap, 36" x 36"	20 / case
	2x drapes, 24" x 24"	
	1x polylined fenestrated drape, 12" x 19"	
	1x neonatal tourniquet	
	2x tape measures, 24"	
	1x pkg. pvp swabsticks (3/pkg)	
	2x saline prep pads	
	2x 12cc syringes	
	1x iris scissor, 3.5"	
	1x eye dressing forcep half curved	
	4x dressing sponges, 2" x 2"	

goldtrace

The Gold-plated Fetal Scalp Electrode

A message from the heart

Goldtrace® Fetal Spiral Electrode (FSE) is used when fetal heart rate monitoring during labor requires continuous and reliable internal registration. Goldtrace FSE is made of materials whose properties ensure the best-possible-quality signal for CTG.

Goldtrace FSE comes in a flexible protective guide tube that keeps the spiral tip inside the tube during application, protecting the user, mother and baby. The flexibility of the guide tube makes it possible to form according to desired anatomical curvature.

Specifications


- Gold-plated spiral electrode
- Single use only
- Includes disposable applicator
- Not made with natural rubber latex
- Sterile product

Article number

N-CNS000004 Goldtrace Fetal Spiral Electrode in boxes of 75 single packed pcs

Materials

Spiral needle	316 Stainless Steel with 10-25µm Gold
Reference plate	316 Stainless Steel with 10-25µm Gold
Hub	PE
Connector	PE and brass, gold-plated
Wires	Copper with PVC jacket
Drive tube	PE
Guide tube	PE
Drive tube grip	PE



FECG Legplate Adapter for Philips

The Connection Between Mother and Child

The FECG Legplate Adapter for Philips provides a signal path between the Fetal Spiral Electrode (FSE), Skin Electrode, and Fetal Monitor when continuously monitoring fetal heart rate during labor and delivery procedures.



Use a FECG Legplate that is compatible with your Fetal Monitor and connect only compatible FSE to the FECG Legplate.

The FECG Legplate Adapter for Philips is compatible with:

- Philips Avalon Monitor

Article Number

N-ACC300120-R0 FECG Legplate Adapter for Philips

neohelp™

A sterile heat loss prevention suit to deliver life-saving thermal care

Immediately after a premature baby is delivered by c-section, the priority is to place them in a warm and sterile environment. neohelp™ provides this essential protection so other essential Golden Hour care and treatment, including delayed cord clamping, can take place.^(2,3)

Integrated adjustable hood

Decreases heat loss through radiation:

- More efficient than a stockinette cap allowing air to pass through the material
- The toggle allows you to adjust to the baby's head

Double layer of soft clear polyethylene

Decreases heat loss through convection and evaporation:

- Creates a warm and humid environment, mimicking the incubator effect
- The thin inner allows for excellent skin contact
- Creates a barrier against drafts
- Allows passage of radiant heat from an additional warming device (if used)⁽⁴⁾
- Transparent material allows vital observations⁽⁵⁾

Pre-shaped foam support


Decreases heat loss through conduction:

- Helps to maintain an open airway by raising the shoulders
- Stabilises the baby's position
- Provides thermal care during transportation
- Provides comfort

Central VELCRO® opening

Provides optimum seal to ensure heat conservation:

- Quick and easy to place around the baby
- Allows full access to the baby's body
- Designed for easy placement of monitoring equipment, IV, umbilical catheters and carrying out Golden Hour care



Why is thermal care important?

Thermal care is vital for a preterm infant because they may have unbalanced skin-surface to weight ratio, very little or no capacity to generate heat (brown adipose tissue), inadequate stores of subcutaneous (insulating) fat and immature epidermal barrier. ^(5,6)

At this vital time, neohelp prevents heat loss through its double layer of soft, clear polyethylene, integrated adjustable hood and VELCRO® seal.

<p>For every 1°C decrease: Risk of sepsis increases by 11% Risk of death increases by 28%⁽⁴⁾</p>	<p>In the first 10-20 minutes, without any protection, temperature can fall by 2-4°C⁽⁷⁾</p>
--	--

Why is delayed cord clamping important?

Provided the baby can be kept warm and does not need immediate resuscitation, the Resuscitation Council UK (RCUK) recommends delayed cord clamping (DCC) for at least 60 seconds whilst breathing is established. ⁽⁸⁾ Delayed cord clamping has been shown to reduce the relative risk of:

- Intraventricular haemorrhage by 41% (RR 0.59, 95% CI 0.41 to 0.85)⁽⁹⁾
- Necrotising enterocolitis by 38% (RR 0.62, 95% CI 0.43 to 0.90)⁽⁹⁾

In addition, this procedure increases circulation of blood volume after birth and an improvement in cardiovascular stability, reducing the need for a blood transfusion.

Ordering codes

Item #	Baby weight	Description	Dimensions	Units/Box
37.09.14	< 1kg	neohelp™ Small	L. 38 cm x W. 30 cm	10
37.09.15	1-2.5 kg	neohelp™ Medium	L. 44 cm x W. 38 cm	10
37.09.16	> 2.5 kg	neohelp™ Large	L. 50 cm x W. 38 cm	10

Warm Gel® Infant Heel Warmer

Uniform heat to reduce heel sticks

WarmGel® Infant Heel Warmer is a gel-filled, disposable heel warmer designed to provide gentle heat to increase blood flow and enhance infant blood sample collection.

WarmGel® Technology

Effective

- Conforms to infant's shape, keeping warmth where needed
- Even heat distribution increases blood flow

Safe

- Non-toxic food-grade gel won't harm infant's skin or eyes
- Single-patient use minimizes risk of infection

Reliable

- Reaches maximum temperature in less than 60 seconds, maintains optimal 104° F for 10 minutes¹
- Trigger disc minimizes auto-activation during shipment



WarmGel® Versus Liquid

- Conforms to infant's shape
- Provides even heat distribution
- Resists migration, keeping warmth where needed
- Competitively priced compared to liquid-filled products

Item #	Description
24401	Liquid Heel Warmer (100/Box)
20418	WarmGel® Infant Heel Warmer (100/Box)

¹"Neonatal Thermoregulation: Guidelines for Practice," National Association of Neonatal Nurses, 1997.



TransWarmer®

Transport Mattress

Consistent, Safe Warming for Newborns



The TransWarmer® Infant Transport Mattress provides safe and reliable heat for up to two hours when cold stress is a concern. This thermostable, nontoxic gel mattress nests the infant in soothing warmth, while minimizing vibration during transport. Based on a 24 °C/75 °F start temperature, the mattress reaches a peak temperature of 40 °C/104 °F in less than 60 seconds - quickly improving the infant's admission temperature and reducing the risk of hypothermia.



Features and Benefits

- Comfortable, gel-filled mattress cradles baby in therapeutic warmth
- Gently conforms to the patient's shape for better nesting
- Minimizes the risk of hypothermia in cesarean deliveries
- Better heat distribution increases blood flow
- Disposable heat pack minimizes risk of cross contamination and infection
- Patented disc mechanism is easy to activate
- TransWarmer is nontoxic, latex- and PVC-free
- Use in labour and delivery, NICU, nursery, OR, ER, and transport units
- 10" wide x 16" tall

Item #	Description
20421	TransWarmer® Infant Transport Mattress (6/Box)




DISTRIBUTED BY:
ADVANCED SURGI-PHARM INC.
www.advancedsurgipharma.com
Tel: 1-800-661-5432

surfcath™

A catheter for surfactant administration using the LISA method

For babies needing treatment for respiratory distress syndrome (RDS), our new surfcath™ uses the LISA (Less Invasive Surfactant Administration) method to place the catheter. It was specifically designed to improve manoeuvrability during placement whilst also eliminating the need for the commonly used Magill forceps.



2cm soft distal black tip


- 2cm mark shows when surfcath™ is in place
- Soft tip minimises risk of tracheal lesions and prevents kinking
- Pre-curved to follow the airways anatomy

Transparent bendable thermoplastic material

- Transparent material allows visual check on the delivery of surfactant
- Thermosensitive material allows you to curve surfcath™ prior to use meaning no need for Magill forceps
- Semi-rigid material also allows high manoeuvrability, helping to follow the airways anatomy


Less invasive 6Fr gauge

- No obstruction of the airways allowing spontaneous breathing
- Low dead space (0,2 mL)



20cm length with cm markings

- Longer length eases manipulation of the surfactant syringe away from the patients head
- Centimetre markings indicate the inserted length and check that surfcath™ stays in place



Treating respiratory distress syndrome (RDS)

RDS is the most common lung disease in premature infants and is a pulmonary disorder resulting from a surfactant deficiency. In Europe RDS is observed for about 90% of babies born at 24 weeks and for 80% babies born at 28 weeks of gestation.⁽¹⁰⁾

International guidelines recommend

- Preterm infants should be managed without mechanical ventilation where possible
- CPAP (Continuous positive airway pressure) with early rescue surfactant is considered optimal management for babies with RDS
- LISA is the preferred mode of surfactant administration for spontaneously breathing babies on CPAP, provided that clinicians are experienced with this technique.⁽¹⁰⁾


Ordering code

Code	Size	Length	Priming vol	Units/Box
5590.106	6 Fr	20 cm	0.2 ml	10

Scan to watch our video

on how to use surfcath™ surfactant injection catheter

Please note: The link will take you to an external website




endotracheal tubes

Single and double lumen tubes in PVC (DEHP-free) for emergency oral or nasal intubation and routine anesthesia. Double lumen tubes have secondary lumen for surfactant administration or airway monitoring.

- Radiopaque
- 15mm adapter
- Markings every half centimeter

Single lumen tubes are also available as a soft tube designed to minimize tracheal damage during long term intubation.



Plain endotracheal tubes

Product Code		Ext. ø	Int. ø	Length	Units/Case
Standard	Soft tube				
516.20	520.20	3.4 mm	2.0 mm	165 mm	25
516.25	520.25	4.1 mm	2.5 mm	165 mm	25
516.30	520.30	4.6 mm	3.0 mm	165 mm	25
516.35	520.35	5.2 mm	3.5 mm	165 mm	25
516.40	520.40	5.7 mm	4.0 mm	230 mm	20
516.45	520.45	6.2 mm	4.5 mm	230 mm	20

Endotracheal tubes with secondary lumen

Product Code	Tube			Secondary lumen		Units/Case
	Ext. ø	Int. ø	Length	Lumen diam.	Flow rate	
5516.20	3.4 mm	2.0 mm	165 mm	0.40 mm	2.15 ml/min	20
5516.25	4.1 mm	2.5 mm	165 mm	0.55 mm	3.35 ml/min	20
5516.30	4.6 mm	3.0 mm	165 mm	0.65 mm	5.00 ml/min	20
5516.35	5.2 mm	3.5 mm	165 mm	0.80 mm	15.0 ml/min	20
5516.40	5.7 mm	4.0 mm	230 mm	0.95 mm	35.0 ml/min	20

mucus extractors

Mucaid allows for quick and safe suctioning and sampling of the upper respiratory tract in the newborn using mechanical / wall suction. This avoids the dangers of self-contamination inherent in oral suctioning (De Lee method).

Included in the range are Mucaid with a secure threaded cap for safe transportation of sample for testing and Mucaid with a filter to prevent any mucus or amniotic fluid from passing through into the line.

Meconium aspirators allows for suctioning below the vocal cords to clear meconium from the lungs and airways.

- Finger tip control
- Endotracheal tube and endotracheal tube hub
- Rounded, soft tip
- Centimeters every half centimeter



Mucaid



Meconium aspirator

Mucaid: Mucus extractor

Product Code	Volume	Tubing to suction source length	Catheter	Units/Case
542.06	25 ml	24 cm	6 Fr - 33 cm	25
542.08	25 ml	24 cm	8 Fr - 33 cm	25
542.10	25 ml	24 cm	10 Fr - 33 cm	25

Mucaid with filter

Product Code	Volume	Tubing to suction source length	Catheter	Units/Case
80.542.10B	25 ml	24 cm	10 Fr - 33 cm	25

Meconium aspirator

Product Code	Int. Ø	Length	Units/Case
80.516.30D-25	3.0 mm	165 mm	25
80.516.35D-25	3.5 mm	165 mm	25

nasal ventilation

Single ventilation tube

Single silicone nasal tubes designed to reduce the stress related to ventilation and to improve the comfort of premature babies and newborns.

- Open tip
- 15mm adapter



Single ventilation tube

Double tube in silicone for nasal neonatal ventilation

Length of tube on 2 branches: 6 cm, except code 2596.30: 7 cm.

Open and rounded distal tips and centimetric markings.

The proximal end has an ISO standard 15 mm adaptor. This tube is specially designed to avoid tracheal intubation which is very traumatic, and to improve the comfort of premature babies and neonates.

4 sizes are available: Ext. Ø 2.7 mm - 3.4 mm - 4.0 mm - 5.0 mm




Double tube in silicone for nasal neonatal ventilation

Product Code	Int. ø	Ext. ø	Diameter	Length	Units/Case
Single Tube					
2595.15	1.5 mm	2.7 mm	8 Fr	9 cm	50
2595.20	2.0 mm	3.4 mm	10 Fr	9 cm	50
2595.25	2.5 mm	4.0 mm	12 Fr	9 cm	50
2595.30	3.0 mm	5.0 mm	14 Fr	9 cm	50
Double Tube					
2596.15	1.5 mm	2.7 mm	8 Fr	6 cm	50
2596.20	2.0 mm	3.4 mm	10 Fr	6 cm	50
2596.25	2.5 mm	4.0 mm	12 Fr	6 cm	50
2596.30	3.0 mm	5.0 mm	15 Fr	7 cm	50

SAFE SPOT

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE

Ventlab



SafeSpot™
SS3200-Series
Shown with integrated pressure manometer, pop-off valve, inflatable mask and inflatable O₂ reservoir

Safe Spot™

- Single patient use resuscitation/ventilation system designed specifically for infant patients
- Bag features ultra-thin 'spots' to provide an enhanced tactile feel for lung compliance and increased patient safety
- Integrated pressure manometer promotes proper delivery (<20 cm H₂O) to guard against over-pressurization of babies' fragile lungs; the tri-color fields alert to the risk level of delivered pressure
- Accurate pressure delivery helps reduce barotrauma or pneumothoraces
- Pop-off valve (40 cm H₂O and 25 cm H₂O available) enables clinicians to use the pop off level appropriate for a given patient population creating a greater margin of safety for the baby
- Optional Vent Mask in various sizes available
- Select from a variety of components to assemble a custom resuscitation bag to meet your unique patient needs

ITEM#	SIZE	VOLUME	STROKE VOLUME	BODY MASS	PK
SS3200MB	Infant	250 mL	150-170 mL	≤ 10 kg (22 lbs)	10

CUSTOMIZE YOUR MANUAL RESUSCITATOR

Choose from the following components to create a custom resuscitation/ventilation bag to best meet your needs:

SAFETY COMPONENTS



INTEGRATED PRESSURE MANOMETER

Ventlab's innovative color-coded dial clearly denotes pressure delivery:

- Green** - Target pressure level <math>< 20 \text{ cm H}_2\text{O}</math>
- Yellow** - Potentially dangerous pressure level
- Red** - Dangerous pressure level



POP-OFF VALVES

40 cm H₂O & 25 cm H₂O
Allows clinician to use an appropriate pop-off level



FILTERS

A bacterial/viral filter provides effective protection for the healthcare provider during expiration



CO₂ INDICATOR

Choose from STAT-Check II or MaxCap for effective assessment of CO₂ exchange



PEEP VALVES

Adjustable PEEP valves;
5-20 cm/H₂O

SafeSpot™, SS3200-Series

Shown with integrated pressure manometer, pop-off valve, inflatable mask and corrugated tubing



SafeSpot™ bags come standard with an inflatable Vent Mask

OXYGEN RESERVOIRS



INFLATABLE O₂ BAG

Closed system oxygen reservoir delivers higher FDO₂ and provides visual indication of proper airflow



CORRUGATED TUBING

For use with secondary oxygen supply; available in 24" and 48" lengths

CONVENIENCE COMPONENTS



FACE MASKS

Inflatable and non-inflatable face masks available



FLEX NECK

Removable, flexible joint allows corrugated tubing to be joined with the patient valve; also reduces the torque on ET tube

INCA®

Infant Nasal CPAP Assembly

INCA® Infant Nasal CPAP Assembly is the most respected Continuous Positive Airway Pressure (CPAP) device among clinicians. INCA affords all the benefits of CPAP/Bubble CPAP for newborns in your care.

Features enhanced precision fit for neonates:

- Five anatomically proportional nasal cannula sizes for neonates
- Clear, extra-soft tubing provides custom shaping and maximal stability
- Chin strap is available to help prevent air leakage around the mouth
- Available in Complete and Replacement Sets



INCA COMPLETE SET CA CE

Consists of a nasal cannula, two lengths of Stay-Flex tubing, two tubing adapters, one in-line "T" pressure adapter, a knit cap with hook/loop strap, two foam blocks and a hook/loop chin strap. One INCA Sizing Gauge is included in every box of five sets.




Item #	Description
44-2707	Set with 7.5 Fr Cannula (5/Box)
44-2709	Set with 9 Fr Cannula (5/Box)
44-2710	Set with 10.5 Fr Cannula (5/Box)
44-2712	Set with 12 Fr Cannula (5/Box)
44-2715	Set with 15 Fr Cannula (5/Box)

INCA ACCESSORY PACK CE

Includes five knit caps with hook/loop strap, five chin straps with ten hook fasteners (not attached) and ten foam blocks.

1 Pack = 5 Hats and 10 Foam Blocks



Item #	Description
44-2001	X-Small Pack (5 XS Caps, 10 Foam Blocks)
44-2002	Small Pack (5 S Caps, 10 Foam Blocks)
44-2003	Medium Pack (5 M Caps, 10 Foam Blocks)
44-2004	Large Pack (5 L Caps, 10 Foam Blocks)

INCA REPLACEMENT SET CA CE

Includes the nasal cannula, two lengths of Stay-Flex tubing, two tubing adapters and one in-line pressure "T" adapter.



Item #	Description
44-0707	Set with 7.5 Fr Cannula (5/Box)
44-0709	Set with 9 Fr Cannula (5/Box)
44-0710	Set with 10.5 Fr Cannula (5/Box)
44-0712	Set with 12 Fr Cannula (5/Box)
44-0715	Set with 15 Fr Cannula (5/Box)

INCA NASAL SIZING GAUGE CA CE

For appropriate nasal prong fitting, use the INCA Sizing Gauge. Wipe gauge with sterile 70 percent isopropyl alcohol prior to each use. Hold the gauge up to the nasal area to determine appropriate cannula size.



Item #	Description
025-001	INCA Nasal Sizing Gauge (1/Pack)

NEO-fit™

Neonatal Endotracheal Tube Grip CA CE

Facilitates Endotracheal Tube Management

The NEO-fit provides quick, secure endotracheal tube stabilization. Designed to firmly grip a Neonatal endotracheal tube, NEO-fit utilizes metal clips to mechanically grip the soft material of the tube. The device is positioned on baby's face using a hydrocolloid adhesive that is hypoallergenic and Hook / Loop strap that allows for rapid release and readjustment of the tube. MRI Conditional.



Item #	Description
42-2540	NEO-Fit Neonatal Endotracheal Tube Grip Set (20/Box) Each Set includes: (1) plastic tube grip, with hook/loop strap, (2) foam foot pads, (2) foam "lollipop" tapes. Non-sterile.

NEO-prep™

Skin Barrier CE

A must-have product for your CPAP patients. Protect their tender tissues during this critical therapy. Neo-prep Skin Barriers are precut and non-sterile.



Item #	Description
44-8000	Precut Hypoallergenic Skin Barrier, non-sterile, small (10/Box)
44-8001	Precut Hypoallergenic Skin Barrier, non-sterile, large (10/Box)

Os Finder™

Cervical Os and Canal Dilator

The Os Finder naturally seeks and gently opens the cervical Os. Made of an extremely smooth and malleable Teflon-like material, this instrument is ideal for use in intracervical dilation or uterine sounding, or for intrauterine procedures.

DEVICE DESCRIPTION

The sterile, single use Os Finder™ naturally seeks and gently opens the cervical Os. It is made of extremely smooth and malleable teflon-like material. The dilator is 21.5 cm long with a taper from 2 to 4 mm. CAUTION: Federal restricts this device to sale by or on the order of physician.

CONTRAINDICATIONS

This instrument should not be used under any of the following conditions:

- Pregnancy or suspicion of pregnancy
- Acute or active inflammation of the cervix, endometrium, fallopian tube, ovary or peritoneum
- In the case of cervicitis, endometritis, tubo-ovarian inflammatory disease or pelvic inflammatory disease
- Invasive cancer that is visible on examination

PRECAUTIONS

Do not use the Os Finder if the package or instrument are damaged.

INDICATIONS FOR USE

This instrument is ideal for use in intracervical dilation, uterine sounding, or for intrauterine procedures.

DIRECTIONS FOR USE

The Os Finder is used to find the external Os by simply pushing the narrow end gently against the cervix until it finds its way. Further insertion of the flexible tip will follow the natural path of the canal until it reaches and goes through the internal Os. Dispose of in accordance with all applicable Federal, State and local Medical/Hazardous waste practices.




Item #	Description
1176	Os Finder (50/Box)

CE
Ver.2021-04-22

Tampa Infusion Catheter

For sonohysterography - the simple, cost-effective choice for ultrasound examination of the uterus.

- Unique introducer sheath (33 cm long) facilitates insertion into cervical Os with minimal trauma
- 7 cm marker provides optimal placement of catheter tip in uterus
- Balloonless configuration is ideal for postmenopausal patients with stenotic cervical canals, as well as for many nulliparous patients
- Small catheter size allows easy placement next to ultrasonography probe



Item #	Description
61-2005	Tampa Catheter (20/Box) 5 French, clear, flexible catheter with introducer sheath, 33 cm long, sterile



Ver.2021-04-22

DISTRIBUTED BY:
ADVANCED SURGI-PHARM INC.
www.advancedsurgipharma.com
Tel.1-800-661-5432



H/S Catheters


For Enhanced Visualization of Uterine Pathology

H/S CATHETER SET CA CE

These easy-to-use devices provide enhanced visualization. Designed for use in sonohysterography and hysterosalpingography, the catheter can be used with oil- or saline-based media depending on the procedure requirements.

Features and Benefits

- Available with integrated stylet for easy insertion
- Unique placement sheath eliminates the need for a tenaculum
- Comes in both 5 and 7 French
- Ideal for office use
- Latex-Free construction



Pre-Saline Infusion

Post-Saline Infusion

Item #	Description
61-5005	H/S Catheter, Kraton, 5 Fr (10/Box)
61-5007	H/S Catheter, Kraton, 7 Fr (10/Box)
61-0020	H/S Stylet Catheter Stiffener (10/Box)
61-3005	H/S Catheter with integrated Stylet 5FR (10/Box)
61-3007	H/S Catheter with integrated Stylet 7FR (10/Box)

Ver.2021-04-22



DISTRIBUTED BY:
 ADVANCED SURGI-PHARM INC.
www.advancedsurgipharm.com
 Tel:1-800-661-5432

H/S Elliptosphere™ Catheter Set and Tray

The Safe and Simple Answer for Both
Sonohysterography and Hysterosalpingography



The patented Elliptosphere H/S catheter offers several advantages over the traditional balloon catheters including the ability to gradually inflate the balloon for placement in any size cervical canal, the ability to see the entire lower uterine segment during an SIS procedure, and greater patient comfort.

The Elliptosphere balloon contour and variable size permit secure placement in the cervical canal or uterine cavity for optimal visualization. The balloon's unique design prevents retrograde leakage of distension media, which facilitates peak uterine distention. A small catheter size enables placement next to ultrasound probe for minimal patient discomfort. The H/S Elliptosphere Catheter can be used with oil- or saline-based medium. Available in two versions – with and without an integrated stylet – the H/S Elliptosphere features latex-free construction.



Item #	Description
61-4005	H/S Elliptosphere Catheter, 5 Fr (10/Box)
61-4005S	H/S Elliptosphere Catheter Set with Integrated Stylet (10/Box)



Comfort Belly Bands

Abdominal Fetal Monitoring Band

A lighter, cooler abdominal fetal monitoring band, Belly Bands will be appreciated by your patients as they hold vital tocos and transducers firmly in place. Belly Bands are available in a wide range of widths to comfortably cover any patient abdomen.



Features and Benefits

- Breathable material helps wick away moisture from the mother's skin
- Lighter and cooler to ensure patient comfort
- Flexible precut lengths cover the belly and secure transducer placement
- Available in 30-foot rolls or individually wrapped packages of precut 14" lengths
- Choose from 10-inch to 18-inch widths
- Economical, easy-to-use product

Item #	Description
13010	10" Width Belly Bands on Roll (10" x 10 yds) (1 Roll)
13014	14" Width Belly Bands on Roll (14" x 10 yds) (1 Roll)
13015	14" Width Belly Bands - Precut (14" x 14") (25/Box)



LAMINARIA

MedGyn's Laminaria is all-natural cervical disposable and sterile dilator made from a sea-grown plant containing no synthetic materials. Laminaria provides gentle, gradual dilation.

Features & Benefits

MedGyn's disposable and sterile Laminaria is all-natural cervical dilator made from a sea-grown plant containing no synthetic materials. Safely used for over a century, Laminaria provides gentle, gradual dilation. Once inserted into cervical canal, Laminaria absorbs liquid, swelling several times its initial dried diameter, gently dilating the cervix.



McMahon Laminaria Inserter

MedGyn Laminaria Inserter

Indication for use:

Laminaria should be considered whenever cervical dilation and/or softening of the cervix is desired. Some examples are: cervical stenosis, related to dysmenorrhea, considered a possible cause of infertility, resulting from cauterization or conization, placement and removal of intrauterine devices, induction of labor, radium placement, curettage and drainage of uterine cavity, endometrial biopsy, and suction cannula aspiration.

Standard, approx. length 60mm

021003 Laminaria 3mm Small

leadercath

Arterial LeaderCath (PE)

Although it can be inserted in a variety of arteries, the most common sites for the placement of this catheter are the radial and the femoral artery.

The Seldinger technique makes arterial puncture easier particularly in hypotensive patients with a weak peripheral pulse and when arteries are atherosclerosed or narrow.

The product is supplied in a rigid blister pack containing:

- 1 introducer needle
- 1 straight guidewire
- 1 transparent and radiopaque catheter in polyethylene with winged hub



Product Code	Catheter						Puncture Needle				Guidewire		Unit/Box
	Int. Ø mm	Ext. Ø mm	Length cm	Gauge G	Ø Fr	Flow rate ml/min	Int. Ø mm	Ext. Ø	Gauge G	Length mm	Ø mm	Length mm	
115.090	0.6	0.9	8	20	3	24	0.6	0.9	20	38	0.53	20	20
115.11	0.8	1.2	10	18	4	24	0.8	1.1	19	54	0.71	30	20

lectrospiral®

Coiled extension tube

Lectrospiral® is a flexible coiled extension tube that provides patients with freedom of movement.

Flexible

- Coiled tubing expands and contracts providing patients with greater freedom of movement than straight extension lines.
- Ideal to use with pediatric, home care and oncology patients.

Safe

- Polyethylene material allows for compatibility with commonly infused drugs and reduces adsorption, absorption or permeation that may occur with PVC.
- Extension tubing minimizes potential of trauma at IV site.



Product Code	Length	Inner Diameter	Outer Diameter	Priming Volume	Units/Case
1155.60	100 cm	1.0 mm	2.0 mm	1.1 ml	15
1155.65	150 cm	1.0 mm	2.0 mm	1.5 ml	15
1155.70	200 cm	1.0 mm	2.0 mm	2.0 ml	15
1155.80	300 cm	1.0 mm	2.0 mm	2.8 ml	15
1155.90	400 cm	1.0 mm	2.0 mm	3.7 ml	15
1159.60	100 cm	1.5 mm	2.5 mm	1.9 ml	15
1159.65	150 cm	1.5 mm	2.5 mm	2.7 ml	15
1159.70	200 cm	1.5 mm	2.5 mm	3.6 ml	15
1159.80	300 cm	1.5 mm	2.5 mm	5.2 ml	15
1159.90	400 cm	1.5 mm	2.5 mm	6.9 ml	15

Ver.2021-05-20

DISTRIBUTED BY:
 ADVANCED SURGI-PHARM INC.
www.advancedsurgipharm.com
 Tel.1-800-661-5432




lectrocath

A range of polyethylene extension lines, suitable for arterial pressure monitoring and venous medication delivery.

The advantages include:

- a Polyethylene material which allows for compatibility with commonly used drugs and reduces adsorption, absorption or permeation that may occur with PVC
- Latex-Free
- DEHP-Free
- Pyrogen-Free



Feature	Benefit
PE tubing	Drug-compatible
Colour-coded clips	For arterial or venous identification

Product Code	ID - OD	Length	Capacity	Units/Box
1155.01	1.0 - 2.0 mm	11 cm	0.4 ml	50
1155.03	1.0 - 2.0 mm	30 cm	0.6 ml	50
1155.05	1.0 - 2.0 mm	50 cm	0.7 ml	50
1155.10	1.0 - 2.0 mm	100 cm	1.2 ml	50
1155.15	1.0 - 2.0 mm	150 cm	1.6 ml	50
1155.20	1.0 - 2.0 mm	200 cm	2.0 ml	50
1159.01	1.5 - 2.5 mm	11 cm	0.5 ml	50
1159.03	1.5 - 2.5 mm	30 cm	0.8 ml	50
1159.05	1.5 - 2.5 mm	50 cm	1.1 ml	50
1159.10	1.5 - 2.5 mm	100 cm	2.0 ml	50
1159.15	1.5 - 2.5 mm	150 cm	2.8 ml	50
1159.20	1.5 - 2.5 mm	200 cm	3.6 ml	50

biocap

Needle-free drug aspiration; for multiple vial access

- Contains Latex: No
- Contains DEHP: No
- Contains biological or animal-based product: No
- Non-pyrogenic: Yes

Product Code	Cap Ø	Units/Box
821.21	14 mm	50
821.22	20 mm	50



bionector

Bionector is a multi-purpose, closed, needle-free connection device.

When Bionector is in place, you can infuse, inject sample and change the I.V. tubing without the I.V. circuit to the atmosphere. Bionector's membrane automatically opens the fluid pathway only when a male Luer is disconnected, the membrane automatically seals the fluid pathway.

- Bionector preserves constantly a closed system
- 7 days or 360 connections
- Straight fluid pathway
- Neutral displacement of fluid
- Minimal dead space
- Can be used with power injectors (CT-scan).
Pressure resistance: 350 psi
- MRI conditional
- Contains Latex: No
- Contains DEHP: No
- Contains biological or animal-based product: No
- Non-pyrogenic: Yes




Product Code	Flow rate	Priming volume	Equivalent Gauge G	Number of accesses	Max. duration of use	Length	Diameter	Units/Box
896.038	105 ml/min	0.03 ml	-	360	7 days	25.7 mm	10 mm	50

Ver.2021-04-22

IV Accessories

Ramp of stopcocks (1 way)

Product Code	Description	Extension tube			Lateral accesses	Units/Box
		Ext. Ø	Int. Ø	Length	number	
873.24	Ramp of stopcocks (1 way)	4.0mm	2.5 mm	200 cm	4	15



Ramp of Stopcocks (1 way): female Luer-lock/female Luer-lock

Stopcock

Product Code	Description	Units/Box
872.10	Lockable one-way stopcock	50



Lockable one-way stopcock

Connectors

Product Code	Description	Units/Box
892.00	Female Luer-lock to female Luer-lock connector	30
893.00	Male Luer Lock to Male Luer Lock Connector	30



Female Luer Lock Connector



Male Luer-lock to male Luer-lock connector

Extension tubes (single lumen)

Product Code	Description	Extension tube			Units/Box
		Ext. Ø	Int. Ø	Length	
220.050	Lectroflex extension tube	4.0 mm	2.5 mm	50 cm	25
221.100	Male Luer-lock to male Luer-lock extension tube	4.0 mm	2.5 mm	100 cm	25



Lectroflex



Male Luer-lock to male Luer-lock extension tube

Rapid Infusion

INFU-SURG® PRESSURE INFUSION BAG

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE

Ethox



Three-way stopcock ensures precise control of pressure while the large oval bulb allows for quick bladder inflation, all with ONE hand



Solution bag hook enables easy bag changing without having to remove the Infu-Surg from the IV pole; keeps fluid bag in proper position while in use



360° view of the color-coded gauge facilitates accurate pressure monitoring while the safety pressure relief valve prevents over-inflation

- Ideal for intravenous and intra-arterial infusion treatment including A-line pressure monitoring
- Pressure relief valve prevents over inflation (330 mm Hg pressure relief)
- Large, oval shaped bulb allows for quick and easy inflation of bladder
- Single-handed inflation and deflation design makes it easy to use and requires minimal training
- Dedicated fluid bag hook enables efficient loading and unloading of the system without having to remove the infusion bag from the IV pole
- Reinforced hanging loop allows for quick and easy hanging on IV pole
- Color-coded gauge makes for accurate pressure monitoring (150-300 mm Hg)
- Three-way stopcock ensures precise control of pressure
- Piston-style gauge allows for 360° visibility

ITEM#	DESCRIPTION	VOLUME	PK
4005H	Infu-Surg® with Hook	500 mL	5
4010H	Infu-Surg® with Hook	1000 mL	5
4030	Infu-Surg®	3000 mL	5

Ver.2021-04-22

DISTRIBUTED BY:
 ADVANCED SURGI-PHARM INC.
 www.advancedsurgipharm.com
 Tel.1-800-661-5432



Rapid Infusion

INFU-SURG® CLARITY PRESSURE INFUSION BAG

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE

Ethox



Clear panel enables full visibility of solution bag to monitor fluid level and easily scan patient barcode



- Transparent panel for clear visibility of solution bag
- Ideal for intravenous and intra-arterial infusion treatment including A-line pressure monitoring
- Dedicated hook enables quick change of solution bag while providing proper fluid bag positioning in pocket
- Pressure relief valve prevents over-inflation (330 mm Hg)
- Color-coded gauge provides accurate pressure monitoring (150-300 mm Hg)
- Single-handed inflation and deflation with oval-shaped bulb makes it easy to use and requires minimal training
- Welded panel and hook construction are more durable than sewn designs to withstand repeated use
- Reinforced hanging loop allows for quick and easy hanging on IV pole
- Three-way stopcock ensures precise control of pressure
- Piston-style gauge allows for 360° visibility



Three-way stopcock ensures precise control of pressure while the large oval bulb allows for quick bladder inflation, all with ONE hand



Solution bag hook enables easy bag changing without having to remove the Infu-Surg from the IV pole; keeps fluid bag in proper position while in use



360° view of the color-coded gauge facilitates accurate pressure monitoring while the safety pressure relief valve prevents over-inflation

ITEM	DESCRIPTION	VOLUME	PK
4050	Infu-Surg® Clarity	500 mL	5
4100	Infu-Surg® Clarity	1000 mL	25

Ver.2021-04-22



DISTRIBUTED BY:
ADVANCED SURGI-PHARM INC.
www.advancedsurgipharm.com
Tel.1-800-661-5432

ASEPT® Catheter and Insertion Kits

The ASEPT® Pleural and Peritoneal Drainage Systems are designed for long-term, intermittent drainage of symptomatic, recurrent fluid accumulation.



- | | | | |
|---|------------------------------------|------------------------|----------------------|
| 1 ASEPT® Catheter | 6 16G x 1.0" (25.4mm) Needle | 11 Foam Catheter Pad | 16 10mL Syringe |
| 2 ASEPT® Drainage Line | 7 25G x 1.0" (25.4mm) Needle | 12 Introducer Needle | 17 J-Tip Guidewire |
| 3 5-in-1 Drain Line Adapter | 8 22G x 1.5" (38.1mm) Needle | 13 Tearaway Introducer | 18 Tunneler |
| 4 Chloraprep® 10.5mL Vial | 9 Tegaderm® Self Adhesive Dressing | 14 8F (2.67mm) Dilator | 19 Forceps |
| 5 Synthetic Suture, 2-0 Straight and Synthetic Suture, 3-0 Curved | 10 Gauze Pads, 4"x4" | 15 12F (4.0mm) Dilator | 20 Safety Scalpel |
| | | | 21 Fenestrated Drape |

The primary component of the ASEPT® Drainage system is the catheter, which has an aseptic, self-sealing valve attached that provides a secure and reliable connection. The valve can be connected to an ASEPT® vacuum drainage bottle, drainage bag, or a stand-alone drainage line to connect to other fluid collection equipment.

Features:

- ▶ Secure connection
- ▶ Cap-less and Aseptic design
- ▶ Needle-free and swabable
- ▶ Easy to replace
- ▶ Self-sealing

ASEPT® Catheter Portfolio

Long-Term Intermittent Drainage Solution

The pfm ASEPT™ Catheter System is designed to provide controlled, long-term, intermittent drainage for patients with pleural effusions and malignant ascites – supporting symptom relief in hospital or outpatient settings.

Here are a few key features of the PFM Catheter 15Fr

- Designed for long-term, intermittent drainage of pleural effusions (M7001-CAN) and malignant ascites (M7008-CAN).
- Aseptic, needle-free, swab-able valve helps maintain a secure and hygienic catheter connection.
- Carefully engineered to allow periodic access to the pleural or peritoneal cavity for symptom relief.
- Included in the kit with most items generally required for catheter implantation, streamlining the procedure.
- Works with compatible drainage lines, bottles, or bags for flexible clinical or outpatient use.

Single Catheter Platform for Pleural & Peritoneal Use

Designed for procedural efficiency and flexible patient management.

15Fr Silicone Catheter

Smooth, flexible design for patient comfort and reliable, long-term placement.

M7001-CAN - 66 cm

Designed for pleural effusions.

M7008-CAN - 71 cm

Designed for malignant ascites



Aseptic Needle-Free Swabable Valve

Helps maintain a secure, hygienic connection and reduces infection risk.

Engineered for Intermittent Drainage

Allows periodic access to the pleural or peritoneal cavity for symptom relief.

Complete Implantation Kit

Includes most items generally required for safe and efficient catheter placement.

ASEPT® Roadmap

Integrated Pathway for Long-Term Intermittent Drainage

1 Patient Consultation & Catheter Placement

Patients are scheduled and seen by their surgeon for ASEPT® catheter insertion either Pleural (M7001-CAN) or Peritoneal (M7008-CAN). Patients equipped with an IPC catheter require ongoing drainage as directed by their healthcare provider.




2 pfm Universal Valve Replacement

Many patients have IPC catheters with brand-specific valves that are only compatible with certain bottles. This creates challenges during urgent or emergency drainage situations. pfm offers a universal Luer-Lok valve that can be replaced bedside, allowing compatibility with any Luer-Lok connection (including wall suction), and most importantly, enabling use with pfm drainage bottles for immediate relief.



3 Clinical Drainage - Jet-Vac

Following insertion, patients may require periodic drainage in clinic. pfm Jet-Vac vacuum bottles are available in 600 mL, 1000 mL, and 1500 mL volumes. Options include attached drain lines, valve-only configurations, and sampling ports for patients requiring additional testing.



4 ASEPT® Patient Homecare

For patients draining at home, pfm offers ready-to-use, patient-friendly drainage kits. These include personal-use bottles and all necessary dressing materials in a convenient carrying case.

This solution is ideal for both individual patient use and home healthcare teams conducting home visits.



Designed for long-term, intermittent drainage of symptomatic fluid accumulation.

ET-VA Drainage Bottles

Breakdown

- ▶ Available in three sizes:
 - ▶ 600 mL
 - ▶ 1,000 mL
 - ▶ 1,500 mL
- ▶ Available in two configurations:
 - ▶ Drain line featuring sampling port for easy access to contents of bottle via luer-lock syringe:
 - Connection is made at catheter level
 - Drain line cannot be removed or changed
 - ▶ Luer-lock valve:
 - Connection is made at bottle level
 - Fluid samples can also be drawn with luer-lock syringe connected directly to valve
 - Internal Sterility



ASEPT[®]

Breakdown

- ▶ Available in two sizes:
 - ▶ 600 mL
 - ▶ 1,000 mL
- ▶ Only available with drain line can be used to regulate flow, should patient feel any discomfort during drainage procedure.

Differences between ASEPT[®] and Jet-Vac Drainage Systems



Sampling Port on Jet-Vac allows for safe and easy fluid sample collection.

Pinch Clamp on ASEPT[®] Bottle allows for ease of use with ailing patients.

JET-VAC®

Jet-Vac Centesis Catheters, Jet-Vac Procedure Trays and Jet-Vac Bottles provide a complete drainage system for physicians, nurses, and hospitals. The Jet-Vac drainage bottle is a disposable vacuum device used to aspirate, remove, or sample body fluids.

Jet-Vac Centesis Catheters

Large circular side holes designed to increase fluid drainage

Tapered catheter tip

- Provides a smooth transition
- Helps minimize patient discomfort

HDPE material

- Reduces catheter collapsing
- Promotes a smooth catheter insertion

Clear needle hub with catheter lock

- Increases fluid visibility
- Ergonomic design provides a secure grip
- Prevents needle/catheter separation

Plastic vent plug

- Allows for air to exit
- Provides added safety for physician

Echogenic needle

- Ensures needle visibility during ultrasound placement



Ordering Information

Jet-Vac Centesis Catheters

Item No.	Catheter Size	Length	Needle Gauge	Box Qty
M7511	4F	7 cm	20G	5/Box
M7512	4F	10 cm	20G	5/Box
M7513	5F	7 cm	19G	5/Box
M7514	5F	10 cm	19G	5/Box
M7515	5F	15cm	19G	5/Box

ASEPT® Accessories

pfm medical, inc. is proud to offer a full line of accessories ensuring that the physician or patient will have everything he or she needs to complete the drainage procedure. Our accessory portfolio contains the following items:

ASEPT® Replacement Valve



ASEPT® Drainage Line Set



ASEPT® Procedure Pack



ASEPT® Permanent Drainage Catheter



Universal Connection to Drainage Bottle



Ordering Information

ASEPT® Drainage Bottles and Bags - Patient Home Care

Item No.	Description	Box Qty
M7002	ASEPT® 600mL Drainage Kit (Includes vacuum bottle and Procedure Pack)	10 each
M7000	ASEPT® 1000mL Drainage Kit (Includes vacuum bottle and Procedure Pack)	10 each
M7006	ASEPT® 600mL Select Drainage Kit (Includes low vacuum bottle and Procedure Pack)	10 each
M7007	ASEPT® 1000mL Select Drainage Kit (Includes low vacuum bottle and Procedure Pack)	10 each
M7060	ASEPT® 2000mL Drainage Starter Kit (Includes gravity bag with drainage line - for Peritoneal drainage only)	10 each
M7058	ASEPT® 2000mL Drainage Bag (Includes gravity bag with drainage line - for Peritoneal drainage only)	10 each
M7080	ASEPT® 600mL Drainage Starter Kit (Includes vacuum bottle and Procedure Pack in patient box)	5 each
M7070	ASEPT® 1000mL Bottle (Includes vacuum bottle - does not include Procedure Pack)	10 each
M7072	ASEPT® 600mL Bottle (Includes vacuum bottle - does not include Procedure Pack)	10 each
M7081	ASEPT® 1000mL Drainage Starter Kit (Includes vacuum bottle and Procedure Pack in patient box)	5 each

Ordering Information

ASEPT® Insertion Kits

Item No.	Description	Box Qty
M7001-CAN	ASEPT® Pleural Drainage System (Includes ASEPT® catheter and insertion kit)	1 each
M7008-CAN	ASEPT® Peritoneal Drainage System (Includes ASEPT® catheter and insertion kit)	1 each
M7009	ASEPT® Supplemental Insertion Kit (Includes 7F dilators & 16F valved peelable introducers)	10 each

Jet-Vac Bottles - Clinical Use

Item No.	Item	Description	Box Qty
M7090	Jet-Vac 1000	1,000 mL Evacuated Drainage Bottle in Valve-Only Configuration	10/Case
M7092	Jet-Vac 600	600 mL Evacuated Drainage Bottle in Valve-Only Configuration	10/Case
M7093	Jet-Vac 1500	1,500 mL Evacuated Drainage Bottle in Valve-Only Configuration	10/Case
M7094	Jet-Vac 1000D	1,000 mL Evacuated Drainage Bottle with Drain Line	10/Case
M7095	Jet-Vac 1500D	1,500 mL Evacuated Drainage Bottle with Drain Line	10/Case
M7096	Jet-Vac 600D	600 mL Evacuated Drainage Bottle with Drain Line	10/Case

Accessories

Item No.	Description	Box Qty
M7004	Catheter Replacement Valve	5/Box
M7005	Drainage Line Set (Includes Drainage Line with End Caps, and 5-in-1 Adapter)	10/Case
M7040	Procedure Pack (Traditional Dressing Change)	10/Case
M7097	Jet-Vac Drainage Line Extension (two female luer connectors)	10/Case
M7098	Jet-Vac Large Bore Drainage Line	5/Box
M7099	Jet-Vac Drainage Line	10/Case
M7555	2000cc Drainage Bag	10/Case



Bronchitrac® “L” Set

There is no need to twist the catheter or patient's body. Simply advance the Bronchitrac® “L” through a properly placed endotracheal or tracheostomy tube. The dual-plane, patented curvature of the catheter ensures consistent entry into the left main-stem bronchus.

The time savings and reduced patient trauma provided by the Bronchitrac® “L” make it the obvious choice when suctioning the left-bronchial airway.

- Ideal for suctioning of post-cardiac surgical, neuro and critical care patients
- Can be used for deep-lobe suctioning
- Completely latex-free



Item #	Description
41-7814	21" long, 14 Fr Bronchial Suction Catheter, Slide-cover Connector and (2) Vinyl Gloves, Sterile (20/Box)

TRAKE-fit™ Collar Set

Tracheostomy Tube Holder For Maximum Patient Comfort

- Available in single- or double-collar sets
- Installs in seconds, remains secure
- Maximum comfort with cotton-lined neckband
- No adhesive, no narrow tape
- Completely adjustable
- Narrow tabs available for pediatric and small adult use
- Third-party reimbursable



Item #	Description
43-1850	Double Collar Set (25/Box) Each Set includes: (2) neckbands, (2) wide hook/loop tabs
43-1851	Single Collar Set (25/Box) Each Set includes: (1) neckband, (2) wide hook/loop tabs
43-1862	Double Collar/Narrow Tab Set (25/Box) Each Set includes: (2) neckbands, (2) narrow hook/loop tabs



COMFIT™ Endotracheal Tube Holder

The Tapeless Way to Secure an Endotracheal Tube



COMFIT™ is a reliable, tapeless solution used to secure endotracheal tubes and better manage patient airways. Physicians appreciate COMFIT'S compliance with Advanced Cardiac Life Support (ACLS) Guidelines that recommend endotracheal tube stabilization devices as a substitute for tape. For patients, COMFIT™ provides a comfortable fit without the skin irritation that tape can cause.

Features and Benefits

- Significantly decreases the incidence of self-extubation¹
- Wide cotton-lined neckband minimizes skin irritation, providing maximum patient comfort
- No tape on the tube or the patient
- Completely adjustable
- Minimal plastic loop around the endotracheal tube allows access to the oral cavity
- A low-cost device that requires infrequent changing
- Latex-free product
- Complies with American Heart Association's ACLS Guidelines recommending the use of commercially available endotracheal tube holders instead of tape
- Plastic double securing loop, ½- x 19-inch neckband



1 McAndrew J, RR. Reduction of self extubation using Comfit™ tracheal tube holder. Respiratory Care. November 1994; 39 (11).

Item #	Description
42-1820	Tapeless Endotracheal Tube Holder (25/Box) Each set includes: (1) plastic tube grip, with hook/loop strap, foam foot pads, (2) foam "lollipop" tapes




DISTRIBUTED BY:
ADVANCED SURGI-PHARM INC.
www.advancedsurgipharma.com
Tel.1-800-661-5432



Esophageal Balloon Catheter Set

For Intrathoracic Pressure Monitoring

- Ease of Insertion - The catheter comes with a new, easy-to-remove PTFE-coated stylet to aid the insertion procedure
- Comfort - The 5 French catheter is soft and pliable, offering maximum patient comfort over extended time periods
- Ease of Positioning - The new proximal Y-connector allows for precise balloon positioning with the stylet in place
- Depth Markings - The catheter is radio opaque with numerical depth markings to facilitate positioning
- Ease of Connection - The catheter has a standard female Luer connector
- Several Applications - Ideal for monitoring ventilator and sleep apnea patients and calculating lung compliance



Item #	Description
47-9005	5 Fr Catheter, latex-free Balloon, Stylet, 36" Pressure Monitoring Extension, 3-way Stopcock (5/Box)

CA  CE

Gastric Lavage

TUM-E-VAC[®] GASTRIC LAVAGE KIT

LATEX FREE | SINGLE USE | DISPOSABLE

Ethox



- Pre-assembled kits are quick and easy to use, reducing the time for lavage procedures
- Promotes quick, efficient lavage of stomach contents
- Closed system helps prevent exposure to stomach contents and charcoal solutions
- Available with activated charcoal suspension (with or without Sorbitol)
- Twist-Lok™ expressing key aids in charcoal administration
- Easy-to-fill top port
- Available with Universal Solution Spike (for use with pre-filled solution containers)
- No clean up or disinfecting procedures necessary
- All kits come standard with Bite Block

WITHOUT CHARCOAL

ITEM#	DESCRIPTION	PK
2054	Tum-E-Vac + 3L Solution Bag + 32Fr Lavage Tubing	10
2055	Tum-E-Vac + 3L Solution Bag	10
2075	Tum-E-Vac + Universal Solution Spike	10

WITH ACTIVATED CHARCOAL


ITEM#	DESCRIPTION	PK
2057	Tum-E-Vac + 32Fr Lavage Tubing + (2) 25g Activated Charcoal	10
2058	Tum-E-Vac + 32Fr Lavage Tubing + (2) 25g Activated Charcoal with Sorbitol	10
2059	Tum-E-Vac + 32Fr Lavage Tubing + (1) 25g Activated Charcoal with Sorbitol + (1) 25g Activated Charcoal	10
2081	Tum-E-Vac + Universal Solution Spike + (2) 25g Activated Charcoal with Sorbitol	10

LAVAGE TUBING

LATEX FREE | SINGLE USE | DISPOSABLE

Ethox

Radiopaque strip confirms proper tube placement



Reference marks at 18", 20" and 22" from distal end




- Open bullet tip design with multiple eyelets to create free-flowing movement with low resistance
- Soft clear tube with smooth tip
- Suction connector on proximal end
- Radiopaque stripe
- All tubing is 48" long with reference marks at 18", 20" and 22" from distal end
- Individually packaged

ITEM#	SIZE	DESCRIPTION	PK
2124	24Fr	48" Lavage Tubing	10
2128	28Fr	48" Lavage Tubing	10
2132	32Fr	48" Lavage Tubing	10
2136	36Fr	48" Lavage Tubing	10

Irrijet®

Splatter-Gard™

Simple, one-hand, thumb-controlled pressure provides an appropriately forceful stream (8 psi) directly to the desired area (or a gentler stream if needed). The soft flexible tip has an inside diameter equivalent to an 18-19 gauge needle.



- Splatter-Gard™ prevents contamination from body secretions – new concave design provides maximum splash-back protection
- Soft, flexible plastic tip eliminates needle-stick for user, minimizes trauma to patient
- Irrigant is delivered with force equivalent up to 8 psi
- Latex-free
- Sterile


Item #	Description
55-0004	Irrijet GP (5/Box)
55-0504	Multi-Use Irrijet (5/Box)
55-2004	Single-Use Irrijet Z (5/Box)
55-2104	Single-Use Irrijet (5/Box)
55-2134	Single-Use Irrijet (5/Box)
55-6504	Irrijet DS (5/Box)
59-5014	Infuset GP (5/Box)



Irrijet GP

55-0004


Spring-loaded syringe connected to a double check-valve, through a side arm and 6-ft. length of tubing to a spike for entry into a bag of irrigating solution. Double check-valve prevents retrograde flow. Alternate movement of the syringe plunger delivers irrigating solution as needed. Syringe Capacity: 12cc.



Multi-Use Irrijet

55-0504


Same basic assembly as 55-0004, plus one package of five individually sealed jet tips.



Single-Use Irrijet Z

55-2004

Same basic assembly as 55-0004, plus one Zerowet Splashshield®.




Single-Use Irrijet

55-2104

Same basic assembly as 55-0004, plus includes permanent flexible plastic jet tip and Splatter-Gard™.

55-2134

Same basic assembly as 55-2104, but with 35 cc capacity syringe.




Wound Irrigation

COMBIPORT® + COMBIGUARD II™

LATEX FREE | STERILE | SINGLE PATIENT USE | DISPOSABLE

Ethox



The Combiguard II™ attaches to a syringe (or Combiport) to assist in wound cleansing

Fluid flows through angled port, directed away from cleansed area and toward the outlet ports

Working together, a saline bag uses a Combiport® attached to a Combiguard II™ to effectively cleanse the wound without splash-back on the caregiver



COMBIPORT® Wound Irrigation Device

COMBIGUARD II™ Irrigation Splash Guard

COMBIPORT® WOUND IRRIGATION DEVICE

- Attaches to a standard saline bag to administer a high volume of fluid in a short period of time
- Irrigation helps remove foreign bodies, reduce pathogen numbers and moisten wound tissue
- Patented spike attaches to standard 500 mL or 1000 mL saline bag
- For use only with Combiguard Irrigation Splash Guard (not sold separately)
- Removable cap for easy stop/start during cleansing
- Easy to assemble and use

COMBIGUARD II™ IRRIGATION SPLASH GUARD

- Egg-shaped cup covers wound site to prevent irrigation splash-back, helping to protect clinician
- Dual outlet ports enable contaminated fluids to safely escape
- Fluid delivery inlet delivers a stream equivalent to an 18 g needle
- Angled delivery inlet directs fluid away from cleansed area toward outlet ports
- Rounded edges provide a smooth surface against the skin
- Fits on 30-40 cc luer-tipped syringe
- Can be combined with Combiport Wound Irrigation Device

STERILE

ITEM#	DESCRIPTION	PK
10222	Combiguard II™ Irrigation Splash Guard	50
10100	Combiguard II™ + Combiport®	50

NON-STERILE

ITEM#	DESCRIPTION	PK
10222NS	Combiguard II™ Irrigation Splash Guard Non-Sterile	500

SMOOTH MOVER



PRODUCT DESCRIPTION

An economical ultra-smooth polyethylene patented board especially designed to furnish expert patient transference with minimal trauma to the patient. The Smooth Mover is capable of safely moving extremely large patients with minimal staffing requirements. It is radiotranslucent; it can be left in place during x-ray procedures. Is lightweight yet durable. Perfect for daily patient moving needs: EMS or hospital, nursing home or clinic applications. Easy on ...easy off. Yes, even one person can move a large patient!

Optional Accessory:

Pull Straps... loops straps slip through hand-holes enabling easy board pick-up on hard surfaces tables (x-ray, etc.) and minimize reach required when transferring patients from bed to stretcher to table during routine moving.

PRODUCT ORDERING INFORMATION

Hospital Shifter Model SM3


Is widely used for patient shifting, bed to stretcher to table, etc..
 Dimensions: 3/16" x 22" x 72"
 Weight: 9 lbs.

Ambulance Patient Shifter Model SM40

Is used as an aid in moving patients from vehicle to stretcher or backboard.
 Dimensions: 3/16" x 18" x 48"
 Weight: 5½ lbs.

ADJUSTABLE PATIENT POSITIONING DEVICES

Patented, inflatable, adjustable OR Comfort™ Patient Positioning Devices provide patient stability and comfort, as well as optimal surgical positioning. Devices can be adjusted throughout a procedure without having to manually lift or move the patient.



SHOULDER-FLOAT™ ADJUSTABLE AXILLARY SUPPORT DEVICE

- Relieves patient's shoulder discomfort during surgical procedures in the lateral decubitus position. Minimizes pressure on shoulder and head to help reduce post surgical pain.
- Useful during total hip replacement, thoracotomy, shoulder arthroscopy and kidney surgery. It can benefit mildly sedated patients undergoing epidural anesthesia and patients undergoing lengthy procedures.

Shoulder-Float Device			
PART #	PACKAGING DESCRIPTION	PART #	PACKAGING DESCRIPTION
SF05	5 per BX	SF10	10 per BX



PELVIC-TILT™ ADJUSTABLE DISPLACEMENT DEVICE

- Designed to lift a patient's hip between 20-40 degrees during surgical procedures. Includes a non-inflatable flap to place underneath the patient to prevent a lateral shifting of the device.
- Useful in obstetrics to avoid compression of the vena cava and/or aorta during cesarean sections and deliveries in the supine position. Practical during kidney surgery and surgery on the side of the leg or ankle.

Pelvic-Tilt Device			
PART #	PACKAGING DESCRIPTION	PART #	PACKAGING DESCRIPTION
PT05	5 per BX	PT10	10 per BX



DELGADO™ INFLATABLE POST

- Slides over the stabilizer post of the operating room table to provide an adjustable "pillow" between the post and the patient.
- Maximizes the stability of the patient and distributes the pressure, minimizing the likelihood of tissue damage or discomfort after surgery.

Delgado™ Inflatable Post			
PART #	PACKAGING DESCRIPTION	PART #	PACKAGING DESCRIPTION
DG05	5 per BX	DG10	10 per BX

REQUEST YOUR COPY OF THE COMPLETE IMPERIAL SURGICAL CATALOGUE TODAY!



THE IMPERIAL SURGICAL PORTFOLIO INCLUDES:

- Warming Cabinets
- IV Poles
- Carts & Cabinets
- Tables
- Maternity Bassinets
- Stools & Chairs
- Imperial Ware, Stands & Receptacles
- Specialized Equipment
- Sinks & Sink Systems



Learn more and view the complete digital catalogue at

www.imperialsurgical.com



Ver. 2025-10-14



2165 46th Ave., Lachine, Quebec H8T 2P1
Phone: +1.514.631.7988 | Toll Free: 800.661.5432 | Fax: +1.514.631.9083
Email: info@surgmed.com | Website: www.advancedsurgipharm.com