

NEW!

BATRIK

# TAG-IT™ | SCOPE

Batrik Scope Tags are tear-resistant, double-sided tags that will prove indispensable throughout the reprocessing journey. Our tags make it easy to record important information, including processing dates, with a user-friendly feature that allows you to select the month, day, and date. Crafted from waterproof synthetic paper, these tags feature a self-looping design, making it effortless to attach them to the scope without the need for additional fastening materials. Elevate your reprocessing protocols with Batrik!

## KEY FEATURES

- The paper, tie wraps, and smudge-proof ink are designed to withstand high temperatures of up to 273.2°F (134°C) and are compatible with Gas plasma sterilization (STERRAD®100 NX Sterilization System).
- Tags include the month, day, and date selection on one side, making it convenient to mark the processing date accurately



## COLOR-CODED AND CUSTOMIZABLE

- Available in seven colors
- The tags facilitate color coding for easy identification

## SELF LOOPING

- The self-looping design allows for easy attachment to the scope
- Eliminating the need for additional fastening materials

## SINGLE-USE AND LATEX FREE

- Batrik Scope Tags are designed for single use and are latex-free
- Promoting hygiene and reducing the risk of cross-contamination

## PRODUCT SELECTION

ITEM #	COLOR	SIZE	PACKAGING
TI-ST05	Yellow	11.8" x 1.55" / 33.8 cm x 4 cm	100 / PKG
TI-ST02	Blue	11.8" x 1.55" / 33.8 cm x 4 cm	100 / PKG
TI-ST03	Green	11.8" x 1.55" / 33.8 cm x 4 cm	100 / PKG
TI-ST07	Orange	11.8" x 1.55" / 33.8 cm x 4 cm	100 / PKG
TI-ST06	Purple	11.8" x 1.55" / 33.8 cm x 4 cm	100 / PKG
TI-ST04	Red	11.8" x 1.55" / 33.8 cm x 4 cm	100 / PKG
TI-ST01	White	11.8" x 1.55" / 33.8 cm x 4 cm	100 / PKG

All information tags are single-use. PACKAGING: 1 pack of 100 tags and 1 permanent marker.

