


Your First ™

ENZYMATIC BEDSIDE SPONGE KITS

RightSwipe™

Introducing RightSwipe™, our Enzymatic Bedside Sponge Kits now feature sponges treated with an enhanced formula for a more effective pre-clean or surface clean solution. This unique three-enzyme blend ensures confident and efficient cleaning, making RightSwipe™ the perfect match for your endoscope care needs. The tubular sponge design remains ideal for all endoscopes, providing versatility alongside the advanced benefits of our enzymatic solution.



FLAT SPONGE



TUBULAR SPONGE



FREESTANDING POUCH



PLASTIC BAG

PRODUCT SELECTION

ITEM #	DESCRIPTION	PACKAGING
17-111	PRE-TREATED TUBULAR SPONGE Add 100-500 ml of water to activate detergent solution. Comes in plastic bag.	25 / Box (4 Boxes / Case)
17-211	PRE-TREATED FLAT SPONGE Add 100-500 ml of water to activate detergent solution. Comes in plastic bag.	50 / Box (4 Boxes / Case)
17-122	PRE-TREATED TUBULAR SPONGE Comes in a freestanding pouch to use as a basin. Simply add up to 500 ml of water.	15 / Box (4 Boxes / Case)
17-222	FLAT SPONGE Comes in a freestanding pouch to use as a basin. Simply add up to 500 ml of water.	25 / Box (4 Boxes / Case)
17-133	READY TO USE TUBULAR SPONGE Pre-treated in 500 ml of an enzymatic solution in freestanding pouch.	24 / Case
17-233	READY TO USE FLAT SPONGE Comes bathing in 500 ml of enzymatic solution in freestanding pouch.	24 / Case
17-144	READY TO USE DUAL PURPOSE TUBULAR SPONGE Comes in a freestanding pouch soaked in 100 ml of enzymatic solution.	12 / Box (4 Boxes / Case)
17-244	READY TO USE DUAL PURPOSE FLAT SPONGE Comes in a freestanding pouch soaked in 100 ml of enzymatic solution.	12 / Box (4 Boxes / Case)
17-155	READY TO USE TUBULAR SPONGE Comes bathing in 220 ml of enzymatic solution in freestanding pouch.	24 / Case
17-255	READY TO USE FLAT SPONGE Comes bathing in 220 ml of enzymatic solution in freestanding pouch.	24 / Case