




# COMPLETE PRODUCT CATALOGUE

Established in 1997. Proudly Canadian. 🇨🇦  
**Sold in over 85 countries worldwide.**



Proudly Canadian and established in 1997, we specialize in the manufacturing of enzymatic solutions offered in different formats as well as the fabrication of silicone and plastic products. Our products are sold in over 85 countries worldwide and used primarily in surgery, endoscopy, and central sterilization areas.

Our in-house controlled environment production facility in Montreal, Canada adheres to strictly enforced quality and regulatory control systems. We are registered as a device manufacturing establishment with the United States FDA and we are in full compliance with CE and Health Canada Regulations. We continually invest in upgrading our plant, machinery, and processes to ensure high quality products at high volumes. Most products are available off-the-shelf.

At Batrik, we are committed to innovation and providing our customers with the finest quality products backed by world-class customer care.

## SURGICAL PRODUCTS & GENERAL ACCESSORIES

04	SUTURE CLAMP TAGS	
05	VASCULAR SILICONE TIES	
06	CLAMP JAW COVERS	
07	G.O.L.F.F.™ ANTI-FOG SOLUTION	
08	SKINS™ <a href="#">Surgical MULTI-TIP Markers</a>	
09	PANDAPATCH™ EYE SAFETY TAPE	NEW!
11	MAGNETIC MAT	
12	PIN GUARDS	
13	CAUTERY TIP CLEANER	
14	SAFE-SOAK™ SPAC & MATS	
15	BLUE CLOUD-GEL™ ANTI-FATIGUE MATS	
16	BLUE-CLOUD NIMBUS™ ANTI-FATIGUE MATS	
17	STACKABLE STEP STOOLS	
18	TOURNY™+ TOURNIQUETS	
20	SHOOZCOVERS™ SOFT LINE™	
22	SHOOZCOVERS™ COATED LINE™	NEW!
24	SHOOZCOVERS™ TOUGH LINE™	
26	SHOOZCOVERS™ GRIPPEEZ™	

## ENDOSCOPY PRODUCTS

28	BITE BLOCKS <a href="#">Endoscopy Bo2™ / Bitestrong™</a>
29	ENDOSCOPY BITESTRONG™ BITE BLOCKS
30	COUVRETTE™ VALVE COVERS
31	DAYBRUSH™ CHANNEL BRUSHES
32	ARC™ CHANNEL BRUSHES
35	VALVE / CONTROL BRUSHES
36	V-BRUSHES


## INFECTION CONTROL & STERILE PROCESSING


37	TAG-IT™	NEW!
42	SILICONE PIN MAT	
43	CORNER POCKETS™ & SOFT POCKETS™	
44	LILA'S INSTRUMENT STAIN REMOVER	
45	INSTRUMENT PROTECTORS	
48	PrAmLi Power <sup>3</sup> LUBE™ INSTRUMENT LUBRICANT	
49	S.T.I.C.S.	
50	ID TAPE	
51	TAPE-B-GONE™ TAPE REMOVER	


52 PADLOCKS SECURITY TAGS


56 SAVEFACE™ DRAPED SHIELDS


58 SCRUB PAD


59 **Your First** ™ ENZYMATIC & NON-ENZYMATIC CLEANERS


60 **Your First** ™ ENVIRO-SACHET™


61 **Your First** ™ FOAMING SPRAY

62 **Your First** ™ F1RST WIPES

63 **Your First** ™ DRY SPONGES

64 **Your First** ™ RIGHTSWIPE™ ENZYMATIC SPONGES  
BEDSIDE SPONGE KITS **NEW!**

65 **Your First** ™ RIGHTSWIPE™ NON-ENZYMATIC  
DETERGENT BEDSIDE SPONGE KITS **NEW!**

66 **Your First** ™ TRANSPORT BAGS FOR ENDOSCOPE  
AND MEDICAL INSTRUMENTS

## BRUSHES

67	ARC™ WIRE STYLE BRUSHES	
71	DOUBLE ARC™ WIRE BRUSHES	
	<a href="#">Asymmetrics / Symmetric</a>	
74	IN-BRUSH™ ANTIMICROBIAL TWISTED WIRE BRUSHES	
76	STAINLESS STEEL ARC™ MICRO BRUSHES	
77	INSTRUMENT CLEANING FLAT BRUSHES	
78	PRECISION BRUSH	
79	TRIFECTA™ BRUSH	NEW!
80	ALL PURPOSE REUSABLE BRUSHES	
81	ACETABULAR REAMER BRUSH	
82	ORTHOPEDIC RASPS CLEANING BRUSH	
83	BOTTLE & DRAIN CLEANING BRUSHES	
84	AUTOClave CLEANING BRUSH KITS	

## IMPERIAL WARE, STANDS & RECEPTACLES



85	BOWL STANDS & RECEPTACLES
86	GRADUATED MEASURING
87	WASH BASINS & BOWLS
88	INSTRUMENT TRAYS
89	CATHETER & INSTRUMENT TRAYS
91	URINALS & BEDPANS
93	JARS & COVERS
94	PETRI DISHES
95	EMESIS BASINS
96	CUPS, TUMBLERS, AND PITCHERS
98	ANCILLARY ITEMS

# BATRIK SUTURE CLAMP TAGS

Non-toxic, radiopaque, and latex free, our Suture Clamp Tags are designed to protect and safely hold sutures without any damage.

Additionally, its foam base, affixed to a sterile area with the clamps, allows for easy insertion.

- Packed in a convenient foam base with easy peel-off adhesive backing, even while wearing gloves
- Available in sterile and non-sterile
- Available in 2, 3, or 5 pairs
- Pre-counted



OUR "EASY PEEL CORNER" IS A FEATURE UNIQUE TO BATRIK

## PRODUCT SELECTION

ITEM # STERILE	ITEM# NON-STERILE	COLOR	DIMENSIONS	PACKAGING
<b>11-1002</b>	<b>911-1002</b>	Yellow	Standard ø (0.06" x 0.5" / 0.15 cm x 1.27 cm)	2 pairs / foam base
<b>11-1003</b>	<b>911-1003</b>	Yellow	Standard ø (0.06" x 0.5" L / 0.15 cm x 1.27 cm)	3 pairs / foam base
<b>11-1005</b>	<b>911-1005</b>	Yellow	Standard ø (0.06" x 0.5" L / 0.15 cm x 1.27 cm)	5 pairs / foam base
<b>11-1025</b>	<b>911-1025</b>	Assorted	Standard ø (0.06" x 0.5" L / 0.15 cm x 1.27 cm)	5 pairs / foam base
<b>11-1045</b>	<b>911-1045</b>	White	Standard ø (0.06" x 0.5" L / 0.15 cm x 1.27 cm)	5 pairs / foam base
<b>11-1065</b>	<b>911-1065</b>	Blue	Standard ø (0.06" x 0.5" L / 0.15 cm x 1.27 cm)	5 pairs / foam base
<b>11-1085</b>	<b>911-1085</b>	Red	Standard ø (0.06" x 0.5" L / 0.15 cm x 1.27 cm)	5 pairs / foam base
<b>11-3002</b>	<b>911-3002</b>	Pink	Large ø (0.09" x 0.75" L / 0.23 cm x 1.91 cm)	2 pairs / foam base
<b>11-3003</b>	<b>911-3003</b>	Pink	Large ø (0.09" x 0.75" L / 0.23 cm x 1.91 cm)	3 pairs / foam base
<b>11-3005</b>	<b>911-3005</b>	Pink	Large ø (0.09" x 0.75" L / 0.23 cm x 1.91 cm)	5 pairs / foam base


STERILE PACKAGING: 1 Foam Base / Pouch – 5 Pouches / Box

NON-STERILE PACKAGING: Sold per foam base

BATRIK

# VASCULAR SILICONE TIES

AVAILABLE IN 4 SIZES  
AND 4 COLORS




## PRODUCT SELECTION

ITEM #	DIMENSIONS	COLOR
10-1000	0.075" X 16" (0.2 cm x 40 cm)	Mini-Red
10-1001	0.075" X 16" (0.2 cm x 40 cm)	Mini-Blue
10-1002	0.075" X 16" (0.2 cm x 40 cm)	Mini-Yellow
10-1003	0.075" X 16" (0.2 cm x 40 cm)	Mini-White
10-1010	0.120" x 16" (0.3 cm x 40 cm)	Maxi-Red
10-1011	0.120" x 16" (0.3 cm x 40 cm)	Maxi-Blue
10-1012	0.120" x 16" (0.3 cm x 40 cm)	Maxi-Yellow
10-1013	0.120" x 16" (0.3 cm x 40 cm)	Maxi-White
10-1101	0.021" x 16" (0.5 cm x 40 cm)	Extra-Mini-Blue
10-1111	0.200" x 22" (0.5 cm x 56 cm)	Extra-Maxi-Blue
10-1113	0.200" x 22" (0.5 cm x 56 cm)	Extra-Maxi-White
10-1800	0.075" x 18" / 0.2 cm x 46 cm	Mini-Red
10-1801	0.075" x 18" / 0.2 cm x 46 cm	Mini-Blue
10-1802	0.075" x 18" / 0.2 cm x 46 cm	Mini-Yellow
10-1803	0.075" x 18" / 0.2 cm x 46 cm	Mini-White
10-1810	0.12" x 18" / 0.3 cm x 46 cm	Maxi-Red
10-1811	0.12" x 18" / 0.3 cm x 46 cm	Maxi-Blue
10-1812	0.12" x 18" / 0.3 cm x 46 cm	Maxi-Yellow
10-1813	0.12" x 18" / 0.3 cm x 46 cm	Maxi-White

**PACKAGING:** 2 ties per sterile pouch – 10 pouches per box

Our Vascular Silicone Ties permit safe and smooth retraction, occlusion and identification of vessels during surgical procedures.

Unlike clamps and retractors which can pinch the vessels flat, or even crack them, the BATRIK elliptical shape allows for the gentle handling of arteries, veins, nerves, tendons, spermatic cords, ureters, or whenever delicate retraction / occlusion is required.



- Manufactured of medical grade silicone
- Do not twist on themselves when unpacked
- Easy to separate, eliminating the need the need for
- Sterile, radiopaque, and latex free cardboard inserts

## PRODUCT SELECTION

ITEM #	DIMENSIONS	COLOR
10-1831	0.04" x 18" / 0.1 cm x 46 cm	Extra-Mini-Blue
10-1841	0.20" x 18" / 0.5 cm x 46 cm	Extra-Maxi-Blue
10-1843	0.20" x 18" / 0.5 cm x 46 cm	Extra-Maxi-White
10-3010	0.075" x 32" (0.2 cm x 81 cm)	Maxi-Red
10-3202	0.075" x 32" (0.2 cm x 81 cm)	Mini-Yellow
10-3200	0.075" x 32" (0.2 cm x 81 cm)	Mini-Red
10-3201	0.075" x 32" (0.2 cm x 81 cm)	Mini-Blue
10-3202	0.075" x 32" (0.2 cm x 81 cm)	Mini-Yellow
10-3203	0.075" x 32" (0.2 cm x 81 cm)	Mini-White
10-3210	0.120" x 32" (0.3 cm x 81 cm)	Maxi-Red
10-3211	0.120" x 32" (0.3 cm x 81 cm)	Maxi-Blue
10-3212	0.120" x 32" (0.3 cm x 81 cm)	Maxi-Yellow
10-3213	0.120" x 32" (0.3 cm x 81 cm)	Maxi-White

**PACKAGING:** 2 ties per sterile pouch – 10 pouches per box

# BATRIK CLAMP JAW COVERS

Manufactured from medical grade silicone, our Clamp Jaw Covers reduce trauma to vessels requiring occlusion, particularly when a hemostatic clamp is utilized. The thread-like indents allow for gentle, constant and non-slip tension.

GOOD TO KNOW: OUR CLAMP JAW COVERS CONVERT ALMOST ANY INSTRUMENT TO ATRAUMATIC



- Gamma sterilized
- Available in 3 colors and 3 widths to fit most clamp sizes
- Radiopaque and latex free

## PRODUCT SELECTION

ITEM # STERILE	ITEM # NON-STERILE	DIMENSIONS	COLOR
<b>12-1000</b>	<b>912-1000</b>	0.1" X 5" / 0.25 cm X 12.7 cm	Red
<b>12-1001</b>	<b>912-1001</b>	0.2" X 5" / 0.51 cm X 12.7 cm	White
<b>12-1002</b>	<b>912-1002</b>	0.3" X 5" / 0.76 cm X 12.7 cm	Blue

STERILE PACKAGING: 2 clamps per sterile pouch - 10 pouches per box

NON-STERILE PACKAGING: 20 clamps per bag (also available 100 units / bag)

# BATRIK G.O.L.F.F.™ ANTI-FOG SOLUTION

Our G.O.L.F.F.™ Anti-Fog Solution solves the fogging problem often encountered during Laparoscopic surgical procedures.

- Prevents the occurrence of fog on lenses
- Ensures clearer vision and avoids unnecessary interruptions during surgery
- Set includes an easy-squeeze 6 mL bottle and an adhesive radiopaque sponge for secure placement
- OEM (available)
- Sterile
- Latex free



## PRODUCT SELECTION

ITEM #	DESCRIPTION	PACKAGING	CASE DIMENSIONS
<b>G.O.L.F.F.</b>	Anti-Fog Solution & Adhesive Sponge set	20 / Box	20" x 16"x 18"

# BATRIK SKINS™ SURGICAL MULTI-TIP SKIN MARKERS

Our SKINS™ Surgical Multi-Tip Skin Markers are designed with a beveled tip to enable the creation of fine or bold tracing lines by holding the pen at different angles.



- Available with or without 6 labels and a flexible 6" calibrated (in/mm) ruler
- Available in non-sterile for pre-operative use, or sterile for use in the operating room
- Surgical prep resistant
- Made of lasting gentian violet ink that is non-toxic and non-tattooing
- Latex free



## PRODUCT SELECTION

ITEM #	TYPE	DESCRIPTION	PACKAGING
<b>S-3000</b>	Sterile	Multi-Tip fine marker with ruler and label set	50 / Box
<b>S-3001</b>	Sterile	Multi-Tip fine marker	50 / Box
<b>S-30005</b>	Sterile	Multi-Tip fine marker with ruler and label set (Volume Packaging)	500 / Box
<b>S-30015</b>	Sterile	Multi-Tip fine marker (Volume Packaging)	500 / Box
<b>S-93000</b>	Non-Sterile	Multi-Tip fine marker with ruler and label set	1000 / Box
<b>S-93001</b>	Non-Sterile	Multi-Tip fine marker	1000 / Box
<b>S-94000</b>	Non-Sterile	Ruler	100 / Bag
<b>S-94001</b>	Non-Sterile	Ruler (Bulk)	20,000 / Box
<b>S-95000</b>	Non-Sterile	Set of 3 labels	100 / Bag

SAVE by ordering our SKINS™ MARKERS in volume packaging. Cost effective and reduces waste.

**NEW!**

## Introducing PandaPatch™: Innovative Adhesive Eye Protection for Safer Patient Care

PandaPatch™ offers an effective alternative to traditional tape for protecting patients' eyes during surgery. These adhesive eye protectors are designed to reduce the risk of corneal injuries and ophthalmic complications for patients under general anesthesia.

With PandaPatch™, you ensure superior patient outcome by preventing avoidable eye damage associated with the use of standard surgical tape.

### KEY FEATURES:

- **Ergonomic Design:** Contours comfortably around the eyes.
- **Post-Surgery Friendly:** Particularly effective for children, the elderly, and individuals with sensitive or delicate skin.
- **Latex-Free & Hypoallergenic Adhesive:** Made from medical-grade material, gentle on the skin and designed to minimize irritation.
- **Non-Adhesive Colored Tabs:** Allowing quick and hassle free application and removal, facilitating effortless monitoring of pupil dilation during surgery.



**70-04ADU**  
Adult Sensitive Eye Protector

### HERE'S WHY IT MATTERS:

- During prolonged surgery under general anesthesia, **up to 59%<sup>1</sup>** of patients may not fully close their eyes, risking corneal abrasions.
- **PandaPatch™** gentle adhesive eye protectors ensure eyelids remain safely and comfortably closed, minimizing the risks typically associated with the use of standard surgical tape.



#### REFERENCES


1. Anson, Jonathan, M.D. "Perioperative Corneal Abrasions: Etiology, Prevention, and Management." Clinical Updates, Penn State University College of Medicine, 2021.




PANDAPATCH™ RANGE:

Three adhesive eye protectors to accomodate different age groups and skin sensitivities. Providing optimal protection for every patient.

- Adult Extra Sensitive Adhesive (Blue Tabs)
- Adult Adhesive (Red Tabs)
- Pediatric Adhesive (Red Tabs)



70-04ADU  
Adult Adhesive Eye Protector




70-04PED  
Pediatric Adhesive Eye Protector

PRODUCT SELECTION

ITEM #	DESCRIPTION	DIMENSIONS	PACKAGING
70-02ADU	ADULT EXTRA SENSITIVE ADHESIVE EYE PROTECTORS (BLUE)	2.15" x 1.35"	50 / Box
70-04ADU	ADULT ADHESIVE EYE PROTECTORS (RED)	2.15" x 1.35"	50 / Box
70-04PED	PEDIATRIC ADHESIVE EYE PROTECTORS (RED)	1.85" x 0.93"	50 / Box

# BATRIK MAGNETIC MAT

Our Magnetic Mat provides a sterile, neutral field for the safe placement and hands-free transfer of instruments. It is latex free, reusable, and autoclavable (up to 134° C).



## PRODUCT SELECTION


ITEM #	DESCRIPTION	DIMENSIONS	PACKAGING
<b>MAG-02-1</b>	Magnetic Mat	12" x 15" / 30 cm x 38 cm	1 Unit
<b>MAG-02-5</b>	Magnetic Mat (Volume Packaging)	12" x 15" / 30 cm x 38 cm	5 Units

\*Independent laboratory tested to validate pre-vacuum steam sterilization efficacy for reprocessing of medical devices in a healthcare setting

# BATRIK PIN GUARDS



- Safely cover the sharp ends of your K-Wire and Steinmann Pins
- Color coded by wire size which eliminates the risk of selecting of the wrong size
- A removable sleeve is provided to avoid injury and facilitate application (sterile only)
- Non-sterile model is available with or without plastic sleeves
- Packaging designed to reduce waste
- Latex free



## PRODUCT SELECTION

ITEM # STERILE PIN GUARDS + SLEEVE	ITEM # NON-STERILE PIN GUARDS + SLEEVE	ITEM # NON-STERILE PIN GUARDS (NO SLEEVE)	COLOR	DIMENSIONS (WIRE SIZE)
<b>14-1028</b>	<b>914-1028 - SLV</b>	<b>914-1028</b>	Yellow	0.028" / 0.71 mm
<b>14-1035</b>	<b>914-1035 - SLV</b>	<b>914-1035</b>	Light Blue	0.035" / 0.89 mm
<b>14-1045</b>	<b>914-1045 - SLV</b>	<b>914-1045</b>	White	0.045" / 1.14 mm
<b>14-1054</b>	<b>914-1054 - SLV</b>	<b>914-1054</b>	Dark Red	0.054" / 1.40 mm
<b>14-1062</b>	<b>914-1062 - SLV</b>	<b>914-1062</b>	Green	0.062" / 1.57 mm
<b>14-1078</b>	<b>914-1078 - SLV</b>	<b>914-1078</b>	Pink	0.078" / 1.98 mm
<b>14-1093</b>	<b>914-1093 - SLV</b>	<b>914-1093</b>	Grey	0.093" / 2.36 mm
<b>14-1109</b>	<b>914-1109 - SLV</b>	<b>914-1109</b>	Beige	0.109" / 2.77 mm
<b>14-1125</b>	<b>914-1125 - SLV</b>	<b>914-1125</b>	Purple	0.125" / 3.17 mm
<b>14-1141</b>	<b>914-1141- SLV</b>	<b>914-1141</b>	Brown	0.141" / 3.58 mm
<b>14-1156</b>	<b>914-1156 - SLV</b>	<b>914-1156</b>	Dark Blue	0.156" / 3.96 mm
<b>14-1188</b>	<b>914-1188 - SLV</b>	<b>914-1188</b>	Red	0.188" / 4.77 mm

STERILE PACKAGING: 80 units/box (2 units per sterile pouch - 40 pouches per box)

NON-STERILE PACKAGING: 100 units/box

## BATRIK CAUTERY TIP CLEANER

Our Cautery Tip Cleaner assists in removing material from cautery tips for both mono-polar and bi-polar electrodes.

- Foam pad with integrated abrasive material
- Adheres to disposable and non-disposable drapes
- Radiopaque
- Available exclusively in Non-Sterile
- Latex free



### PRODUCT SELECTION

ITEM #	DESCRIPTION	DIMENSIONS	PACKAGING
60-10000	Cautery Tip Cleaner Non-Sterile	1.8" x 1.8" / 4.5 cm x 4.5 cm	2,200 /Box

## BATRIK SAFE-SOAK™ SPAC

Our Safe-Soak™ SPAC floor suction disk quickly removes fluid waste from the operating room floor through the application of suction.



- Highly durable: it won't be damaged if stepped on
- Easily maneuverable by foot around operating room floor
- Latex free, disposable, and non-sterile
- The bottom opening permits fluid collection through a highly absorbent material that comes into contact with its surface
- The hose is securely affixed on its side to reduce risk of catching equipment
- Its highly absorbent base material allows for quick fluid wicking and removal when suction is applied

## BATRIK SAFE-SOAK™ ABSORBENT FLOOR MATS

Our Safe-Soak™ Absorbent Floor Mats are designed to absorb a variety of fluids from the operating room floor quickly and efficiently, ensuring the floor beneath remains slip-free, safe, and dry during surgical procedures.

- Highly absorbent
- Facilitates quick and easy clean up
- Latex free, single-use, disposable, and non-sterile
- Conforms to uneven surfaces
- Skid-resistant backing (blue absorbent mat only)



### PRODUCT SELECTION

ITEM #	DESCRIPTION	ABSORBENCY	SIZE	PACKAGING
<b>SS-20-1000</b>	Safe-Soak™ / White extra absorbent mat	4.5 L / mat	28" x 48" (71 cm x 122 cm)	20 / Box
<b>SS-20-2000</b>	Safe-Soak™ / Blue absorbent mat with skid-resistant backer	2.5 L / mat	28" x 40" (71 cm x 102 cm)	20 / Box
<b>SS-20-3000</b>	Safe-Soak™ / Yellow absorbent mat	3 L / mat	30" x 52" (76 cm x 132 cm)	15 / Box
<b>SS-21-1000</b>	Safe-Soak™ SPAC	N/A	8.625" (22 cm diameter)	12 / Box

# BATRIK BLUE CLOUD-GEL™ ANTI-FATIGUE MATS


Our Blue Cloud-Gel™ Anti-Fatigue Mats have been scientifically engineered for use in the O.R. using patented gel technology and premium foam to provide the ultimate comfort and support for surgeons, nurses, and technicians particularly when movement is restricted for prolonged durations.

- Proprietary 3/4" thickness Blue Cloud-Gel™ foam technology will not "bottom-out" like conventional foam mats
- Conforms to and supports feet promoting efficient circulation, reducing physical discomfort and fatigue
- Promotes proper posture while diminishing back, knee, and foot discomfort
- Latex free, water resistant, and nonabsorbent
- Reusable, easy to clean, and Cart Washer compatible (contact us for more details)
- Excellent friction properties when damp or dry surgical booties are worn
- Completely sealed with an RF weld for excellent protection against moisture absorption
- US Standards: CALL 117 and NFPA 260 Compliant
- 1 Year Warranty

DID YOU KNOW? STUDIES HAVE DEMONSTRATED THAT ANTI-FATIGUE MATS IMPROVE PERFORMANCE, ENERGY, AND POST-OPERATIVE DISCOMFORT.

## PRODUCT SELECTION

ITEM #	DIMENSIONS	PACKAGING
<b>800-2401</b>	18" X 24" / 46 cm x 61 cm	1 / Box
<b>800-2405</b>	18" X 24" / 46 cm x 61 cm	5 / Box
<b>800-3201</b>	20" X 32" / 51 cm x 81 cm	1 / Box
<b>800-3203</b>	20" X 32" / 51 cm x 81 cm	3 / Box
<b>800-4801</b>	20" X 48" / 51 cm x 122 cm	1 / Box
<b>800-4802</b>	20" X 48" / 51 cm x 122 cm	2 / Box
<b>800-7201</b>	30" X 72" / 76 cm x 183 cm (Special Order)	1 / Box



18" X 24"

20" X 32"


20" X 48"

# BATRIK BLUE-CLOUD NIMBUS™ ANTI-FATIGUE MATS

Our Blue-Cloud Nimbus™ Anti-Fatigue Mats are designed for general hospital use. Made with proprietary high-density “energy-return” polyurethane materials, they provide unmatched comfort, relief, and ergonomic support.

- Our Energy-Return properties will help reduce leg, foot, lower back discomfort, and work-related injuries
- Longer lasting support as 3/4” thickness proprietary foam blend will not bottom out like most other foam based anti-fatigue mats
- Antimicrobial properties blended into the mat provides added protection
- Highly water resistant surface will not absorb water or dirt
- Lightweight and easy to clean
- Textured design on surface and high-traction base ensures greater safety. (US Standard NFSI 101-C)
- Beveled edging reduces tripping hazards and anti-slip backing enhances safety
- Durable design is ideal for commercial applications, yet still ergonomically and aesthetically pleasing for home or office environments
- Flammability/Burn resistant per CFR 16 Parts 1630/1631
- Made with eco-friendly materials, is non-toxic, and latex free
- 5 Year Warranty (non-industrial usage), 4 Year Warranty (industrial usage)

GOOD TO KNOW: IF WITHIN  
30 DAYS YOU ARE NOT 100%  
SATISFIED, CONTACT YOUR  
DEALER FOR A REFUND  
(EXCLUDES TRANSPORT,  
HANDLING, AND TAXES)



Available in  
Taupe or Black


## PRODUCT SELECTION

ITEM #	DIMENSIONS INCHES / CM	PACKAGING
<b>BLACK COLOR</b>		
<b>809-3201</b>	20" X 32" / 51 cm x 81 cm	1 / Box
<b>809-3204</b>	20" X 32" / 51 cm x 81 cm	4 / Box
<b>809-4801</b>	20" X 48" / 51 cm x 122 cm	1 / Box
<b>809-4804</b>	20" X 48" / 51 cm x 122 cm	4 / Box
<b>809-3601</b>	24" X 36" / 61 cm x 91 cm	1 / Box
<b>809-7201</b>	20" X 72" / 51 cm x 183 cm	1 / Box
<b>809-6001</b>	36" X 60" / 91 cm x 152 cm	1 / Box
<b>TAUPE COLOR</b>		
<b>810-3201</b>	20" X 32" / 51 cm x 81 cm	1 / Box
<b>810-3204</b>	20" X 32" / 51 cm x 81 cm	4 / Box
<b>810-4801</b>	20" X 48" / 51 cm x 122 cm	1 / Box
<b>810-4804</b>	20" X 48" / 51 cm x 122 cm	4 / Box
<b>810-3601</b>	24" X 36" / 61 cm x 91 cm	1 / Box
<b>810-7201</b>	20" X 72" / 51 cm x 183 cm	1 / Box
<b>810-6001</b>	36" X 60" / 91 cm x 152 cm	1 / Box

**BATRIK**

# STACKABLE STEP STOOLS

Made of lightweight cast aluminium and protected with an antimicrobial powder coating, our patented Stackable Step Stools have been designed with infection prevention as a prime focus.



## PRODUCT FEATURES

- Weight: 8.3 lbs (3.8 kg)
- Dimensions: 17.4" L x 13" W x 4.7" H (44.2 cm L x 33 cm W x 11.9 cm H)
- Transport Cart Available
- Weight Capacity: 300 lb (136 kg)
- Warranty: 7 years (step stool only)

## ANTIMICROBIAL COATING:

Prevents the spread of hazardous microorganisms and suppresses the growth of mold, bacteria, yeast, and mildew

## ANTI-SLIP SURFACE:

Provides enhanced stability.

## RUBBER-FREE CORNERS:

Avoids fluid collection.

## WASHABLE CART:

Model OR-36364 (without mat) can be sterilized.

## HYGIENIC DESIGN:

No welds or large hollow openings to trap debris or fluids

## PRODUCT SELECTION

ITEM #	DESCRIPTION	COLOR	CASE DIMENSIONS
<b>OR-36363</b>	Stackable Step Stool with mat	Black	19" x 14" x 6"
<b>OR-36363-02</b>	Stackable Step Stool with mat	Dark Blue	19" x 14" x 6"
<b>OR-36363-03</b>	Stackable Step Stool with mat	Green	19" x 14" x 6"
<b>OR-36363-07</b>	Stackable Step Stool with mat	Orange	19" x 14" x 6"
<b>OR-36364</b>	Stackable Step Stool without mat	N/A	19" x 14" x 6"


**PACKAGING:** Sold by Unit

## REPLACEMENT PARTS

ITEM #	DESCRIPTION	PACKAGING	CASE DIMENSIONS
<b>OR-36363-1A</b>	Replacement mat for OR-36363	1 Unit	15" x 18" x 5" (Ziplock)
<b>OR-36363-1B</b>	Replacement Foot Pegs	1 Unit	Ziplock Based on Quantity ordered

## OPTIONAL ACCESSORIES

ITEM #	DESCRIPTION	PACKAGING	CASE DIMENSIONS
<b>OR-3693</b>	Transport cart	1 Unit	34" x 17" x 7"
<b>809-1705</b>	Reusable Anti-Fatigue Mat with adhesive strips to secure it in place	5 / Box	12" x 12" x 4"
<b>899-1710</b>	Single-Use Anti-Fatigue Mat with adhesive strips to secure it in place	10 / Box	12" x 12" x 8"


**OR-3693**

Ver. 2024-11-29

BATRIK

# TOURNY™ + TOURNIQUETS

Our Tourniquets are comfortable and easy to use on all patients, including the elderly, undernourished, children, and newborns.



## AVAILABLE COLORS



- Available in different sizes including compact sizes, which are convenient for phlebotomy trays
- Innovative latex free material has a more comfortable texture for patients, puts less stress on skin, and safeguards hospital staff and patients from latex allergens

## SINGLE-USE BENEFITS

- Reduced risk of cross-contamination (patient to patient or patient to healthcare provider)
- Reduced risk of spread of nosocomial infections
- May improve your facility's compliance with OSHA's 1910-1030 (d)(1) "Blood-borne Pathogens Standard"



FOR YOUR CONVENIENCE:  
 25 UNITS PER ROLL  
 COME IN A DISPENSER BOX  
 FOR EASY ACCESS



## PRODUCT SELECTION

ITEM #	DESCRIPTION	COLOR	DIMENSIONS	PACKAGING
<b>K-002</b>	TOURNY™+ Non-Sterile	Blue	1" x 18" / 2.5 cm x 45 cm	1 / Bag
<b>K-004</b>	TOURNY™+ Non-Sterile	Blue	1" x 18" / 2.5 cm x 45 cm	100 / Box
<b>K-009</b>	TOURNY™+ Non-Sterile	Blue	0.200" x 22" / 0.5 cm x 56 cm	20 / Box
<b>K-010</b>	TOURNY™+ Non-Sterile	White	0.200" x 22" / 0.5 cm x 56 cm	20 / Box
<b>K-400</b>	TOURNY™+ Non-Sterile	Blue	1" x 18" / 2.5 cm x 45 cm	16 rolls 400 / Box
<b>K-400PB</b>	TOURNY™+ Non-Sterile (Pebbled)	Light Blue	1" x 18" / 2.5 cm x 45 cm	16 rolls 400 / Box
<b>K-100</b>	TOURNY™+ Non-Sterile	Assorted (3 rolls of each) Blue, Orange, Green, and Pink	1" x 18" / 2.5 cm x 45 cm	12 rolls 300 / Box
<b>K-254</b>	TOURNY™+ Non-Sterile (Pebbled)	Light Blue	1" x 18" / 2.5 cm x 45 cm	250 / Box
<b>K-2-18RL</b>	TOURNY™+ Non-Sterile	Blue	1" x 18" / 2.5 cm x 45 cm	1 Roll
<b>K-2-18PBRL</b>	TOURNY™+ Non-Sterile (Pebbled)	Light Blue	1" x 18" / 2.5 cm x 45 cm	1 Roll
<b>K-3-18RL</b>	TOURNY™+ Non-Sterile	Green	1" x 18" / 2.5 cm x 45 cm	1 Roll
<b>K-6-18RL</b>	TOURNY™+ Non-Sterile	Pink	1" x 18" / 2.5 cm x 45 cm	1 Roll
<b>K-7-18RL</b>	TOURNY™+ Non-Sterile	Orange	1" x 18" / 2.5 cm x 45 cm	1 Roll
<b>K-2-29RL</b>	Long TOURNY™+ Non-Sterile	Blue	1" x 29.5" / 2.5 cm x 75 cm	1 Roll

# SHOOZ<sup>COVERS</sup>™ | **SOFT** LINE™


Batrik ShoozCovers™ Soft Line™ is manufactured in Canada from high quality 35 g/m<sup>2</sup> polypropylene fabric and 100% latex free elastics for performance and durability you can rely on. Our anti-skid tread reduces the risk of slipping and our smart design delivers easy on-and-off convenience and great fit.

## SHOOZ-REPELLENT™ TECHNOLOGY:

The **Soft Line™** integrates Shooz-Repellent™ technology which creates a protective barrier that repels fluids, delivering basic incidental fluid contact repellent. Shoe covers are not waterproof or water-resistant.

## ANTI-SKID TREAD:

Our anti-skid tread reduces the risk of slipping with **triple ultrasonic welds** that increases tear strength and delivers easy on-and-off convenience.



MADE IN CANADA 



## OUR SOFT LINE™ IS PERFECT FOR:



MEDICAL  
ESTABLISHMENTS



PHARMACEUTICAL  
LABS



GMP  
SITES



ELECTRONIC  
MANUFACTURING

LOLLY BLUE | **SOFT LINE™**


## WHY YOU'LL LOVE THE SOFT LINE™:



- Both shoe covers are manufactured in Canada with top-quality 35 g/m<sup>2</sup> polypropylene fabric and 100% latex-free USA-made elastics!
- This line is particle count tested and complies with IEST-RP-CC003.4.
- Other colors available upon request.

UN XL

AVAILABLE IN TWO VERSATILE SIZES

BUBBLE GUM PINK | **SOFT LINE™**

## SUPPORT A CAUSE, WEAR PINK:



- A portion of our Bubble Gum Pink shoe cover sales will be donated to charity to help further breast cancer research. Join us in making a meaningful contribution by choosing Pink ShoozCovers™!

UN

AVAILABLE IN UNIVERSAL



## PRODUCT SELECTION:

ITEM #	DESCRIPTION	COLOR	SIZE	PACKAGING	CASE DIMENSIONS
<b>SC-1210-PB</b>	ShoozCovers™ - Soft Line™ <i>polypropylene non-woven fabric</i>	Lolly Blue	Universal Up to shoe size 12 ~ 42 cm x 16.5 cm / 16.5" x 6.5"	Case of 1,000 (10 bags of 100 units)	20"x 16"x 16"
<b>SC-1216-PDB</b>	ShoozCovers™ - Soft Line™ <i>polypropylene non-woven fabric</i>	Lolly Blue	Universal Shoe size 12 ~ 42 cm x 16.5 cm / 16.5" x 6.5"	Case of 600 (6 Dispenser Boxes of 100 units)	16.25"x 8.875"x 18.75"
<b>SC-1216-PB</b>	ShoozCovers™ - Soft Line™ <i>polypropylene non-woven fabric</i>	Lolly Blue	Universal Shoe size 12 ~ 42 cm x 16.5 cm / 16.5" x 6.5"	Case of 600 (6 bags of 100 units)	19.5"x 19.5"x 14.25"
<b>SC-1226-PDB</b>	ShoozCovers™ - Soft Line™ <i>polypropylene non-woven fabric</i>	Lolly Blue	XL Shoe size 12-15 ~ 47 cm x 19 cm / 18.5" x 7.5"	Case of 600 (6 Dispenser Boxes of 100 units)	17.25"x 9.375"x 21.75"
<b>SC-1220-PB</b>	ShoozCovers™ - Soft Line™ <i>polypropylene non-woven fabric</i>	Lolly Blue	XL Shoe size 12-15 ~ 47 cm x 19 cm / 18.5" x 7.5"	Case of 1,000 (10 bags of 100 units)	22"x 16"x 16"
<b>SC-1616-PDB</b>	ShoozCovers™ - Soft Line™ <i>polypropylene non-woven fabric</i>	Bubble Gum Pink	Universal Shoe size 12-15 ~ 42 cm x 16.5 cm / 16.5" x 6.5"	Case of 600 (6 Dispenser Boxes of 100 units)	16.25"x 8.875"x 18.75"
<b>SC-1616-PB</b>	ShoozCovers™ - Soft Line™ <i>polypropylene non-woven fabric</i>	Bubble Gum Pink	Universal Shoe size 12-15 ~ 42 cm x 16.6 cm / 16.5" x 6.5"	Case of 600 (6 bags of 100 units)	19.5"x 19.5"x 14.25"

ALTERNATIVE PACKAGING: Add "-PDB" to end of code for a dispenser box instead of bags.

<sup>1</sup>Many suppliers claim their products are latex free, however this often only relates to the fabric and not the elastics used. We source our ShoozCovers™ elastics from the USA to ensure that our components conform fully as a 100% latex free product.


# SHOOZCOVERS™ | COATED LINE™

## STAY DRY:

Our Batrik ShoozCovers™ Coated Line™ is made from a ~55 g/m<sup>2</sup> blend of **polypropylene** and **polyethylene**, providing **high-quality water resistance**.<sup>1</sup> Ideal for environments where **maximum protection against fluids** is required.

## TRIPLE WELDED TECHNOLOGY:

These shoe covers are reinforced with **triple ultrasonic welds**, a technique that joins materials together without the need for adhesives or stitching. This technique **increase tear strength and reduce linting**, for long lasting, resilient covers.



MADE IN CANADA 



## OUR COATED LINE™ IS PERFECT FOR:



MEDICAL  
ESTABLISHMENTS



PHARMACEUTICAL  
LABS



GMP SITES



ELECTRONIC  
MANUFACTURING

<sup>1</sup>OUR SHOOZCOVERS™ COATED LINE™ IS RESISTANT AGAINST ALL NON PETROLEUM BASED LIQUIDS.

## MARSHMALLOW WHITE | COATED LINE™

### WHY YOU'LL LOVE THE COATED LINE™:

- **Superior protection using our unique blend of fabrics** which protects against not only water but all non-petroleum-based liquids, making them the ideal choice for environments where spills and leaks are a concern
- **Elastics around the ankle create a strong and secure fit** which is especially important when working in hazardous environments where a misplaced shoe cover could lead to contamination or injury
- **Smart design** with a large opening, providing a better fit and easy on-and-off convenience
- **Other colors available upon request**



UN XL

AVAILABLE IN TWO VERSATILE SIZES

### PRODUCT SELECTION:

ITEM #	DESCRIPTION	COLOR	SIZE	PACKAGING	CASE DIMENSIONS
<b>SC-2112-NPB</b>	ShoozCovers™ - Coated Line™ <i>Polypropylene + polyethylene non-woven fabric, waterproof</i>	Marshmallow White	Universal Up to shoe size 12 ~ 42 cm x 17.5 cm / 16.5" x 6.88"	Case of 300 (2 bags of 150 units)	16" x 13.375" x 10"
<b>SC-2112-NPDB</b>	ShoozCovers™ - Coated Line™ <i>Polypropylene + polyethylene non-woven fabric, waterproof</i>	Marshmallow White	Universal Up to shoe size 12 ~ 42 cm x 17.5 cm / 16.5" x 6.88"	Case of 300 (2 boxes of 150 units)	16.62" x 14.875" x 8.75"
<b>SC-2122-NPB</b>	ShoozCovers™ - Coated Line™ <i>Polypropylene + polyethylene non-woven fabric, waterproof</i>	Marshmallow White	XL Shoe size 12-15 ~ 47 cm x 19 cm / 18.5" x 6.88"	Case of 240 (2 bags of 120 units)	16" x 13.375" x 10"
<b>SC-2122-NPDB</b>	ShoozCovers™ - Coated Line™ <i>Polypropylene + polyethylene non-woven fabric, waterproof</i>	Marshmallow White	XL Shoe size 12-15 ~ 47 cm x 19 cm / 18.5" x 6.88"	Case of 240 (2 boxes of 120 units)	16.62" x 14.875" x 8.75"

ALTERNATIVE PACKAGING: Add "-PDB" to end of code for a dispenser box instead of bags

# SHOOZCOVERS™ | TOUGH LINE™

## EXTRA-LARGE FOR EXTRA-TOUGH WORK:


Our Batrik ShoozCovers™ Tough Line™ comes in an **extra-large** size that fits over a **men's size 15 steel-toed boot**, delivering optimal full coverage and protection while you work.<sup>1</sup>

## ENHANCED DURABILITY:

These shoe covers are formulated with our **thickest fabric yet!** Made of **~55 g/m² polypropylene fabric**, these highly durable covers provide added protection, increase tear strength and reduce linting.



SHOE SIZES 12-15 (45-49.5 EU)  
~ 47 CM X 19 CM / 18.5" X 7.5"



MADE IN CANADA 🇨🇦



## OUR TOUGH LINE™ IS PERFECT FOR:



HOMES OR  
COMMERCIAL  
BUILDINGS



SHIELDING  
AGAINST MUD  
& GREASE



MEDICAL  
ESTABLISHMENTS



PHARMACEUTICAL  
LABS

<sup>1</sup> SHOOZCOVERS™ ARE MADE LONGER (42 & 47CM) FOR ACCURATE SHOE SIZE SELECTION.

## BLACK LIQUORICE | TOUGH LINE™

### WHY YOU'LL LOVE THE TOUGH LINE™:

- Triple ultrasonic welds that are long lasting to reduce waste
- High quality made in the USA latex-free elastics are triple reinforced to contract the shoe cover around your foot and deliver a strong and secure fit<sup>2</sup>
- Shooz-Repellent™ technology delivers basic incidental fluid contact repellency<sup>3</sup>
- Other colors available upon request



XL

AVAILABLE IN EXTRA LARGE

### ANY TIME, ANYWHERE:



RESIDENTIAL



COMMERCIAL



INDUSTRIAL

### PRODUCT SELECTION:

ITEM #	DESCRIPTION	COLOR	SIZE	PACKAGING	CASE DIMENSIONS
<b>SC-1922-PB</b>	ShoozCovers™ · Tough Line™ Polypropylene + polyethylene non-woven fabric	Black Liquorice	XL Shoe size 12-15 ~ 47 cm x 19 cm / 18.5" x 7.5"	Case of 200 (2 bags of 100 units)	16" x 13.375" x 10"

ALTERNATIVE PACKAGING: Add "-PDB" to end of code for a dispenser box instead of bags

<sup>2</sup> Many suppliers claim their products are latex free, however this often only relates to the fabric and not the elastics used. We source our ShoozCovers™ elastics from the USA to ensure that our components conform fully as a 100% latex free product. <sup>3</sup> Shooz-Repellent™ technology provides basic incidental fluid contact repellent. It is not waterproof nor water resistant.

# GRIPPEEZ SHOOZ™



## RELIABLE, DURABLE, AND TRULY ANTI-SLIP SHOE COVERS

ShoozCovers™ Grippeez provide exceptional protection and **enhanced traction** for professionals working in the most demanding, slippery environments. Featuring waterproof protection, our shoe covers ensure safety and hygiene without compromising on quality. They are the ideal solution for tough jobs.

MADE IN CANADA 



## PERFECT FOR:



AEROSPACE



FOOD  
PROCESSING  
INDUSTRY



OIL, GAS AND  
MINING SITES



CONSTRUCTION &  
RENOVATION SITES

<sup>1</sup>Industrial strength elastics

GRIPPEEZSHOOZ™ ARE A PERFECT FIT FOR EVERYONE:

- Fast, Reliable Fit:** Stretch material ensures a secure, easy-on fit. Coating provides protection against slippery surfaces.
- Enhanced Safety:** Anti-slip rubberized coating provides reliable grip on all surfaces, helping to prevent sliding in critical areas.
- Comes in an **extra-large size** that fits over a **men's size 15 steel-toed boot**.
- Available in two colors: **Gecko Navy** and **Rhino Grey**.

RHINO GREY



GECKO NAVY



PRODUCT SELECTION:

ITEM #	DESCRIPTION	COLOR	SIZE	PACKAGING	CASE DIMENSIONS
SC-29312-NPDB	Grippez ShoozCovers™ Polypropylene + polyethylene non-woven fabric, waterproof	Rhino Grey	Universal Up to shoe size 12 ~ 42 cm x 17.5 cm / 16.5" x 6.88"	Case of 300 (2 dispenser boxes of 150 units)	16.62" x 14.875" x 8.75"
SC-29312-NPB	Grippez ShoozCovers™ Polypropylene + polyethylene non-woven fabric, waterproof	Rhino Grey	Universal Up to shoe size 12 ~ 42 cm x 17.5 cm / 16.5" x 6.88"	Case of 300 (2 bags of 150 units)	16" x 13.375" x 10"
SC-29322-NPDB	Grippez ShoozCovers™ Polypropylene + polyethylene non-woven fabric, waterproof	Rhino Grey	XL Shoe size 12-15 ~ 47 cm x 19 cm / 18.5" x 6.88"	Case of 240 (2 dispenser boxes of 120 units)	16.62" x 14.875" x 8.75"
SC-29322-NPB	Grippez ShoozCovers™ Polypropylene + polyethylene non-woven fabric, waterproof	Rhino Grey	XL Shoe size 12-15 ~ 47 cm x 19 cm / 18.5" x 6.88"	Case of 240 (2 bags of 120 units)	16" x 13.375" x 10"
SC-22212-NPDB	Grippez ShoozCovers™ Polypropylene + polyethylene non-woven fabric, waterproof	Gecko Navy	Universal Up to shoe size 12 ~ 42 cm x 17.5 cm / 16.5" x 6.88"	Case of 300 (2 dispenser boxes of 150 units)	16.62" x 14.875" x 8.75"
SC-22212-NP	Grippez ShoozCovers™ Polypropylene + polyethylene non-woven fabric, waterproof	Gecko Navy	Universal Up to shoe size 12 ~ 42 cm x 17.5 cm / 16.5" x 6.88"	Case of 300 (2 bags of 150 units)	16" x 13.375" x 10"
SC-22222-NPDB	Grippez ShoozCovers™ Polypropylene + polyethylene non-woven fabric, waterproof	Gecko Navy	XL Shoe size 12-15 ~ 47 cm x 19 cm / 18.5" x 6.88"	Case of 240 (2 dispenser boxes of 120 units)	16.62" x 14.875" x 8.75"
SC-22222-NPB	Grippez ShoozCovers™ Polypropylene + polyethylene non-woven fabric, waterproof	Gecko Navy	XL Shoe size 12-15 ~ 47 cm x 19 cm / 18.5" x 6.88"	Case of 240 (2 bags of 120 units)	16" x 13.375" x 10"

ALTERNATIVE PACKAGING: Add "-PDB" to end of code for a dispenser box instead of bags


# BATRIK ENDOSCOPY BO2™ BITE BLOCKS WITH NASAL OXYGEN DELIVERY

Our BO2™ Bite Blocks significantly enhance patient safety by providing supplemental oxygenation during gastrointestinal endoscopy, preventing hypoxia in sedated patients, and ensuring superior patient comfort.

- Saves time and money by eliminating the need for separate nasal prongs and oxygen masks
- Designed to be free of sharp edges to ensure patient comfort and safety
- Available exclusively in Single-Use
- Optional retention strap is available to ensure a secure fit
- The Bite Blocks and optional retention straps are all latex free, lead free and DEHP free

**BO2**  
(FRONT SIDE)

**BO2**  
(BACK SIDE)



## PRODUCT SELECTION

ITEM #	TYPE	DESCRIPTION	PACKAGING	CASE DIMENSION
<b>BO2</b>	Single-Use	18mm / 54FR, includes strap	100 / Box	10" x 10" x 10"
<b>BO2-N</b>	Single-Use	18mm / 54FR, without strap	100 / Box	10" x 10" x 10"

# BATRIK ENDOSCOPY BITESTRONG™ BITE BLOCKS

## ADULT & PEDIATRIC BITE BLOCKS

Our Bitestrong™ Bite Blocks provide a safe passage through the patient's mouth during TEE, protecting the patient's teeth and endoscope.

Our molding process ensures a smooth finish with no sharp edges (flashing) that can often irritate a patient's mouth. All Bite Blocks are individually wrapped in sealed, easy to open PE bags and are latex free, lead free, and DEHP free. Additional straps are available, if required

### SINGLE USE

- Available without dental retention

### REUSABLE

- Designed for repeated usage and will withstand commonly practiced chemical and steam sterilization cycles (134°C / 275°F).
- Available with or without an adjustable retention strap

## PRODUCT SELECTION

ITEM #	TYPE	DESCRIPTION	PACKAGING	CASE DIMENSIONS
<b>B1</b>	Single-Use <a href="#">Pediatric</a>	16 mm / 48FR, includes strap	100 / Box	10" x 10" x 10"
<b>B1-25</b>	Single-Use <a href="#">Pediatric</a>	16 mm / 48FR, includes strap	25 / Box	10" x 7" x 7"
<b>B1-N</b>	Single-Use <a href="#">Pediatric</a>	16 mm / 48FR, without strap	100 / Box	10" x 10" x 10"
<b>B3</b>	Single-Use	20 mm / 60FR, includes strap	100 / Box	10" x 10" x 10"
<b>B3-N</b>	Single-Use	20 mm / 60FR, without strap	100 / Box	10" x 10" x 10"
<b>B3-SA</b>	Single-Use	20 mm / 60FR, with strap attached to one side	100 / Box	10" x 10" x 10"
<b>B3-NDR</b>	Single-Use	20 mm / 60FR, no dental retention rim	100 / Box	10" x 10" x 10"
<b>B3-NDR-N</b>	Single-Use	Same as B3-N, no dental retention rim, without strap.	100 / Box	10" x 10" x 10"
<b>B4</b>	Single-Use	19 mm / 60FR	100 / Box	10" x 10" x 10"
<b>B4T</b>	Single-Use	19 mm / 60FR, includes strap	100 / Box	10" x 10" x 10"
<b>B4T-N</b>	Single-Use	19 mm / 60FR, without strap	100 / Box	10" x 10" x 10"
<b>B6</b>	Reusable	19 mm / 60FR	100 / Box	10" x 10" x 10"
<b>B6T</b>	Reusable	19 mm / 60FR, includes strap	100 / Box	10" x 10" x 10"
<b>B6T-N</b>	Reusable	19 mm / 60FR, without strap	100 / Box	10" x 10" x 10"
<b>B6-25</b>	Reusable	19 mm / 60FR	25 / Box	10" x 7" x 7"
<b>B6T-25</b>	Reusable	19 mm / 60FR, includes strap	25 / Box	10" x 7" x 7"
<b>B6T-N-25</b>	Reusable	19 mm / 60FR, without strap	25 / Box	10" x 7" x 7"
<b>BS</b>	Single-Use	Bite-Block straps 21.7" x 0.79" / 55 cm x 2 cm	100 / Box	10" x 13" x 1.5"

**BS****B1****B3****B4****B6**

BATRIK

# COUVRETTE™ VALVE COVERS

Our Couvrette™ Valve Covers eliminate the need for manual cleaning and reprocessing of reusable valves. Our single use Valve Covers effectively reduce cross contamination between patients and also reduce fluid leakage from the parts during procedures.

- Easily and securely fit to all valves
- Reduce the potential for errors
- Single-Use
- Non-Sterile
- Latex free



16-0002

16-0004

TIP: REUSABLE VALVE COVERS ARE NOTORIOUSLY DIFFICULT TO CLEAN.

SAVE TIME AND REDUCE RISK WITH COUVRETTE™.

COMING SOON!



16-0022

16-0044

## PRODUCT SELECTION

ITEM #	DESCRIPTION	COLOR	PACKAGING
16-0002	Biopsy Valve Cover - Olympus®	Blue	100 / Box
16-0004	Biopsy Valve Cover - Pentax®	Red	100 / Box
16-0022	Biopsy Valve Cover - Olympus® & Fujinon®	Light Blue	100 / Box <a href="#">COMING SOON!</a>
16-0044	Biopsy Valve Cover - Pentax®	Pink	100 / Box <a href="#">COMING SOON!</a>


NOTE: Contact us to verify compatibility with your scope

# BATRIK DAYBRUSH™ CHANNEL BRUSHES

## REUSABLE

Our Reusable DayBrush™ Channel Brushes can be used for cleaning endoscopes for an entire day. In fact, it's been proven to effectively clean up to 10 scopes before the need to discard.\*

- Soft tip ensures easy manipulation within the channel with no risk of damage
- The back-end is sealed and softened to prevent fluid from getting trapped inside the tube. The softening process also prevents damage if the cleaning brush is used in reverse.
- Can be steam sterilized.
- Bristle volume is 20mm in length – nearly double most brushes on the market. A large brush provides more surface contact during the cleaning process.
- Significant per use savings when compared to other reusable or disposable channel cleaning brushes.



20 mm / 0.79"

*(\*) Discard after cleaning 10 scopes or sooner if any brush degradation is noticed. After ten applications, the brush may have the appearance that it can continue to clean effectively, but we discourage prolonged use due to risk of cross contamination. Actual number of uses will vary based on handling and processing.*

## PRODUCT SELECTION

### 5 MM CYLINDRICAL BRUSHES – SHEATH Ø 1.75 MM

ITEM #	COLOUR	LENGTH	BRUSH LENGTH	BRUSH (Ø)	SHEATH (Ø)	BRUSH ENDS
<b>30-18050</b>	White	70" / 180 cm	0.79" / 20 mm	0.20" / 5 mm	0.07" / 1.75 mm	Single
<b>30-18050-T</b>	Yellow	70" / 180 cm	0.79" / 20 mm	0.20" / 5 mm	0.07" / 1.75 mm	Single
<b>30-18050-2T</b>	Yellow	70" / 180 cm	0.79" / 20 mm	0.20" / 5 mm	0.07" / 1.75 mm	Double
<b>30-23050</b>	White	90" / 230 cm	0.79" / 20 mm	0.20" / 5 mm	0.07" / 1.75 mm	Single
<b>30-23050-T</b>	Yellow	90" / 230 cm	0.79" / 20 mm	0.20" / 5 mm	0.07" / 1.75 mm	Single
<b>30-23050-2T</b>	Yellow	90" / 230 cm	0.79" / 20 mm	0.20" / 5 mm	0.07" / 1.75 mm	Double

PACKAGING : 50 / BOX

### OPTIONAL KIT ADD-ONS

ADD THE SUFFIX:

"F" to receive a kit with our Flat Sponge (Packaging: Box of 40). Applicable on certain models only, call for details.


"K" to receive a kit with the single use Valve / Control Cleaning Brush 5/11 Double Brushes (item # 27-57787). Applicable on all models.

"FK" to receive both the sponge and valve cleaning brushes.

# BATRIK ARC™ CHANNEL BRUSHES SINGLE-USE

BATRIK offers one of the largest selection of Endoscopy Channel Cleaning Brushes in the industry. Our brushes are available in a variety of lengths, materials, and styles to accommodate the widest range of endoscopic devices used in the healthcare industry today.

No matter what your requirements, you'll be able to find the brush that suits your specific needs. If you don't see the brush you require, simply call us. We have over 50 different brushes available and can custom fabricate to your precise specifications.



## OPTIONAL KIT ADD-ONS


ADD THE SUFFIX:

"F" to receive a kit with our Flat Sponge (Packaging: Box of 40). Applicable on certain models only, call for details.

"K" to receive a kit with the single use Valve / Control Cleaning Brush 5/11 Double Brushes (item # 27-57787). Applicable on all models.

"FK" to receive both the sponge and valve cleaning brushes.

## SERIES 311 & 333



### PRODUCT SELECTION

#### 5 MM CYLINDRICAL BRUSHES – SHEATH Ø 1.85 MM - PUNCHY PINK POLYPROPYLENE

ITEM #	LENGTH	BRUSH LENGTH	BRUSH (Ø)	SHEATH (Ø)	BRUSH ENDS
<b>311-10050</b>	39" / 100 cm	0.94" / 24 mm	0.20" / 5 mm	0.073" / 1.85 mm	Single
<b>311-18050</b>	70" / 180 cm	0.94" / 24 mm	0.20" / 5 mm	0.073" / 1.85 mm	Single
<b>311-23050</b>	90" / 230 cm	0.94" / 24 mm	0.20" / 5 mm	0.073" / 1.85 mm	Single
<b>311-23050-2</b>	90" / 230 cm	0.94" / 24 mm	0.20" / 5 mm	0.073" / 1.85 mm	Double


PACKAGING: 50 / BOX

#### 5/6 MM TAPERED BRUSHES – SHEATH Ø 1.85 MM - PUNCHY PINK POLYPROPYLENE

ITEM #	LENGTH	BRUSH LENGTH	BRUSH (Ø)	SHEATH (Ø)	BRUSH ENDS
<b>333-23056</b>	90" / 230 cm	0.94" / 24 mm	0.20", 0.24" / 5, 6 mm	0.073" / 1.85 mm	Single
<b>333-23056-2</b>	90" / 230 cm	0.94" / 24 mm	0.20", 0.24" / 5, 6 mm	0.073" / 1.85 mm	Double

PACKAGING: 50 / Box

## SERIES 300 & 322



### PRODUCT SELECTION

#### 5 MM CYLINDRICAL BRUSHES – SHEATH Ø 1.75 MM - METALLIC GREY DELRIN

ITEM #	LENGTH	BRUSH LENGTH	BRUSH (Ø)	SHEATH (Ø)	BRUSH ENDS
<b>300-18050</b>	70" / 180 cm	0.79" / 20 mm	0.20" / 5 mm	0.07" / 1.75 mm	Single
<b>300-23050</b>	90" / 230 cm	0.79" / 20 mm	0.20" / 5 mm	0.07" / 1.75 mm	Single
<b>300-23050-2</b>	90" / 230 cm	0.79" / 20 mm	0.20" / 5 mm	0.07" / 1.75 mm	Double


PACKAGING: 50 / Box

#### 5/6 MM TAPERED BRUSHES – SHEATH Ø 1.75 MM - METALLIC GREY DELRIN

ITEM #	LENGTH	BRUSH LENGTH	BRUSH (Ø)	SHEATH (Ø)	BRUSH ENDS
<b>322-23056</b>	90" / 230 cm	0.79" / 20 mm	0.20", 0.24" / 5, 6 mm	0.07" / 1.75 mm	Single
<b>322-23056-2</b>	90" / 230 cm	0.79" / 20 mm	0.20", 0.24" / 5, 6 mm	0.07" / 1.75 mm	Double

PACKAGING: 50 / Box

## SERIES 344




### PRODUCT SELECTION

#### 3 MM CYLINDRICAL BRUSHES – SHEATH Ø 1.55 MM - METALLIC PURPLE DELRIN

ITEM #	LENGTH	BRUSH LENGTH	BRUSH (Ø)	SHEATH (Ø)	BRUSH ENDS
<b>344-10030</b>	39" / 100 cm	0.79" / 20 mm	0.12" / 3 mm	0.06" / 1.55 mm	Single
<b>344-18030</b>	70" / 180 cm	0.79" / 20 mm	0.12" / 3 mm	0.06" / 1.55 mm	Single
<b>344-23030</b>	90" / 230 cm	0.79" / 20 mm	0.12" / 3 mm	0.06" / 1.55 mm	Single

PACKAGING: 50 / Box

## SERIES 35



### PRODUCT SELECTION

#### 2.5 MM CYLINDRICAL BRUSHES – SHEATH Ø 1.00 MM - YELLOW DELRIN

ITEM #	LENGTH	BRUSH LENGTH	BRUSH (Ø)	SHEATH (Ø)	BRUSH ENDS
<b>35-12025</b>	47.2" / 120 cm	0.59" / 15 mm	0.098" / 2.5 mm	0.039" / 1.00 mm	Single
<b>35-23025</b>	90" / 230 cm	0.59" / 15 mm	0.098" / 2.5 mm	0.039" / 1.00 mm	Single
<b>35-23025-2</b>	90" / 230 cm	0.59" / 15 mm	0.098" / 2.5 mm	0.039" / 1.00 mm	Double

PACKAGING: 50 / Box



#### 2 MM CYLINDRICAL BRUSHES – SHEATH Ø 1.00 MM - YELLOW DELRIN

ITEM #	LENGTH	BRUSH LENGTH	BRUSH (Ø)	SHEATH (Ø)	BRUSH ENDS
<b>35-12020</b>	47.2" / 120 cm	0.59" / 15 mm	0.079" / 2.0 mm	0.039" / 1.00 mm	Single
<b>35-18020</b>	70.8" / 180 cm	0.59" / 15 mm	0.079" / 2.0 mm	0.039" / 1.00 mm	Single
<b>35-23020</b>	90" / 230 cm	0.59" / 15 mm	0.079" / 2.0 mm	0.039" / 1.00 mm	Single
<b>35-23020-2</b>	90" / 230 cm	0.59" / 15 mm	0.079" / 2.0 mm	0.039" / 1.00 mm	Double

PACKAGING: 50 / Box

BATRIK

# VALVE / CONTROL BRUSHES


**27-05197****27-57787****40-18511****36-10010****27-15311**

## PRODUCT SELECTION

ITEM #	DESCRIPTION	LENGTH	BRUSH 1 LENGTH	BRUSH 2 LENGTH	BRUSH 1 (Ø)	BRUSH 2 (Ø)
<b>27-05197</b>	Double Trumpet 5/7 White Color	5.00" / 13 cm	0.20" / 5 mm	1.00" / 25 mm	0.20" / 5 mm	0.28" / 7 mm
PACKAGING: Bag of 3 (single-use recommended). Volume Packaging: Bag of 50 use ITEM # 27-05197-50						
<b>27-57787</b>	Double. Universal 5/11 Yellow Color	6.00" / 15 cm	1.18" / 30 mm	0.79" / 20 mm	0.43" / 11 mm	0.20" / 5 mm
PACKAGING: Bag of 3 (single-use recommended). Volume Packaging: Bag of 50 use ITEM # 27-57787-50						
<b>40-18511</b>	Double Universal 5/11 Extra Long Blue Color	8.07" / 20.5 cm	1.57" / 40 mm	0.79" / 20 mm	0.43" / 11 mm	0.20" / 5 mm
PACKAGING: Bag of 3 (reusable and can be steam sterilized). Volume Packaging: Bag of 50 use ITEM # 40-18511-50						
<b>36-10010</b>	Single end large Red Color	4.3" / 11 cm	1.57" / 40 mm	N/A	0.39" / 10 mm	N/A
PACKAGING: Bag of 3 (single-use recommended). Volume Packaging: Bag of 50 use ITEM # 36-10010-50						
<b>27-15311</b>	Single end large Blue Color	4" / 10.3 cm	1.57" / 40 mm	N/A	0.39" / 10 mm	N/A
PACKAGING: Bag of 3 (single-use recommended). Volume Packaging: Bag of 50 use ITEM # 27-15311-50						

## BATRIK V-BRUSHES EXTERNAL CLEANING BRUSHES

Our External Cleaning Brushes for endoscopes / optical instruments provide excellent cleaning capabilities. With their unique U-shape bristles and innovative design, our brushes are the ultimate tool for cleaning the exterior of your endoscopic instruments. Free from sharp or metallic ends that can scratch or cause damage to valuable equipment, both single-use and reusable models are latex free and feature ergonomic handles.



### PRODUCT SELECTION

ITEM #	DESCRIPTION	BRUSH LENGTH / (Ø)	HANDLE LENGTH / (Ø)	TOTAL LENGTH	PACKAGING
<b>36-20005</b>	Single-Use External Cleaning Brush	70 mm / 25 mm	80 mm / 3 mm	150 mm	50 / Box
<b>36-20010</b>	Reusable External Cleaning Brush	90 mm / 25 mm	120 mm / 12 mm	210 mm	3 / Bag

NEW!

BATRIK


# TAG-IT™ | NOTES

## WHEN SIMPLY BROKEN DOESN'T CUT IT—TAG IT, TELL IT, OWN IT!

Tag-It™ Notes offer the freedom to **identify and document instructions for instruments, trays, and carts** before and after **sterilization**. Tag-It™ Notes are **tear-resistant, withstand extreme chemicals, high temperature** and ensures **reliable communication**.

### WHY YOU'LL LOVE THEM:

- Tag-It™ Notes gives you the **freedom** to clearly document important information **without the clutter**.
- Robust and easy-to-use **tie wraps** designed to **securely attach** to instruments (134°C & chemical tested).
- Includes a **fast-drying, heat-resistant permanent marker** for writing your important notes.
- Choose from other **customizable cards** for personalized needs or **pre-made cards** for fast, convenient use.
- When you need more than a "repair" or "sharpen" seal to identify an issue. Instead, **Tag it, Tell it, Own it.**



**WATERPROOF**  
The synthetic paper withstands spills and liquid exposure without any damage.

**ECO-FRIENDLY**  
Made from sustainable materials.

**HEAT RESISTANT**  
Can resist high temperatures of up to 273.2°F (134°C). Compatible with Gas Plasma Sterilization (STERAD®100 NX Sterilization).

**DURABLE & TEAR-PROOF**  
Handles rigorous cleaning, tough conditions, and resists tearing.

### PRODUCT SELECTION

ITEM #	COLOR	SIZE	PACKAGING
TI-NOTE-R0	Yellow	NOTES	3" x 5" / 8 cm x 13cm
			100 / PKG


All Tag-It™ are single-use. PACKAGING: 1 pack of 100 tags, 1 permanent marker, and 100 ties.

# TAG-IT™

Tag-It™ Cards are an effortless way to add information with either pre-printed or blank formats. These **durable** tags streamline your workflow, **withstand tough conditions**, and feature **(optional) indicator dots** for quick verification. They help keep **instruments, trays, and carts** organized, ensuring efficiency before and after sterilization/decontamination.

## KEY PRODUCT FEATURES:

- **Special Permanent Marker:** The ink will not run or appear blurry after exposure to chemicals or high temperature (max. 260°C/500°F) sterilization. *(Regular pens or markers can be used only if the tags are not going to be exposed to any sterilization/decontamination).*
- **Tie Wraps:** Assist in securing the tags to almost any instrument or container. *(Tie wraps have been tested for chemical and steam sterilization).*
- **Our Tag-It™:** Will not shrink, absorb water/chemicals, or lose its intended integrity after chemical or steam sterilization.



**(Optional)**  
**All Tags Available with Indicator Dots**  
Provides visual confirmation that the tag was processed.



CLICK BELOW TO  
WATCH OUR VIDEO



## PRODUCT SELECTION:

ITEM #	COLOR	DESCRIPTION	SIZE	PACKAGING
<b>TI-3502</b>	Blue	VENDOR	3" x 5" / 8 cm x 13cm	100 / PKG
<b>TI-3503</b>	Green	LOANER	3" x 5" / 8 cm x 13cm	100 / PKG
<b>TI-3504</b>	Red	PRIORITY	3" x 5" / 8 cm x 13cm	100 / PKG
<b>TI-3505</b>	Yellow	RE-RUN	3" x 5" / 8 cm x 13cm	100 / PKG
<b>TI-3506</b>	Pink	MISSING	3" x 5" / 8 cm x 13cm	100 / PKG
<b>TI-3507</b>	Orange	BROKEN	3" x 5" / 8 cm x 13cm	100 / PKG
<b>TI-3509</b>	Black	BLANK	3" x 5" / 8 cm x 13cm	100 / PKG
<b>TI-35B1</b>	White	BLANK	3" x 5" / 8 cm x 13cm	100 / PKG
<b>TI-35B4</b>	Red	BROKEN	3" x 5" / 8 cm x 13cm	100 / PKG
<b>TI-35C22</b>	Light Blue	VENDOR CHECK IN + OUT	3" x 5" / 8 cm x 13cm	100 / PKG
<b>TI-35L3</b>	Green	LOANER CHECK IN + OUT	3" x 5" / 8 cm x 13cm	100 / PKG
<b>TI-35T2</b>	Blue	QUICK TURNOVER	3" x 5" / 8 cm x 13cm	100 / PKG
<b>TI-PM09</b>	Black	FINE PERMANENT MARKERS	1.5" Ø x 5.5" / 3.81cm Ø x 12.7 cm	36 / Bag

PACKAGING: 1 pack of 100 tags, 100 tie wraps, and 1 permanent marker.

OPTIONAL | PROCESSING INDICATOR DOTS: Add suffix "ST" to end of code for steam dots, "PS" for plasma dots.

Available in other languages or customize your own.

## PRODUCT SELECTION: (customization also available)

### VENDOR / BLUE

ITEM #	SIZE
TI-3502	3" x 5" (8 cm x 13cm)



**VENDOR**

Company Name \_\_\_\_\_

Tray Name \_\_\_\_\_

Date of Surgery \_\_\_\_\_

WWW.BATRIK.COM TAG-IT™ TI-3502

### NOTES / YELLOW

ITEM #	SIZE
TI-NOTE-R0	3" x 5" (8 cm x 13cm)



**Notes**

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_


\_\_\_\_\_

\_\_\_\_\_

WWW.BATRIK.COM TAG-IT™ TI-NOTE-R0

### LOANER / GREEN

ITEM #	SIZE
TI-3503	3" x 5" (8 cm x 13cm)



**LOANER**

Loaned From \_\_\_\_\_


Loaned To \_\_\_\_\_

By \_\_\_\_\_ Date \_\_\_\_\_

WWW.BATRIK.COM TAG-IT™ TI-3503

### BLANK / WHITE

ITEM #	SIZE
TI-35B1	3" x 5" (8 cm x 13cm)



\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

WWW.BATRIK.COM TAG-IT™ TI-35B1

### PRIORITY / RED

ITEM #	SIZE
TI-3504	3" x 5" (8 cm x 13cm)



**PRIORITY**

Room # \_\_\_\_\_

Time Needed \_\_\_\_\_

Initials \_\_\_\_\_

WWW.BATRIK.COM TAG-IT™ TI-3504

### BROKEN / RED

ITEM #	SIZE
TI-35B4	3" x 5" (8 cm x 13cm)



**BROKEN DO NOT USE**

DATE: \_\_\_\_\_ PHONE: \_\_\_\_\_

DEPARTMENT/CONTACT: \_\_\_\_\_

ITEM: \_\_\_\_\_

CRITICAL: (☠) Yes No

INSTRUMENT SET: \_\_\_\_\_

ASSEMBLED BY: \_\_\_\_\_

WWW.BATRIK.COM TAG-IT™ TI-35B4

### RE-RUN / YELLOW

ITEM #	SIZE
TI-3505	3" x 5" (8 cm x 13cm)



**RE-RUN**

Room # \_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

WWW.BATRIK.COM TAG-IT™ TI-3505

### VENDOR IN + OUT / LIGHT BLUE

ITEM #
TI-35C22

SIZE
5" x 3" (13 cm x 8 cm)



**VENDOR**

Check In

Tray Received From: \_\_\_\_\_

Tray Name: \_\_\_\_\_

Received By: \_\_\_\_\_

Date: \_\_\_\_\_ Surgery Date: \_\_\_\_\_

☐ CPD Notified Prior to Receipt of Trays

☐ Received >24 Hours Before Case

☐ Inventory List Available

☐ Inventory Check Complete

☐ All Instruments in Good Condition

☐ Containers in Good Condition

☐ All Trays Labeled (Patient & Surgeon)

☐ All Trays <25 pounds

☐ Written Instructions For Use Available

☐ Disassembled/Decontaminated

WWW.BATRIK.COM TAG-IT™ TI-35C22

**VENDOR**

Check Out

☐ Disassembled/Decontaminated

☐ Inventory Check Complete

☐ All Items Complete and Accounted For

Returned By: \_\_\_\_\_

Received By: \_\_\_\_\_

Date: \_\_\_\_\_

Notes: \_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

WWW.BATRIK.COM TAG-IT™ TI-35C22

### MISSING / PINK

ITEM #	SIZE
TI-3506	3" x 5" (8 cm x 13cm)



**MISSING**

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_


\_\_\_\_\_

WWW.BATRIK.COM TAG-IT™ TI-3506

### LOANER IN + OUT / GREEN

ITEM #
TI-35L3

SIZE
5" x 3" (13 cm x 8 cm)



**LOANER**

Check In

Tray Received From: \_\_\_\_\_

Tray Name: \_\_\_\_\_

Received By: \_\_\_\_\_

Date: \_\_\_\_\_ Surgery Date: \_\_\_\_\_

☐ CPD Notified Prior to Receipt of Trays

☐ Received >24 Hours Before Case

☐ Inventory List Available

☐ Inventory Check Complete

☐ All Instruments in Good Condition

☐ Containers in Good Condition

☐ All Trays Labeled (Patient & Surgeon)

☐ All Trays <25 pounds

☐ Written Instructions For Use Available

☐ Disassembled/Decontaminated

WWW.BATRIK.COM TAG-IT™ TI-35L3

**LOANER**

Check Out

☐ Disassembled/Decontaminated

☐ Inventory Check Complete

☐ All Items Complete and Accounted For

Returned By: \_\_\_\_\_

Received By: \_\_\_\_\_

Date: \_\_\_\_\_

Notes: \_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

WWW.BATRIK.COM TAG-IT™ TI-35L3

### BROKEN/ ORANGE

ITEM #	SIZE
TI-3507	3" x 5" (8 cm x 13cm)



**BROKEN DO NOT USE**

DATE: \_\_\_\_\_ PHONE: \_\_\_\_\_

DEPARTMENT/CONTACT: \_\_\_\_\_


ITEM: \_\_\_\_\_

CRITICAL: (☠) Yes No

WWW.BATRIK.COM TAG-IT™ TI-3507

### BLANK / BLACK

ITEM #	SIZE
TI-3509	3" x 5" (8 cm x 13cm)



\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_


\_\_\_\_\_

WWW.BATRIK.COM TAG-IT™ TI-3509

### QUICK TURNOVER

### QUICK TURNOVER/ BLUE

ITEM #	SIZE
TI-35T2	3" x 5" (8 cm x 13cm)



**QUICK TURNOVER**

Tech Initials: \_\_\_\_\_

Room Number: \_\_\_\_\_

Item Number: \_\_\_\_\_

Time Needed: \_\_\_\_\_

☐ Decontam/LD/LD/HLD

☐ Immediate Use Sterilization

☐ Terminal Sterilization in SPD

WWW.BATRIK.COM TAG-IT™ TI-35T2

All information tags are single-use. PACKAGING: 1 pack of 100 tags, 100 tie wraps, and 1 permanent marker.

NEW!

BATRIK

## TAG-IT™ | SCOPE

Batrik Scope Tags are **tear-resistant, double-sided** tags that will prove indispensable throughout the reprocessing journey. Our tags make it easy to **record important information**, including processing dates, with a user-friendly feature that allows you to select the month, day, and date. Crafted from **waterproof synthetic paper**, these tags feature a **self-looping design**, making it effortless to attach them to the scope without the need for additional fastening materials. Elevate your reprocessing protocols with Batrik!

### KEY FEATURES

- The paper, tie wraps, and smudge-proof ink are designed to withstand high temperatures of up to 273.2°F (134°C) and are compatible with Gas plasma sterilization (STERRAD®100 NX Sterilization System).
- Tags include the month, day, and date selection on one side, making it convenient to mark the processing date accurately



### COLOR-CODED AND CUSTOMIZABLE

- Available in seven colors
- The tags facilitate color coding for easy identification

### SELF LOOPING

- The self-looping design allows for easy attachment to the scope
- Eliminating the need for additional fastening materials



### SINGLE-USE AND LATEX FREE

- Batrik Scope Tags are designed for single use and are latex-free
- Promoting hygiene and reducing the risk of cross-contamination

### PRODUCT SELECTION

ITEM #	COLOR	SIZE	PACKAGING
TI-ST05	Yellow	11.8" x 1.55" / 33.8 cm x 4 cm	100 / PKG
TI-ST02	Blue	11.8" x 1.55" / 33.8 cm x 4 cm	100 / PKG
TI-ST03	Green	11.8" x 1.55" / 33.8 cm x 4 cm	100 / PKG
TI-ST07	Orange	11.8" x 1.55" / 33.8 cm x 4 cm	100 / PKG
TI-ST06	Purple	11.8" x 1.55" / 33.8 cm x 4 cm	100 / PKG
TI-ST04	Red	11.8" x 1.55" / 33.8 cm x 4 cm	100 / PKG
TI-ST01	White	11.8" x 1.55" / 33.8 cm x 4 cm	100 / PKG



All information tags are single-use. PACKAGING: 1 pack of 100 tags and 1 permanent marker.



**NEW!****BATRIK**

# TAG-IT™ | SCOPE REPROCESSING

Batrik Scope Reprocessing Tags are durable and can withstand handling during reprocessing **without tearing**, ensuring that information remains intact. These tags are vital in recording important information such as **date to reprocess on, scope model and serial number, sterilizer type, meeting sterilizer parameter conditions, alcohol flush, and hang date**. These tags have been meticulously designed to elevate **operational efficiency** and contribute to the delivery of the highest standard of patient care.

**TI-RST****TI-MPB**

\*Suitable for various information tagging purposes

## LOGGING IMPORTANT INFORMATION

- Easy for staff to indicate when the scope was processed
- Key scope, technician, and cycle-related information easily reported and identified.

## SMART DESIGN

- The self-looping design allows for easy attachment to the scope
- No requirement for extra fastening materials
- Secure fastening as the loop is tightened around the instrument
- Features an integrated barcode area



## PRODUCT SELECTION



ITEM #	DESCRIPTION	DIMENSIONS	PACKAGING
<b>TI-RST</b>	REPROCESSING SCOPE TAG	3" x 10.5" / 8 cm x 27 cm	100 / PKG
<b>TI-MPB</b>	MULTI-PURPOSE BLANK TAG	3" x 10.5" / 8 cm x 27 cm	100 / PKG

All information tags are single-use. **PACKAGING:** 1 pack of 100 tags and 1 permanent marker.


# BATRIK SILICONE PIN MAT

Our Silicone Pin Mat cushions, protects, and organizes delicate surgical instruments for sterilization. The elevated nubs provide greater air-flow ventilation and drying, significantly reducing residual moisture collection.

- Reusable
- Autoclavable
- Latex free



**FRONT SIDE**  
Pins keep instruments securely in place



**BACK SIDE**  
Numerous holes and elevated nubs provide superior venting and reduce water collection

## PRODUCT SELECTION

ITEM #	DIMENSIONS	PACKAGING
<b>13-2001</b>	13" x 18" (33 cm x 45.72 cm) 11.97 mm high from base of mat to end of pin 15.53 mm high from base of elevated nub to end of pin	1 unit
<b>13-2005</b>	13" x 18" (33 cm x 45.72 cm) 11.97 mm high from base of mat to end of pin 15.53 mm high from base of elevated nub to end of pin	5 units

GOOD TO KNOW:  
CUSTOM SIZES  
ARE AVAILABLE  
UPON REQUEST

## BATRIK CORNER POCKETS™


With three models to choose from, our Corner Pockets™ portfolio is a truly comprehensive line of products designed to prevent tearing associated with sharp corners or feet on surgical trays.



18-CP02



18-CP07




18-CP06

COLOUR	BLUE	ORANGE	PINK
USAGE	SINGLE-USE	REUSABLE <small>withstands repeated sterilization up to 134°C</small>	REUSABLE <small>withstands repeated sterilization up to 134°C</small>
DIMENSIONS	2.8" x 4" (7.3 cm x 10.3 cm)	2.8" x 4" (7.3 cm x 10.3 cm)	3.4" x 4.7" (8.5 cm x 12 cm)
LATEX-FREE	✓	✓	✓
NO MOISTURE RETENTION	✓	✓	✓
FABRICATION	Flexible, silicone-like plastic	Flexible, silicone-like plastic	Our most flexible plastic for easiest placement and most secure protection

## SOFT POCKETS™

Our Soft Pockets™ Foam Corner Protectors are made from high density foam to effectively cushion and protect the corners of surgical trays to prevent tearing. They are especially useful with heavier instrument sets.



### PRODUCT SELECTION

ITEM #	DESCRIPTION	PACKAGING	CASE DIMENSIONS
18-CP02**	Corner Pockets™ Single-Use, Blue	240 / Bag	80 / Box (24 Boxes / Pallet)
18-EVAL2**	Corner Pockets™ Single-Use, Blue	24 / Bag	-
18-CP07	Corner Pockets™ Reusable, Orange	80 / Bag	-
18-CORNER	Corner Pockets™ Reusable, Orange	320 / Box (4 Bags of 80 pcs)	-
18-EVAL7	Corner Pockets™ Reusable, Orange	24 / Bag	-
18-CP06	Corner Pockets™ Reusable, Pink	80 / Bag	-
18-CPF-02	Soft Pockets™ Single-Use Foam Corner Protectors	80 / Box (24 Boxes / Pallet)	-

\*\* Although the single use Corner Pockets™ appear similar to the reusable models they are blended differently and have only been validated as single-use. Re-using them may compromise sterility.

# LILA'S INSTRUMENT STAIN REMOVER

Lila's special blend of naturally sourced powders provides outstanding removal of stubborn stains.

Convenient-to-use formula, you can achieve spotless results every time!

It can be used on surgical/dental instruments, autoclaves, ultrasonic cleaners, stainless steel countertops, trays, bowls and more.

Specifically designed for use on stainless steel, chrome, nickel, brass, copper, and silver.



LSR-8501


## EASE TO USE AS 1,2,3...

1. Dampen a cloth or grab a Lila's Low Linting Dura-Towel and dab into the jar.
2. Gently rub the area that is stained.<sup>1</sup>
3. Rinse thoroughly for a flawless shine.<sup>2</sup>

## A little Goes a long Way!

Even great to tackle stubborn stains on steel pots and pans in your kitchen with ease.

Lila's - If you have used Spectra-Scrub or Integra's instrument stain remover then you will love Lila's!



60-46003



17-DT1213C

## LILA'S SCRUB PAD

- Say goodbye to extra stubborn stains with Lila's Scrub Pad without leaving any metal residue.
- Simply add a small amount of Lila's Stain Remover with one of our Scrub Pad's for very large areas.
- The gentle yet effective scrub is the perfect choice for Ophthalmic Instruments.<sup>1</sup>

## LILA'S LOW-LINTING DURA-TOWEL

- Lila's Dura-Towel has long lasting strength and has a minimal lint composition.
- These single-use towels are clean room certified and are specifically designed for instrument wiping.
- Its ideal dimensions and unique attributes guarantee optimal performance, setting it apart from ordinary paper towels.

## PRODUCT SELECTION

ITEM #	DESCRIPTION	DIMENSIONS	PACKAGING
LSR-8501	LILA'S INSTRUMENT STAIN REMOVER 3 oz (85g) plastic jar with a convenient screw-on lid	Ø 2.8" / 70 mm	1 Jar / Bag
LSR-8512	LILA'S INSTRUMENT STAIN REMOVER 3 oz (85g) plastic jar with a convenient screw-on lid	Ø 2.8" / 70 mm	12 Jars / Box
17-DT1213C	LILA'S LOW LINTING DURA-TOWEL 1080 towels	12" x 13" / 31 cm x 33 cm Folded: 3" x 3.25" / 7.5 cm x 8.25 cm	90 / Bag 12 Bags / Box
60-46003	LILA'S SCRUB PAD	4" x 6" / 10 cm x 15 cm	3 / Bag
LSR-BUNDLE	LILA'S STAIN REMOVER + LILA'S DURA-TOWEL 3 oz (85g) plastic jar + 90 towels	Ø 2.8" / 70 MM 12" x 13" / 31 cm x 33 cm	1 Jar / Bag 90 / Bag




1 : Although Lila's is completely safe, it is not recommended for use on ophthalmic instruments. Use Lila's Scrub Pad on its own.

2: The stain remover should not be mixed with bleach and is for external use only to remove instrument stains.

**BATRIK****INSTRUMENT PROTECTORS**



Our Instrument Protectors protect the tips and edges of instruments to prevent damage during handling and sterilization. Available in a wide array of color coded sizes for almost every application, including osteotome, retractor, double, and perforated.

- Steam and gas permeable
- Perforated protectors reduce vacuum seal and moisture build-up
- Validated for pre-vacuum sterilization\*
- Multiple color-coded sizes
- Latex Free
- Non-Sterile

**REGULAR / CLEAR****PERFORATED / CLEAR****REGULAR / TINTED****PERFORATED / TINTED****REGULAR / SOLID COLOR****PERFORATED / SOLID COLOR**

\*Independent laboratory tested to validate pre-vacuum steam sterilization efficacy in a healthcare setting


## PRODUCT SELECTION

	ITEM # REGULAR	ITEM # PERFORATED	COLOR	DIMENSIONS
 RC-07  PC-07	<b>REGULAR &amp; PERFORATED / CLEAR</b>			
	RC-01	PC-01	Clear	0.062" x 0.75" (1.57 mm x 19 mm)
	RC-02	PC-02	Clear	0.075" x 0.75" (2 mm x 19 mm)
	RC-03	PC-03	Clear	0.109" x 0.75" (2.77 mm x 19 mm)
	RC-04	PC-04	Clear	0.125" x 1" (3.18 mm x 25.4 mm)
	RC-05	PC-05	Clear	0.187" x 1" (4.75 mm x 25.4 mm)
	RC-06	PC-06	Clear	0.406" x 0.75" (10.3 mm x 19 mm)
	RC-07	PC-07	Clear	0.062" x 0.375" x 1" (1.57 mm x 9.53 mm x 25.4 mm)
	RC-08	PC-08	Clear	0.078" x 0.625" x 1" (2 mm x 15.875 mm x 25.4 mm)
	RC-09	PC-09	Clear	0.125" x 1" x 1" (3.18 mm x 25.4 mm x 25.4 mm)
	RC-10	PC-10	Clear	Assorted
	PACKAGING: Bag of 100			
 RT-07  PT-07	<b>REGULAR &amp; PERFORATED / TINTED</b>			
	RT-01	PT-01	Tinted White	0.062" x 0.75" (1.57 mm x 19 mm)
	RT-02	PT-02	Tinted Blue	0.075" x 0.75" (2 mm x 19 mm)
	RT-03	PT-03	Tinted Green	0.109" x 0.75" (2.77 mm x 19 mm)
	RT-04	PT-04	Tinted Red	0.125" x 1" (3.18 mm x 25.4 mm)
	RT-05	PT-05	Tinted Yellow	0.187" x 1" (4.75 mm x 25.4 mm)
	RT-06	PT-06	Tinted Maroon	0.406" x 0.75" (10.3 mm x 19 mm)
	RT-07	PT-07	Tinted Orange	0.062" x 0.375" x 1" (1.57 mm x 9.53 mm x 25.4 mm)
	RT-08	PT-08	Tinted Brown	0.078" x 0.625" x 1" (2 mm x 15.875 mm x 25.4 mm)
	RT-09	PT-09	Tinted Black	0.125" x 1" x 1" (3.18 mm x 25.4 mm x 25.4 mm)
	RT-10	PT-10	Tinted Assorted	Assorted
	PACKAGING: Bag of 100			
 RS-07  PS-07	<b>REGULAR &amp; PERFORATED / SOLID</b>			
	RS-01	PS-01	Solid White	0.062" x 0.75" (1.57 mm x 19 mm)
	RS-02	PS-02	Solid Blue	0.075" x 0.75" (2 mm x 19 mm)
	RS-03	PS-03	Solid Green	0.109" x 0.75" (2.77 mm x 19 mm)
	RS-04	PS-04	Solid Red	0.125" x 1" (3.18 mm x 25.4 mm)
	RS-05	PS-05	Solid Yellow	0.187" x 1" (4.75 mm x 25.4 mm)
	RS-06	PS-06	Solid Maroon	0.406" x 0.75" (10.3 mm x 19 mm)
	RS-07	PS-07	Solid Orange	0.062" x 0.375" x 1" (1.57 mm x 9.53 mm x 25.4 mm)
	RS-08	PS-08	Solid Brown	0.078" x 0.625" x 1" (2 mm x 15.875 mm x 25.4 mm)
	RS-09	PS-09	Solid Black	0.125" x 1" x 1" (3.18 mm x 25.4 mm x 25.4 mm)
	RS-10	PS-10	Solid Assorted	Assorted
	PACKAGING: Bag of 100			

ITEM #	DIMENSIONS
<b>OP-11</b>	0.25" x 0.75" (6.35 mm x 19.05 mm)
<b>OP-12</b>	0.375" x 1" (9.53 mm x 25.4 mm)
<b>OP-13</b>	0.5" x 1" (12.7 mm x 25.4 mm)
<b>OP-14</b>	0.75" x 1" (19.05 mm x 25.4 mm)
<b>OP-15</b>	1" x 1" (25.4 mm x 25.4 mm)
<b>OP-16</b>	1.25" x 1.25" (31.75 mm x 31.75 mm)
<b>OP-17</b>	1.5" x 1.5" (38.1 mm x 38.1 mm)
<b>OP-18</b>	1.75" x 1.25" (44.45 mm x 31.75 mm)
<b>OP-20</b>	Assorted


PACKAGING: Bag of 50

### OSTEOTOME PROTECTORS



ITEM #	DIMENSIONS	PACKAGING
<b>RGS-01</b>	1.5" (38 mm)	55 / Bag
<b>RGS-02</b>	2" (51 mm)	45 / Bag
<b>RGS-03</b>	3" (76.2 mm)	40 / Bag


### SINGLE-TIP RETRACTOR PROTECTORS



ITEM #	COLOR	SIZE
<b>RGM-01</b>	Pink	Small
<b>RGM-02</b>	Clear	Small
<b>RGM-03</b>	Green	Medium
<b>RGM-04</b>	Clear	Medium

PACKAGING: Bag of 50


### WEITLANER RETRACTOR PROTECTORS



ITEM #	COLOR	DIMENSIONS
<b>DG-01</b>	Clear	0.0625" x 0.75" (1.59 mm x 19.05 mm)
<b>DG-02</b>	Blue	0.078" x 0.75" (1.98 mm x 19.05 mm)
<b>DG-03</b>	Green	0.109" x 0.75" (2.78 mm x 19.05 mm)
<b>DG-04</b>	Red	0.125" x 1" (3.18 mm x 25.4 mm)
<b>DG-05</b>	Yellow	0.187" x 1" (4.76 mm x 25.4 mm)
<b>DG-07</b>	Orange	0.062" x 0.375" x 1" (1.59 mm x 9.53 mm x 25.4 mm)
<b>DG-10</b>	Assorted	Assorted

PACKAGING: Bag of 50

### PERFORATED DOUBLE PROTECTORS



# BATRIK PrAmLi Power<sup>3</sup> LUBE™ Instrument Lubricant

Specially formulated to lubricate surgical instruments, hinges, box locks, and general equipment during autoclaving. PrAmLi Power<sup>3</sup> LUBE™ is uniquely formulated with an all natural corrosion inhibitor.

- Forms a protective coating which can extend the instrument's, useful life
- Bacteriostatic
- Neutral pH 6.5 - 7.5
- Prevents corrosion

**TIP:**  
FOR MILK BATH SOLUTION,  
DILUTE 1 PART SOLUTION  
IN 6 PARTS WATER

ITEM #	DIMENSIONS
OP-11	0.25" x 0.75" (6.35 mm x 19.05 mm)
OP-12	0.375" x 1" (9.53 mm x 25.4 mm)
OP-13	0.5" x 1" (12.7 mm x 25.4 mm)
OP-14	0.75" x 1" (19.05 mm x 25.4 mm)
OP-15	1" x 1" (25.4 mm x 25.4 mm)
OP-16	1.25" x 1.25" (31.75 mm x 31.75 mm)
OP-17	1.5" x 1.5" (38.1 mm x 38.1 mm)
OP-18	1.75" x 1.25" (44.45 mm x 31.75 mm)
OP-20	Assorted

PACKAGING: Bag of 50




## PRODUCT SELECTION

ITEM #	DILUTION LUBRICANT / WATER	PACKAGING	CASE DIMENSIONS
SS-21-6001	1 to 4 ml / L (0.125 to 0.5 oz / Gallon)	1 x 3.8 L (1 Gallon)	6" x 6" x 13"
SS-21-6004	1 to 4 ml / L (0.125 to 0.5 oz / Gallon)	4 x 3.8 L (4 Gallons)	13" x 13" x 14"
SS-21-6025	1 to 4 ml / L (0.125 to 0.5 oz / Gallon)	2 x 9.5 L (2 x 2.5 Gallons)	14.5" x 9.25" x 15"
SS-21-6032	No dilution required	1 x 946 ml (32 oz) Spray bottle	10" x 4"
SS-21-6012	No dilution required	12 X 946 ml (32 oz) Spray Bottles	18" x 12" x 11"

# BATRIK S.T.I.C.S.

Our Simple Tracking Identification Color System, or S.T.I.C.S., is an easy to use, latex free, tool that allows you to quickly identify your instruments by department and/or procedure. The tape adheres securely and every tape sheet comes with a protective cover to ensure it lasts as long as you need.



Each 8.5" x 11" (216 mm x 279 mm) sheet contains the following adhesive strips:

- 1" (25.4 mm) x 220
- 1.5" (38.1 mm) x 44
- 2" (50.8 mm) x 44
- Total tape per sheet 374" (9500 mm)

## PRODUCT SELECTION

SOLID COLORS	
ITEM #	COLOR
T-850-2-01	White
T-850-2-02	Royal Blue
T-850-2-03	Kelly Green
T-850-2-04	Red
T-850-2-05	Yellow
T-850-2-06	Purple
T-850-2-07	Orange
T-850-2-08	Brown
T-850-2-09	Black
T-850-2-10	Gray
T-850-2-11	Lavender
T-850-2-12	Beige
T-850-2-13	Light Blue
T-850-2-14	Lime Green
T-850-2-15	Pink
T-850-2-17	Plum
T-850-2-18	Mint Green

BLACK STRIPES	
ITEM #	COLOR
T-850-3-01	Black/White
T-850-3-03	Black/Kelly Green
T-850-3-04	Black/Red
T-850-3-05	Black/Yellow
T-850-3-07	Black/Orange
T-850-3-10	Black/Gray
T-850-3-11	Black/Lavender
T-850-3-12	Black/Beige
T-850-3-13	Black/Light Blue
T-850-3-14	Black/Lime Green
T-850-3-15	Black/Pink
T-850-3-17	Black/Plum

TRI-COLOR	
ITEM #	COLOR
T-850-5-51	Red/White/Blue
T-850-5-54	Red/Black/Green

SPECIALTY	
ITEM #	COLOR
T-850-7-04	Yellow/Red


WHITE STRIPES	
ITEM #	COLOR
T-850-4-02	White/Royal Blue
T-850-4-03	White/Kelly Green
T-850-4-04	White/Red
T-850-4-05	White/Yellow
T-850-4-07	White/Orange
T-850-4-08	White/Brown
T-850-4-11	White/Lavender
T-850-4-13	White/Light Blue
T-850-4-14	White/Lime Green
T-850-4-15	White/Pink
T-850-4-17	White/Plum

FOUR-COLOR	
ITEM #	COLOR
T-850-2-31	Red/Yellow/ Purple/Orange
T-850-2-32	Blue/White/ Black/Green

# BATRIK ID TAPE



## AVAILABLE WIDTHS



A simple solution for marking surgical instruments with the aid of color-coding, our ID Tape is durable and will not chip, flake, or peel. It can be used in autoclaves, withstanding repeated exposure to sterilization. Independent validation tests demonstrated the efficiency of sterilization of marked instruments, even underneath the BATRIK ID Tape.



Quickly identifies instruments by department or procedure



Easy to apply



Withstands repeated exposure to sterilization



Has more per tape roll (300" / 7620 mm)



Available in 7 colors and 2 widths



Conveniently packaged in a clear box



Non-Sterile

## PRODUCT SELECTION

ITEM # 0.125" / 318 MM	ITEM# 0.25" / 635 MM	COLOR
T-125-01	T-250-01	White
T-125-02	T-250-02	Blue
T-125-03	T-250-03	Green
T-125-04	T-250-04	Red
T-125-05	T-250-05	Yellow
T-125-07	T-250-07	Orange
T-125-09	T-250-09	Black

**BATRIK**

# TAPE-B-GONE™ TAPE REMOVER

**Tape-B-Gone™** is a uniquely formulated and ready-to-use proprietary formula that easily removes adhesives from medical equipment and instruments. Unlike many other brands, Tape-B-Gone™ does not contain highly dangerous or caustic chemicals, is not flammable or combustible, and is safe for use on metal, glass, plastic, and rubber without any special equipment or ventilation.



## PRODUCT SELECTION

ITEM #	DESCRIPTION	PACKAGING
<b>T-08RMV</b>	8 oz (237 ml) bottle	1 Unit
<b>T-32RMV</b>	8 oz (237 ml) bottle	4 / Pack
<b>T-SAMP-RMV</b>	0.5 oz (15 ml) bottle	1 Unit
<b>T-96RMVSC-R0</b>	8oz (237 ml) bottle Screwcap*	12 / Box

\*Samples will be available in Batrik by April 21st

BATRIK


# PADLOCKS SECURITY TAGS

Our Padlock Security Tags for instrument containers and general hospital applications are specifically designed for the healthcare market. Our material selection, fabrication process, and engineering ensure consistent quality and high temperature compliance with the various steam sterilization temperatures in worldwide hospital applications.


Available with optional reactive ink indicator dots that provide visual sterility indication. Light blue dots turn brown/dark brown when steam sterilized.

## AVAILABLE IN THREE FORMATS


### 700-SERIES



### 800-SERIES



### 900-SERIES



Made of polypropylene



Latex and DEHP Free



Will not melt or soften with any established steam sterilization protocol




Fully open design (45°) facilitates insertion into instrument containers



Easy twist breaking design



Large base for easier handling



Available in different size



Custom Colors and OEM available upon request

## PRINTED & CUSTOM NUMBERED PADLOCKS

Available in 700, 800 & 900 SERIES

### PRODUCT SELECTION

COLOR	PRINTED WORD
<b>PRINTED PADLOCKS</b>	
WHITE	CLEAN
BLUE	STERILE
RED	REPAIR
YELLOW	SHARPEN
YELLOW	LOANER
BLACK	DISCARD
<b>CUSTOM NUMBERED</b>	
ALL COLORS	up to 8 numbers




**CUSTOM NUMBERED**  
Custom printed padlocks (up to 8 characters)



**PRINTED**  
6 different printed labels available

NOTE: Printed & Custom Numbered Padlocks are available in 700, 800, and 900 Series and can be ordered in 100, 200, and 500 unit quantities. Please contact us for ordering details.

INDICATOR  
DOTS

## REGULAR



## 700-SERIES PADLOCKS

The 700 series padlocks are perfect for "Aesculap" containers and those with a short, large opening.

## PRODUCT SELECTION

## REGULAR PADLOCKS

COLOR	100 UNITS	200 UNITS	500 UNITS
WHITE	15-701	15-711	15-710
BLUE	15-702	15-712	15-720
GREEN	15-703	15-713	15-730
RED	15-704	15-714	15-740
YELLOW	15-705	15-715	15-750
PINK	15-706	15-716	15-760
ORANGE	15-707	15-717	15-770
BLACK	15-709	15-719	15-790
LIGHT BLUE	15-722	15-732	15-742

## WITH INDICATORS DOTS

COLOR	100 UNITS	200 UNITS	500 UNITS
WHITE	15-701-IND	15-711-IND	15-710-IND
BLUE	15-702-IND	15-712-IND	15-720-IND
GREEN	15-703-IND	15-713-IND	15-730-IND
RED	15-704-IND	15-714-IND	15-740-IND
YELLOW	15-705-IND	15-715-IND	15-750-IND
PINK	15-706-IND	15-716-IND	15-760-IND
ORANGE	15-707-IND	15-717-IND	15-770-IND
BLACK	15-709-IND	15-719-IND	15-790-IND
LIGHT BLUE	15-722-IND	15-732-IND	15-742-IND

Double locking hooks can be inserted in a large guide hole



The **grooves** on the padlock facilitate handling



## INDICATOR DOTS



## REGULAR



# 800-SERIES PADLOCKS

The 800 series padlocks have a longer arm (32 mm) and have a single, thinner hook than the 700 series. They are perfect for almost any type of container, even those with long and thin channels.

## PRODUCT SELECTION

### REGULAR PADLOCKS

COLOR	100 UNITS	200 UNITS	500 UNITS
WHITE	15-801	15-811	15-810
BLUE	15-802	15-812	15-820
GREEN	15-803	15-813	15-830
RED	15-804	15-814	15-840
YELLOW	15-805	15-815	15-850
PINK	15-806	15-816	15-860
ORANGE	15-807	15-817	15-870
BLACK	15-809	15-819	15-890
LIGHT BLUE	15-822	15-832	15-842

### WITH INDICATORS DOTS

COLOR	100 UNITS	200 UNITS	500 UNITS
WHITE	15-801-IND	15-811-IND	15-810-IND
BLUE	15-802-IND	15-812-IND	15-820-IND
GREEN	15-803-IND	15-813-IND	15-830-IND
RED	15-804-IND	15-814-IND	15-840-IND
YELLOW	15-805-IND	15-815-IND	15-850-IND
PINK	15-806-IND	15-816-IND	15-860-IND
ORANGE	15-807-IND	15-817-IND	15-870-IND
BLACK	15-809-IND	15-819-IND	15-890-IND
LIGHT BLUE	15-822-IND	15-832-IND	15-842-IND

**Single locking hook**  
can be inserted  
in a large guide hole



The locking system  
allows for **easy closure**



INDICATOR  
DOTS

## REGULAR



## 900-SERIES PADLOCKS

The 900 series padlocks are identical to the 800 series, but with a larger rectangular base which allows for the placement of a 20 mm x 30 mm tracking label.

## PRODUCT SELECTION

## REGULAR PADLOCKS

COLOR	100 UNITS	200 UNITS	500 UNITS
WHITE	15-901	15-911	15-910
BLUE	15-902	15-912	15-920
GREEN	15-903	15-913	15-930
RED	15-904	15-914	15-940
YELLOW	15-905	15-915	15-950
PINK	15-906	15-916	15-960
ORANGE	15-907	15-917	15-970
BLACK	15-909	15-919	15-990
LIGHT BLUE	15-922	15-932	15-942

## WITH INDICATORS DOTS

COLOR	100 UNITS	200 UNITS	500 UNITS
WHITE	15-901-IND	15-911-IND	15-910-IND
BLUE	15-902-IND	15-912-IND	15-920-IND
GREEN	15-903-IND	15-913-IND	15-930-IND
RED	15-904-IND	15-914-IND	15-940-IND
YELLOW	15-905-IND	15-915-IND	15-950-IND
PINK	15-906-IND	15-916-IND	15-960-IND
ORANGE	15-907-IND	15-917-IND	15-970-IND
BLACK	15-909-IND	15-919-IND	15-990-IND
LIGHT BLUE	15-922-IND	15-932-IND	15-942-IND

**Single locking hook**  
for use with smaller  
diameters



Ideal for adding labels/  
indicators indicators or  
markings to the large  
rectangular base



# SAVEFACE™ DRAPED SHIELDS

**REDUCE PERSONAL INJURY**

**SUPERIOR INSTRUMENT DECONTAMINATION PROTECTION**

## SAVEFACE DRAPED-FACE SHIELD

In the Sterile Processing Department (SPD) Central Sterile Supply Department (CSSD), personal safety is paramount.

Our Draped Face Shield provides essential protection during high-risk instrument decontamination.

## EFFECTIVE PROTECTION IS NON-NEGOTIABLE

SAVEFACE Draped dynamically reduces harmful exposure to caustic chemicals as well as the risk of cross-contamination. Effective face protection is non-negotiable.

## PPE CATEGORY 2 CERTIFIED

SAVAFACE DRAPED is tested to protect. Wear it with confidence.



UKAS PRODUCT  
CERTIFIED



PATENT PENDING



INCLUDES CLOSING SYSTEM TO  
REDUCE RISK OF FLUID CONTACT  
FROM BELOW THE SHIELD

## SAVEFACE: KEY FEATURES

### ANSI/AAMI PB70 LEVEL 3 COMPLIANT DRAPE:

- Provides a high level of fluid penetration protection.
- Tuck into clothing or pinch together.

### EXTRA-LONG 9"/23CM WRAP-AROUND SHIELD:

- Ensures full coverage, reducing the risk of exposure from different angles.

### COMFORTABLE:

- The soft, self-shaping foam headpiece molds comfortably to your forehead, allowing for extended wear without discomfort.

### OPTICAL CLARITY WITH ANTI-FOG AND ANTI-GLARE

ASSEMBLED IN CANADA

## DON'T MISS THE OPPORTUNITY TO EXPLORE OUR BROADER RANGE OF PROTECTIVE SHIELDS.

Visit our website [www.surgmed.com](http://www.surgmed.com) or consult our detailed literature to find the shield that best fits your needs.

### VERY IMPORTANT

#### IS YOUR PRESENT FACE SHIELD CATEGORY 2 CERTIFIED?

Batrik face shields are independently certified to comply with the PPE Category 2 standard. Many facilities are now requesting compliance to this standard.

SHIELD LENGTH  
9" / 23 CM

DRAPE LENGTH  
6" / 15 CM



**SF-8240**

**SF-240**




## PRODUCT SELECTION

ITEM #	DESCRIPTION	DIMENSIONS	TOTAL LENGTH	PACKAGING	CASE DIMENSIONS
<b>SF-240</b>	BATRIK SAVEFACE™ no drape	Shield Length: 9" / 23cm Shield width: 13" / 33cm	9" (23cm)	240 / Case (8 Boxes of 30 Pcs)	26.5" x 21" x 15"
<b>SF-8240</b>	BATRIK SAVEFACE™ with ANSI/AAMI PB70 Level 3 compliant protection drape	Shield Length: 9" / 23cm Shield width: 13" / 33cm Shield width: 6" / 15cm	58 cm / 23"	240 / Case (8 Boxes of 30 Pcs)	26.5" x 21" x 15"

# BATRIK SCRUB PAD

Our Scrub Pad is ideal for pre-cleaning your surgical instruments without leaving any metal residue that often occurs when using steel wool or metal sponges. The abrasiveness is perfectly balanced to provide you with efficient cleaning without damaging your instruments. Its fine fibers allow the user to polish, clean, and debur all types of material surfaces.

- Abrasive Material: Aluminum Oxide
- Grit: 400
- Grade: Very Fine
- Style: Flexible
- Pad material: Non-Woven
- Color: Maroon
- Type: Hand Pad



## PRODUCT SELECTION

ITEM #	DESCRIPTION	DIMENSIONS	PACKAGING
<b>60-46050</b>	Scrub Pad	4" x 6" / 10 cm x 15 cm	150 / Box (50 bags of 3 units each)
<b>60-46003</b>	Scrub Pad	4" x 6" / 10 cm x 15 cm	3 / Bag

Your First ™

# ENZYMATIC & NON-ENZYMATIC CLEANERS

Our Enzymatic and Non-Enzymatic Cleaners are specially formulated to remove organic material from medical equipment and instrumentation.

Our **Enzymatic Cleaner** formats are formulated with PrAmLi Power3 Enzymes™, our proprietary unique blend of enzymes, and a surfactant to efficiently remove all proteins, blood, fatty tissue, and more. They also remove stains from surgical scrubs, linens, and medical clothing.

Our **Non-Enzymatic Cleaner** formats are formulated with a free-rinsing, neutral pH detergent blend for effective cleaning.



Both formats are powerful solutions that break down all bioburden for effective surface and instrument pre-cleaning\*

## PRODUCT SELECTION

ENZYMATIC CLEANER	NON-ENZYMATIC CLEANER	DILUTION DETERGENT / WATER	PACKAGING	CASE DIMENSIONS
<b>LOW FOAM</b>				
<b>SS-21-1001</b>	<b>SS-721-1001</b>	4 to 8 ml / L (0.5 to 1 oz / Gal)	1 x 3.8 L (1 Gal)	6" x 6" x 13"
<b>SS-21-1004</b>	<b>SS-721-1004</b>	4 to 8 ml / L (0.5 to 1 oz / Gal)	4 x 3.8 L (4 Gal)	13" x 13" x 14"
<b>MILD FOAM</b>				
<b>SS-21-2001</b>	<b>SS-721-2001</b>	4 to 8 ml / L (0.5 to 1 oz / Gal)	1 x 3.8 L (1 Gal)	6" x 6" x 13"
<b>SS-21-2004</b>	<b>SS-721-2004</b>	4 to 8 ml / L (0.5 to 1 oz / Gal)	4 x 3.8 L (4 Gal)	13" x 13" x 14"
<b>SUPER CONCENTRATE</b>				
<b>SS-21-3001</b>	<b>SS-721-3001</b>	1 to 2 ml / L (0.125 to 0.25 oz / Gal)	1 x 3.8 L (1 Gal)	6" x 6" x 13"
<b>SS-21-3004</b>	<b>SS-721-3004</b>	1 to 2 ml / L (0.125 to 0.25 oz / Gal)	4 x 3.8 L (4 Gal)	13" x 13" x 14"

### GOOD TO KNOW:

OUR LOW FOAM & SUPER CONCENTRATE DETERGENTS (ENZYMATIC & NON-ENZYMATIC) MAY BE USED IN WASHERS / ULTRASONIC CLEANERS.

WE RECOMMEND CONFIRMING SUGGESTED DILUTION RATIO AND WATER TEMPERATURE WITH THE MANUFACTURER OF YOUR WASHER / ULTRASONIC CLEANER AS MACHINES / SETTINGS VARY.


NOTE: Our Low Foam & Super Concentrate Detergent formulas are available in a 2.5 Gallon format based on a minimum order quantity. Contact us for more information.

\*Caution: We recommend consulting your hospital pre-cleaning/cleaning protocols for various equipment/instrument cleaning detergent preferences.

Your First ™

# ENVIRO-SACHET™

## ENZYMATIC GEL & NON-ENZYMATIC CLEANER

Your First ™ Enviro-Sachet™, our unique cleaning solutions now available in Enzymatic & Non-Enzymatic formats.



- Our Enzymatic Gel format is formulated with the same triple action **PrAmLi Power<sup>3</sup> Enzymes™** that you know and trust.
- Our Non-Enzymatic Cleaner format is formulated with a free-rinsing, neutral pH detergent blend.

Both formats are powerful solutions that break down all bioburden for effective surface and instrument pre-cleaning.

### PRODUCT SELECTION

ENZYMATIC GEL	NON-ENZYMATIC CLEANER	DESCRIPTION	PACKAGING	CASE DIMENSIONS
<b>17-166</b>	<b>717-166</b>	1 Dry Tubular Sponge & Enviro-Sachet™ set	25 / Box 4 Boxes / Case	21" x 14" x 13.5"
<b>17-266</b>	<b>717-266</b>	1 Dry Flat Sponge & Enviro-Sachet™ set	50 / Box 4 Boxes / Case	21" x 14" x 13.5"
<b>17-177</b>	<b>717-177</b>	1 Dry Tubular Sponge, Enviro-Sachet™ & Stainless Steel Bowl kit	1 Kit	7" x 6" x 3"
<b>17-277</b>	<b>717-277</b>	1 Dry Flat Sponge, Enviro-Sachet™ & Stainless Steel Bowl kit	1 Kit	7" x 6" x 3"
<b>17-EZSAMP</b>	<b>717-EZSAMP</b>	SAMPLE KIT INCLUDES: 3 Dry Flat Sponges with Enviro-Sachet™ sets 1 Dry Tubular Sponge with Enviro-Sachet™ set 1 Reusable Stainless Steel Bowl	1 Kit	7" x 6" x 3"
<b>H-EZ500</b>		Reusable Stainless Steel Bowl 4.3" x 5.3" x 3" / 11 cm x 13.5 cm x 7.5 cm	1 Unit	7" x 6" x 3"

Your First ™

# FOAMING SPRAY

## READY-TO-USE ENZYME FOAMING SPRAY

Our Foaming Spray is ideal for instrument pre-cleaning as the thick foam spray enzymatically breaks down all bioburden to allow for effective cleaning and sterilization.

Formulated with **PrAmLi Power<sup>3</sup> Enzymes™**, a unique blend of three powerful enzymes along with a surfactant to provide the ultimate in cleaning efficacy



## INSTRUCTIONS



## PRODUCT SELECTION

ITEM #	PACKAGING	CASE DIMENSIONS
<b>SS-21-4001</b>	1 x 946 ml (32 oz) Spray Bottle	10" x 4"
<b>SS-21-4012</b>	12 x 946 ml (32 oz) Spray Bottles	18" x 12" x 11"

Important: In addition to our enzymatic formulation, our spray bottle contains as much as 25% more solution when compared to the competition.

Your First ™

## F1RST WIPES

Our cleanroom wipes are designed to do the first wipe-down of an endoscope after a procedure. They are made from hydroentangled non-woven material composed of 45% polyester and 55% cellulose. These wipes are manufactured and packaged in a cleanroom environment.



- Can be used for a cleanroom environment of ISO Class 5 through 8
- Production process guarantees extremely low particle generation and chemical extractables
- Compatible with virtually any solvent or solution for cleaning and wiping (such as **PrAmLi Power<sup>3</sup> Enzymes™**)
- Double-bag packaging
- Available in two versatile sizes


## PRODUCT SELECTION

ITEM #	DESCRIPTION	DIMENSIONS	PACKAGING	CASE DIMENSIONS
<b>17-A0909</b>	Non-woven Cleanroom Wipes	9" x 9" (23 cm x 23 cm)	300 / Bag 14 Bags / Case	18.5" x 10" x 25.25"
<b>17-A1212</b>	Non-woven Cleanroom Wipes	12" x 12" (31 cm x 31 cm)	150 / Bag 10 Bags / Case	13" x 13" x 19"

Your First ™

# DRY SPONGES

Our sponges are specifically fabricated for use with all medical instrumentation. Unlike most competitor sponges, our sponges are white to better evidence the effectiveness of the cleaning process. Dark colored sponges easily hide contaminants, but with our white sponges you can visually evaluate the pre-clean efficacy.



OUR WHITE SPONGES MAY APPEAR SLIGHTLY YELLOWED. THIS IS ONLY DUE TO UV EXPOSURE, IS NOT A REFLECTION OF ITS AGE.

## PRODUCT SELECTION

ITEM #	DESCRIPTION	DIMENSIONS	PACKAGING
<b>17-100</b>	DRY TUBULAR SPONGE <a href="#">White Dispenser Box</a>	2.0" ø x 4.25" (2.7 cm ø x 10.8 cm)	25 / Box (4 Boxes / Case)
<b>17-101</b>	DRY TUBULAR SPONGE <a href="#">Bulk</a>	2.0" ø x 4.25" (2.7 cm ø x 10.8 cm)	560 / Case
<b>17-200</b>	DRY FLAT SPONGE <a href="#">White Dispenser Box</a>	4" x 5" x 0.5" (10.2 cm x 12.7 cm x 1.3 cm)	50 / Box (4 Boxes / Case)
<b>17-201</b>	DRY FLAT SPONGE <a href="#">Bulk</a>	4" x 5" x 0.5" (10.2 cm x 12.7 cm x 1.3 cm)	1,000 / Case

Your First ™

## ENZYMATIC BESIDE SPONGE KIT

### RightSwipe™

Introducing RightSwipe™, our newly launched Enzymatic Bedside Sponge Kits now feature sponges treated with an enhanced formula for a more effective pre-clean or surface clean solution. Our improved three-enzyme blend ensures thorough and safe cleaning and a simplified application process, making RightSwipe™ the perfect match for your endoscope care needs. The tubular sponge design is compatible with diverse endoscope models, providing versatility alongside the superior benefits of our enzymatic solution.



FLAT SPONGE



TUBULAR SPONGE



FREESTANDING POUCH



ENVIRO-SACHET™  
NON-ENZYMATIC CLEANER



DISPENSER BOX

### PRODUCT SELECTION

ITEM #	DESCRIPTION	PACKAGING	CASE DIMENSIONS
<b>17-111</b>	PRE-TREATED TUBULAR SPONGE Add 100-500 ml of water to activate detergent solution. Comes in plastic bag <sup>1</sup>	25 / Box (4 Boxes / Case)	21" x 14" x 13.5"
<b>17-211</b>	PRE-TREATED FLAT SPONGE Add 100-500 ml of water to activate detergent solution. Comes in plastic bag <sup>1</sup>	50 / Box (4 Boxes / Case)	21" X 14" X 13.5"
<b>17-122</b>	PRE-TREATED TUBULAR SPONGE Comes in a freestanding pouch to use as a basin. Simply add 500 ml of water to the pouch.	15 / Box (4 Boxes / Case)	21" X 14" X 13.5"
<b>17-222</b>	PRE-TREATED FLAT SPONGE Comes in a freestanding pouch to use as a basin. Simply add 500 ml of water to the pouch.	25 / Box (4 Boxes / Case)	21" X 14" X 13.5"
<b>17-133</b>	READY TO USE TUBULAR SPONGE Comes pre-treated in 500 ml of an enzymatic solution in freestanding pouch.	24 / Case	14" X 12" X 9"
<b>17-233</b>	READY TO USE FLAT SPONGE Comes pre-treated in 500 ml of an enzymatic solution in freestanding pouch.	24 / Case	14" X 12" X 9"
<b>17-144</b>	READY TO USE DUAL PURPOSE TUBULAR SPONGE Comes pre-treated in 100 ml of an enzymatic solution in freestanding pouch.	12 / Box (4 Boxes / Case)	21" X 14" X 13.5"
<b>17-244</b>	READY TO USE DUAL PURPOSE FLAT SPONGE Comes pre-treated in 100 ml of an enzymatic solution in freestanding pouch.	12 / Box (4 Boxes / Case)	21" X 14" X 13.5"
<b>17-155</b>	READY TO USE TUBULAR SPONGE Comes pre-treated in 220 ml of an enzymatic solution in freestanding pouch.	24 / Case	14" x 12" x 9"
<b>17-255</b>	READY TO USE FLAT SPONGE Comes pre-treated in 220 ml of an enzymatic solution in freestanding pouch.	24 / Case	14" X 12" X 9"
<b>17-277</b>	ENZYMATIC GEL Comes with 1 Dry Flat Sponge, Tubular Sponge, Enviro-Sachet™ & Stainless Steel Bowl kit	1 Kit	1 Kit


<sup>1</sup> Our stainless steel bowl can be sold separately and used as a basin for easy and convenient sponge activation.

Your First ™

# NON-ENZYMATIC BESIDE SPONGE KIT

## RightSwipe™

Presenting our enhanced Non-Enzymatic Bedside Sponge Kits. Each kit includes sponges pre-treated with a free-rinsing, neutral pH detergent blend from RightSwipe™, ensuring an effective pre-cleaning and surface cleaning solution for medical instruments. This kit stands as Your First™ step in the medical instrument cleaning process, showcasing the advanced benefits for superior instrument care. Our enhanced bedside kits will provide you with a streamlined workflow.



FLAT SPONGE



TUBULAR SPONGE



FREESTANDING POUCH

ENVIRO-SACHET™  
NON-ENZYMATIC CLEANER

DISPENSER BOX

## PRODUCT SELECTION

ITEM #	DESCRIPTION	PACKAGING	CASE DIMENSIONS
<b>717-122</b>	PRE-TREATED TUBULAR SPONGE Comes in a freestanding pouch to use as a basin. Simply add 500 ml of water to the pouch. <sup>1</sup>	15 / Box (4 Boxes / Case)	21" X 14" X 13.5"
<b>717-222</b>	PRE-TREATED FLAT SPONGE Comes in a freestanding pouch to use as a basin. Simply add 500 ml of water to the pouch. <sup>1</sup>	25 / Box (4 Boxes / Case)	21" X 14" X 13.5"
<b>717-133</b>	READY TO USE TUBULAR SPONGE Comes pre-treated in 500 ml of non-enzymatic solution in freestanding pouch.	24 / Case	14" X 12" X 9"
<b>717-233</b>	READY TO USE FLAT SPONGE Comes pre-treated in 500 ml of non-enzymatic solution in freestanding pouch.	24 / Case	14" X 12" X 9"
<b>717-144</b>	READY TO USE DUAL PURPOSE TUBULAR SPONGE Comes pre-treated in 100 ml of non-enzymatic solution in freestanding pouch.	12 / Box (4 Boxes / Case)	21" X 14" X 13.5"
<b>717-244</b>	READY TO USE DUAL PURPOSE FLAT SPONGE Comes pre-treated in 100 ml of non-enzymatic solution in freestanding pouch.	12 / Box (4 Boxes / Case)	21" X 14" X 13.5"
<b>717-155</b>	READY TO USE TUBULAR SPONGE Comes pre-treated in 220 ml of non-enzymatic solution in freestanding pouch.	24 / Case	14" x 12" x 9"
<b>717-255</b>	READY TO USE FLAT SPONGE Comes pre-treated in 220 ml of non-enzymatic solution in freestanding pouch.	24 / Case	14" X 12" X 9"
<b>717-277</b>	NON-ENZYMATIC GEL Comes with 1 Dry Flat Sponge, Tubular Sponge, Enviro-Sachet™ & Stainless Steel Bowl kit	1 Kit	1 Kit


<sup>1</sup> Our stainless steel bowl can be sold separately and used as a basin for easy and convenient sponge activation.

Your First ™

# TRANSPORT BAGS FOR ENDOSCOPE AND MEDICAL INSTRUMENTS

Thoughtfully designed to eliminate touch-points during endoscope transport from point-of-care to the reprocessing area, and back, focused on safety of patients and medical staff.

- Made of recyclable LDPE material known to be quite flexible, highly resilient, and tough
- Latex and DEHP-free
- Excellent resistance to dilute and concentrated acids, alcohols, bases and esters\*
- Good resistance to aldehydes, ketones, and vegetable-based oils\*
- Can withstand high temperatures of up to 90 °C (194°F)
- Tamper-evident seal for added safety, also aids in keeping the instrument moist to avoid drying of bio-debris
- Clearly marked bags with color-coded visuals facilitates immediate identification of clean vs. used instruments
- Users can write important information directly on the bag in spaces marked with visual icons
- Available in two convenient sizes


GOOD TO KNOW: USING **Your First ™** SEALABLE TRANSPORT BAGS REDUCES USER ERROR, PREVENTS CROSS-CONTAMINATION, AND PROMOTES PATIENT AND MEDICAL STAFF SAFETY.

## PRODUCT SELECTION


ITEM #	DESCRIPTION	DIMENSIONS	PACKAGING	CASE DIMENSIONS
<b>17-HAP24</b>	Clean Instruments Bag	24" x 22" / 61.0 cm x 56.0 cm	200 / Case	18" x 14" x 5"
<b>17-SAD24</b>	Soiled Instruments Bag	24" x 22" / 61.0 cm x 56.0 cm	200 / Case	18" x 14" x 5"
<b>17-HAP30</b>	Clean Instruments Bag	30.5" x 26" / 77.0 cm x 66.0 cm	200 / Case	21" x 19" x 7"
<b>17-SAD30</b>	Soiled Instruments Bag	30.5" x 26" / 77.0 cm x 66.0 cm	200 / Case	21" x 19" x 7"

\*Reference on file

### CLEAN INSTRUMENTS BAG



### SOILED INSTRUMENTS BAG



MADE IN CANADA 

BATRIK

# ARC™ WIRE STYLE BRUSHES

## SINGLE-USE AND REUSABLE

You've got options with our ARC™ Wire Brushes. Our single-use and reusable brushes are available in both flexible and rigid varieties and come in a multitude of sizes. ARC™ Wire Brushes are guaranteed to last longer than twisted wire and preserve the integrity of your valuable instruments.

### BRUSH SELECTION TIPS:

#### INSTRUMENT SPECIFICATIONS

##### What is the maximum length of the instrument(s) to be cleaned?

Add at least 2" (5cm) to that length, so you have an adequate length to work with.

##### What is the inner diameter of the instrument?

Add at most 0.125" (3mm) to determine the diameter of the bristle to ensure effective cleaning. The bristle diameter has to slightly exceed the inner tube diameter. Bristles that are too large will bend too much and not clean effectively.

#### BRUSH SELECTION

Simply find the section in our catalogue that best corresponds to the overall length and diameter of your instrument. As a rule, if your present brush is 8" (20.32cm) in length and we only offer 10" (25.4cm) length with your required diameter, that is entirely acceptable.


In many cases, our ARC™ Wire Brushes will have more bristles on the brush. As an example, many 5mm (0.19") diameter brushes are only 1" (2.5cm) in length. At BATRIK, we strived to provide the greatest amount of bristles to ensure the most effective cleaning.

#### CUSTOMIZABLE WIRE BRUSHES


A unique aspect of single-use ARC™ Wire Brushes is that they can be customized to any length:

- Up to 2000mm for Wire Brushes on plastic semi-rigid sheaths, rigid brushes from 5 to 20mm in diameter.
- No limit in length for Wire Brushes on flexible plastic sheaths, flexible brushes from 2-5mm in diameter

## RIGID



## FLEXIBLE



## 12" / 30 CM SINGLE-END WIRE BRUSH - DELRIN PLASTIC

	ITEM # SINGLE-USE	ITEM # REUSABLE	COLOR	TOTAL LENGTH	BRUSH (Ø)	BRUSH LENGTH	SHEATH (Ø)
<b>FLEXIBLE</b> BRUSH & SHEATH	<b>41-31502-S</b>	<b>41-31502</b>	White	12.401" / 31.5 cm	0.078" / 2 mm	0.59" / 15 mm	0.39" / 1 mm
	<b>41-32503-S</b>	<b>41-32503</b>	Yellow	12.795" / 32.5 cm	0.118" / 3 mm	0.98" / 25 mm	0.39" / 1 mm
	<b>411-32504-S</b>	<b>411-32504</b>	Violet	12.795" / 32.5 cm	0.157" / 4 mm	0.98" / 25 mm	0.61" / 1.55 mm
	<b>411-32505-S</b>	<b>411-32505</b>	Grey	12.795" / 32.5 cm	0.19" / 5 mm	0.98" / 25 mm	0.68" / 1.75 mm
<b>RIGID</b> BRUSH & SHEATH	<b>41-35005-S</b>	<b>41-35005</b>	White	13.779" / 35 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
	<b>41-35006-S</b>	<b>41-35006</b>	Green	13.779" / 35 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
	<b>41-35008-S</b>	<b>41-35008</b>	Blue	13.779" / 35 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
	<b>41-35010-S</b>	<b>41-35010</b>	Red	13.779" / 35 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
	<b>41-35012-S</b>	<b>41-35012</b>	Yellow	13.779" / 35 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
	<b>41-35015-S</b>	<b>41-35015</b>	Green	13.779" / 35 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
	<b>41-35020-S</b>	<b>41-35020</b>	Blue	13.779" / 35 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING Single-Use: Bag of 20 / Reusable: Bag of 3

## 18" / 45 CM SINGLE-END WIRE BRUSH - DELRIN PLASTIC

	ITEM # SINGLE-USE	ITEM # REUSABLE	COLOR	TOTAL LENGTH	BRUSH (Ø)	BRUSH LENGTH	SHEATH (Ø)
<b>FLEXIBLE</b> BRUSH & SHEATH	<b>41-46502-S</b>	<b>41-46502</b>	White	18.307" / 46.5 cm	0.078" / 2 mm	0.59" / 15 mm	0.39" / 1 mm
	<b>41-47503-S</b>	<b>41-47503</b>	Yellow	18.7" / 47.5 cm	0.118" / 3 mm	0.98" / 25 mm	0.39" / 1 mm
	<b>41-47505-S</b>	<b>N/A</b>	White	18.7" / 47.5 cm	0.19" / 5 mm	0.98" / 25 mm	0.68" / 1.75 mm
	<b>411-47504-S</b>	<b>411-47504</b>	Violet	18.7" / 47.5 cm	0.157" / 4 mm	0.98" / 25 mm	0.61" / 1.55 mm
	<b>411-47505-S</b>	<b>411-47505</b>	Grey	18.7" / 47.5 cm	0.19" / 5 mm	0.98" / 25 mm	0.68" / 1.75 mm
<b>RIGID</b> BRUSH & SHEATH	<b>41-50005-S</b>	<b>41-50005</b>	White	19.685" / 50 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
	<b>41-50006-S</b>	<b>41-50006</b>	Green	19.685" / 50 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
	<b>41-50008-S</b>	<b>41-50008</b>	Blue	19.685" / 50 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
	<b>41-50010-S</b>	<b>41-50010</b>	Red	19.685" / 50 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
	<b>41-50012-S</b>	<b>41-50012</b>	Yellow	19.685" / 50 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
	<b>41-50015-S</b>	<b>41-50015</b>	Green	19.685" / 50 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
	<b>41-50020-S</b>	<b>41-50020</b>	Blue	19.685" / 50 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Single-Use: Bag of 20 / Reusable: Bag of 3

**24" / 60 CM SINGLE-END WIRE BRUSH - DELRIN PLASTIC**

	ITEM # SINGLE-USE	ITEM # REUSABLE	COLOR	TOTAL LENGTH	BRUSH (Ø)	BRUSH LENGTH	SHEATH (Ø)
<b>FLEXIBLE</b> BRUSH & SHEATH	<b>41-61502-S</b>	<b>41-61502</b>	White	24.212" / 61.5 cm	0.078" / 2 mm	0.59" / 15 mm	0.39" / 1mm
	<b>41-62503-S</b>	<b>41-62503</b>	Yellow	24.6" / 62.5 cm	0.118" / 3 mm	0.98" / 25 mm	0.39" / 1mm
	<b>41-62505-S</b>	<b>N/A</b>	White	24.6" / 62.5 cm	0.19" / 5 mm	0.98" / 25 mm	0.68" / 1.75mm
	<b>411-62504-S</b>	<b>411-62504</b>	Violet	24.6" / 62.5 cm	0.157" / 4 mm	0.98" / 25 mm	0.61" / 1.55mm
	<b>411-62505-S</b>	<b>411-62505</b>	Grey	24.6" / 62.5 cm	0.19" / 5 mm	0.98" / 25 mm	0.68" / 1.75mm
<b>RIGID</b> BRUSH & SHEATH	<b>41-65005-S</b>	<b>41-65005</b>	White	25.59" / 65 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3mm
	<b>41-65006-S</b>	<b>41-65006</b>	Green	25.59" / 65 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3mm
	<b>41-65008-S</b>	<b>41-65008</b>	Blue	25.59" / 65 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
	<b>41-65010-S</b>	<b>41-65010</b>	Red	25.59" / 65 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
	<b>41-65012-S</b>	<b>41-65012</b>	Yellow	25.59" / 65 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
	<b>41-65015-S</b>	<b>41-65015</b>	Green	25.59" / 65 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
	<b>41-65020-S</b>	<b>41-65020</b>	Blue	25.59" / 65 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm



PACKAGING: Single-Use: Bag of 20 / Reusable: Bag of 3

**30" / 75 CM SINGLE-END WIRE BRUSH - DELRIN PLASTIC**

	ITEM # SINGLE-USE	ITEM # REUSABLE	COLOR	TOTAL LENGTH	BRUSH (Ø)	BRUSH LENGTH	SHEATH (Ø)
<b>FLEXIBLE</b> BRUSH & SHEATH	<b>41-76502-S</b>	<b>41-76502</b>	White	30.118" / 76.5 cm	0.078" / 2 mm	0.59" / 15 mm	0.39" / 1 mm
	<b>41-77503-S</b>	<b>41-77503</b>	Yellow	30.511" / 77.5 cm	0.118" / 3 mm	0.98" / 25 mm	0.39" / 1 mm
	<b>411-77504-S</b>	<b>411-77504</b>	Violet	30.511" / 77.5 cm	0.157" / 4 mm	0.98" / 25 mm	0.61" / 1.55 mm
	<b>411-77505-S</b>	<b>411-77505</b>	Grey	30.511" / 77.5 cm	0.19" / 5 mm	0.98" / 25 mm	0.68" / 1.75 mm
<b>RIGID</b> BRUSH & SHEATH	<b>41-80005-S</b>	<b>41-80005</b>	White	31.496" / 80 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
	<b>41-80006-S</b>	<b>41-80006</b>	Green	31.496" / 80 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
	<b>41-80008-S</b>	<b>41-80008</b>	Blue	31.496" / 80 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
	<b>41-80010-S</b>	<b>41-80010</b>	Red	31.496" / 80 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
	<b>41-80012-S</b>	<b>41-80012</b>	Yellow	31.496" / 80 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
	<b>41-80015-S</b>	<b>41-80015</b>	Green	31.496" / 80 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
	<b>41-80020-S</b>	<b>41-80020</b>	Blue	31.496" / 80 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm


PACKAGING: Single-Use: Bag of 20 / Reusable: Bag of 3

## BRUSHES



## FLEXIBLE

Flexible brush diameter (Ø) ranging from 2 mm to 5 mm



## RIGID

Rigid brush diameter (Ø) ranging from 5 mm to 20 mm

### BRUSH SELECTION SYSTEM (Actual Size)

2mm 0.078"	3 mm 0.118"	4 mm 0.157"	5 mm 0.196"	6 mm 0.236"	8 mm 0.314"	10 mm 0.394"	12 mm 0.472"	15 mm 0.591"	20 mm 0.787"

# BATRIK REUSABLE DOUBLE ARC™ WIRE BRUSHES

## ASYMMETRIC

### 14" / 35 CM ASYMMETRIC DOUBLE-END WIRE BRUSHES

ITEM #	COLOR	TOTAL LENGTH	BRUSH 1 (Ø)	BRUSH 2 (Ø)	BRUSH LENGTH	SHEATH (Ø)
40-35056	White	13.779" / 35 cm	0.19" / 5 mm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
40-35112	Yellow	13.779" / 35 cm	0.39" / 10 mm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
40-35152	Green	13.779" / 35 cm	0.59" / 15 mm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm
40-35810	Blue	13.779" / 35 cm	0.31" / 8 mm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm

PACKAGING: Bag of 3

### 20" / 50 CM ASYMMETRIC DOUBLE-END WIRE BRUSHES

ITEM #	COLOR	TOTAL LENGTH	BRUSH 1 (Ø)	BRUSH 2 (Ø)	BRUSH LENGTH	SHEATH (Ø)
40-50056	White	19.685" / 50 cm	0.19" / 5 mm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
40-50112	Yellow	19.685" / 50 cm	0.39" / 10 mm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
40-50152	Green	19.685" / 50 cm	0.59" / 15 mm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm
40-50810	Blue	19.685" / 50 cm	0.31" / 8 mm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm

PACKAGING: Bag of 3

### 26" / 65 CM ASYMMETRIC DOUBLE-END WIRE BRUSHES

ITEM #	COLOR	TOTAL LENGTH	BRUSH 1 (Ø)	BRUSH 2 (Ø)	BRUSH LENGTH	SHEATH (Ø)
40-65056	White	25.59" / 65 cm	0.19" / 5 mm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
40-65112	Yellow	25.59" / 65 cm	0.39" / 10 mm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
40-65152	Green	25.59" / 65 cm	0.59" / 15 mm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm
40-65810	Blue	25.59" / 65 cm	0.31" / 8 mm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm

PACKAGING: Bag of 3

### 31" / 80 CM ASYMMETRIC DOUBLE-END WIRE BRUSHES

ITEM #	COLOR	TOTAL LENGTH	BRUSH 1 (Ø)	BRUSH 2 (Ø)	BRUSH LENGTH	SHEATH (Ø)
40-80056	White	31.496" / 80 cm	0.19" / 5 mm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
40-80112	Yellow	31.496" / 80 cm	0.39" / 10 mm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
40-80152	Green	31.496" / 80 cm	0.59" / 15 mm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm
40-80810	Blue	31.496" / 80 cm	0.31" / 8 mm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm

PACKAGING: Bag of 3

**BATRIK**

# REUSABLE DOUBLE ARC™ WIRE BRUSHES

## SYMMETRIC

### 14" / 35 CM SYMMETRIC DOUBLE-END WIRE BRUSHES

ITEM #	COLOR	TOTAL LENGTH	BRUSH 1 (Ø)	BRUSH LENGTH	SHEATH (Ø)
<b>40-35005</b>	White	13.779" / 35 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
<b>40-35006</b>	Green	13.779" / 35 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
<b>40-35008</b>	Blue	13.779" / 35 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
<b>40-35010</b>	Red	13.779" / 35 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
<b>40-35012</b>	Yellow	13.779" / 35 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
<b>40-35015</b>	Green	13.779" / 35 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
<b>40-35020</b>	Blue	13.779" / 35 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Bag of 3

### 20" / 50 CM SYMMETRIC DOUBLE-END WIRE BRUSHES

ITEM #	COLOR	TOTAL LENGTH	BRUSH 1 (Ø)	BRUSH LENGTH	SHEATH (Ø)
<b>40-50005</b>	White	19.685" / 50 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
<b>40-50006</b>	Green	19.685" / 50 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
<b>40-50008</b>	Blue	19.685" / 50 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
<b>40-50010</b>	Red	19.685" / 50 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
<b>40-50012</b>	Yellow	19.685" / 50 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
<b>40-50015</b>	Green	19.685" / 50 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
<b>40-50020</b>	Blue	19.685" / 50 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Bag of 3

### 26" / 65 CM SYMMETRIC DOUBLE-END WIRE BRUSHES

ITEM #	COLOR	TOTAL LENGTH	BRUSH 1 (Ø)	BRUSH LENGTH	SHEATH (Ø)
<b>40-65005</b>	White	25.59" / 65 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
<b>40-65006</b>	Green	25.59" / 65 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
<b>40-65008</b>	Blue	25.59" / 65 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
<b>40-65010</b>	Red	25.59" / 65 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
<b>40-65012</b>	Yellow	25.59" / 65 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
<b>40-65015</b>	Green	25.59" / 65 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
<b>40-65020</b>	Blue	25.59" / 65 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Bag of 3

### 31" / 80 CM SYMMETRIC DOUBLE-END WIRE BRUSHES

ITEM #	COLOR	TOTAL LENGTH	BRUSH 1 (Ø)	BRUSH LENGTH	SHEATH (Ø)
<b>40-80005</b>	White	31.496" / 80 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
<b>40-80006</b>	Green	31.496" / 80 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
<b>40-80008</b>	Blue	31.496" / 80 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
<b>40-80010</b>	Red	31.496" / 80 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
<b>40-80012</b>	Yellow	31.496" / 80 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
<b>40-80015</b>	Green	31.496" / 80 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
<b>40-80020</b>	Blue	31.496" / 80 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Bag of 3

## ASYMMETRIC BRUSHES

**5 MM**

**6 MM**



**8 MM**

**10 MM**



**10 MM**

**12.5 MM**



**15 MM**

**20 MM**




BATRIK

# REUSABLE IN-BRUSH™ ANTIMICROBIAL TWISTED WIRE BRUSHES

Our In-Brush™ Antimicrobial Reusable Twisted Wire Brushes are designed for use with Frazier/Baron suction tubes. The twisted wires have a fine/medium diameter yet are rigid enough to ensure excellent cleaning capabilities.

Our wide assortment of Nylon Twisted Wire Brushes are ideal for large bore, uniquely shaped, and endoscopic instruments.

GOOD TO KNOW: OUR TWISTED WIRE BRUSHES FEATURE A PROTECTIVE BLUE TIP FOR ADDED SAFETY.



## 10" / 250 MM FRAZIER /BARON SUCTION TUBE BRUSHES

ITEM #	SIZE	BRUSH LENGTH	BRUSH (Ø)
<b>21-10039</b>	3 FR	0.5" / 12.7 mm	0.039" / 1 mm
<b>21-10067</b>	5 FR	0.5" / 12.7 mm	0.065" / 1.67 mm
<b>20-10059</b>	5 FR	1" / 25.4 mm	0.059" / 1.50 mm
<b>20-10079</b>	6 FR	1" / 25.4 mm	0.079" / 2 mm
<b>20-10090</b>	7 FR	1" / 25.4 mm	0.090" / 2.3 mm
<b>20-10106</b>	8 FR	1" / 25.4 mm	0.106" / 2.7 mm
<b>20-10118</b>	9 FR	1" / 25.4 mm	0.118" / 3 mm
<b>20-10129</b>	10 FR	1.5" / 38 mm	0.129" / 3.3 mm
<b>20-10146</b>	11 FR	1.5" / 38 mm	0.0146" / 3.7 mm
<b>20-10157</b>	12 FR	1.5" / 38 mm	0.157" / 4 mm

PACKAGING: Bag of 3.  
Also available in bags of 50, add the suffix "-50".

## 12" / 305 MM NYLON TWISTED WIRE BRUSHES

ITEM #	BRUSH LENGTH	BRUSH (Ø)
<b>22-12078</b>	2" / 50.8 mm	0.078" / 2 mm
<b>22-12118</b>	2" / 50.8 mm	0.118" / 3 mm
<b>22-12157</b>	3" / 76 mm	0.157" / 4 mm
<b>22-12197</b>	3" / 76 mm	0.197" / 5 mm
<b>22-12236</b>	3" / 76 mm	0.236" / 6 mm
<b>22-12393</b>	3" / 76 mm	0.394" / 10 mm

PACKAGING: Bag of 3.  
Also available in bags of 50, add the suffix "- 50".

### 16" / 406 MM NYLON TWISTED WIRE BRUSHES

ITEM #	BRUSH LENGTH	BRUSH (Ø)
<b>23-16197</b>	2.5" / 63.5 mm	0.197" / 5 mm
<b>23-16236</b>	2.5" / 63.5 mm	0.236" / 6 mm
<b>23-16276</b>	3" / 76 mm	0.276" / 7 mm
<b>23-16394</b>	3" / 76 mm	0.394" / 10 mm
<b>23-16472</b>	3" / 76 mm	0.472" / 12 mm
<b>23-16591</b>	3" / 76 mm	0.591" / 15 mm
<b>23-16787</b>	3" / 76 mm	0.787" / 20 mm

PACKAGING: Bag of 3.

Also available in bags of 50, add the suffix "-50".

### 18" / 457 MM NYLON TWISTED WIRE BRUSHES

ITEM #	BRUSH LENGTH	BRUSH (Ø)
<b>19-18010</b>	0.47" / 12 mm	0.039" / 1 mm
<b>19-18015</b>	0.5" / 12 mm	0.06" / 1.6 mm

PACKAGING: Bag of 3.


Also available in bags of 50, add the suffix "-50".

### 24" / 610 MM NYLON TWISTED WIRE BRUSHES

ITEM #	BRUSH LENGTH	BRUSH (Ø)
<b>24-24098</b>	2.5" / 63.5 mm	0.098" / 2.5 mm
<b>24-24197</b>	2.5" / 63.5 mm	0.197" / 5 mm
<b>24-24276</b>	3" / 76 mm	0.276" / 7 mm
<b>24-24394</b>	3" / 76 mm	0.394" / 10 mm
<b>24-24472</b>	3" / 76 mm	0.472" / 12 mm

PACKAGING: Bag of 3.

Also available in bags of 50, add the suffix "-50".



### BRONCHO BRUSHES

ITEM #	DESCRIPTION	BRUSH LENGTH	BRUSH (Ø) INCHES / MM
<b>26-06175</b>	Small	6" / 152 mm	1.75" / 44.5 mm
<b>26-08325</b>	Large	8" / 203 mm	3.25" / 82.6 mm

PACKAGING: Bag of 3.

Also available in bags of 50, add the suffix "-50".

### TRACHEOSTOMY BRUSHES

ITEM #	DESCRIPTION	TOTAL LENGTH	BRUSH LENGTH	BRUSH (Ø)
<b>25-08500</b>	Large 6-12.3 3 / 4"	8" / 203 mm	3.75" / 95 mm	0.5" / 12.7 mm
<b>25-37550</b>	Small 000-5.1 7 / 8"	4" / 93 mm	1.88" / 48 mm	0.5" / 12.7 mm

PACKAGING: Bag of 3. Also available in bags of 50, add the suffix "-50".

# BATRIK STAINLESS STEEL ARC™ MICRO BRUSHES

Our Stainless Steel ARC™ Micro Brushes are some of the smallest brushes in the industry. They are designed to clean channels with an inner diameter ranging from 0.5 mm to 1.40 mm.

- Flexible stainless steel sheath will not bend like the regular metal twisted wire brushes
- Polyamide bristles prevents or slows bacterial growth between sterilizations
- Reusable\*
- Available in two lengths

TIP: USING THE PULL-THROUGH  
TECHNIQUE IS HIGHLY  
RECOMMENDED WITH THESE  
BRUSHES. SIMPLY INSERT THE  
BLUNT, NON-BRUSH END FIRST.



## PRODUCT SELECTION

ITEM #	TOTAL LENGTH	BRUSH SIZE (Ø)	BRUSH LENGTH	SHEATH (Ø)	PACKAGING
<b>35-03015</b>	12" / 30 cm	0.06" / 1.5 mm	0.6" / 15 mm	0.02" / 0.40 mm	6 / Tube
<b>35-04515</b>	18" / 45 cm	0.06" / 1.5 mm	0.6" / 15 mm	0.02" / 0.40 mm	6 / Tube

\*Although all materials used in the fabrication of these brushes can be re-processed repeatedly, due to the delicate nature of these brushes re-use is conditional on their handling and application.

# BATRIK

# INSTRUMENT CLEANING

# FLAT BRUSHES

Our Instrument Cleaning Flat Brushes are specifically designed for hospital use in the pre-clean and preparation of instruments.

- Long ergonomic handle permits ease-of-use
- Extra long length of brush head (60 [±5] mm in Single-Use and 75 [±5] mm in Reusable)
- Transparent brush fibers allow for quick verification of its cleanliness
- The plastic handle will not release any toxic components
- Available in 6 models (4 reusable/autoclavable and 2 single-use)



## FLEXIBLE

This brush is designed for large instruments that require a more gentle cleaning application such as optics, plastics, or hard to reach spot.



## EXTRA FLEXIBLE

The very long, flexible bristles eliminate any stress during sanitation, and makes this the ideal brush for fragile instruments as well as both pliable and rigid endoscopes. The length of the bristles assist in cleaning areas that would otherwise be difficult to reach, such as the accessory channel or tip deflection of pliable endoscopes.



## RIGID

The rigidity of this brush, along with the length of the bristles and handle, improves and simplifies instrument cleaning.



## EXTRA RIGID

Instrument safe polyamide fibers offer steel-like rigidity and allows deep cleaning of soiled instruments, particularly in grooved or uneven surfaces.



50-40919



50-40921



50-40920



50-40918



50-40922



50-40917

## PRODUCT SELECTION

### DELTRIN PLASTIC

ITEM #	BRUSH TYPE	DESCRIPTION	COLOR	BRUSH LENGTH	BRUSH SIZE	PACKAGING
50-40917	Extra Rigid	Reusable Autoclavable Flat Brush	Blue	220 mm	75 mm x 20 mm x 10 mm	Bag of 3
50-40918	Rigid	Reusable Autoclavable Flat Brush	Green	220 mm	75 mm x 20 mm x 15 mm	Bag of 3
50-40919	Flexible	Reusable Autoclavable Flat Brush	Yellow	220 mm	75 mm x 20 mm x 15 mm	Bag of 3
50-40920	Extra Flexible	Reusable Autoclavable Flat Brush	Red	220 mm	75 mm x 20 mm x 25 mm	Bag of 3
50-40921	Flexible	Single-Use Flat Brush	Orange	220 mm	60 mm x 20 mm x 15 mm	Bag of 20
50-40922	Rigid	Single-Use Flat Brush	Light Blue	220 mm	60 mm x 20 mm x 15 mm	Bag of 20

Reusable models are also available in bags of 50. Add the suffix "- 50".

BATRIK

# REUSABLE PRECISION BRUSH

## INSTRUMENT CLEANING

Our Reusable Precision Instrument Cleaning Brush is the ideal brush for the quick removal of difficult debris. It cleans externally as well as deep into instrument jaws. Although the firm and durable bristles feel like stainless steel, there is absolutely no risk of instrument damage or personal injury.

- Bristles have an extra firm density and superior durability
- Polyamide (TYNEX®) bristles
- Steamable thermoplastic ergonomic handle
- Steam / chemical sterilizable
- Latex free



GOOD TO KNOW: ELIMINATES  
THE NEED FOR THE STAINLESS  
STEEL TOOTHBRUSH-STYLE BRUSH,  
WHICH CAN DAMAGE VALUABLE  
INSTRUMENTS!

## PRODUCT SELECTION

### BLUE DELRIN PLASTIC

ITEM #	DESCRIPTION	TOTAL LENGTH
50-40800	Reusable Precision External Cleaning Brush	6.29" / 160 mm

PACKAGING: Bag of 3



NEW!


# TRIFECTA™ BRUSH

## INSTRUMENT CLEANING

Introducing the **Trifecta™ Brush**, a versatile **3-in-1 tool** designed to meet all your instrument cleaning needs. This reusable instrument brush is available in two styles: one with all-nylon bristles for gentle cleaning and another with nylon bristles on two sides and **stainless steel bristles** on the third, ideal for tackling tougher buildup and stubborn cleaning tasks. This innovative design ensures effective, thorough maintenance for a wide range of instruments.

### TRIFECTA™ BRUSH - NYLON

- 1 NYLON BRISTLES**  
Gentle, transparent nylon bristles provide an effective yet safe cleaning solution. This side allows you to verify cleanliness easily, without risking damage to delicate surfaces. Ideal for a light, non-invasive clean.
- 2 PRECISION TIP**  
Reach those hard-to-clean spots with the precision tip. This small, firm-bristle tip is designed for tight, hard-to-reach areas, such as instrument jaws and box joints. It ensures a deep clean where it's most needed, without the risk of damage or injury to delicate parts.
- 3 FIRM NYLON BRISTLES**  
For tougher residues, the firm nylon side delivers the cleaning power you need. With bristles that feel similar to stainless steel but are safe for your instruments, this side easily handles stubborn dirt and buildup.



### TRIFECTA™ BRUSH - NYLON & STAINLESS STEEL



29-23012

28-23012

### PRODUCT SELECTION

ITEM #	COLOR	DESCRIPTION	LENGTH TOTAL	BRISTLE LENGTH	BRISTLES 1 Ø	BRISTLES 2 Ø	BRISTLES 3 Ø
28-23012	Blue	Triple Nylon Brush	8.66" / 220 mm	1.20" x 1.20" x 0.47" / 30.5 x 30.5 x 12.0 mm	0.47" / 12.0 mm	0.47" / 12.0 mm	0.315" / 8.0 mm
29-23012	Blue	Triple Nylon / Stainless Steel	8.66" / 220 mm	1.20" x 1.20" x 0.47" / 30.5 x 30.5 x 12.0 mm	0.47" / 12.0 mm	0.47" / 12.0 mm	0.47" / 12.0 mm

PACKAGING: Bag of 3 - For a bag of 50, add "-50" following the item code.

BATRIK

# ALL PURPOSE REUSABLE BRUSHES

## Plastic Handle Brushes

We offer a wide selection of brushes for instrument cleaning applications. These brushes are ideal for large bore, uniquely shaped, and endoscopic instruments.

### In-Brush™ Plastic Handle

Available in regular toothbrush style (nylon or stainless) and double-ended brushes.



28-72550



28-69394



29-72550



29-69394



29-70050



50-40911

## PRODUCT SELECTION


ITEM #	DESCRIPTION	LENGTH TOTAL	BRUSH LENGTH	BRUSH HEIGHT
<b>28-72550</b>	Plastic Single End Nylon Reusable Handle Brush	7.25" / 184 mm	1.5" / 38 mm	0.5" / 12.7 mm
<b>29-72550</b>	Plastic Single End Stainless Steel Reusable Handle Brush	7.25" / 184 mm	1.5" / 38 mm	0.5" / 12.7 mm
<b>29-70050</b>	Plastic Single End Brass Reusable Handle Brush	7.25" / 184 mm	1.5" / 38 mm	0.5" / 12.7 mm
<b>28-69394</b>	Plastic Double End Nylon Reusable Handle Brush	6.9" / 175 mm	1.38" / 35 mm 1" / 25 mm	0.39" / 10 mm 0.2" / 5 mm
<b>29-69394</b>	Plastic Double End Stainless Steel Reusable Handle Brush	6.9" / 175 mm	1.38" / 35 mm 1" / 25 mm	0.39" / 10 mm 0.2" / 5 mm
<b>50-40911</b>	Hand Scrub Nylon Brush Reusable autoclavable	4.09" / 104 mm	3.74" / 95 mm	1.1" / 28 mm

**PACKAGING:** Bag of 3 - For a bag of 50, add "-50" following the item code.

**BATRIK**

# ACETABULAR REAMER BRUSH

Our double-ended cleaning brush is manufactured exclusively for cleaning acetabular bone reamers. The "Pig Tail" rounded style circular head brush is designed to effectively clean inside the reamer, while the rigid straight brush on the reverse side is used for cleaning individual cutting holes. Special fusion of the brush heads to the hand piece along with its high quality antimicrobial medical grade bristles ensures a durable, reusable, and highly effective cleaning tool.

**STRAIGHT BRUSH****ROUND BRUSH**

## PRODUCT SELECTION

### GREEN DELRIN PLASTIC

ITEM #	DESCRIPTION	HANDLE LENGTH	STRAIGHT BRUSH	ROUND BRUSH HEAD
<b>50-ACETAB</b>	Acetabular Reamer Cleaning Brush	4.72" / 120 mm ø 0.47" / 12 mm	1.18" / 30 mm ø 0.3125" / 8 mm	2" / 50 mm ø 1.75" / 45 mm

PACKAGING: Bag of 3


NOTE: Specifications subject to small variances

BATRIK

# REUSABLE ORTHOPEDIC RASPS CLEANING BRUSH

Our Reusable Orthopedic Rasps Cleaning Brush is used to clean orthopedic rasps and other similar instruments used in orthopedic surgery.

- Long lasting and can withstand many sterilization cycles
- Ergonomic steamable thermoplastic handle
- Antibacterial polyamide bristles
- Reusable
- Latex free



## PRODUCT DESCRIPTION


ITEM #	TOTAL LENGTH	HANDLE LENGTH	BRUSH SIZE	PACKAGING
50-RASP	9" / 23 cm	5.5" / 14 cm	0.78" x 2.7" x 3.5" (2 cm x 7 cm x 9 cm)	3 / Bag

BATRIK

# BOTTLE & DRAIN CLEANING BRUSHES

Designed specifically for hospital applications, our Bottle and Drain Cleaning Brushes feature long and flexible bristles to maximize the effectiveness of large-surface cleaning.

- High quality polyamide clear bristles
- Easy to identify debris
- Autoclavable
- Latex free



## PRODUCT SELECTION

### BLUE DELRIN PLASTIC

ITEM #	DESCRIPTION	HANDLE LENGTH	BRUSH LENGTH	BRUSH (Ø)	TOTAL LENGTH
50-10295	Reusable Jar / Bottle Brush	120 mm ± 2 mm	120 mm ± 5 mm	7 mm ± 2 mm	13.19" / 355 mm
50-10224	Sterilizer Drain Cleaning Brush	N/A	4" / 10.1 cm	2" / 5.08 cm	24" / 609 mm

PACKAGING: Bag of 3

BATRIK

# AUTOCLAVE CLEANING BRUSH KITS



## ROUND BRUSH

- Features 3 holes for adjustable pole placement
- Adjustable brush placement (parallel or perpendicular) and longer brush length provides easy access to hard-to-reach corners
- Made with FDA approved polyester bristles

## HAND BRUSH

- Ergonomic handle for better hand comfort
- Longer, FDA approved polyester bristles
- Easy to clean (no cavity)

## PRODUCT SELECTION

### REGULAR AUTOCLAVE CLEANING BRUSH KIT

ITEM # (KIT)	DESCRIPTION	ITEM # (UNIT)	DESCRIPTION	DIMENSIONS
50-CLAVEREG	Autoclave REGULAR Brush Cleaning Kit	50-35900	1 Hand Brush	Length: 3.5" x 8" / 8.89 cm x 20.32 cm
		50-55800	1 Rounded Corner Brush (Head Brush)	Head size: 5.5 x 10" / 14 x 25.4 cm
		50-EXTPOLE	1 Anodized Aluminum Extension Pole	2-4 ft / 60-122 cm

### SUPERIOR AUTOCLAVE CLEANING BRUSH KIT

ITEM # (KIT)	DESCRIPTION	ITEM # (UNIT)	DESCRIPTION	DIMENSIONS	CASE DIMENSIONS
50-CLAVESUP	Autoclave SUPERIOR Brush Cleaning Kit	50-35900	1 Hand Brush	Length: 3.5" x 8" / 11.5 cm x 20.32cm	
		50-55800	1 Rounded Corner Brush (Head Brush)	Head size: 5.5 x 10" / 14 x 25.4 cm	
		50-EXTPOLE	1 Anodized Aluminum Extension Pole	2-4 ft / 60-122 cm	
		SF-240	2 SAVEFACE™ Face Shields	9" / 23 cm	CASE OF 240 / 26.5" x 21" x 15"

## REPLACEMENT COMPONENTS



### HAND BRUSH

50-35900

Packaging: 1



### BRUSH HEAD

50-55800

Packaging: 1



### EXTENSION POLE

50-EXTPOLE

Packaging: 1



### SAVEFACE™ FACE SHIELDS SF-240

Packaging: 100 / Box

# BOWL STANDS & RECEPTACLES

Our Bowl Stands and Receptacles are designed for convenience, mobility, and durability. Available in five versatile models that accommodate most standard sizes of stainless steel bowls and pails. Fully welded stainless steel construction throughout with a #4 polished finish.

All bowl stands are backed by an outstanding **7 year warranty** on all parts (see warranty for details).

## KEY FEATURES

- Seamless welded tubular frame with lower cross-bracing provides maximum rigidity and extended product life
- Optional lower stainless steel shelf is available for all configurations
- Non-marking conductive swivel design 2" (51 mm) casters

## PRODUCT SELECTION

ITEM #	DESCRIPTION
<b>OR-36310</b>	<b>SPONGE RECEPTACLES: FLOOR UNITS</b> Mobile Kick Bucket with Qt. Pail and Perimeter Bumper Receptacle: 14.5" Ø x 8" H / 368 mm Ø x 203 mm H Receptacle and bucket: 14.5" Ø x 15" H / 368 mm Ø x 381 mm H Bucket: 12" Ø x 11" H / 305 mm Ø x 279 mm H
<b>OR-3632</b>	<b>MOBILE SPONGE RECEPTACLES</b> with 10 Qt. Oval Basin and Perimeter Bumper
<b>OR-3631-TB</b>	<b>SPONGE RECEPTACLES: TABLE HEIGHT</b> Mobile Table Height Receptacle with 20 Qt. Pail 15" L x 16" W x 38" H (381 mm x 406 mm x 965 mm)
<b>OR-3651</b>	<b>SPONGE RECEPTACLES: TABLE HEIGHT</b> Single Bowl Stand with 7 Qt. Basin 15" L x 16" W x 34" H (381 mm x 406 mm x 864 mm)
<b>OR-3652</b>	<b>SPONGE RECEPTACLES: TABLE HEIGHT</b> Double Bowl Stand with Two 7 Qt. Basins 33" L x 16" W x 34" H (838 mm x 406 mm x 864 mm)



THE LOWER SHELF IS AN OPTIONAL FEATURE THAT MUST BE ADDED TO THE ITEM IF DESIRED

## ACCESSORIES

ITEM #	DESCRIPTION
<b>OR-3631-PS</b>	12 QT. STAINLESS STEEL PAIL for OR-36310
<b>OR-3631-TB-PS</b>	20 QT. STAINLESS STEEL PAIL for OR-3631-TB
<b>OR-3632-BS</b>	10 QT. OVAL STAINLESS STEEL BASIN for OR 3632
<b>OR-3651-BS</b>	7 QT. STAINLESS STEEL BASIN for OR-3651 and OR-3652

Add suffix -LS to any Item # to add a lower shelf

# GRADUATED MEASURING

## PRODUCT SELECTION

ITEM #	DESCRIPTION (Ø x H)
H-0201	300 CC (10 OZ) 83 mm x 114 mm
H-0202	500CC (16 OZ) 96 mm x 124 mm
H-0203	1000 CC (32 OZ) 128 mm x 134 mm
H-0204	2000 CC (64 OZ) 151 mm x 172 mm



H-0201

H-0203

# GALLIPOTS

## PRODUCT SELECTION

ITEM #	DESCRIPTION (Ø x H)
H-0361	GALLIPOT, 2 OZ 65 mm x 35 mm
H-0362	GALLIPOT, 4 OZ 77 mm x 45 mm
H-0363	GALLIPOT, 6 OZ 83 mm x 50 mm
H-0364	GALLIPOT, 10 OZ 95 mm x 55 mm
H-0365	GALLIPOT, 14 OZ 112 mm x 58 mm
H-0366	GALLIPOT, 8 OZ 90 mm x 50 mm
H-0367	GALLIPOT, 12 OZ 103 mm x 60 mm



H-0361

H-0363

All Imperial Ware products are backed by our 1 year warranty.

# WASH BASINS & BOWLS



H-0321

## PRODUCT SELECTION

ITEM #	DESCRIPTION (Ø x H)
H-1228	WASH BASIN, SHALLOW TYPE, 12.5" 320 mm x 80 mm
H-1229	WASH BASIN, DEEP TYPE, 12.5" 320 mm x 105 mm
H-1230	WASH BASIN, DEEP TYPE, 13.375" 340 mm x 100 mm
H-1231	WASH BASIN, DEEP TYPE, 14.4" 365 mm x 115 mm
H-0301	BOWL, 10 CM, 3/8 QT, 14 OZ 108 mm x 50 mm
H-0302	BOWL, 12 CM, 5/8 QT, 20 OZ 133 mm x 56 mm
H-0303	BOWL, 14 CM, 1 QT, 32 OZ 145 mm x 60 mm
H-0304	BOWL, 16 CM, 1 3/8 QT, 44 OZ 167 mm x 70 mm
H-0305	BOWL, 18 CM, 1 7/8 QT, 60 OZ 185 mm x 70 mm
H-0321	SOLUTION BOWL, 4 QT, 128 OZ 330 mm x 84 mm
H-0322	SOLUTION BOWL, 5 1/8 QT, 164 OZ 308 mm x 110 mm
H-0323	SOLUTION BOWL, 7 QT, 224 OZ Ø 355 mm x 115 mm



H-0302

All Imperial Ware products are backed by our 1 year warranty.

# FOOD SERVICE CONTAINERS



**H-0802**

## PRODUCT SELECTION

ITEM #	DESCRIPTION (L x W x H)
<b>H-0801</b>	FOOD SERVICE CONTAINER, FULL SIZE 530 mm x 320 mm x 65 mm
<b>H-0801P</b>	PERFORATED FOOD SERVICE CONTAINER, FULL SIZE 530 mm x 320 mm x 65 mm
<b>H-0802</b>	FOOD SERVICE CONTAINER, FULL SIZE 530 mm x 320 mm x 100mm
<b>H-0802P</b>	PERFORATED FOOD SERVICE CONTAINER, FULL SIZE 530 mm x 320 mm x 100mm
<b>H-0803</b>	COVER FOR FOOD SERVICE CONTAINER, FULL SIZE 533 mm x 325 mm x 35 mm
<b>H-0807</b>	FOOD SERVICE CONTAINER, 1/2 SIZE 321 mm x 262 mm x 65 mm
<b>H-0807P</b>	PERFORATED FOOD SERVICE CONTAINER, 1/2 SIZE 321 mm x 262 mm x 65 mm
<b>H-0808</b>	FOOD SERVICE CONTAINER, 1/2 SIZE 321 mm x 262 mm x 100 mm
<b>H-0808P</b>	PERFORATED FOOD SERVICE CONTAINER, 1/2 SIZE 321 mm x 262 mm x 100 mm
<b>H-0809</b>	COVER FOR FOOD SERVICE CONTAINER, 1/2 SIZE 324 mm x 265 mm x 32 mm



**H-0801P**




**H-0808P**

All Imperial Ware products are backed by our **1 year warranty**.

# INSTRUMENT TRAYS

## shallow



**H-0557**



**H-0553**



**H-0558**

## PRODUCT SELECTION

ITEM #	DESCRIPTION (L x W x H)
<b>H-0553</b>	INSTRUMENT TRAY 210 mm x 110 mm x 8 mm
<b>H-0554</b>	SHALLOW RECTANGULAR TRAY 217 mm x 128 mm x 12 mm
<b>H-0556</b>	SHALLOW RECTANGULAR TRAY 225 mm x 162 mm x 14 mm
<b>H-0557</b>	SHALLOW RECTANGULAR TRAY 285 mm x 185 mm x 17 mm
<b>H-0558</b>	SHALLOW RECTANGULAR TRAY 338 mm x 245 mm x 17 mm
<b>H-0559</b>	SHALLOW RECTANGULAR TRAY 390 mm x 271 mm x 18 mm
<b>H-0559P</b>	PERFORATED SHALLOW TRAY 390 mm x 271 mm x 18 mm
<b>H-0560</b>	SHALLOW RECTANGULAR TRAY 440 mm x 295 mm x 18 mm
<b>H-0560P</b>	PERFORATED SHALLOW TRAY 440 mm x 295 mm x 18 mm
<b>H-0561</b>	SHALLOW RECTANGULAR TRAY 491 mm x 322 mm x 18 mm

All Imperial Ware products are backed by our **1 year warranty**.

# INSTRUMENT TRAYS

## deep & perforated



H-0563

## PRODUCT SELECTION

ITEM #	DESCRIPTION (L x W x H)
H-0562	DEEP RECTANGULAR TRAY 292 mm x 195 mm x 50 mm
H-0563	DEEP RECTANGULAR TRAY 335 mm x 246 mm x 55 mm
H-0563P	PERFORATED RECTANGULAR TRAY 335 mm x 246 mm x 55 mm
H-0564	DEEP RECTANGULAR TRAY 387 mm x 272 mm x 60 mm
H-0564P	PERFORATED DEEP RECTANGULAR TRAY 387 mm x 272 mm x 60 mm
H-0565	DEEP RECTANGULAR TRAY 445 mm x 281 mm x 97 mm
H-0565P	PERFORATED DEEP RECTANGULAR TRAY 445 mm x 281 mm x 97 mm
H-0566	DEEP RECTANGULAR TRAY 501 mm x 336 mm x 101 mm



H-0565



H-0562

All Imperial Ware products are backed by our 1 year warranty.

# CATHETER & INSTRUMENT TRAYS



**H-0153**




**H-0156**



**H-0159**

## PRODUCT SELECTION

ITEM #	DESCRIPTION
<b>H-0151</b>	CATHETER & INSTRUMENT TRAY (NO COVER)
<b>H-0152</b>	COVER FOR H-0151
<b>H-0153</b>	CATHETER & INSTRUMENT TRAY SET (WITH COVER, RAISED HANDLE) 219 mm L x 81 mm W x 35 mm H
<b>H-0154</b>	CATHETER & INSTRUMENT TRAY (NO COVER)
<b>H-0155</b>	COVER FOR H-0154
<b>H-0156</b>	CATHETER & INSTRUMENT TRAY SET (WITH COVER, RAISED HANDLE) 231 mm L x 132 mm W x 49 mm H
<b>H-0157</b>	CATHETER & INSTRUMENT TRAY (NO COVER)
<b>H-0158</b>	COVER FOR H-0157
<b>H-0159</b>	CATHETER & INSTRUMENT TRAY SET (WITH COVER, RAISED HANDLE) 320 mm L x 95 mm W x 50 mm H



All Imperial Ware products are backed by our **1 year warranty**.

# CATHETER & INSTRUMENT TRAYS

## PRODUCT SELECTION

ITEM #	DESCRIPTION
<b>H-0160</b>	CATHETER & INSTRUMENT TRAY (NO COVER)
<b>H-0160P</b>	PERFORATED CATHETER & INSTRUMENT TRAY (NO COVER)
<b>H-0161</b>	COVER FOR H-0160 AND H-0160P
<b>H-0162</b>	CATHETER & INSTRUMENT TRAY SET (WITH COVER, FLUSH HANDLE) 310 mm L x 196 mm W x 50 mm H
<b>H-0163</b>	CATHETER & INSTRUMENT TRAY (NO COVER)
<b>H-0163P</b>	PERFORATED CATHETER & INSTRUMENT TRAY (NO COVER)
<b>H-0164</b>	COVER FOR H-0163 AND H-0163P
<b>H-0165</b>	CATHETER & INSTRUMENT TRAY SET (WITH COVER, FLUSH HANDLE) 310 mm L x 196 mm W x 65 mm H
<b>H-0166</b>	CATHETER & INSTRUMENT TRAY (NO COVER)
<b>H-0167</b>	COVER FOR H-0166
<b>H-0168</b>	CATHETER & INSTRUMENT TRAY SET (WITH COVER, RAISED HANDLE) 441 mm L x 116 mm W x 45 mm H



**H-0162**



**H-0165**



**H-0168**

All Imperial Ware products are backed by our 1 year warranty.

# URINALS & BEDPANS



**H-0102**



**H-0108**



**H-0109**

## PRODUCT SELECTION

ITEM #	DESCRIPTION
<b>H-0101</b>	ADULT BEDPAN WITHOUT COVER 355 mm x 290 mm x 50 mm
<b>H-0102</b>	ADULT BEDPAN WITH COVER 355 mm x 165 mm
<b>H-0105</b>	FRACTURE BEDPAN (LONG COVER) 310 mm x 245 mm
<b>H-0106</b>	FRACTURE BEDPAN (SHORT COVER) 355 mm x 290 mm x 50
<b>H-0107</b>	MALE URINAL (VERTICAL TYPE) 110 mm x 230 mm
<b>H-0108</b>	MALE URINAL (HORIZONTAL TYPE) 110 mm x 230 mm
<b>H-0109</b>	FEMALE URINAL 124 mm x 125 mm x 216 mm
<b>H-0111</b>	ADULT BEDPAN WITH 2 HANDLES, RIGHT & LEFT 110 mm x 230 mm
<b>H-0112</b>	ADULT BEDPAN WITH 2 HANDLES, FRONT & BACK 110 mm x 230 mm

All Imperial Ware products are backed by our **1 year warranty**.

# JARS & COVERS

## PRODUCT SELECTION

ITEM #	DESCRIPTION (Ø x H)
H-0391	NEEDLE OR OINTMENT JAR WITH COVER (RAISED HANDLE) 76 mm x 63 mm
H-0392	NEEDLE OR OINTMENT JAR WITH COVER (RAISED HANDLE) 95 mm x 65 mm
H-0393	DRESSING JAR WITH COVER (FLUSH HANDLE) 90 mm x 110 mm
H-0394	DRESSING JAR WITH COVER (FLUSH HANDLE) 102 mm x 80 mm
H-0395	DRESSING JAR WITH COVER (FLUSH HANDLE) 102 mm x 130 mm
H-0396	DRESSING JAR WITH COVER (FLUSH HANDLE) 123 mm x 130 mm
H-0397	DRESSING JAR WITH COVER (FLUSH HANDLE) 123 mm x 166 mm
H-0398	DRESSING JAR WITH COVER (FLUSH HANDLE) 180 mm x 180 mm
H-0399	DRESSING JAR WITH COVER 62 mm x 82 mm
H-0400	DRESSING JAR WITH COVER (FLUSH HANDLE) 82 mm x 100 mm
H-1214	FORCEP JAR WITH HANDLE, 60 MM 65 mm x 100 mm
H-1215	FORCEP JAR WITH HANDLE, 75 MM 82 mm x 123 mm
H-1216	FORCEP JAR WITH HANDLE, 90 MM 98 mm x 145 mm
H-1217	FORCEP JAR WITH HANDLE, 100 MM 105 mm x 170 mm
H-1218	FORCEP JAR, 50 MM 55 mm x 80 mm



H-0398



H-1215

All Imperial Ware products are backed by our 1 year warranty.

# PETRI DISHES

H-1223

H-1224

H-1226

H-1227



H-1225

## PRODUCT SELECTION

ITEM #	DESCRIPTION (Ø x H)
H-1223	PETRI DISH WITH COVER 60 mm x 20 mm
H-1224	PETRI DISH WITH COVER 75 mm x 20 mm
H-1225	PETRI DISH WITH COVER 90 mm x 17 mm
H-1226	PETRI DISH WITH COVER 90 mm x 20 mm
H-1227	PETRI DISH WITH COVER 120 mm x 20 mm

All Imperial Ware products are backed by our 1 year warranty.

# EMESIS BASINS



H-0262

## PRODUCT SELECTION

ITEM #	DESCRIPTION (Ø x W x H)
H-0261	EMESIS BASIN, 6" WITHOUT COVER 172 mm x 96 mm x 34 mm
H-0262	COVER FOR H-0261 175 mm x 100 mm x 7 mm
H-0264	EMESIS BASIN, 8" WITHOUT COVER 208 mm x 98 mm x 36 mm
H-0265	COVER FOR H-0264 210 mm x 100 mm x 7 mm
H-0267	EMESIS BASIN, 10" WITHOUT COVER 256 mm x 119 mm x 45 mm
H-0268	COVER FOR H-0267 260 mm x 125 mm x 7 mm
H-0270	EMESIS BASIN, 12" WITHOUT COVER 310 mm x 152 mm x 55 mm
H-0271	COVER FOR H-0270 315 mm x 155 mm x 10 mm



H-0261

All Imperial Ware products are backed by our 1 year warranty.

# CUPS, TUMBLERS, AND PITCHERS

## PRODUCT SELECTION

ITEM #	DESCRIPTION (Ø x H)
H-0341	MEDICINE CUP, 1 OZ 48 mm x 30 mm
H-0342	MEDICINE CUP, 2 OZ 52 mm x 50 mm
H-0351	IODINE CUP, 2 OZ 60 mm x 30 mm
H-0352	IODINE CUP, 6 OZ 83 mm x 55 mm
H-0353	IODINE CUP, 14 OZ 120 mm x 60 mm
H-0381	TUMBLER, 5 OZ 65 mm x 75 mm
H-0382	TUMBLER, 7 OZ 77 mm x 80 mm
H-1233	TUMBLER, 10 OZ 80 mm x 98 mm
H-0383	TUMBLER, 12 OZ 83 mm x 115 mm
H-0502	SOLUTION PITCHER, 1 QT 100 mm x 150 mm
H-0503	SOLUTION PITCHER, 2 QT 130 mm x 165 mm
H-0505	BELL-SHAPED PITCHER, 0.9 LT 150 mm x 170 mm
H-0506	BELL-SHAPED PITCHER, 1.8 LT 160 mm x 205 mm
H-0507	BELL-SHAPED PITCHER, 2.7 LT 195 mm x 226 mm



H-0342



H-0352



H-0381



H-0502

All Imperial Ware products are backed by our 1 year warranty.

# ANCILLARY ITEMS



H-0241



H-0233

## PRODUCT SELECTION

ITEM #	DESCRIPTION (Ø x H)
H-0211	THERMOMETER / FORCEP JAR 27 mm x 102 mm
H-0212	THERMOMETER / FORCEP JAR 55 mm x 120 mm
H-0213	THERMOMETER / FORCEP JAR 55 mm x 190 mm
H-0232	TONGUE DEPRESSOR 22 mm x 170 mm
H-0233	TONGUE DEPRESSOR 23 mm x 190 mm
H-0241	TABLET COUNTING DISH L 215 mm x W 150 mm x H 135 mm



H-0212



H-0213

All Imperial Ware products are backed by our 1 year warranty.

100

[illegible]

# SURGMED

GROUP

Proudly Canadian, **The Surgmed™ Group of Companies** began operations in 1992 as a single distribution company known as **Advanced Surgi-Pharm Inc.** The company quickly became known as an industry leader and specialty distributor of value-added surgical products in the Canadian hospital market. Five years later, **Surgmed expanded upon its success with the launch of Batrik Medical Manufacturing Inc.** in 1997 and the subsequent acquisition of **Imperial Surgical Ltd.** one year later in 1998.


**Together, our three unique companies serve the Canadian and global healthcare markets with high quality and competitively priced products, all backed by world-class customer care.**



Founded in 1992, Advanced Surgi-Pharm Inc. is the distribution arm of Surgmed™ in Canada, representing leading US, international and domestic manufacturers for a wide range of specialty products with unique benefits for applications in the Operating Room, Women's Health, Labor & Delivery, ICU, and Sterile Processing areas. Equipped with a direct sales team, ASP provides effective nationwide coverage as the Canadian distributor of Batrik, Vygon, Rocket Medical, Cooper Surgical, SunMed, and other industry leading products including a wide-ranging portfolio of top-rated PPE products.

Established in 1997 and sold in 80 countries globally, Batrik Medical Manufacturing Inc. specializes in the manufacturing of endoscopy accessories, enzymatic solutions, infection control, central sterilization, and surgical devices. Batrik's in-house controlled environment production facility adheres to strictly enforced quality and regulatory control systems. They are registered as a device manufacturing establishment with the United States FDA and are in full compliance with CE and Health Canada Regulations. All Batrik products are rigorously tested, liability insured, and meet/exceed high international quality standards.

Established in 1935, Imperial Surgical Ltd. is one of Canada's oldest hospital capital equipment manufacturers, well-known for their stainless steel healthcare equipment for Sterile Processing, Operating Room, Examination and Treatment, Patient Handling, and Mortuary/ Pathology areas. Their growing portfolio of products also has applications in Infection Control & Prevention, an important area of focus for their future.



SURGMED GROUP

2165 46th Ave., Lachine, QC, Canada, H8T 2P1

Phone: +1.514.631.7988 | Toll Free: 800.661.5432 | Fax.: +1.514.631.9083

Email: [info@surgmed.com](mailto:info@surgmed.com) | Website: [www.batrik.com](http://www.batrik.com)