


BATRIK

BRIGHTBRUSH

THE RIGHT BRUSH FOR EVERY JOB

Tackle even the toughest cleaning jobs with ease! Our Color-Coded Twisted Wire Brushes are designed for precision, convenience, and superior performance—ensuring spotless results every time.



Color-Coded Bead

Each brush features a color-coded bead that indicates the bristle diameter and size, making it easy to select the right tool quickly and efficiently.



Hanging Loop

The brush includes a loop at the end for convenient hanging.

PRODUCT SELECTION

ITEM #	WIRE (Ø) (mm)	BEAD (Ø) (mm)	BEAD COLOR	DIMENSIONS Length x Brush (Ø) x Brush Length
12" Color-coded Twisted Wire Brushes				
22-12078B	0.7	4	Red	12" x .078" (2M)x 58mm
22-12118B	0.7	4	White	12" x .118" (3M)x 58mm
22-12157B	0.8	4	Blue	12" x .157" (4M)x 76mm
22-12197B	0.8	4	Green	12" x .197" (5M)x 76mm
22-12236B	0.8	4	Black	12" x .236" (6M)x 76mm
16" Color-coded Twisted Wire Brushes				
23-16197B	1.0	4	Green	16" X .197" (5M)x 63.5mm
23-16276B	1.0	4	Yellow	16" X .276" (7M)x 76mm
23-16394B	1.0	4	Pink	16" X .394" (10M)x 76mm
18" Color-coded Twisted Wire Brushes				
19-18078B	0.6	4	Red	18" x .078" (2M)x 63.5mm
19-18118B	0.8	4	White	18" x .118" (3M)x 63.5mm
19-18157B	0.8	4	Blue	18" x .157" (4M)x 63.5mm
19-18197B	1.0	4	Green	18" x .197" (5M)x 63.5mm
19-18236B	1.0	4	Black	18" x .236" (6M)x 76mm
24" Color-coded Twisted Wire Brushes				
24-24098B	0.8	4	Orange	24" X .098" (2.5M)x 63.5mm
24-24197B	0.8	4	Green	24" X .197" (5M)x 63.5mm
24-24276B	1.0	4	Yellow	24" X .276" (7M)x 82.6mm
24-24394B	1.0	4	Pink	24" X .394" (10M)x 82.6mm

PACKAGING: Bag of 3. Also available in bags of 50, add the suffix "-50".