

BATRIK

MEDICAL MANUFACTURING INC.

The **GLOBAL LEADER** in
SINGLE-SOURCE
BRUSH
SOLUTIONS

for **CSSD** and **HEALTHCARE**

PRODUCT BROCHURE



Established in 1997. Proudly Canadian. 🍁
Sold in over 85 countries worldwide.



Batrik Medical Manufacturing Inc., founded in 1997, is a globally recognized manufacturer with a strong presence in more than 85 countries through a network of reputable distributors. The company specializes in the production of Sterile Processing solutions, Surgical Accessories and Endoscopy.

Headquartered in Montreal (Lachine), Canada, **BATRIK** Medical upholds rigorous quality and regulatory standards. Our products undergo comprehensive testing, are liability-insured, and consistently exceed international quality requirements.



Batrik offers the world's most comprehensive brush portfolio, designed to support every critical step of cleaning and reprocessing in sterile processing departments and endoscopy units. Trusted in over 85 countries, Batrik delivers proven quality, innovative design, and dependable performance backed by rigorous regulatory compliance.

WHY CHOOSE BATRIK BRUSHES




Broadest range of styles, sizes and diameters, leading to full single-source convenience




Quality-controlled manufacturing



Medical-grade materials



Solutions for all CSSD and endoscopy products



OEM options available

Customer & Distributor Support

Our team provides full support to help end-users and distributors select the right brushes for their markets, including:

- Product selection guidance
- Technical specifications and compatibility information
- Samples for evaluation
- Marketing assets and educational materials
- Pricing, MOQ, and lead-time information

We are committed to helping our partners grow with a complete, single-source brush solution.

WHY PARTNER WITH US



QUALITY EXCELLENCE

Certifications

ISO 13485:2016, CE, FDA, MDSAP and EU MDR.

Quality Control

Rigorous testing and assurance.

Innovation

Ongoing R&D to meet evolving needs.




SUSTAINABILITY & ECO-FRIENDLY PRACTICES

Eco-Friendly Initiatives:

Reducing packaging waste, using sustainable fixtures, and embracing digital solutions.

Sustainability Commitment:

Hydro-powered manufacturing which minimizes our ecological footprint.



COMMUNITY ENGAGEMENT

Support for Local Events:

Participation in the Crohn's and Colitis Montreal Gutsy Walk.

Ongoing Donations:

Supporting the Missing Children's Network with every Imperial Surgical Kerri Bear Bassinet sold.



OUR COMPETITIVE ADVANTAGE

Over **25 years** of trusted expertise.

In-house design, manufacturing, and packaging.

Custom OEM/private label programs.

Sterile Processing

ARC™ Wire Style Brushes - Twisted Wire Brush With Plastic Handle	06
Double ARC™ Wire Brushes - Twisted Wire Brush With Plastic Handle Mixed Size / Same Size	10
In-Brush™ Antimicrobial Twisted Wire Brushes	13
BrightBrush - Color Beaded Twisted Wire Brushes NEW!	15
ComfyClean NEW!	16
ComfyClean Mini NEW!	17
ComfyClean Mini Max NEW!	18
Instrument Cleaning Flat Brush	19
Flat Brush Mini NEW!	20
Reusable Precision Brush	21
Trifecta™ Brush	22
Acetabular Reamer Brush	23
Orthopedic Rasps Cleaning Brush	24
All Purpose Reusable Brushes	25
Stainless Steel ARC™ Micro Brushes	26
Bottle & Drain Cleaning Brushes	27
Autoclave Cleaning Brush Kits	28

Endoscopy Solutions

ARC™ Channel Brushes	29
Daybrush™ Channel Brushes	32
Valve / Control Brushes	33
V Brushes	34

BATRIK

ARC™ Wire Style Brushes

SINGLE-USE and REUSABLE

You've got options with our **ARC™** Wire Brushes. Our single-use and reusable brushes are available in both flexible and rigid varieties and come in a multitude of sizes. **ARC™** Wire Brushes are guaranteed to last longer than twisted wire and preserve the integrity of your valuable instruments.

BRUSH SELECTION TIPS:

INSTRUMENT SPECIFICATIONS

What is the maximum length of the instrument(s) to be cleaned?

Add at least 2" (5 cm) to that length, so you have an adequate length to work with.

What is the inner diameter of the instrument?

Add at most 0.125" (3 mm) to determine the diameter of the bristle to ensure effective cleaning. The bristle diameter has to slightly exceed the inner tube diameter. Bristles that are too large will bend too much and not clean effectively.

BRUSH SELECTION

Simply find the section in our catalogue that best corresponds to the overall length and diameter of your instrument. As a rule, if your present brush is 8" (20.32 cm) in length and we only offer 10" (25.4 cm) length with your required diameter, that is entirely acceptable.

In many cases, our **ARC™** Wire Brushes will have more bristles on the brush. As an example, many 5 mm (0.19") diameter brushes are only 1" (2.5 cm) in length. At **BATRIK**, we strived to provide the greatest amount of bristles to ensure the most effective cleaning.

Too small

Bristles will not make sufficient contact and surfaces may go untouched.



Too large

Applies unnecessary pressure on bristles and may cause brush to get stuck in device.



Just right


Bristles make sufficient contact with the instrument.



RIGID



FLEXIBLE



CUSTOMIZABLE WIRE BRUSHES

A unique aspect of single-use **ARC™** Wire Brushes is that they can be customized to any length:

- Up to 2000 mm for Wire Brushes on plastic semi-rigid sheaths with brushes from 5 to 20 mm in diameter.
- Flexible plastic sheaths with brushes from 2-5 mm in diameter.

12" / 30 CM SINGLE-END WIRE BRUSH

	ITEM # SINGLE-USE	ITEM # REUSABLE	COLOR	TOTAL LENGTH	BRUSH (Ø)	BRUSH LENGTH	SHEATH (Ø)
FLEXIBLE BRUSH & SHEATH	41-31502-S	41-31502	White	12.401" / 31.5 cm	0.078" / 2 mm	0.59" / 15 mm	0.39" / 1 mm
	41-32503-S	41-32503	Yellow	12.795" / 32.5 cm	0.118" / 3 mm	0.98" / 25 mm	0.39" / 1 mm
	411-32504-S	411-32504	Violet	12.795" / 32.5 cm	0.157" / 4 mm	0.98" / 25 mm	0.61" / 1.55 mm
	411-32505-S	411-32505	Grey	12.795" / 32.5 cm	0.19" / 5 mm	0.98" / 25 mm	0.68" / 1.75 mm
RIGID BRUSH & SHEATH	41-35005-S	41-35005	White	13.779" / 35 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
	41-35006-S	41-35006	Green	13.779" / 35 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
	41-35008-S	41-35008	Blue	13.779" / 35 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
	41-35010-S	41-35010	Red	13.779" / 35 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
	41-35012-S	41-35012	Yellow	13.779" / 35 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
	41-35015-S	41-35015	Green	13.779" / 35 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
	41-35020-S	41-35020	Blue	13.779" / 35 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Single-Use: Bag of 20 / Reusable: Bag of 3

18" / 45 CM SINGLE-END WIRE BRUSH

	ITEM # SINGLE-USE	ITEM # REUSABLE	COLOR	TOTAL LENGTH	BRUSH (Ø)	BRUSH LENGTH	SHEATH (Ø)
FLEXIBLE BRUSH & SHEATH	41-46502-S	41-46502	White	18.307" / 46.5 cm	0.078" / 2 mm	0.59" / 15 mm	0.39" / 1 mm
	41-47503-S	41-47503	Yellow	18.7" / 47.5 cm	0.118" / 3 mm	0.98" / 25 mm	0.39" / 1 mm
	411-47504-S	411-47504	Violet	18.7" / 47.5 cm	0.157" / 4 mm	0.98" / 25 mm	0.61" / 1.55 mm
	411-47505-S	411-47505	Grey	18.7" / 47.5 cm	0.19" / 5 mm	0.98" / 25 mm	0.68" / 1.75 mm
RIGID BRUSH & SHEATH	41-50005-S	41-50005	White	19.685" / 50 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
	41-50006-S	41-50006	Green	19.685" / 50 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
	41-50008-S	41-50008	Blue	19.685" / 50 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
	41-50010-S	41-50010	Red	19.685" / 50 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
	41-50012-S	41-50012	Yellow	19.685" / 50 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
	41-50015-S	41-50015	Green	19.685" / 50 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
	41-50020-S	41-50020	Blue	19.685" / 50 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Single-Use: Bag of 20 / Reusable: Bag of 3

24" / 60 CM SINGLE-END WIRE BRUSH


	ITEM # SINGLE-USE	ITEM # REUSABLE	COLOR	TOTAL LENGTH	BRUSH (Ø)	BRUSH LENGTH	SHEATH (Ø)
FLEXIBLE BRUSH & SHEATH	41-61502-S	41-61502	White	24.212" / 61.5 cm	0.078" / 2 mm	0.59" / 15 mm	0.39" / 1 mm
	41-62503-S	41-62503	Yellow	24.6" / 62.5 cm	0.118" / 3 mm	0.98" / 25 mm	0.39" / 1 mm
	411-62504-S	411-62504	Violet	24.6" / 62.5 cm	0.157" / 4 mm	0.98" / 25 mm	0.61" / 1.55 mm
	411-62505-S	411-62505	Grey	24.6" / 62.5 cm	0.19" / 5 mm	0.98" / 25 mm	0.68" / 1.75 mm
RIGID BRUSH & SHEATH	41-65005-S	41-65005	White	25.59" / 65 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
	41-65006-S	41-65006	Green	25.59" / 65 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
	41-65008-S	41-65008	Blue	25.59" / 65 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
	41-65010-S	41-65010	Red	25.59" / 65 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
	41-65012-S	41-65012	Yellow	25.59" / 65 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
	41-65015-S	41-65015	Green	25.59" / 65 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
	41-65020-S	41-65020	Blue	25.59" / 65 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Single-Use: Bag of 20 / Reusable: Bag of 3

30" / 75 CM SINGLE-END WIRE BRUSH


	ITEM # SINGLE-USE	ITEM # REUSABLE	COLOR	TOTAL LENGTH	BRUSH (Ø)	BRUSH LENGTH	SHEATH (Ø)
FLEXIBLE BRUSH & SHEATH	41-76502-S	41-76502	White	30.118" / 76.5 cm	0.078" / 2 mm	0.59" / 15 mm	0.39" / 1 mm
	41-77503-S	41-77503	Yellow	30.511" / 77.5 cm	0.118" / 3 mm	0.98" / 25 mm	0.39" / 1 mm
	411-77504-S	411-77504	Violet	30.511" / 77.5 cm	0.157" / 4 mm	0.98" / 25 mm	0.61" / 1.55 mm
	411-77505-S	411-77505	Grey	30.511" / 77.5 cm	0.19" / 5 mm	0.98" / 25 mm	0.68" / 1.75 mm
RIGID BRUSH & SHEATH	41-80005-S	41-80005	White	31.496" / 80 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
	41-80006-S	41-80006	Green	31.496" / 80 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
	41-80008-S	41-80008	Blue	31.496" / 80 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
	41-80010-S	41-80010	Red	31.496" / 80 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
	41-80012-S	41-80012	Yellow	31.496" / 80 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
	41-80015-S	41-80015	Green	31.496" / 80 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
	41-80020-S	41-80020	Blue	31.496" / 80 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Single-Use: Bag of 20 / Reusable: Bag of 3



FLEXIBLE


Flexible brush diameter (Ø) ranging from 2 mm to 5 mm



RIGID

Rigid brush diameter (Ø) ranging from 5 mm to 20 mm

BRUSH SELECTION SYSTEM (Actual Size)



BATRIK

Reusable Double Arc™ Wire Brushes

MIXED SIZE BRUSHES

14" / 35 CM MIXED SIZE DOUBLE-END WIRE BRUSHES

ITEM #	COLOR	TOTAL LENGTH	BRUSH 1 (Ø)	BRUSH 2 (Ø)	BRUSH LENGTH	SHEATH (Ø)
40-35056	White	13.779" / 35 cm	0.19" / 5 mm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
40-35810	Blue	13.779" / 35 cm	0.31" / 8 mm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
40-35112	Yellow	13.779" / 35 cm	0.39" / 10 mm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
40-35152	Green	13.779" / 35 cm	0.59" / 15 mm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Bag of 3

20" / 50 CM MIXED SIZE DOUBLE-END WIRE BRUSHES

ITEM #	COLOR	TOTAL LENGTH	BRUSH 1 (Ø)	BRUSH 2 (Ø)	BRUSH LENGTH	SHEATH (Ø)
40-50056	White	19.685" / 50 cm	0.19" / 5 mm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
40-50810	Blue	19.685" / 50 cm	0.31" / 8 mm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
40-50112	Yellow	19.685" / 50 cm	0.39" / 10 mm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
40-50152	Green	19.685" / 50 cm	0.59" / 15 mm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Bag of 3

26" / 65 CM MIXED SIZE DOUBLE-END WIRE BRUSHES

ITEM #	COLOR	TOTAL LENGTH	BRUSH 1 (Ø)	BRUSH 2 (Ø)	BRUSH LENGTH	SHEATH (Ø)
40-65056	White	25.59" / 65 cm	0.19" / 5 mm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
40-65810	Blue	25.59" / 65 cm	0.31" / 8 mm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
40-65112	Yellow	25.59" / 65 cm	0.39" / 10 mm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
40-65152	Green	25.59" / 65 cm	0.59" / 15 mm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Bag of 3

31" / 80 CM MIXED SIZE DOUBLE-END WIRE BRUSHES

ITEM #	COLOR	TOTAL LENGTH	BRUSH 1 (Ø)	BRUSH 2 (Ø)	BRUSH LENGTH	SHEATH (Ø)
40-80056	White	31.496" / 80 cm	0.19" / 5 mm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
40-80810	Blue	31.496" / 80 cm	0.31" / 8 mm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
40-80112	Yellow	31.496" / 80 cm	0.39" / 10 mm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
40-80152	Green	31.496" / 80 cm	0.59" / 15 mm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Bag of 3

BATRIK

Reusable Double Arc™ Wire Brushes

SAME SIZE BRUSHES

14" / 35 CM SAME SIZE DOUBLE-END WIRE BRUSHES

ITEM #	COLOR	TOTAL LENGTH	BRUSH 1 (Ø)	BRUSH LENGTH	SHEATH (Ø)
40-35005	White	13.779" / 35 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
40-35006	Green	13.779" / 35 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
40-35008	Blue	13.779" / 35 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
40-35010	Red	13.779" / 35 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
40-35012	Yellow	13.779" / 35 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
40-35015	Green	13.779" / 35 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
40-35020	Blue	13.779" / 35 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Bag of 3

20" / 50 CM SAME SIZE DOUBLE-END WIRE BRUSHES

ITEM #	COLOR	TOTAL LENGTH	BRUSH 1 (Ø)	BRUSH LENGTH	SHEATH (Ø)
40-50005	White	19.685" / 50 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
40-50006	Green	19.685" / 50 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
40-50008	Blue	19.685" / 50 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
40-50010	Red	19.685" / 50 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
40-50012	Yellow	19.685" / 50 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
40-50015	Green	19.685" / 50 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
40-50020	Blue	19.685" / 50 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Bag of 3

26" / 65 CM SAME SIZE DOUBLE-END WIRE BRUSHES

ITEM #	COLOR	TOTAL LENGTH	BRUSH 1 (Ø)	BRUSH LENGTH	SHEATH (Ø)
40-65005	White	25.59" / 65 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
40-65006	Green	25.59" / 65 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
40-65008	Blue	25.59" / 65 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
40-65010	Red	25.59" / 65 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
40-65012	Yellow	25.59" / 65 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
40-65015	Green	25.59" / 65 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
40-65020	Blue	25.59" / 65 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Bag of 3

31" / 80 CM SAME SIZE DOUBLE-END WIRE BRUSHES

ITEM #	COLOR	TOTAL LENGTH	BRUSH 1 (Ø)	BRUSH LENGTH	SHEATH (Ø)
40-80005	White	31.496" / 80 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
40-80006	Green	31.496" / 80 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
40-80008	Blue	31.496" / 80 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
40-80010	Red	31.496" / 80 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
40-65112	Yellow	31.496" / 80 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
40-80015	Green	31.496" / 80 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
40-80020	Blue	31.496" / 80 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Bag of 3

MIXED SIZE BRUSHES

5 MM

6 MM



8 MM

10 MM



10 MM

12.5 MM



15 MM

20 MM



BATRIK


In-Brush™ Antimicrobial Twisted Wire Brushes

REUSABLE

Our In-Brush™ Antimicrobial Reusable Twisted Wire Brushes are designed for use with Frazier/Baron suction tubes. The twisted wires have a fine/medium diameter yet are rigid enough to ensure excellent cleaning capabilities.

Our wide assortment of nylon twisted wire brushes are ideal for large bore, uniquely shaped, and endoscopic instruments.

GOOD TO KNOW: OUR TWISTED WIRE BRUSHES FEATURE A PROTECTIVE BLUE TIP FOR ADDED SAFETY.



10" / 250 MM FRAZIER /BARON SUCTION TUBE BRUSHES

ITEM #	SIZE	BRUSH LENGTH	BRUSH (Ø)
21-10039	3 FR	0.5" / 12.7 mm	0.039" / 1 mm
21-10067	5 FR	0.5" / 12.7 mm	0.065" / 1.67 mm
20-10059	5 FR	1" / 25.4 mm	0.059" / 1.50 mm
20-10079	6 FR	1" / 25.4 mm	0.079" / 2 mm
20-10090	7 FR	1" / 25.4 mm	0.090" / 2.3 mm
20-10106	8 FR	1" / 25.4 mm	0.106" / 2.7 mm
20-10118	9 FR	1" / 25.4 mm	0.118" / 3 mm
20-10129	10 FR	1.5" / 38 mm	0.129" / 3.3 mm
20-10146	11 FR	1.5" / 38 mm	0.0146" / 3.7 mm
20-10157	12 FR	1.5" / 38 mm	0.157" / 4 mm

PACKAGING: Bag of 3.

Also available in bags of 50, add the suffix "-50".

12" / 305 MM NYLON TWISTED WIRE BRUSHES

ITEM #	BRUSH LENGTH	BRUSH (Ø)
22-12078	2.28" / 58 mm	0.078" / 2 mm
22-12118	2.28" / 58 mm	0.118" / 3 mm
22-12157	3" / 76 mm	0.157" / 4 mm
22-12197	3" / 76 mm	0.197" / 5 mm
22-12236	3" / 76 mm	0.236" / 6 mm
22-12393	3" / 76 mm	0.394" / 10 mm

PACKAGING: Bag of 3.

Also available in bags of 50, add the suffix "- 50".

16" / 406 MM NYLON TWISTED WIRE BRUSHES

ITEM #	BRUSH LENGTH	BRUSH (Ø)
23-16197	2.5" / 63.5 mm	0.197" / 5 mm
23-16236	2.5" / 63.5 mm	0.236" / 6 mm
23-16276	3" / 76 mm	0.276" / 7 mm
23-16394	3" / 76 mm	0.394" / 10 mm
23-16472	3" / 76 mm	0.472" / 12 mm
23-16591	3" / 76 mm	0.591" / 15 mm
23-16787	3" / 76 mm	0.787" / 20 mm

PACKAGING: Bag of 3.

Also available in bags of 50, add the suffix "-50".

18" / 457 MM NYLON TWISTED WIRE BRUSHES

ITEM #	BRUSH LENGTH	BRUSH (Ø)
19-18010	0.47" / 12 mm	0.039" / 1 mm
19-18015	0.5" / 12 mm	0.06" / 1.6 mm

PACKAGING: Bag of 3.

Also available in bags of 50, add the suffix "-50".

24" / 610 MM NYLON TWISTED WIRE BRUSHES

ITEM #	BRUSH LENGTH	BRUSH (Ø)
24-24098	2.5" / 63.5 mm	0.098" / 2.5 mm
24-24197	2.5" / 63.5 mm	0.197" / 5 mm
24-24276	3.25" / 82.6 mm	0.276" / 7 mm
24-24394	3.25" / 82.6 mm	0.394" / 10 mm
24-24472	3.25" / 82.6 mm	0.472" / 12 mm

PACKAGING: Bag of 3.

Also available in bags of 50, add the suffix "-50".



BRONCHO BRUSHES

ITEM #	DESCRIPTION	BRUSH LENGTH	BRUSH (Ø) INCHES / MM
26-06175	Small	6" / 152 mm	1.75" / 44.5 mm
26-08325	Large	8" / 203 mm	3.25" / 82.6 mm

PACKAGING: Bag of 3.

Also available in bags of 50, add the suffix "-50".

TRACHEOSTOMY BRUSHES

ITEM #	DESCRIPTION	TOTAL LENGTH	BRUSH LENGTH	BRUSH (Ø)
25-08500	Large 6-12.3 3 / 4"	8" / 203 mm	3.75" / 95 mm	0.5" / 12.7 mm
25-37550	Small 000-5.1 7 / 8"	4" / 93 mm	1.88" / 48 mm	0.5" / 12.7 mm

PACKAGING: Bag of 3. Also available in bags of 50, add the suffix "-50".

BATRIK


BRIGHTBRUSH ^{NEW!}

THE RIGHT BRUSH FOR EVERY JOB

Tackle even the toughest cleaning jobs with ease! Our Color-Coded Twisted Wire Brushes are designed for precision, convenience, and superior performance—ensuring spotless results every time.



PRODUCT SELECTION



Color-Coded Bead

Each brush features a color-coded bead that indicates the bristle diameter and size, making it easy to select the right tool quickly and efficiently.



Hanging Loop

The brush includes a loop at the end for convenient hanging.

ITEM #	WIRE (Ø) (mm)	BEAD (Ø) (mm)	BEAD COLOR	DIMENSIONS Length x Brush (Ø) x Brush Length
12" COLOR-CODED TWISTED WIRE BRUSHES				
22-12078B	0.7 mm	4 mm	Red	12" x .078" (2M)x 58 mm
22-12118B	0.7 mm	4 mm	White	12" x .118" (3M)x 58 mm
22-12157B	0.8 mm	4 mm	Blue	12" x .157" (4M)x 76 mm
22-12197B	0.8 mm	4 mm	Green	12" x .197" (5M)x 76 mm
22-12236B	0.8 mm	4 mm	Black	12" x .236" (6M)x 76 mm
16" COLOR-CODED TWISTED WIRE BRUSHES				
23-16197B	1.0 mm	4 mm	Green	16" X .197" (5M)x 63.5 mm
23-16276B	1.0 mm	4 mm	Yellow	16" X .276" (7M)x 76 mm
23-16394B	1.0 mm	4 mm	Pink	16" X .394" (10M)x 76 mm
18" COLOR-CODED TWISTED WIRE BRUSHES				
19-18078B	0.6 mm	4 mm	Red	18" x .078" (2M)x 63.5 mm
19-18118B	0.8 mm	4 mm	White	18" x .118" (3M)x 63.5 mm
19-18157B	0.8 mm	4 mm	Blue	18" x .157" (4M)x 63.5 mm
19-18197B	1.0 mm	4 mm	Green	18" x .197" (5M)x 63.5 mm
19-18236B	1.0 mm	4 mm	Black	18" x .236" (6M)x 76 mm
24" COLOR-CODED TWISTED WIRE BRUSHES				
24-24098B	0.8 mm	4 mm	Orange	24" X .098" (2.5M)x 63.5 mm
24-24197B	0.8 mm	4 mm	Green	24" X .197" (5M)x 63.5 mm
24-24276B	1.0 mm	4 mm	Yellow	24" X .276" (7M)x 82.6 mm
24-24394B	1.0 mm	4 mm	Pink	24" X .394" (10M)x 82.6 mm

PACKAGING: Bag of 3. Also available in bags of 50, add the suffix "-50".

BATRIK NEW!

COMFYCLEAN

Elevate Your Cleaning Power

Engineered for precision, comfort, and durability. Our new **COMFYCLEAN** range- featuring angled handles for a more natural, controlled grip. This design reduces hand fatigue during repetitive tasks and improves precision when cleaning intricate instruments. Made from long-lasting materials, this reusable brush not only boosts cleaning performance but also supports hospital sustainability goals by minimizing single-use waste and lowering overall costs.

KEY PRODUCT FEATURES

Four Bristle Hardness Levels for Every Cleaning Need

Select the optimal bristle stiffness to match the instrument type and contamination level, ensuring both thorough cleaning and equipment protection:

- Extra-Rigid Bristles**
 Nearly as rigid as steel, these bristles are ideal for removing stubborn debris from grooved or heavily soiled instruments.
- Rigid Bristles**
 Offers a balanced combination of stiffness and length for effective cleaning of standard instruments.
- Extra-Flexible Bristles**
 Designed for delicate instruments, including optical devices and both flexible and rigid endoscopes. Long fibers enable safe cleaning of sensitive or narrow areas without causing damage.
- Flexible Bristles**
 Ideal for cleaning fragile instruments, such as those made of optical or plastic materials, as well as tools with hollow or narrow areas. The semi-rigid fibers offer just enough flexibility to navigate tight spaces while remaining gentle enough to prevent damage.

Simplified Color-Coding System

A clear, visual color-coding system makes it easy to distinguish at a glance, reducing training time, minimizing risk of error, and improving overall workflow efficiency in fast-paced environments.

Transparent brush fibers allow for quick verification of its cleanliness.



50-40717
 EXTRA-RIGID
 Bristles Length 10 mm




50-40720
 EXTRA-FLEXIBLE
 Bristles Length 25 mm



50-40718
 RIGID
 Bristles Length 15 mm



50-40719
 FLEXIBLE
 Bristles Length 15 mm



PRODUCT SELECTION

ITEM #	DESCRIPTION	COLOR	BRUSH LENGTH	BRISTLE SIZE	PACKAGING
50-40717	Reusable ComfyClean Ergonomic Brush - Extra Rigid	Blue	222 mm	75 mm x 18 mm x 10 mm	Bag of 3
50-40718	Reusable ComfyClean Ergonomic Brush - Rigid	Green	222 mm	75 mm x 18 mm x 15 mm	Bag of 3
50-40719	Reusable ComfyClean Ergonomic Brush - Flexible	Yellow	222 mm	75 mm x 18 mm x 15 mm	Bag of 3
50-40720	Reusable ComfyClean Ergonomic Brush - Extra Flexible	Red	222 mm	75 mm x 18 mm x 25 mm	Bag of 3

Reusable models are also available in bags of 50. Add the suffix “ - 50”

BATRIK **NEW!** COMFYCLEAN MINI

Compact. Precise. Reusable.

The **COMFYCLEAN** MINI offers the same high-performance cleaning as the original **COMFYCLEAN**, in a smaller, more agile format. Its angled handle ensures a comfortable and controlled grip. Ideal for detailed cleaning in confined areas, this durable, reusable brush supports both precision and sustainability.

KEY PRODUCT FEATURES

Two Bristle Hardness Levels for Every Cleaning Need.

Select the optimal bristle stiffness to match the instrument type and contamination level, ensuring both thorough cleaning and equipment protection:

- **Extra-Rigid Bristles**


Nearly as rigid as steel, these bristles are ideal for removing stubborn debris from grooved or heavily soiled instruments.

- **Extra-Flexible Bristles**

Designed for delicate instruments, including optical devices and both flexible and rigid endoscopes. Long fibers enable safe cleaning of sensitive or narrow areas without causing damage.

Simplified Color-Coding System

A clear, visual color-coding system makes it easy to distinguish at a glance, reducing training time, minimizing risk of error, and improving overall workflow efficiency in fast-paced environments.



50-40417
EXTRA-RIGID
Bristles Length 10 mm



50-40420
EXTRA-FLEXIBLE
Bristles Length 25 mm

PRODUCT SELECTION

ITEM #	DESCRIPTION	COLOR	BRUSH LENGTH	BRISTLE SIZE	PACKAGING
50-40417	Reusable ComfyClean Ergonomic Mini Brush - Extra-Rigid	Blue	187 mm	39 mm x 10 mm x 10 mm	Bag of 3
50-40420	Reusable ComfyClean Ergonomic Mini Brush - Extra-Flexible	Red	187 mm	39 mm x 10 mm x 25 mm	Bag of 3

Reusable models are also available in bags of 50. Add the suffix " - 50"

BATRIK **COMFYCLEAN** MINI MAX **NEW!**

Maximum Power. Minimal Size.

The **COMFYCLEAN** MINI MAX is a compact, reusable brush designed for maximum cleaning power in the most hard-to-reach areas. As part of the **COMFYCLEAN** family, it features the same ultra-comfortable ergonomic grip, with a reinforced tip and high-density fibers for enhanced scrubbing precision where it matters most.

KEY PRODUCT FEATURES

Four Bristle Hardness Levels for Every Cleaning Need

Select the optimal bristle stiffness to match the instrument type and contamination level, ensuring both thorough cleaning and equipment protection:

- Extra-Rigid Bristles**
 Nearly as rigid as steel, these bristles are ideal for removing stubborn debris from grooved or heavily soiled instruments.
- Rigid Bristles**
 Offers a balanced combination of stiffness and length for effective cleaning of standard instruments.
- Extra-Flexible Bristles**
 Designed for delicate instruments, including optical devices and both flexible and rigid endoscopes. Long fibers enable safe cleaning of sensitive or narrow areas without causing damage.
- Flexible Bristles**
 Ideal for cleaning fragile instruments, such as those made of optical or plastic materials, as well as tools with hollow or narrow areas. The semi-rigid fibers offer just enough flexibility to navigate tight spaces while remaining gentle enough to prevent damage.

Simplified Color-Coding System

A clear, visual color-coding system makes it easy to distinguish at a glance, reducing training time, minimizing risk of error, and improving overall workflow efficiency in fast-paced environments.



50-40317
 EXTRA-RIGID
 Bristles Length 10 mm




50-40320
 EXTRA-FLEXIBLE
 Bristles Length 25 mm



50-40318
 RIGID
 Bristles Length 15 mm



50-40319
 FLEXIBLE
 Bristles Length 15 mm



PRODUCT SELECTION

ITEM #	DESCRIPTION	COLOR	BRUSH LENGTH	BRISTLE SIZE	PACKAGING
50-40317	Reusable ComfyClean Mini Max Ergonomic Brush - Extra Rigid	Blue	195 mm	55 mm x 10 mm x 10 mm	Bag of 3
50-40318	Reusable ComfyClean Mini Max Ergonomic Brush - Rigid	Green	195 mm	55 mm x 10 mm x 15 mm	Bag of 3
50-40319	Reusable ComfyClean Mini Max Ergonomic Brush - Flexible	Yellow	195 mm	55 mm x 10 mm x 15 mm	Bag of 3
50-40320	Reusable ComfyClean Mini Max Ergonomic Brush - Extra Flexible	Red	195 mm	55 mm x 10 mm x 25 mm	Bag of 3

Reusable models are also available in bags of 50. Add the suffix "- 50"

BATRIK

Instrument Cleaning Flat Brushes

Our Instrument Cleaning Flat Brushes are specifically designed for hospital use in the pre-clean and preparation of instruments.

- Long ergonomic handle permits ease-of-use
- Extra long length of brush head (60 [±5] mm in Single-Use and 75 [±5] mm in Reusable)
- Transparent brush fibers allow for quick verification of its cleanliness
- The plastic handle will not release any toxic components
- Available in 6 models (4 reusable/autoclavable and 2 single-use)



FLEXIBLE

This brush is designed for large instruments that require a more gentle cleaning application such as optics, plastics, or hard to reach spot.



EXTRA FLEXIBLE

The very long, flexible bristles eliminate any stress during sanitation, and makes this the ideal brush for fragile instruments as well as both pliable and rigid endoscopes. The length of the bristles assist in cleaning areas that would otherwise be difficult to reach, such as the accessory channel or tip deflection of pliable endoscopes.



RIGID

The rigidity of this brush, along with the length of the bristles and handle, improves and simplifies instrument cleaning.



EXTRA RIGID

Instrument safe polyamide fibers offer steel-like rigidity and allows deep cleaning of soiled instruments, particularly in grooved or uneven surfaces.



50-40919



50-40921



50-40920



50-40918



50-40922



50-40917

PRODUCT SELECTION

ITEM #	BRUSH TYPE	DESCRIPTION	COLOR	BRUSH LENGTH	BRUSH SIZE	PACKAGING
50-40917	Extra Rigid	Reusable Autoclavable Flat Brush	Blue	220 mm	75 mm x 20 mm x 10 mm	Bag of 3
50-40918	Rigid	Reusable Autoclavable Flat Brush	Green	220 mm	75 mm x 20 mm x 15 mm	Bag of 3
50-40919	Flexible	Reusable Autoclavable Flat Brush	Yellow	220 mm	75 mm x 20 mm x 15 mm	Bag of 3
50-40920	Extra Flexible	Reusable Autoclavable Flat Brush	Red	220 mm	75 mm x 20 mm x 25 mm	Bag of 3
50-40921	Flexible	Single-Use Flat Brush	Orange	220 mm	60 mm x 20 mm x 15 mm	Bag of 20
50-40922	Rigid	Single-Use Flat Brush	Light Blue	220 mm	60 mm x 20 mm x 15 mm	Bag of 20

Reusable models are also available in bags of 50. Add the suffix "- 50".

BATRIK

NEW!

Flatbrush Mini

Need to clean smaller instruments with precision?

This compact, reusable brush delivers the trusted performance of our signature Flatbrush in a reduced format, perfect for detailed, high-precision work. Built from long-lasting materials, this brush not only enhances cleaning performance but also supports hospital sustainability goals by reducing single-use waste and associated costs.

KEY PRODUCT FEATURES

Two Bristle Hardness Levels for Every Cleaning Need

Select the optimal bristle stiffness to match the instrument type and contamination level, ensuring both thorough cleaning and equipment protection:

- Extra-Rigid Fibers**
 Nearly as rigid as steel, these bristles are ideal for removing stubborn debris from grooved or heavily soiled instruments.
- Extra-Flexible Fibers**
 Designed for delicate instruments, including optical devices and both flexible and rigid endoscopes. Long fibers enable safe cleaning of sensitive or narrow areas without causing damage.

Simplified Color-Coding System

A clear, visual color-coding system makes it easy to distinguish at a glance, reducing training time, minimizing risk of error, and improving overall workflow efficiency in fast-paced environments.

Transparent brush fibers allow for quick verification of its cleanliness.

50-40617
EXTRA-RIGID



EXTRA-RIGID
reusable cleaning brush
Fiber Length 10 mm



EXTRA-FLEXIBLE
reusable cleaning brush
Fiber Length 25 mm

BRISTLED PART
L x l (mm)
38 x 9



50-40620
EXTRA-FLEXIBLE



PRODUCT SELECTION

ITEM #	DESCRIPTION	COLOR	BRUSH LENGTH	BRISTLE SIZE	PACKAGING
50-40617	Reusable Mini Flat Brush - Extra Rigid	Blue	185 mm	38 mm x 9 mm x 10 mm	Bag of 3
50-40620	Reusable Mini Flat Brush - Extra Flexible	Red	185 mm	38 mm x 9 mm x 25 mm	Bag of 3

Reusable models are also available in bags of 50. Add the suffix "-50"


BATRIK

Reusable Precision Brush

INSTRUMENT CLEANING

Our Reusable Precision Instrument Cleaning Brush is the ideal brush for the quick removal of difficult debris. It cleans externally as well as deep into instrument jaws. Although the firm and durable bristles feel like stainless steel, there is absolutely no risk of instrument damage or personal injury.

- Bristles have an extra firm density and superior durability
- Polyamide (TYNEX®) bristles
- Steamable thermoplastic ergonomic handle
- Steam / chemical sterilizable
- Latex free



GOOD TO KNOW: ELIMINATES THE NEED FOR THE STAINLESS STEEL TOOTHBRUSH-STYLE BRUSH, WHICH CAN DAMAGE VALUABLE INSTRUMENTS!

PRODUCT SELECTION

ITEM #	DESCRIPTION	TOTAL LENGTH
50-40800	Reusable Precision External Cleaning Brush	6.29" / 160 mm

PACKAGING: Bag of 3

BATRIK

Trifecta™ Brush

INSTRUMENT CLEANING

Introducing the **Trifecta™** Brush, a versatile 3-in-1 tool designed to meet all your instrument cleaning needs. This reusable instrument brush is available in two styles: one with all-nylon bristles for gentle cleaning and another with nylon bristles on two sides and stainless steel bristles on the third, ideal for tackling tougher buildup and stubborn cleaning tasks. This innovative design ensures effective, thorough maintenance for a wide range of instruments.

TRIFECTA™ BRUSH - NYLON

1 NYLON BRISTLES

Gentle, transparent nylon bristles provide an effective yet safe cleaning solution. This side allows you to verify cleanliness easily, without risking damage to delicate surfaces. Ideal for a light, non-invasive clean.

2 PRECISION TIP

Reach those hard-to-clean spots with the precision tip. This small, firm-bristle tip is designed for tight, hard-to-reach areas, such as instrument jaws and box joints. It ensures a deep clean where it's most needed, without the risk of damage or injury to delicate parts.

3 FIRM NYLON BRISTLES

For tougher residues, the firm nylon side delivers the cleaning power you need. With bristles that feel similar to stainless steel but are safe for your instruments, this side easily handles stubborn dirt and buildup.



TRIFECTA™ BRUSH - NYLON & STAINLESS STEEL



→ STAINLESS STEEL BRISTLES ARE MORE DURABLE AND EFFECTIVE FOR CLEANING STUBBORN BUILDUP.

PRODUCT SELECTION

ITEM #	COLOR	DESCRIPTION	LENGTH TOTAL	BRISTLE LENGTH	BRISTLES 1 Ø	BRISTLES 2 Ø	BRISTLES 3 Ø
28-23012	Blue	Triple Nylon Brush	8.66" / 220 mm	1.20" x 1.20" x 0.47" / 30.5 x 30.5 x 12.0 mm	0.47" / 12.0 mm	0.47" / 12.0 mm	0.315" / 8.0 mm
29-23012	Blue	Triple Nylon / Stainless Steel	8.66" / 220 mm	1.20" x 1.20" x 0.47" / 30.5 x 30.5 x 12.0 mm	0.47" / 12.0 mm	0.47" / 12.0 mm	0.47" / 12.0 mm

PACKAGING: Bag of 3 - For a bag of 50, add "-50" following the item code.

BATRIK

Acetabular Reamer Brush

Our double-ended cleaning brush is manufactured exclusively for cleaning acetabular bone reamers. The “Pig Tail” rounded style circular head brush is designed to effectively clean inside the reamer, while the rigid straight brush on the reverse side is used for cleaning individual cutting holes. Special fusion of the brush heads to the hand piece along with its high quality antimicrobial medical grade bristles ensures a durable, reusable, and highly effective cleaning tool.

**STRAIGHT BRUSH****ROUND BRUSH**

PRODUCT SELECTION

ITEM #	DESCRIPTION	HANDLE LENGTH	STRAIGHT BRUSH	ROUND BRUSH HEAD	PACKAGING
50-ACETAB	Acetabular Reamer Cleaning Brush	4.72" / 120 mm Ø 0.47" / 12 mm	1.18" / 30 mm Ø 0.3125" / 8 mm	2" / 50 mm Ø 1.75" / 45 mm	Bag of 3


NOTE: Specifications subject to small variances

BATRIK

Reusable Orthopedic Rasps Cleaning Brush

Our Reusable Orthopedic Rasps Cleaning Brush is used to clean orthopedic rasps and other similar instruments used in orthopedic surgery.

- Long lasting and can withstand many sterilization cycles
- Ergonomic steamable thermoplastic handle
- Antibacterial polyamide bristles
- Reusable
- Latex free



PRODUCT SELECTION

ITEM #	TOTAL LENGTH	HANDLE LENGTH	BRUSH SIZE	PACKAGING
50-RASP	9" / 23 cm	5.5" / 14 cm	0.78" x 2.7" x 3.5" 2 cm x 7 cm x 9 cm	Bag of 3

BATRIK

All Purpose Reusable Brushes

Plastic Handle Brushes

We offer a wide selection of brushes for instrument cleaning applications. These brushes are ideal for large bore, uniquely shaped, and endoscopic instruments.

In-Brush™ Plastic Handle

Available in regular toothbrush style (nylon or stainless) and double-ended brushes.



PRODUCT SELECTION

ITEM #	DESCRIPTION	COLOR	BRUSH LENGTH	BRISTLE SIZE	PACKAGING
28-72550	Plastic Single End Nylon Reusable Handle Brush	Black	184 mm	38 mm x 12.7 mm	Bag of 3
29-72550	Plastic Single End Stainless Steel Reusable Handle Brush	Black	184 mm	38 mm x 12.7 mm	Bag of 3
29-70050	Plastic Single End Brass Reusable Handle Brush	Black	184 mm	38 mm x 12.7 mm	Bag of 3
28-69394	Plastic Double End Nylon Reusable Handle Brush	Black	175 mm	35 mm x 10 mm 25 mm x 5 mm	Bag of 3
29-69394	Plastic Double End Stainless Steel Reusable Handle Brush	Black	175 mm	35 mm x 10 mm 25 mm x 5 mm	Bag of 3
50-40911	Hand Scrub Nylon Brush Reusable autoclavable	White	104 mm	95 mm x 28 mm	Bag of 3

Reusable models are also available in bags of 50. Add the suffix "- 50"

NEW!

Engineered for precision and durability, this reusable dual-ended is designed to tackle even the most complex instruments. Its ultra-fine tip reaches into the smallest, hardest-to-access areas—delivering exceptional hygiene with every use. Built from long-lasting materials, this brush not only enhances cleaning performance but also supports hospital sustainability goals by reducing single-use waste and associated costs.



PRODUCT SELECTION

ITEM #	DESCRIPTION	COLOR	BRUSH LENGTH	BRISTLE SIZE 1	BRISTLE SIZE 2	PACKAGING
50-40517	Reusable Extra Rigid Flat Brush Dual Ended	Blue	183 mm	38 mm x 9 mm x 10 mm	28 mm x 3 mm x 8 mm	Bag of 3
50-40518	Reusable Rigid Flat Brush Dual Ended	Green	183 mm	38 mm x 9 mm x 15 mm	28 mm x 3 mm x 13 mm	Bag of 3
50-40520	Reusable Extra Flexible Flat Brush Dual Ended	Red	183 mm	38 mm x 9 mm x 25 mm	28 mm x 3 mm x 23 mm	Bag of 3

Reusable models are also available in bags of 50. Add the suffix "- 50"

BATRIK

Stainless Steel ARC™ Micro Brushes

Our **Stainless Steel ARC™** Micro Brushes are some of the smallest brushes in the industry. They are designed to clean channels with an inner diameter ranging from 0.5 mm to 1.40 mm.

- Flexible stainless steel sheath will not bend like the regular metal twisted wire brushes
- Polyamide bristles prevents or slows bacterial growth between sterilizations
- Reusable*
- Available in two lengths

TIP: USING THE PULL-THROUGH TECHNIQUE IS HIGHLY RECOMMENDED WITH THESE BRUSHES. SIMPLY INSERT THE BLUNT, NON-BRUSH END FIRST.



PRODUCT SELECTION

ITEM #	TOTAL LENGTH	BRUSH SIZE (Ø)	BRUSH LENGTH	SHEATH (Ø)	PACKAGING
35-03015	12" / 30 cm	0.06" / 1.5 mm	0.6" / 15 mm	0.02" / 0.40 mm	6 / Tube
35-04515	18" / 45 cm	0.06" / 1.5 mm	0.6" / 15 mm	0.02" / 0.40 mm	6 / Tube

*Although all materials used in the fabrication of these brushes can be re-processed repeatedly, due to the delicate nature of these brushes re-use is conditional on their handling and application. Re-use it at the user's discretion.

BATRIK

Bottle & Drain Cleaning Brushes

REUSABLE INSTRUMENT CLEANING BRUSHES


Designed specifically for hospital applications, our Reusable Instrument Cleaning Brushes feature long and flexible bristles to maximize the effectiveness of large-surface cleaning.

- High quality polyamide clear bristles
- Easy to identify debris
- Autoclavable
- Latex free

REUSABLE JAR/BOTTLE BRUSH


 Ergonomic
 Handle

STERILIZER DRAIN BRUSH


 Blue Tip
 for added safety

 Heavy Duty
 Stainless Wire

PRODUCT SELECTION

BLUE DELRIN PLASTIC

ITEM #	DESCRIPTION	HANDLE LENGTH	BRUSH LENGTH	BRUSH (Ø)	TOTAL LENGTH	PACKAGING
50-10295	Reusable Jar / Bottle Brush	120 mm ± 2 mm	120 mm ± 5 mm	70 mm ± 2 mm	13.19" / 355 mm	Bag of 3
50-10224	Sterilizer Drain Cleaning Brush	N/A	4" / 10.1 cm	2" / 5.08 cm	24" / 609 mm	Bag of 3

BATRIK

Autoclave Cleaning Brush Kits



ROUND BRUSH

- Features 3 holes for adjustable pole placement
- Adjustable brush placement (parallel or perpendicular) and longer brush length provides easy access to hard-to-reach corners
- Made with FDA approved polyester bristles

HAND BRUSH

- Ergonomic handle for better hand comfort
- Longer, FDA approved polyester bristles
- Easy to clean (no cavity)

PRODUCT SELECTION

REGULAR AUTOCLAVE CLEANING BRUSH KIT

ITEM # (KIT)	DESCRIPTION	ITEM # (UNIT)	DESCRIPTION	DIMENSIONS
50-CLAVEREG	Autoclave REGULAR Brush Cleaning Kit	50-35900	1 Hand Brush	Length: 3.5" x 8.3" / 9 cm x 21 cm
		50-55800	1 Rounded Corner Brush (Head Brush)	Head size: 5.5 x 9.8" / 14 x 25 cm
		50-EXTPOLE	1 Anodized Aluminum Extension Pole	2-4 ft / 68-120 cm

SUPERIOR AUTOCLAVE CLEANING BRUSH KIT

ITEM # (KIT)	DESCRIPTION	ITEM # (UNIT)	DESCRIPTION	DIMENSIONS
50-CLAVESUP	Autoclave SUPERIOR Brush Cleaning Kit	50-35900	1 Hand Brush	Length: 3.5" x 8.3" / 9 cm x 21 cm
		50-55800	1 Rounded Corner Brush (Head Brush)	Head size: 5.5 x 9.8" / 14 x 25 cm
		50-EXTPOLE	1 Anodized Aluminum Extension Pole	2-4 ft / 68-120 cm
		SF-100	2 SAVEFACE™ Face Shields	9" / 23 cm

REPLACEMENT COMPONENTS



HAND BRUSH
50-35900

Packaging: 1



BRUSH HEAD
50-55800

Packaging: 1



EXTENSION POLE
50-EXTPOLE

Packaging: 1



SAVEFACE™
FACE SHIELDS SF-100


Packaging: 100 / Box

BATRIK

ARC™ Channel Brushes

BATRIK offers one of the largest selection of Endoscopy Channel Cleaning Brushes in the industry. Our brushes are available in a variety of lengths, materials, and styles to accommodate the widest range of endoscopic devices used in the healthcare industry today.

No matter what your requirements, you'll be able to find the brush that suits your specific needs. If you don't see the brush you require, simply call us. We have over 50 different brushes available and can custom fabricate to your precise specifications.



OPTIONAL KIT ADD-ONS

ADD THE SUFFIX:

"F" to receive a kit with our Flat Sponge. Applicable on certain models only, call for details. (Box of 40)

"K" to receive a kit with the single use Valve / Control Cleaning Brush 5/11 Double Brushes (item # 27-57787). Applicable on all models. (Box of 50)

"FK" to receive both the sponge and valve cleaning brushes. (Box of 40)

SERIES 333 & 311



PRODUCT SELECTION

5 MM CYLINDRICAL BRUSHES – SHEATH Ø 1.85 MM - PUNCHY PINK POLYPROPYLENE

ITEM #	LENGTH	BRUSH LENGTH	BRUSH (Ø)	SHEATH (Ø)	BRUSH ENDS	PACKAGING
311-23050	90" / 230 cm	0.94" / 24 mm	0.20" / 5 mm	0.073" / 1.85 mm	Single	Box of 50
311-23050-2	90" / 230 cm	0.94" / 24 mm	0.20" / 5 mm	0.073" / 1.85 mm	Double	Box of 50
311-18050	70" / 180 cm	0.94" / 24 mm	0.20" / 5 mm	0.073" / 1.85 mm	Single	Box of 50
311-10050	39" / 100 cm	0.94" / 24 mm	0.20" / 5 mm	0.073" / 1.85 mm	Single	Box of 50

5/6 MM TAPERED BRUSHES – SHEATH Ø 1.85 MM - PUNCHY PINK POLYPROPYLENE

ITEM #	LENGTH	BRUSH LENGTH	BRUSH (Ø)	SHEATH (Ø)	BRUSH ENDS	PACKAGING
333-23056	90" / 230 cm	0.79" / 20 mm	0.20", 0.24" / 5, 6 mm	0.073" / 1.85 mm	Single	Box of 50
333-23056-2	90" / 230 cm	0.79" / 20 mm	0.20", 0.24" / 5, 6 mm	0.073" / 1.85 mm	Double	Box of 50

SERIES 300 & 322



PRODUCT SELECTION

5 MM CYLINDRICAL BRUSHES – SHEATH Ø 1.75 MM - METALLIC GREY DELRIN

ITEM #	LENGTH	BRUSH LENGTH	BRUSH (Ø)	SHEATH (Ø)	BRUSH ENDS	PACKAGING
300-18050	70" / 180 cm	0.79" / 20 mm	0.20" / 5 mm	0.070" / 1.75 mm	Single	Box of 50
300-23050	90" / 230 cm	0.79" / 20 mm	0.20" / 5 mm	0.070" / 1.75 mm	Single	Box of 50
300-23050-2	90" / 230 cm	0.79" / 20 mm	0.20" / 5 mm	0.070" / 1.75 mm	Double	Box of 50

5/6 MM TAPERED BRUSHES – SHEATH Ø 1.75 MM - METALLIC GREY DELRIN

ITEM #	LENGTH	BRUSH LENGTH	BRUSH (Ø)	SHEATH (Ø)	BRUSH ENDS	PACKAGING
322-23056	90" / 230 cm	0.79" / 20 mm	0.20", 0.24" / 5, 6 mm	0.070" / 1.75 mm	Single	Box of 50
322-23056-2	90" / 230 cm	0.79" / 20 mm	0.20", 0.24" / 5, 6 mm	0.070" / 1.75 mm	Single	Box of 50

SERIES 344




PRODUCT SELECTION

3 MM CYLINDRICAL BRUSHES – SHEATH Ø 1.55 MM - METALLIC PURPLE DELRIN

ITEM #	LENGTH	BRUSH LENGTH	BRUSH (Ø)	SHEATH (Ø)	BRUSH ENDS	PACKAGING
344-23030	90" / 230 cm	0.79" / 20 mm	0.12" / 3 mm	0.060" / 1.55 mm	Single	Box of 50
344-18030	70" / 180 cm	0.79" / 20 mm	0.12" / 3 mm	0.060" / 1.55 mm	Single	Box of 50
344-10030	39" / 100 cm	0.79" / 20 mm	0.12" / 3 mm	0.060" / 1.55 mm	Single	Box of 50

SERIES 35



PRODUCT SELECTION

2.5 MM CYLINDRICAL BRUSHES – SHEATH Ø 1.00 MM - YELLOW DELRIN - WITH TIP LIMITED AVAILABILITY

ITEM #	LENGTH	BRUSH LENGTH	BRUSH (Ø)	SHEATH (Ø)	BRUSH ENDS	PACKAGING
35-23025	90" / 230 cm	0.59" / 15 mm	0.098" / 2.5 mm	0.039" / 1.00 mm	Single	Box of 50
35-12025	47.2" / 120 cm	0.59" / 15 mm	0.098" / 2.5 mm	0.039" / 1.00 mm	Single	Box of 50

2 MM CYLINDRICAL BRUSHES – SHEATH Ø 1.00 MM - YELLOW DELRIN - WITHOUT TIP

ITEM #	LENGTH	BRUSH LENGTH	BRUSH (Ø)	SHEATH (Ø)	BRUSH ENDS	PACKAGING
35-23020	90" / 230 cm	0.59" / 15 mm	0.079" / 2.0 mm	0.039" / 1.00 mm	Single	Box of 50
35-12020	47.2" / 120 cm	0.59" / 15 mm	0.079" / 2.0 mm	0.039" / 1.00 mm	Single	Box of 50
35-18020	70.8" / 180 cm	0.59" / 15 mm	0.079" / 2.0 mm	0.039" / 1.00 mm	Single	Box of 50


BATRIK

DAYBRUSH™ Channel Brushes

REUSABLE

Our Reusable **DAYBRUSH™ Channel Brushes** can be used for cleaning endoscopes for an entire day. In fact, it's been proven to effectively clean up to 10 scopes before the need to discard.*

- Soft tip ensures easy manipulation within the channel with no risk of damage
- The back-end is sealed and softened to prevent fluid from getting trapped inside the tube. The softening process also prevents damage if the cleaning brush is used in reverse.
- Can be steam sterilized.
- Bristle volume is 20 mm in length – nearly double most brushes on the market. A large brush provides more surface contact during the cleaning process.
- Significant per use savings when compared to other reusable or disposable channel cleaning brushes.



(*)Discard after cleaning 10 scopes or sooner if any brush degradation is noticed. After ten applications, the brush may have the appearance that it can continue to clean effectively, but we discourage prolonged use due to risk of cross contamination. Actual number of uses will vary based on handling and processing.

PRODUCT SELECTION

5 MM CYLINDRICAL BRUSHES – SHEATH Ø 1.75 MM

ITEM #	COLOUR	LENGTH	BRUSH LENGTH	BRUSH (Ø)	SHEATH (Ø)	BRUSH ENDS	PACKAGING
30-18050	White	70" / 180 cm	0.79" / 20 mm	0.20" / 5 mm	0.07" / 1.75 mm	Single	Box of 50
30-18050-T	Yellow	70" / 180 cm	0.79" / 20 mm	0.20" / 5 mm	0.07" / 1.75 mm	Single	Box of 50
30-18050-2T	Yellow	70" / 180 cm	0.79" / 20 mm	0.20" / 5 mm	0.07" / 1.75 mm	Double	Box of 50
30-23050	White	90" / 230 cm	0.79" / 20 mm	0.20" / 5 mm	0.07" / 1.75 mm	Single	Box of 50
30-23050-T	Yellow	90" / 230 cm	0.79" / 20 mm	0.20" / 5 mm	0.07" / 1.75 mm	Single	Box of 50
30-23050-2T	Yellow	90" / 230 cm	0.79" / 20 mm	0.20" / 5 mm	0.07" / 1.75 mm	Double	Box of 50

OPTIONAL KIT ADD-ONS

ADD THE SUFFIX:

"F" to receive a kit with our Flat Sponge. Applicable on certain models only, call for details. (Box of 40)

"K" to receive a kit with the single use Valve / Control Cleaning Brush 5/11 Double Brushes (item # 27-57787). Applicable on all models. (Box of 50)

"FK" to receive both the sponge and valve cleaning brushes. (Box of 40)

BATRIK

Valve/Control Brushes

LARGEST SELECTION AVAILABLE


27-05197



27-57787



40-18511



36-10010



27-15311



PRODUCT SELECTION


ITEM #	DESCRIPTION	LENGTH	BRUSH 1 LENGTH	BRUSH 2 LENGTH	BRUSH 1 (Ø)	BRUSH 2 (Ø)
27-05197	Double Trumpet 5/7 White Color	5.00" / 13 cm	0.20" / 5 mm	1.00" / 25 mm	0.20" / 5 mm	0.28" / 7 mm
PACKAGING: Bag of 3 (single-use recommended). Volume Packaging: Bag of 50 use ITEM # 27-05197-50						
27-57787	Double Universal 5/11 Yellow Color	6.00" / 15 cm	1.18" / 30 mm	0.79" / 20 mm	0.43" / 11 mm	0.20" / 5 mm
PACKAGING: Bag of 3 (single-use recommended). Volume Packaging: Bag of 50 use ITEM # 27-57787-50						
40-18511	Double Universal 5/11 Extra Long Blue Color	8.07" / 20.5 cm	1.57" / 40 mm	0.79" / 20 mm	0.43" / 11 mm	0.20" / 5 mm
PACKAGING: Bag of 3 (reusable and can be steam sterilized). Volume Packaging: Bag of 50 use ITEM # 40-18511-50						
36-10010	Single end large Red Color	4.3" / 11 cm	1.57" / 40 mm	N/A	0.39" / 10 mm	N/A
PACKAGING: Bag of 3 (single-use recommended). Volume Packaging: Bag of 50 use ITEM # 36-10010-50						
27-15311	Single end large Blue Color	4" / 10.3 cm	1.57" / 40 mm	N/A	0.39" / 10 mm	N/A

BATRIK

V Brushes

EXTERNAL CLEANING BRUSHES

Our External Cleaning Brushes for endoscopes / optical instruments provide excellent cleaning capabilities. With their unique U-shape bristles and innovative design, our brushes are the ultimate tool for cleaning the exterior of your endoscopic instruments. Free from sharp or metallic ends that can scratch or cause damage to valuable equipment, both single-use and reusable models are latex free and feature ergonomic handles.



PRODUCT SELECTION

ITEM #	DESCRIPTION	BRUSH LENGTH	HANDLE LENGTH / (Ø)	TOTAL LENGTH	PACKAGING
36-20005	Single-Use External Cleaning Brush	70 mm / 25 mm	80 mm / 3 mm	150 mm	50 / Box
36-20010	Reusable External Cleaning Brush	90 mm / 25 mm	120 mm / 12 mm	210 mm	3 / Bag

Ordering & Contact

How to Order

Batrik products are available worldwide through our network of authorized distributors.

To become an authorized distributor, please contact:

General Inquiries:

customerservice@surgmed.com

Visit Us Online

www.batrik.com

Explore our expanded product categories, technical sheets, and educational material.