



# Premicath®

A 1Fr PICC line designed specially for newborns under 1kg

Infants born prematurely weighing under 1kg often need the nutritional and therapeutic support that can only be delivered into their tiny vessels via venous access. At just 1Fr, premicath® has been specially designed for the smallest and most fragile babies to deliver optimum accuracy during placement as well as minimising trauma for patients.

Available with choice of introducers including splitting needle and microflash®



## Premicath®

Product Code	Catheter			Introducer		Stylet	Units/Box
	Length	Priming Volume	Flow rate (Pressure 1 bar)	Type	Size		
1261.20G	20 cm	0.09 ml	0.7 ml / min	Splitting Needle	24 G	No	10
1261.203G	20 cm	0.09 ml	0.7 ml / min	Splitting Needle	24 G	Yes	10
1261.205G	20 cm	0.09 ml	0.7 ml / min	Neocath Split	22 G	Yes	10
1261.208G	20 cm	0.09 ml	0.7 ml / min	Microflash	20 G	Yes	10
1261.306G	30 cm	0.11 ml	0.6 ml / min	No Introducer	-	Yes	10
<b>Higher Barium</b>							
1261.203A	20 cm	0.09 ml	0.7 ml / min	Splitting Needle	24 G	Yes	10
1261.306A	30 cm	0.11 ml	0.6 ml / min	No Introducer	-	Yes	10

## Kits & Accessories

Product Code	Description	Units/Box
7371.19	20 G Microflash Introducer	10
1261.20N	24 G Splitting Needle	20
AMS-850TC	V-Cutter Catheter Trimming Tool	20



V-Cutter



Splitting needle



Microflash