


nutriline®

TWINFLO

A 2Fr dual lumen catheter with peelable needle

nutriline® twinflo has been designed specifically for babies who require multiple infusions. Its dual lumen tip design offers enhanced fluid management and avoids incompatible drugs mixing.



Product Code	Catheter			Introducer		Stylet	Units / Case
	Length	Priming Case Volume	Flow rate (pressure 1bar)	Type	Size		
1252.232G	30 cm	2 x 0.2 ml	2 x 1.45 ml / min	No introducer	-	No	10
1252.232M	30 cm	2 x 0.2 ml	2 x 1.45 ml / min	Microflash	20 G	No	10
1252.230G	30 cm	2 x 0.2 ml	2 x 1.45 ml / min	Microflash	20 G	Yes	10

Kits & Accessories

Product Code	Description	Units/Box
7371.19	20 G Microflash Introducer	10
1261.20N	24 G Splitting Needle	20
AMS-850TC	V-Cutter Catheter Trimming Tool	20

