


nutri^{line}[®]

A PICC line with a unique & easy-to-use peelable needle available in 2 Fr

nutri^{line}'s one-piece catheter provides clinicians with a high degree of safety. Available in a choice of lengths to ensure accurate tip placement for most IV access sites.



Product code	Length	Size	Priming volume	Flow rate	Introducer	Stylet	Units/Case
I252.30G	30 cm	2 Fr	0.12 ml	5.0 ml / min	Splitting needle 20 G	No	10
I252.31G	30 cm	2 Fr	0.12 ml	5.0 ml / min	No Introducer	No	10
I252.35G	30 cm	2 Fr	0.12 ml	5.0 ml / min	Microflash 20 G	No	10
I252.030G	30 cm	2 Fr	0.12 ml	5.0 ml / min	Microflash 20 G	Yes	10
I252.03I	30 cm	2 Fr	0.12 ml	5.0 ml / min	No Introducer	Yes	10