INSTRUMENT TABLES

Our Instrument Tables come in six versatile models and are an economical back table for instruments or dressings. Fully welded heavy gauge stainless steel construction polished to a #4 finish provides excellent rigidity and extended product life.

All instrument tables are backed by an outstanding **7 year warranty** on structural and all components (excluding casters).

KEY FEATURES

- Channel support sub-frame for maximum weight capacity and structural rigidity
- Fully welded seamless tubular legs H-frame supports or lower shelf assembly
- Non-marking swivel design casters come in 2" (51 mm) or 3" (76 mm) depending on table size
- All models are 34" (864 mm) high and feature a lower shelf
- Open tubular base design is available upon request (reference # OR-3643)



OR-3641L

PRODUCT SELECTION

ITEM #	DIMENSIONS (L x W x H)
OR-3641A	20" x 16" (506 mm x 406 mm)
OR-3641B	30" x 16" (762 mm x 406 mm)
OR-3641C	34" x 18" (864 mm x 457 mm)
OR-3641D	36" x 20" (914 mm x 508 mm)
OR-3641H*	46" x 24" (1168 mm x 610 mm)
OR-3641L*	58" x 24"(1473 mm x 610 mm)
OR-3641M	60" x 24" x 34" (1524 mm × 610 mm × 864 mm)
OR-3641-48W	48" x 24" x 38" (1219 mm × 610 mm × 965 mm) (No shelf)
OR-3641-30W	72" x 30" x 36" (1829 mm × 762 mm × 914 mm)
OR-3641-24W	72" x 24" x 34" (1829 mm × 610 mm × 864 mm)
OR-3641G	48" x 24" x 34" (1219 mm × 610 mm × 864 mm)

*H and L models supplied with 3" / 76 mm caster assemblies

ACCESSORIES

ITEM #	DESCRIPTION
OR-3641X-R3	THREE-SIDED GUIDE RAIL Assembly for Table Top or Lower Shelf (Specify size)
OR-3641X-R4	FOUR-SIDED GUIDE RAIL Assembly for Table Top or Lower Shelf (Specify size)
OR-3640	DRAWER ASSEMBLY Inside Dimensions: 15"L x 16 3/4" W x 3 3/4" D (381 mm x 426 mm x 95 mm)